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Armando Moretti

ID.10301

This manual is the property of DORSAL®

Reproduction, disclosure and translation into any language, transfer and copying, in any form, even partially, including any magnetic or typewritten media, to third parties without the consent of Dorsal® is prohibited.

Storage of the instruction manual

Keep the instruction manual in a suitable place for its maintenance and, in any case, in the location where the Anatomical Support is used. The document should be kept for the entire lifespan of the Anatomical Support to serve as an adequate reference for any possible needs.

The devices **XPoint, Experience, Human, Ergo, Forest, Sana, Super, Super Tonic, Super Flexi, Superego, and Supernova** are distributed by

Dorsal s.r.l.
via Madonna di Loreto 3
31020 Corbanese (TV) – Italy
C.F e P.iva 00317220267
Ph. +39 0438 933003
info@dorsal.it
www.dorsal.it

General index

1	Terms and definition	50
1.1	Glossary	50
2	Technical characteristics	52
2.1	Device	52
2.2	Materials used	54
2.2.1	Anatomical Supports with perimeter frame in plywood	54
2.2.2	Anatomical Supports with perimeter frame in steel	55
2.3	Technical data related to the motor	56
2.3.1	Version RA – FRA – MRA with wired handset	56
2.3.2	Version RA – FRA – MRA with Bluetooth® remote	57
2.3.3	Version R-4 with Bluetooth® remote	57
2.3.4	Experience 600 version RT – FRT	58
2.3.5	Experience 500 version RT – FRT	58
2.4	Conditions of use and storage of the device	59
3	Warranty conditions	59
4	Intended use	59
5	Device classification	59
6	Warnings and contraindications	60
6.1	Mandatory behaviour	60
6.2	Warnings before use	61
6.3	Warnings during use	62
6.4	Batteries	62
6.5	Residual risks	63
7	Principles of operation	64
8	Standard included	64
9	Installation of the device	65
9.1	Safety precautions for the user	65
9.2	Installation	65
10	Description of the Anatomical Support (motorized versions)	67
11	Handset	68
11.1	Version RA – FRA – MRA with cable control	68
11.2	Version RA – FRA – MRA with Bluetooth® remote control	69
11.2.1	Synchronization for bluetooth® remote to RA – FRA – MRA motor	70
11.2.2	Synchronization for two RA Bluetooth® motorized supports with two remotes and synchronization cable (optional)	70

11.2.3	Synchronization for two RA Bluetooth® motorized Supports with one remote	70
11.2.4	Memories	70
11.3	Version R-4.....	71
11.3.1	Bluetooth® remote synchronization for the R-4	72
11.3.2	Synchronization for two R-4 Bluetooth® motorized Supports with two remotes and synchronization cable (optional).....	72
11.3.3	Synchronization for two Bluetooth® R-4 motorized supports with one remote.....	72
11.3.4	Saving memory positions for Bluetooth® R-4 motorized Support	72
11.3.5	Procedure to restore the Zero Gravity and Anti-snore positions for Bluetooth® R-4 motorized Supports.....	73
11.4	APP Dorsal Bliss from smartphone	74
11.4.1	Procedure to operate RA BT and R-4 motorized supports with Bluetooth® technology from the smartphone app.....	74
11.4.2	Procedure to operate two or more motorized Supports with Bluetooth® technology from the app	74
11.5	Version RT – FRT with cable control.....	75
11.5.1	Synchronization of two motors with a single wired remote control	75
11.5.2	Synchronization of two motors with two wired remote controls	75
11.6	Version RT - FRT with RF remote control.....	76
11.6.1	Installation of RF radiofrequency receiver.....	77
11.6.2	Synchronization of two motors with a single RF remote control.....	77
11.6.3	Synchronization of two motors with two RF remote controls	77
12	Mode of use.....	77
12.1	General instructions	77
12.2	Behaviours to adopt in case of electrical blackout.....	78
12.2.1	Emergency system for RA and R-4 versions	78
12.2.2	Emergency system for Experience 600 RT – FRT versions.....	78
12.2.3	Emergency system for Experience 500 RT – FRT versions.....	79
13	Installation of the independent battery power kit for Experience RT – FRT models	80
13.1	Instructions for correct use of the Experience battery power kit	80
14	Maintenance and cleaning of the device	81
14.1	Periodic maintenance.....	82
14.2	Spare parts.....	82
14.3	Storage instructions.....	82
14.4	Disposal.....	83
15	Labelling and symbols.....	84
15.1	Device	85
16	Declaration of conformity	86

1 Terms and definitions

1.1 Glossary

Anatomical Support:	This refers to the detachable bed frame as defined in technical standards and is identified with the items from the Dorsal® catalogue models XPoint, Experience, Human, Ergo, Forest, Sana, Super, Super Tonic, Super Flexi, Superego, Supernova in the versions NA, FA, FXA, RA, RT, FRT, R-4, FRA, MRA . If a citation refers to only one of the mentioned variants, this will be clearly specified in the text.
Version:	Defines the presence and type of movement for the headboard and footboard.
Single Bed:	The surface of the bed intended for one person to rest. The width of a single anatomical support does not exceed 120 cm.
Double Bed:	The surface of the bed of such dimensions that it allows two people to rest. The width of an anatomical support for double use is greater than 140 cm. A double anatomical support is equivalent to two single beds.
NA Version:	Detachable single bed without joints and without the possibility of adjusting the headboard and footboard.
FA Version:	Detachable double bed without joints and without the possibility of adjusting the headboard and footboard.
FXA Version:	Detachable double bed without joints and without the possibility of adjusting the headboard and footboard. In the Superego model, this version allows for the combination of two different comforts in a single double version.
RA Version:	Detachable single bed, with the possibility of adjusting the headboard and footboard via a motor controlled by electric pulses with maintained action. Available with a wired remote control or Bluetooth® technology (RA BT).
RT Version:	Detachable single bed with the possibility of adjusting the headboard and footboard via an integrated motor controlled by electric pulses with maintained action. Available in 6, 5, or 4-section movement versions, with wired or radio frequency (RF) technology.

- R-4 Version: Detachable single bed with the possibility of adjusting the headboard and footboard via a motor controlled by electric pulses with maintained action. Additionally, this version allows for the independent adjustment of the head and foot angles using two electric actuators controlled by a single remote.
- FRA Version: Detachable double bed with the possibility of independently adjusting each headboard and footboard via two motors controlled by electric pulses with maintained action.
Available in versions with a wired remote or Bluetooth® remote (FRA BT).
- FRT Version: Detachable double bed with the possibility of independently adjusting each headboard and footboard via two integrated motors controlled by electric pulses with maintained action.
Available in 6- and 5-section movement versions with wired or radio frequency (RF) remote.
- MRA Version: Detachable single slat double bed with the possibility of simultaneously adjusting the headboard and footboard via a motor controlled by electric pulses with maintained action.
Available in versions with a wired remote or Bluetooth® technology (MRA BT).

2 Technical characteristics

2.1 Device

Model	Version	Available Widths (cm)	Available Lengths (cm)	Outer Frame	Ergonomically recommended heights (*1) (cm)	Safe working load (kg) (*2)	Person weight (kg) (*2)	Weight of Device (*3) (kg) ±5%
					IEC 60601	IEC 60601	IEC 60601	
XPoint	NA	80, 85, 90, 120	190, 195, 200	Multilayer Beech profile 60x30 mm	30, 35, 40	173	160	18
	FA	140, 160, 165, 170, 180	190, 195, 200		37			
	RA	80, 85, 90, 120	190, 195, 200		33			
	R-4	80, 85, 90, 120	190, 195, 200		40			
Experience	RT	80, 85, 90, 120	190, 195, 200	Multilayer Beech profile 70x25 mm	40	150	130	27
	FRT* Available only in the 600 version with 6 sections and the 500 version with 5 sections	160, 165, 170, 180	190, 195, 200					54
Human	NA	80, 85, 90, 120* * Size not available in the Human Active version	190, 195, 200	Multilayer Beech profile 60x30 mm	30, 35, 40	173	160	16
	FA	140, 160, 165, 170	190, 195, 200		32			
	RA	80, 85, 90, 120* * Size not available in the Human Active version	190, 195, 200		30			
	R-4	80, 85, 90, 120* * Size not available in the Human Active version	190, 195, 200		38			
	FRA	140, 160, 165, 170	190, 195, 200		60			
Ergo	NA	80, 85, 90, 120	190, 195, 200	Multilayer Beech profile 60x30 mm	30, 35, 40	173	160	16
	FA	140, 160, 165, 170	190, 195, 200		32			
	RA	80, 85, 90, 120	190, 195, 200		40			32
Forest	NA	80, 85, 90, 120	190, 195, 200	Multilayer Beech profile 60x25 mm.	30, 35, 40	173	160	13
	FA	140, 160, 165, 170	190, 195, 200		26			
	RA	80, 85, 90, 120	190, 195, 200		40			29

Sana	NA	80, 85, 90, 100, 120	190, 195, 200	Multilayer Beech profile 60x30 mm.	30, 35, 40	173	160	14
	FA	140, 160, 165, 170, 180	190, 195, 200		28			
	RA	80, 85, 90, 100, 120	190, 195, 200		29			
	R-4	80, 85, 90, 120	190, 195, 200		37			
	FRA	140, 160, 165, 170,180	190, 195, 200		58			
Super	NA	80, 85, 90, 100, 120	190, 195, 200	Cold-rolled Steel Section 40x30 mm. Section 50x30 mm in sizes 120, 180, and in the MRA version."	30, 35, 40	200	160	14
	FA	140, 160, 165, 170, 180	190, 195, 200		28			
	RA	80, 85, 90, 120	190, 195, 200		28			
	FRA	140, 160, 165, 170, 180	190, 195, 200		56			
	MRA	140, 160, 165, 170	190, 195, 200		400	320	49	
Super Tonic	NA	80, 85, 90, 100, 120	190, 195, 200	Cold-rolled Steel Section 40x30 mm.	30, 35, 40	200	160	14
	FA	140, 160, 165, 170, 180	190, 195, 200		28			
	RA	80, 85, 90, 100, 120	190, 195, 200		40			29
Super Flexi	NA	80, 85, 90	190, 195, 200	Cold-rolled Steel Section 40x30 mm.	30, 35, 40	173	160	13
	FA	140, 160, 165, 170, 180	190, 195, 200		26			
	RA	80, 85, 90	190, 195, 200		40			27
Superego	NA	80, 85, 90, 120	190, 195, 200	Cold-rolled Steel Section 40x30 mm.	30, 35, 40	200	160	15
	FA	140, 160, 165, 170	190, 195, 200		30			
	RA	80, 85, 90, 120	190, 195, 200		40			31
	FXA	140, 160, 165, 170	190, 195, 200		30, 35, 40			30
Supernova	NA	80, 85, 90, 120	190, 195, 200	Cold-rolled Steel Section 40x30 mm.	30, 35, 40	200	160	15
	FA	140, 160, 165, 170	190, 195, 200		30			
	RA	80, 85, 90, 120	190, 195, 200		40			31

1. Height calculated from the floor to the top edge of the Anatomical Support (mattress excluded). The minimum recommended height is to be considered mandatory in the case of the motorized version. The maximum recommended height is to be considered mandatory in the case of the Anatomical Support (with mounted feet) used independently, meaning not installed within a bed structure. In any case, if the Anatomical Support is used at a height greater than the recommended one, it is the user's responsibility to ensure the presence and compliance of safety features in the adopted bed structure.

2. The indicated weight refers to the load for each individual bed.

3. The measured weight refers to the product without packaging and without mounted feet, measuring 80x190 cm for the NA, RA, R4, and RT versions, and 160x190 cm for the FA, FRA, FRT, MRA, and FXA versions.

2.1 Materials used

2.2.1 Anatomical Supports with perimeter frame in plywood

FRAME

Beech wood plywood made using non-toxic resins.

SLATS

Beech wood plywood made using non-toxic resins.

SLATS HOLDERS

Suspensions in thermoplastic rubber Hytrel and SBS.
Caps with polyethylene.

FEET

Optional upon request. Compatible only with Dorsal models.
Material: solid beech.
Finish: polyurethane varnish.

MOVEMENT OF THE JOINTS VERSION RA – R-4 – FRA

Head-section and foot-section: joints in plywood and steel levers controlled by the push of an electric pulse motor with maintained action commands.

MOVEMENT OF THE JOINTS VERSION RT-FRT

Head-section and foot-section: joints in plywood and thermoplastic plastic levers, controlled by the push of an integrated electric pulse motor with maintained action commands.

2.2.2 Anatomical Supports with perimeter frame in steel

FRAME

Material: cold-rolled steel E195.
Finish: epoxy polyester paint or thermosetting polyester powder.

SLATS

Beech plywood made using non-toxic resins.

SLATS HOLDERS

Suspensions in thermoplastic rubber Hytrel and SBS.
Caps with polyethylene.

FEET

Optional upon request. Compatible only with Dorsal models.

Material: cold-rolled steel E195.

Finish: epoxy polyester paint or thermosetting polyester powder.

MOVEMENT OF THE JOINTS VERSION RA – FRA – MRA

Head-section and foot-section: joints and levers in steel E195 controlled by the push of an electric pulse motor with maintained action commands.

2.3 Technical data related to the motor

2.3.1 Version RA – FRA – MRA with wired handset

Power supply:	100-240 VAC 50/60 Hz
Insulation class:	II – No grounding required
Protection rating:	IP20 – Not waterproof, protected against access with a finger

Motor and remote control: DC, low voltage 24V DC

Standby consumption: less than 0.5 W

Fuse with thermal switch for overload protection

Maximum connection time: 10% with continuous operation of 2/18 minutes

Maximum push when operating the head-section or foot-section separately: 6500 N, equivalent to 662 kg (*)

(*) For Forest, Supernova, and Supernova Plus models: 4500 N, equivalent to 458 kg

Device for lowering the joints to a flat position in case of an electrical blackout (includes 9V alkaline battery).

2.3.2 Version RA – FRA – MRA with Bluetooth® remote

Power supply:	100-240 VAC 50/60 Hz
Insulation class:	II – No grounding required
Protection rating:	IP20 – Not waterproof, protected against access with a finger
Motor and remote control:	DC, low voltage 24V DC
Independent Bluetooth® remote control,	powered by two included AAA batteries
Standby consumption:	less than 0.5 W
Fuse with thermal switch for overload protection	
Maximum connection time:	10% with continuous operation of 2/18 minutes
Power required depending on the load and models:	30-300 W
Device for lowering the joints to a flat position in case of an electrical blackout (includes 9V alkaline battery)	
Courtesy light	
2 memory positions	
Integrated Bluetooth® communication	
Maximum push when operating the head-section or foot-section separately:	6500 N, equivalent to 663 kg

2.3.3 Version R-4 with Bluetooth® remote

Motor:

Power supply:	100-240 VAC 50/60 Hz
Insulation class:	II – No grounding required
Protection rating:	IP20 – Not waterproof, protected against access with a finger
Motor and remote control:	DC, low voltage 24V DC
Independent Bluetooth® remote control,	powered by two included AAA batteries
Standby consumption:	less than 0.5 W
Fuse with thermal switch for overload protection	
Maximum connection time:	10% with continuous operation of 2/18 minutes
Power required depending on the load and models:	30-300 W
Device for lowering the joints to a flat position in case of an electrical blackout (includes 9V alkaline battery)	
Courtesy light	
4 memory positions	
Maximum push when operating the head-section or foot-section separately:	6500 N pari a 663 kg

Monoactuators:

Protection class: II – No grounding required
 Protection rating: IP20 – Not waterproof, protected against access with a finger

Motor: DC, low voltage 24V DC

Maximum push of monoactuator for head-section and foot-section: 4500 N, equivalent to 459 kg

Lowering power of monoactuator for head-section and foot-section: 3000 N, equivalent to 306 kg

Integrated Bluetooth® communication via control unit.

Device for lowering the joints to a flat position in case of an electrical blackout (includes two 9V emergency batteries).

2.3.4 Experience 600 version RT – FRT

Power supply: 100-240 V 50/60 Hz
 Insulation class: II – No grounding required
 Protection rating: IP20 – Not waterproof, protected against access with a finger

Motor and remote control: DC, low voltage 24V DC

Standby consumption: less than 0.5 W

Maximum connection time: 10% with continuous operation of 2/10 minutes

Maximum push when operating the head-section or foot-section separately: 7500 N, equivalent to 459 kg

Mechanical joint lowering system to a flat position in case of an electrical blackout.

2.3.5 Experience 500 version RT – FRT

Power supply: 100-240 V 50/60 Hz
 Insulation class: II – No grounding required
 Protection rating: IP20 – Not waterproof, protected against access with a finger

Motor and remote control: DC, low voltage 24V DC

Standby consumption: less than 0.5 W

Maximum connection time: 10% with continuous operation of 2/10 minutes

Maximum push when operating the head-section or foot-section separately: 7500 N, equivalent to 765 kg

Mechanical joint lowering system to a flat position in case of an electrical blackout.

2.4 Conditions of use and storage of the device

Operating temperature:	from 5 °C to 40 °C
Maximum relative humidity:	max 50 %
Storage and transport temperature:	from -25 °C to 55 °C (up to 70°C if not exceeding 24 hours)

Dorsal Anatomical Supports are designed to operate safely for objects and people, within the limits and conditions specified in this manual, for approximately 12 years. Improper use of the device reduces its expected lifespan.

3 Warranty conditions

The Anatomical Supports are guaranteed for 2 years. Any claims must be reported to the authorized Dorsal® dealer.

The warranty period can be extended free of charge to 5 years (for purchases made within Italy) by completing the form on the website www.garanzia.dorsal.it. Electrical parts are excluded from this extension. For further details, carefully read the warranty certificate included with the product.

4 Intended use

The Anatomical Supports XPoint, Experience, Human, Ergo, Forest, Sana, Super, Super Tonic, Super Flexi, Superego, and Supernova are devices designed to improve sleep quality. They are classified as Medical Devices with the following intended uses:

- improvement of joint comfort with a reduction of rheumatic pain;
- reduction of insomnia and its physical and psycho-emotional effects;
- improvement of posture, particularly of the spine, with a reduction of associated painful symptoms;
- reduction of thermal and mechanical effects on the skin.

The use of these devices is not intended for individuals with musculoskeletal, vascular, or nervous (central or peripheral) system problems that specifically require the use of different types of supports.

5 Device classification

The Anatomical Supports XPoint, Experience, Human, Ergo, Forest, Sana, Super, Super Tonic, Super Flexi, Superego, and Supernova are classified as Class I Medical Devices in accordance with Rule 1 of Annex VIII, in compliance with EU Regulation 2017/745.

6 Warnings and contraindications

6.1 Mandatory behaviour

These instructions for use contain essential information that must be followed during installation, use, and maintenance. For this reason, these instructions must be carefully read by the user. The instructions described here must be followed in addition to general home safety precautions.

To maintain the highest level of safety, the user must operate the device in accordance with the specifications and limits outlined in the user manual.

Failure to observe safety precautions will void the warranty rights and may cause harm to individuals or other risks not specifically described in this instruction manual.

The manufacturer disclaims any responsibility for improper use as indicated and prescribed in this manual.

To ensure the safety of the user, the operators assisting with the device, and to safeguard the longevity of the product, it is essential to follow the information provided by the manufacturer.

It is prohibited to remove labels indicating danger conditions. Refer to the Maintenance section if the labels are damaged and/or illegible.

The Anatomical Support must be used in normal positions for correct rest. It should not be subjected to or used for improper purposes (e.g., jumping).

To avoid tipping risks, the user must take into account the maximum load allowed on the long and short sides, as well as the safe operating load suggested by the manufacturer (see paragraph 2).

The Anatomical Support should not be placed in damp locations or near sources of ignition or heat sources exceeding 100°C.

Do not use in environments with potentially explosive or flammable atmospheres.

If the Anatomical Support is used with feet, these should not exceed the height recommended by the manufacturer (see paragraph 2).

The motorized version of the Anatomical Support must be connected to the electrical network through its own power cable. Before performing this operation, ensure that the installation complies with the regulations in force in your country.

Check that the voltage and frequency of the network match the data indicated on the motor label.

The Anatomical Support is not suitable for outdoor use, should not be placed in damp locations, and should not be exposed to direct sunlight.

6.2 Warnings before use

The handling of the product (packaged and unpackaged) must be carried out under safe conditions and with appropriate handling methods. It must be done using means that do not damage the product. Unloading and handling, if done manually, must be performed by at least two people. For the weight of the product, refer to paragraph 2.1 .

In addition to the Anatomical Support, the optimal bed system includes the use of a mattress and pillow. It is recommended to use a mattress and pillow with high ergonomic adaptability and pressure load distribution capabilities, preferably of the Class I Medical Device type.

To ensure the safety of the user and the operators assisting with the device, it is not recommended to use multiple power strips connected directly to the motor of the motorized version, as they are not compatible with the safe use of the Medical Device.

If the Anatomical Support is used with feet, verify that they are properly screwed to the frame and that they all rest on the ground and are level. It is MANDATORY to use only Dorsal® products.

If the feet are not used, it must be confirmed in advance that the structure where the Anatomical Support is installed is properly sized and adequately supportive.

For the motorized version, check that nothing obstructs or interferes with the movement of the movable sections of the Anatomical Support.

For the motorized version, check that the power cables:

- do not interfere with the movement of the sections;
- are not placed in walkways, and if this is unavoidable, they must be adequately protected by conduits or other types of protection to allow safe passage;
- do not intersect with the power cables of other electrical devices in a way that obstructs the movement of the joints or causes damage.

6.3 Warnings during use

For the motorized versions **RA, RT, FRA, FRT, MRA, and R-4**, the head-section and foot-section movement system is not a continuous-use device. Follow the usage times indicated in the paragraph 2.3.

6.4 Batteries

Attention

- Keep the batteries out of reach of children.
- Batteries must not be thrown into fire, short-circuited, disassembled, or charged. Risk of explosion!
- Batteries contain harmful acids. Replace weak batteries as soon as possible to avoid leakage.
- Never use new and old batteries together, and do not mix different types of batteries.

Battery leakage handling:

- when replacing leaking batteries, wear gloves and chemical-resistant safety glasses!
- avoid contact with bare hands.
- if contact occurs with leaked battery components, carefully wash your hands.
- remove any electrolyte residue with a damp cloth. Wash all clothing that came into contact with the electrolyte.

6.5 Residual risks

The following risks cannot be eliminated:

Risk of crushing

WARNING! Do not insert hands or objects between the moving sections and the bed frame.

The motorized versions of the Anatomical Supports are equipped with an electric motor that brings the joints to a flat position at a very slow speed. However, maximum caution is required during the movement and adjustment of the sections of the Anatomical Support to ensure that no one inserts hands or other body parts, or any foreign objects, between the moving sections. See also paragraph 15.

There are some labels on the head-section and foot-section of the adjustable internal frame. These labels indicate areas of higher risk of crushing and are clearly visible during the section lifting phase.

If an object is caught between the joints during the downward movement, return the movable elements to a position that clears the obstruction. Do not remove the object without first lifting the movable part. Removing the object without first lifting the movable part will result in the part falling.

The user or the operator assisting with the device must pay close attention, ensuring that their hands or those of others are not placed in areas of higher risk.

Risk of falling in case of joint collapse

WARNING! It is prohibited to sit at the ends of the sections when they are raised.

This risk occurs only when a person sits on the head-section or foot-section when they are not in the flat position of the motorized version of the Anatomical Support. There are some labels on the head-section and foot-section of the adjustable internal frame. See also the paragraph 15. These labels indicate areas of higher risk of joint collapse and are clearly visible during the section lifting phase.

It is absolutely forbidden to lie down or sit on the Anatomical Support in the motorized version when the head-section or foot-section is adjusted to a non-flat position.

The collapse of the joints, in addition to significantly increasing the chances of crushing limbs or the user falling from the Anatomical Support, can damage the

functionality of the joints themselves and compromise the safety of the Anatomical Support.

7 Principles of operation

The function of the Anatomical Support is to adequately support the mattress and, together with the pillow, create an ideal Bed System for rest. **XPoint, Experience, Human, Ergo, Forest, Sana, Super, Super Tonic, Super Flexi, Superego, and Supernova** are Anatomical Supports suitable for the environments specified in paragraph 4. They are designed to provide balanced support for the mattress and contribute to pressure load distribution from the user's body.

The versions **NA, FA, FXA, RA, RT, R-4, FRA, FRT, MRA** are available in various sizes, as specified in the paragraph 2.

The motorized versions (**RA, RT, R-4, FRA, FRT, MRA**) are particularly indicated for:

- people with mild or moderate venous and lymphatic stasis in the lower limbs who can benefit from the declined positioning of the lower limbs;
- people with respiratory issues who need an inclined trunk position to aid expectoration;
- people with heart conditions and breathing difficulties due to pulmonary stasis;
- people who, unable to get up independently or with proper assistance, need to be fed in bed and thus benefit from correct trunk and head positioning to facilitate physiological swallowing.

8 Standard included

A. Nr. 1 Anatomical Support

B. Nr. 1 foot-section mattress stopper (only in motorized versions **RA, FRA, MRA, R4, RT, FRT**)

Optional accessories

C. Nr. 4 supportive feet for standing alone (only upon request)

9 Installation of the device

9.1 Safety precautions for the user

The electrical power system of the motorized version must comply with current regulations for injury prevention and safety in domestic environments.

In particular, the electrical safety of the motor is only ensured if it is properly connected to a system that complies with current safety standards.

9.2 Installation


The installation of the Anatomical Supports can be performed by specialized personnel or directly by the user, who must adhere to the following instructions. If there are doubts regarding the installation procedure, completeness of the equipment, or the integrity of the product after unpacking, contact the authorized Dorsal® dealer.

The Anatomical Support is supplied in an appropriate package that varies depending on the models and versions.


- Unpack the Anatomical Support, placing the packaging so that it does not interfere with the installation process.
- Keep packaging elements out of the reach of children and do not dispose of them in the environment. Carefully read the manual containing installation, use, and maintenance instructions for the product.
- For the disposal of packaging materials, refer to paragraph 14.4 .
- Keep the manual in a suitable location, preferably in the place where the support is installed and used.
- For motorized versions, it is recommended to perform a complete movement cycle in advance to verify the correct assembly of the system.

Assembly of the feet on the Anatomical Supports

If provided, screw the feet into the designated slots, ensuring that the screws are fully tightened.

Assembly of feet on plywood supports.


Assembly of feet on steel supports.


Dragging the product with the feet mounted may cause damage to the feet and compromise the integrity of the anchoring points. The handling of the support must be carried out by 2 operators, lifting the Anatomical Support.

Release the head-section and foot-section (motorized versions).


In the RA, RT, R-4, FRA, FRT, MRA versions, release the head-section and foot-section by cutting the plastic straps.

Caution: Perform this operation only once the support is positioned horizontally!


Assembly of the foot-section mattress stopper

In the case of the motorized versions **RA, RT, R-4, FRA, FRT, MRA**, also install the foot-side mattress stopper.

Connect the power cable (motorized versions)

In the RA, RT, R-4, FRA, FRT, MRA versions, connect the power cable to the electrical socket. Also, refer to paragraph 6.1 .


Perform a complete test cycle by fully raising and then fully lowering the head-section and foot-section.

Completion of installation

Place the mattress on top of the Anatomical Support.



10 Description of the Anatomical Support (motorized versions)

Placement of the main components and labels.

11 Handsets


11.1 Version RA – FRA – MRA with cable control

Head-section


Foot-section


Head and Foot sections

1	Head-section raise
2	Head-section lower
3	Foot-section raise
4	Foot-section lower
5	Simultaneous raise of the head-section and foot-section
6	Simultaneous lower of the head-section and foot-section


11.2 Version RA – FRA – MRA with Bluetooth® remote control

Head-section

Foot-section

Head and Foot sections

1	Head-section raise
2	Simultaneous lower of the head-section and foot-section
3	Foot-section raise
4	Foot-section lower
5	Head-section lower
6	Zero Gravity position /memory
7	Torch light on remote
8	Anti-snore position/ memory
9	Courtesy light positioned under the bed
10	Synchronize ON of two RA Bluetooth® motorized supports
11	Synchronize OFF of two RA Bluetooth® motorized supports

11.2.1 Synchronization for Bluetooth® remote to RA – FRA – MRA motor

The Bluetooth® motor is already synchronized with its remote control.

If re-synchronization is needed, follow these steps:

- unplug the power cord and wait for 2 minutes.
- reconnect the power cord to the socket. Within 10 seconds, press buttons 3 and 5 simultaneously for 2 seconds. If synchronization is successful, a double beep will be heard from the motor.

11.2.2 Synchronization for two RA Bluetooth® motorized Supports with two remotes and synchronization cable (optional)

Connect the synchronization cable (optional, not provided with the support) between the two motors, and on one of the remotes, press button 10 "SYNC ON" for 3 seconds. A beep will indicate that the two beds are synchronized, and both remotes can now be used.

Upon the first use of the remotes, the two supports will align to the same position. Even if synchronized, each motor will retain its own Zero Gravity and Anti-snore position in memory.

To disconnect the two supports and make them independent again, press button 11 "SYNC OFF" for 3 seconds on one of the remotes. A beep will indicate that they are no longer synchronized.

11.2.3 Synchronization for two RA Bluetooth® motorized supports with one remote

To pair the second motor with a single Bluetooth® remote, follow the procedure outlined in paragraph 11.2.1 using the remote control of the first motor (do not unplug the power cord of the first motor).

11.2.4 Memories

The Zero Gravity and Anti-snore positions are pre-programmed at the factory. They can always be changed using the following procedure:


- stop the motor in the desired position and press buttons 6 or 8 for 3 seconds to save a different position.
- press buttons 6 or 8 for 1 second to recall the stored positions.

To restore the factory settings, press buttons 2 and 5 for 3 seconds.

11.3 Version R-4

1	Torch light on remote
2	Zero Gravity position /memory
3	Memory position 1
4	Anti-head section raise
5	Anti-snore position/ memory
6	Memory position 2
7	Anti-foot section lower
8	Anti-foot section raise
9	Anti-head section lower
10	Head section raise and lower
11	Courtesy light positioned under the bed
12	Simultaneous raise and descent of the anti-head-section+head-section and the anti-foot-section+foot-section
13	Foot-section raise and lower
14	Synchronize ON/OFF of two R-4 Bluetooth® motorized Supports

11.3.1 Bluetooth® remote synchronization for the R-4 motor

The Bluetooth® motor is already synchronized with its remote control.

If re-synchronization is needed, follow these steps:

- unplug the power cord and wait for 2 minutes.
- reconnect the power cord to the socket. Within 10 seconds, press buttons 8 and 9 simultaneously for 2 seconds. If synchronization is successful, a double beep will be heard from the motor.

11.3.2 Synchronization for two R-4 Bluetooth® motorized Supports with two remotes and synchronization cable (optional)

Connect the synchronization cable (optional, not provided with the support) between the two motors, and on one of the remotes, press button 14 “SYNC ON” for 3 seconds. A beep will indicate that the two beds are synchronized, and both remotes can now be used. Upon the first use of the remotes, the two supports will align to the same position. Even if synchronized, each motor will retain its own Zero Gravity and Anti-snore position in memory. To disconnect the two supports and make them independent again, press button 14 “SYNC OFF” for 3 seconds on one of the remotes. A beep will indicate that they are no longer synchronized.

11.3.3 Synchronization for two Bluetooth® R-4 motorized Supports with one remote

To pair the second motor with a single Bluetooth® remote for R-4, follow the procedure outlined in section 11.3.1 using the remote control of the first motor (do not unplug the power cord of the first motor).

11.3.4 Saving memory positions for R-4 Bluetooth® motorized Support

The Zero Gravity and Anti-snore positions are pre-programmed at the factory. They can always be changed. Simply stop the motor in the desired position and press buttons 2 or 5 for 3 seconds to save a different position. To save the desired positions, press buttons 3 or 6 for 3 seconds, until a confirmation beep is heard. Press buttons 2, 3, 5, or 6 for 1 second to recall the saved positions.


Do not look directly at the LED!

11.3.5 Procedure to restore the Zero Gravity and Anti-snore positions for Bluetooth® R-4 motorized support

Press and hold button 12  until the support returns to a flat position and a beep is heard.

To restore the Zero Gravity position:


- return the support to the position indicated in the following figure, using buttons 4, 7, 8, 9, 10, and 13:

- save the reached position by holding button 2 for 3 seconds, until a beep is heard.

To restore the Anti-snore position:

- return the support to the position indicated in the following figure, using buttons 4 and 9:

- save the reached position by holding button 5 for 3 seconds, until a beep is heard.

11.4 APP Dorsal Bliss from smart phones

It is possible to download the free Dorsal app from Google Play for Android and from the App Store for iOS

ANDROID – Google Play™

IOS – Apple Store®

11.4.1 Procedure to operate RA BT and R-4 motorized supports with Bluetooth® technology from the smartphone app

- Download the 'Dorsal Bliss' app and launch it.
- Ensure that the Bluetooth® function on the smartphone is active.
- Accept all security settings. If not accepted, the app will not work properly.
- Select the virtual remote control for RA BT or R-4.
- Activate the desired functions directly from the touch screen.
- You can operate up to ten different Bluetooth® motorized supports from the Dorsal Bliss app.


11.4.2 Procedure to operate two or more motorized supports with Bluetooth® technology from the app

If there are other Bluetooth® motorized Dorsal supports within a 10-meter range, the app will detect all available devices.


To operate one of the RA BT or R-4 motorized supports, select the "Settings" button within the app and choose the desired device.

Activate the desired functions directly from the smartphone.


11.5 Version RT – FRT with cable control

Head-section

Foot-section

Head and Foot sections

1	Head-section raise
2	Head-section lower
3	Foot-section raise
4	Foot-section lower
5	Simultaneous raise of the head-section and foot-section
6	Simultaneous lower of the head-section and foot-section

11.5.1 Synchronization of two motors with a single wired remote control

Connect both ends of the synchronization cable to the 6-pin ports of the two motors and insert the remote control cable into one of the synchro boxes.


Now both beds can be moved simultaneously.

11.5.2 Synchronization of two motors with two wired remote controls


Connect both ends of the synchronization cable to the 6-pin ports of the two motors and insert the remote control cables into the two synchro boxes.

Now both beds can be moved simultaneously with both remote controls.


11.6 Version with RF remote control RT - FRT

Head-section

Foot-section

Head and Foot sections

1	Head-section raise
2	Head-section lower
3	Foot-section raise
4	Foot-section lower
5	Simultaneous raise of the head-section and foot-section
6	Simultaneous lower of the head-section and foot-section
7	Allows activation of the courtesy light by pressing buttons 5 and 6 simultaneously.

11.6.1 Installation of RF radiofrequency receiver

Connect the radiofrequency receiver to the motor of the support via the dedicated port. The remote control is already synchronized with its receiver.

11.6.2 Synchronization of two motors with a single RF remote control

Connect both ends of the synchronization cable to the 6-pin ports of the two motors and insert the radiofrequency receiver cable into one of the synchro boxes. Now both beds can be moved simultaneously.

11.6.3 Synchronization of two motors with two RF remote controls

Connect both ends of the synchronization cable to the 6-pin ports of the two motors and insert the radiofrequency receiver cables into the two synchro boxes. Now both beds can be moved simultaneously with both remote controls.

12 Mode of use

12.1 General Instructions

The support (motorized versions) can be used safely by keeping in mind all the following information:

- the adjustment can be performed by the user while lying in bed or by an external operator;
- it is essential that the user only access the bed when all sections are in the flat position;
- the movement of the movable sections is possible by pressing the buttons on the supplied remote control (see paragraph 11);
- in case of an electrical blackout, there is a system that allows the joints to return to the flat position (see paragraph 12.2).

NOTE: For 'warnings before and during use' refer to paragraph 6.

12.2 Behaviours to adopt in case of electrical blackout

12.2.1 Emergency system for RA and R-4 versions

The motor is equipped with a system that can return the support to the horizontal position in the event of a sudden power outage. This device cannot be used for the normal movement of the Anatomical Support and should not be considered as an alternative to the electrical supply.

This device, which should only be used in cases of absolute necessity, is powered by a non-rechargeable 9V alkaline battery (3 alkaline 9V batteries for the R-4 version) that Dorsal® provides in combination with the motor.


The installation of the battery that powers the emergency device must be performed by the user in the event of an emergency.

Once the auxiliary battery is installed, press the foot-section lowering command and the head-section lowering command on the remote control to return the sections to the flat position.


The head-section and foot-section must be returned to the flat position separately, first activating one and then the other command.

Although the battery can support several emergency descents, it is recommended to replace it periodically or after a maximum of two years, due to its deterioration.

12.2.2 Emergency system for Experience 600 RT - FRT versions


In case of power failure, the backrest and foot section can be lowered manually.

To achieve the emergency lowering, remove the mattress and pull the red lever of the section to be lowered towards the centre of the bed.


Important! No one should be on the bed.

With the lever pulled, press on the corresponding area until it reaches the horizontal position. Once the power is restored, the system can be repositioned as usual.

12.2.3 Emergency system for Experience 500 RT - FRT versions


In case of power failure, the backrest and foot section can be lowered manually.

Remove the mattress and act on the head-section and foot-section joints.

Manually raise the head-section or foot-section.

Caution! During this procedure, it is necessary to support the end of the head-section or foot-section to prevent it from falling.

Using a screwdriver, remove the steel cable that ends with a small ball from the plastic joint that holds it.

Guide the downward movement of the movable section to the flat position. Repeat the operation on the opposite side of the Anatomical Support.

Important! No one should be on the bed. Once the power is restored, the system can be repositioned as usual.

13 Installation of the independent battery power kit for Experience RT - FRT models.

The independent power battery allows the operation of the Experience motorized supports without the use of the power cable.

Screw the battery onto the wooden spindle at the foot side for the Experience 500 model or at the head side for the Experience 600 model using the screws provided.

The battery should be positioned near the motor to facilitate access to the charging/power cable.


Photo of the Experience 500 model
View from below.

Once the battery is screwed in, follow this procedure:

- disconnect the power cord of the motorized support from the power supply.
- disconnect the power supply from the motor.
- connect the battery to the motor.
- connect the power supply to the battery.
- plug the power supply into the socket and charge the battery.
- once charging is complete, if desired, disconnect the power supply and operate the device freely.

13.1 Instructions for correct use of the Experience battery power kit

The battery should be charged for at least 8 hours with the provided power cable during its first use.

The battery charge lasts about 30 days with one daily usage cycle.

If the battery has sufficient charge, it will continue to operate the motorized support, including the courtesy light, even in case of a blackout.

14 Maintenance and cleaning of the device

The cleaning of surfaces should be done with a dry cloth suitable for dust removal. No specific cleaning products are required for the elements of the Anatomical Support.

In case of extreme necessity (e.g., removal of stains, grease, liquid spills, etc.), use a damp cloth or lightly dampened with mild detergent and immediately dry it.

For any cleaning done under the bed, do not use cleaning tools that cause splashing of water or other liquids.

For motorized Anatomical Supports, cleaning should always be performed with the joints in the flat position and the power turned off.

For motorized Anatomical Supports, if the auxiliary battery has been unused for an extended period, the user should take note of the expiration date on the battery and possibly replace it.

If the labels on the Anatomical Support are so worn that they are no longer visible, request an immediate replacement from the manufacturer.

For sanitizing the product, it is recommended to use products suitable for this purpose.

For motorized Anatomical Supports, it is advisable to periodically check the integrity of all electrical cables.

If the user notices any operational issues, component wear, malfunctions, or breakages, they should request the intervention of authorized personnel, who will handle the repair.

No repairs are allowed by the user for any reason, and it is strictly forbidden to remove or open the motor or the control panel of the motorized version.

Repairs must only be performed by specialized personnel using original spare parts. Interventions made by unauthorized personnel or the use of non-original spare parts voids product compliance (and therefore any warranty on the product).

14.1 Periodic maintenance

Component	Type of check	Frequency	Action
Battery	Expiration date	Every 6 months	Replacement
Handset	Check that each movement corresponds to the expected movement	Every 3 months	Contact support
Motor	Fully raise and lower the head-section and foot-section	Every 3 months	Contact support

14.2 Spare parts

The installation of spare parts is only permitted if performed by specialized personnel.

For the procurement of spare parts, please contact only authorized Dorsal® dealers. For information on the availability of the nearest authorized personnel, please visit the website www.dorsal.it in the "Sales points" section.

14.3 Storage instructions

If not in use, store the product protected by nylon or cardboard in a dry environment, away from light. As these are not fireproof products, the Anatomical Supports should not be placed near ignition sources or heat sources exceeding 100°C.

Also, refer to the information in the paragraph 2 and 6.1 sections.

14.4 Disposal

Nylon, cardboard, and other packaging materials should be kept out of reach of children and sent for recycling. The discarded product must be disposed of through the public collection and recycling system according to the regulations in the user's area, or at an authorized collection centre for the recycling of electrical and electronic waste (WEEE*).

Improper handling of this type of waste can have negative consequences for the environment and human health due to potentially harmful substances commonly found in such waste.

The user's cooperation for the proper disposal of this product will contribute to effective use of natural resources and help avoid administrative penalties under Article 38 and subsequent provisions of Legislative Decree No. 49/2014. For more information on recycling this product, please contact local authorities, the waste collection agency, authorized personnel, or the domestic waste collection service.

* For example (WEEE in mainly in the EU). In the United States, the e-Waste Recycling Act and similar state laws govern the disposal of electronic waste. In China, regulations like the Solid Waste Pollution Prevention Law cover the recycling and disposal of electronic waste. In Japan has the Home Appliance Recycling Law and other related regulations that promote the recycling of electronics. In India, the E-Waste (Management) Rules govern the disposal and recycling of e-waste.

15 Labelling and symbols

	Risk of crushing.
	Do not sit on the head-section and foot-section when not in the flat position.
	Pay attention.
	This symbol on your device indicates that it complies with the requirements of the EU Medical Device Directives 2017/745.
	Indicates the medical device.
	Mandatory behavior.
	Manufacturer information.
	Year of manufacture.
	Obligation to read the manual.
	Indicates the model reference.
	Indicates the website where users can obtain more information about medical products.

	Place in a dry location.
	Temperature range for storage and transport of the product.
	Not sold in sterile condition.
	Do not use if the packaging is damaged.
	WEEE Symbol (Waste Electrical and Electronic Equipment). Recycling symbol. In particular, for the European Union and the European Economic Area, the presence of this symbol on a component indicates that the item cannot be disposed of with household waste, according to the WEEE Directive (2012/19/EC), and/or national laws implementing these directives.

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The label above is located on the head-section side of the Anatomical Support in the central position.

16 Declaration of conformity

EU Declaration of Conformity

The manufacturer: Dorsal Srl,
Via Madonna di Loreto, 3
31020 Corbanese (TV), ITALY

SRN: Not available

Declares, under its own and exclusive responsibility, that the device(s):

CODE	MODEL	ID BD/RDM	Basic UDI-DI
SUPLE	ERGO FA	2215803	805201108SUPLE7V
SUPLE	ERGO NA	2215802	
SUPLE	FOREST FA	2215805	
SUPLE	FOREST NA	2215804	
SUPLE	HUMAN ACTIVE FA	2215817	
SUPLE	HUMAN ACTIVE NA	2215798	
SUPLE	HUMAN BALANCE FA	2215818	
SUPLE	HUMAN BALANCE NA	2215799	
SUPLE	HUMAN FAMILY ACTIVE/STRONG FA	2215844	
SUPLE	HUMAN FAMILY BALANCE/ACTIVE FA	2215837	
SUPLE	HUMAN FAMILY BALANCE/STRONG FA	2215831	
SUPLE	HUMAN STRONG FA	2215801	
SUPLE	HUMAN STRONG NA	2215800	
SUPLE	SANA FA	2215857	
SUPLE	SANA NA	2215852	
SUPLE	XPOINT 200 FA	2215112	
SUPLE	XPOINT 200 NA	2215097	
SUPLE	XPOINT 300 FA	2215101	
SUPLE	XPOINT 300 NA	2215098	
SUPLE	XPOINT 400 FA	2215102	
SUPLE	XPOINT 400 NA	2215099	
SUPAC	SUPER FA	2215807	805201108SUPAC6Q
SUPAC	SUPER FAMILY SUPER/FLEXI FA	2265538	
SUPAC	SUPER FAMILY SUPER/TONIC FA	2265535	
SUPAC	SUPER FAMILY TONIC/FLEXI FA	2265537	

ID.10301

follows...

SUPAC	SUPER FLEXI FA	2265545	805201108SUPAC6Q
SUPAC	SUPER FLEXI NA	2265543	
SUPAC	SUPER NA	2215806	
SUPAC	SUPER TONIC FA	2265532	
SUPAC	SUPER TONIC NA	2265530	
SUPAC	SUPEREGO FA	2215808	
SUPAC	SUPEREGO FIRM FA	2215811	
SUPAC	SUPEREGO FIRM NA	2215810	
SUPAC	SUPEREGO FXA	2215809	
SUPAC	SUPEREGO NA	2215816	
SUPAC	SUPERNOVA FA	2215813	
SUPAC	SUPERNOVA NA	2215812	
SUPAC	SUPERNOVA PLUS FA	2215815	
SUPAC	SUPERNOVA PLUS NA	2215814	

CND/EMDN: V080699 - LETTI MEDICALI – ALTRI

Intended use:

The Anatomical Supports are new devices, designed to improve sleep quality. They are classed as medical devices, and have been designed to be used to:

- improve joint comfort and reduce rheumatic pain;
- alleviate insomnia and its associated physical, psychological and emotional effects;
- improve posture, in particular for the spinal column, with associated alleviation of symptomatic pain;
- reduce thermal and mechanical effects on the skin.

These devices are not suitable for use by individuals suffering from problems associated with the musculo-skeletal, vascular or nervous (both central and peripheral) systems, who need to use specific bed bases of a different nature.

Risk class:

Class I (in accordance with Rule 1, Annex VIII of Regulation (EU) 2017/745)

They complies with the following Union legislative acts:

- Regulation (EU) 2017/745 of the European Parliament and of the Council, of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC of the Council and amended by the Regulation (EU) 2020/561 of 23 April 2020, 2023/607 of 15 March 2023 and 2024/1860 of 13 June 2024

Complies with the following technical/harmonised standards and/or common specifications:

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
The devices are associated with the evaluation procedure provided for by article 52, point 7 of Regulation (EU) 2017/745

Corbanese (TV),
24/02/2025

Signature of legal representative

Armando Moretti

ID.10301

EU Declaration of Conformity

The manufacturer: Dorsal Srl,
Via Madonna di Loreto, 3
31020 Corbanese (TV), ITALY

SRN: Not available

Declares, under its own and exclusive responsibility, that the device(s):

CODE	MODEL	ID BD/RDM	Basic UDI-DI
SUPLE	ERGO RA	2215077	805201108SUPLE7V
SUPLE	ERGO RA BT	2215094	
SUPLE	EXPERIENCE 500 RT 4 SEZIONI	2750002	
SUPLE	EXPERIENCE 500 RT 5 SEZIONI	2750003	
SUPLE	EXPERIENCE 500 FRT 5 SEZIONI	2750800	
SUPLE	EXPERIENCE 600 FRT 6 SEZIONI	2753888	
SUPLE	EXPERIENCE 600 RT 4 SEZIONI	2753960	
SUPLE	EXPERIENCE 600 RT 6 SEZIONI	2753963	
SUPLE	FOREST RA	2215078	
SUPLE	HUMAN ACTIVE FRA	2215070	
SUPLE	HUMAN ACTIVE FRA BT	2215089	
SUPLE	HUMAN ACTIVE RA	2215068	
SUPLE	HUMAN ACTIVE RA BT	2215088	
SUPLE	HUMAN ACTIVE RQUATTRO	2215069	
SUPLE	HUMAN BALANCE FRA	2215073	
SUPLE	HUMAN BALANCE FRA BT	2215091	
SUPLE	HUMAN BALANCE RA	2215071	
SUPLE	HUMAN BALANCE RA BT	2215090	
SUPLE	HUMAN BALANCE RQUATTRO	2215072	
SUPLE	HUMAN FAMILY ACTIVE/STRONG FRA	2215843	
SUPLE	HUMAN FAMILY ACTIVE/STRONG FRA BT	2215840	
SUPLE	HUMAN FAMILY BALANCE/ACTIVE FRA	2215835	
SUPLE	HUMAN FAMILY BALANCE/ACTIVE FRA BT	2215834	
SUPLE	HUMAN FAMILY BALANCE/STRONG FRA	2215823	
SUPLE	HUMAN FAMILY BALANCE/STRONG FRA BT	2215821	

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follows...

Intended use:

The Anatomical Supports are new devices, designed to improve sleep quality. They are classed as medical devices, and have been designed to be used to:

- improve joint comfort and reduce rheumatic pain;
- alleviate insomnia and its associated physical, psychological and emotional effects;
- improve posture, in particular for the spinal column, with associated alleviation of symptomatic pain;
- reduce thermal and mechanical effects on the skin.

These devices are not suitable for use by individuals suffering from problems associated with the musculo-skeletal, vascular or nervous (both central and peripheral) systems, who need to use specific bed bases of a different nature.

Risk class:

Class I (in accordance with Rule 13, Annex VIII of Regulation (EU) 2017/745)

It complies with the following Union legislative acts:

- Regulation (EU) 2017/745 of the European Parliament and of the Council, of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC of the Council and amended by the Regulation (EU) 2020/561 of 23 April 2020, 2023/607 of 15 March 2023 and 2024/1860 of 13 June 2024
- Directive 2006/42/EC of the European Parliament and of the Council, of 17 May 2006 on machinery, and amending Directive 95/16/EC
- Directive 2014/30/EU of the European Parliament and of the Council, of 17 May 2006 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility
- Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment with modification by the delegated Directive (EU) 2015/863 of 31 March 2015

Complies with the following technical/harmonised standards and/or common specifications:

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The device is associated with the evaluation procedure provided for by article 52, point 7 of Regulation (EU) 2017/745

Person authorised to compile the technical file in compliance with Directive 2006/42/EC: Dorsal s.r.l.

Corbanese (TV),
24/02/2025

Signature of legal representative

Armando Moretti

ID.10301

00000000425/03.25

www.dorsal.it **info@dorsal.it**

Dorsal Srl Via Madonna di Loreto, 3 - 31020 Corbanese (TV) Italy
Telefono +39 0438 933003 Fax +39 0438 933023