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Classe I / Class I

Conformemente a / In accordance with

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Firma del legale rappresentante
 Armando Moretti

This manual is a property of Dorsal®

The reproduction, disclosure to third parties and translation into any language, transfer and copying, in any way, even partial, including any magnetic or typed media, without the consent of Dorsal® is prohibited.

Retention of the instruction manual

Keep the instruction manual in a place suitable for its retention and in any case in the place where the mattress is being used. The document must be kept for the entire life of the mattress thereby constituting an appropriate reference for all possible needs.

Zefiro Contour, Therapeutic Contour and Vitalay Azzurro Pillows are distributed by:

Dorsal Srl
via Madonna di Loreto, 3
31020 Corbanese (TV) – Italy
C.F e P.iva 00317220267

Requests for information and technical assistance

To request the intervention of authorised personnel, for any further information on the mattress or on anything not already mentioned in this manual, contact Dorsal® using the following references:

Tel. +39 0438 933003

info@dorsal.it

www.dorsal.it

Index

1	Terms and definitions	3
1.1	Glossary	3
2	Technical characteristics.....	4
2.1	Composition	4
2.2	Conditions of use and storage of the device.....	5
3	Warranty conditions	5
4	Intended use	5
5	Classification of the device	6
6	Warnings and contraindications.....	6
6.1	Mandatory behaviour.....	6
6.2	Warnings before use	7
7	Principles of action	8
8	Equipment.....	9
9	Placing in service of the device.....	9
9.1	Installation.....	9
10	Maintenance and cleaning.....	10
10.1	Washing of the external covering	11
10.2	Repairs and spare parts.....	11
10.3	Storage mode	11
10.4	Disposal	12
11	Labelling and symbols	12
11.1	Device	12
11.2	Covering.....	13
12	Declaration of conformity	14

1 Terms and definitions

1.1 Glossary

Pillow:	Refers to all pillows variants. Any references to individual products will be specifically indicated.
Talalay Latex:	Talalay is latex par excellence: a blend of natural ingredients crafted using an exclusive and complex manufacturing process. The result is a special structure of very open and homogeneous round cell that lends the material its distinctive elasticity and superb breathability. Nonallergenic and antibacterial.
Internal layer:	This refers to the interior part of the pillows.
Tubes®:	A complex system of highly-resilient tubular elements that form a network of ventilation channels, designed to quickly remove humidity. Made with raw materials from renewable resources.
Natur Memory:	Viscoelastic made with raw materials from renewable resources. Reacting to heat it adapts to the shape of the body. Thanks to its open cell structure, it provides enhanced breathability. It's hypoallergenic, antibacterial and odourless.
Soft Care:	Balanced-pressure padding, made with raw materials from renewable resources. The open cell structure makes Soft Care extremely breathable, hypoallergenic, antibacterial, odourless and completely washable in water. Thanks to its softness and adaptability it provides a unique level of comfort.

2 Technical characteristics

Variant	Width (cm)	Lenght(cm)	Height (cm)	Weight (kg)
Zefiro Contour	70	38	9/11	1,5
Therapeutic Contour	70	40	8/10	0,9
	70	40	10/12	1,3
Vitalay Azzurro	70	40	8/10	0,8
	70	40	10/12	1,1

2.1 Composition

Zefiro Contour


Foam based on vegetable oils and water printed with Tubes technology.

Therapeutic Contour

Thermosensitive foam based on vegetable water, with elastic memory.

Vitalay Azzurro

Latex produced with the Talalay method.

2.2 Conditions of use and storage of the device

Operating temperature.....from -30°C to 50°C

Maximum relative humidity.....max 65%

Storage and transport temperature.....from -40°C to 70°C

The product is not fireproof.

The product is not impermeable to liquids.

3 Warranty conditions

The **Zefiro Contour**, **Therapeutic Contour** and **Vitalay Azzurro Pillows** are guaranteed for 2 years. Any complaints can be reported to Dorsal® in the ways and under the terms defined in the Consumer Code Legislative Decree 21 of 21 February 2014 (ITA) and subsequent amendments and additions, in the Directive 2001/95/EC (EU) concerning the general safety of products placed on the market and in the 1999 Directive/34/EC (EU) on liability for damage caused by defective products. For any other details carefully read the warranty certificate attached to the product.

4 Intended use

The function of the pillow is to support the head and the head-neck segment of the user during the rest, and together with the mattress and anatomical support create an ideal bed system.

The **Zefiro Contour**, **Therapeutic Contour** and **Talalay Azzurro Pillows** are new devices, designed to improve sleep quality.

They are configured as medical devices with the following intended uses:

- Improvement of joint comfort with reduction of rheumatic pains;
- Reduction of insomnia and its physical and psycho-emotional effects;
- Improvement of the postural structure in particular of the vertebral column with reduction of the associated painful symptoms;
- Reduction of thermal and mechanical effects on the skin;

The use of these devices is not intended for those people suffering from problems with the musculoskeletal, vascular or nervous system (central or peripheral) who are used specifically for use of another nature.

5 Classification of the device

The **Zefiro Contour**, **Therapeutic Contour** and **Vitalay Azzurro Pillows** have been classified as:

- Class I Medical Devices in accordance with rule 1 according to Annex VIII

In compliance:

- To regulation EU 2017/745.

Articles and registration number at the Ministry of Health:

DORSAL PILLOWS CLASS I MEDICAL DEVICES	REGISTRATION NUMBER TO THE MINISTRY OF HEALTH	CND CLASSIFICATION
ZEFIRO CERVICALE	2215888	V080399 - MEDICAL PILLOWS - OTHERS
THERAPEUTIC CERVICALE	2215891	V080399 - MEDICAL PILLOWS - OTHERS
VITALAY AZZURRO	2215895	V080399 - MEDICAL PILLOWS - OTHERS

6 Warnings and contraindications

6.1 Mandatory behaviour

These instructions for use contain essential information that must be observed during installation, use and maintenance. For these reasons these instructions must be read carefully by the user. The instructions described here must be observed in addition to general home safety precautions.

To maintain the maximum level of safety the user must use the pillow in compliance with the prescriptions and limitations of use of the user manual.

Failure to comply with the safety precautions will void any warranty claims and may result in personal injury or other risks not described in detail in this instruction manual.

The manufacturer declines all responsibility for any use other than what is indicated and prescribed in this manual.

In order to guarantee the safety of the user and to safeguard the product's durability, it is essential to comply with the information provided by the manufacturer.

The pillow must be used in normal positions for correct rest. It must not be strained or used for improper uses.

The pillow is not fireproof and must not be placed near sources of ignition (such as naked flames, cigarettes, flammable anaesthetics, etc.) or near heat sources above 100°C (such as electric or gas heaters).

The use of electric powered thermal blankets is prohibited.

The pillow is not suitable for outdoor use, it should not be placed in damp places and should not be exposed to sunlight.

The from using anatomical supports with adjustable sections without bed base cover.

6.2 Warnings before use

Handling of the product (packaged or not) must be performed in conditions of safety and with suitable gripping methods to avoid damaging the product. Once released from the packaging, to move or position the pillow, always use the two-handed grip.

In addition to the pillow, the optimal bed system involves the use of:

- An anatomical support;
- A mattress.

To guarantee the effectiveness of the pillow, the use of an anatomical support and mattress Dorsal® is recommended. However, the use of an anatomical support and

a mattress with high capacity for ergonomic adaptation and distribution of pressure loads, preferably of the Class I Medical Device type, is in any case advisable.

The pillow is permeable to water and body fluids: the presence of wet spots on its surface can result in the potential proliferation of pathogenic micro-organisms. It is therefore absolutely necessary to protect it with pillow cases, which are permeable to air, impermeable to liquids and must be of soft non-synthetic material, clean and adherent in order to minimise friction with the skin of persons.

The pillow must be placed on a mattress of adequate length. The maximum width of the mattress must be greater than that of the pillow.

7 Principles of action

The bed has the function of allowing optimal rest of the body facilitating the recovery of physical and mental functions. If we consider that the sleep requirement of an average life is approximately 192,000 hours equal to around 22 years and that sleep disorders involve the entire psycho-physical/relational sphere of the individual, the importance of the quality of the environment in which we rest is evident. The lack of sleep is related to the decreased speed and correctness in the response to external stimuli therefore making driving and work activity potentially dangerous, as well as negatively interfering with social relationships and mood.

A study conducted in the USA notes that around 50% of adults suffer from sleep disorders. The reasons for the great interest and the huge investments in the evaluation and the study of increasingly sophisticated "systems" for the qualitative improvement of sleep are therefore easily understandable. We must then consider that the bed becomes the "hospitalisation" of the person when they are sick, and it is precisely in these cases that the advantages of correct "support" of the body are amplified.

In Western tradition the "bed" refers to the following components

- "Bed Base" or anatomical support;
- Mattress;
- Pillow;
- Sheets and blankets.

The pillow supports the entire structure of your head and spine. It can be made from different materials, ranging from natural ones such as down or Talalay latex, to materials such as breathable viscoelastic or breathable foam based on natural oils.

Being in direct contact with your body, your pillow plays a major role in ensuring the overall quality and effectiveness of your bed system, and is capable of regulating the distribution of pressure points, soft tissue tension, and maintaining the correct position of cervical and dorsal sections the spine while you sleep. For this reason, the contribution an anatomical pillow makes can be considered fundamental in determining both sleep quality, preventing disorders of the musculoskeletal system and secondary complications which may occur when lying down.

The Dorsal® **Zefiro Contour**, **Therapeutic Contour** and **Vitalay Azzurro** Pillows are the result of our close partnership with the Italian Association of Physiotherapists, which aims to GUARANTEE clinical efficacy by determining the optimum mechanical and ergonomic performance necessary to ensure they fulfill their intended purpose.

8 Equipment

- Nr. 1 Pillow
- Nr. 1 User manual and 1 guarantee sheet

9 Placing in service of the device

9.1 Installation

The pillow can be installed by a specialist dealer or directly by the user who must follow the instructions below.

In case of uncertainty concerning the methods of installation, the completeness of the equipment or the integrity of the product after unpacking, contact the Dorsal® technical personnel.

- If The pillow is protected by a special bag. Once the pillow has been unpacked, depositing the packaging in a way that does not interfere with the installation operations.
- Do not leave packaging items within reach of children and do not dispose of them in the environment.
- The product can be used immediately without any prior washing.
- Carefully read the manual containing installation, use and maintenance instructions for the product.
- Keep the manual in a suitable place and in any case in the place where the pillow is installed and being used.
- Place the pillow on the mattress, taking care to handle it with a two-handed grip and making sure that the flat part of the pillow leaning to the mattress.
- Search the curvature of the pillow that is most suitable for the user's cervical structure.
- The pillow must be placed on a mattress of the correct width, not smaller than the pillow width.

10 Maintenance and cleaning

Maintenance can be carried out by the operator by following the instructions given by the manufacturer using the symbols on the product and the instructions contained in this manual.

It is strictly forbidden to wash the inner layer of the pillow, if it must be washed, contact the manufacturer. For any sterilisation of the product it is advisable to contact exclusively a firm specialised in the sanitisation of products of the health/hospital sector.

Remove any dust that is deposited on the surface of the pillow every week, using a soft brush.

In case of localised stains never use chemicals or solvents and do not use steam machines, never use the vacuum cleaner or beat the pillow violently.

10.1 Washing of the external covering

- The removable covering can be washed under the following conditions:
 - in water and with the limitations shown in the labelling sewn on the inside of the covering.
 - dry cleaning at a specialist laundry (not at self-service laundries).
- In the specific case of washing in water:
 - the covering must only be washed by itself, avoiding combining it with other linen.
 - spinning suitable for delicate items must be carried out.
 - ironing is not permitted.
 - The adoption of washing temperatures lower than the maximum limit indicated by the manufacturer or the preference for dry cleaning with respect to washing in water will ensure cleaning of the covering and less wear of the fabrics.

Washing of the external covering is recommended every 6 months of continuous use.

If, due to repeated washing, the information contained in the label sewn to the covering should not be legible, refer to the symbols shown in the Labelling and symbols section before proceeding with new washes.

10.2 Repairs and spare parts

If the users finds anomalies with the product, they must request only the intervention of a specialist dealer who will examine the problem. For the supply of spare parts, contact an authorised Dorsal® products dealer only. For information on the nearest authorised dealer, contact Dorsal® directly.

10.3 Storage mode

If not being use, store the pillow inside its original packaging, in a horizontal position, on a flat surface, in a dry and protected from light.

10.4 Disposal

The discarded product must not be dispersed in the environment but must be given to a public waste collection and recycling system in accordance with the procedures and regulations in force in the relevant municipality or Country.



11 Labelling and symbols

11.1 Device

	This symbol on your device indicates that it complies with the requirements of the Medical Device Directives EU 2017/745.
	It indicates the medical device
	It informs the user of mandatory behaviour
	Manufacturer identification data
	Year of manufacture
	It informs on the obligation to read the manual
	It indicates the model reference
	It indicates the website from which patients can obtain further information on medical products
	Place in a dry place, protected from moisture
	Temperature required for - storage and transportation of the product (shown on the packaging) - the use (shown on the product)
	Humidity envisaged for - storage and transportation of the product (shown on the packaging) - the use (shown on the product)
	Keep away from sunlight

	Not sold in sterile packaging
	Do not use if the packaging is open or damaged
	Presence of natural latex in the device (see the technical characteristics for the relevant models)

11.2 Covering

	Wash with a special gentle cycle at a maximum temperature of 30°C. • reduced mechanical action. • rinse at a gradually decreasing temperature. • reduced centrifugation.
	Wash with a special gentle cycle at a maximum temperature of 40°C. • reduced mechanical action. • rinse at gradually decreasing temperature. • reduced centrifugation.
	Wash with a special gentle cycle at a maximum temperature of 60°C. • reduced mechanical action. • rinse at gradually decreasing temperature. • reduced centrifugation.
	The fabric cannot be treated with bleaches.
	Do not iron.
	Do not use a dryer.
	The garment can be cleaned with perchlorethylene, hydrocarbons or solvents R113 and R11, respecting the temperature and humidity limits. No added water or stressful mechanical pressures.
	Hand wash.

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under the sole responsibility issue this

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declaration of conformity at Regulation (EU) 2017/745

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Dorsal Srl, Via Madonna di Loreto 3 31020 Corbanese (TV), ITALIA
Dorsal Srl Via Madonna di Loreto, 3 - 31020 Corbanese (TV) Italy