Natural Cycles

User Manual / Instructions for Use

Please read this manual before you start using Natural Cycles.

Intended Use/Intended Purpose

Natural Cycles is a stand-alone software application, intended for women 18 years and older, to monitor their fertility. Natural Cycles can be used for preventing a pregnancy (contraception) or planning a pregnancy (conception).

Intended User

The intended target group is women of reproductive age. The pediatric population <18 years of age is excluded from the target group. The user shall be capable of reading this User Manual and operating the Natural Cycles application.

Intended Use Environment

The device is intended for use in non-medical settings.

Contraindications

There are no known contraindications for use of Natural Cycles.

Side-Effects

The product has no known direct side effects on the human body.

Natural Cycles may not be right for you if

- You have a medical condition where pregnancy would be associated with a significant risk to the mother or the

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fetus. Using Natural Cycles does not guarantee 100% that you will not get pregnant. If pregnancy presents a significant risk, you should talk to your doctor about the best contraceptive option for you.

- You are currently taking hormonal birth control or undergoing hormonal treatment that inhibits ovulation. Natural Cycles will mainly provide red days if you do not ovulate. You can explore the app in Demo Mode until you discontinue your hormonal birth control or treatment.

Warnings

1. No method of contraception is 100% effective. Natural Cycles does not guarantee that you cannot get pregnant, it only informs you whether and when you can become pregnant based on the information that you enter into the application.
2. Even with using the app perfectly, you can still have an unintended pregnancy. Natural Cycles is 93% effective under typical use, which means that 7 women out of 100 get pregnant during 1 year of use. With using the app perfectly, i.e. if you never have unprotected intercourse on red days, Natural Cycles is 98% effective, which means that 2 women out of 100 get pregnant during 1 year of use.
3. On red days, you must either abstain or use protection, such as condoms or another form of barrier protection, to prevent pregnancy. If you have not been using a non-hormonal method of contraception before, make sure to have a method on hand to use during red days.
4. Our data show that women transitioning to Natural Cycles from hormonal contraception have a higher risk of becoming pregnant due to generally not yet being accustomed to using other forms of protection, such as a condom or another form of barrier protection.
5. If you are discontinuing hormonal contraception, you may first experience a withdrawal bleed, which is much lighter than your true period. You will normally get your true period a few weeks later. Withdrawal bleeding

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should not be entered as “period”. If you are unsure if you are experiencing withdrawal bleeding or a period, please consult your doctor.

6. Natural Cycles does not protect against sexually transmitted infections (STIs). Use a condom to protect against STIs.
7. Always check your fertility status for the day. Be aware that fertility predictions for upcoming days are only predictions and may change in the future. If you are unsure of your fertility status and can't access the app, abstain or use protection in the meantime.
8. If you are exploring the demo version of the app, this is for informational or instructional purposes only and does not display real data. Do not use the information displayed in Demo Mode for contraceptive purposes.
9. If you have taken the emergency contraceptive pill, you should indicate it in the app. You can then continue measuring as normal and the algorithm will take the additional hormones from the pill into account. Note that the bleeding you may experience a few days after taking emergency contraception is not the same as your regular period and should not be entered into the app. If your period is more than one week late after it normally would start you should take a pregnancy test.

Precautions

- As a contraceptive, Natural Cycles may be less suitable for you if you have irregular menstrual cycles (i.e. cycles with length less than 21 days or greater than 35 days) and/or fluctuating temperatures as predicting fertility is more difficult in these circumstances. Natural Cycles' effectiveness as a contraceptive is not affected by irregular cycles, but you will experience an increased number of red days, which can reduce your satisfaction with Natural Cycles. This can be the case:
 - **If you discontinue hormonal contraception.** It can take several cycles after discontinuing hormonal contraception before becoming more regular again. You should expect more red days during these

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cycles and therefore expect to use protection or abstain more frequently.

- **If you have medical conditions that lead to irregular cycles such as Polyendocrine Metabolic Ovarian Syndrome (PMOS) or thyroid-related conditions.**
- **If you experience symptoms of menopause.**
- After giving birth, we recommend that you start using Natural Cycles again when you get your first period back, regardless of whether you are breastfeeding or not. If you wish to start using Natural Cycles before, you should expect to have only red days until your first ovulation is confirmed. How long it may take to get your first period back can vary significantly from person to person, and it is generally longer for nursing mothers. Please note that you are at risk of becoming pregnant even before your first period comes back, and you should use some form of contraception if you wish to prevent a new pregnancy. If you are fully breastfeeding and you have given birth less than 6 months ago, LAM (Lactational Amenorrhea Method) can serve as an effective method of contraception postpartum.
- When you start using Natural Cycles, it can take some time to get to know your unique cycle, resulting in an increased number of red days during your first 1-3 cycles. This is to ensure that the app is effective from the first day of use, and once Natural Cycles has detected your ovulation you are likely to receive more green days. It is important to note that you can rely on the fertility status indicated by the app from the first day of use.
- Always make sure to update your software app, so that you have the latest version. Natural Cycles recommends that you set your smartphone to update the app automatically.
- Always make sure to update your software app, so that you have the latest version. Using the latest version helps keep your app safe from cybersecurity threats and vulnerabilities. Your app is updated via your app store. Natural Cycles recommends setting your smartphone to update the app automatically.
- Always follow the directions in 'Natural Cycles and Measuring Device' to increase the accuracy of ovulation detection.

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Checklist

To get started with Natural Cycles you will need the following:

1. A Natural Cycles account

For access via the app, download the app on your Android, iPhone or iPad on Google Play or the App Store. Follow the registration steps to create your personal account. Make sure the device you use to access Natural Cycles is connected to the internet. Your Natural Cycles account is personal and should not be used by anyone except you.

2. Measuring Device

To get started, you need a compatible wearable device or an oral basal thermometer. Compatible wearable devices are: NC° Band, an Oura Ring (Gen2, Gen3, Gen4, or Gen 5), an Apple Watch with temperature sensors (Series 8, Series 9, Series 10, Series 11, Ultra, Ultra 2, Ultra 3, or SE3), a Garmin wearable (Enduro™ 3, fēnix® 8, fēnix® 8 Pro, tactix® 8, quantix® 8, quantix® 8 Pro, Forerunner® 570, Forerunner® 970, Venu® 3, Venu® 3S, Venu® 4, Venu® X1, CIRQA Band), or a WHOOP (MG or 5.0). If you choose to measure with a thermometer, you must use a basal thermometer, which is more sensitive than regular fever thermometers and shows two decimals (for example, 36.72°C / 98.11°F).

Note: Do not use the wearable temperature data shown in the app to make medical decisions (for example, diagnosing fever).

Note: If you have a condition that affects nocturnal skin temperature, using a wearable measuring device as an input source to Natural Cycles is not recommended.

Note: If you need to change your measuring device for any reason, you should start using your new device at the start of your period.

3. Ovulation tests (optional)

Ovulation (LH) tests are urine test strips which detect the surge of luteinizing hormone that occurs 1-2 days before ovulation. It is optional to use LH tests with Natural Cycles,

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but they can help to increase the accuracy of detecting ovulation and can, therefore, increase the number of green days you will see if you are preventing a pregnancy. Using LH tests will not affect the effectiveness of Natural Cycles as a contraceptive. If you are planning a pregnancy, we recommend using LH tests as they can help to find your most fertile days. You can purchase LH tests at shop.naturalcycles.com or from within the app itself.

Natural Cycles and Measuring Device

How to measure

Follow the user manual provided with your measuring device to either measure overnight (for wearables) or in the morning (for basal body thermometers).

When to measure - wearable device

If you are using a wearable device for overnight measuring, your temperature syncs according to the wearable device manufacturer's instructions. There may be steps required to sync your temperature from the wearable to the NC° app – see in-app messaging for more details. Aim to measure at least five days a week.

When to measure - basal body thermometer

Measure your temperature once a day, on as many days as you can, so the algorithm can quickly learn your unique cycle. Refer to your thermometer manufacturer's instructions for information on how to use your thermometer. Aim to measure at least five days a week. One tip to help you remember is to place your thermometer on top of your phone before going to bed.

When to exclude temperature data

Feel hungover / Feel sick / Slept differently (>2 hours more or less than usual – only applicable if using a thermometer)
Medication can also affect your temperature. If you take regular medication or antipyretics (fever-reducing substances), it's worth talking with your doctor to understand how they might impact your temperature.

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
2. Nospiediet pogu tikai tad, kad termometrs ir mutē un

Natural Cycles products

Natural Cycles offers five modes for our users to choose from depending on their fertility goal:

- NC° Birth Control - for preventing pregnancy
- NC° Plan Pregnancy - for increasing chances of conception
- NC° Follow Pregnancy - for monitoring the course of a pregnancy
- NC° Postpartum - for supporting recovery after childbirth
- NC° Perimenopause - for decoding and navigating the menopausal shift

These five modes can be chosen by the user prior to signing up and then within the app, if the user would like to switch between modes.

Contraception: NC° Birth Control

How it works

Your daily morning routine with Natural Cycles is as follows:

- Measure or sync your data when you wake up, as described in the 'Natural Cycles and Measuring Device' section.
- Add or sync your temperature and whether you have your period in the app. Tap the temperature field on the Today tab to add, exclude, or edit a temperature (for thermometer users only). To log a period, LH test, sex, or other trackers, tap the plus button ("+" icon). Once you add or update your data, the app will recalculate your status for the day.
- The app will show predicted and confirmed ovulation, as well as predicted period days. The app also shows fertility status predictions for upcoming days and months. Keep in mind that predictions can change and should not be used as the final result, which is why you should check your status for the day – every day.

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Green days indicate you are not fertile and can have sex without protection. Green days are indicated by a green outlined circle with the words "Not fertile" inside.

Red = Use protection

Red days indicate you are likely to be fertile and at risk of becoming pregnant. Use protection (such as condoms or another form of barrier protection) or abstain in order to prevent a pregnancy. Red days are indicated by a red outlined circle with the words "Use protection" inside.

The app will indicate your ovulation as well as when you have entered and are expected to have your period. The application also shows predictions for upcoming days and months. The predictions may change and should not be used as the final result, which is why you should check your status for the day - every day.

Summary of clinical studies - how effective is Natural Cycles for contraception?

Natural Cycles is 93% effective under typical use. Under perfect use, i.e. if you never have unprotected intercourse on red days, Natural Cycles is 98% effective.

Clinical studies have been conducted in order to evaluate the effectiveness of Natural Cycles for contraception. The current version of the algorithm (v.3) has been investigated on 15,570 women (on average 29 years old). At the end of the study, the following commonly used failure rates were determined:

- The app has a method failure rate of 0.6, which is a measurement of how often the app incorrectly displays a green day when the woman is actually fertile and gets pregnant after having unprotected intercourse on this green day. This means that 0.6 out of 100 women who use the app for one year get pregnant due to this type of failure.
- The app has a perfect use Pearl Index of 1, which means that 1 out of 100 women who use the app for

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- Jums ir pagiras

one year and who get pregnant will do so either because:

- They had unprotected intercourse on a green day that was falsely attributed as non-fertile (i.e., method failure); or
- They had protected intercourse on a red day, but the chosen method of contraception failed.
- The app has a typical use Pearl Index of 6.5, which means that in total 6.5 women out of 100 get pregnant during one year of use due to all possible reasons (e.g. falsely attributed green days, having unprotected intercourse on red days, and failure of the contraceptive method used on red days).

Summary of clinical data from 15,570 women on the effectiveness of Natural Cycles of the current algorithm version.

Algorithm Version	v.3
Study Date Range	Sept 2017 - Apr 2018
# women	15,570
Exposure time (Women-years)	7,353
# pregnancies [worst-case]	475 [584]
Typical use PI (95% confidence interval) [worst-case]	6.5 (5.9-7.1) [7.9]
Method failure rate (95% confidence interval)	0.6 (0.4-0.8)

Effectiveness of Natural Cycles for two subgroups - women who used hormonal contraception within 60 days prior to using the app and women who did not use hormonal contraception within 12 months prior to using the app.

Subgroup	Typical Use PI (95% confidence interval)
Recent Hormonal Contraception 3779 women	8.6 % (7.2-10.0)
No Hormonal Contraception	5.0 %

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8412 women

(4.3-5.7)

You can find an overview of the effectiveness of different contraceptive methods in the Appendix of this manual. This will help you to understand the expected effectiveness of all forms of contraception.

Instructions if you are discontinuing hormonal contraception

Type of Contraception	When can you start Natural Cycles?
The Pill (Combined, Medium or Mini)	We recommend that you finish all the active pills of your packet. You can start using Natural Cycles as soon as you finish. You will have a high number of red days to begin with but this will improve over time.
Hormonal IUD (Intrauterine Device)	You can start using Natural Cycles the day after it has been removed.
Contraceptive Implant	You can start using Natural Cycles the day after it has been removed.
Hormonal Patch, Hormonal Contraceptive Ring	You can start using Natural Cycles the day after it has been removed. You will have a high number of red days to begin with but this will improve over time.
Contraceptive Injection	You can start using Natural Cycles when the effect of the injection has worn off. This takes 8-13 weeks depending on which injection you had.

If you are discontinuing hormonal contraception you may first experience a withdrawal bleed, which is much lighter than your true period. You will normally get your true period a few weeks later. Withdrawal bleeding should **not** be entered in the app as “period”.

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Note: a copper IUD does not contain hormones and does therefore not inhibit ovulation. You can thus use Natural Cycles alongside a copper IUD to get insights on your body and menstrual cycle without having to use a barrier method on red days.

If you become pregnant on Natural Cycles

If your period is overdue and your temperature does not drop, you may be pregnant and the app will encourage you to take a pregnancy test to confirm the pregnancy. If the test is positive, please talk to your doctor and indicate it in the app.

NC° Plan Pregnancy

When you are planning a pregnancy, select Natural Cycles° Plan Pregnancy while signing up in the app. If you are already signed into the app, you can go to the settings page and choose NC° Plan Pregnancy. In this mode, your fertility is displayed as a scale, so you can identify the days you are most likely to become pregnant. The same morning routine as in NC° Birth Control mode applies here as well - measure, add or sync data and check your fertility status.

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When you're planning a pregnancy, we recommend taking LH tests since the occurrence of LH coincides with your most fertile days. If your period is overdue and your temperature does not drop, the app will encourage you to take a pregnancy test to confirm the pregnancy.

NC° Follow Pregnancy

Once you are pregnant, the app can be used as a pregnancy calendar and you can follow your weekly development and the development of your baby.

NC° Postpartum

After you have given birth, you have the option to enter NC° Postpartum mode. NC° Postpartum guides you through the physical and mental recovery of childbirth with resources tailored to help support your overall health and well-being during this temporary transitional phase. While taking your temperature is not required, measuring is encouraged as it

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helps our algorithm detect when your fertility returns so you can better prepare. Until the NC° app detects your ovulation, it will give you Red Days. You can exit the NC° Postpartum mode at any time.

NC° Perimenopause

As you approach the menopause transition, you have the option to enter NC° Perimenopause mode. When users enter this mode, the NC° Menopause Algorithm will analyze data to help users discover what stage of perimenopause they are most likely to be in. NC° Perimenopause offers daily check-ins for users to track their symptoms, bleeding, temperatures, and other data so it can provide daily insights based on the user's input. Cycle changes and trends can be seen in this mode and data can be downloaded in PDF form to share with doctors. NC° Perimenopause can be accessed by users regardless of their birth control needs.

Demo Mode

In Demo Mode you can see how the app could look with a few cycles of data, but you cannot add any of your own data. This is for demonstration purposes only, demo mode cannot provide you with any fertility information.

If you are currently on hormonal contraception you will access the app in Demo Mode. Once you stop using hormonal contraception you can then easily access the full version of the app. The app will now function for contraception or for planning a pregnancy.

Problems/ Troubleshooting

If you cannot access the application, the application is malfunctioning, or if you suspect something is wrong with your current fertility status, please follow the steps below:

- Use protection, such as a condom or another form of barrier protection, until the error has been fixed.
- Please contact our service desk via help.naturalcycles.com and our trained support agents will assist you.
- Serious incidents involving this device should be

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reported both to Natural Cycles and to the competent authority in your country. A 'serious incident' means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat;

Application Overview

Your fertility status for the day is shown through different views within the app.

Calendar tab

This tab provides the overview of today's fertility and gives insights into the data logged in the past through various views as described below.

Today view

On this screen, you will get all of the info you need in one glance. A screen to enter your temperature pops up automatically when you open your app in the morning. Once you've entered your data, your fertility status will appear as a colored circle. Below the circle, you will see the weekly predictions and additional information. You can add or update your input data, such as temperature, period, or LH tests, at any time by tapping the '+' icon.

Month view

Tap the "Month" tab in the calendar view to access the monthly view. Here you can see your predicted red and green days for the month, as well as when you are predicted to ovulate and have your period. These may change over time, so you should only rely on the information provided in the "Today" view for contraceptive or pregnancy planning purposes. You can easily access your past data by tapping on any past day.

Graph view

The graph is a visualization of your temperature curve which will continue to develop as you measure throughout your cycle. Information about ovulation, period, and fertile

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days is also visible. Swipe from left to right to view your past cycles, and zoom in and out to compare them to each other. You can also access the compare mode to compare your period and cycle length, temperature variation, ovulation window and day, and tracker graph.

Top menu

Under the top menu, you can access and manage all the information regarding your personal profile, the app and privacy settings, and your account including subscriptions. Here you can also manage whether or not you will get reminders about when to use protection, take an LH test, expect PMS, measure your temperature and check your breasts. Note that you must enable Natural Cycles to send you push notifications in your device settings to receive these.

Offline mode

When you are offline you can continue to add your temperature and other data to the app. When offline, the app does not calculate your fertility status. Instead, the app will show an orange outlined circle indicating your status as "Use protection". Use protection (such as condoms or another form of barrier protection) or abstain until you are back online and the app has calculated today's fertility status.

Platform description

Natural Cycles Application Version: 5.0.0 and onwards

Medical Device version: E


Devices that may be used to access the application:

- A mobile phone running Android OS (version 11 or above) or iOS (version 17 or above), 250 kbps or faster internet.

Accessing the Instructions for Use

The Instructions for Use can be found online on the Natural Cycles website (<https://naturalcycles.app/IFU>) and in your app under the "Regulatory" page.

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If you would like a paper copy of the Instructions for Use we will send you one at no cost. Contact customer support with your request.

Cease using Natural Cycles

When you sign up for Natural Cycles, whether it's for a Monthly or Yearly subscription, you will need to cancel it if you don't want your subscription to be automatically renewed for the next renewal period. Cancellation of your subscription must be done at least 24 hours before the subscription is scheduled to be renewed. Please refer to help.naturalcycles.com for a full description on how to cancel your subscription successfully. Please refer to our Privacy policy for further information on how Natural Cycles processes your data.

Security

Natural Cycles maintains a high level of security and data privacy in the app. The app is continuously monitored for security events and data breaches. If a security event or data breach occurs, it will be addressed promptly according to established processes, and you will receive the information you need for your continued safe use of the app. Despite these measures, some cybersecurity risks remain.

Unauthorized access to your account

If your login credentials are compromised, an unauthorized person may gain access to your account and personal data. To reduce this risk, you should use a strong password and keep your login details confidential.

Phishing and fraudulent communication

Attackers may try to obtain your login details through emails or messages that appear to come from Natural Cycles. Natural Cycles will not ask for your password unless you have started a password reset or support request. You should not share your login details in response to unsolicited messages.

Access from shared or compromised devices

If you use your account on shared or unsecured devices, your data may be accessed by others. You can reduce this risk by logging out after use and ensuring your device is secure. If you suspect unauthorized access, you should change your password immediately and contact support to terminate active sessions.

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Device Name:

Natural Cycles

We declare that this product meets all applicable requirements of the Regulation (EU) 2017/745 for medical devices (MDR) and bears the mark CE0123.

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● No, never	● Ne pēdējo 12 mēnešu laikā ● Nē, nekad
Using Natural Cycles while still on hormonal birth control When you use hormonal birth control (such as the Pill, an implant, or a vaginal ring), you don't ovulate. This means that Natural Cycles won't be able to track your menstrual cycle. While you are using hormonal birth control, we recommend you explore the demo version of the app. ● Use in Demo Mode	Lietojot Natural Cycles, vēl lietojot hormonālo dzimstības kontroli Lietojot hormonālos kontracepcijas līdzekļus (piemēram, tabletes, implantu vai maksts gredzenu), jums nav ovulācijas. Tas nozīmē, ka Natural Cycles nevarēs izsekot jūsu menstruālajam ciklam. Kamēr lietojat hormonālo kontracepciju, iesakām izpētīt lietotnes demonstrācijas versiju. ● izmantot demonstrācijas režīmā
Approximately how long does your cycle usually last? A cycle is the time span between one period and the next.	Cik apmēram parasti ilgst jūsu cikls? Cikls ir laika posms starp vienu un nākamo periodu.
How much does the length of your cycle vary? Understanding the regularity of your cycles helps the algorithm get to know you better. ● 0-4 Days ● 5-9 Days ● 10+ Days ● I don't know	Cik daudz svārstās jūsu cikla ilgums? Izpratne par jūsu ciklu regularitāti palīdz algoritmam jūs labāk iepazīt. ● 0-4 dienas ● 5-9 dienas ● 10+ dienas ● Es nezinu
For how long did you use hormonal contraception? ● 0 - 1 Years	Cik ilgi lietojāt hormonālo kontracepciju? ● 0 - 1 gadi

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Using Natural Cycles while breastfeeding For the first 4–6 months after giving birth, many women who are exclusively breastfeeding may not ovulate. Natural Cycles will detect your first ovulation 2 weeks before your first menstruation, but until then you will only see red days in the app. Even though you are not ovulating, it is still important to use contraception or abstain from intercourse on red days. ● I understand ● Learn more	Natural Cycles lietošana zīdīšanas laikā Pirmos 4–6 mēnešus pēc dzemdībām daudzām sievietēm, kuras baro bērnu ar krūti, var nebūt ovulācijas. Natural Cycles noteiks jūsu pirmo ovulāciju 2 nedēļas pirms pirmajām menstruācijām, bet līdz tam lietotnē redzēsiet tikai sarkanās dienas. Pat ja jums nav ovulācijas, sarkanajās dienās joprojām ir svarīgi lietot kontracepciju vai atturēties no dzimumakta. ● Es saprotu ● Uzziniet vairāk
--	--

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General

- Confirmed ovulation
- Ovulation
- Today
- Cycle day
- Monday
- Tuesday
- Wednesday
- Thursday
- Friday
- Saturday
- Sunday
- Predictions - do not take as result
- Exclude temperature
- Sick
- Slept differently
- Hungover
- Bleeding
- Period
- Spotting
- Sex
- Protected
- Unprotected
- None
- LH test
- Positive
- Negative
- Compare
- Pregnancy test
- No longer pregnant
- More
- Less
- Save
- Done
- Clear
- Skip
- Regulatory
- Version

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		change is not necessary for safety reasons.	
--	--	---	--