








Target Identification Methods in Phenotypic Drug Discovery

Chemical Probe-Based

Label-Free Methods

Indirect Methods

METHOD	STRENGTHS	LIMITATIONS	INDUSTRY APPLICATION	EXAMPLE REFERENCE
Affinity-Based Pulldown (e.g., kinobeads, SILAC, click, iTRAQ)	<ul style="list-style-type: none"> Probes can be reversible <ul style="list-style-type: none"> Relatively high probability of success 	<ul style="list-style-type: none"> Requires probe linker modification Potent reversible probe needed ($K_D = \text{nM}$) Requires cell lysis 	 <p>KYOWA KIRIN, NAMPT</p>	Yamaguchi, D. et al. <i>Sci. Rep.</i> 2019 , <i>9</i> , 7742.
Photoaffinity Labeling (e.g., diazirines, competitive labeling)	<ul style="list-style-type: none"> Irreversible binding, can provide binding site information Can be done in living cells and for insoluble proteins 	<ul style="list-style-type: none"> Requires photoreactive group and purification tag group <ul style="list-style-type: none"> Relatively lower probability of success 	 <p>EISAI, SF3b</p>	Kotake, Y. et al. <i>Nat. Chem. Biol.</i> 2007 , <i>3</i> , 570.
Activity-Based Protein Profiling (e.g., sulfonyl fluorides, beta-lactones)	<ul style="list-style-type: none"> Irreversible binding, can provide binding site information Can be done in living cells and for insoluble proteins 	<ul style="list-style-type: none"> Requires reactive probe for target class Limited to specific classes of enzymes 	 <p>PFIZER, ALDH2+</p>	Huang, Z. et al. <i>ACS Chem. Biol.</i> 2019 , <i>14</i> , 192.
Thermal Stability Profiling (CETSA® Explore / TPP / PISA)	<ul style="list-style-type: none"> No modification of compound needed Reversible molecules can be used Can be done in living cells <ul style="list-style-type: none"> High protein coverage among label-free methods 	<ul style="list-style-type: none"> Does not provide binding site information May have lower sensitivity than probe-based methods Requires soluble target 	 <p>MERCK, DPP9</p>	Moore, K. P. et al. <i>ACS Chem. Biol.</i> 2022 , <i>17</i> , 2595.
Proteolytic or Chemical Stability Profiling (e.g., PP, DARTS, LiP-MS, SPROX)	<ul style="list-style-type: none"> No modification of compound needed Reversible molecules can be used Provides alternative methods of detection 	<ul style="list-style-type: none"> Can theoretically localize binding site to peptide regions Cannot be conducted in live cells May be limited in protein coverage (e.g., SPROX requires methionines) 	 <p>BAYER, FUNGAL CK1</p>	Piazza, I. et al. <i>Nat. Commun.</i> 2020 , <i>11</i> , 4200.
Resistance Screening (e.g. anti-infectives, cytotoxics)	<ul style="list-style-type: none"> No modification of compound needed Reversible molecules can be used Provides relatively strong evidence for target 	<ul style="list-style-type: none"> Only applicable in certain indications (e.g., anti-infectives, cancer) Can be time and biology resource-intensive 	 <p>NOVARTIS, DENV NS4B</p>	Moquin, S. A. et al. <i>Sci. Transl. Med.</i> 2021 , <i>13</i> , issue 579.
Chemogenomic Profiling (e.g. HIP-HOP, CRISPR mutagenesis scanning)	<ul style="list-style-type: none"> No modification of compound needed Reversible molecules can be used Unbiased for target abundance 	<ul style="list-style-type: none"> Only applicable in certain settings (e.g., antifungals, cancer) Can be time and biology resource-intensive 	 <p>NOVARTIS, SEC14P</p>	Pries, V. et al. <i>Cell Chem. Biol.</i> 2018 , <i>25</i> , 279.