



HoLEP/ThuLEP

HoLEP/ThuLEP



Cystoscopes

Part Number	Diameter	View	Length
STS-2040-302BA	4.0 mm	30°	302 mm

26 Fr. Laser Bayonet-Lock Sheaths

Part Number	Description	Size
STS-URO-340-20-SET	Rotatable continuous flow sheath rotatable	26 Fr.
Consisting of:		
STS-URO-209-06	26Fr. outer sheath	
STS-URO-209-26	Rotatable inner sheath with ceramic tip	
STS-URO-320-10	Obturator	

Working Elements

Part Number	Description
STS-URO-047-08	Laser-working element for 0.8mm fiber laser
STS-URO-047-10	Laser-working element for 1.2mm fiber laser

Nephroscopes

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Length
STS-5026-250-90HA	19Fr.	90°	0°	15Fr.	250 mm
STS-5026-250-45HE	19Fr.	45°	15°	15Fr.	250 mm

Laser Adapters

Part Number	Description
STS-NEP-420-25QL	Adapter to use with nephroscope and resectoscope sheath



3020 NW 82nd Avenue
Miami, Florida USA 33122

sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System