




## Mechanical Lithotripsy

# Mechanical Lithotripsy



## Cystoscopes

| Part Number    | Diameter | View | Length |
|----------------|----------|------|--------|
| STS-2040-302BA | 4.0 mm   | 30°  | 302 mm |

## 4mm Cystoscopy Sheaths

| Part Number        | Color Code | Channel sizes    | Size   |
|--------------------|------------|------------------|--------|
| STS-URO-001-25-SET | White      | 1x10Fr. / 2x8Fr. | 25 Fr. |

## Adaptors

| Part Number    | Description                                                                        |
|----------------|------------------------------------------------------------------------------------|
| STS-URO-035-25 | Adaptor for optical instrument in resectoscopy sheaths                             |
| STS-URO-035-40 | Adapter for optical stone punch for use with resectoscope outer sheath             |
| STS-URO-082-90 | Adaptor for stone crushing forceps to be used with 26Fr. resectoscopy outer sheath |

## Resectoscopes

| Part Number         | Description                                   | Size   |
|---------------------|-----------------------------------------------|--------|
| STS-URO-340-20-SET* | Resectoscopy continuous flow rotatable sheath | 26 Fr. |
| Consisting of:      |                                               |        |
| STS-URO-209-06      | Outer sheath                                  |        |
| STS-URO-209-26      | Rotatable inner sheath with ceramic tip       |        |
| STS-URO-320-10      | Obturator                                     |        |

\* 28Fr. monopolar resectoscope available, and upon request, higher flow outer sheath

## Optical Stone Crushing Forceps

| Part Number    | Description                                                                                                    |
|----------------|----------------------------------------------------------------------------------------------------------------|
| STS-URO-082-00 | Optical Stone Crushing Forceps to be used with 25Fr. Sheath or with an adaptor and a 26Fr. Resectoscopy sheath |

# Mechanical Lithotripsy



## Cystoscopes

| Part Number    | Diameter | View | Length |
|----------------|----------|------|--------|
| STS-2040-302BA | 4.0 mm   | 30°  | 302 mm |

## Optical Stone Crushing Forceps

| Part Number    | Description                                                    |
|----------------|----------------------------------------------------------------|
| STS-URO-076-QL | Cystoscopy mechanical lithotripter for hard stones, quick-lock |



3020 NW 82nd Avenue  
Miami, Florida USA 33122

[sales@endoscopia.com](mailto:sales@endoscopia.com) | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399  
**ISO 13485:2016 – Certified Quality Management System**