


See the difference



Arthroscopy
Instrument Product Guide

Index

Linear Line® Instruments	1
Duckbee Basket Punch	3
Ducklee Basket Punch	5
Basket Punch Slim Style	7
Basket Punch Oval	8
Basket Punch Scoop	9
Basket Punch Meniscal Elevator	11
Hook Punch	12
Posterior Basket Punch	13
Mega Biter	14
Basket Punch Unicorn	15
Rotary Basket Punch	16
Scissor Punch	17
Hook Scissors	18
Rotary Hook Scissors	19
Rongeur Raptor	20
Suture Cutter	21
Grasping Forceps	22
Drill Guide System	24
One-Piece Instruments	25
Chondro Pick.....	27
Simple 7 Instrument Set	28
Flexible 15 Instrument Set	29
Knee Suture Passer	30
Shoulder Suture Passer	31
Modulock®	32
Arthroscopes	33
Sheaths Sets	35
Accessories	36




The Strauss Linear Line[®]

When performance is on the line, go Linear.

The Strauss Linear Line of knee arthroscopy instruments is a model of thoughtful design and manufacture. Built for comfort and performance, no hinges or pins ensures high durability without compromising strength or cutting precision. The aerospace-grade aluminum handle reduced overall instrument weight by 60%; enhancing maneuverability while helping to lessen grip fatigue. The entire line features the lowest tip profile of open shaft instruments – with open-shafts that are easier to clean and sterilize.





Straight

The primary style for a direct approach.



Upbiter

An angled tip increases access to the posterior meniscus.



Upswept

An angled shaft follows the contour of the condyle and keeps the cutting surface parallel to the tibial plateau.



Curved left and curved right

Angled to reach the anterior and posterior meniscus.



Classic Handle

Scissors, basket punches and graspers are equipped with traditional classic handle.



Slide-Lock Handle

All graspers are equipped with this infinite-position locking-handle design.



Rotatable Handle

Available on rotary basket punches and scissors.

Duckbee Basket Punch

provides maximum tissue resection with its waved tooth and wide bite.



Duckbee Basket Punch, small



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	2.50mm	4.15mm	STS-ART-10-0001
Upswept shaft	2.52 mm	2.50mm	4.15mm	STS-ART-10-0002



Duckbee Basket Punch, small, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	2.50mm	4.15mm	STS-ART-10-0003
Left curved shaft	2.52 mm	2.50mm	4.15mm	STS-ART-10-0004
Right curved shaft	2.52 mm	2.50mm	4.15mm	STS-ART-10-0005



Duckbee Basket Punch, medium



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0006
Left curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0007
Right curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0008
Upswept shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0009



Duckbee Basket Punch, aggressor medium



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Upswept shaft	2.52 mm	3.20mm	4.15mm	STS-ART-10-0010



Duckbee Basket Punch, medium, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0011
Left curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0012
Right curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0013
Upswept shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0014



Duckbee Basket Punch, large



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	4.50mm	6.15mm	STS-ART-10-0015



Duckbee Basket Punch, large, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.52 mm	4.50mm	6.15mm	STS-ART-10-0016



Duckbee Basket Punch, upbiter, diagonal cut



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Left curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0017
Right curved shaft	2.52 mm	3.20mm	5.05mm	STS-ART-10-0018

Ducklee Basket Punch

combines a wide bite with the thinnest possible tip to provide optimum maneuverability under the condyle.



Ducklee Basket Punch, small



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93 mm	2.50mm	4.15mm	STS-ART-10-0019
Left curved shaft	1.93 mm	2.50mm	4.15mm	STS-ART-10-0020
Right curved shaft	1.93 mm	2.50mm	4.15mm	STS-ART-10-0021
Upswept shaft	1.93 mm	2.50mm	4.15mm	STS-ART-10-0022



Ducklee Basket Punch, small, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93 mm	2.50mm	4.15mm	STS-ART-10-0023



Ducklee Basket Punch, medium



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0024
Left curved shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0025
Right curved shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0026
Upswept shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0027



Ducklee Basket Punch, medium, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0028
Left curved shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0029
Right curved shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0030
Upswept shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0031



Ducklee Basket Punch, large



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93mm	4.50mm	6.15mm	STS-ART-10-0032
Upswept shaft	1.93mm	4.50mm	6.15mm	STS-ART-10-0033
Right curved shaft	1.93mm	4.50mm	6.15mm	STS-ART-10-0034
Left curved shaft	1.93mm	4.50mm	6.15mm	STS-ART-10-0035



Ducklee Basket Punch, large, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93mm	4.50mm	6.15mm	STS-ART-10-0036



Ducklee Punch Elevator



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0037

Basket Punch Slim Style

offers a small aggressive tip for excellent maneuverability and multiple plane cutting capability in the tightest areas.



Basket Punch, slim style



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.90mm	1.80mm	2.89mm	STS-ART-10-0038
Upswept shaft	1.90mm	1.80mm	2.89mm	STS-ART-10-0039



Basket Punch, slim style



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft, 45° angled left jaw	1.50mm	1.80mm	2.89mm	STS-ART-10-0040
Straight shaft, 45° angled right jaw	1.50mm	1.80mm	2.89mm	STS-ART-10-0041



Hook Punch, slim style, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.90mm	1.80mm	2.89mm	STS-ART-10-0042
Left curved shaft	1.90mm	1.80mm	2.89mm	STS-ART-10-0043
Right curved shaft	1.90mm	1.67mm	2.89mm	STS-ART-10-0044

Basket Punch Oval

gives the lowest profile tip with the largest bite for the most versatility.



Basket Punch, oval



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.65mm	4.39mm	6.57mm	STS-ART-10-0045
Left curved shaft	1.65mm	4.39mm	6.57mm	STS-ART-10-0046
Right curved shaft	1.65mm	4.39mm	6.57mm	STS-ART-10-0047
Upswept shaft	1.65mm	4.39mm	6.57mm	STS-ART-10-0048



Basket Punch, oval, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.65mm	4.40mm	6.57mm	STS-ART-10-0049
Left curved shaft	1.65mm	4.40mm	6.57mm	STS-ART-10-0050
Right curved shaft	1.65mm	4.40mm	6.57mm	STS-ART-10-0051
Upswept shaft	1.65mm	4.40mm	6.57mm	STS-ART-10-0052



Basket Punch, half oval



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft, jaw right	1.93mm	3.20mm	5.00mm	STS-ART-10-0053
Straight shaft, jaw left	1.93mm	3.20mm	5.00mm	STS-ART-10-0054
15° shaft curved right, jaw right	1.93mm	3.20mm	5.00mm	STS-ART-10-0055
15° shaft curved left, jaw left	1.93mm	3.20mm	5.00mm	STS-ART-10-0056

Basket Punch Scoop

serves different purposes due to their. different sizes.



Basket Punch Scoop



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
1.0mm, straight shaft	1.10mm	1.80mm	3.25mm	STS-ART-10-0057
1.0mm, upswept shaft	1.10mm	1.80mm	3.25mm	STS-ART-10-0058
1.3mm, shaft curved left	1.13mm	1.80mm	3.25mm	STS-ART-10-0059
1.3mm, shaft curved right	1.13mm	1.80mm	3.25mm	STS-ART-10-0060



Basket Punch Scoop, 1.3mm, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.60mm	1.80mm	3.27mm	STS-ART-10-0061
Left curved shaft	1.60mm	1.80mm	3.27mm	STS-ART-10-0062
Right curved shaft	1.60mm	1.80mm	3.27mm	STS-ART-10-0063



Basket Punch Scoop



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Upswept shaft	1.30mm	1.80mm	3.50mm	STS-ART-10-0064
Straight shaft	1.60mm	1.80mm	3.50mm	STS-ART-10-0065



Basket Punch Scoop, downbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.30mm	1.80mm	3.27mm	STS-ART-10-0066



Basket Punch Scoop



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Upswept shaft	1.60mm	1.80mm	3.88mm	STS-ART-10-0067
Straight shaft	2.50mm	1.80mm	3.58mm	STS-ART-10-0068
Straight shaft	3.00mm	1.80mm	3.76mm	STS-ART-10-0069



Basket Punch Scoop, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.28mm	1.70mm	3.88mm	STS-ART-10-0070



Basket Punch, short cut



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0071
15° upswept shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0072

Basket Punch Meniscal Elevator

uses the extended tip of the lower jaw to slide easily under a tight meniscus.



Basket Punch Meniscal Elevator



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
1.3mm straight shaft	1.90mm	1.80mm	3.20mm	STS-ART-10-0073



Basket Punch Meniscal Elevator, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
1.3mm straight shaft	1.90mm	1.80mm	3.20mm	STS-ART-10-0074

Hook Punch

enables this instrument to get into the posterior horn for flush trimming of the meniscus.



Hook Punch, Blunt Nose 1.0mm



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0075
Upswept shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0076
Left curved shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0077
Right curved shaft	1.83mm	1.80mm	3.32mm	STS-ART-10-0078

Posterior Basket Punch

utilizes the short jaw and gradually tapered shaft to allow the tip to open wide for cutting in tight posterior compartments.



Posterior Basket Punch



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0079
Left curved shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0080
Right curved shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0081
Upswept shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0082



Posterior Basket Punch, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0083
Left curved shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0084
Right curved shaft	2.46mm	2.20mm	4.00mm	STS-ART-10-0085

Mega Biter

has a wide jaw to enable a highly efficient tissue resection.



Mega Biter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	1.93mm	5.20mm	7.30mm	STS-ART-10-0086
Upswept shaft	1.93mm	5.20mm	7.30mm	STS-ART-10-0087

Basket Punch Unicorn

delivers backbiting capability to offer the most direct access to the anterior meniscus.



Basket Punch Unicorn



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft, left cut	3.93mm	2.30mm	5.58mm	STS-ART-10-0088
Straight shaft, right cut	3.93mm	2.30mm	5.58mm	STS-ART-10-0089



Basket Punch Unicorn



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft, middle cut	3.50mm	3.00mm	5.00mm	STS-ART-10-0090

Rotary Basket Punch

allows successive bites, ideal for balancing the meniscal rim, this punch offers an effective angle for the anterior meniscus.



Rotary Basket Punch



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
2.2mm shaft straight, jaw 90°, angled left	1.60mm	2.20mm	4.10mm	STS-ART-20-0001
2.2mm shaft straight, jaw 90°, angled right	1.60mm	2.20mm	4.10mm	STS-ART-20-0002
3.4mm shaft straight, jaw 90°, angled left	1.60mm	3.40mm	5.30mm	STS-ART-20-0003
3.4mm shaft straight, jaw 90°, angled right	1.60mm	3.40mm	5.30mm	STS-ART-20-0004

Scissor Punch

enables small bites to ensure increased visibility while consistently providing a clean shear cut.



Scissors Punch



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0091
Left curved shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0092
Right curved shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0093



Scissors Punch, upbiter



Description	Tip Profile	Bite Width	Tip Width	Classic Handle
Straight shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0094
Left curved shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0095
Right curved shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0096
Upswept shaft	2.20mm	1.10mm	2.80mm	STS-ART-10-0097

Hook Scissors

leverages its hook teeth to prevent the tissue from slipping and delivers a clean cut.

 Hook Scissors		
Description		Classic Handle
Straight shaft		STS-ART-10-0098
Shaft 12° curved left		STS-ART-10-0099
Shaft 12° curved right		STS-ART-10-0100

Rotary Hook Scissors

offers a unique sideways cutting action and ability to fit into tight areas making it ideal for manipulation of all aspects of the meniscus.



Rotary Hook Scissors



Description	Classic Handle
Straight shaft, jaw 20° angled left	STS-ART-20-0005
Straight shaft, jaw 20° angled right	STS-ART-20-0006



Rotary Scissors, serrated



Description	Classic Handle
Straight shaft, jaw 20° angled left	STS-ART-20-0007
Straight shaft, jaw 20° angled right	STS-ART-20-0008



Rotary Hook Scissors



Description	Classic Handle
Straight shaft, jaw 60° angled left	STS-ART-20-0009
Straight shaft, jaw 60° angled right	STS-ART-20-0010

Rongeur Raptor

enables cutting and removal of tissue particularly useful in ACL surgery.



Rounger Raptor



Description		Classic Handle
Straight shaft		STS-ART-10-0101

Suture Cutter

enables suture cutting during meniscal surgery.



Meniscal Suture Cutter

Description	Classic Handle
Meniscal suture cutter	STS-ART-10-0102

Grasping Forceps

offers various versions to provide all of the possibilities for grasping tissue in each joint area.



Grasping Forceps 1x2 teeth



Description	Slide-Lock Handle
Straight shaft	STS-ART-30-0001



Grasping Forceps Cupped, serrated



Description	Slide-Lock Handle
Straight shaft, jaw wide 2.7mm	STS-ART-30-0002
Straight shaft, jaw wide 3.4mm	STS-ART-30-0003



Grasping Forceps Alligator, upbiter



Description	Slide-Lock Handle
Straight shaft, jaw wide 2.7mm	STS-ART-30-0004
Upswept shaft, jaw wide 2.7mm	STS-ART-30-0005



Grasping Forceps Alligator



Description	Slide-Lock Handle
Straight shaft, jaw wide 2.7mm	STS-ART-30-0006
Upswept shaft, jaw wide 2.7mm	STS-ART-30-0007



Grasping Forceps Alligator



Description	Classic Handle
Straight shaft, jaw wide 3.2mm	STS-ART-30-0008
Upswept shaft, jaw wide 3.2mm	STS-ART-30-0009



Tissue Grasping Forceps



Description	Classic Handle
Straight shaft, jaw wide 2.7mm	STS-ART-30-0004
Upswept shaft, jaw wide 2.7mm	STS-ART-30-0005



Knot Pusher, closed




Description	Working Length	Classic Handle
Straight shaft	125mm	STS-ART-30-0006

Drill Guide System

The femoral ACL aiming hook adapts to the anatomy of each patient and offers the user:


- Flexibility
- Precision
- Control

For the anterior cruciate ligament reconstruction, the angle of the aimer and the bullet can be set precisely. The length of the tibial drill channel can be read from the bullet before drilling and adjusted as necessary.



	Description	Part Number
	ACL / PCL Handle	STS-ART-1000-H
	ACL / PCL Bullet	STS-ART-1025-A
	ACL Tibial aimer	STS-ART-1052-T
	ACL Femoral (LP) aimer	STS-ART-1050-F
	Multiuse aimer	STS-ART-1054-M

One-Piece Instruments




	Description	Probes	Silicone Handle
	Hook Probe pointed, shaft straight	-	STS-ART-40-0002
	Hook Probe, graduated, shaft straight	3mm	STS-ART-40-0003
	Hook Probe, graduated, shaft curved right 30°	-	STS-ART-40-0004
	Hook Probe, graduated, shaft curved left 30°	-	STS-ART-40-0005
	Hook Probe, graduated, shaft straight	5mm	STS-ART-40-0006

	Description	Probes	Silicone Handle
	Spoon 30°, shaft straight	3mm	STS-ART-40-0007

	Description	Probes	Silicone Handle
	Ring curette 30°, shaft straight	3mm	STS-ART-40-0008
	Ring curette 15°, shaft straight	5mm	STS-ART-40-0009
	Ring curette 30°, shaft straight	7mm	STS-ART-40-0010

One-Piece Instruments



	Description	Probes	Silicone Handle
	Knife 45°, retrograde, shaft straight	5mm	STS-ART-40-0011
	Bajonet Knife, shaft straight	-	STS-ART-40-0012
	Banana Knife, serrated, shaft	-	STS-ART-40-0013
	Banana Knife, smooth, shaft straight	-	STS-ART-40-0014
	Scalpel Knife, shaft straight	-	STS-ART-40-0015
	Knife, serrated, shaft straight	3mm	STS-ART-40-0016
	Sickle Knife, shaft straight	-	STS-ART-40-0017
	Hook Knife 90°, shaft straight		STS-ART-40-0018
	Description	Probes	Silicone Handle
	Meniscotome Smillie, shaft straight	3mm	STS-ART-40-0019
	Meniscotome Smillie, shaft straight	5mm	STS-ART-40-0020
	Meniscotome Smillie, shaft straight	7mm	STS-ART-40-0021
	Description	Probes	Silicone Handle
	Cartilage Rasp, fine, shaft straight	-	STS-ART-40-0022
	Cartilage Rasp, rough, shaft straight	-	STS-ART-40-0023

Chondro Pick

The Chondro Pick was specifically designed for cartilage treatment defects of micro fracture on large and medium joints.

The instruments are available in 5 different tips.



	Description	Length	Ø	Silicone Handle
	Chondro Pick, straight	125mm	6mm	STS-ART-40-0024
	Chondro Pick, angled tip 25°	125mm	6mm	STS-ART-40-0025
	Chondro Pick, angled tip 40°	125mm	6mm	STS-ART-40-0026
	Chondro Pick, angled tip 60°	125mm	6mm	STS-ART-40-0027
	Chondro Pick, angled tip 90°	125mm	6mm	STS-ART-40-0028



Description	Silicone Handle
Mallet for Picks	STS-ART-40-0029



Simple 7 Instruments Set


Instrument Type	Description	Tip Profile	Bite Width	Tip Width	Part Number
Ducklee Basket Punch, medium	Straight shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0024
Ducklee Basket Punch, medium	Shaft curved left	1.93mm	3.20mm	5.00mm	STS-ART-10-0025
Ducklee Basket Punch, medium	Shaft curved right	1.93mm	3.20mm	5.00mm	STS-ART-10-0026
Ducklee Basket Punch, medium	Shaft upswept	1.93mm	3.20mm	5.00mm	STS-ART-10-0027
Ducklee Basket Punch, medium, upbiter	Straight shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0028
Hook Punch, slim style, upbiter	Shaft straight	1.90mm	1.80mm	2.89mm	STS-ART-10-0042
Grasping Forceps Alligator, 3.4mm	Shaft straight	N/A	N/A	N/A	STS-ART-30-0008



Flexible 15 Instruments Set

Instrument Type	Description	Tip Profile	Bite Width	Tip Width	Part Number
Ducklee Basket Punch, medium	Straight shaft	1.93 mm	3.20mm	5.00mm	STS-ART-10-0024
Ducklee Basket Punch, medium	Shaft curved left	1.93mm	3.20mm	5.00mm	STS-ART-10-0025
Ducklee Basket Punch, medium	Shaft curved right	1.93mm	3.20mm	5.00mm	STS-ART-10-0026
Ducklee Basket Punch, medium	Shaft upswept	1.93mm	3.20mm	5.00mm	STS-ART-10-0027
Ducklee Basket Punch, medium, upbiter	Straight shaft	1.93mm	3.20mm	5.00mm	STS-ART-10-0028
Hook Punch, slim style, upbiter	Shaft straight	1.90mm	1.80mm	2.89mm	STS-ART-10-0042
Basket Punch Scoop	1.0mm, straight shaft	1.10mm	1.80mm	3.25mm	STS-ART-10-0057
Posterior Basket Punch, upbiter	Shaft straight	2.46mm	2.20mm	4.00mm	STS-ART-10-0083
Basket Punch unicorn	Shaft straight, left cut	3.93mm	2.30mm	5.58mm	STS-ART-10-0088
Basket Punch unicorn	Shaft straight, right cut	3.93mm	2.30mm	5.58mm	STS-ART-10-0089
Scissors Punch	Shaft straight	2.20mm	1.10mm	2.80mm	STS-ART-10-0091
Rotary Basket Punch	2.2mm, shaft straight, jaw 90°, angled left	1.60mm	2.20mm	4.10mm	STS-ART-20-0001
Rotary Basket Punch	2.2mm, shaft straight, jaw 90 degree, angled right	1.60	2.20	4.10	STS-ART-20-0002
Grasping Forceps 1 x 2 teeth	Shaft straight	N/A	N/A	N/A	STS-ART-30-0001
Grasping Forceps Alligator, 3.4mm	Shaft straight	N/A	N/A	N/A	STS-ART-30-0008

Knee Suture Passer



Knee Suture Passer


Description	Length	Classic Handle
Suture passer with open middle/straight jaw opening, without ratchet.	167mm	STS-ART-200-167SP

Knee Suture Passer Needles



Description	Working Length	Diameter	Classic Handle
Non-sterile packed single use nitinol needle	167mm	0.30mm	STS-ART-230-167NN
Sterile packed single use nitinol needle, boxes w/ 6 units	167mm	0.30mm	STS-ART-231-167NN

Shoulder Suture Passer



Shoulder Suture Passer

Description	Length	Bite Depth	Working Opening	Diameter	Part Number
Suture Passer with middle/straight jaw opening, without ratchet	148mm	22 mm	8 mm	4.5 mm	STS-ART-400-170SP
Suture Passer with middle/straight jaw opening, with ratchet	148mm	22 mm	8 mm	4.5 mm	STS-ART-401-170SP
Suture Passer with sideward jaw opening, without ratchet	148mm	22 mm	8 mm	4.5 mm	STS-ART-402-170SP
Suture Passer with sideward jaw opening, with ratchet	148mm	22 mm	8 mm	4.5 mm	STS-ART-403-170SP

Shoulder Suture Passer Needles



Description	Working Length	Diameter	Part Number
Non-sterile packed single use nitinol needle	170mm	0.30mm	STS-ART-430-170NN
Sterile packed single use nitinol needle, boxes w/ 6 units	170mm	0.30mm	STS-ART-431-170NN

Modulock®

The Strauss ModuLock Universal Arthroscopic Cannula System

The Strauss ModuLock Universal Arthroscopic Cannula System has been engineered to provide optimal joint insertion for hip arthroscopic repair procedures. The versatility and easy portal exchange offered by the Strauss ModuLock System helps improve surgery time and patient outcomes.



Cannulas

Description	Size	Color Code	Part Number
Regular flow cannula	4.5mm		STS-ART-110-45
Regular flow cannula	5.0mm		STS-ART-110-50
Regular flow cannula	5.5mm		STS-ART-110-55

Obturator

Description	Size	Color Code	Part Number
Peek cannulated obturator	4.5mm		STS-ART-112-45
Peek cannulated obturator	5.0mm		STS-ART-112-50
Peek cannulated obturator	5.5mm		STS-ART-112-55

Bridge

Description	Part Number
Universal rotatable snap lock bridge	STS-ART-115-00
Inflow cannula adapter with stopcock	STS-ART-IC

HD Arthroscopes



General Information:


- HD Image
- Color Coded
- Compatible with different manufacturers
- Autoclavable for 18min at 134°C / 273°F

- The optical system is calculated and aligned with the best possible resolution, thus suited for use with HD resolution cameras.
- The dimensions and optical parameters allows the scope to be compatible with other manufactures.
- Scopes are subjected to endurance tests under maximum load to ensure material compatibility.

- Outstanding performance (picture quality, insertion, irrigation)
- Safe and reliable (robust, laser welded, easy to clean)
- Durability extends life time, resulting in low cost of ownership.
- Optimized optical system for very bright images with rich contrast.


Color Coding

Color coding allows for fast and easy identification of the scope's field of view orientation.



Angle	Letter	Color
0°	A	Green
30°	B	Red
45°	F	Black
70°	C	Yellow

HD Arthroscopes




For Knee, Shoulder and Hip

Part Number	Diameter	View	Length	Compatible
STS-1040-175AA	4.0 mm	0°	175 mm	Karl Storz®
STS-1040-175BA	4.0 mm	30°	175 mm	Karl Storz®
STS-1040-175CA	4.0 mm	70°	175 mm	Karl Storz®
STS-1040-175FA	4.0 mm	45°	175 mm	Karl Storz®

For Small Joints

Part Number	Diameter	View	Length	Compatible
STS-1019-060AA	1.9 mm	0°	60 mm	Karl Storz®
STS-1019-060BA	1.9 mm	30°	60 mm	Karl Storz®
STS-7027-110ASA	2.7mm	0°	110mm	Karl Storz®
STS-7027-110BSA	2.7mm	30°	110mm	Karl Storz®
STS-7027-110CSA	2.7mm	70°	110mm	Karl Storz®
STS-1027-187ASA	2.7 mm	0°	187 mm	Karl Storz®
STS-1027-187BSA	2.7 mm	30°	187 mm	Karl Storz®
STS-1027-187CSA	2.7 mm	70°	187 mm	Karl Storz®

Arthroscope Sheaths



Sheath Sets

All arthroscopy sheath sets includes a rotatable single or dual stopcock sheath and obturator (blunt, conical blunt or triangular obturator).

Part Number	Diameter	Length	Lock System	Degree Scope Use	Compatible
STS-ART-190-60	2.8 mm	60 mm	Bayonet lock	30°	Karl Storz®
STS-ART-201-27	3.8 mm	110 mm	Bayonet lock	0°, 30° & 70°	Karl Storz®
STS-ART-113-28	3.8 mm	187 mm	Bayonet lock	0°, 30° & 70°	Karl Storz®
STS-ART-113-28K	3.8 mm	187 mm	Quick lock	0°, 30° & 70°	Karl Storz®
STS-ART-116-40	4.0 mm	175 mm	Quick lock	0°, 30°, 45° & 70°	Karl Storz®
STS-ART-116-40-HF-SET	4.0 mm	175 mm	Quick Lock	0°, 30° & 70°	Karl Storz®
STS-ART-116-40/Y	4.0 mm	175 mm	Bayonet lock	0°, 30°, 45° & 70°	Karl Storz®
STS-ART-115-40D	4.0 mm	175 mm	Bayonet lock	0°, 30°, 45° & 70°	Dyonics®

Accessories

Arthroscope Replacement Accessories

Part Number	Description
Light Cable Adaptor	Karl Storz® connection light cable adaptor
Light Cable Adaptor	Olympus® / Richard Wolf® connection light cable adaptor



Sterilization Tray

Part Number	Description
STS-STB-H28-140	Sterilization tray for arthroscopy forceps instruments, holds up to 12 forceps
STS-STB-H26	Sterilization tray for arthroscope and accessories 260mm x 120mm x 50mm



Fiber Optic Light Cables

Part Number	Description
STS-LC-KK-10GR	For Karl Storz Light Source to Karl Storz Endoscopes
STS-LC-KO-10GR	For Olympus Light Source to Storz/Olympus Endoscopes
STS-LC-STR-10GR	For Stryker Light Source to Karl Storz Endoscopes



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System