

See the difference



NEPHROLOGY

Instrument Product Guide

Index


Nephroscopes 2

13.5Fr. Mini-Nephroscopy 3

26Fr. Nephroscopy 4

Nephroscopy Instruments 5

Nephroscopes



Mini-Nephroscopes

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Length
STS-5135-220HA	13.5Fr.	45°	8°	6Fr.	220 mm
STS-5135-320HA	13.5Fr.	45°	8°	6Fr.	320 mm
STS-URO-5000-BRGA-002	13.5Fr. mini-nephroscope bridge				
STS-5190-220HA	19Fr.	45°	8°	9Fr.	220mm

Nephroscopes

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Flushing Channel	Length
STS-5229-250AA	22Fr.	90°	15°	14Fr.	ø1.5mm	250 mm
STS-5269-250AA	26Fr.	90°	0°	ø4.7mm	Two of ø1.5mm	250 mm
STS-5269-250HA	26Fr.	90°	6°	ø5mm	Two of ø1.5mm	250 mm

13.5Fr. Mini-Nephroscopy



Nephroscope

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Length
STS-5135-220HA	13.5Fr.	45°	8°	6Fr.	220 mm

Nephroscope Sheath

Part Number	Description	Size
STS-NEP-420-19	1 fixed stopcock sheath	16.5/19.5Fr.

Applicator for Nephroscope Sheath

Part Number	Description	Size
STS-NEP-420-14-SET	Applicator for nephroscope sheath	19.5Fr.

Dilator

Part Number	Description	Size
STS-NEP-420-16	Metal dilator with central hole for guide wire. For use with 16.5/19.5Fr. nephroscope sheath.	16.5Fr. X 300mm

26Fr. Nephroscopy



Nephroscope

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Flushing Channel	Length
STS-5269-250AA	26Fr.	90°	0°	ø4.7mm	Two of ø1.5mm	250 mm

Nephroscope Sheath

Part Number	Description	Size
STS-NEP-526-250	1 rotatable stopcock sheath	26Fr. X 250mm

Obturator

Part Number	Description	Size
STS-NEP-521-250	Centrally hollow obturator for 26Fr. sheath	250mm

Dilator Set

Part Number	Description
STS-NEP-500-DIL-SET	Metal dilator tubes with 9, 12,15, 18, 21, 24, 27, 30Fr., 1 rigid and 1 flexible guide wire

Flexible Guide Wire

Part Number	Description
STS-NEP-500-202	Hollow flexible guide wire


Nephroscopy Instruments



Grasping Forceps

For use with the 26Fr. Nephrology Set

Part Number	Description	Size
STS-NEP-500-142	Grasping forceps with fenestrated jaw for softer stones	3.5mm x 350mm
STS-NEP-500-143	Grasping forceps with serrated jaw for larger stones	3.5mm x 350mm



Spring Handle Trident

For use with the 26Fr. Nephrology Set

Part Number	Description	Size
STS-NEP-500-140	Trident forceps for large stones and fragments	3.5mm x 380mm



+1 (305) 436-0599
sales@straussurgical.com
www.straussurgical.com