


See the difference



Meniscectomy

Instrument Product Guide

Index

Arthroscopes & Sheaths 3

Essential Instruments 5

Arthroscopes & Sheaths

Arthroscopes & Sheath


Our high definition endoscopes use a complex proximal and distal sapphire window system that is scratch free, fogs less and prevents corrosion and fluid invasion.

Our low-maintenance non coated stainless steel body is resistant to corrosion, rust and staining.

We include three light cable adapters, allowing compatibility with ACMI®, Richard Wolf® and Karl Storz®.


Our scope body and sheath are laser-welded tightly to ensure that all parts are properly secured together.

Our scratch-resistant insulated synthetic ocular objective and distal sapphire-window prevents internal glass corrosion to ensure a clear, precisely defined image.



Strauss Surgical HD Arthroscopes

Part Number	Diameter	View	Length	Compatible
STS-1040-175BA	4.0 mm	30°	175 mm	Karl Storz®
STS-1040-175CA	4.0 mm	70°	175 mm	Karl Storz®



High Flow

Regular Flow

Strauss Surgical Sheaths

Part Number	Description	Diameter	Length	Lock System	Compatible
STS-ART-116-40	Regular flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-116-40-HF	High flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40-HF	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®

Essential Instruments

Essential Instruments



Part Number	Description	Tip Profile	Bite Width	Tip Width
STS-ART-10-0053	Basket punch, half oval, shaft straight, jaw right	1.93mm	3.20mm	5.00mm
STS-ART-10-0054	Basket punch, half oval, shaft straight, jaw left	1.93mm	3.20mm	5.00mm



Part Number	Description	Tip Profile	Bite Width	Tip Width
STS-ART-10-0091	Scissor punch, straight shaft	2.20mm	1.10mm	2.80mm



Part Number	Description	Tip Profile	Bite Width	Tip Width
STS-ART-10-0042	Hook punch, slim style upbiter, straight shaft	1.90mm	1.80mm	2.89mm




Part Number	Description	Tip Profile	Bite Width	Tip Width
STS-ART-10-0024	Duckly basket punch, medium, straight shaft	1.93 mm	3.20mm	5.00mm
STS-ART-10-0025	Duckly basket punch, medium, curved shaft left	1.93 mm	3.20mm	5.00mm
STS-ART-10-0026	Duckly basket punch, medium, curved shaft right	1.93 mm	3.20mm	5.00mm



Part Number	Description	Jaw Width
STS-ART-30-0008	Alligator grasping forceps, straight shaft	3.2mm

Essential Instruments



Part Number	Description	Probes
STS-ART-40-0006	Hook probe, graduated, shaft straight	5mm



Part Number	Description
STS-ART-40-0023	Cartilage rasp, rough, shaft straight



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System