



strauss®
SURGICAL

See the difference



Rotator Cuff Repair

Instrument Product Guide



Index

Arthroscopes & Sheaths	2
Suture Passers	4
Essential Instruments	6
Complementary Instruments	8

Arthroscopes & Sheaths

Arthroscopes & Sheath

Our high definition endoscopes use a complex proximal and distal sapphire window system that is scratch free, fogs less and prevents corrosion and fluid invasion.

Our low-maintenance non coated stainless steel body is resistant to corrosion, rust and staining.

We include three light cable adapters, allowing compatibility with ACMI®, Richard Wolf® and Karl Storz®.


Our scope body and sheath are laser-welded tightly to ensure that all parts are properly secured together.

Our scratch-resistant insulated synthetic ocular objective and distal sapphire-window prevents internal glass corrosion to ensure a clear, precisely defined image.



Strauss Surgical HD Arthroscopes

Part Number	Diameter	View	Length	Compatible
STS-1040-175BA	4.0 mm	30°	175 mm	Karl Storz®
STS-1040-175CA	4.0 mm	70°	175 mm	Karl Storz®




Strauss Surgical Sheaths

Part Number	Description	Diameter	Length	Lock System	Compatible
STS-ART-116-40	Regular flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-116-40-HF	High flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40-HF	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®

Suture Passers

Suture Passers



Strauss Surgical Suture Passers

Part Number	Description	Length	Diameter
STS-ART-401-170SP	Middle/straight jaw opening, with ratchet	170mm	4.5 mm
STS-ART-403-170SP	Sideward jaw opening, with ratchet	170mm	4.5 mm



Strauss Surgical Suture Passer Needles

Part Number	Description	Length	Diameter
STS-ART-430-170NN	Non-sterile packed single use nitinol needle	232mm	0.30mm
STS-ART-431-170NN	Sterile packed single use nitinol needle, boxes w/ 6 units	232mm	0.30mm

Essential Instruments

Essential Instruments



Part Number	Description
TAG-Pb338	T.A.G. Medical calibrated probe



Part Number	Description
TAG-301200	T.A.G. Knot manipulator full loop black



Part Number	Description
TAG-403300	T.A.G. Suture scissors




Part Number	Description
TAG-221200	T.A.G. Sliding suture cutter



Part Number	Description
TAG-601200	T.A.G. Suture manipulator

Complementary Instruments

Complementary Instruments



T.A.G. Suture Leaders

	Part Number	Description
a.	TAG-232012	T.A.G. Suture Leader 70° degree right green
b.	TAG-232022	T.A.G. Suture Leader 70° degree left red
c.	TAG-232013	T.A.G. Suture Leader 30° degree right green
d.	TAG-232023	T.A.G. Suture Leader 30° degree left red
e.	TAG-232011	T.A.G. Suture Leader 45° degree right green
f.	TAG-232021	T.A.G. Suture Leader 45° degree left red



Complementary Instruments



T.A.G. Penetrating Graspers

Part Number	Description
TAG-398300	T.A.G. Penetrating grasper straight
TAG-398331	T.A.G. Penetrating grasper 35° degree up
TAG-398321	T.A.G. Penetrating grasper 45° degree left
TAG-398311	T.A.G. Penetrating grasper 45° degree right



T.A.G. Rasp

Part Number	Description
TAG-901131	T.A.G. Rasp 20° Degree Up Bend Blue
TAG-901141	T.A.G. Rasp 20° Degree Down Bend

Complementary Instruments



T.A.G. Tissue Liberator

Part Number	Description
TAG-902100	T.A.G. Tissue liberator blade up black
TAG-902141	T.A.G. Tissue liberator blade down black
TAG-902200	T.A.G. Tissue liberator atraumatic up black
TAG-902241	T.A.G. Tissue liberator atraumatic down black



T.A.G. Clever Hook

Part Number	Description
TAG-328332	T.A.G. Clever hook right green
TAG-328342	T.A.G. Clever hook left red



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System