




Laser Resection

26Fr. Laser Resection



Cystoscopes

Part Number	Diameter	View	Length
STS-2040-302BA	4.0mm	30°	302mm

26 Fr. Laser Bayonet-Lock Sheaths

Part Number	Description	Size
STS-URO-340-20-SET	Rotatable continuous flow sheath rotatable	26Fr.
Consisting of:		
STS-URO-209-06	26Fr. outer sheath	
STS-URO-209-26	Rotatable inner sheath with ceramic tip	
STS-URO-320-10	Obturator	

Visual Obturator

Part Number	Description	Size
STS-URO-235-24	Laser working insert, passive, for use with laser probes up to 2.3mm/7Fr.	24/26Fr.


Working Elements

Part Number	Description
STS-URO-047-08	Laser-working element for 0.8mm fiber laser
STS-URO-047-10	Laser-working element for 1.2mm fiber laser

Laser Adapters

Part Number	Description
STS-URO-047-26	Laser adapter
STS-URO-047-82	Tuohy Borts adapter

22Fr. Laser Resection



Cystoscopes

Part Number	Diameter	View	Length
STS-2029-302BA	2.9mm	30°	302mm

22 Fr. Laser Bayonet-Lock Sheaths

Part Number	Description	Size
STS-HYS-310-22-SET	Rotatable continuous flow sheath rotatable	22Fr.
Consisting of:		
STS-HYS-309-01	22Fr. outer sheath	
STS-HYS-315-21	19Fr. rotatable inner sheath with ceramic tip	
STS-HYS-310-19	Obturator	


Working Elements

Part Number	Description
STS-HYS-046-08	Passive laser-working element for up 0.8mm fiber laser
STS-HYS-046-12	Passive laser-working element for up 1.2mm fiber laser

Laser Adapters

Part Number	Description
STS-URO-047-26	Laser adapter
STS-URO-047-82	Tuohy Borts adapter

18.5Fr. Laser Resection



Cystoscopes

Part Number	Diameter	View	Length
STS-2029-302BA	2.9mm	30°	302mm

18.5 Fr. Laser Bayonet-Lock Sheaths

Part Number	Description	Size
STS-HYS-175-18-SET	Rotatable continuous flow sheath rotatable	22Fr.
Consisting of:		
STS-HYS-175-23-BP	18.5Fr. obturator	
STS-HYS-175-24-BP	17.5/18.5Fr. outer sheath	
STS-HYS-175-26-BP	17.5/18.5Fr. rotatable inner sheath with ceramic tip	

Working Elements

Part Number	Description
STS-HYS-175-08	Passive laser-working element for up 0.8mm fiber laser

Laser Adapters

Part Number	Description
STS-URO-047-26	Laser adapter
STS-URO-047-82	Tuohy-borst laser adapter, for ø 06.-1.4mm laser probes



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System