




strauss®
SURGICAL

See the difference



Hip Arthroscopy for FAI

Instrument Product Guide



Index

Arthroscopes & Sheaths	4
Instruments	6

Arthroscopes & Sheaths

Arthroscopes & Sheath

Our high definition endoscopes use a complex proximal and distal sapphire window system that is scratch free, fogs less and prevents corrosion and fluid invasion.

Our low-maintenance non coated stainless steel body is resistant to corrosion, rust and staining.

We include three light cable adapters, allowing compatibility with ACMI®, Richard Wolf® and Karl Storz®.

Our scope body and sheath are laser-welded tightly to ensure that all parts are properly secured together.



Our scratch-resistant insulated synthetic ocular objective and distal sapphire-window prevents internal glass corrosion to ensure a clear, precisely defined image.

Strauss Surgical HD Arthroscopes

Part Number	Diameter	View	Length	Compatible
STS-1027-187BSA	2.7mm	30°	187mm	Karl Storz®
STS-1027-187CSA	2.7mm	70°	187mm	Karl Storz®
STS-1040-175BA	4.0 mm	30°	175 mm	Karl Storz®



High Flow



Regular Flow



Strauss Surgical Sheaths

Part Number	Description	Diameter	Length	Lock System	Compatible
STS-ART-113-28	Regular flow sheath	2.7mm	187mm	Bayonet lock	Karl Storz®
STS-ART-216-28	Blunt obturator	2.7mm	187mm	Bayonet lock	Karl Storz®
STS-ART-116-40	Regular flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-116-40-HF	High flow sheath	4.0 mm	175 mm	Quick lock	Karl Storz®
STS-ART-216-40-HF	Blunt obturator	4.0 mm	175 mm	Quick lock	Karl Storz®


Instruments

Instruments




T.A.G. Knot Pusher

Part Number	Description
TAG-301400	T.A.G. Knot pusher



T.A.G. Stitchers

Part Number	Description
TAG-328422	T.A.G. Medical sixter left
TAG-328412	T.A.G. Medical sixter right



T.A.G. Suture Passer

Part Number	Description
TAG-399400	T.A.G. Medical suture passer straight
TAG-399431	T.A.G. Medical suture passer 35° up

Instruments



T.A.G. Slotted Cannula & Cannulated Switching Stick

Part Number	Description
TAG-754420	T.A.G Cannulated switching stick 4.0mm
TAG-754480	T.A.G. Cannulated switching stick 5.5mm
TAG-754470	T.A.G. Cannulated switching stick 6.0mm

Part Number	Description
TAG-754490	T.A.G. Slotted cannula 5.5mm
TAG-754410	T.A.G. Slotted cannula 6.0mm



T.A.G. Manipulator

Part Number	Description
TAG-601400	T.A.G Medical suture manipulator

Instruments



T.A.G. Grasper

Part Number	Description
TAG-391400	T.A.G. Medical grasper straight 5mm



T.A.G. Suture Cutter

Part Number	Description
TAG-231400	T.A.G. Medical suture cutter; trim cord



T.A.G. Probe

Part Number	Description
TAG-138415	T.A.G. Calibrated probe, curved
TAG-138404	T.A.G. Calibrated probe 90°

Instruments



T.A.G. Bitter

Part Number	Description
TAG-727400	T.A.G Medical bitter straight 3mm
TAG-727431	T.A.G. Medical bitter 3mm, 15° up



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System