

See the difference




Nephroscopy

Instrument Product Guide


Index

Nephroscopes	3
13.5Fr. Mini-Nephroscopy	5
26Fr. Nephroscopy	8



Nephroscopes

Nephroscopes



Mini-Nephroscopes

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Length
STS-5135-220HA*	13.5Fr.	45°	8°	6Fr.	220mm
STS-5135-320HA*	13.5Fr.	45°	8°	6Fr.	320mm
STS-URO-5000-BRGA-002	13.5Fr. mini-nephroscope bridge				

*Nephroscope bridge sold separately.



Nephroscopes

Part Number	Diameter of Sheath	Eye Piece Angle	Direction of View	Channel	Length
STS-5022-180-45AA	22Fr.	45°	0°	15.3Fr.	180mm
STS-5047-200AA	26Fr.	45°	0°	15Fr.	200mm
STS-5269-250AA	26Fr.	90°	0°	15Fr.	250mm
STS-5269-250HA	26Fr.	90°	6°	15Fr.	250mm
STS-5229-250HA	26Fr.	90°	15°	15Fr.	250mm

Mini-Nephroscopy

13.5Fr. Mini-Nephroscopy



Nephroscope

Part Number	Diameter	Eye Piece Angle	Direction of View	Channel	Length
STS-5135-220HA-SET*	13.5Fr.	45°	8°	6Fr.	220mm

*SET includes the Nephroscope and bridge.



Nephroscope Sheath

Part Number	Description	Size
STS-NEP-420-19	1 fixed stopcock sheath	16.5/19.5Fr.



Applicator for Nephroscope Sheath


Part Number	Description	Size
STS-NEP-420-14-SET	Applicator for nephroscope sheath	19.5Fr.

13.5Fr. Mini-Nephroscopy



Dilator

Part Number	Description	Size
STS-NEP-420-16	Metal dilator with central hole for guide wire. For use with 16.5/19.5Fr. nephroscope sheath.	16.5Fr. x 300mm



Mini-Nephroscopy Forceps

Part Number	Description	Size
STS-URO-BIO-50-40F	Flexible biopsy forceps, oval spoon, double action	5Fr. x 40cm
STS-URO-GRA-50-40AF	Flexible grasping forceps, alligator jaw, double action	5Fr. x 40cm
STS-HYS-SCI-50-40BF	Scissor pointed, single action jaws	5Fr. x 40cm

26Fr. Nephroscopy

26Fr. Nephroscopy



Nephroscopes

Part Number	Diameter of Sheath	Eye Piece Angle	Direction of View	Channel	Length
STS-5047-200AA	26Fr.	45°	0°	15Fr.	200mm
STS-5269-250AA	26Fr.	90°	0°	15Fr.	250mm



Nephroscopy Sheaths


Part Number	Description	Size	Length
STS-NEP-526-200	Rotatable 1 stopcock sheath	26Fr.	200mm
STS-NEP-521-200	Centrally hollow obturator	26Fr.	200mm
STS-NEP-526-250	Rotatable 1 stopcock sheath	26Fr.	250mm
STS-NEP-521-250	Centrally hollow obturator	26Fr.	250mm



Metal Dilator Set


Part Number	Description
STS-NEP-500-DIL-SET	Metal dilator set containing 8 dilators (9, 12, 15, 18, 21, 24, 27 and 30Fr.) with one rigid hollow and one flexible hollow guide wire

26Fr. Nephroscopy



10Fr. Grasping Forceps For use with the 26Fr. Nephroscopy Set

Part Number	Description	Size
STS-NEP-500-142	Stone grasping forceps with fenestrated jaws	35cm
STS-NEP-500-143	Stone grasping forceps with serrated jaws	35cm



Spring Handle 10Fr. Forceps For use with the 26Fr. Nephroscopy Set

Part Number	Description	Size
STS-NEP-500-140	Trident forceps	38cm



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 436-0599 | Fax: +1 (305) 436-0399
ISO 13485:2016 – Certified Quality Management System