



ACE SOUTHERN
Surgical Solutions

PRODUCT CATALOG

2022-2023



WELCOME.

HERE ARE A FEW TIPS ON FINDING THE PRODUCT YOU'RE LOOKING FOR



PRODUCT CATALOG

TABLE OF CONTENTS

SECTION QUICK LINKS

JUMP DIRECTLY TO A SECTION

Or you can click the INDEX link and find the page number of a specific product.

SEARCH BY NAME OR ITEM NUMBER

Click FIND and enter the name of the product or item number you are looking for.

PHARMACEUTICALS	1
ANESTHESIA SUPPLIES EMERGENCY RESPONSE	2
BONE GRAFTING	3
COSMETIC AND RESTORATIVE DENTISTRY	4
DISPOSABLES	5
HANDPIECES AND BURS	6
IMPRESSION MATERIAL	7
INFECTION CONTROL	8
INSTRUMENTS	9
IV ADMINISTRATION	10
PATIENT DIAGNOSTIC	11
STERILIZATION	12
SUCTION / IRRIGATION	13
SYRINGE AND NEEDLE	14
UNIFORMS	15
WOUND CLOSURE / SUTURE	16
WOUND MANAGEMENT	17
X-RAY	18
MISCELLANEOUS	19
INDEX	20

AN ALLIANCE

Built for Your Success

TRUST



ACE SOUTHERN

ACE SOUTHERN: A Solid Foundation You Can Trust



- **Dr. Carchidi, DDS** founds ACE Surgical Supply Company in 1966
- New company headquarters in Brockton, MA
- First ACE branded product (IV Extension tube)

- **Mike Stamps**, a pharmacist, founded Southern Anesthesia and Surgical (SAS) in 1984 for Oral and Maxillofacial Surgeons practicing in-office sedation
- ACE expands campus to accommodate ISO8 Clean Room and achieves ISO9001 Certification

- SAS establishes AAOMS Advantage partnership providing benefits to AAOMS members
- ACE launches a full regenerative product portfolio
- Henry Schein forms partnerships with both ACE and Southern Anesthesia



ACE Surgical Supply and Southern Anesthesia and Surgical have joined forces as **ACE SOUTHERN**, a partnership that simply made sense given their decades long history and similar core values. Each company, established by experts in their field, identified gaps in the dental surgical field and set forth to support and empower the surgical dental community.

Today, with 90+ years of combined industry experience **ACE SOUTHERN** will continue the legacy each company built individually as a leading and trusted partner of dental surgical solutions.

With the new **ACE SOUTHERN** you now have access to:

- > Best-in-class dental surgical supplier collaboration for more comprehensive product offerings
- > Superior support dedicated to exceeding your expectations across all channels: in-person, by telephone or through a new, easy to navigate website
- > Enhanced manufacturing processes and warehouses to ensure on time, every time product delivery
- > Rewarded loyalty with easy to reach levels with the Alliance Program (ASAP)
- > Time and money saving Replenish Program for automated product fulfillment



- SAS opens new offices and distribution center in West Columbia, SC
- SAS forms partnership with AAP, Perio Perks to benefit AAP members
- SAS drives e222 online ordering system of schedule II narcotics
- ACE is first dental manufacturer to obtain American Association of Tissue Bank (AATB) accreditation
- SAS develops Rx Alerts to circumvent difficulties from shortages and backorders and opens new distribution facilities in Nevada and Oklahoma
- ACE Surgical and Southern Anesthesia form an alliance to become **ACE SOUTHERN**

Learn more about ACE SOUTHERN by visiting ACESOUTHERN.COM today!



TABLE OF CONTENTS

PHARMACEUTICALS	1
ANESTHESIA SUPPLIES EMERGENCY RESPONSE	2
BONE GRAFTING	3
COSMETIC AND RESTORATIVE DENTISTRY	4
DISPOSABLES	5
HANDPIECES AND BURS	6
IMPRESSION MATERIAL	7
INFECTION CONTROL	8
INSTRUMENTS	9
IV ADMINISTRATION	10
PATIENT DIAGNOSTIC	11
STERILIZATION	12
SUCTION / IRRIGATION	13
SYRINGE AND NEEDLE	14
UNIFORMS	15
WOUND CLOSURE / SUTURE	16
WOUND MANAGEMENT	17
X-RAY	18
MISCELLANEOUS	19
INDEX	20

COLOR-CODED SECTION REFERENCE

MAKES FINDING WHAT YOU NEED A BREEZE.

Simply identify the section you want and its color, match it to the colors on the outside edges of the pages, and flip right to your section!

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

Introducing the ACE SOUTHERN Alliance Program (ASAP)

When you create an alliance with **ACE SOUTHERN** you have a committed partner in your success. Part of that commitment is to reward you **for your loyalty with significant savings on every dollar you spend.**

As an ASAP Qualified Customer, you will benefit from savings and special promotions across all the products **ACE SOUTHERN** offers to help your practice prosper.

	Premier \$15K Annual Spend	Premier Plus \$25K Annual Spend*
5% savings on all Rx, Supplies and Equipment	●	
10% savings on all Instruments	●	●
15-20% savings on Bone Grafting Materials and Implants	●	●
10% savings on all Rx, Supplies and Equipment		●

**AAOMS/AAP members only*

Become an ASAP Member by Enrolling Today • acesouthern.com/ASAP

I'm ready to buy.
**Which Price
Tier Column
Do I Use?**

Regular Price

Most customers fall into this category. It's not "suggested retail", but reflects the strong buying power of **ACE SOUTHERN**, which we pass along to you as **everyday low prices.**

Member Price

In partnership with the **AAOMS** and **AAP** societies, their members are given a special 5% discount on general merchandise, including Rx, supplies and equipment.

Premier Price

As shown above, **enrolled members** who have met the required annual spend of \$15,000 are qualified to receive the **Premier Discount.**

Premier+ Price

As a special reward for our **enrolled society members**, once they have met the required annual spend of \$25,000 they will receive the **Premier+ Discount.**



Avoid having your form returned! We're here to help!!!

Call our DEA Compliance Team **1.888.222.3722, prompt 5.**

DEA 222 Form Sample (New Single-Page)

Once completed, make a copy for your records and mail original to us.

PURCHASER INFORMATION Dr. John Doe 123 Anywhere Street Anywhere, US 12345			REGISTRATION INFORMATION REGISTRATION #: REGISTERED AS: SCHEDULES: ORDER FORM NUMBER: DATE ISSUED:			FILL IN SUPPLIER DEA NUMBER: ① <input type="text"/> PART 2: TO BE FILLED IN BY PURCHASER ② ACE SOUTHERN Surgical Solutions BUSINESS NAME ③ One Southern Court STREET ADDRESS ④ West Columbia, SC 29169 CITY, STATE, ZIP CODE					
PART 1: TO BE FILLED IN BY PURCHASER ① Print or Type Name and Title ⑩ <i>John Doe, M.D.</i> Signature of Requesting Official (must be authorized to sign order form) Date				PART 5: TO BE FILLED IN BY PURCHASER ⑤ Today's Date		PART 3: ALTERNATE SUPPLIER IDENTIFICATION - to be filled in by first supplier (name in part 2) if order is endorsed to another supplier to fill. ALTERNATE DEA# <input type="text"/> Signature - by first supplier Official authorized to execute on behalf of supplier Date					
ITEM	NO. OF PACKAGES ⑥	PACKAGE SIZE ⑦	NAME OF ITEM ⑧	NUMBER REC'D	DATE REC'D	PART 4: TO BE FILLED IN BY SUPPLIER				NUMBER SHIPPED	DATE SHIPPED
1	2	10 x 2ml	Fentanyl amps								
2	4	20ml	Fentanyl vial								
3	3	30ml	Demerol 50mg/ml								
4	1	10/box	Morphine 2mg/ml 1ml Carpujects								
5											
6											
7											
8											
9											
10											
11											
12											
13											
14											
15											
16											
17											
18											
19											
20											
4 ⑨	← LAST LINE COMPLETED (MUST BE 20 OR LESS)										

For information on the e222 ELECTRONIC FORM and online ordering, please...

SEE PAGE 246

Requirements For Properly Completed 222 Forms:

DEA requires that your 222 form address be the same as the address on your current DEA Certificate. DO NOT fill out Supplier DEA Number in Part 2, or anything in Part 3. This information will be completed by ACE SOUTHERN after your order has shipped.

- ① **Fill in Supplier DEA Number:** Please call 888.222.3722, prompt 5 for assistance
- ② **Business Name:** ACE SOUTHERN Surgical Solutions
- ③ **Street Address:** One Southern Court
- ④ **City, State, Zip Code:** West Columbia, SC 29169
- ⑤ **Date:** Today's Date (the date you are filling out the form)
- ⑥ **Number of Packages:** The quantity of the drug being ordered
- ⑦ **Package Size:** The size of the drug being ordered (ie. 10x2ml, 20ml)
- ⑧ **Name of Item:** The name and description/strength of the drug being ordered (ie. Fentanyl amps, Fentanyl vial, Demerol 50mg/ml)
- ⑨ **Last Line Completed:** This number should correspond to the Line Item No. of the last line on which a product was entered
- ⑩ **Signature of Physician or Power of Attorney***
*If the signature is anyone other than the Registered Physician, we must have a copy of the Power of Attorney on file.
- ⑪ **Printed name and title of person signing the form**

Mistake anywhere on this form? You will need to VOID your entire form, keep on file for your records, and start with a new one. Please do not write over mistakes in an attempt to correct. This is considered an alteration. We cannot accept forms with alterations or errors.

PHARMACEUTICALS

1 Brand/Generic Pharmaceutical Comparison Quick Reference

Brand Name	Generic Comparison	Brand Name	Generic Comparison
Achromycin	Tetracycline HCl	Lasix	Furosemide
Adenocard	Adenosine	Levophed	Norepinephrine
Adrenalin Chloride	Epinephrine HCl	Marcaine HCl	Bupivacaine HCl
Alfenta	Alfentanil HCl	Motrin	Ibuprofen
Alupent	Metaproterenol Sulfate	Mycitracin	Neomycin, Polymyxin-B Sulfate & Bacitracin
Ambien	Zolpidem Tartrate	Narcan	Naloxone HCl
Amoxil	Amoxicillin	Neo-Synephrine	Phenylephrine HCl
Anaprox	Naproxen Sodium	Neosporin	Neomycin, Polymyxin-B Sulfate & Bacitracin
Ancef/Kefzol	Cefazolin Sodium	Nitrostat	Nitroglycerin
Anectine/Quelicin	Succinylcholine Chloride	Normodyne/Trandate	Labetalol HCl
Aristocort	Triamcinolone	Nubain	Nalbuphine HCl
Ativan	Lorazepam	Numorphan	Oxymorphone HCl
Benadryl	Diphenhydramine	Pavulon	Pancuronium Bromide
Betadine	Povidone-Iodine	Pentothal	Thiopental Sodium
Brethine	Terbutaline Sulfate	Pfizerpen	Penicillin G Potassium
Brevibloc	Esmolol HCl	Phenergan	Promethazine HCl
Carbocaine/Polocaine	Mepivacaine HCl	Pontocaine HCl	Tetracaine HCl
Cardizem	Diltiazem HCl	Precedex	Dexmedetomidine HCl
Cipro	Ciprofloxacin	Procardia	Nifedipine
Cleocin	Clindamycin	Pronestyl	Procainamide HCl
Compazine	Prochlorperazine	Prostigmin	Neostigmine
Cordarone	Amiodarone	Proventil	Albuterol
Cortef	Hydrocortisone	Quelicin	Succinylcholine Chloride
Cortisporin	Neomycin, Poly B, Hydrocortisone	Recombivax HB	Hepatitis B Vaccine
Dalmane	Flurazepam HCl	Reglan	Metoclopramide HCl
Dantrium	Dantrolene Sodium	Revonto	Dantrolene Sodium
Decadron	Dexamethasone	Robinul	Glycopyrrolate
Demerol HCl	Meperidine HCl	Rocephin	Ceftriaxone Sodium
Dilantin	Phenytoin	Romazicon	Flumazenil
Dilaudid	Hydromorphone HCl	Sublimaze	Fentanyl
Diprivan	Propofol	Sufenta	Sufentanil Citrate
Duramorph	preservative free Morphine Sulfate	Sumycin	Tetracycline HCl
Engerix-B	Hepatitis B Vaccine	Suprane	Desflurane
Erythrocin Stearate	Erythromycin Stearate	Tagamet	Cimetidine
Forane	Isoflurane	Toradol	Ketorolac Tromethamine
Garamycin	Gentamicin	Trandate/Normodyne	Labetalol HCl
Halcion	Triazolam	Trimox	Amoxicillin
Havrix	Hepatitis B Vaccine	Ultane	Sevoflurane
Hespan	Hetastarch	Ultiva	Remifentanil
Inapsine	Droperidol	Valium	Diazepam
Inderal	Propranolol HCl	Versed	Midazolam HCl
Isoptin	Verapamil HCl	Vistaril	Hydroxyzine Pamoate
Keflex	Cephalexin	Vitamin B-12	Cyanocobalamin
Kefzol	Cefazolin Sodium	Xanax	Alprazolam
Kenacort	Triamcinolone	Xylocaine	Lidocaine HCl
Ketalar	Ketamine HCl	Zemuron	Rocuronium Bromide
Klonopin	Clonazepam	Zofran	Ondansetron HCl
Kenalog	Triamcinolone Acetonide		
Lanoxin	Digoxin		

SCHEDULE II PHARMACEUTICALS

Cocaine Topical Solution	4%	4ml vial	each	\$575.79	\$547.00	\$547.00	\$518.21
Cocaine Topical Solution	4%	10ml vial	each	\$598.99	\$569.04	\$569.04	\$539.09
Demerol	50mg/ml	30ml vial	each	\$116.29	\$110.48	\$110.48	\$104.66
Demerol	25mg	1ml needleless luer lock carpuject	10/box	\$77.49	\$73.62	\$73.62	\$69.74
Demerol	50mg	1ml needleless luer lock carpuject	10/box	\$81.29	\$77.23	\$77.23	\$73.16
Demerol	100mg	1ml needleless luer lock carpuject	10/box	\$87.29	\$82.93	\$82.93	\$78.56
Duramorph	1mg/ml	10ml ampule	10/box	\$185.49	\$176.22	\$176.22	\$166.94
Fentanyl	10 x 2ml	ampule	10/box	\$29.99	\$28.49	\$28.49	\$26.99
Fentanyl	10 x 5ml	ampule	10/box	\$35.99	\$34.19	\$34.19	\$32.39
Fentanyl	0.05mg/ml	2ml vial	25/box	\$71.99	\$68.39	\$68.39	\$64.79
Fentanyl	0.05mg/ml	5ml vial	25/box	\$74.99	\$71.24	\$71.24	\$67.49
Fentanyl	0.05mg/ml	10ml vial	25/box	\$137.79	\$130.90	\$130.90	\$124.01
Fentanyl	0.05mg/ml	20ml vial	each	\$20.99	\$19.94	\$19.94	\$18.89
Fentanyl	0.05mg/ml	20ml vial	25/box	\$242.29	\$230.18	\$230.18	\$218.06
Fentanyl	0.05mg/ml	50ml vial	each	\$44.99	\$42.74	\$42.74	\$40.49
Fentanyl	0.05mg/ml	50ml vial	25/box	\$599.99	\$569.99	\$569.99	\$539.99
Hydromorphone	2mg/ml	20ml vial	each	\$37.29	\$35.43	\$35.43	\$33.56
Morphine	10mg/ml	1ml vial	25/box	\$92.79	\$88.15	\$88.15	\$83.51
Morphine	2mg/ml	1ml needleless luer lock carpuject	10/box	\$31.29	\$29.73	\$29.73	\$28.16
Morphine	4mg/ml	1ml needleless luer lock carpuject	10/box	\$31.29	\$29.73	\$29.73	\$28.16
Morphine	8mg/ml	1ml needleless luer lock carpuject	10/box	\$31.29	\$29.73	\$29.73	\$28.16
Morphine	10mg/ml	1ml needleless luer lock carpuject	10/box	\$31.29	\$29.73	\$29.73	\$28.16
Ultiva	1mg/3ml	3ml vial	10/box	\$875.79	\$832.00	\$832.00	\$788.21
Ultiva	2mg/5ml	5ml vial	10/box	\$1,634.49	\$1,552.77	\$1,552.77	\$1,471.04

A properly completed DEA 222 Form (see example on page 3) AND a copy of your current Federal DEA Registration Certificate are required for the purchase of all Class II Pharmaceuticals. If you have questions or would like assistance in completing your form, please call our DEA Compliance Team at 1.888.222.3722.

Compliance Centered



Our NABP Drug Distribution Accreditation ensures that the pharmaceuticals we ship are 100% authentic, protecting you from counterfeit drugs. And, our turnkey services will help keep you compliant.



PHARMACEUTICALS

1 SCHEDULE III, IV & V PHARMACEUTICALS

Alprazolam	0.25mg	tablets	100/bottle	\$10.99	\$10.44	\$10.44	\$9.89
Alprazolam	0.5mg	tablets	100/bottle	\$13.49	\$12.82	\$12.82	\$12.14
Alprazolam	1mg	tablets	100/bottle	\$16.49	\$15.67	\$15.67	\$14.84
Alprazolam	2mg	tablets	100/bottle	\$29.99	\$28.49	\$28.49	\$26.99
Brevital	500mg	50ml vial	each	\$115.29	\$109.53	\$109.53	\$103.76
Buprenex	0.3mg	1ml ampule	5/box	\$149.99	\$142.49	\$142.49	\$134.99
Byfavo	20mg		10/box	\$557.29	\$529.43	\$529.43	\$501.56
Clonazepam	0.5mg	tablets	100/bottle	\$9.49	\$9.02	\$9.02	\$8.54
Clonazepam	1mg	tablets	100/bottle	\$14.99	\$14.24	\$14.24	\$13.49
Clonazepam	2mg	tablets	100/bottle	\$21.79	\$20.70	\$20.70	\$19.61
Diazepam	5mg/ml	10ml vial	each	\$61.29	\$58.23	\$58.23	\$55.16
Diazepam	5mg/ml	10ml vial	10/box	\$589.79	\$560.30	\$560.30	\$530.81
Diazepam	2mg	tablets	100/bottle	\$4.79	\$4.55	\$4.55	\$4.31
Diazepam	5mg	tablets	100/bottle	\$7.29	\$6.93	\$6.93	\$6.56
Diazepam	5mg	tablets	500/bottle	\$25.29	\$24.03	\$24.03	\$22.76
Diazepam	10mg	tablets	100/bottle	\$7.79	\$7.40	\$7.40	\$7.01
Diazepam	10mg	tablets	500/bottle	\$28.29	\$26.88	\$26.88	\$25.46
Diazepam Oral Solution	5mg/ml	30ml	each	\$54.99	\$52.24	\$52.24	\$49.49
Diazepam Oral Solution	5mg/5ml	500ml	each	\$97.99	\$93.09	\$93.09	\$88.19
Diprivan	10mg/ml	20ml vial	10/box	\$106.29	\$100.98	\$100.98	\$95.66
Diprivan	10mg/ml	100ml vial	10/box	\$495.99	\$471.19	\$471.19	\$446.39
Ephedrine	50mg/ml	1ml sd vial	each	\$65.29	\$62.03	\$62.03	\$58.76
Ephedrine	50mg/ml	1ml sd vial	25/box	\$1,511.29	\$1,435.73	\$1,435.73	\$1,360.16
Halcion	0.25mg	tablets	10/bottle	\$834.79	\$793.05	\$793.05	\$751.31
Ketalar	10mg/ml	20ml md vial	each	\$33.99	\$32.29	\$32.29	\$30.59
Ketalar	10mg/ml	20ml md vial	10/box	\$333.29	\$316.63	\$316.63	\$299.96
Ketalar	50mg/ml	10ml md vial	10/box	\$213.79	\$203.10	\$203.10	\$192.41
Ketalar	100mg/ml	5ml md vial	10/box	\$412.29	\$391.68	\$391.68	\$371.06
Ketamine	10mg/ml	20ml md vial	10/box	\$309.99	\$294.49	\$294.49	\$278.99
Ketamine	50mg/ml	10ml md vial	10/box	\$102.99	\$97.84	\$97.84	\$92.69
Ketamine	100mg/ml	5ml md vial	10/box	\$129.99	\$123.49	\$123.49	\$116.99
Ketamine	100mg/ml	10ml md vial	10/box	\$224.99	\$213.74	\$213.74	\$202.49
Lorazepam	2mg/ml	1ml carpuject	10/box	\$53.79	\$51.10	\$51.10	\$48.41
Lorazepam	2mg/ml	1ml vial	10/box	\$45.79	\$43.50	\$43.50	\$41.21
Lorazepam	2mg/ml	1ml vial	25/box	\$102.49	\$97.37	\$97.37	\$92.24
Lorazepam	2mg/ml	10ml vial	10/box	\$258.79	\$245.85	\$245.85	\$232.91
Lorazepam	1mg	tablets	100/bottle	\$19.99	\$18.99	\$18.99	\$17.99
Lorazepam	2mg	tablets	100/bottle	\$35.99	\$34.19	\$34.19	\$32.39
Midazolam	1mg/ml	2ml vial	25/box	\$39.79	\$37.80	\$37.80	\$35.81
Midazolam	1mg/ml	5ml ft vial	10/box	\$48.49	\$46.07	\$46.07	\$43.64
Midazolam	1mg/ml	10ml vial	each	\$10.29	\$9.78	\$9.78	\$9.26
Midazolam	1mg/ml	10ml ft vial	10/box	\$49.29	\$46.82	\$46.82	\$44.36

continued on following page

A copy of your current Federal DEA Registration Certificate is required for the purchase of controlled substances.

SCHEDULE III, IV & V PHARMACEUTICALS - CONTINUED

1

Midazolam	5mg/ml	1ml ft vial	10/box	\$20.99	\$19.94	\$19.94	\$18.89
Midazolam	5mg/ml	2ml vial	10/box	\$37.99	\$36.09	\$36.09	\$34.19
Midazolam	5mg/ml	5ml ft vial	each	\$9.29	\$8.83	\$8.83	\$8.36
Midazolam	5mg/ml	5ml ft vial	10/box	\$87.79	\$83.40	\$83.40	\$79.01
Midazolam	5mg/ml	10ml vial	each	\$13.29	\$12.63	\$12.63	\$11.96
Midazolam	5mg/ml	10ml vial	10/box	\$102.99	\$97.84	\$97.84	\$92.69
Midazolam Syrup	2mg/ml	118ml bottle	each	\$99.99	\$94.99	\$94.99	\$89.99
Nalbuphine	10mg/ml	10ml ft vial	each	\$54.79	\$52.05	\$52.05	\$49.31
Nalbuphine	20mg/ml	1ml ampule	10/box	\$117.79	\$111.90	\$111.90	\$106.01
Nalbuphine	20mg/ml	10ml ft vial	each	\$79.49	\$75.52	\$75.52	\$71.54
Propofol	10mg/ml	20ml vial	5/box	\$124.29	\$118.08	\$118.08	\$111.86
Propofol	10mg/ml	50ml vial	20/box	\$299.29	\$284.33	\$284.33	\$269.36
Propofol	10mg/ml	100ml vial	10/box	\$344.99	\$327.74	\$327.74	\$310.49
Temazepam	15mg	capsules	100/bottle	\$17.99	\$17.09	\$17.09	\$16.19
Temazepam	30mg	capsules	500/bottle	\$154.29	\$146.58	\$146.58	\$138.86
Triazolam	0.125mg	tablets	100/bottle	\$315.29	\$299.53	\$299.53	\$283.76
Triazolam	0.25mg	tablets	100/bottle	\$234.29	\$222.58	\$222.58	\$210.86
Zolpidem Tartrate	10mg	tablets	100/bottle	\$42.99	\$40.84	\$40.84	\$38.69

A copy of your current Federal DEA Registration Certificate is required for the purchase of controlled substances.

EMERGENCY BOX PREFILLED SYRINGES

Atropine Sulfate	0.1mg/ml	10ml syringe	each	\$19.99	\$18.99	\$18.99	\$17.99
Calcium Chloride	10%	10ml syringe	each	\$23.79	\$22.60	\$22.60	\$21.41
Dextrose	25%	10ml needleless syringe	each	\$19.29	\$18.33	\$18.33	\$17.36
Dextrose	50%	50ml needleless syringe	each	\$19.29	\$18.33	\$18.33	\$17.36
Epinephrine	1:10,000	10ml syringe	10/box	\$143.49	\$136.32	\$136.32	\$129.14
Epinephrine	0.15mg	auto injector	2/box	\$514.99	\$489.24	\$489.24	\$463.49
Epinephrine	0.3mg	auto injector	2/box	\$514.99	\$489.24	\$489.24	\$463.49
Epipen Junior	0.15mg	auto injector	2/box	\$799.99	\$759.99	\$759.99	\$719.99
Epipen	0.3mg	auto injector	2/box	\$799.99	\$759.99	\$759.99	\$719.99
Furosemide	10mg/ml	10ml needleless syringe	each	\$35.79	\$34.00	\$34.00	\$32.21
Glucagon Emergency Kit	1mg	1ml syringe	kit	\$407.79	\$387.40	\$387.40	\$367.01
Lidocaine	1%	5ml needleless syringe	each	\$13.29	\$12.63	\$12.63	\$11.96
Lidocaine	2%	5ml needleless syringe	each	\$8.29	\$7.88	\$7.88	\$7.46
Naloxone	1mg/ml	2ml syringe	each	\$51.49	\$48.92	\$48.92	\$46.34
Sodium Bicarbonate	8.4%	10ml syringe	each	\$25.99	\$24.69	\$24.69	\$23.39
Sodium Bicarbonate	8.4%	50ml syringe	10/box	\$154.79	\$147.05	\$147.05	\$139.31

PHARMACEUTICALS

1

MISCELLANEOUS EMERGENCY BOX PHARMACEUTICALS

Adenosine	6mg	2ml vial	each	\$24.79	\$23.55	\$23.55	\$22.31
Amiodarone	50mg/ml	3ml vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Atrovent HFA	17mcg	12.9gm inhaler	each	\$645.99	\$613.69	\$613.69	\$581.39
Charcoal Powder, Activated		1 oz unit dose	each	\$36.99	\$35.14	\$35.14	\$33.29
Clonidine	0.1mg	tablets	100/bottle	\$12.99	\$12.34	\$12.34	\$11.69
Diazepam*	5mg/ml	10ml vial	each	\$54.79	\$52.05	\$52.05	\$49.31
Digoxin	0.25mg/ml	2ml ampule	each	\$10.99	\$10.44	\$10.44	\$9.89
Diltiazem	5mg/ml	5ml sd vial	each	\$10.79	\$10.25	\$10.25	\$9.71
Diphenhydramine	25mg	capsules	100/bottle	\$5.79	\$5.50	\$5.50	\$5.21
Diphenhydramine	50mg	capsules	100/bottle	\$6.79	\$6.45	\$6.45	\$6.11
Diphenhydramine	50mg/ml	1ml sd vial	each	\$3.99	\$3.79	\$3.79	\$3.59
Ephedrine*	50mg/ml	1ml sd vial	each	\$65.29	\$62.03	\$62.03	\$58.76
Epinephrine	1:1000	1ml ampule	each	\$19.29	\$18.33	\$18.33	\$17.36
Flumazenil	0.1mg/ml	5ml vial	each	\$23.99	\$22.79	\$22.79	\$21.59
Flumazenil	0.1mg/ml	10ml md vial	each	\$36.29	\$34.48	\$34.48	\$32.66
Furosemide	10mg/ml	2ml sd vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Furosemide	10mg/ml	4ml sd vial	each	\$9.99	\$9.49	\$9.49	\$8.99
Furosemide	10mg/ml	10ml sd vial	each	\$10.49	\$9.97	\$9.97	\$9.44
Insta-Glucose		31gm tube	each	\$19.29	\$18.33	\$18.33	\$17.36
Magnesium Sulfate	50%	2ml vial	each	\$3.79	\$3.60	\$3.60	\$3.41
Naloxone	0.4mg/ml	1ml vial	each	\$33.79	\$32.10	\$32.10	\$30.41
Nitrolingual Spray	0.4mg	4.9gm pump	each	\$404.79	\$384.55	\$384.55	\$364.31
Nitrolingual Spray	0.4mg	12gm pump	each	\$439.99	\$417.99	\$417.99	\$395.99
Nitrostat	0.4mg	tablets	25/bottle	\$22.29	\$21.18	\$21.18	\$20.06
Nitrostat	0.4mg	tablets	100/bottle	\$89.79	\$85.30	\$85.30	\$80.81
Phenylephrine	10mg/ml	1ml sd vial	each	\$14.99	\$14.24	\$14.24	\$13.49
ProAir HFA	90mcg	8.5gm inhaler	each	\$98.29	\$93.38	\$93.38	\$88.46
Procainamide	500mg/ml	2ml vial	each	\$133.49	\$126.82	\$126.82	\$120.14
Propranolol	1mg/ml	1ml vial	each	\$13.79	\$13.10	\$13.10	\$12.41
Revonto	20mg	vial	6/box	\$537.99	\$511.09	\$511.09	\$484.19
Solu-Cortef	250mg	2ml act-o-vial	each	\$48.99	\$46.54	\$46.54	\$44.09
Ventolin HFA	90mcg	8gm inhaler	each	\$46.79	\$44.45	\$44.45	\$42.11
Ventolin HFA	90mcg	18gm inhaler	each	\$79.79	\$75.80	\$75.80	\$71.81
Verapamil	2.5mg/ml	2ml vial	each	\$42.79	\$40.65	\$40.65	\$38.51

*A copy of your current Federal DEA Registration Certificate is required for the purchase of controlled substances.

MISCELLANEOUS PHARMACEUTICALS

Adenosine	6mg	2ml vial	each	\$24.79	\$23.55	\$23.55	\$22.31
Adenosine	6mg	2ml vial	10/box	\$235.79	\$224.00	\$224.00	\$212.21
Adrenalin	1mg/ml	1ml sd vial	25/box	\$568.79	\$540.35	\$540.35	\$511.91
Albuterol Inhaler • see ProAir Inhaler above or page 14							
Aminocaproic Acid	250mg/ml	20ml vial	each	\$16.99	\$16.14	\$16.14	\$15.29
Aminocaproic Acid	250mg/ml	20ml vial	25/box	\$319.29	\$303.33	\$303.33	\$287.36

continued on following page

MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Aminophylline	250mg	10ml vial	each	\$25.79	\$24.50	\$24.50	\$23.21
Aminophylline	250mg	10ml vial	25/box	\$527.79	\$501.40	\$501.40	\$475.01
Aminophylline	500mg	20ml vial	each	\$26.29	\$24.98	\$24.98	\$23.66
Aminophylline	500mg	20ml vial	25/box	\$603.79	\$573.60	\$573.60	\$543.41
Amiodarone	50mg/ml	3ml vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Amiodarone	50mg/ml	3ml vial	25/box	\$162.79	\$154.65	\$154.65	\$146.51
Amoxicillin	250mg	capsules	100/bottle	\$16.49	\$15.67	\$15.67	\$14.84
Amoxicillin	500mg	capsules	100/bottle	\$22.49	\$21.37	\$21.37	\$20.24
Amoxicillin	500mg	capsules	500/bottle	\$86.79	\$82.45	\$82.45	\$78.11
Ampicillin	1gm	vial	each	\$10.99	\$10.44	\$10.44	\$9.89
Ampicillin	1gm	vial	10/box	\$69.99	\$66.49	\$66.49	\$62.99
Ampicillin	2gm	vial	each	\$16.49	\$15.67	\$15.67	\$14.84
Ampicillin	2gm	vial	10/box	\$153.99	\$146.29	\$146.29	\$138.59
Aspirin, Chewable	81mg	chewable tablets	36/bottle	\$5.29	\$5.03	\$5.03	\$4.76
Aspirin	325mg tab		100/bottle	\$5.49	\$5.22	\$5.22	\$4.94
Atropine Sulfate	0.1mg/ml	10ml syringe	each	\$19.99	\$18.99	\$18.99	\$17.99
Atropine Sulfate	0.1mg/ml	10ml syringe	10/box	\$163.79	\$155.60	\$155.60	\$147.41
Atropine Sulfate	0.4mg/ml	20ml vial	each	\$53.49	\$50.82	\$50.82	\$48.14
Atropine Sulfate	0.4mg/ml	20ml vial	10/box	\$467.79	\$444.40	\$444.40	\$421.01
Atrovent HFA	17mcg	12.9gm inhaler	each	\$645.99	\$613.69	\$613.69	\$581.39
Azithromycin	250mg	unit of use tablets	3x6 tablets	\$74.29	\$70.58	\$70.58	\$66.86
Benadryl Ultra	25mg	tablets	24/box	\$6.79	\$6.45	\$6.45	\$6.11
Botox Cosmetic	50 units/vial	5ml vial	each	\$469.99	\$446.49	\$446.49	\$422.99
Botox Cosmetic	100 units/vial	10ml vial	each	\$799.99	\$759.99	\$759.99	\$719.99
Botox Therapeutic	100 units/vial	10ml vial	each	\$814.99	\$774.24	\$774.24	\$733.49
Brevibloc • see Esmolol HCl - page 11							
Bupivacaine	0.5%	50ml md vial	each	\$4.49	\$4.27	\$4.27	\$4.04
Calcium Chloride	10%	10mg/ml syringe	each	\$23.79	\$22.60	\$22.60	\$21.41
Calcium Chloride	10%	10mg/ml syringe	10/box	\$225.29	\$214.03	\$214.03	\$202.76
Calcium Gluconate	10%	10ml vial	each	\$18.29	\$17.38	\$17.38	\$16.46
Calcium Gluconate	10%	10ml vial	25/box	\$435.49	\$413.72	\$413.72	\$391.94
Cefazolin	1gm	vial	each	\$3.99	\$3.79	\$3.79	\$3.59
Cefazolin	1gm	vial	25/box	\$87.99	\$83.59	\$83.59	\$79.19
Ceftriaxone	250mg	vial	each	\$5.29	\$5.03	\$5.03	\$4.76
Ceftriaxone	250mg	vial	10/box	\$45.99	\$43.69	\$43.69	\$41.39
Ceftriaxone	1gm	vial	10/box	\$75.29	\$71.53	\$71.53	\$67.76
Celestone Soluspan	6mg/ml	5ml vial	each	\$56.29	\$53.48	\$53.48	\$50.66
Cephalexin	500mg	capsules	100/bottle	\$43.99	\$41.79	\$41.79	\$39.59
Cephalexin	500mg	capsules	500/bottle	\$129.99	\$123.49	\$123.49	\$116.99
Charcoal Powder, Activated		1 oz unit dose	each	\$36.99	\$35.14	\$35.14	\$33.29
Chlorhexidine Oral Rinse	0.12%	480ml bottle	each	\$6.29	\$5.98	\$5.98	\$5.66
Chlorpromazine	25mg/ml	1ml ampule	each	\$48.29	\$45.88	\$45.88	\$43.46
Ciprofloxacin	250mg	tablets	100/bottle	\$105.99	\$100.69	\$100.69	\$95.39

continued on following page

PHARMACEUTICALS

1 MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Ciprofloxacin	500mg	tablets	100/bottle	\$199.99	\$189.99	\$189.99	\$179.99
Ciprofloxacin	750mg	tablets	50/bottle	\$125.99	\$119.69	\$119.69	\$113.39
Citric Acid Granules	500g	16 oz	each	\$76.99	\$73.14	\$73.14	\$69.29
Cleocin	300mg	2ml sd vial	each	\$8.99	\$8.54	\$8.54	\$8.09
Cleocin	600mg	4ml sd vial	25/box	\$298.49	\$283.57	\$283.57	\$268.64
Clindamycin	150mg	capsules	100/bottle	\$45.99	\$43.69	\$43.69	\$41.39
Clindamycin	300mg	capsules	100/bottle	\$138.99	\$132.04	\$132.04	\$125.09
Clindamycin	300mg	2ml ft vial	each	\$5.49	\$5.22	\$5.22	\$4.94
Clindamycin	300mg	2ml ft vial	25/box	\$129.99	\$123.49	\$123.49	\$116.99
Clindamycin	600mg	4ml ft vial	each	\$8.49	\$8.07	\$8.07	\$7.64
Clindamycin	600mg	4ml ft vial	25/box	\$192.49	\$182.87	\$182.87	\$173.24
Clindamycin	900mg	6ml sd vial	each	\$9.99	\$9.49	\$9.49	\$8.99
Clindamycin	900mg	6ml ft vial	25/box	\$223.79	\$212.60	\$212.60	\$201.41
Clonidine	0.1mg	tablets	100/bottle	\$12.99	\$12.34	\$12.34	\$11.69
Cyanocobalamin	1000mcg/ml	1ml md vial	25/box	\$315.99	\$300.19	\$300.19	\$284.39
Cyklokapron	100mg/ml	10ml vial	10/box	\$985.99	\$936.69	\$936.69	\$887.39
Dantrolene • see Revonto - pages 8, 15							
Debacterol Canker Sore Relief	swabs		12/box	\$73.79	\$70.10	\$70.10	\$66.41
Depo-Medrol	20mg/ml	5ml vial	each	\$45.79	\$43.50	\$43.50	\$41.21
Depo-Medrol	40mg/ml	1ml sd vial	each	\$18.49	\$17.57	\$17.57	\$16.64
Depo-Medrol	40mg/ml	5ml vial	each	\$70.29	\$66.78	\$66.78	\$63.26
Depo-Medrol	40mg/ml	10ml vial	each	\$141.29	\$134.23	\$134.23	\$127.16
Depo-Medrol	80mg/ml	1ml vial	each	\$30.29	\$28.78	\$28.78	\$27.26
Depo-Medrol	80mg/ml	5ml vial	each	\$139.99	\$132.99	\$132.99	\$125.99
Dexamethasone	4mg/ml	1ml vial	25/box	\$35.99	\$34.19	\$34.19	\$32.39
Dexamethasone	4mg/ml	5ml md vial	each	\$6.49	\$6.17	\$6.17	\$5.84
Dexamethasone	4mg/ml	5ml md vial	25/box	\$129.99	\$123.49	\$123.49	\$116.99
Dexamethasone	4mg/ml	30ml md vial	each	\$27.49	\$26.12	\$26.12	\$24.74
Dexamethasone	4mg/ml	30ml md vial	25/box	\$549.99	\$522.49	\$522.49	\$494.99
Dexamethasone	10mg/ml	10ml vial	each	\$19.99	\$18.99	\$18.99	\$17.99
Dexamethasone	10mg/ml	10ml vial	10/box	\$185.59	\$176.31	\$176.31	\$167.03
Dexamethasone	0.75mg	tablets	100/bottle	\$16.49	\$15.67	\$15.67	\$14.84
Dexamethasone	4mg	tablets	100/bottle	\$33.29	\$31.63	\$31.63	\$29.96
Dexmedetomidine	100mcg/ml	2ml vial	25/box	\$464.29	\$441.08	\$441.08	\$417.86
Dextrose	25%	10ml needleless syringe	each	\$19.29	\$18.33	\$18.33	\$17.36
Dextrose	25%	10ml needleless syringe	10/box	\$155.79	\$148.00	\$148.00	\$140.21
Dextrose	50%	50ml needleless syringe	each	\$19.29	\$18.33	\$18.33	\$17.36
Dextrose	50%	50ml needleless syringe	10/box	\$136.49	\$129.67	\$129.67	\$122.84
Dextrose	50%	50ml vial	each	\$8.79	\$8.35	\$8.35	\$7.91
Dextrose	50%	50ml vial	25/box	\$217.79	\$206.90	\$206.90	\$196.01
Digoxin	0.25mg/ml	2ml ampule	each	\$10.99	\$10.44	\$10.44	\$9.89

continued on following page

MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Digoxin	0.25mg/ml	2ml ampule	25/box	\$255.79	\$243.00	\$243.00	\$230.21
Diltiazem	5mg/ml	5ml sd vial	each	\$10.79	\$10.25	\$10.25	\$9.71
Diltiazem	5mg/ml	5ml sd vial	10/box	\$92.29	\$87.68	\$87.68	\$83.06
Diphenhydramine	25mg	capsules	100/bottle	\$5.79	\$5.50	\$5.50	\$5.21
Diphenhydramine	50mg	capsules	100/bottle	\$6.79	\$6.45	\$6.45	\$6.11
Diphenhydramine	50mg/ml	1ml sd vial	each	\$3.99	\$3.79	\$3.79	\$3.59
Diphenhydramine	50mg/ml	1ml sd vial	25/box	\$54.99	\$52.24	\$52.24	\$49.49
Diphenhydramine	50mg/ml	10ml vial	each	\$19.99	\$18.99	\$18.99	\$17.99
Dobutamine	250mg	20ml vial	10/box	\$97.79	\$92.90	\$92.90	\$88.01
Dopamine	200mg/5ml	vial	each	\$10.99	\$10.44	\$10.44	\$9.89
Dopamine	200mg/5ml	vial	25/box	\$256.99	\$244.14	\$244.14	\$231.29
Dopamine	400mg/10ml	ft vial	each	\$9.29	\$8.83	\$8.83	\$8.36
Dopamine	400mg/10ml	ft vial	25/box	\$220.29	\$209.28	\$209.28	\$198.26
Dopram	20mg/ml	20ml md vial	6/box	\$456.29	\$433.48	\$433.48	\$410.66
Doxycycline	100mg	tablets	50/bottle	\$59.99	\$56.99	\$56.99	\$53.99
Doxycycline	100mg	20ml vial	10/box	\$421.29	\$400.23	\$400.23	\$379.16
EDTA Solution	17%	120ml bottle	each	\$18.99	\$18.04	\$18.04	\$17.09
EDTA Solution	17%	480ml bottle	each	\$65.29	\$62.03	\$62.03	\$58.76
Engerix B	20mcg	1ml sd vial	10/box	\$942.99	\$895.84	\$895.84	\$848.69
Epinephrine	1:10,000	10ml needleless syringe	10/box	\$143.49	\$136.32	\$136.32	\$129.14
Epinephrine	1:1000	1ml ampule	each	\$19.29	\$18.33	\$18.33	\$17.36
Epinephrine	0.15mg	auto injector	2/box	\$514.99	\$489.24	\$489.24	\$463.49
Epinephrine	0.30mg	auto injector	2/box	\$514.99	\$489.24	\$489.24	\$463.49
Epipen Junior	0.15mg	auto injector	2/box	\$799.99	\$759.99	\$759.99	\$719.99
Epipen	0.30mg	auto injector	2/box	\$799.99	\$759.99	\$759.99	\$719.99
Erythrocin	500mg	vial	10/box	\$1,162.79	\$1,104.65	\$1,104.65	\$1,046.51
Erythromycin	500mg	tablets	100/bottle	\$1,909.29	\$1,813.83	\$1,813.83	\$1,718.36
Esmolol HCl	10mg/ml	10ml sd vial	each	\$35.99	\$34.19	\$34.19	\$32.39
Esmolol HCl	10mg/ml	10ml sd vial	10/box	\$334.49	\$317.77	\$317.77	\$301.04
Eugenol		1 oz bottle	each	\$8.79	\$8.35	\$8.35	\$7.91
Eugenol		4 oz bottle	each	\$25.99	\$24.69	\$24.69	\$23.39
Flumazenil	0.1mg/ml	5ml vial	each	\$23.99	\$22.79	\$22.79	\$21.59
Flumazenil	0.1mg/ml	5ml vial	10/box	\$199.99	\$189.99	\$189.99	\$179.99
Flumazenil	0.1mg/ml	10ml md vial	each	\$36.29	\$34.48	\$34.48	\$32.66
Flumazenil	0.1mg/ml	10ml md vial	10/box	\$289.99	\$275.49	\$275.49	\$260.99
Furosemide	10mg/ml	2ml sd vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Furosemide	10mg/ml	2ml sd vial	25/box	\$164.99	\$156.74	\$156.74	\$148.49
Furosemide	10mg/ml	4ml sd vial	each	\$9.99	\$9.49	\$9.49	\$8.99
Furosemide	10mg/ml	4ml sd vial	25/box	\$233.49	\$221.82	\$221.82	\$210.14
Furosemide	10mg/ml	4ml needleless syringe	10/box	\$160.79	\$152.75	\$152.75	\$144.71
Furosemide	10mg/ml	10ml sd vial	each	\$10.49	\$9.97	\$9.97	\$9.44
Furosemide	10mg/ml	10ml sd vial	25/box	\$189.79	\$180.30	\$180.30	\$170.81
Furosemide	10mg/ml	10ml needleless syringe	each	\$35.79	\$34.00	\$34.00	\$32.21
Furosemide	10mg/ml	10ml needleless syringe	10/box	\$294.99	\$280.24	\$280.24	\$265.49

continued on following page

PHARMACEUTICALS

1

MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Gentamicin	40mg/ml	2ml vial	each	\$7.49	\$7.12	\$7.12	\$6.74
Gentamicin	40mg/ml	2ml vial	25/box	\$130.79	\$124.25	\$124.25	\$117.71
Gentamicin	40mg/ml	20ml md vial	each	\$44.99	\$42.74	\$42.74	\$40.49
Gentamicin	40mg/ml	20ml md vial	25/box	\$735.29	\$698.53	\$698.53	\$661.76
GlucaGen Diagnostic Kit	1mg w/1ml diluent		kit	\$255.99	\$243.19	\$243.19	\$230.39
Glucagon Emergency Kit	1mg	1ml syringe	kit	\$407.79	\$387.40	\$387.40	\$367.01
Glycopyrrolate	0.2mg/ml	1ml vial	each	\$26.99	\$25.64	\$25.64	\$24.29
Glycopyrrolate	0.2mg/ml	1ml vial	25/box	\$583.99	\$554.79	\$554.79	\$525.59
Glycopyrrolate	0.2mg/ml	2ml sd vial	each	\$43.29	\$41.13	\$41.13	\$38.96
Glycopyrrolate	0.2mg/ml	2ml sd vial	25/box	\$795.29	\$755.53	\$755.53	\$715.76
Glycopyrrolate	0.2mg/ml	5ml md vial	each	\$68.49	\$65.07	\$65.07	\$61.64
Glycopyrrolate	0.2mg/ml	20ml vial	each	\$209.99	\$199.49	\$199.49	\$188.99
Glycopyrrolate	0.2mg/ml	20ml vial	10/box	\$1,909.99	\$1,814.49	\$1,814.49	\$1,718.99
Hydralazine	20mg/ml	1ml sd vial	each	\$74.79	\$71.05	\$71.05	\$67.31
Hydralazine	20mg/ml	1ml sd vial	25/box	\$899.99	\$854.99	\$854.99	\$809.99
Hydrogen Peroxide	3%	pint	each	\$3.29	\$3.13	\$3.13	\$2.96
Hydroxyzine HCl	10mg/5ml	16 oz syrup	each	\$48.29	\$45.88	\$45.88	\$43.46
Hydroxyzine HCl	10mg	tablets	100/bottle	\$79.99	\$75.99	\$75.99	\$71.99
Hydroxyzine HCl	25mg	tablets	100/bottle	\$83.49	\$79.32	\$79.32	\$75.14
Hydroxyzine HCl	50mg	tablets	100/bottle	\$104.99	\$99.74	\$99.74	\$94.49
Hylenex	150u/ml	1ml vial	4/box	\$318.29	\$302.38	\$302.38	\$286.46
Ibuprofen	200mg	tablets	100/bottle	\$5.29	\$5.03	\$5.03	\$4.76
Ibuprofen	400mg	tablets	100/bottle	\$13.99	\$13.29	\$13.29	\$12.59
Ibuprofen	400mg	tablets	500/bottle	\$51.79	\$49.20	\$49.20	\$46.61
Ibuprofen	600mg	tablets	100/bottle	\$16.99	\$16.14	\$16.14	\$15.29
Ibuprofen	600mg	tablets	500/bottle	\$52.99	\$50.34	\$50.34	\$47.69
Ibuprofen	800mg	tablets	100/bottle	\$19.99	\$18.99	\$18.99	\$17.99
Ibuprofen	800mg	tablets	500/bottle	\$62.79	\$59.65	\$59.65	\$56.51
Insta-Glucose	31gm	tube	each	\$6.99	\$6.64	\$6.64	\$6.29
Isoflurane	100ml	bottle	each	\$27.49	\$26.12	\$26.12	\$24.74
Isoflurane	250ml	bottle	6/case	\$368.79	\$350.35	\$350.35	\$331.91
Kenalog-10	10mg/ml	5ml md vial	each	\$20.29	\$19.28	\$19.28	\$18.26
Kenalog-40	40mg/ml	1ml vial	each	\$15.29	\$14.53	\$14.53	\$13.76
Kenalog-40	40mg/ml	5ml md vial	each	\$63.99	\$60.79	\$60.79	\$57.59
Kenalog-40	40mg/ml	10ml md vial	each	\$106.99	\$101.64	\$101.64	\$96.29
Ketorolac	15mg/ml	1ml sd vial	25/box	\$123.79	\$117.60	\$117.60	\$111.41
Ketorolac	30mg/ml	1ml sd vial	each	\$4.79	\$4.55	\$4.55	\$4.31
Ketorolac	30mg/ml	1ml sd vial	25/box	\$71.79	\$68.20	\$68.20	\$64.61
Ketorolac	60mg/2ml	2ml sd vial	each	\$4.99	\$4.74	\$4.74	\$4.49
Ketorolac	60mg/2ml	2ml sd vial	25/box	\$98.99	\$94.04	\$94.04	\$89.09
Labetalol	5mg/ml	20ml md vial	each	\$10.29	\$9.78	\$9.78	\$9.26
Levophed	1mg/ml	4ml ampule	10/box	\$485.79	\$461.50	\$461.50	\$437.21
Lidocaine	0.5%	50ml ft vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Lidocaine	0.5%	50ml vial	25/box	\$109.29	\$103.83	\$103.83	\$98.36
Lidocaine	1%	5ml ampule	25/box	\$59.29	\$56.33	\$56.33	\$53.36

continued on following page

MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Lidocaine	1%	5ml syringe	each	\$13.29	\$12.63	\$12.63	\$11.96
Lidocaine	1%	20ml md vial	each	\$3.99	\$3.79	\$3.79	\$3.59
Lidocaine	1%	20ml md vial	25/box	\$43.49	\$41.32	\$41.32	\$39.14
Lidocaine	1%	30ml vial	25/box	\$78.29	\$74.38	\$74.38	\$70.46
Lidocaine	1%	50ml vial	each	\$4.99	\$4.74	\$4.74	\$4.49
Lidocaine	1% w/epi 1:100,000	20ml md vial	each	\$4.99	\$4.74	\$4.74	\$4.49
Lidocaine	1% w/epi 1:100,000	30ml md vial	each	\$7.99	\$7.59	\$7.59	\$7.19
Lidocaine	1% w/epi 1:100,000	30ml md vial	25/box	\$143.29	\$136.13	\$136.13	\$128.96
Lidocaine	1% w/epi 1:100,000	50ml md vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Lidocaine	1% w/epi 1:100,000	50ml md vial	25/box	\$129.29	\$122.83	\$122.83	\$116.36
Lidocaine	2%	5ml syringe	each	\$7.79	\$7.40	\$7.40	\$7.01
Lidocaine	2%	5ml needleless syringe	each	\$8.29	\$7.88	\$7.88	\$7.46
Lidocaine	2%	5ml needleless syringe	10/box	\$76.49	\$72.67	\$72.67	\$68.84
Lidocaine	2%	5ml sd vial	each	\$5.29	\$5.03	\$5.03	\$4.76
Lidocaine	2%	5ml sd vial	25/box	\$106.79	\$101.45	\$101.45	\$96.11
Lidocaine	2%	20ml md vial	each	\$3.49	\$3.32	\$3.32	\$3.14
Lidocaine	2%	20ml md vial	25/box	\$64.99	\$61.74	\$61.74	\$58.49
Lidocaine	2%	50ml md vial	each	\$5.79	\$5.50	\$5.50	\$5.21
Lidocaine	2%	50ml md vial	25/box	\$102.79	\$97.65	\$97.65	\$92.51
Lidocaine	2% w/epi 1:100,000	20ml md vial	each	\$5.79	\$5.50	\$5.50	\$5.21
Lidocaine	2% w/epi 1:100,000	20ml md vial	25/box	\$130.79	\$124.25	\$124.25	\$117.71
Lidocaine	2% w/epi 1:100,000	30ml md vial	each	\$8.49	\$8.07	\$8.07	\$7.64
Lidocaine	2% w/epi 1:100,000	50ml md vial	each	\$5.99	\$5.69	\$5.69	\$5.39
Lidocaine	2% w/epi 1:100,000	50ml md vial	25/box	\$123.29	\$117.13	\$117.13	\$110.96
Lincocin	300mg/ml	2ml vial	each	\$65.99	\$62.69	\$62.69	\$59.39
Lincocin	300mg/ml	10ml md vial	each	\$244.79	\$232.55	\$232.55	\$220.31
Lubricating Jelly	2.7gm	foil packs (unit dose)	144/box	\$9.99	\$9.49	\$9.49	\$8.99
Magnesium Sulfate	50%	2ml vial	each	\$3.79	\$3.60	\$3.60	\$3.41
Magnesium Sulfate	50%	2ml vial	25/box	\$65.29	\$62.03	\$62.03	\$58.76
Magnesium Sulfate	50%	10ml syringe	each	\$36.99	\$35.14	\$35.14	\$33.29
Marcaine	0.25%	30ml sd vial	10/box	\$60.99	\$57.94	\$57.94	\$54.89
Marcaine	0.25%	50ml md vial	each	\$8.99	\$8.54	\$8.54	\$8.09
Marcaine	0.25% w/epi 1:200,000	50ml md vial	each	\$14.79	\$14.05	\$14.05	\$13.31
Marcaine	0.5% mpf	10ml sd vial	each	\$6.99	\$6.64	\$6.64	\$6.29
Marcaine	0.5% mpf	10ml sd vial	10/box	\$56.29	\$53.48	\$53.48	\$50.66
Marcaine	0.5% mpf	30ml sd vial	each	\$7.99	\$7.59	\$7.59	\$7.19
Marcaine	0.5% mpf	30ml sd vial	10/box	\$65.49	\$62.22	\$62.22	\$58.94
Marcaine	0.5%	50ml md vial	each	\$10.29	\$9.78	\$9.78	\$9.26
Marcaine	0.5% w/epi 1:200,000	50ml md vial	each	\$21.29	\$20.23	\$20.23	\$19.16
Marcaine	0.75% mpf	10ml sd vial	10/box	\$89.49	\$85.02	\$85.02	\$80.54
Marcaine	0.75% mpf	30ml sd vial	10/box	\$116.29	\$110.48	\$110.48	\$104.66
Metoclopramide	5mg/ml	2ml sd vial	each	\$2.99	\$2.84	\$2.84	\$2.69
Metoclopramide	5mg/ml	2ml sd vial	25/box	\$60.99	\$57.94	\$57.94	\$54.89
Metoprolol	1mg/ml	5ml vial	10/box	\$184.99	\$175.74	\$175.74	\$166.49

continued on following page

PHARMACEUTICALS

1 MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Naloxone	0.4mg/ml	1ml vial	each	\$33.79	\$32.10	\$32.10	\$30.41
Naloxone	0.4mg/ml	1ml vial	10/box	\$287.99	\$273.59	\$273.59	\$259.19
Naloxone	0.4mg/ml	10ml md vial	each	\$166.99	\$158.64	\$158.64	\$150.29
Naloxone	1mg/ml	2ml syringe	each	\$51.49	\$48.92	\$48.92	\$46.34
Naloxone	1mg/ml	2ml syringe	10/box	\$41.99	\$39.89	\$39.89	\$37.79
Neostigmine	1mg/ml	10ml vial	each	\$60.49	\$57.47	\$57.47	\$54.44
Neostigmine	1mg/ml	10ml vial	10/box	\$574.99	\$546.24	\$546.24	\$517.49
Nitrolingual Spray	0.4mg	4.9gm pump spray	each	\$404.79	\$384.55	\$384.55	\$364.31
Nitrolingual Spray	0.4mg	12gm pump spray	each	\$439.99	\$417.99	\$417.99	\$395.99
Nitrostat	0.4mg	tablets	25/bottle	\$22.29	\$21.18	\$21.18	\$20.06
Nitrostat	0.4mg	tablets	100/bottle	\$89.79	\$85.30	\$85.30	\$80.81
Ondansetron	2mg/ml	2ml sd vial	each	\$2.79	\$2.65	\$2.65	\$2.51
Ondansetron	2mg/ml	2ml sd vial	25/box	\$11.49	\$10.92	\$10.92	\$10.34
Ondansetron	2mg/ml	20ml md vial	each	\$17.29	\$16.43	\$16.43	\$15.56
Ondansetron ODT	4mg	tablets (unit dose)	30/count	\$14.99	\$14.24	\$14.24	\$13.49
Penicillin G Potassium	5mmu	vial	each	\$56.99	\$54.14	\$54.14	\$51.29
Penicillin G Potassium	5mmu	vial	10/box	\$498.99	\$474.04	\$474.04	\$449.09
Penicillin VK	250mg	tablets	100/bottle	\$48.99	\$46.54	\$46.54	\$44.09
Penicillin VK	500mg	tablets	100/bottle	\$55.29	\$52.53	\$52.53	\$49.76
Peridex Oral Rinse		pint	each	\$15.49	\$14.72	\$14.72	\$13.94
Peridex Oral Rinse		64 oz bottle	each	\$41.49	\$39.42	\$39.42	\$37.34
Petrolatum Jelly		4 oz tube	each	\$1.49	\$1.42	\$1.42	\$1.34
Phenergan	25mg/ml	1ml ampule	25/box	\$113.79	\$108.10	\$108.10	\$102.41
Phenergan	25mg/ml	1ml vial	25/box	\$113.79	\$108.10	\$108.10	\$102.41
Phenylephrine	10mg/ml	1ml sd vial	each	\$14.99	\$14.24	\$14.24	\$13.49
Phenylephrine	10mg/ml	1ml sd vial	25/box	\$304.49	\$289.27	\$289.27	\$274.04
Phenytoin	50mg/ml	2ml vial	each	\$3.49	\$3.32	\$3.32	\$3.14
Phenytoin	50mg/ml	2ml vial	25/box	\$83.79	\$79.60	\$79.60	\$75.41
Phenytoin	50mg/ml	5ml vial	each	\$4.29	\$4.08	\$4.08	\$3.86
Phenytoin	50mg/ml	5ml vial	25/box	\$99.49	\$94.52	\$94.52	\$89.54
Physostigmine	1mg/ml	2ml ampule	each	\$127.99	\$121.59	\$121.59	\$115.19
Physostigmine	1mg/ml	2ml ampule	10/box	\$1,207.79	\$1,147.40	\$1,147.40	\$1,087.01
Precedex	100mcg/ml	2ml ft vial	each	\$110.49	\$104.97	\$104.97	\$99.44
Precedex	100mcg/ml	2ml ft vial	25/box	\$2,594.79	\$2,465.05	\$2,465.05	\$2,335.31
ProAir HFA	90mcg/inh	8.5gm inhaler	each	\$98.29	\$93.38	\$93.38	\$88.46
Procainamide	100mg/ml	10ml md vial	each	\$170.79	\$162.25	\$162.25	\$153.71
Procainamide	500mg/ml	2ml vial	each	\$133.49	\$126.82	\$126.82	\$120.14
Prochlorperazine	5mg/ml	2ml vial	10/box	\$305.99	\$290.69	\$290.69	\$275.39
Promethazine	25mg	tablets	100/bottle	\$51.49	\$48.92	\$48.92	\$46.34
Promethazine	25mg/ml	1ml ampule	each	\$4.79	\$4.55	\$4.55	\$4.31
Promethazine	25mg/ml	1ml vial	each	\$4.79	\$4.55	\$4.55	\$4.31
Promethazine	25mg/ml	1ml vial	25/box	\$109.29	\$103.83	\$103.83	\$98.36
Promethazine	50mg/ml	1ml ampule	each	\$6.49	\$6.17	\$6.17	\$5.84
Promethazine	50mg/ml	1ml ampule	25/box	\$155.99	\$148.19	\$148.19	\$140.39

continued on following page

MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Propranolol	1mg/ml	1ml vial	each	\$13.79	\$13.10	\$13.10	\$12.41
Propranolol	1mg/ml	1ml vial	10/box	\$99.99	\$94.99	\$94.99	\$89.99
Quelicin	20mg/ml	10ml vial	25/box	\$880.49	\$836.47	\$836.47	\$792.44
Revonto	20mg	vial	6/box	\$537.99	\$511.09	\$511.09	\$484.19
Rocephin • see Ceftriaxone - page 9							
Rocuronium	10mg/ml	5ml vial	each	\$39.99	\$37.99	\$37.99	\$35.99
Rocuronium	10mg/ml	5ml vial	10/box	\$155.99	\$148.19	\$148.19	\$140.39
Rocuronium	10mg/ml	10ml vial	each	\$49.99	\$47.49	\$47.49	\$44.99
Rocuronium	10mg/ml	10ml vial	10/box	\$191.79	\$182.20	\$182.20	\$172.61
Romazicon • see Flumazenil - pages 8, 11							
Scopolamine Transdermal	1.5mg	patch	10/box	\$312.29	\$296.68	\$296.68	\$281.06
Sevoflurane		250ml bottle	each	\$199.99	\$189.99	\$189.99	\$179.99
Silver Nitrate Applicators, Non-Sterile			100/box	\$37.79	\$35.90	\$35.90	\$34.01
Sodium Bicarbonate	7.5%	50ml syringe	each	\$21.49	\$20.42	\$20.42	\$19.34
Sodium Bicarbonate	8.4%	10ml syringe	each	\$25.99	\$24.69	\$24.69	\$23.39
Sodium Bicarbonate	8.4%	10ml syringe	10/box	\$189.99	\$180.49	\$180.49	\$170.99
Sodium Bicarbonate	8.4%	50ml sd vial	each	\$38.99	\$37.04	\$37.04	\$35.09
Sodium Bicarbonate	8.4%	50ml sd vial	25/box	\$353.29	\$335.63	\$335.63	\$317.96
Sodium Bicarbonate	8.4%	50ml syringe	10/box	\$154.79	\$147.05	\$147.05	\$139.31
Sodium Chloride	0.9%	10ml sd vial	each	\$2.79	\$2.65	\$2.65	\$2.51
Sodium Chloride	0.9%	10ml sd vial	25/box	\$46.29	\$43.98	\$43.98	\$41.66
Sodium Chloride	0.9%	20ml sd vial	each	\$2.99	\$2.84	\$2.84	\$2.69
Sodium Chloride	0.9%	20ml sd vial	25/box	\$58.49	\$55.57	\$55.57	\$52.64
Sodium Chloride	0.9%	50ml sd vial	each	\$5.79	\$5.50	\$5.50	\$5.21
Sodium Chloride	0.9%	50ml sd vial	25/box	\$95.29	\$90.53	\$90.53	\$85.76
Sodium Chloride - Bacteriostatic	0.9%	30ml vial	each	\$2.79	\$2.65	\$2.65	\$2.51
Sodium Chloride - Bacteriostatic	0.9%	30ml vial	25/box	\$51.49	\$48.92	\$48.92	\$46.34
Sodium Chloride IV Flush Solution	0.9%	10ml syringe	60/box	\$49.99	\$47.49	\$47.49	\$44.99
Solu-Cortef	100mg	2ml act-o-vial	each	\$22.99	\$21.84	\$21.84	\$20.69
Solu-Cortef	250mg	2ml act-o-vial	each	\$48.99	\$46.54	\$46.54	\$44.09
Solu-Medrol	500mg	vial	each	\$97.99	\$93.09	\$93.09	\$88.19
Solu-Medrol	1gm w/o diluent	vial	each	\$62.99	\$59.84	\$59.84	\$56.69
Solu-Medrol	40mg	1ml act-o-vial	each	\$9.29	\$8.83	\$8.83	\$8.36
Solu-Medrol	40mg	1ml act-o-vial	25/box	\$221.29	\$210.23	\$210.23	\$199.16
Solu-Medrol	125mg	2ml act-o-vial	each	\$16.29	\$15.48	\$15.48	\$14.66
Solu-Medrol	125mg	2ml act-o-vial	25/box	\$383.79	\$364.60	\$364.60	\$345.41
Solu-Medrol	500mg	4ml act-o-vial	each	\$55.99	\$53.19	\$53.19	\$50.39
Solu-Medrol	1gm	8ml act-o-vial	each	\$92.29	\$87.68	\$87.68	\$83.06
Succinylcholine	20mg/ml	10ml md vial	each	\$43.79	\$41.60	\$41.60	\$39.41
Suprane		240ml bottle	each	\$283.79	\$269.60	\$269.60	\$255.41
Surgilube • see Lubricating Jelly - page 13							
Tamiflu	75mg	unit of use capsules	10/pack	\$69.99	\$66.49	\$66.49	\$62.99
Terbutaline	1mg/ml	1ml vial	each	\$32.49	\$30.87	\$30.87	\$29.24
Terbutaline	1mg/ml	1ml vial	10/box	\$83.99	\$79.79	\$79.79	\$75.59
Thrombin Topical	5000 IU vial	w/5ml diluent	each	\$128.79	\$122.35	\$122.35	\$115.91

continued on following page

PHARMACEUTICALS

1 MISCELLANEOUS PHARMACEUTICALS - CONTINUED

Tigan	100mg/ml	2ml sd vial	each	\$69.99	\$66.49	\$66.49	\$62.99
Tigan	100mg/ml	2ml sd vial	25/box	\$1,443.79	\$1,371.60	\$1,371.60	\$1,299.41
Tranexamic Acid	100mg/ml	10ml sd vial	10/box	\$451.49	\$428.92	\$428.92	\$406.34
Tylenol Extra-Strength	500mg	caplets	100/bottle	\$17.99	\$17.09	\$17.09	\$16.19
Tylenol Extra-Strength	500mg	caplets (dose packet)	50 x 2 count	\$17.49	\$16.62	\$16.62	\$15.74
Ultane		250ml bottle	6/case	\$336.79	\$319.95	\$319.95	\$303.11
Unasyn	1.5gm	vial	10/box	\$123.99	\$117.79	\$117.79	\$111.59
Unasyn	3gm	vial	10/box	\$210.99	\$200.44	\$200.44	\$189.89
Vancomycin	500mg	vial	each	\$12.79	\$12.15	\$12.15	\$11.51
Vancomycin	500mg	vial	10/box	\$113.99	\$108.29	\$108.29	\$102.59
Vasostriect	20u/ml	1ml md vial	each	\$336.29	\$319.48	\$319.48	\$302.66
Ventolin HFA	90mcg	8gm inhaler	each	\$46.79	\$44.45	\$44.45	\$42.11
Ventolin HFA	90mcg	18gm inhaler	each	\$79.79	\$75.80	\$75.80	\$71.81
Verapamil	2.5mg/ml	2ml vial	each	\$42.79	\$40.65	\$40.65	\$38.51
Verapamil	2.5mg/ml	2ml vial	25/box	\$1,022.49	\$971.37	\$971.37	\$920.24
Verapamil	2.5mg/ml	4ml needleless syringe	each	\$51.29	\$48.73	\$48.73	\$46.16
Water For Injection - Bacteriostatic		30ml md vial	each	\$3.99	\$3.79	\$3.79	\$3.59
Water For Injection - Bacteriostatic		30ml md vial	25/box	\$85.79	\$81.50	\$81.50	\$77.21
Water For Injection - Sterile		10ml sd vial	each	\$2.99	\$2.84	\$2.84	\$2.69
Water For Injection - Sterile		10ml sd vial	25/box	\$44.49	\$42.27	\$42.27	\$40.04
Water For Injection - Sterile		20ml sd vial	each	\$3.49	\$3.32	\$3.32	\$3.14
Water For Injection - Sterile		20ml sd vial	25/box	\$65.29	\$62.03	\$62.03	\$58.76
Water For Injection - Sterile		50ml sd vial	each	\$4.29	\$4.08	\$4.08	\$3.86
Water For Injection - Sterile		50ml sd vial	25/box	\$82.29	\$78.18	\$78.18	\$74.06
Water For Injection - Sterile		100ml sd vial	each	\$7.99	\$7.59	\$7.59	\$7.19
Water For Injection - Sterile		100ml sd vial	25/box	\$156.29	\$148.48	\$148.48	\$140.66
Zithromax Z-Pack	250mg	unit of use tablets	3 x 6 tablets	\$59.99	\$56.99	\$56.99	\$53.99

OPHTHALMIC PRODUCTS

Ointments

Lacri-Lube		3.5gm tube	each	\$17.29	\$16.43	\$16.43	\$15.56
Neo/Poly/Bac		3.5gm tube	each	\$68.99	\$65.54	\$65.54	\$62.09
Neo/Poly/Dex		3.5gm tube	each	\$40.79	\$38.75	\$38.75	\$36.71

Solutions

Balanced Salt Solution		15ml bottle	6/pack	\$93.49	\$88.82	\$88.82	\$84.14
Betadine Prep Solution		30ml bottle	each	\$97.99	\$93.09	\$93.09	\$88.19
Eye Wash, Irrigating		118ml w/cup	each	\$4.79	\$4.55	\$4.55	\$4.31

OINTMENTS/CREAMS/LOTIONS

Aquaphor Healing Ointment		1.75 oz tube	each	\$7.99	\$7.59	\$7.59	\$7.19
Aquaphor Healing Ointment		14 oz jar	each	\$18.49	\$17.57	\$17.57	\$16.64
Bacitracin Ointment		1/32 oz foil packs	144/box	\$44.79	\$42.55	\$42.55	\$40.31
Bacitracin Ointment		28gm tube	each	\$3.99	\$3.79	\$3.79	\$3.59
Bacitracin Zinc Ointment		4 oz tube	each	\$19.79	\$18.80	\$18.80	\$17.81

continued on following page

OINTMENTS/CREAMS/LOTIONS - CONTINUED

Bacitracin Zinc Ointment		1 lb jar	each	\$59.29	\$56.33	\$56.33	\$53.36
Cleocin T Topical Lotion	1%	60ml bottle	each	\$197.29	\$187.43	\$187.43	\$177.56
Fluocinonide Cream	0.05%	15gm tube	each	\$46.99	\$44.64	\$44.64	\$42.29
Fluocinonide Cream	0.05%	60gm tube	each	\$146.49	\$139.17	\$139.17	\$131.84
Fluocinonide Ointment	0.05%	15gm tube	each	\$92.79	\$88.15	\$88.15	\$83.51
Hydrocortisone Ointment	1%	1 oz tube	each	\$14.99	\$14.24	\$14.24	\$13.49
Neosporin Ointment		0.5 oz tube	each	\$6.29	\$5.98	\$5.98	\$5.66
Neosporin Ointment		1 oz tube	each	\$9.99	\$9.49	\$9.49	\$8.99
Neosporin Ointment		0.03 oz foil packs (unit dose)	144/box	\$66.49	\$63.17	\$63.17	\$59.84
Nitro-Bid Ointment	2%	30gm tube	each	\$52.49	\$49.87	\$49.87	\$47.24
Nitro-Bid Ointment	2%	60gm tube	each	\$104.79	\$99.55	\$99.55	\$94.31
Silvadene Cream	1%	50gm jar	each	\$23.29	\$22.13	\$22.13	\$20.96
Triple Antibiotic Ointment		28gm tube	each	\$12.79	\$12.15	\$12.15	\$11.51

DENTAL ANESTHETIC CROSS-REFERENCE GUIDE

Generic Name	Cook-Waite (Septodont)		Dentsply		Septodont	
	Reorder No.	Product Name	Reorder No.	Product Name	Reorder No.	Product Name
Articaine Hydrochloride 4% w/Epinephrine 1:100,000 Injection					A1400	Septocaine
Articaine Hydrochloride 4% w/Epinephrine 1:200,000 Injection					A1200	Septocaine
Bupivacaine Hydrochloride 0.5% and Epinephrine 1:200,000 Injection	99184	Marcaine				
Lidocaine Hydrochloride 2% and Epinephrine 1:50,000 Injection	99169	Lidocaine			A1050	Lignospan forte
Lidocaine Hydrochloride 2% and Epinephrine 1:100,000 Injection	99167	Lidocaine	017616	Xylocaine	109560	Lignospan standard
Mepivacaine Hydrochloride 3% Injection	99171	Carbocaine	044116	Polocaine	A1003	Scandonest plain
Prilocaine 4% Hydrochloride Injection			046616	Citanest Plain		
Prilocaine 4% and Epinephrine 1:200,000 Injection			048816	Citanest Forte		

COOK-WAITE DENTAL ANESTHETICS

Lidocaine HCl 2% and Epinephrine 1:100,000 Injection

99167	Lidocaine Hydrochloride and Epinephrine Injection, USP	50 cartridges/box	\$43.99	\$41.79	\$41.79	\$39.59
-------	--	-------------------	---------	---------	---------	---------

Lidocaine HCl 2% and Epinephrine 1:50,000 Injection

99169	Lidocaine Hydrochloride and Epinephrine Injection, USP	50 cartridges/box	\$36.49	\$34.67	\$34.67	\$32.84
-------	--	-------------------	---------	---------	---------	---------

Carbocaine® 3% Injection

99171	Mepivacaine Hydrochloride Injection, USP	50 cartridges/box	\$59.49	\$56.52	\$56.52	\$53.54
-------	--	-------------------	---------	---------	---------	---------

Marcaine® 0.5% w/Epinephrine 1:200,000(as Bitartrate) Injection

99184	Bupivacaine Hydrochloride and Epinephrine Injection, USP	50 cartridges/box	\$81.49	\$77.42	\$77.42	\$73.34
-------	--	-------------------	---------	---------	---------	---------

DENTSPLY ANESTHETICS

2% Xylocaine® DENTAL w/Epinephrine 1:100,000

017616	Lidocaine HCl and Epinephrine Injection, USP	50 cartridges/box	\$52.29	\$49.68	\$49.68	\$47.06
--------	--	-------------------	---------	---------	---------	---------

3% Polocaine® DENTAL

044116	Mepivacaine HCl Injection, USP	50 cartridges/box	\$63.29	\$60.13	\$60.13	\$56.96
--------	--------------------------------	-------------------	---------	---------	---------	---------

4% Citanest® Plain DENTAL

046616	Prilocaine HCl Injection, USP	50 cartridges/box	\$62.49	\$59.37	\$59.37	\$56.24
--------	-------------------------------	-------------------	---------	---------	---------	---------

4% Citanest® Forte DENTAL w/Epinephrine 1:200,000

048816	Prilocaine and Epinephrine Injection, USP	50 cartridges/box	\$60.99	\$57.94	\$57.94	\$54.89
--------	---	-------------------	---------	---------	---------	---------



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

PHARMACEUTICALS

1



SEPTODONT ANESTHETICS

- Septodont anesthetics do not contain methylparaben preservatives
- The convenient blister packs protect cartridges from breakage by eliminating glass-to-glass contact

Lignospan® standard

109560	Lidocaine Hydrochloride 2% and Epinephrine 1:100,000 Injection, USP	50 cartridges/box	\$49.49	\$47.02	\$47.02	\$44.54
--------	---	-------------------	---------	---------	---------	---------

Scandonest® 3% plain

A1003	Mepivacaine Hydrochloride 3% Injection, USP	50 cartridges/box	\$60.79	\$57.75	\$57.75	\$54.71
-------	---	-------------------	---------	---------	---------	---------

Septocaine® w/Epinephrine 1:200,000

A1200	Articaine Hydrochloride 4% with Epinephrine 1:200,000 Injection	50 cartridges/box	\$71.79	\$68.20	\$68.20	\$64.61
-------	---	-------------------	---------	---------	---------	---------

Septocaine® w/Epinephrine 1:100,000

A1400	Articaine Hydrochloride 4% with Epinephrine 1:100,000 Injection	50 cartridges/box	\$71.79	\$68.20	\$68.20	\$64.61
-------	---	-------------------	---------	---------	---------	---------



KEYSTONE INDUSTRIES GELATO TOPICAL ANESTHETICS

- Has no bitter aftertaste
- Fast acting with no systemic absorption
- Contains 20% Benzocaine for effective pain relief
- Provides temporary relief of pain during procedures including local anesthetic injections, periodontal curettage, impression taking, scaling, intra-oral radiographs, root planing and prophylaxis

0302219	Strawberry	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302319	Cherry	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302819	Mint	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302419	Bubble gum	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302119	Mango	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302019	Raspberry	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0303519	Banana	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0301319	Grape	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0302519	Pina colada	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84
0301719	Pina colada - dye free	1 oz jar	each	\$6.49	\$6.17	\$6.17	\$5.84



CROSSTEX GUMNUMB® TOPICAL ANESTHETICS

- Contains 20% Benzocaine
- Rapid onset within 30 seconds and a short duration with virtually no systemic absorption
- No bitter aftertaste

UTBG	Bubble gum	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11
UTCH	Cherry	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11
UTG	Grape	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11
UTM	Mint	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11
UTPC	Pina colada	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11
UTS	Strawberry	1 oz jar	each	\$6.79	\$6.45	\$6.45	\$6.11

LOOK IN THE SYRINGE AND NEEDLE SECTION FOR ASPIRATING SYRINGES (PAGE 282) AND DENTAL NEEDLES (PAGE 280)

PHARMACEUTICALS

1

**BEUTLICH PHARMACEUTICALS
HURRICAIN[®] TOPICAL
ANESTHETICS**

- 20% Benzocaine Oral Anesthetic
- Safe - virtually no systemic absorption
- Works fast - rapid onset of 20 seconds
- No artificial colors



Spray

- Disposable extension tubes eliminate cross contamination
- Environmentally safe - no CFC's

2 oz Spray with 1 Disposable Extension Tube

067902	each	\$45.79	\$43.50	\$43.50	\$41.21
--------	------	---------	---------	---------	---------

2 oz Spray Kit with 200 Disposable Extension Tubes

067960	kit	\$55.79	\$53.00	\$53.00	\$50.21
--------	-----	---------	---------	---------	---------

Gel

- Eliminates pain and discomfort



087131	Wild cherry	1 oz Jar	each	\$10.99	\$10.44	\$10.44	\$9.89
088631	Pina colada	1 oz Jar	each	\$10.99	\$10.44	\$10.44	\$9.89
029331	Watermelon	1 oz Jar	each	\$10.99	\$10.44	\$10.44	\$9.89
101631	Strawberry	1 oz Jar	each	\$10.99	\$10.44	\$10.44	\$9.89

Liquid

- Perfect for hard to reach areas
- Can be used as a rinse for the back of the throat
- Ideal for irrigating into the sulcus during dental procedures such as scaling and root planing



188631	Pina colada	1 oz Jar	each	\$11.29	\$10.73	\$10.73	\$10.16
--------	-------------	----------	------	---------	---------	---------	---------

Snap-n-Go[™] Swabs

- Eliminates cross contamination
- Each individually wrapped swab applicator contains 0.005 fl oz (0.15ml) of HurriCaine Topical Anesthetic Liquid
- Wild cherry flavor



056972	72 swabs/box	\$45.99	\$43.69	\$43.69	\$41.39
--------	--------------	---------	---------	---------	---------

CETYLITE CETACAINE[®] TOPICAL ANESTHETICS

- Use for all accessible mucous membranes, except the eyes
- Cetacaine is the first topical anesthetic with three active ingredients (Benzocaine 14%, Butamben 2.0%, Tetracaine Hydrochloride 2.0%)
- Rapid onset within 30 seconds - Lasts 30-60 minutes
- Cetacaine controls pain, eases discomfort and suppresses the gag reflex, resulting in patient cooperation



Spray

- Stainless steel spray cannula provides practitioners fingertip control
- Tropical flavoring

020103	20g Bottle	each	\$174.49	\$165.77	\$165.77	\$157.04
--------	------------	------	----------	----------	----------	----------

Gel

The pump-top jar controls the amount of Cetacaine dispensed. The protective lid easily comes off and on to keep the surface clean and intact.



- Ointment applies evenly and consistently, reacts with body temperature to melt and absorb quickly into tissue and has lubricating qualities
- Strawberry flavor

0217	32g Jar	each	\$109.29	\$103.83	\$103.83	\$98.36
------	---------	------	----------	----------	----------	---------

Liquid

Well-suited for scaling and root planing or other procedures requiring site specific topical anesthesia. Apply into sulcus using Cetylite's 0.4ml syringes and microcapillary tips. Draw only the amount needed with the unique luer-lock dispenser cap. 0.4ml maximum dose.



0211	30g Bottle	each	\$299.49	\$284.52	\$284.52	\$269.54
------	------------	------	----------	----------	----------	----------

DENTSPLY INTERNATIONAL ORAQIX[®]

**(Lidocaine and Prilocaine Periodontal Gel)
2.5%/2.5%**

Oraqix is the only FDA-approved, needle-free anesthetic for scaling and root planing.



- Needle-free, increases patient comfort during anesthesia administration
- Safe*, low risk of allergic reaction and not contraindicated for use on inflamed tissue
- Easy to use, applied as a liquid, but sets to a gel, so it stays in place for effective absorption on desired treatment area
- Fast onset, anesthetizing effects in just 30 seconds, allowing treatment to begin quickly
- 20-minute duration, returns patients to normal sensitivity in about 20 minutes, so they can resume normal activity quickly
- High maximum dosage, up to 5 cartridges per appointment
- Site-specific, so multiple areas of the mouth can be treated in one visit, enabling full-mouth procedures

**Note: There are possible risks associated with the use of Oraqix, which are described in the Package Insert. These risks should be reviewed prior to using Oraqix.*

011020	1.7 gram	20/box	\$149.49	\$142.02	\$142.02	\$134.54
--------	----------	--------	----------	----------	----------	----------

PHARMACEUTICALS

1

SYNERA®

(Lidocaine 70mg/Tetracaine 70mg Topical Patch)

Fast, effective needle-stick pain prevention in a simple-to-use patch

- Effective in just 20-30 minutes
- Advanced warming technology

SYNERA is indicated for use on intact skin to provide local dermal analgesia for superficial venous access and superficial dermatological procedures.

GAL000210	10/box	\$257.29	\$244.43	\$244.43	\$231.56
-----------	--------	----------	----------	----------	----------

LIDOCAINE 2.5% AND PRILOCAINE 2.5% CREAM



5gm Tube and Dressing

035756	5 Tubes and 12 dressings/box	box	\$64.29	\$61.08	\$61.08	\$57.86
--------	------------------------------	-----	---------	---------	---------	---------

30gm Tube

066730		each	\$41.99	\$39.89	\$39.89	\$37.79
--------	--	------	---------	---------	---------	---------

LMX 4® (LIDOCAINE 4%) TOPICAL ANESTHETIC CREAM



5gm Tube and Tegaderm Dressing

088207	5 Tubes and 10 dressings/box	box	\$69.99	\$66.49	\$66.49	\$62.99
--------	------------------------------	-----	---------	---------	---------	---------

30gm Tube

088230		each	\$61.49	\$58.42	\$58.42	\$55.34
--------	--	------	---------	---------	---------	---------

SEPTODONT 5% LIDOCAINE OINTMENT

- Amide-based topical anesthetic designed to improve patient's comfort during dental procedures
- Available in a pleasant tasting mint flavor
- Mess-free, value sized jar



01T0160	50g Jar (1-3/4 oz)	each	\$29.29	\$27.83	\$27.83	\$26.36
---------	--------------------	------	---------	---------	---------	---------

MISCELLANEOUS LIDOCAINE TOPICAL ANESTHETICS



Lidocaine Jelly

2%	30ml Tube	each	\$128.79	\$122.35	\$122.35	\$115.91
----	-----------	------	----------	----------	----------	----------

Lidocaine Viscous Solution

2%	100ml Bottle	each	\$21.79	\$20.70	\$20.70	\$19.61
----	--------------	------	---------	---------	---------	---------

Lidocaine Topical Solution

4%	50ml Bottle	each	\$37.79	\$35.90	\$35.90	\$34.01
----	-------------	------	---------	---------	---------	---------

GEBAUER PAIN EASE INSTANT TOPICAL ANESTHETIC



Pain Ease is an instant topical anesthetic skin refrigerant approved to temporarily control the pain associated with needle procedures and minor surgical procedures. It can be applied to minor open wounds and intact oral mucous membranes.

- Non-drug, non-flammable
- 3.9 oz bottle

Mist Spray

GB000802		each	\$45.99	\$43.69	\$43.69	\$41.39
----------	--	------	---------	---------	---------	---------

Medium Stream Spray

000803		each	\$45.99	\$43.69	\$43.69	\$41.39
--------	--	------	---------	---------	---------	---------

GEBAUER ETHYL CHLORIDE INSTANT TOPICAL ANESTHETIC



- Proven topical skin refrigerant
- Temporary numbing lasts up to one minute
- Non-returnable item
- 3.9 oz bottle

Fine Stream Spray

000104BX4	4/box	\$215.99	\$205.19	\$205.19	\$194.39
-----------	-------	----------	----------	----------	----------

000104	12/case	\$539.99	\$512.99	\$512.99	\$485.99
--------	---------	----------	----------	----------	----------

Medium Stream Spray

000103BX4	4/box	\$204.49	\$194.27	\$194.27	\$184.04
-----------	-------	----------	----------	----------	----------

000103	12/case	\$499.99	\$474.99	\$474.99	\$449.99
--------	---------	----------	----------	----------	----------

GEBAUER SPRAY AND STRETCH TOPICAL ANESTHETIC VAPOCOOLANT SPRAY



Gebauer's Spray and Stretch® is an instant topical anesthetic skin refrigerant that effectively manages myofascial pain and trigger point release when used in conjunction with the spray and stretch technique.

- Non-drug, non-flammable
- Decreases pain sensation to allow the muscle to be stretched
- Helps inactivate trigger points, relieve muscle spasms and reduce referred pain
- 3.9 oz bottle

Fine Stream Spray

000404		each	\$49.99	\$47.49	\$47.49	\$44.99
--------	--	------	---------	---------	---------	---------

COLTENE WHALEDENT® HYGENIC ENDO-ICE®

Endo-Ice® spray is applied to a cotton pellet for vitality check

- Much colder than ice or ethyl chloride at -26.2°C
- Includes directional spray extension
- Fresh spearmint scent
- No CFCs



6 oz Spray Can w/Applicator Tip

H05032		each	\$31.29	\$29.73	\$29.73	\$28.16
--------	--	------	---------	---------	---------	---------

Replacement Straws

60022994	100/pack	\$12.49	\$11.87	\$11.87	\$11.24
----------	----------	---------	---------	---------	---------

SAFELY STORE YOUR CONTROLLED SUBSTANCES!

BRANDT® DIGITAL NARCOTIC SAFE



- Heavy duty steel welded box in gray hammer-tone
- Door is suspended on a full piano hinge
- Provides protection for controlled substances
- Elock - All metal construction, motorized locking engagement
- Easy administration with smartphone app
- Access with keypad or free lock app
- Elock "tiles" show locked or unlocked state
- Battery meter
- Multi-use pass code, e.g. company employees
- View list of active codes and add/delete users at any time
- Data logs keep record of the last 600 activities
- Rechargeable lithium battery
- Bluetooth data transmission secured with AES encryption to prevent hacking
- 10" x 15" x 6.5" with two adjustable shelves

80000E	each	\$570.29	\$541.78	\$541.78	\$513.26
--------	------	----------	----------	----------	----------

WATERLOO HEALTHCARE NARCOTIC CABINET



Double Door and Double Lock

- Doors suspended on full-length piano hinges
- Includes two adjustable shelves
- Can be set into the wall or surface mounted
- 15"H x 12"W x 8"D

MPH01WNCS2	White	each	\$350.79	\$333.25	\$333.25	\$315.71
------------	-------	------	----------	----------	----------	----------

GRAHAM-FIELD NARCOTIC SAFE



Double Door and Double Lock

- Double-door protection for narcotics
- A different key for each door
- Doors suspended on full-length piano hinges
- Keys cannot be removed when lock is in open position
- Made of high quality steel with welded construction
- Can be set into the wall or surface mounted
- ⚠️ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov
- 14.5"H X 9.5"W X 4" D

GF3007	each	\$335.79	\$319.00	\$319.00	\$302.21
--------	------	----------	----------	----------	----------

WATERLOO HEALTHCARE MEDICAL CARTS



Features Include:

- White plastic top with ergonomic cutaway front
- Standard drawer handles
- 5" Premium easy roll sealed bearing casters with precision bearing swivels
- Thread guards minimize sutures and other materials from collecting around axles
- Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability/protection
- Single piece construction resists scuffs/scratches
- Dual ergonomic handles
- Smooth ball bearing slides with auto-close mechanism
- Single pull-out shelf provides additional work surface
- Pre-threaded holes to assist in mounting some accessories
- Steel drawers
- Overall dimensions: 23"W x 24.5"D x 36"H

3 Drawer: 2 - 3"H and 2 - 9"H Drawers, Dimensions 16.5"W x 16.5"D

MPH01W3	Key lock, light gray	each	\$1,358.49	\$1,290.57	\$1,290.57	\$1,222.64
---------	----------------------	------	------------	------------	------------	------------

4 Drawer: 2 - 3"H, 1 - 6"H, and 1 - 9"H Drawers, Dimensions 16.5"W x 16.5"D

MPH01W4	Key lock, light gray	each	\$1,461.49	\$1,388.42	\$1,388.42	\$1,315.34
---------	----------------------	------	------------	------------	------------	------------

HEALTH CARE LOGISTICS® NARCOTIC CABINET



Double Door and Double Lock

- Double walled outer door with a cam lock and panel inner door with a separate cam lock
- Two cam locks are keyed differently
- Two keys per lock are included: An original key is needed to make duplicates
- Mount to the wall by two holes in the backside or set on or under countertops
- Made of 20 gauge steel; Beige enamel finish
- Includes one adjustable shelf
- Inside dimensions without shelf: 11-3/4"W x 13-1/2"H x 5-3/4"D

3704	each	\$303.49	\$288.32	\$288.32	\$273.14
------	------	----------	----------	----------	----------

PHARMACEUTICALS

1

HEALTH CARE LOGISTICS® RED PULL-TIGHT SEALS



Numbered

7820	100/box	\$26.29	\$24.98	\$24.98	\$23.66
------	---------	---------	---------	---------	---------

Without Numbers

7821	100/box	\$17.79	\$16.90	\$16.90	\$16.01
------	---------	---------	---------	---------	---------

HEALTH CARE LOGISTICS® EXEMPT NARCOTIC LOG BOOK, SALES AND PURCHASES

- 8-3/8" length x 5-3/8" height
- Light blue paper cover log book
- Total of 100 pages to record on (50 sheets front and back)
- Information is the same on all pages, is sectioned off and has lines to record the following: Date, Preparation Name, Amount, Signature of Purchaser, Address and Disp. By



11654	each	\$12.49	\$11.87	\$11.87	\$11.24
-------	------	---------	---------	---------	---------

LMA® MAD NASAL™ INTRANASAL MUCOSAL ATOMIZATION DEVICE

When faced with a daunting variety of tasks, you shouldn't have to worry about needle stick injuries. MAD Nasal can come to the rescue delivering safe, painless, and rapidly-effective treatment with minimal resource utilization.

With no danger of needle sticks, no need to sterilize, and no cumbersome and time-consuming IV setup, this intranasal mucosal atomization device will deliver a mist of atomized medication that is absorbed directly into your patient's bloodstream, so you spend more time caring for your patient, and less time starting a painful intravenous line.

- Latex-free
- With luer lock connection

Minimal Resource Utilization

- Nasal drug administration is quick and easy, eliminating IV set-up time
- No sterile technique required

Controlled Administration

- Exact dosing, exact volume
- Effective in both children and adults
- Titratable to effect
- Atomizes in any position

Rapidly Effective

- Atomized nasal medications absorb directly into bloodstream, avoiding first-pass metabolism
- Atomized nasal medications absorb directly into the brain and cerebrospinal fluid via olfactory mucosa to nose-brain pathway
- Achieves medication levels comparable to injections

Safe and Painless

- No needle, no shot, no pain, and no risk of needle stick



Without Syringe

MAD300	25/box	\$208.79	\$198.35	\$198.35	\$187.91
--------	--------	----------	----------	----------	----------

With 3ml Syringe

MAD100	25/box	\$237.79	\$225.90	\$225.90	\$214.01
--------	--------	----------	----------	----------	----------



LMA MADGIC™ MALLEABLE TOPICAL ANESTHETIC DEVICE

LMA MADgic is a single-use atomization device. Its unique design has a flexible cannula with memory and a tip that atomizes drugs into a fine mist of particles 30-100 microns in size.

Topical anesthetic delivered to the mucosal surface of the larynx may encourage greater hemodynamic stability and reduced airway stimulation in patients:

- Reduce heart rate vs pre-induction values
- Reduced cardiovascular response
- Easier intubation vs saline
- Fewer attempts for successful insertion of laryngeal mask



With 3ml Syringe

MAD600	25/box	\$363.79	\$345.60	\$345.60	\$327.41
--------	--------	----------	----------	----------	----------

PHILIPS RESPIRONICS OPTICHAMBER DIAMOND HOLDING CHAMBER

For use with inhaled steroids, as well as long-acting and rescue pressurized valved metered dose inhaler (pMDI) medications, OptiChamber Diamond's intuitive design enhances medication delivery and compliance for patients of all ages.



OptiChamber Diamond

1077478	each	\$22.29	\$21.18	\$21.18	\$20.06
---------	------	---------	---------	---------	---------

OptiChamber Diamond w/Medium LiteTouch Mask

1079826	each	\$22.99	\$21.84	\$21.84	\$20.69
---------	------	---------	---------	---------	---------

HEALTH CARE LOGISTICS® DISPOSABLE AMPULE BREAKER

- Reduces the chances of finger cuts and cross contamination.
- Ampules can be safely opened by snapping the top of the ampule breaker to the side; thus breaking the top of the ampule off at the neck of the bottle. The top of the ampule can then be discarded with the Ampule Breaker.



10273	100/box	\$18.99	\$18.04	\$18.04	\$17.09
-------	---------	---------	---------	---------	---------

BD® TWINPAK™ DUAL CANNULA DEVICE

Simplifying Needleless IV Access

A universal IV access device, offering twice the performance of traditional methods. It features two cannulas in one shield; one red hub cannula for syringe filling and one clear, blunt plastic cannula for safe split septum access. The TwinPak reduces the number of individual pieces, thereby saving time and reducing costs.



303390	100/box	\$51.49	\$48.92	\$48.92	\$46.34
--------	---------	---------	---------	---------	---------

TwinPak™ is also available with syringes. Call customer service for details.

MONOJECT™ 10 CC ORAL MEDICATION SYRINGE

- With special tip and ribbed cap
- Clear polypropylene
- Single-use



907102	100/box	\$51.29	\$48.73	\$48.73	\$46.16
--------	---------	---------	---------	---------	---------

HOSPIRA CARPUJECT® SYRINGE SYSTEM

- Easy loading cartridges and one-handed disposal




204902	each	\$3.99	\$3.79	\$3.79	\$3.59
--------	------	--------	--------	--------	--------

Q.I. MEDICAL, INC. SAFEFILL™ VIAL ADAPTERS

Safely fill syringes without needles

- Fits all multi-dose vials
- Use with other needleless systems
- Eliminates needle sticks when filling syringes
- SF200S locks onto vial for improved aseptic technique
- Large lumen speeds solution transfer
- More features than expensive dispensing pins and vial adapters



Universal; For Use on Single or Multi-Dose Drug Vials

SF100		50/case	\$71.99	\$68.39	\$68.39	\$64.79
-------	--	---------	---------	---------	---------	---------

For Use on Vials with 1/2" (13mm) Tops

SF200S	Locks onto vial	50/case	\$89.99	\$85.49	\$85.49	\$80.99
--------	-----------------	---------	---------	---------	---------	---------

B|BRAUN DISPENSING PINS - DEHP AND LATEX-FREE

Standard Spike - Luer Slip Connector


- For preparing and dispensing diluent or additive from multi-dose rubber-stoppered vials
- Utilizes a bacterial retentive air-venting filter



412013		each	\$2.99	\$2.84	\$2.84	\$2.69
--------	--	------	--------	--------	--------	--------

MINI-SPIKE® IV Additive-Proximal Luer Lock Connector


- For preparing and dispensing diluent or additive from multi-dose rubber-stoppered vials
- Utilizes a bacterial retentive air-venting filter



BB412012		50/case	\$132.29	\$125.68	\$125.68	\$119.06
----------	--	---------	----------	----------	----------	----------

MINI-SPIKE® Dispensing Pin with ULTRASITE®


- For rubber-stoppered vials
- Two-way valve for aspiration and injection
- Utilizes a bacterial retentive air-venting filter



412026		50/case	\$216.49	\$205.67	\$205.67	\$194.84
--------	--	---------	----------	----------	----------	----------

MINI-SPIKE® Dispensing Pin with ULTRASITE® Valve and Security Clip


- For attachment to standard 20mm rubber-stoppered vials
- Two-way valve for aspiration and injection
- Utilizes a bacterial retentive air-venting filter



413504		50/case	\$177.49	\$168.62	\$168.62	\$159.74
--------	--	---------	----------	----------	----------	----------

MICRO PIN® Single Use Dispensing Pin with Luer Lock Connector

- For withdrawal or injection of medication from rubber-stoppered vials
- Non-vented



415019		100/case	\$47.79	\$45.40	\$45.40	\$43.01
--------	--	----------	---------	---------	---------	---------

PHARMACEUTICALS

1

ACE MEDICATION LABELS

All medication containers must be labeled to ensure proper identification.


- Inexpensive
- Self-adhering
- Color-coded
- Variety of uses: syringes, med cups, and med bottles
- 700/roll



9150011	Anectine, fluorescent red	each	\$13.49	\$12.82	\$12.82	\$12.14
9150030	Atropine, green	each	\$13.49	\$12.82	\$12.82	\$12.14
9150003	Brevital®, yellow	each	\$13.49	\$12.82	\$12.82	\$12.14
9150041	Decadron®, white	each	\$13.49	\$12.82	\$12.82	\$12.14
9150024	Demerol, blue	each	\$13.49	\$12.82	\$12.82	\$12.14
9150014	Droperidol, salmon	each	\$13.49	\$12.82	\$12.82	\$12.14
9150015	Ephedrine, violet	each	\$13.49	\$12.82	\$12.82	\$12.14
9150110	Fentanyl, tan	each	\$13.49	\$12.82	\$12.82	\$12.14
9150006	Ketamine, yellow	each	\$13.49	\$12.82	\$12.82	\$12.14
9150009	Lidocaine®, gray	each	\$13.49	\$12.82	\$12.82	\$12.14
9150099	Marcaine®, gray	each	\$13.49	\$12.82	\$12.82	\$12.14
9150795	Midazolam, orange	each	\$13.49	\$12.82	\$12.82	\$12.14
9150027	Morphine, blue	each	\$13.49	\$12.82	\$12.82	\$12.14
9150016	Nalbuphine, white/black	each	\$13.49	\$12.82	\$12.82	\$12.14
9150036	Narcan, white	each	\$13.49	\$12.82	\$12.82	\$12.14
9150044	Nubain, blue	each	\$13.49	\$12.82	\$12.82	\$12.14
9150018	Ondansetron, white	each	\$13.49	\$12.82	\$12.82	\$12.14
9151010	Propofol, yellow	each	\$13.49	\$12.82	\$12.82	\$12.14
9150108	Robinul, green	each	\$13.49	\$12.82	\$12.82	\$12.14
9150017	Romazicon, white/red	each	\$13.49	\$12.82	\$12.82	\$12.14
9150025	Sublimaze, blue	each	\$13.49	\$12.82	\$12.82	\$12.14
9150013	Succinylcholine, fluorescent red	each	\$13.49	\$12.82	\$12.82	\$12.14
9150029	Valium, orange	each	\$13.49	\$12.82	\$12.82	\$12.14
9150104	Versed, orange	each	\$13.49	\$12.82	\$12.82	\$12.14
9150008	Xylocaine, gray	each	\$13.49	\$12.82	\$12.82	\$12.14

ACE PLAIN LABELS

- Labels are 1/4" x 1-1/4"
- 700/roll



9150062	White	each	\$13.49	\$12.82	\$12.82	\$12.14
9150063	Light blue	each	\$13.49	\$12.82	\$12.82	\$12.14
9150066	Gray	each	\$13.49	\$12.82	\$12.82	\$12.14
9150039	Pumpkin	each	\$13.49	\$12.82	\$12.82	\$12.14
9150065	Fluorescent pink	each	\$13.49	\$12.82	\$12.82	\$12.14
9150068	Fluorescent red	each	\$13.49	\$12.82	\$12.82	\$12.14
9150069	Fluorescent green	each	\$13.49	\$12.82	\$12.82	\$12.14


UNITED LAB LABEL LARGE LABEL DISPENSER



- For widths up to 8-3/4"
- Holds 10 rolls

ULTLD1	each	\$171.99	\$163.39	\$163.39	\$154.79
--------	------	----------	----------	----------	----------

ACE IDENTIFICATION LABELS



"Allergic To" Labels

- Shown actual size
- 7/8" x 3-1/2", fluorescent pink
- 250/roll

9150083	each	\$17.49	\$16.62	\$16.62	\$15.74
---------	------	---------	---------	---------	---------



Biohazard Labels

- Shown actual size
- 2" x 2", orange/black

0030041	25/pack	\$8.49	\$8.07	\$8.07	\$7.64
---------	---------	--------	--------	--------	--------

PHARMACEUTICALS

1

PDC HEALTHCARE ANESTHESIA DRUG SYRINGE TAPE LABELS

- User applied drug labels are available for accurate and reliable labeling of all drug syringes
- Each 1/2" x 1-1/2" tape label is color coded by drug class and a perforation is provided between each label

AN1	Atropine	Green
AN41	Brevital 10mg/ml, Date	Yellow
AN46	Demerol	Powder blue
AN132	Dexamethasone	White
AN155	Diprivan/Propofol	Yellow
AN7	Fentanyl	Blue
AN9	Glycopyrrolate	Green
AN59	Ketalar	Yellow
AN11	Lidocaine	Gray
AN120	Midazolam	Orange
ANS14	Naloxone	White/blue
AN55	Ondansetron	White
AN27	Propofol	Yellow
AN88	Robinul	Green
ANS19	Romazicon	White/red
AN20	Succinylcholine	Fluorescent red
AN101	Valium	Orange
AN149	Versed	Orange

333/roll each **\$7.49** **\$7.12** **\$7.12** **\$6.74**

ANWC	Clindamycin	White/black
ANWK	Ketorolac	White
ANW	Nalbuphine	White/black

500/roll each **\$16.29** **\$15.48** **\$15.48** **\$14.66**

LAN60	Ketamine	Yellow
-------	----------	--------

600/roll each **\$15.29** **\$14.53** **\$14.53** **\$13.76**

**Additional labels available, call for more information.*

PDC HEALTHCARE PLASTIC MULTIPLE ROLL LABEL DISPENSER

- For dispensing multiple rolls
- Includes flip-top cover and bottom suction cups
- Dispenser size: 9-7/8"W x 3-3/8"H x 4-1/2"D



For Rolls with 1" Core

TM500P each **\$54.29** **\$51.58** **\$51.58** **\$48.86**

COVIDIEN PHARMASAFETY™ WASTE CONTAINERS



- Developed exclusively to meet state regulations for pharmaceutical waste. White with blue lids, they are easily distinguished from other disposal containers. Each has a leak-resistant gasket and absorbent pad to help contain liquid contents.
- Tamper evident label helps deter tampering during final disposal.
- Secure, hinged lid with temporary and final closure accepts larger, bulkier waste; closes temporarily between uses; locks for final disposal.

With Lid and Absorbent Pad

8820	2 Gallon with gasketed hinged lid	each	\$11.49	\$10.92	\$10.92	\$10.34
8836SA	3 Gallon with counterbalance lid	each	\$13.49	\$12.82	\$12.82	\$12.14

BECAUSE WE JUST CAN'T FIT EVERYTHING INTO THIS CATALOG!


SHOP ONLINE AT **ACESOUTHERN.COM**

ACE SOUTHERN
Surgical Solutions



Tier Pricing Details see page 2

Regular Price **Member Price** **Premier Price** **Premier+ Price**



KING VISION™ ABLADE™ VIDEO LARYNGOSCOPE

The King Vision combines the convenience of a durable, reusable, video display with an affordable disposable blade. Own the Airway™ at a price that's affordable for all of your intubations, difficult and routine.

Affordable

- The King Vision aBlade is an excellent solution to solve the issue of limited availability of advanced visualization solutions in the OR. The King Vision is designed to visualize the vocal cords indirectly, which can help minimize the need for tissue manipulation and improve the ability to see anterior airways.
- With King Vision, you can improve patient care by having high performance visualization capabilities for every intubation at an affordable price.

Portable

- King Vision is lightweight, self-contained, battery operated and water resistant. King Vision can be used continuously for about 80 minutes, making it highly portable.
- The aBlades are individually packaged. Each blade is disposable, eliminating concerns of cross contamination. King Vision offers the choice of two blade types depending on preferred tube delivery technique.

Durable

- King Vision and aBlades have been designed to be your primary tool for all intubations.
- Video output capability is compatible with external monitor and recording devices.
- The robust, full-color, non-glare display can resist repeated cleaning and comes with a 1-year warranty. The aBlade video adapter is reusable, designed with materials that can withstand repeated use.

King Vision Display Kit

aBlade Kit Contains: Display, Video Adapter, Case, Batteries, 4 Channeled aBlades and 2 Standard Size 3 aBlades

KVLBKIT3	King Vision™ aBlade™	kit	\$1,799.99	\$1,709.99	\$1,709.99	\$1,619.99
aBlade Video Adapters						
KLVA12	Size 1 and 2	each	\$855.79	\$813.00	\$813.00	\$770.21
KLVA34	Size 3	each	\$580.29	\$551.28	\$551.28	\$522.26
Standard aBlades						
KVLAB1	Size1	20/case	\$399.99	\$379.99	\$379.99	\$359.99
KVLAB2	Size 2	20/case	\$399.99	\$379.99	\$379.99	\$359.99
KVLAB3	Size 3	each	\$20.99	\$19.94	\$19.94	\$18.89
Channeled aBlades						
KVLAB2C	Size 2c	20/case	\$399.99	\$379.99	\$379.99	\$359.99
KVLAB3C	Size 3c	20/case	\$399.99	\$379.99	\$379.99	\$359.99

Size Recommendations:

For infants and children, the small video adapter 1-2 and associated blades are indicated:

- aBlade size 1 (infant patients <4 years)*
- aBlade size 2 and 2C (child patients 1-10 years)*

For larger children and adults, the large video adapter 3-4 and associated blades are indicated:

- aBlade size 3 and 3C (children >5 years and adults)*

* Size ranges are based on teeth-to-vallecula data for the given ages, tracheal tube sizing recommendations and actual experience. A medical professional must evaluate on a patient-by-patient basis.

AVANTE INTEGRA SP VSO2 PORTABLE ANESTHESIA MACHINE

The Gold Standard in Portable Anesthesia

The Avante Integra SP VSO2 is a fully-functional, compact and portable anesthesia system in a lightweight package. The Integra SP VSO2 anesthesia system is compact enough for the smallest surgery or imaging suite, and is fully portable.

Easily transport the Avante Integra SP VSO2 by using the retractable handle, located conveniently on top of the unit. Side-located gas ports allow clinicians to safely and efficiently set up the machine.

In keeping with Avante's commitment to quality, all Integra SP anesthesia systems include a two-year manufacturer's warranty. The Avante Integra SP VSO2 is available with a hardshell carrying case, complete with casters and a retractable handle. A convenient mobile trolley with a sturdy built-in shelf is also available.

Features

- Small and compact unit ideal for office-based anesthesia
- Lightweight machine weighs only 26 pounds!
- Three gas machine (O₂, N₂O and Air)
- Mechanical anti-hypoxic device and air/N₂O interlock
- Retractable top handle for easy transport
- Conveniently located side gas ports
- Selectatec®-compatible vaporizer mounting system
- Optional carrying case with retractable handle and casters
- Optional mobile trolley with shelf
- Holds 1 vaporizer
- Specifically designed to accommodate low flow
- Auxiliary O₂ flowmeter port
- Recessed flush valve
- 2-year warranty



00101T2RS	Avante Integra SP VSO2	each	\$13,571.49	\$12,892.92	\$12,892.92	\$12,214.34
-----------	------------------------	------	-------------	-------------	-------------	-------------

Accessories

15000CASE	SP VSO2 Carrying Case	each	\$1,029.49	\$978.02	\$978.02	\$926.54
00101TAC	SP VSO2 Mobile Cart	each	\$808.99	\$768.54	\$768.54	\$728.09

**A NEW NAME, AN ALLIANCE,
WITH A PROVEN HISTORY
THAT YOU CAN TRUST.**

VISIT ACESOUTHERN.COM TODAY



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



FLEXICARE® BRITEPRO DISPOSABLE FIBEROPTIC HANDLE SYSTEM

BritePro is a single patient use Fiberoptic handle system which is comprised of a disposable handle and a reusable light source utilizing the latest in LED technology, providing clinicians with a complete system.

- Significantly brighter light than conventional blades or halogen handles
- Long bulb life - Never change a bulb again; LED bulbs last for the life of the handle with no replacement
- The main body of the handle incorporates a green textured grip which provides excellent control and feel even when wearing gloves
- LED bulb uses less energy and therefore batteries last more than 10 times longer
- A freely rotating metal hinge allows the blade to be inserted easily while providing a positive locking mechanism when moved into the closed position
- A white LED light provides clear visualization of the vocal chords

BriteBlade Pro, Macintosh


040711U	Size 1	each	\$8.49	\$8.49	\$7.64	\$7.64
040712U	Size 2	each	\$8.49	\$8.49	\$7.64	\$7.64
040713U	Size 3	each	\$8.49	\$8.49	\$7.64	\$7.64
040714U	Size 4	each	\$8.49	\$8.49	\$7.64	\$7.64

BriteBlade Pro, Miller

040721U	Size 1	each	\$8.49	\$8.49	\$7.64	\$7.64
040722U	Size 2	each	\$8.49	\$8.49	\$7.64	\$7.64
040723U	Size 3	each	\$8.49	\$8.49	\$7.64	\$7.64
040724U	Size 4	each	\$8.49	\$8.49	\$7.64	\$7.64

Accessories

040791U	Handle	each	\$13.29	\$13.29	\$11.96	\$11.96
040792U	Light source	each	\$80.79	\$80.79	\$72.71	\$72.71



FLEXICARE® FIBEROPTIC LARYNGOSCOPE HANDLES AND BULB

All handles are bi-polar, allowing the laryngoscope to function irrespective of which way the batteries are inserted. The handles can also be opened from both ends (top and bottom) in order to remove swelled batteries with ease.

- Rough Finish
- Batteries not included: Standard = two C batteries, Mini = two AA batteries

040812U	Mini size	each	\$88.29	\$88.29	\$79.46	\$79.46
040811U	Standard size	each	\$88.29	\$88.29	\$79.46	\$79.46
041601U	Halogen bulb	each	\$23.99	\$22.79	\$22.79	\$21.59

FLEXICARE® FIBEROPTIC LARYNGOSCOPE BLADES

The Flexicare range of Fiberoptic laryngoscopes are manufactured without an inner sleeve, thus producing a lightweight handle. The halogen bulb produces a cold, yet bright illumination, providing an intensity of light required during examination and intubation.

- High quality stainless steel
- Fully autoclavable

Macintosh Fiberoptic Blade

040831U	Size 1	each	\$87.29	\$87.29	\$78.56	\$78.56
040832U	Size 2	each	\$87.29	\$87.29	\$78.56	\$78.56
040833U	Size 3	each	\$87.29	\$87.29	\$78.56	\$78.56
040834U	Size 4	each	\$87.29	\$87.29	\$78.56	\$78.56

Miller Fiberoptic Blade

040851U	Size 1	each	\$85.79	\$85.79	\$77.21	\$77.21
040852U	Size 2	each	\$85.79	\$85.79	\$77.21	\$77.21
040853U	Size 3	each	\$85.79	\$85.79	\$77.21	\$77.21
040854U	Size 4	each	\$85.79	\$85.79	\$77.21	\$77.21

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



WELCH ALLYN® STANDARD FIBEROPTIC LARYNGOSCOPE HANDLES

- Single HPX lamp cartridge design used across all sizes of 2.5v handles.
- Redesigned to reduce total weight, Welch Allyn's handles are up to 40% lighter than our previous design, which means increased responsiveness and maneuverability.
- Knurled texture ensures a secure grip.
- Batteries Included: Penlight = two AA batteries
Medium = two C batteries

60814	Penlight	each	\$107.29	\$107.29	\$96.56	\$96.56
60813	Medium	each	\$107.29	\$107.29	\$96.56	\$96.56

Halogen Replacement Lamp

06000U	2.5V	each	\$25.99	\$24.69	\$24.69	\$23.39
--------	------	------	--------------------	--------------------	--------------------	--------------------

WELCH ALLYN® STANDARD FITTING FIBEROPTIC LARYNGOSCOPE BLADES

- Removable, repolishable light pipes make these the longest-lasting fiberoptic laryngoscope blades in the world.
- Single-piece stainless steel blades ensure blade integrity (no soldered joints).
- Fits Welch Allyn Fiberoptic Handle

Macintosh Fiberoptic Blades

69062	Size 2	each	\$149.29	\$149.29	\$134.36	\$134.36
69063	Size 3	each	\$149.29	\$149.29	\$134.36	\$134.36
69064	Size 4	each	\$149.29	\$149.29	\$134.36	\$134.36

Miller Fiberoptic Blades

68062	Size 2	each	\$149.29	\$149.29	\$134.36	\$134.36
68063	Size 3	each	\$149.29	\$149.29	\$134.36	\$134.36



WELCH ALLYN® FIBEROPTIC LARYNGOSCOPE SET

- Batteries included
- 5-year warranty

Welch Allyn Fiberoptic Laryngoscope Set Includes:

Medium and Penlight Handles, Mac Blades (1,2,3, and 4), and Case

69696	kit	\$902.29	\$902.29	\$812.06	\$812.06
-------	-----	---------------------	---------------------	---------------------	---------------------



FLEXICARE® FIBEROPTIC LARYNGOSCOPE BOX SETS

- Standard size handle, knurled
- Batteries not included: two C batteries



FLEXICARE® Box Sets Includes:

Standard Handle, Macintosh Blades (see below), and Case

040860U	1, 2 & 3	kit	\$414.49	\$414.49	\$373.04	\$373.04
040862U	2, 3 & 4	kit	\$414.49	\$414.49	\$373.04	\$373.04
040861U	1, 2, 3 & 4	kit	\$515.79	\$515.79	\$464.21	\$464.21

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



RÜSCH STANDARD/CONVENTIONAL CHROME PLATED BRASS HANDLES

- Non-sterile
- Batteries not included: Small Penlight = two AA batteries
Standard = two C batteries

8620	Small/Penlight	each	\$30.79	\$30.79	\$27.71	\$27.71
8621	Med/Standard	each	\$30.79	\$30.79	\$27.71	\$27.71

RÜSCH LARYNGOSCOPE BLADES

Standard/Conventional

- Non-sterile
- Stainless steel
- Lamp supplied with blade

Macintosh Style Laryngoscope Blades

601100	1 Infant, 66mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
602200	2 Child, 85mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
603300	3 Med. adult 110mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
603350	3.5 Med. adult 121mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
604400	4 Large adult, 140mm*	each	\$49.99	\$49.99	\$44.99	\$44.99

Miller Style Laryngoscope Blades

615050	0 Premature 53mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
616100	1 Infant, 80mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
617200	2 Child, 130mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
618300	3 Med. adult, 173mm*	each	\$49.99	\$49.99	\$44.99	\$44.99
619400	4 Large adult 180mm*	each	\$49.99	\$49.99	\$44.99	\$44.99

* From inside base to tip

RÜSCH FROSTED REPLACEMENT BULBS

- Non-sterile

Fits: Macintosh 1 & 2; Miller 00, 0, 1, 2 & 3; Wisconsin 1

623100	Small	6/box	\$24.29	\$23.08	\$23.08	\$21.86
--------	-------	-------	---------	---------	---------	---------

Fits: Macintosh 3, 3.5, & 4; Miller 4; Wisconsin 2, 3 & 4

622300	Large	6/box	\$24.29	\$23.08	\$23.08	\$21.86
--------	-------	-------	---------	---------	---------	---------



FLEXICARE® STANDARD REUSABLE LARYNGOSCOPE HANDLES & BULBS

All handles are bipolar, allowing the laryngoscope to function irrespective of which way the batteries are inserted. The handles can also be opened from both ends (top and bottom) in order to remove swelled batteries with ease.

- Knurled Finish Handles
- Batteries not included: Small Penlight = two AA batteries
Standard = two C batteries

040022U	Small/Penlight	each	\$37.79	\$37.79	\$34.01	\$34.01
040021U	Med/Standard	each	\$37.79	\$37.79	\$34.01	\$34.01



Kryton Bulbs

041401U	Small 00, 0, 1	each	\$12.79	\$12.15	\$12.15	\$11.51
041402U	Standard 2, 3, 4	each	\$12.79	\$12.15	\$12.15	\$11.51

FLEXICARE® STANDARD LARYNGOSCOPE BLADES

The standard range of laryngoscopes are manufactured from high quality stainless steel. The light bulb in the blade can be easily changed and the distance from the bulb to the blade tip will protect the patient from accidental burns, while providing a high level of illumination.

Macintosh Blades


040110U	Size 1	each	\$35.49	\$35.49	\$31.94	\$31.94
040111U	Size 2	each	\$35.49	\$35.49	\$31.94	\$31.94
040112U	Size 3	each	\$35.49	\$35.49	\$31.94	\$31.94
040113U	Size 4	each	\$35.49	\$35.49	\$31.94	\$31.94

Miller Blades

040130U	Size 0	each	\$35.49	\$35.49	\$31.94	\$31.94
040131U	Size 1	each	\$35.49	\$35.49	\$31.94	\$31.94
040132U	Size 2	each	\$35.49	\$35.49	\$31.94	\$31.94
040133U	Size 3	each	\$35.49	\$35.49	\$31.94	\$31.94
040134U	Size 4	each	\$35.49	\$35.49	\$31.94	\$31.94

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



WELCH ALLYN® STANDARD LARYNGOSCOPE HANDLES

- Universal design allows the handle to accept any ASTM standard blade.
- Crafted from durable chrome-plated brass; knurled finish for sure grip.
- Batteries Included: Medium = two C batteries

60300	Medium	each	\$75.79	\$75.79	\$68.21	\$68.21
-------	--------	------	---------	---------	---------	---------

WELCH ALLYN® STANDARD LARYNGOSCOPE BLADES

- Blades easily convert from lamp to Fiberoptic illumination, making upgrading economical (Mac, E. Mac and Miller only).
- Repairable, allowing replacement of damaged electrical components for a fraction of the cost of a replacement blade (Mac, E. Mac and Miller only).
- Single-piece stainless steel blades ensure blade integrity (no soldered joints) (Mac, E. Mac and Miller only).

Macintosh Blades - Curved Blade

69041	Size 1	each	\$140.49	\$140.49	\$126.44	\$126.44
69042	Size 2	each	\$140.49	\$140.49	\$126.44	\$126.44
69043	Size 3	each	\$140.49	\$140.49	\$126.44	\$126.44
69044	Size 4	each	\$140.49	\$140.49	\$126.44	\$126.44

Miller Blades - Straight Blade

68041	Size 1	each	\$140.49	\$140.49	\$126.44	\$126.44
68042	Size 2	each	\$140.49	\$140.49	\$126.44	\$126.44
68043	Size 3	each	\$140.49	\$140.49	\$126.44	\$126.44

WELCH ALLYN® REPLACEMENT LAMPS



Fits: Blade size 0, 1 & 2

04700U	Small	each	\$13.49	\$12.82	\$12.82	\$12.14
--------	-------	------	---------	---------	---------	---------

Fits: Blade size 3 & 4, Miller size 2

04800U	Large	each	\$13.49	\$12.82	\$12.82	\$12.14
--------	-------	------	---------	---------	---------	---------

SUNMED CONVENTIONAL LARYNGOSCOPE KIT



Macintosh American Profile Kit

- Blades constructed of 303/304 surgical stainless steel
- High impact plastic case for ease of transport
- Latex-free
- Batteries not included: two C batteries

Macintosh American Profile Kit Includes:

Medium Handle, 2 Extra Lamps, Macintosh Blades (1,2,3, and 4), and Case

5505244	kit	\$157.29	\$157.29	\$141.56	\$141.56
---------	-----	----------	----------	----------	----------

SUNMED CONVENTIONAL LARYNGOSCOPE KIT



Macintosh/Miller American Profile Kit

- Blades constructed of 303/304 surgical stainless steel
- High impact plastic case for ease of transport
- Latex-free
- Batteries not included: two C batteries

Macintosh/Miller American Profile Kit Includes: Medium Handle, 2 Extra Lamps, Macintosh Blades (2,3, and 4), Miller Blades (0,1,2, And 3), and Case

5506247	kit	\$230.79	\$230.79	\$207.71	\$207.71
---------	-----	----------	----------	----------	----------

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



**SMITHS MEDICAL PORTEX™
AIRCARE® ENDOTRACHEAL TUBES**

Ideal for short-term intubations, AIRCARE® tubes come in a wide selection of sizes, each with a standard 15mm ISO connector.

- Made of transparent PVC that softens at body temperature
- Radiopaque blue line x-ray visualization
- Oral or nasal intubation
- Sterile
- Sold individually
- Magill curve
- Latex-free
- Murphy eye
- Single-use

AIRCARE® Cuffed	each	\$3.79	\$3.60	\$3.60	\$3.41
AIRCARE® Uncuffed	each	\$2.99	\$2.84	\$2.84	\$2.69
AIRCARE® Cuffed w/ Preloaded Stylet	each	\$6.29	\$5.98	\$5.98	\$5.66
AIRCARE® Uncuffed w/ Preloaded Stylet	each	\$5.79	\$5.50	\$5.50	\$5.21

Size	AIRCARE® Cuffed	AIRCARE® Uncuffed	AIRCARE® Cuffed w/ Preloaded Stylet	AIRCARE® Uncuffed w/ Preloaded Stylet
3.0mm	100100030	100101030	100102030	100103030
3.5mm	100100035	100101035	100102035	100103035
4.0mm	100100040	100101040	100102040	100103040
4.5mm	100100045	100101045	100102045	100103045
5.0mm	100100050	100101050	100102050	100103050
5.5mm	100100055	100101055	100102055	-
6.0mm	100100060	100101060	100102060	-
6.5mm	100100065	100101065	100102065	-
7.0mm	100100070	100101070	100102070	-
7.5mm	100100075	100101075	100102075	-
8.0mm	100100080	100101080	100102080	-
8.5mm	100100085	100101085	100102085	-
9.0mm	100100090	100101090	100102090	-
10.0mm	100100010	100101010	100102010	-

BE THE OFFICE HERO.

Save Time and Money with...



REPLENISH

Our Replenish program allows you to set up scheduled deliveries for the products your practice uses on a regular basis. **Ensure that your practice never runs out of essential items that keep your office running smoothly.**

- **SHOP** from hundreds of eligible Replenish products
- **CHOOSE** your delivery frequency; every week, every month, you're in control
- **SAVE** an additional 3% on all Replenish-eligible products


See page 155 for complete details or visit acesouthern.com/replenish

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

**FLEXICARE®
ENDOTRACHEAL TUBES**

Magill Uncuffed	each	\$4.79	\$4.55	\$4.55	\$4.31
Murphy Uncuffed	each	\$4.79	\$4.55	\$4.55	\$4.55
Preformed, Oral Uncuffed	each	\$6.29	\$5.98	\$5.98	\$5.66
Preformed, Nasal Uncuffed	each	\$6.29	\$5.98	\$5.98	\$5.66

Size	Magill Uncuffed	Murphy Uncuffed	Preformed, Oral Uncuffed	Preformed Nasal Uncuffed
2.5mm	-	038961025U	-	-
3.0mm	038961030MU	038961030U	-	038964030U
3.5mm	038961035MU	-	-	-
4.0mm	038961040MU	038961040U	-	038964040U
4.5mm	038961045MU	038961045U	-	038964045U
5.0mm	038961050MU	038961050U	038962050U	038964050U
5.5mm	-	038961055U	038962055U	038964055U
6.0mm	-	038961060U	038962060U	038964060U
6.5mm	-	-	-	038964065U
7.0mm	-	-	038962070U	-



**SUNMED PREFORMED AGT RAE NASAL
ENDOTRACHEAL TUBES**

- Designed to direct tube over patient's forehead, reducing pressure on nares
- Allows circuit to be positioned out of surgical field; ideal for oral and maxillofacial surgery applications
- Polished murphy eye
- 15mm male fitting included
- Radiopaque strip
- Smooth beveled tip for atraumatic introduction
- Latex-free

For Cuffed Only:

- High volume, low pressure cuff
- Valve with sensitive pilot balloon indicating cuff status

Preformed AGT RAE, Nasal Cuffed	10/box	\$77.29	\$73.43	\$73.43	\$69.56
Preformed AGT RAE, Nasal Uncuffed	10/box	\$77.29	\$73.43	\$73.43	\$69.56

Size	Preformed AGT RAE Nasal Cuffed	Preformed AGT RAE Nasal Uncuffed
3.0mm	-	1732230
3.5mm	-	1732235
4.0mm	1732340	1732240
4.5mm	1732345	1732245
5.0mm	1732350	1732250
5.5mm	1732355	1732255
6.0mm	1732360	-
6.5mm	1732365	-
7.0mm	1732370	-
7.5mm	1732375	-
8.0mm	1732380	-



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



Cuffed, Murphy Eye-Nasal



Uncuffed, 2 Eyes-Nasal



Cuffed, Murphy Eye



Uncuffed, 2 Eyes



Cuffed, Murphy Eye



Uncuffed, Murphy Eye



Cuffed, Magill



Uncuffed, Magill

RÜSCH® ENDOTRACHEAL TUBES

- Made from thermo sensitive Ruschelit® PVC with clear depth markings and black positioning ring
- All tubes have atraumatic hooded tips and smooth edges
- Premium high volume/low pressure cuff and 15mm connectors are standard features
- Sterile, latex-free
- Single-use

AGT Nasal Preformed, Cuffed, Murphy Eye-Nasal	each	\$9.99	\$9.49	\$9.49	\$8.99
AGT Nasal Preformed, Uncuffed, 2 Eyes-Nasal	10/box	\$45.49	\$43.22	\$43.22	\$40.94
AGT Oral Preformed, Cuffed, Murphy Eye	each	\$8.79	\$8.35	\$8.35	\$7.91
AGT Oral Preformed, Uncuffed, 2 Eyes	each	\$6.29	\$5.98	\$5.98	\$5.66
Safety Clear® Plus, Cuffed, Murphy Eye	each	\$4.49	\$4.27	\$4.27	\$4.04
Safety Clear® Plus, Uncuffed, Murphy Eye	each	\$3.99	\$3.79	\$3.79	\$3.59
Safety Clear® Plus, Cuffed, Magill Eye	10/box	\$30.79	\$29.25	\$29.25	\$27.71
Safety Clear® Plus, Uncuffed, Magill Eye	each	\$4.79	\$4.55	\$4.55	\$4.31

Size	AGT Nasal Preformed Cuffed, Murphy Eye-Nasal	AGT Nasal Preformed Uncuffed, 2 Eyes-Nasal	AGT Oral Preformed Cuffed, Murphy Eye	AGT Oral Preformed Uncuffed, 2 Eyes	Safety Clear® Plus Cuffed, Murphy Eye	Safety Clear® Plus Uncuffed, Murphy Eye	Safety Clear® Plus Cuffed Magill	Safety Clear® Plus Uncuffed Magill
2.5mm	-	-	-	-	-	10038225	-	-
3.0mm	-	100181030	-	10018030	-	10038230	-	10038030
3.5mm	111781035	-	11178035	-	-	10038235	-	10038035
4.0mm	111781040	100181040	11178040	10018040	-	10038240	-	10038040
4.5mm	111781045	100181045	11178045	10018045	-	10038245	-	-
5.0mm	111781050	100181050	11178050	10018050	11208250	10038250	11208050	10038050
5.5mm	111781055	100181055	11178055	10018055	11208255	10038255	11208055	10038055
6.0mm	111781060	100181060	11178060	-	11208260	10038260	11208060	10038060
6.5mm	111781065	100181065	11178065	-	11208265	10038265	11208065	10038065
7.0mm	111781070	-	11178070	10018070	11208270	-	11208070	-
7.5mm	111781075	-	11178075	-	11208275	-	11208075	-
8.0mm	111781080	-	11178080	-	11208280	-	11208080	-
8.5mm	-	-	11178085	-	11208285	-	-	-
9.0mm	-	-	-	-	11208290	-	-	-
9.5mm	-	-	-	-	11208295	-	-	-
10.0mm	-	-	-	-	11208210	-	-	-

SALTER LABS® PARKER FLEX-TIP® ENDOTRACHEAL TUBES

Salter's endotracheal tubes are clinically preferred due to their innovative features. The clampless double-lumen valve system simplifies endo-bronchial intubation. The Parker Flex-Tip glides down the airway by flexing when it encounters airway structures. Flex tips are clinically proven to have faster intubations with higher first-pass success rates. Preloaded stylets improve audit compliance.

Oral/Nasal High Volume Low Pressure Cuffed	each	\$2.79	\$2.65	\$2.65	\$2.51
Oral/Nasal Uncuffed	each	\$3.49	\$3.32	\$3.32	\$3.14
Preformed, Nasal, High Volume Low Pressure Cuffed	*5/box	\$23.79	\$22.60	\$22.60	\$21.41
Preformed, Nasal, High Volume Low Pressure Cuffed	each	\$7.29	\$6.93	\$6.93	\$6.56
Preformed, Oral, High Volume Low Pressure Cuffed	each	\$7.29	\$6.93	\$6.93	\$6.56



Size	Oral/Nasal High Volume Low Pressure Cuffed	Oral/Nasal Uncuffed	Preformed, Nasal, High Volume Low Pressure Cuffed	Preformed, Oral, High Volume Low Pressure Cuffed
3.0mm	-	HPFU030	-	-
3.5mm	-	HPFU035	-	-
4.0mm	HPFHV40	HPFU040	HPFNC40 *	-
4.5mm	HPFHV45	HPFU045	HPFNC45 *	-
5.0mm	HPFHV50	HPFU050	HPFNC50 *	HPFOC50
5.5mm	HPFHV55	HPFU055	HPFNC55 *	HPFOC55
6.0mm	HPFHV60	HPFU060	HPFNC60	HPFOC60
6.5mm	HPFHV65	HPFU065	HPFNC65	HPFOC65
7.0mm	HPFHV70	HPFU070	HPFNC70	HPFOC70
7.5mm	HPFHV75	-	HPFNC75	HPFOC75
8.0mm	HPFHV80	-	-	HPFOC80
8.5mm	HPFHV85	-	-	-
9.0mm	HPFHV90	-	-	-
10.0mm	-	-	-	-



* Indicates that this item is sold 5 per box, not individually

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

RÜSCH® SLICK® DISPOSABLE ENDOTRACHEAL STYLETS

- Sureloc® Stop helps maintain desired length for patient safety
- Sterile, single-use
- Nylon-coated wire, malleable
- Handle forming capability for better grip

RS500	32cm, 6 Fr (tube sizes 2.0–3.5)	each	\$7.49	\$7.12	\$7.12	\$6.74
RS750	42cm, 8 Fr (tube sizes 4.0–6.5)	each	\$7.49	\$7.12	\$7.12	\$6.74
1000R	45.5cm, 10 Fr (tube sizes 7.0–10)	each	\$7.49	\$7.12	\$7.12	\$6.74

SUNMED STYLETS

9020701	Infant, malleable stainless, 6 Fr, 11.5"	each	\$22.99	\$21.84	\$21.84	\$20.69
9020410	Aluminum stylet, 10 Fr, 16"	each	\$7.79	\$7.40	\$7.40	\$7.01
9021005	Stilette stop for #9020410	each	\$4.79	\$4.55	\$4.55	\$4.31

ANESTHESIA ASSOCIATES STYLETS

00478	Copper & suction tube stylet, for intubation and suction, 14" (36cm)	each	\$29.99	\$28.49	\$28.49	\$26.99
-------	--	------	---------	---------	---------	---------



HUDSON ENDOTRACHEAL TUBE HOLDER

- Keeps an endotracheal tube in place without applying tape to patient's face.
- Tube may be attached to holder using either tape or safety pin
- Soft foam cushion protects patient's lips and face, cloth tie for securing holder & safety pin to tube
- 100% latex-free
- Individually packaged

1065	Endotracheal tube holder	each	\$5.79	\$5.50	\$5.50	\$5.21
------	--------------------------	------	--------	--------	--------	--------



WESTMED CO2 EASY™ CARBON DIOXIDE DETECTOR

- Built in bacterial filter protects both patient and clinician
- 3-year shelf life
- Designed for fast, easy airway assessment
- Verifies proper endotracheal tube placement
- Helps prevent unwarranted intubations
- Confirms cardiac output during CPR
- Expansive 2" diameter visual indicator makes rapid breath-to-breath response easy to see
- Anti-fogging agent reduces susceptibility to moisture and humidity, while ensuring continuous lens clarity



562135	Pediatric 1-15kg	each	\$12.79	\$12.15	\$12.15	\$11.51
562134	Adult, over 15kg	each	\$12.79	\$12.15	\$12.15	\$11.51

SUNMED TOOTH PROTECTOR

- Helps prevent broken or chipped teeth during incubation, bronchoscopy, laryngoscopy, tonsliectomy, and tooth extraction
- Minimizes tooth injury
- Latex-Free • Single Use • Disposable

9304103	Adult	50/box	\$103.29	\$98.13	\$98.13	\$92.96
---------	-------	--------	----------	---------	---------	---------



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

ACE ENDOTRACHEAL TUBE CONNECTORS

Straight - No port

0032186	4.0mm	each	\$25.99	\$24.69	\$24.69	\$23.39
0031005	5.0mm	each	\$25.99	\$24.69	\$24.69	\$23.39
003106	6.0mm	each	\$25.99	\$24.69	\$24.69	\$23.39
003107	7.0mm	each	\$25.99	\$24.69	\$24.69	\$23.39
003109	9.0mm	each	\$25.99	\$24.69	\$24.69	\$23.39

Right Angle - No port

ANE003125	5.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
0031206	6.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
003127	7.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44

Acute Angle - No port

003134	4.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
0031305	5.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
0031306	6.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
0031306	7.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
003138	8.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44
003139	9.0mm	each	\$40.49	\$38.47	\$38.47	\$36.44



SURGILUBE LUBRICANT

Fougera

020537	4-1/2 oz. tube	each	\$6.29	\$5.98	\$5.98	\$5.66
--------	----------------	------	--------	--------	--------	--------

Henry Schein

904953	3gm Foil Pack	144/box	\$27.49	\$26.12	\$26.12	\$24.74
--------	---------------	---------	---------	---------	---------	---------

SUNMED INTUBATION KIT WITH DISPOSABLE LARYNGOSCOPES




- Disposable
- Latex-Free

Kit Includes:

- 8 oropharyngeal color coded guedel airways (sizes 40mm - 100mm)
- 1 penlight chrome plated laryngoscope handle
- 1 medium chrome plated laryngoscope handle
- 5 miller laryngoscope blades (sizes 0-4)
- 4 macintosh laryngoscope blades (sizes 1-4)
- 4 extra laryngoscope lamps (2 small, 2 large)
- 1 adult magill forceps
- 1 pediatric magill forceps
- 2 "C" batteries
- 2 "AA" batteries
- 4 packets of lubricating jelly
- 2 endotracheal uncuffed tubes (sizes 2.5, 4mm)
- 8 endotracheal cuffed tubes (sizes 5-10mm)
- 4 nasopharyngeal airways (sizes 20, 24, 28, 32Fr)
- 1 10ml syringe
- 1 endotracheal tube introducer (bougie) with coude tip
- 1 adult stylette
- 1 nylon airway management case

8105094	SunMed intubation kit	\$447.79	\$425.40	\$425.40	\$403.01
---------	-----------------------	----------	----------	----------	----------



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



SALTER LABS® PARKER LARYNGEAL MASK DEVICE

Disposable Silicone Laryngeal Masks

- Gentle to the airway.
- The cuff's matt finish is designed to reduce the risk of the mask slipping out of place.
- Disposable mask eliminates costly reprocessing and the risk of cross-contamination.
- The smooth contouring of the 100% medical-grade silicone device enables easy insertion and is designed to produce an effective seal.
- The unique color-coded and printed pilot balloon allows for fast and correct size identification.
- The clear 15mm connector allows potential blockages at that site to be seen.
- Latex-free
- DEHP-free

HPLMD105	Infant, size 1	under 5kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD155	Infant, size 1.5	5kg-10kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD205	Infant, size 2	10kg-20kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD255	Child, size 2.5	20kg-30kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD305	Child, size 3	30kg-50kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD405	Adult, size 4	50kg-70kg	each	\$13.79	\$13.10	\$13.10	\$12.41
HPLMD505	Adult, size 5	70kg-100kg	each	\$13.79	\$13.10	\$13.10	\$12.41

SALTER LABS® PARKER REINFORCED LARYNGEAL MASK DEVICE

Disposable Silicone Laryngeal Masks

- Gentle to the airway.
- Spiral-wound reinforcing wire in the Parker flexible silicone tube reduces the risk of kinking.
- The smooth contouring of the laryngeal mask device enables easy insertion and is designed to produce an effective seal.
- The cuff's matt finish is designed to reduce the risk of the mask slipping out of place.
- Disposable mask eliminates costly reprocessing and the risk of cross-contamination.
- The unique color-coded and printed pilot balloon allows for fast and correct size identification.
- The clear 15mm connector allows potential blockages at that site to be seen.
- Latex-free
- DEHP-free

HPLMDR205	Infant, size 2	10kg-20kg	each	\$25.29	\$24.03	\$24.03	\$22.76
HPLMDR255	Child, size 2.5	20kg-30kg	each	\$25.29	\$24.03	\$24.03	\$22.76
HPLMDR305	Child, size 3	30kg-50kg	each	\$25.29	\$24.03	\$24.03	\$22.76
HPLMDR405	Adult, size 4	50kg-70kg	each	\$25.29	\$24.03	\$24.03	\$22.76
HPLMDR505	Adult, size 5	70kg-100kg	each	\$25.29	\$24.03	\$24.03	\$22.76

AMBU® AURAONCE™ LARYNGEAL MASK AIRWAYS

The Ambu® AuraOnce™ provides the safety and convenience of a sterile, single-use product, yet is remarkably cost-effective to use while the Auro40™ gives you the same features in an autoclavable version.

Benefits

- Unique curved airway tube replicates human anatomy
- Reinforced tip ensures correct placement
- Extra-soft cuff provides better seal with less internal pressure
- No sharp ridges that can scratch airway walls
- Cuff, mask, and airway molded in single unit
- Disposable
- Full range of sizes

Individual Sterile Peel-packs

321300000U	Sm. adult, size 3	30kg-50kg	each	\$34.49	\$32.77	\$32.77	\$31.04
321400000U	Adult, size 4	50kg-70kg	each	\$34.49	\$32.77	\$32.77	\$31.04
321500000U	Lg. adult, size 5	70kg-100kg	each	\$34.49	\$32.77	\$32.77	\$31.04



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

AMBU® AURA-I™ LARYNGEAL MASK AIRWAYS

The Ambu® Aura-I™ is the single-use supraglottic airway (SGA) that utilizes advanced design for superior performance during routine use. Industry leading safety features provide users with the tools needed for the unexpected.

- Easy and correct placement the first time. The airway tube is preformed to follow the anatomy of the human airway.
- Enhanced patient safety - The unique bite block helps prevent occlusion and ensure a stable position in the mouth throughout the procedure.
- Extra soft cuff ensures the best possible seal with least possible mucosal pressure.
- Convenient depth marks for monitoring correct position.
- Disposable
- Latex-free

329250000U	Child, size 2.5	20kg-30kg	each	\$32.29	\$30.68	\$30.68	\$29.06
329300000U	Sm. adult, size 3	30kg-50kg	each	\$32.29	\$30.68	\$30.68	\$29.06
329400000U	Adult, size 4	50kg-70kg	each	\$32.29	\$30.68	\$30.68	\$29.06
329500000U	Lg. adult, size 5	70kg-100kg	each	\$32.29	\$30.68	\$30.68	\$29.06



FLEXICARE® LARYSEAL™ BLUE LARYNGEAL MASK AIRWAYS

The first silicone single patient use laryngeal mask airway introduced to provide the benefits of the well-accepted reusable LarySeal™ with the added benefits of single use, protecting against cross contamination.

- Medical grade silicone
- Single-use
- Cost effective
- Crush Resistant - The main tube wall thickness reduces the risk of crushing.
- Color Coded Pilot Balloon - Ease of identification of sizes by color coded pilot balloons.
- Less Risk of Slipping or Dislocation - The silicone matt surface finish of the cuff provides a good adhesion to the oropharyngeal area and reduces the risk of movement and dislocation of the cuff.
- Universal Connection - 15mm M connector provides a universal connection.

03894230U	Child, size 3	30kg-50kg	each	\$15.29	\$14.53	\$14.53	\$13.76
03894240U	Adult, size 4	50kg-70kg	each	\$15.29	\$14.53	\$14.53	\$13.76
03894250U	Lg. adult, size 5	70kg-100kg	each	\$15.29	\$14.53	\$14.53	\$13.76



FLEXICARE® LARYSEAL™ CLEAR LARYNGEAL MASK AIRWAYS

LarySeal™ Clear is a single patient use PVC Laryngeal Mask Airway engineered using extra soft matt textured materials to replicate as closely as possible that of silicone, providing anesthetists the same feel as LarySeal™ Blue and LarySeal™ Multiple.

- Medical grade PVC
- Single-use
- Crush Resistant - The main tube wall thickness reduces the risk of crushing.
- Universal Connection - 15mm M connector provides a universal connection.

03894330U	Sm. adult, size 3	30kg-50kg	each	\$13.79	\$13.10	\$13.10	\$12.41
03894340U	Adult, size 4	50kg-70kg	each	\$13.79	\$13.10	\$13.10	\$12.41
03894350U	Lg. adult, size 5	70kg-100kg	each	\$13.79	\$13.10	\$13.10	\$12.41



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**FLEXICARE® LARYSEAL™
FLEXI LARYNGEAL MASK AIRWAYS**

Specially designed for use in dental and other head and neck surgeries. The construction of the main tube allows it to be positioned well away from the surgical field, improving surgical access without compressing the seal. The LarySeal™ Flexi main tube is reinforced to avoid kinking and eliminate the risk of airway tube occlusion.

- Reinforced - The main tube is reinforced with kink-free wire to eliminate the risk of airway tube occlusion
- Universal Connection - 15mm M connector provides a universal connection

03894530U	Child, size 3	30kg-50kg	each	\$23.99	\$22.79	\$22.79	\$21.59
03894540U	Adult, size 4	50kg-70kg	each	\$23.99	\$22.79	\$22.79	\$21.59
03894550U	Lg. adult, size 5	70kg-100kg	each	\$23.99	\$22.79	\$22.79	\$21.59

**FLEXICARE® LARYSEAL™
MULTIPLE LARYNGEAL MASK AIRWAYS**

A silicone reusable laryngeal mask airway designed to be sterilized and used up to a recommended 40 times. The LarySeal™ Multiple's cuff adapts to the contour of the oropharyngeal area to provide a secure seal.

- Medical grade silicone
- Autoclavable
- Crush Resistant - The main tube wall thickness reduces the risk of crushing.
- Color Coded Pilot Balloon - Ease of identification of sizes by color coded pilot balloons.
- Secure Seal for Airway Management - Highly elastic silicone material allows the cuff of the mask to form itself to the contour of the oropharyngeal area and provide a seal.
- Universal Connection - 15mm M connector provides a universal connection.

03894130U	Child, size 3	30kg-50kg	each	\$199.99	\$189.99	\$189.99	\$179.99
03894140U	Adult, size 4	50kg-70kg	each	\$199.99	\$189.99	\$189.99	\$179.99
03894150U	Lg. adult, size 5	70kg-100kg	each	\$199.99	\$189.99	\$189.99	\$179.99

**FLEXICARE® LARYSEAL™
COLOR CODED CUFF INFLATOR**

The LarySeal Cuff Inflator has clearly marked increments, each one specifically color coded to match the size color coding of the LarySeal pilot balloon, enabling easier and quicker identification of the correct measure for inflating the cuff of the LarySeal.

03894001U	Color coded cuff inflator	each	\$37.79	\$35.90	\$35.90	\$34.01
-----------	---------------------------	------	---------	---------	---------	---------



GET YOUR ORDER FASTER.

TO YOU
IN 2



Did you know 99% of ACE SOUTHERN

customers receive their order within 1-2 business days of shipment? ACE SOUTHERN has four order fulfillment facilities throughout the United States to ensure your order is **TO YOU IN 2***.

*See page 261 for complete details or visit acesouthern.com/shipping-information

ACE SOUTHERN

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

INTERSURGICAL I-GEL® SUPRAGLOTTIC AIRWAYS

The second generation supraglottic airway with a non-inflatable cuff for use in anesthesia.

Quick, easy and reliable to insert, i-gel is ideal for use as a routine airway in anesthesia as it provides high seal pressure and reduced trauma. It also incorporates a gastric channel to provide additional protection against aspiration and has the versatility to be applicable for use during difficult airway management, as a rescue device and a conduit for intubation (with fiberoptic guidance). Single-use.

8215	Infant, size 1.5	5kg-12kg	each	\$29.79	\$28.30	\$28.30	\$26.81
8202	Small pediatric, size 2	10kg-25kg	each	\$29.79	\$28.30	\$28.30	\$26.81
8225	Large pediatric, size 2.5	25kg-35kg	each	\$29.79	\$28.30	\$28.30	\$26.81
8203	Small adult, size 3	30kg-60kg	each	\$29.79	\$28.30	\$28.30	\$26.81
8204	Medium adult, size 4	50kg-90kg	each	\$29.79	\$28.30	\$28.30	\$26.81
8205	Large adult, size 5	90kg and over	each	\$29.79	\$28.30	\$28.30	\$26.81



2

COVIDIEN MALLINCKRODT COMBITUBE™ ESOPHAGEAL TRACHEAL DOUBLE LUMEN AIRWAYS

When a difficult airway must be established, time and accuracy are critical. You can be assured of rapid intubation and effective ventilation with the COMBITUBE airway. Unlike traditional endotracheal tubes, the COMBITUBE airway is designed to establish a non-surgical patient airway when placed into either the trachea or esophagus.

- Unique double-lumen design allows for rapid establishment of an effective airway through either esophageal or tracheal placement
- Blind placement eliminates the need for a laryngoscope
- Pharyngeal balloon inflates to hold device firmly in place and helps prevent the escape of gas through the nose or mouth
- Full-length lumen allows for suctioning of gastric contents with no interruption of patient ventilation should the COMBITUBE airway be placed in the esophagus
- Esophageal cuff inflates to seal the esophagus so that gas does not enter the stomach and gastric contents are not aspirated
- Appropriate for pre-hospital, surgical and emergency use

518537	For patients 4' to 5'6" tall	37 Fr	each	\$67.99	\$64.59	\$64.59	\$61.19
518541	For patients 5' tall and above	41 Fr	each	\$67.99	\$64.59	\$64.59	\$61.19



AMBU® KING LTS-D™ LARYNGEAL TUBES

The Ambu King LTS-D laryngeal tube is a disposable, simple to use alternative airway device that provides superior patient ventilation. The King LTS-D allows the passage of the gastric tube through a separate channel.

Designed with a straightened, beveled distal tip that assists in directing the airway posterior to the larynx and into the upper esophagus. Due to this unique configuration, there is minimal risk of the device entering the trachea.

A single port inflates both proximal and distal cuffs. The proximal cuff stabilizes the King LTS-D and seals the oropharynx. The distal cuff blocks entry into the esophagus, reducing the possibility of gastric insufflation. Ventilatory seal pressures of 30cm H₂O or more are achievable with the King LTS-D. Multiple distal ventilatory openings and bilateral ventilation eyelets facilitate air flow.

Key Benefits

- Seals in the esophagus and oropharynx to provide positive pressure ventilation
- Single inflation port
- Phthalate and latex-free material
- Drain tube for gastric and suction catheters
- Color-coded full range of sizes

KLTS420	Pediatric, size 0 - transparent	under 5kg	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS421	Pediatric, size 1 - white	5kg-12kg	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS422	Pediatric, size 2 - green	12kg-25kg	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS4225	Pediatric, size 2.5 - orange	25kg-35kg	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS423	Adult, size 3 - yellow	4' - 5' height	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS424	Adult, size 4 - red	5' - 6' height	each	\$43.49	\$41.32	\$41.32	\$39.14
KLTS425	Adult, size 5 - purple	over 6' height	each	\$43.49	\$41.32	\$41.32	\$39.14



Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**HUDSON RCI®
DUAL-CHANNEL AIRWAYS**

- Two-channel, Berman design
- Semi-rigid, non-toxic polyethylene



1152	#4, 40mm	each	\$1.49	\$1.42	\$1.42	\$1.34
1151	#6, 60mm	each	\$1.49	\$1.42	\$1.42	\$1.34
1150	#8, 80mm	each	\$1.49	\$1.42	\$1.42	\$1.34
1149	#9, 90mm	each	\$1.49	\$1.42	\$1.42	\$1.34
1148	#10, 100mm	each	\$1.49	\$1.42	\$1.42	\$1.34

**RÜSCH®
BERMAN AIRWAYS**

- Latex-free
- Single-use
- Non-sterile
- mm length printed on airway



121801	#4, 40mm, pink	each	\$0.79	\$0.75	\$0.75	\$0.71
121850	#5, 50mm, mint	each	\$0.79	\$0.75	\$0.75	\$0.71
121802	#6, 60mm, black	each	\$0.79	\$0.75	\$0.75	\$0.71
121870	#7, 70mm, white	each	\$0.79	\$0.75	\$0.75	\$0.71
121803	#8, 80mm, green	each	\$0.79	\$0.75	\$0.75	\$0.71
121804	#9, 90mm, yellow	each	\$0.79	\$0.75	\$0.75	\$0.71
121805	#10, 100mm, red	each	\$0.79	\$0.75	\$0.75	\$0.71
121806	#11, 110mm, orange	each	\$0.79	\$0.75	\$0.75	\$0.71

Emergency Pack

Includes: 1 of Each Size, Sizes 40-110mm

396274	pack	\$14.99	\$14.24	\$14.24	\$13.49
--------	------	---------	---------	---------	---------

**HUDSON RCI®
CATH-GUIDE® AIRWAYS**

- Flexible vinyl with a rigid bite block
- Guedel style with three internal channels



1171	55mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1170	60mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1169	70mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1168	80mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1164	90mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1167	100mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1166	110mm	each	\$1.99	\$1.89	\$1.89	\$1.79
1165	120mm	each	\$1.99	\$1.89	\$1.89	\$1.79

Emergency Pack

Includes: 1 of Each Size, Sizes 55mm - 120mm

1173	pack	\$53.49	\$50.82	\$50.82	\$48.14
------	------	---------	---------	---------	---------

**RÜSCH® LATEX-FREE
ROBERTAZZI AIRWAYS**

- Mediprene® construction
- Flexible for maximum patient comfort
- Soft material for enhanced nasal passageway comfort
- Rounded tip allows for gentle insertion
- Trumpet design for secure placement
- Individually packaged, sterile



123120	5.0mm, 20 Fr	light green	each	\$4.79	\$4.55	\$4.55	\$4.31
123122	5.5mm, 22 Fr	light green	each	\$4.79	\$4.55	\$4.55	\$4.31
123124	6.0mm, 24 Fr	black	each	\$4.79	\$4.55	\$4.55	\$4.31
123126	6.5mm, 26 Fr	black	each	\$4.79	\$4.55	\$4.55	\$4.31
123128	7.0mm, 28 Fr	dark green	each	\$4.79	\$4.55	\$4.55	\$4.31
123130	7.5mm, 30 Fr	dark green	each	\$4.79	\$4.55	\$4.55	\$4.31
123132	8.0mm, 32 Fr	orange	each	\$4.79	\$4.55	\$4.55	\$4.31
123134	8.5mm, 34 Fr	orange	each	\$4.79	\$4.55	\$4.55	\$4.31
123136	9.0mm, 36 Fr	red	each	\$4.79	\$4.55	\$4.55	\$4.31
			10/box	\$43.99	\$41.79	\$41.79	\$39.59

Emergency Pack

Includes: 1 of Each, Sizes 22-32 Fr, and Foil Packages of Lubricating Jelly

12312232	pack	\$52.29	\$49.68	\$49.68	\$47.06
----------	------	---------	---------	---------	---------

**FLEXICARE®
GUEDEL AIRWAYS**

- Integrated bite block
- Smooth bevelled tip
- Sterile



03892902U	#5, 50mm	blue	\$0.99	\$0.94	\$0.94	\$0.89
03892900U	#6, 60mm	black	\$0.99	\$0.94	\$0.94	\$0.89
03891910U	#7, 70mm	white	\$0.99	\$0.94	\$0.94	\$0.89
03891920U	#8, 80mm	green	\$0.99	\$0.94	\$0.94	\$0.89
03891930U	#9, 90mm	yellow *	\$0.99	\$0.94	\$0.94	\$0.89
03891940U	#10, 100mm	red	\$0.99	\$0.94	\$0.94	\$0.89

* formerly orange

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

RÜSCH® PVC NASAL AIRWAYS

- Soft PVC construction
- Inner striations for smooth suction catheter insertion
- Beveled tip for comfortable, gentle insertion
- Individually packaged, sterile



Item #	Size	Fr	Unit	Reg Price	Member Price	Premier Price	Premier+ Price
123318	4.0mm	18	each	\$5.99	\$5.69	\$5.69	\$5.39
123320	5.0mm	20	each	\$5.99	\$5.69	\$5.69	\$5.39
123322	5.5mm	22	each	\$5.99	\$5.69	\$5.69	\$5.39
123324	6.0mm	24	each	\$5.99	\$5.69	\$5.69	\$5.39
123326	6.5mm	26	each	\$5.99	\$5.69	\$5.69	\$5.39
123328	7.0mm	28	10/box	\$50.49	\$47.97	\$47.97	\$45.44
123330	7.5mm	30	each	\$5.99	\$5.69	\$5.69	\$5.39
123332	8.0mm	32	each	\$5.99	\$5.69	\$5.69	\$5.39
			10/box	\$50.49	\$47.97	\$47.97	\$45.44

RÜSCH® NASOPHARYNGEAL AIRWAYS

- Adjustable flange
- Sterile, disposable
- Latex-free
- Soft elastic finish
- Single-use



Item #	Size	Fr	Unit	Reg Price	Member Price	Premier Price	Premier+ Price
18542012	12 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542014	14 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542016	16 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542018	18 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542020	20 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542022	22 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542024	24 Fr		5/box	\$80.79	\$76.75	\$76.75	\$72.71
18542028	28 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542030	30 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81
18542032	32 Fr		5/box	\$80.79	\$76.75	\$76.75	\$72.71
18542034	34 Fr		5/box	\$80.79	\$76.75	\$76.75	\$72.71
18542036	36 Fr		each	\$19.79	\$18.80	\$18.80	\$17.81

FLEXICARE® NASOSAFE® NASOPHARYNGEAL AIRWAYS

- Swivel Safety Grip - The NasoSafe swivel grip protects the nasopharyngeal airway from loss in the nasal passage
- Reduced Risk of Injury During Intubations - The rounded edge of the tip combined with angle cut design helps to reduce the risk of injury during intubation
- Size Color-Coding - The swivel safety grip is color coded to match suction catheter size color-coding
- Patient Comfort - The swivel safety grip can be positioned so that it does not block the free nostril
- Reduced Risk of Patient Reaction - NasoSafe is made of medical grade PVC and is implant tested
- Ideal when oral passage is not accessible
- Allows access for nasotracheal suctioning



Item #	OD	ID (Fr)	Color	Unit	Reg Price	Member Price	Premier Price	Premier+ Price
03895004U	4.0mm	20	light green	each	\$6.29	\$5.98	\$5.98	\$5.66
03895045U	4.5mm	22	light green	each	\$6.29	\$5.98	\$5.98	\$5.66
03895005U	5.0mm	24	black	each	\$6.29	\$5.98	\$5.98	\$5.66
03895055U	5.5mm	26	black	each	\$6.29	\$5.98	\$5.98	\$5.66
03895006U	6.0mm	28	dark green	each	\$6.29	\$5.98	\$5.98	\$5.66
03895065U	6.5mm	28	dark green	each	\$6.29	\$5.98	\$5.98	\$5.66
03895007U	7.0mm	30	orange	each	\$6.29	\$5.98	\$5.98	\$5.66
03895075U	7.5mm	32	orange	each	\$6.29	\$5.98	\$5.98	\$5.66
03895008U	8.0mm	34	red	each	\$6.29	\$5.98	\$5.98	\$5.66
03895085U	8.5mm	36	red	each	\$6.29	\$5.98	\$5.98	\$5.66
03895009U	9.0mm	38	yellow	each	\$6.29	\$5.98	\$5.98	\$5.66
				10/box	\$59.99	\$56.99	\$56.99	\$53.99

FLEXICARE® NASOCLEAR™ NASOPHARYNGEAL AIRWAYS

NasoClear is a simple and effective airway adjunct designed for easy insertion and maximum patient compliance.

- Flared Nasal Opening - The large flared opening prevents the nasopharyngeal airway from physically traveling down the nasal passage but also allows for comfortable positioning without obstructing the other nostril
- Clear Construction - Provides visual confirmation before insertion that there are no blockages
- Reduced Risk of Patient Reaction - NasoClear is made of medical grade PVC and is implant tested
- Allows Access - Ideal when the oral passage is not accessible; It also allows for nasotracheal suctioning
- Reduced Risk of Injury During Intubation - Rounded bevelled tip reduces mucosal wall trauma and guides the airway along the floor of the nasal passage



Item #	OD	ID (Fr)	Unit	Reg Price	Member Price	Premier Price	Premier+ Price
03895004CU	4.0mm	20	each	\$3.79	\$3.60	\$3.60	\$3.41
03895045CU	4.5mm	22	each	\$3.79	\$3.60	\$3.60	\$3.41
03895005CU	5.0mm	24	each	\$3.79	\$3.60	\$3.60	\$3.41
03895055CU	5.5mm	26	each	\$3.79	\$3.60	\$3.60	\$3.41
03895006CU	6.0mm	28	each	\$3.79	\$3.60	\$3.60	\$3.41
03895065CU	6.5mm	28	each	\$3.79	\$3.60	\$3.60	\$3.41
03895007CU	7.0mm	30	each	\$3.79	\$3.60	\$3.60	\$3.41
03895075CU	7.5mm	32	each	\$3.79	\$3.60	\$3.60	\$3.41
03895008CU	8.0mm	34	each	\$3.79	\$3.60	\$3.60	\$3.41
03895085CU	8.5mm	36	each	\$3.79	\$3.60	\$3.60	\$3.41
03895009CU	9.0mm	38	each	\$3.79	\$3.60	\$3.60	\$3.41
			10/box	\$36.49	\$34.67	\$34.67	\$32.84

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

SMITHS MEDICAL PORTEX™ SINGLE-USE, DISPOSABLE ANESTHESIA BREATHING CIRCUITS

All Portex circuit systems deliver exceptional value by providing the configuration that meets your specific needs. With several standard circuit configurations and a variety of components in either adult or pediatric sizes, regular or expandable tubing, as well as adult and pediatric single-limb circuits, Portex has a circuit that fits your requirements.



40" Pediatric Corrugated Circuit

Gas sampling elbow, straight wye, breathing filter, 2L breathing bag

372007	20/case	\$131.79	\$125.20	\$125.20	\$118.61
--------	---------	----------	----------	----------	----------

60" Pediatric Corrugated Circuit with Gas Sampling Connectors

Gas sampling elbow, straight wye, 1L breathing bag

372904	15/case	\$95.99	\$91.19	\$91.19	\$86.39
--------	---------	---------	---------	---------	---------

40" - 100" Pediatric Expandable Circuit

Gas sampling elbow, straight wye, and 1L (latex-free) breathing bag

492804NL	20/case	\$215.29	\$204.53	\$204.53	\$193.76
----------	---------	----------	----------	----------	----------

33" - 87" Adult Expandable Circuit

Gas sampling elbow, straight wye, breathing filter, 3L (latex-free) breathing bag

650904	each	\$8.49	\$8.07	\$8.07	\$7.64
	20/case	\$155.79	\$148.00	\$148.00	\$140.21

40" Adult Circuits with Gas Sampling Connectors

Gas sampling elbow, straight wye, breathing filter, 3L (latex-free) breathing bag

670904	20/case	\$132.79	\$126.15	\$126.15	\$119.51
--------	---------	----------	----------	----------	----------

60" Adult Circuits with Gas Sampling Connectors

Gas sampling elbow, straight wye, breathing filter, 3L breathing bag

3718806	each	\$6.29	\$5.98	\$5.98	\$5.66
	20/case	\$88.79	\$84.35	\$84.35	\$79.91

FLEXICARE® ANESTHESIA BREATHING CIRCUITS



- Secure connections are welded together to prevent leakage and maintain connection
- High quality tubing made of stronger material to minimize risk of leakage with excellent drapability
- Soft, tactile breathing bag provides excellent feel with a secure connection to the machine
- Soft cushion anesthesia masks provide superior seal and patient comfort

96" Pediatric Circuit

Gas sampling elbow, breathing filter, 1L (latex-free) breathing bag

03802612U	20/box	\$140.49	\$133.47	\$133.47	\$126.44
-----------	--------	----------	----------	----------	----------

60" Adult Circuit

Gas sampling elbow, fixed wye, BV filter, 3L (latex-free) breathing bag

03801601U	20/box	\$145.99	\$138.69	\$138.69	\$131.39
-----------	--------	----------	----------	----------	----------

96" Adult Circuit

Gas sampling elbow, breathing filter, 3L (latex-free) breathing bag

03801612U	20/box	\$154.79	\$147.05	\$147.05	\$139.31
-----------	--------	----------	----------	----------	----------

SUNMED LATEX-FREE DISPOSABLE BREATHING BAGS



- Constructed of lightweight neoprene material, assuring optimum sensitivity to patients' breathing
- Individually inspected and packaged
- 22mm bushing

4104105	0.5 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104110	1 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104115	1.5 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104120	2 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104130	3 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104140	4 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89
4104150	5 Liter	each	\$20.99	\$19.94	\$19.94	\$18.89

PORTER BREATHING BAGS



- Latex-free
- Push on

41003NL	3 Liter breathing bag	each	\$38.29	\$36.38	\$36.38	\$34.46
---------	-----------------------	------	---------	---------	---------	---------

ANESTHESIA ASSOCIATES BREATHING BAG



- Black, with bushings
- Latex-free
- Fluted

E198	3 Liter breathing bag	each	\$49.29	\$46.83	\$46.83	\$44.36
LF12450	5 Liter breathing bag	each	\$90.79	\$86.25	\$86.25	\$81.71

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



TELEFLEX CORRUGATED BREATHING TUBES

- Can only be steam sterilized up to 270°
- 22mm female bushings

22675DP	22" Black	each	\$30.79	\$29.25	\$29.25	\$27.71
---------	-----------	------	---------	---------	---------	---------



ANESTHESIA ASSOCIATES BREATHING CIRCUIT ADAPTERS

- Flexible connector - 15mm male / 15mm female
- Fit patient end to facilitate positioning of tracheal tubes at any convenient angle
- Flexible metal adapters – chrome plated

00484		each	\$47.49	\$45.12	\$45.12	\$42.74
-------	--	------	---------	---------	---------	---------

SMITHS MEDICAL PORTEX™ FLEXTUBE



- 6.5" long
- 15mm I.D. one end; 22mm I.D. other end

002837		50/case	\$125.79	\$119.50	\$119.50	\$113.21
--------	--	---------	----------	----------	----------	----------

SMITHS MEDICAL PORTEX™ GAS SAMPLING ELBOW CONNECTOR, CLEAR, WITH CAP



- Clear patient connectors allow for easy monitoring during procedures
- Connectors have a 15mm I.D./22mm O.D patient fitting and a 15mm circuit fitting

3523804		50/case	\$55.99	\$53.19	\$53.19	\$50.39
---------	--	---------	---------	---------	---------	---------

SMITHS MEDICAL PORTEX™ BACTERIAL/VIRAL BREATHING FILTER



- Filter efficiency > 99.9%
- Latex-free

002862		40/case	\$47.29	\$44.93	\$44.93	\$42.56
--------	--	---------	---------	---------	---------	---------

MEDLINE STRAIGHT BV FILTER



- Designed to meet needs of respiratory and anesthesia applications
- Exceptional bacterial and viral filtration efficiency helps protect patients, anesthesia machines and mechanical ventilators
- Translucent housing for easy viewing of vomitus or other possible obstructions in system
- Latex-free

DYNJAABV1		50/case	\$66.79	\$63.45	\$63.45	\$60.11
-----------	--	---------	---------	---------	---------	---------

TELEFLEX® GIBECK HUMID-VENT® HEAT/MOISTURE EXCHANGER



- Transparent housing provides excellent visibility
- Very low resistance to flow
- 2 port (for tidal volumes of 150-1500ml)

13312		each	\$6.79	\$6.45	\$6.45	\$6.11
-------	--	------	--------	--------	--------	--------

SMITHS MEDICAL PORTEX™ THERMOVENT® 600 HEAT/MOISTURE EXCHANGER



Portex HMEs provide efficient heat and moisture recovery. All HMEs exhibit low resistance to flow and high moisture output which helps to minimize the loss of heat and moisture in ventilated or spontaneously breathing patients. All HMEs and filters have ISO standard tapered ends to assure a consistent leak-free seal.

- All Portex HMEs and filters are latex-free.
- For Tidal Volumes up to 600ml

580011		each	\$3.29	\$3.13	\$3.13	\$2.96
--------	--	------	--------	--------	--------	--------



SUNMED MASK HARNESS

- Webbed head portion for ventilation and pressure distribution over large area
- Conductive

105422	Child	each	\$17.79	\$16.90	\$16.90	\$16.01
--------	-------	------	---------	---------	---------	---------

105423	Adult	each	\$16.79	\$15.95	\$15.95	\$15.11
--------	-------	------	---------	---------	---------	---------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

SMITHS MEDICAL PORTEX™ DISPOSABLE ANESTHESIA BREATHING MASKS (DABM)

- Single-use
- Malleable cone conforms to patient facial features
- Clear, pre-inflated, anatomically shaped cushions for gas tight, low-pressure seals
- Color-coded hook rings on selected masks
- Inflation valve on Premium Soft Plus™ masks provide a custom fit



Premium™ No Hook Ring

5052	Infant	each	\$3.49	\$3.32	\$3.32	\$3.14
------	--------	------	--------	--------	--------	--------

Green Hook Ring

5070	Toddler	each	\$3.49	\$3.32	\$3.32	\$3.14
------	---------	------	--------	--------	--------	--------

Premium Soft Plus™

- With extra-soft cushion and inflation valve

Clear Hook Ring

005448	Child	each	\$3.79	\$3.60	\$3.60	\$3.41
--------	-------	------	--------	--------	--------	--------

Yellow Hook Ring

005455	Small adult	each	\$3.79	\$3.60	\$3.60	\$3.41
--------	-------------	------	--------	--------	--------	--------

Darl Blue Hook Ring

005445	Medium adult	each	\$3.79	\$3.60	\$3.60	\$3.41
--------	--------------	------	--------	--------	--------	--------

Light Blue Hook Ring

005465	Large adult	each	\$3.79	\$3.60	\$3.60	\$3.41
--------	-------------	------	--------	--------	--------	--------

FLEXICARE® VENTICAIRE® INFLATABLE ANESTHESIA FACE MASKS

- Inflatible Cushion - The seal can be enhanced by inflating/deflating the cushion to suit the patient's anatomical contour
- Ease of Observation - Clear face piece allows visual observation of patient status
- Soft Cushion - Smooth material used for the cushion allows for a tight seal while providing patient comfort
- Size Color-Coding - The hooks are color-coded to provide mask size identification



03853440SFPU	Size 0, Neonatal	each	\$4.79	\$4.55	\$4.55	\$4.31
03852430SFPU	Size 1, Small child	each	\$4.79	\$4.55	\$4.55	\$4.31
03852420SFPU	Size 2, Medium child	each	\$4.79	\$4.55	\$4.55	\$4.31
03852415SFPU	Size 3, Large child	each	\$4.79	\$4.55	\$4.55	\$4.31
03851410SFPU	Size 4, Adult	each	\$4.79	\$4.55	\$4.55	\$4.31
03851400SFPU	Size 5, Large adult	each	\$4.79	\$4.55	\$4.55	\$4.31

ALLIED HEALTHCARE LSP VINYL MASK

- Anatomically designed clear vinyl mask
- Pliable single-piece construction
- Cup automatically positions mask on face, creating positive seal



L595060020	Child	each	\$17.49	\$16.62	\$16.62	\$15.74
L595060050	Adult	each	\$20.79	\$19.75	\$19.75	\$18.71

AMBU® TRANSPARENT FACE MASKS

- Blue silicone
- Autoclavable up to 138°
- Latex-free



000315000	Size 0, Infant	each	\$43.49	\$41.32	\$41.32	\$39.14
000316000	Size 2, Child	each	\$43.49	\$41.32	\$41.32	\$39.14
000312000	Size 3-4, Small adult	each	\$43.49	\$41.32	\$41.32	\$39.14
000317000	Size 5, Adult	each	\$43.49	\$41.32	\$41.32	\$39.14

AMBU® SILICONE FACE MASKS

- Conforms to anatomical features of the face to ensure a good fit
- Clear plastic body to monitor lip coloration, bleeding or vomitus
- Soft inflatable cuff ensures comfortable, leak-resistant seal
- Totally autoclavable • Latex-free



000251002	Size 1, Infant	each	\$25.99	\$24.69	\$24.69	\$23.39
000251004	Size 3, Child	each	\$34.99	\$33.24	\$33.24	\$31.49
000251005	Size 4, Adult	each	\$37.79	\$35.90	\$35.90	\$34.01
000251006	Size 5, Large adult	each	\$37.79	\$35.90	\$35.90	\$34.01

TELEFLEX ANATOMICAL FACE MASKS

- 22mm fitting
- Antistatic/conductive rubber with neoprene coating
- With universal harness retention ring



3101102	Size 2, Child	each	\$62.29	\$59.18	\$59.18	\$56.06
3101103	Size 3, Small adult	each	\$62.29	\$59.18	\$59.18	\$56.06
3101104	Size 4, Medium adult	each	\$62.29	\$59.18	\$59.18	\$56.06
3101105	Size 5, Large adult	each	\$62.29	\$59.18	\$59.18	\$56.06

**HUDSON RCI®
DUAL SPRAY METERED DOSE
INHALER ADAPTER**



- Bidirectional nozzle creates a chamber effect
- 22mm OD x 22mm ID connections
- Extended universal port
- Dual-spray nozzle directs plumes in both directions simultaneously
- Tethered cap

1759	Inhaler adapter	each	\$5.79	\$5.50	\$5.50	\$5.21
------	-----------------	------	--------	--------	--------	--------

**MERCURY MEDICAL®
HANK'S POP-OFF VALVE**



- Conductive neoprene (latex-free)
- Lightweight
- Autoclavable (Caution: separate components before sterilizing)
- 15mm male on machine end
- Patient end is nominal 22mm male with taper and shoulder for secure fit in masks and 15mm female for tracheal connectors
- Spring-loaded, adjustable pop-off valve

1900452P	Pop-off valve	each	\$220.79	\$209.75	\$209.75	\$198.71
----------	---------------	------	----------	----------	----------	----------

**GENERAL
ANESTHETIC SERVICE
SEVOFLURANE
KEY FILLER**



GABASV	Key filler	each	\$140.29	\$133.28	\$133.28	\$126.26
--------	------------	------	----------	----------	----------	----------

AIR GAS SODA LIME



12003000	Bag - approx 3lbs	\$16.49	\$15.67	\$15.67	\$14.84
12004000	PRE-PAK canister - approx 3lbs	\$20.79	\$19.75	\$19.75	\$18.71

**ALLOGRAFTS.
XENOGRAFTS.
ALLOPLASTS.**

**THE NEW HOME OF
BONE REGENERATION
STARTS ON PAGE 70
OF THIS CATALOG!**

**OR VISIT US TODAY AT
ACESOUTHERN.COM**



ACE SOUTHERN

LOOKING FOR MOISTURE TRAPS?

Moisture Traps are located in the PATIENT DIAGNOSTIC section.

EDAN PAGE 225	CRITICARE PAGE 227	NONIN PAGE 232
--------------------------	-------------------------------	---------------------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

O₂ only

SALTER LABS® CONVENTIONAL STYLE NASAL CANNULAS

Salter cannulas provide real patient comfort and convenience. Unique features set them apart, while providing patients with many facets they quickly find unique and satisfying.

You can expect Salter cannulas to provide:

- Curved, tapered nasal prongs
- Comfortably shaped facepieces
- Lightweight feel; not made with natural rubber latex
- Twist-resisting slide bolo
- Quality and wearability
- Feature non-flared tips with adjustable tubing
- For low flow volumes up to 6 LPM

1053750	7' Adult with 3-channel safety supply tubing	50/case	\$33.49	\$31.82	\$31.82	\$30.14
1056750	7' Adult with smooth bore supply tubing	50/case	\$33.49	\$31.82	\$31.82	\$30.14

O₂ only

SALTER LABS® 16SOFT NASAL CANNULAS

Designed to deliver comfort with soft material in both the tubing and facepiece for a soft-to-the-touch feel throughout.

The softness is everywhere. Curved, comfortable prongs fit gently into patients' nostrils without pressing too hard upon the upper lip. Soft, pliable, stretchable tubing curves around the cheeks and ears, lessening any pressure. And the Salter 16SOFT features a totally redesigned slide bolo to resist twisting and turning.

- With 3-channel safety supply tubing for low flow volumes up to 6 LPM

16SOFTPED750	7' Pediatric	50/case	\$49.49	\$47.02	\$47.02	\$44.54
16SOFT750	7' Adult	50/case	\$32.79	\$31.15	\$31.15	\$29.51

O₂ only

SMITHS MEDICAL PORTEX™ FIRST BREATH™ OVER-THE-EAR NASAL CANNULAS


- Latex-free

7' Curved

001289	Flared tip	50/case	\$27.99	\$26.59	\$26.59	\$25.19
001283	Non-flared tip	50/case	\$33.99	\$32.29	\$32.29	\$30.59

7' Straight

001285	Flared tip	50/case	\$27.99	\$26.59	\$26.59	\$25.19
--------	------------	---------	---------	---------	---------	---------



O₂ only

HUDSON RCI® NASAL CANNULAS

- Over-the-ear design provides secure fit
- Soft transparent vinyl
- Star Lumen (kink/crush-resistant) or standard tubing
- Flared or nonflared nasal tips
- Latex-free

7' of Star Lumen tubing

1103	Nonflared nasal tips	50/case	\$47.29	\$44.93	\$44.93	\$42.56
1104	Flared nasal tips	50/case	\$47.29	\$44.93	\$44.93	\$42.56

Cannula Only - No tubing

1108	Flared nasal tips	50/case	\$43.29	\$41.13	\$41.13	\$38.96
1109	Nonflared nasal tips	50/case	\$43.29	\$41.13	\$41.13	\$38.96

Adjustable Elastic Holds Cannula

- Adjustable elastic band for comfortable, snug fit below the ears
- Soft transparent vinyl
- 7' of Standard tubing

TEL1100	Nonflared tips	each	\$4.79	\$4.55	\$4.55	\$4.31
---------	----------------	------	--------	--------	--------	--------




O₂ only

FLEXICARE® ADULT NASAL CANNULAS

- Two O₂ ports deliver oxygen to both nostrils – important for patients with a blocked passage way or congestion
- Soft-tip cannula provides superior patient comfort
- 7' Non-kinking tube ensures consistent accuracy and O₂ delivery
- Latex-free • Adult

03210022U	Adult	50/box	\$35.99	\$34.19	\$34.19	\$32.39
03210101U	Soft Tip Curved	50/box	\$41.99	\$39.89	\$39.89	\$37.79



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

TOP SELLING CANNULA SELECTION GUIDE See pages 48-50 for additional product details (tubing characteristics, length, etc.)

BLUE = Male Luer PINK = Female Luer BLACK = NA

	Portex	Salter Lab		FlexiCare		Teleflex Hudson RCI	
O₂ only	001283 001285 001289	1053750 16SOFTPED750	1056750 16SOFT750	03210022U 03210101U		1103 1108	1104 1109
CO₂ only		4000725	4000F725	03210108U			
Dual O₂ and CO₂	001292 001293	47067010 47077725 47037725 4707101025	4706F7025 4707F7725 4703F7725 4706F7010	03210125U 03210127U 03212125U	03210126U 03210128U	RHET43U RHET44U RHET45U	TEL2844 2845

CO₂ only
SALTER LABS® NASAL CO₂ SAMPLING CANNULAS

- End tidal sampling from both nares
- Curved, tapered nasal prongs
- Anatomically curved to fit upper lip
- Delivers accurate, quantitative reading
- Doesn't interfere with patient observation
- Latex-free



7' Adult

4000725	Male luer	25/case	\$51.79	\$49.20	\$49.20	\$46.61
4000F725	Female luer	25/case	\$51.79	\$49.20	\$49.20	\$46.61

CO₂ only
FLEXICARE® ADULT NASAL CO₂ MONITORING CANNULAS

- Two O₂ ports deliver oxygen to both nostrils – important for patients with a blocked passage way or congestion
- Soft-tip cannula provides superior patient comfort
- 7' Non-kinking tube ensures consistent accuracy and O₂ delivery
- Latex-free • Adult



03210108U	CO ₂ Monitoring, Male Luer	25/box	\$47.49	\$45.12	\$45.12	\$42.74
-----------	---------------------------------------	--------	---------	---------	---------	---------

Dual O₂ and CO₂
SALTER LABS® ETCO₂ SAMPLING CANNULAS W/O₂ DELIVERY

- Salter's special End Tidal CO₂ Sampling Cannula permits an undiluted end tidal gas sample, even with simultaneous insufflation of oxygen or gaseous analgesia; Designed for optimum performance with maximum comfort
- Divided design delivers O₂ through one nare and samples CO₂ through the other for more accurate end tidal gas monitoring
- Tapered nasal prongs and curved facepiece provide anatomical fit and comfort
- Light green tubing supports clinical efficiency
- With Dual Port Salter Eyes® (safety apertures help to reduce possible occlusions)
- Latex-free



7' Adult with 2" CO₂ Pigtail

47067010	Male luer	10/case	\$45.79	\$43.50	\$43.50	\$41.21
4706F7010	Female luer	10/case	\$45.79	\$43.50	\$43.50	\$41.21

7' Adult with 2" CO₂ Pigtail

4706F7025	Female luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84
-----------	-------------	---------	----------	----------	----------	---------

7' Pediatric with 7' CO₂ Line

47037725	Male luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84
4703F7725	Female luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84

7' Adult with 7' CO₂ Line

47077725	Male luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84
4707F7725	Female luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84

10' Adult with 10' CO₂ Line

4707101025	Male luer	25/case	\$106.49	\$101.17	\$101.17	\$95.84
------------	-----------	---------	----------	----------	----------	---------

Dual O₂ and CO₂
HUDSON RCI® TELEFLEX® SOFTECH® PLUS ETCO₂ SAMPLING CANNULAS

- Adult cannula with 7' oxygen supply tubing and 7' CO₂ sampling line
- Uses a divided design to deliver oxygen from one nare and sample carbon dioxide from the other, permitting an undiluted end tidal gas sample
- Features softer nasal prongs and lariat for greater patient comfort
- Tubing management clips help keep cables organized
- Smaller nare design allows for enhanced patient comfort
- Non-DEHP, not made with natural rubber latex




2845	Male luer	50/case	\$136.99	\$130.14	\$130.14	\$123.29
TEL2844	Female luer	50/case	\$136.99	\$130.14	\$130.14	\$123.29

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE



2



**HUDSON RCI®
SOFTECH BI-FLO®
NASAL CANNULAS**

- The Softech BI-FLO cannula provides the capability of delivering supplemental oxygen while simultaneously sampling expired O₂.
- The BI-FLO Cannula patented design allows delivery of oxygen to both nasal prongs and sampling of expired gases from both nares, minimizing complications caused by deviated septa and occlusion.
- Available in convenient adult configurations for use with any sidestream O₂ monitor to track respiratory rate and trend expired O₂.
- Latex-Free

7' oxygen supply tubing with CO₂ sampling line and Luer connector

RHET43U	3" CO ₂ , female luer	25/box	\$130.79	\$124.25	\$124.25	\$117.71
RHET44U	7' CO ₂ , male luer	25/box	\$141.29	\$134.23	\$134.23	\$127.16
RHET45U	7' CO ₂ , female luer	25/box	\$141.29	\$134.23	\$134.23	\$127.16



**SMITHS MEDICAL PORTEX™
NASAL CANNULAS W/ETCO₂
MONITORING AND O₂ DELIVERY**

- 7' star lumen oxygen tubing
- Anesthesia circuit oxygen adapter
- 10' gas monitoring line
- Latex-free



001292	Male luer	50/case	\$174.99	\$166.24	\$166.24	\$157.49
001293	Female luer	50/case	\$174.99	\$166.24	\$166.24	\$157.49



**FLEXICARE®
DUAL CO₂/O₂
NASAL CANNULAS**



- Two end tidal CO₂ sampling nares provide excellent sampling accuracy and waveforms
- Two O₂ ports deliver to both nostrils – important for blocked passage way/congestion
- Soft-tip cannula provides superior patient comfort
- Non-kinking tube ensures consistent accuracy and O₂ delivery
- Latex-free

03210125U	Male luer, Adult 7'	25/box	\$105.99	\$100.69	\$100.69	\$95.39
03210126U	Female luer, Adult 7'	25/box	\$105.99	\$100.69	\$100.69	\$95.39
03210127U	Male luer, Adult 10'	25/box	\$110.79	\$105.25	\$105.25	\$99.71
03210128U	Female luer, Adult 10'	25/box	\$110.79	\$105.25	\$105.25	\$99.71
03212125U	Male luer, Pedi 10'	25/box	\$110.79	\$105.25	\$105.25	\$99.71

**PORTEX™ PVC GAS SAMPLING LINES
FOR END TIDAL CO₂ MONITORING**

- Each gas sampling line fitting is designed to fit standard luer connections
- 0.05" I.D.
- Latex-free



3450800	10' Male/Female Luer	50/case	\$85.79	\$81.50	\$81.50	\$77.21
---------	----------------------	---------	---------	---------	---------	---------

**SALTER LABS® ETCO₂
SAMPLING LINES**



- Sample lines compatible with most systems and modalities
- Maximum versatility for connection to most monitors
- Easy to install
- Luer-lok® ends provide positive attachment
- Latex-free

4507725	7' Male/Male	25/case	\$41.29	\$39.23	\$39.23	\$37.16
9300850	8' Male/Male	50/case	\$78.79	\$74.85	\$74.85	\$70.91
45101025	10' Male/Male	25/case	\$50.49	\$47.97	\$47.97	\$45.44
44101025	10' Male/Female	25/case	\$50.49	\$47.97	\$47.97	\$45.44
45151525	15' Male/Male	25/case	\$50.49	\$47.97	\$47.97	\$45.44

**HUDSON RCI®
GAS SAMPLING LINES**

- PVC sampling lines for measurement of end-tidal CO₂
- Coextruded polyethylene sampling lines for precision measurement of gases
- All gas sampling lines include a male luer connector on one end for connection to the patient circuit gas sampling port
- 0.05" I.D., 10 ft.



1233H	Male/Male Luers	50/case	\$156.79	\$148.95	\$148.95	\$141.11
-------	-----------------	---------	----------	----------	----------	----------

**HUDSON RCI®
OXYGEN SUPPLY TUBING**

Used in oxygen administration, this soft and flexible vinyl tubing comes with tubing connectors that fit most oxygen source adapters. Individually packaged.



1115	7' Tubing	each	\$1.29	\$1.23	\$1.23	\$1.16
		50/case	\$56.29	\$53.48	\$53.48	\$50.66

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**HUDSON RCI®
NIPPLE AND NUT ADAPTER**

- Hose adapter for standard DISS oxygen outlets



2555	Adapter	each	\$2.99	\$2.84	\$2.84	\$2.69
------	---------	------	--------	--------	--------	--------

**HUDSON RCI®
TUBING ADAPTER**



1423	22mm ID x 5-7mm tapered nipple	50/case	\$54.99	\$52.24	\$52.24	\$49.49
1846	15mm x 6mm	50/box	\$26.49	\$25.17	\$25.17	\$23.84

**SEDATION RESOURCE
METAL CANNULA ADAPTER**



SDR101001	15mm OD	each	\$18.79	\$17.85	\$17.85	\$16.91
SDR101003	22mm OD	each	\$44.79	\$42.55	\$42.55	\$40.31

**HUDSON RCI® MEDIUM
CONCENTRATION MASKS**

- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 7' oxygen supply tubing



1040	Adult	each	\$3.29	\$3.13	\$3.13	\$2.96
		50/case	\$152.29	\$144.68	\$144.68	\$137.06
1035	Pediatric	each	\$2.49	\$2.37	\$2.37	\$2.24
		50/case	\$98.29	\$93.38	\$93.38	\$88.46

**HUDSON RCI® HIGH
CONCENTRATION MASKS**



1009	Adult	each	\$4.79	\$4.55	\$4.55	\$4.31
		50/case	\$98.29	\$93.38	\$93.38	\$88.46

**FLEXICARE® MEDIUM
CONCENTRATION MASKS**


- The soft and pliable mask molds easily to the patient's face; the feather rim provides maximum comfort to the patient
- All masks are supplied with knitted latex-free straps
- 7' oxygen tubing



03210002U	Adult	each	\$1.49	\$1.42	\$1.42	\$1.34
		50/case	\$60.29	\$57.28	\$57.28	\$54.26

**SMITHS MEDICAL
ADAPTERS AND
CONNECTORS**

- Single-use
- Non-sterile
- Plastic



Straight Connector

717670	4.5mm I.D. 6mm O.D.	10/pack	\$38.29	\$36.38	\$36.38	\$34.46
--------	---------------------	---------	---------	---------	---------	---------

Connectors

715630	3mm I.D. 4mm O.D.	10/pack	\$34.99	\$33.24	\$33.24	\$31.49
715634	5mm I.D. 6mm O.D.	10/pack	\$34.99	\$33.24	\$33.24	\$31.49

**CARDINAL HEALTH
MEDI-VAC® CONNECTORS**

- Lightweight autoclavable plastic connectors come in a variety of configurations and are transparent to allow easy observation of flow. Packaged in a dispenser carton.



5-in-1 Straight Connector

Fits five sizes of tubing. Connect same size tubes or use as a reducer for any combination of 3/16" to 7/16" tubing

350A	5-in-1 connector	50/box	\$89.49	\$85.02	\$85.02	\$80.54
------	------------------	--------	---------	---------	---------	---------

1/4" Straight Connector

Used to join two sections of 1/4" tubing



353A	Straight connector	50/box	\$95.29	\$90.53	\$90.53	\$85.76
------	--------------------	--------	---------	---------	---------	---------

**HUDSON RCI® NON-REBREATHING
MASKS W/SAFETY VENT**

- Clear, soft vinyl for patient comfort and visual assessment
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape
- Safety vent allows for entrainment of room air
- Complete with 7' oxygen supply tubing



1059	Adult	each	\$3.49	\$3.32	\$3.32	\$3.14
		50/case	\$129.79	\$123.30	\$123.30	\$116.81
HU1058	Pediatric	each	\$3.79	\$3.60	\$3.60	\$3.41
		50/case	\$139.29	\$132.33	\$132.33	\$125.36

**FLEXICARE®
NON REBREATHING MASKS**

- Star lumen tubing - ensures that the flow of oxygen is maintained even if the tube is accidentally kinked
- Totally latex-free - eliminates the risk of allergic reaction to latex
- 7' oxygen tubing



03210072NCU	Adult w/nose clip	each	\$9.79	\$9.30	\$9.30	\$8.81
		50/case	\$73.79	\$70.10	\$70.10	\$66.41

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



**THE AXESS MASK:
NEW LOOK OF
N₂O/O₂ SEDATION**

**ACCUTRON INC. AXESS™
LOW PROFILE
SINGLE-USE NASAL MASKS**

Patient comfort and clinical efficiency - a perfect combination!

- Low profile - expands clinical field of view and access to patient's mouth
- Translucent colors and clear unscented mask - shows patient properly breathing through their nose (mask fogs upon exhalation)
- Contoured mask design - provides a good seal minimizing ambient nitrous oxide and gas flow into patient's eyes
- Fits snugly without tape or uncomfortable cannula protruding into patient's nostril
- Two scents (bubblegum and mint) and unscented for scent-sensitive patients - engages patients in mask selection
- Single-use - minimizes cross-contamination; saves staff time
- Not made with natural rubber latex

Masks and CO₂ Adapter - 24/box

530379	Small, unscented	\$148.29	\$140.88	\$140.88	\$133.46
5303712	Small, birthday bubblegum	\$148.29	\$140.88	\$140.88	\$133.46
5303716	Small, fresh mint	\$148.29	\$140.88	\$140.88	\$133.46
530359	Medium, unscented	\$148.29	\$140.88	\$140.88	\$133.46
5303512	Medium, birthday bubblegum	\$148.29	\$140.88	\$140.88	\$133.46
5303516	Medium, fresh mint	\$148.29	\$140.88	\$140.88	\$133.46
530349	Large, unscented	\$148.29	\$140.88	\$140.88	\$133.46
5303412	Large, birthday bubblegum	\$148.29	\$140.88	\$140.88	\$133.46
5303416	Large, fresh mint	\$148.29	\$140.88	\$140.88	\$133.46

Low Profile Nasal Hood CO₂ Adapter* - 24/box

53055	\$63.79	\$60.60	\$60.60	\$57.41
-------	---------	---------	---------	---------

*Allows connection of a CO₂ sample line to monitor EtCO₂



**ACCUTRON INC. AXESS™
SCAVENGING CIRCUITS**

For RFS™ (Remote Flow System) and Accu-Vac™ Vacuum Controller

Includes: 6-foot long low-profile scavenging circuit (to be used with RFS™ and Accu-Vac™ Vacuum Controllers - not provided), complimentary low profile nasal mask assortment, adapters for connecting most non-Accutron vacuum controllers, instructions for use and reprocessing.

52006	Vacuum controller	each	\$169.29	\$160.83	\$160.83	\$152.36
-------	-------------------	------	----------	----------	----------	----------

For In-line Vacuum Controller

Includes: 6-foot long low profile scavenging circuit with in-line vacuum controller, complimentary low profile nasal mask assortment, adapters for connecting most non-Accutron vacuum controllers, instructions for use and reprocessing.

52005	Vacuum controller	each	\$184.79	\$175.55	\$175.55	\$166.31
-------	-------------------	------	----------	----------	----------	----------

**Axess™ Scavenging
Circuit Extension**

32" autoclavable circuit tubing with fittings that allow extension of Axess™ Scavenging Circuit



53004	Circuit extension	each	\$29.99	\$28.49	\$28.49	\$26.99
-------	-------------------	------	---------	---------	---------	---------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**PORTER SILHOUETTE
LOW PROFILE NITROUS OXIDE MASK**

A single patient use nasal mask and breathing circuit designed with lightweight tubing and a comfortable and conforming mask. It's unique design offers advantages over traditional nasal hoods

- Unobstructed access to the oral cavity
- Adhesive strip secures mask in place
- A more predictable nitrous experience
- Easily retrofits to most flowmeters
- Disposable mask/circuit improves infection control
- Enhanced scavenging efficiency
- Four sizes for a perfect fit on all patients reduced N₂O exposure



Porter SILHOUETTE Masks

SILPED012	Pediatric	12/pack	\$143.99	\$136.79	\$136.79	\$129.59
SILSM12	Adult small	12/pack	\$143.99	\$136.79	\$136.79	\$129.59
SILMED12	Medium	12/pack	\$143.99	\$136.79	\$136.79	\$129.59
SILLG12	Large	12/pack	\$143.99	\$136.79	\$136.79	\$129.59
SILVAR	Variety pack (3x of each)	12/pack	\$143.99	\$136.79	\$136.79	\$129.59

Porter SILHOUETTE Accessories

SILSIZER4	SILHOUETTE sizer masks color coded, 1 for each size, autoclavable	4/pack	\$62.99	\$59.84	\$59.84	\$56.69
B5603005	Cannula adapter	5/pack	\$12.29	\$11.68	\$11.68	\$11.06

SILHOUETTE to Flowmeter Connection Kit

Includes: cannula adapter, breathing bag cap, 6.5' fresh gas/vacuum tubing, 4 color coded sizer masks, and clip

SILCONNKIT	Flowmeter connection kit	kit	\$52.99	\$50.34	\$50.34	\$47.69
------------	--------------------------	-----	---------	---------	---------	---------



SILSIZER



SILCONKIT

**MATRX DYNAMITE SINGLE-USE
SCENTED NITROUS SCAVENGING NASAL HOODS**

- Scented hoods come with a colorful sticker
- Single-use reduces risk of cross infection
- Put your patient at ease
- Latex-free

SCENT ▼	SIZE ►	Small	Medium	Large	Assorted
Plain		91316500	91316507	91316514	-
Bubble Gum Scent		91316496	91316503	91316510	-
Strawberry Scent		91316497	91316504	91316511	-
Vanilla Scent		91316499	-	91316513	-
Assorted Scents		91316501	91316508	91316515	91316516

Various Sizes	Scented hoods	12/bag	\$49.99	\$47.49	\$47.49	\$44.99
---------------	---------------	--------	---------	---------	---------	---------



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

ACCUTRON INC. SCAVENGING HUB



The transparent hub for the Scavenging Circuit from Accutron permits visual monitoring of patient respiration during dental procedures. In addition, the hub's low-profile design offers an expanded field of operation. The circuit's precision vacuum-control assembly allows for vacuum flow at the ADA's recommended rate of 45 liters per minute.

22099	Scavenging hub	each	\$52.49	\$49.87	\$49.87	\$47.24
-------	----------------	------	---------	---------	---------	---------

ACCUTRON INC. PIP+™ SINGLE PATIENT USE N₂O NASAL HOODS*

Accutron's Personal Inhaler Plus+ single-use nasal hoods provide your patients with reliable and comfortable nitrous oxide delivery while minimizing the potential for cross contamination.



- 8 soothing scents (including unscented for patients who prefer it)
- Saves valuable staff time
- Non-latex material minimizes potential for allergic reactions
- Thermoplastic elastomer reduces concern for nasal hood collapse

Small

3301712	Birthday bubblegum	24/box	\$109.29	\$103.83	\$103.83	\$98.36
3301717	Groovy grape	24/box	\$109.29	\$103.83	\$103.83	\$98.36
3301711	Sassy strawberry	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330179	Unscented	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330171	Variety pack 1†	24/box	\$109.29	\$103.83	\$103.83	\$98.36

Medium

3301612	Birthday bubblegum	24/box	\$109.29	\$103.83	\$103.83	\$98.36
3301611	Sassy strawberry	24/box	\$109.29	\$103.83	\$103.83	\$98.36
3301614	Vanilla	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330169	Unscented	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330161	Variety pack 1†	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330162	Variety pack 2†	24/box	\$109.29	\$103.83	\$103.83	\$98.36

Large

330159	Unscented	24/box	\$109.29	\$103.83	\$103.83	\$98.36
330153	Variety pack 3†	24/box	\$109.29	\$103.83	\$103.83	\$98.36

- † Variety Pack 1 Contents: orange, strawberry, bubblegum
 Variety Pack 2 Contents: piña colada, vanilla, orange
 Variety Pack 3 Contents: mint, peach, vanilla

* Adapters available for all brands of scavenging systems. Identify your brand and call Accutron at 1-800-531-2221 for your free adapter.

ACCUTRON INC. CLEARVIEW™ CLASSIC SINGLE-USE NASAL MASKS*



- Hood-in-hood design draws excess gas into vacuum reducing ambient nitrous oxide in operatory
- Clear outer hood shows patient "breathing through nose" respiration (condensation appears in hood)
- Soft inner hood creates desirable nasal/facial seal
- Low profile expands clinical field of view
- Single-use minimizes cross-contamination; saves staff time
- Latex-free

Pediatric

3303712	Birthday bubblegum	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303717	Groovy grape	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303710	Outlaw orange	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303711	Sassy strawberry	12/box	\$76.49	\$72.67	\$72.67	\$68.84
330379	Unscented grey	12/box	\$76.49	\$72.67	\$72.67	\$68.84

Adult

3303512	Birthday bubblegum	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303514	French vanilla	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303516	Fresh mint	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303510	Outlaw orange	12/box	\$76.49	\$72.67	\$72.67	\$68.84
330359	Unscented grey	12/box	\$76.49	\$72.67	\$72.67	\$68.84

Large Adult

3303414	French vanilla	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303416	Fresh mint	12/box	\$76.49	\$72.67	\$72.67	\$68.84
3303410	Outlaw orange	12/box	\$76.49	\$72.67	\$72.67	\$68.84
330349	Unscented grey	12/box	\$76.49	\$72.67	\$72.67	\$68.84

ACCUTRON INC. CLEARVIEW™ ADAPTER

Allows connection of the ClearView single-use nasal hoods to a Matrix scavenging circuit that has one smooth tube on each side and a cone with two connections.



35360	Adapter	each	\$61.29	\$58.23	\$58.23	\$55.16
-------	---------	------	---------	---------	---------	---------

ACCUTRON INC. CLEARVIEW™ NASAL MASK AND CAPNOGRAPHY BUNDLE

Single-use nasal masks and capnograph fittings. Adapter allows CO₂ capnograph monitor sample line connection for the purpose of monitoring patient end tidal CO₂ (ETCO₂).



Pediatric

3303712CAP	Unscented grey	12/box	\$96.79	\$91.95	\$91.95	\$87.11
330379CAP	Birthday bubblegum	12/box	\$96.79	\$91.95	\$91.95	\$87.11

Adult

3303514CAP	French vanilla	12/box	\$96.79	\$91.95	\$91.95	\$87.11
330359CAP	Unscented grey	12/box	\$96.79	\$91.95	\$91.95	\$87.11

Large Adult

3303416CAP	Fresh mint	12/box	\$96.79	\$91.95	\$91.95	\$87.11
330349CAP	Unscented grey	12/box	\$96.79	\$91.95	\$91.95	\$87.11

Adapters Only

43055		12/box	\$45.29	\$43.03	\$43.03	\$40.76
-------	--	--------	---------	---------	---------	---------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**FLEXICARE®
SINGLE-USE
NASAL HOODS****



- The feathered edges provide a comfortable fit while also maintaining a secure seal
- Latex-free

Small						
03555210U	Unscented	20/box	\$67.29	\$63.93	\$63.93	\$60.56
03555210GU	Grape	20/box	\$67.29	\$63.93	\$63.93	\$60.56
03555210SU	Strawberry	20/box	\$67.29	\$63.93	\$63.93	\$60.56

Medium						
03555220U	Unscented	20/box	\$73.99	\$70.29	\$70.29	\$66.59
03555220VU	Vanilla	20/box	\$73.99	\$70.29	\$70.29	\$66.59
03555220SU	Strawberry	20/box	\$73.99	\$70.29	\$70.29	\$66.59

Large						
03555230U	Unscented	20/box	\$75.79	\$72.00	\$72.00	\$68.21
03555230VU	Vanilla	20/box	\$75.79	\$72.00	\$72.00	\$68.21
03555230MU	Mint	20/box	\$75.79	\$72.00	\$72.00	\$68.21

**FLEXICARE® DUAL
NASAL HOODS****



The unique Dual Nasal Hood with LuerSafe provides optional CO2 monitoring during dental procedures that utilize nitrous oxide or other forms of sedation.

By tracking the patient's expired CO2, clinicians can evaluate the ventilation condition of the patient during sedation and make the judgment to assess and treat the patient.

Designed to provide desired gas to both nostrils and also sample CO2, at the same time. Quality CO2 waveforms are achieved as the patient's expiratory flow will divert fresh gas away from the inhalation channel.

Feature Flexicare's unique cap-less self-sealing LuerSafe port, which will only accept the correct corresponding color-coded LuerSafe 6' monitoring line to be fitted to open the valve and provide uninterrupted flow to sample CO2 levels.

- Plain scent
- Single-use
- Latex-free

03555310U	Small	20/box	\$100.49	\$95.47	\$95.47	\$90.44
03555320U	Medium	20/box	\$107.29	\$101.93	\$101.93	\$96.56
03555330U	Large	20/box	\$108.79	\$103.35	\$103.35	\$97.91

**ACCUTRON INC.
MULTI-USE NASAL MASKS***



- Latex-free

33020	Small	each	\$102.99	\$97.84	\$97.84	\$92.69
33019	Medium	each	\$121.29	\$115.23	\$115.23	\$109.16
33018	Large	each	\$130.29	\$123.78	\$123.78	\$117.26

**WHITE GAPS NASAL
HOOD LINERS**



WGPE	Pediatric	50/box	\$323.29	\$307.13	\$307.13	\$290.96
WGAD	Adult	50/box	\$334.29	\$317.58	\$317.58	\$300.86

**ACE PRE-PUNCHED MATRIX
AUTOCLAVABLE
NASAL HOODS**



- 3 Pre-punched autoclavable nasal hoods per pack
- 25 Luer Lock adapters with caps

0039620	Small, blue	3/pack	\$264.79	\$251.55	\$251.55	\$238.31
0039640	Medium, blue	3/pack	\$264.79	\$251.55	\$251.55	\$238.31
0039680	Large, blue	3/pack	\$264.79	\$251.55	\$251.55	\$238.31

ACE NASAL HOOD PUNCH



- For Preparing Nasal Hoods For ETCO2 (Capnography) Monitoring
- Easily punch through rubber nasal hoods
- Hole diameter set for a snug Luer Lock fit
- Manufactured from German stainless steel

9600456	Nasal hood punch	each	\$104.99	\$104.99	\$94.49	\$94.49
---------	------------------	------	----------	----------	---------	---------

**QOSINA NASAL HOOD
LUER LOCK ADAPTERS
WITH CAPS**



- 25 Luer Lock adapters with caps only (hoods sold separately)


S66079	25/pack	\$29.49	\$28.02	\$28.02	\$26.54
--------	---------	---------	---------	---------	---------

* Adapters available for all brands of scavenging systems. Identify your brand and call Accutron at 1-800-531-2221 for your free adapter.

** Adapters available for all brands of scavenging systems. Identify your brand and call Flexicare at 1-800-985-3314 for your free adapter.

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



MATRIX AUTOCLAVABLE NITROUS SCAVENGER REPLACEMENT PARTS (CURRENT)

91515094	A	Small hood, blue, autoclavable*	each	\$70.29	\$66.78	\$66.78	\$63.26
91515095	B	Medium hood, blue, autoclavable*	each	\$70.29	\$66.78	\$66.78	\$63.26
91515096	C	Large hood, blue, autoclavable*	each	\$84.29	\$80.08	\$80.08	\$75.86
91515202	D	Vacuum tubing, 5 feet long	each	\$105.99	\$100.69	\$100.69	\$95.39
91515203	E	Gas tubing with slide	each	\$113.79	\$108.10	\$108.10	\$102.41
91515201	F	Hood/hose connector	2/pack	\$24.99	\$23.74	\$23.74	\$22.49
91515207	G	Cage valve	each	\$39.49	\$37.52	\$37.52	\$35.54
91515162	H	Non-breathing valve (not autoclavable)	5/pack	\$15.99	\$15.19	\$15.19	\$14.39
91515208	I	Matrix scavenger cone	each	\$50.99	\$48.44	\$48.44	\$45.89

MATRIX AUTOCLAVABLE NITROUS SCAVENGER REPLACEMENT PARTS (OLDER)

91515152	K	Scavenging cone assembly	each	\$86.79	\$82.45	\$82.45	\$78.11
91515160	L	1-½" Connector tubing	each	\$12.29	\$11.68	\$11.68	\$11.06
91515153	M	Y-piece	each	\$16.49	\$15.67	\$15.67	\$14.84
S506007A	N	Hood/hose connector	each	\$17.49	\$16.62	\$16.62	\$15.74
91515154	O	Vacuum tubing, 5 feet long	each	\$90.29	\$85.78	\$85.78	\$81.26
91515155	P	Fresh gas tubing	each	\$90.29	\$85.78	\$85.78	\$81.26
91515161	Q	Slide	each	\$21.99	\$20.89	\$20.89	\$19.79
91515158	R	Vacuum nipple (small orifice)	each	\$28.49	\$27.07	\$27.07	\$25.64
91515156	S	Fresh gas connector	each	\$15.49	\$14.72	\$14.72	\$13.94
91515159	T	Shut-off valve	each	\$65.49	\$62.22	\$62.22	\$58.94
91515205	U	Vacuum adapter	each	\$14.49	\$13.77	\$13.77	\$13.04
91515151	V	22mm external fresh gas adapter	each	\$42.49	\$40.37	\$40.37	\$38.24



MATRIX UNIVERSAL CONVERSION KIT


Use Matrix DynoMite Hoods on Porter/Brown and Accutron Circuits

- In one easy step, you can convert your existing nitrous breathing circuit to use the popular Matrix DynoMite Nasal Hoods
- Latex-free

91515142	Conversion kit	each	\$74.79	\$74.79	\$74.79	\$74.79
----------	----------------	------	----------------	----------------	----------------	----------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



PORTER NITROUS SCAVENGER

Complete Assemblies

- Includes In-Line Control Kit and 3 Liners for Nasal Hoods
- All components are latex-free

51551	Porter nitrous scavenger, adult	each	\$608.49	\$578.07	\$578.07	\$547.64
-------	---------------------------------	------	----------	----------	----------	----------

REPLACEMENT PARTS FOR PORTER NITROUS SCAVENGER

5054A	A	Adult Hood Kit, Beige, w/3 Liners, autoclavable	each	\$270.49	\$256.97	\$256.97	\$243.44
PID50541	B	Adult Hood Liners, autoclavable	3/pack	\$49.99	\$47.49	\$47.49	\$44.99
5053AD12	-	Adult Hood Liners, disposable	12/pack	\$39.49	\$37.52	\$37.52	\$35.54
5054B	A	Pedo Hood Kit, Beige, w/3 Liners, autoclavable	each	\$270.49	\$256.97	\$256.97	\$243.44
50542	B	Pedo Hood Liners, autoclavable	3/pack	\$58.79	\$55.85	\$55.85	\$52.91
5053PD12	-	Pedo Hood Liners, disposable	12/bag	\$39.49	\$37.52	\$37.52	\$35.54
5061	C	Mask to Tubing Connectors, autoclavable	2/pack	\$27.49	\$26.12	\$26.12	\$24.74
5056	D	Coaxial Inhaler Tubing, 30", autoclavable	each	\$197.29	\$187.43	\$187.43	\$177.56
5057	E	Tension Slide Adjuster, autoclavable	each	\$23.49	\$22.32	\$22.32	\$21.14
5058	F	Fresh Gas Y-Connector, autoclavable	each	\$21.99	\$20.89	\$20.89	\$19.79
PID50603	G	3' Fresh Gas Corrugated Tubing, non-autoclavable	each	\$40.49	\$38.47	\$38.47	\$36.44
50606	G	6' Fresh Gas Corrugated Tubing, non-autoclavable	each	\$61.99	\$58.89	\$58.89	\$55.79
5062	H	Vacuum Line Y-Connector, autoclavable	each	\$24.99	\$23.74	\$23.74	\$22.49
5059	I	Vacuum Hose, 8', autoclavable	each	\$42.49	\$40.37	\$40.37	\$38.24
5501RK	J	In-line Vacuum Kit, non-autoclavable	each	\$96.79	\$91.95	\$91.95	\$87.11



Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



MATRIX MDM FLOWMETER HEAD

The MDM Features Include:

- **Min. 30% O₂ Concentration** – The MDM provides a minimum 30% oxygen concentration at all times. If the oxygen supply is interrupted or inadvertently cut off, the MDM will discontinue the flow of nitrous oxide. If the oxygen supply becomes exhausted, both gases will diminish proportionately.
- **In Ratio Fail-Safe** – Should the oxygen supply be exhausted both gases will diminish in the same proportion, maintaining the selected ratio as the gas pressure falls.
- **Automatic Air Intake Valve** – In the event the nitrous oxide and oxygen contents of the breathing bag are depleted, the air intake valve automatically opens to provide room air to the patient.
- **Non-Rebreathing Patient Circuit** – A check valve between the bag and the patient protects against the rebreathing of exhaled gases.
- **Oxygen Flush** – When the oxygen flush button is depressed, the valve opens and supplies 20 LPM of oxygen immediately.
- **Precision Hand-Calibrated Tapered Glass Flow Tubes** – Feature accurate scales from 1-10 liters permanently fused onto the glass surface for durability and accuracy.
- **Auxiliary Oxygen Outlet** – An outlet is provided for the attachment of the optional Machine Quick Connect and "3-in-1" Demand Valve Resuscitator.
- **Warranty:** 1 year from date of delivery. All rubber and plastic parts and accessories are warranted for 90 days from purchase.

MDM Specifications:

Overall Dimensions: 5"W x 9-3/4"H x 12-1/2"D
 Flush Mount..... 6-1/4" x 12-3/4" x 6"
 Weight: 7 lbs.
 Flush Mount..... 14 lbs.
 Mixed Gas Outlet Diameter:..... 22mm O.D. nominal
 Oxygen Inlet:..... Male DISS
 Nitrous Oxide Inlet:..... Male DISS
 Oxygen Mix Ratio Range: 30%-100%
 Flowmeter Range (each): 1-10 LPM



Shown with Assembly, not included.

94500150	MDM Flowmeter Head	each	\$2,723.29	\$2,587.13	\$2,587.13	\$2,450.96
----------	--------------------	------	------------	------------	------------	------------

MOUNTS AND ACCESSORIES FOR FLOWMETER HEADS

Flowmeter Mounts

92505044	Top-mount bracket, for wall arms	each	\$170.99	\$162.44	\$162.44	\$153.89
92505045	Bottom-mount stud	each	\$65.49	\$62.22	\$62.22	\$58.94

Mobile Stand

2042	5-leg base, adjusts from 34-46"	each	\$411.79	\$391.20	\$391.20	\$370.61
------	---------------------------------	------	----------	----------	----------	----------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

MATRIX FLOW ASSEMBLIES

MDM Flowmeter Head

94500150	A	each	\$2,723.29	\$2,587.13	\$2,587.13	\$2,450.96
----------	----------	------	------------	------------	------------	------------

Matrix products

E198	B	Breathing Bag with Bushing, 3 Liter	each	\$49.29	\$46.83	\$46.83	\$44.36
94525192	C	Y-Selector Valve	each	\$492.79	\$468.15	\$468.15	\$443.51
22698DP	D	42" Corrugated Breathing Tube, black with 22mm female bushings (can only be steam sterilized up to 270°)	each	\$34.29	\$32.58	\$32.58	\$30.86
91515194	G	Scavenging Nasal Inhaler Assembly, Pediatric	each	\$429.79	\$408.30	\$408.30	\$386.81
91515193	G	Scavenging Nasal Inhaler Assembly, Medium	each	\$429.79	\$408.30	\$408.30	\$386.81
91515192	G	Scavenging Nasal Inhaler Assembly, Adult	each	\$429.79	\$408.30	\$408.30	\$386.81

ACE products

0032276	E	Universal Female Adapter, 15mm	each	\$23.49	\$22.32	\$22.32	\$21.14
---------	----------	--------------------------------	------	---------	---------	---------	---------


Anesthesia Associates products

00233P	F1	Hank's Elbow with Pop-Off	each	\$153.29	\$145.63	\$145.63	\$137.96
00422	F2	Non-rebreathing Valve-Positive Press Thumb Button	each	\$418.79	\$397.85	\$397.85	\$376.91
00233	F3	Mask Elbow Valve, 22mm male/15mm female fitting	each	\$38.79	\$36.85	\$36.85	\$34.91


H Full Face Masks - see page 46

2


DIRECTIONAL FLOW ASSEMBLY



NON-SCAVENGING ASSEMBLY



STANDARD SCAVENGER ASSEMBLIES



ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



PORTER MXR® 3000 OS FLOWMETER

The MXR-3000 has been the most trusted and recognized flowmeter featuring Porter precision glass indicator tubes and ergonomic single motion controls.

Features include on/off switch, oxygen flush button, and oxygen quick connect. The Porter MXR can be equipped with either a standard breathing circuit with In-Line vacuum control or a breathing circuit with the patented Automatic Vacuum Switch (AVS) control. Mounting options include: mobile cart, portable "E" cylinder stand.

- Includes oral surgery thread bag T with adapter
- Warranty: 1 year from date of delivery. All rubber and plastic parts and accessories are warranted for 90 days from purchase.

3000 OS	MXR 3000 OS	\$2,556.49	\$2,428.67	\$2,428.67	\$2,300.84
PID2020	MXR 3000 OS Wall Mount	\$348.79	\$331.35	\$331.35	\$313.91




PORTER MXR® 3000 OS FLOWMETER ORAL SURGERY SYSTEM WITH MOBILE STAND


Oral Surgery System Includes:

- MXR-3000 OS Flowmeter with oral surgery threaded bag-tee adapter
- 3-Liter Breathing Bag
- Porter Adult Inhaler Scavenger Assembly
- Directional Y-Valve
- Non-Rebreathing Valve
- Corrugated Tube
- Oral Surgery Star-base Stand, telescopes the top of the MXR-3000 from 46" to 55"

3444 OS	MXR 3444 OS System	\$5,207.29	\$4,946.93	\$4,946.93	\$4,686.56
---------	--------------------	------------	------------	------------	------------



PORTER FLOW ACCESSORIES



Porter Scavenging Nasal Inhaler Complete Assembly

51551	Adult	each	\$608.49	\$578.07	\$578.07	\$547.64
-------	-------	------	----------	----------	----------	----------

Anesthesia Associates Breathing Valves

00422	Non-rebreathing valve-positive press thumb button	each	\$418.79	\$397.85	\$397.85	\$376.91
00233	Mask elbow with 22mm male and 15mm female fitting	each	\$38.79	\$36.85	\$36.85	\$34.91

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

PORTER FLOWMETER Y-SELECTORS

- Provides simple selection of nasal or face mask breathing circuit



2210K	Screw-on Y-selector	each	\$581.29	\$552.23	\$552.23	\$523.16
-------	---------------------	------	----------	----------	----------	----------

OXYGEN AND NITROUS HOSE ASSEMBLIES


OXYGEN - Green	Length	Fitting 1	Fitting 2						
1031100	A	36"	DISS female	DISS female	each	\$57.29	\$54.43	\$54.43	\$51.56
5231220	B	6'	DISS female	DISS female	each	\$44.79	\$42.55	\$42.55	\$40.31
5231222	C	6'	DISS female	Puritan male	each	\$83.79	\$79.60	\$79.60	\$75.41
5231234	D	6'	DISS female	Ohio/MatrX male	each	\$73.79	\$70.10	\$70.10	\$66.41
5231227	E	6'	DISS female	Porter O ₂ DISS QC	each	\$132.29	\$125.68	\$125.68	\$119.06

NITROUS OXIDE - Blue

1031090	F	36"	DISS female	DISS female	each	\$57.29	\$54.43	\$54.43	\$51.56
5231221	G	6'	DISS female	DISS female	each	\$62.99	\$59.84	\$59.84	\$56.69
5231223	H	6'	DISS female	Puritan male	each	\$78.49	\$74.57	\$74.57	\$70.64
5231235	I	6'	DISS female	Ohio/MatrX male	each	\$87.99	\$83.59	\$83.59	\$79.19

COMBO OXYGEN and NITROUS OXIDE

8065	J	5'	DISS female	Porter dual stem	each	\$257.79	\$244.90	\$244.90	\$232.01
------	----------	----	-------------	------------------	------	----------	----------	----------	----------



WESTERN MEDICA STANDARD, PIPED SUPPLY REGULATORS

Fits H cylinder, without hose

34M1540PG	O ₂ regulator	each	\$414.49	\$393.77	\$393.77	\$373.04
-----------	--------------------------	------	----------	----------	----------	----------

PORTER TRIPLE OUTLET

Use with single (5231227) or dual (8065) Porter hose assemblies

62563	Triple outlet O ₂ /N ₂ O/vac	each	\$534.79	\$508.05	\$508.05	\$481.31
-------	--	------	----------	----------	----------	----------



		00129	00132	0032261		
ADAPTERS						
00129	Non-rebreathing valve (N.R.B.) with positive pressure thumb button and adjustable release valve	each	\$48.29	\$45.88	\$45.88	\$43.46
00132	Hank's elbow	each	\$45.79	\$43.50	\$43.50	\$41.21
0032261	Universal female adapter, fits 15mm/20mm tubes	each	\$24.99	\$23.74	\$23.74	\$22.49

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



**ALLIED HEALTHCARE
B&F OXYGEN PRESSURE REGULATORS
WITH CONSTANT FLOW CONTROLLERS**

B&F Oxygen Regulators offer an array of options, such as an easy to read flow rate window and a dial-in click-style flow controller... making operation convenient for the variety of users. B&F oxygen regulators are available with or without a hose barb outlet... with or without a DISS outlet... and with a one-year manufacturer's warranty.

- Uni-body construction designed to strengthen the regulator and prevent breakage
- Nickel-plated brass body for safe use in oxygen service
- Two body designs: yoke inlet for 'D' or 'E' size cylinders, nut and nipple inlet for large cylinders
- Hose barb or DISS outlet to provide secure connection of oxygen supply line to patient
- Compact 1-1/2" diameter low-profile contents gauge with gauge protector
- Serrated dial-in flow control knob for ease of flow rate selection
- Large easy to read flow rate numbers
- Filter: 50 micron sintered bronze
- Safety relief pressure: 75 to 140 psi
- Supply pressure: 500 to 2200 psi
- Outlet pressure: High flow outlets: 40 to 60 psi at 500 to 2200 psi inlet pressure
- Flow capacity: High flow outlets: 100 lpm min. single or combined flow w/500 psi inlet pressure
- Constant flow outlet selections:
0-25 lpm: 0, 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, 25
- 1-year warranty

Single Gauge Regulator

21025	With 2 check valves and constant flow outlet	\$119.49	\$113.52	\$113.52	\$107.54
21030	With constant flow outlet	\$96.79	\$91.95	\$91.95	\$87.11

**ALLIED HEALTHCARE
DUAL-GAUGE OXYGEN
PRESSURE REGULATORS
WITH CONSTANT
FLOW CONTROLLERS**



Chemtron Dual Gauge Regulator

- Pin index yoke inlet connection
- Flow rate: 1-8 LPM
- Barbed outlet • For 'D' and 'E' size oxygen cylinders

32262058	each	\$213.79	\$203.10	\$203.10	\$192.41
----------	------	----------	----------	----------	----------

**MEDLINE
BEST-VALUED OXYGEN
REGULATOR, 0-8 LPM**



- This low cost regulator is fitted with a brass sleeve around key parts
- Laser etched to eliminate flaking and rub-off
- Test method for evaluating the ignition
- Sensitivity and fault tolerance of oxygen
- For all medical and emergency applications
- Connector type - CGA870

HCS8708M	each	\$42.99	\$40.84	\$40.84	\$38.69
----------	------	---------	---------	---------	---------



**ALLIED HEALTHCARE
CHEMTRON OXYGEN FLOWMETERS**

The Chemtron Flowmeters offered here are dual scale 1-15 LPM flowmeters. Dual scale flowmeters allow for .5 LPM increments for precise flow requirements between 1 and 5 LPM, and 1 LPM increments between 5 and 15 LPM settings for general use. All units provide a minimum 65 LPM in full open position marked flood.

CHEMTRON OXYGEN Flowmeters

1500103	Schrader O ₂ male flowmeter	\$44.79	\$42.55	\$42.55	\$40.31
---------	--	---------	---------	---------	---------

**MADA MEDICAL
OXYGEN FLOWMETERS**

**50 psi, 1/8" NPT
Female Back Inlet**

1701	1-15 LPM	each	\$46.29	\$43.98	\$43.98	\$41.66
17018	1/2 -8 LPM	each	\$43.79	\$41.60	\$41.60	\$39.41

**BAY CORPORATION
FLOWMETER TO OUTLET ADAPTERS**

- All have 1/8" NPT

SH252S	Schrader O ₂ male flowmeter	\$53.79	\$51.10	\$51.10	\$48.41
12426	DISS female flowmeter	\$19.79	\$18.80	\$18.80	\$17.81
PQ1242	Puritan male flowmeter	\$15.49	\$14.72	\$14.72	\$13.94
CM1242	Chemtron/NCG Flowmeter	\$29.49	\$28.02	\$28.02	\$26.54
OH1242	Ohio/Matrix flowmeter	\$30.49	\$28.97	\$28.97	\$27.44

ACE FLOWMETER TO OUTLET ADAPTERS

- 1/8" NPT

5231232	Porter flowmeter	\$131.29	\$124.73	\$124.73	\$118.16
---------	------------------	----------	----------	----------	----------

**HUDSON RCI®
NIPPLE AND NUT ADAPTER**

- Hose adapter for standard DISS oxygen outlets

2555	Adapter	each	\$2.99	\$2.84	\$2.84	\$2.69
------	---------	------	--------	--------	--------	--------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



ALLIED HEALTHCARE OXYGEN-POWERED ASPIRATOR

- Operates in conjunction with portable oxygen cylinder or central oxygen source
- Inlet pressure 40-90 psi with O₂ flow – 30 LPM
- Inlet fitting O₂ DISS male

Includes:

- 18 Fr suction catheter without control valve
- feature to clear catheter of entrapped material (hose not included)

L146020	Oxygen-powered aspirator	each	\$233.79	\$222.10	\$222.10	\$210.41
---------	--------------------------	------	----------	----------	----------	----------



ALLIED HEALTHCARE RESUSCITATOR VALVES

Demand Valve Only

L06305R	40 LPM setting	each	\$313.99	\$298.29	\$298.29	\$282.59
L063050	120 LPM setting	each	\$457.79	\$434.90	\$434.90	\$412.01



PORTABLE UNIT ACCESSORIES

Allied Healthcare Cylinder Handwheel

- Lightweight plastic to ensure a firm comfortable grip
- Built-in hose fitting wrench
- Tie-down chain included to prevent loss
- For use with 'D' and 'E' size cylinders

L002865050	Cylinder handwheel	each	\$16.49	\$16.49	\$14.84	\$14.84
------------	--------------------	------	---------	---------	---------	---------



Teleflex Cylinder Wrenches

5080	Die-cast zinc	each	\$7.79	\$7.79	\$7.01	\$7.01
------	---------------	------	--------	--------	--------	--------



MADA MEDICAL OXYGEN CARTS

Demand Valve Resuscitator Kit on Cart

1631E Contents:

- Aluminum cylinder, empty (1602E)
- 50 psi high flow regulator (1635-1)
- Demand valve with adult mask and 6' hose (1534MA)
- Cart (370C)

1631E	each	\$731.79	\$695.20	\$695.20	\$658.61
-------	------	----------	----------	----------	----------

Demand Valve Resuscitator Kit with Inhalator on Cart

1630A15E Contents:

- Aluminum cylinder, empty (1602E)
- High-low flow regulator (R1835GBHF-25)
- Demand valve with adult mask and 6' hose (1534MA)
- Adult oxygen mask with 7' tubing, single-use
- Cart (370C)

1634E	each	\$806.29	\$765.98	\$765.98	\$725.66
-------	------	----------	----------	----------	----------

Oxygen Cylinder on Cart

1630A15E Contents:

- Aluminum cylinder, empty (1602E)
- Adjustable flow regulator, 2-15 LPM
- Nasal cannula, single-use (1335-15)
- Cart (370C)

1630A15E	each	\$216.99	\$206.14	\$206.14	\$195.29
----------	------	----------	----------	----------	----------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

RÜSCH® DISPOSABLE RESUSCITATOR BAGS

- Lightly textured PVC bag for sure grip and comfort
- Flexible, transparent PVC cushion face mask ensures a tight seal and patient comfort
- Transparent re-breathing valve assembly
- 360 degree swivel connection
- Single patient use
- Latex-free



Adult Medium Air Cushion Mask

157100300	With bag reservoir	each	\$31.79	\$30.20	\$30.20	\$28.61
157100100	With tube reservoir	each	\$38.79	\$36.85	\$36.85	\$34.91

Pediatric Small Air Cushion Mask

157300100	With 45cm H2O pressure relief valve	each	\$32.79	\$31.15	\$31.15	\$29.51
-----------	-------------------------------------	------	----------------	----------------	----------------	----------------

Infant Sealflex™ Infant Mask

157400100	With pressure relief valve	each	\$46.29	\$43.98	\$43.98	\$41.66
-----------	----------------------------	------	----------------	----------------	----------------	----------------

PORTEX® 1ST RESPONSE™ ADULT MANUAL RESUSCITATOR

- Collapsible volume bag design for ease of storage
- Horizontally-oriented, low dead-space patient valve assembly
- Patient connector with integral swivel improves range of motion
- 30mm exhalation port accepts 8501 PEEP valve
- Long tethered cap for clean, convenient storage of bag after use
- Enhanced SureGrip™ textured surface helps you maintain a secure grasp



8500	Adult Manual Resuscitator	each	\$17.49	\$16.62	\$16.62	\$15.74
------	---------------------------	------	----------------	----------------	----------------	----------------

MERCURY MEDICAL® ADULT CPR-2 BAGS WITH OPTIONAL COLOR-CODED MANOMETER

The CPR-2 is a full line of resuscitators manufactured with new thermoplastic bag material providing superior lung compliance "feel" with exceptional responsiveness. The CPR-2 now offers clinicians a choice of the same unparalleled quality performance and accessories as the regular synthetic rubber CPR bag product line.

- Color Coded Manometer, accuracy is ± 3cm H₂O between 0-15 cm H₂O
- Monitors both Airway and PEEP Pressure
- Fast re-expansion/high breath-per-minute capability
- Latex-free medical grade thermoplastic bag material



1056002	Adult bag	each	\$25.29	\$24.03	\$24.03	\$22.76
1056009	Adult bag and PEEP valve	each	\$31.49	\$29.92	\$29.92	\$28.34
1056102	Child bag	each	\$28.79	\$27.35	\$27.35	\$25.91

DYNAREX DISPOSABLE RESUSCITATOR BAGS

- MPR/intake valve ensures one-way air-flow
- Swivel port rotates 360° for optimal freedom of movement
- Translucent components allow for improved visualization
- Disposable, single-patient use
- Not made with natural rubber latex or DEHP



36000	Adult	each	\$15.29	\$14.53	\$14.53	\$13.76
36003	Pediatric	each	\$15.29	\$14.53	\$14.53	\$13.76

SMITHS MEDICAL INTERTECH FIRST RESPONSE MANUAL RESUSCITATOR

- Disposable
- High-volume reservoir tube
- Needs oxygen tank with ability to produce flow rates of only 5-10 LPM
- Transparent valve and bag
- Includes clear air cushioned adult mask



SMM8506	Manual resuscitator	each	\$27.29	\$25.93	\$25.93	\$24.56
---------	---------------------	------	----------------	----------------	----------------	----------------

LSP BAG MASK RESUSCITATOR

- Disposable
- Transparent valve
- Oxygen enrichment assembly enables the resuscitator to deliver 100% oxygen
- Includes: adult clear mask, supply tubing and reservoir bag
- Latex-Free



L670040	Bag mask resuscitator	each	\$20.79	\$19.75	\$18.71	\$18.71
---------	-----------------------	------	----------------	----------------	----------------	----------------

AMBU® SPUR, SPUR II RESUSCITATORS

- Latex-Free • Disposable
- Transparent valve can be disassembled and reassembled instantly for obstruction removal
- Support strap helps maintain a firm grip and uniform compression rate
- Includes: - Manual Resuscitator with Oxygen Reservoir permanently attached, 7' Oxygen Tubing and Storage Bag



520211000	Adult, size 5	each	\$25.99	\$24.69	\$24.69	\$23.39
530213000	Pediatric	each	\$18.99	\$18.04	\$18.04	\$17.09

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**HUDSON LIFESAVER®
DISPOSABLE RESUSCITATION BAG**



- Foldable bag collapses for compact storage.
- Designed to meet resuscitation bag standard (ISO 10651-4)
- Unique right-angled oxygen tubing connector swivels 360°
- Includes Star Lumen® Tubing and universal oxygen connector
- All configurations include swivel flow diverter for PEEP valve connection
- Device not made with natural rubber latex
- Manometer option allows for convenient pressure monitoring
- Available in bag reservoir and 40" (1 m) tube reservoir option
- Individually packaged in convenient, tamper-resistant, drawstring bag.

5387	Disposable	each	\$35.99	\$34.19	\$34.19	\$32.39
------	------------	------	---------	---------	---------	---------

**HUDSON LIFESAVER®
RESUSCITATION BAG**



- Autoclavable (except reservoir bag)
- Easy to assemble
- Collapsible for compact convenient storage
- Nonbreathing valve offers low resistance to both inspiratory and expiratory gas flows
- Adapter positioned between mask and bag may be used for pressure or CO₂ monitoring
- Includes mask, reservoir bag and reservoir valve assembly

5345	Adult	each	\$374.79	\$356.05	\$356.05	\$337.31
------	-------	------	----------	----------	----------	----------

**AMBU®
MARK IV RESUSCITATOR
MASK WITH OXYGEN RESERVOIR**



- Double wall bag provides built-in pressure limitation for patient safety
- Unique single-shutter valve system for reliable functionality
- Integrated handle for user comfort and uniform compression
- Autoclavable at 134°: Entire Mark IV resuscitator, including the closed reservoir
- Pressure-limiting valve with override cap on Mark IV Baby for added infant protection
- Reservoir allows delivery of oxygen concentration up to 100% by minimizing escape of gas into atmosphere

304006000	Mark IV, adult	each	\$413.29	\$392.63	\$392.63	\$371.96
-----------	----------------	------	----------	----------	----------	----------

**ALLIED HEALTHCARE
LSP BAG MASK
RESUSCITATOR**



- Autoclavable
- Silicone resuscitator
- Tubing not included
- Includes clear Tru-Fit mask
- May be used with oxygen enrichment assembly to deliver up to 100% oxygen concentration

L238220	Adult	each	\$202.79	\$192.65	\$192.65	\$182.51
L238210	Child	each	\$189.79	\$180.30	\$180.30	\$170.81

**ALLIED HEALTHCARE
OXYGEN ENRICHMENT
ASSEMBLY**



- Simple, effective means for delivery of high oxygen concentration
- The oxygen adapter is inserted into the 22mm inlet port on the valve
- With oxygen input of 15 LPM at 15 cycles per minute the bag mask will deliver up to 100% oxygen
- Bag mask not included
- To be used with L238220 and L238220 Bag Mask Resuscitator

L003482	Oxygen enrichment	each	\$36.99	\$35.14	\$35.14	\$33.29
---------	-------------------	------	---------	---------	---------	---------

**COVIDIEN ARGYLE™ SUCTION
CATHETERS WITH CHIMNEY VALVE**



- Chimney style valve has a large opening on the control vent that minimizes residual suction
- Catheters are available in two tip styles - the smooth atraumatic DeLee tip (5 Fr - 8 Fr) with staggered eyes and the angled Whistle tip (10 Fr - 18 Fr)
- Latex-free

30500	5 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
30888	8 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31000	10 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31200	12 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31400	14 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31600	16 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31800	18 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

**LAERDAL V•VAC™
MANUAL SUCTION SYSTEM**



One-hand operation - no batteries required.
The V•Vac unit works the instant you squeeze it!

- Adjustment for high and low-pressure suctioning
- Efficient: Maximum suction exceeds 380mmHg
- Peak air flow over 70 liters per minute
- Versatile catheter connection
- Disposable one-piece cartridge provides a new, clean pump, valve, and suction tube
- Reusable handle is immersible for cleaning

Starter Kit includes

- Handle, 2 cartridges, short suction catheter w/adaptor tip, double male connector and directions for use

985000	Starter kit	kit	\$268.99	\$255.54	\$255.54	\$242.09
985001	Replacement cartridge	each	\$54.79	\$52.05	\$52.05	\$49.31

**PRECISION MEDICAL
EASYGO VAC ASPIRATOR
PORTABLE SUCTION UNIT**



- High performance unit delivers a vacuum range of 2-21 inHg
- Battery allows 1.5 hours of run time
- Uses any standard 800 cc canister
- A dual pump protection system features two safeguards – a float shut-off in the canister to protect overflow and an In-line Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump
- 2-year warranty on the unit and a 6-month warranty on the battery

Suction Unit includes

- Suction unit, 6' patient suction tubing, 800 cc canister, carrying bag, AC hospital grade power cord, filter assembly (filter, tubing and elbow), battery charger/power supply

PM65HG		kit	\$618.99	\$588.04	\$588.04	\$557.09
--------	--	-----	----------	----------	----------	----------

**GOMCO OPTIVAC®
AC/DC PORTABLE
SUCTION UNIT**



- 3 hours run time at maximum vacuum on a fully-charged battery
- Long-life internal rechargeable battery eliminates the need for frequent replacement
- Sealed water-resistant unit is easy to clean and minimizes potential contamination or damage to the compressor and battery
- Angled control panel for easy viewing
- Impact-resistant and durable
- 3-year warranty on the unit, 1-year on the battery

Suction Unit includes

- 1500ml Disposable Canister, Disposable Suction Tubing and 3 Disposable Hydrophobic Bacteria Filters

G180		kit	\$797.99	\$758.09	\$758.09	\$718.19
------	--	-----	----------	----------	----------	----------


1500ml Disposable Collection Canister, Stem Inlet

20080016		each	\$8.29	\$7.88	\$7.88	\$7.46
		16/case	\$116.49	\$110.67	\$110.67	\$104.84

Disposable Hydrophobic/Bacteria Filter, 1/4" Hose Barb x 1/4" Hose Barb

985001		3/box	\$54.79	\$52.05	\$52.05	\$49.31
--------	--	-------	---------	---------	---------	---------

PROFEX HEAD CRADLE



- Comfortably stabilizes head in any position
- Used for lateral, supine or prone position
- Vinyl/nylon covering
- Reusable
- Depth=8", Width=9", Height=4-1/2"

0187B	Cradle	each	\$105.49	\$100.22	\$100.22	\$94.94
-------	--------	------	----------	----------	----------	---------

**PROFEX DONUT™
FOR SURGERY**



- For face, ear or brain surgery
- ProResilient™
- Foam covered with black conductive sheeting
- Reusable – 9", Black

0196BB	Donut, 9", black	each	\$59.29	\$56.33	\$56.33	\$53.36
--------	------------------	------	---------	---------	---------	---------

**ACE VELCRO
ARM AND BODY STRAPS**



2" Velcro Arm Restraint, Black

0032012	2" x 12"	each	\$12.29	\$11.68	\$11.68	\$11.06
0032020	2" x 20"	each	\$13.49	\$12.82	\$12.82	\$12.14
0032030	2" x 30"	each	\$15.49	\$14.72	\$14.72	\$13.94



2" Velcro Strap with Buckle, Black

0032112	2" x 12"	each	\$13.49	\$12.82	\$12.82	\$12.14
0032116	2" x 16"	each	\$14.49	\$13.77	\$13.77	\$13.04
0032120	2" x 20"	each	\$15.49	\$14.72	\$14.72	\$13.94
0032128	2" x 28"	each	\$16.49	\$15.67	\$15.67	\$14.84



2" Body Velcro Strap, Black

0030062	2" x 6'	each	\$28.49	\$27.07	\$27.07	\$25.64
0030082	2" x 8'	each	\$30.49	\$28.97	\$28.97	\$27.44



1" Velcro Arm Restraint, Black

0031012	1" x 12"	each	\$11.29	\$10.73	\$10.73	\$10.16
0031024	1" x 24"	each	\$13.49	\$12.82	\$12.82	\$12.14
0031028	1" x 28"	each	\$14.49	\$13.77	\$13.77	\$13.04



1" Velcro Arm Restraint, Black

0031112	1" x 12"	each	\$7.79	\$7.40	\$7.40	\$7.01
0031124	1" x 24"	each	\$12.29	\$11.68	\$11.68	\$11.06
0031128	1" x 28"	each	\$13.49	\$12.82	\$12.82	\$12.14

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

RUSCH QUICKTRACH®

Allows quick and safe access for ventilation in the presence of acute respiratory distress with upper airway obstruction. The kit consists of a pre-assembled emergency cricothyrotomy unit with a 10ml syringe, padded neck strap and a flexible connecting tube.



- Conical plastic cannula with fixation flange and 15mm connector
- Removable stainless steel cricothyrotomy needle
- Removable safety stopper
- Preassembled and ready to use
- Single-use
- Latex-free, sterile
- Includes: 1 QuickTrach - Sterile with Stopper, 1 Syringe, 1 Connecting Tube with 15mm Adapter, 1 Cushion Neckband

120900040	Adult, 4.0mm	each	\$232.29	\$220.68	\$220.68	\$209.06
120900020	Child, 2.0mm	each	\$232.29	\$220.68	\$220.68	\$209.06



PORTEX® CRICOTHYROTOMY KIT - PCK

Emergency Airway Access

In difficult environments, where emergency medical care has to be administered efficiently, with the minimum of delay, you need equipment that you can trust, that is quick and safe to use. The kit incorporates a Veress® needle design that confirms entry into the trachea and indicates any subsequent contact with the posterior tracheal wall. The 6mm bore Cricothyrotomy tube enables spontaneous breathing and the Portex® SOFT SEAL® cuff secures the airway. The kit has been designed to meet emergency requirements. The device is pre-assembled and its packaging compact and robust - ideal for inclusion in "trauma bags".

PCK - The Components

- Veress® needle design indicates entry to the trachea and any subsequent contact with the posterior tracheal wall
- 6mm bore Cricothyrotomy tube with SOFT SEAL® cuff enables spontaneous breathing
- 15mm connector
- Designed to meet emergency requirements
- Pre-assembled kit contains all components required for the procedure
- Includes scalpel, tube holder, nylon sutures, Thermovent® T Heat & Moisture Exchanger

PCK - The Package

- Robust and compact packaging - tested to withstand extreme weather conditions
- Components presented in "ready" state and in sequence of use
- Ideal for inclusion in trauma bags
- Supplied in sterile outer bag

100465060	Cricothyrotomy Kit	each	\$246.99	\$234.64	\$234.64	\$222.29
-----------	--------------------	------	----------	----------	----------	----------

WAISMED LTD NIO™ ADULT AUTOMATIC INTRAOSSEOUS ACCESS DEVICE

The NIO™ Adult is an automatic intraosseous access device packaged for safe, quick, and easy vascular access in 10 seconds or less. No batteries, assembly, or extra parts required.



- Sterile, single-use device—eliminates the risk of cross-contamination
- Safe, no exposed needle
- Patient-friendly and user-friendly
- No assembly required
- No external power source or battery required
- Pocket-sized and lightweight—3.5oz
- Eliminates the risk of CLABSI
- NIO fixation
- IV extension set with 90 degree connector, and 3-way stopcock
- Secure IV strap
- Alcohol prep pads
- Patient bracelet
- 5-year shelf life

660115	Adult, 15ga	each	\$171.49	\$162.92	\$162.92	\$154.34
--------	-------------	------	----------	----------	----------	----------

WAISMED LTD NIO™ PEDIATRIC AUTOMATIC INTRAOSSEOUS ACCESS DEVICE

The NIO™ Pediatric is an automatic intraosseous access device packaged for safe, quick, and easy vascular access in 10 seconds or less. The NIO™ Pediatric is intended for pediatric patients ages 3-12 years. Its needle depth may be changed according to the age of the patient by adjusting the needle stabilizer located at the base of the device. In addition, the needle stabilizer has protruding arrows that assist in correct anatomical placement of the device, using the tibia tuberosity as a landmark. The NIO™ Pediatric is compact and requires no assembly making it ideal for use in both pre-hospital and hospital settings.



- Intended for use in the pediatric population ages 3-12 years
- Sterile, single-use device
- Safe, no exposed needle
- Adjustable needle depth: 14mm for ages 3-9
18mm for ages 9-12
- Needle stabilizer
- No assembly required
- No external power source or battery required
- Pocket-sized and lightweight
- NIO fixation
- IV extension set with 90 degree connector, and 3-way stopcock
- Secure IV strap
- Alcohol prep pads
- Patient bracelet
- 5-year shelf life

660118	Pediatric	each	\$171.49	\$162.92	\$162.92	\$154.34
--------	-----------	------	----------	----------	----------	----------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2

HUDSON RCI® CPR BOARD

Aids in the administration of CPR by helping to maintain the position of the patient and to absorb the compression of the procedure. High compression strength material is lightweight and portable.



Measures 23-1/4" x 17-1/4". Weighs 4.5 lbs. Latex-free.

1178	CPR Board	each	\$102.29	\$97.18	\$97.18	\$92.06
------	-----------	------	----------	---------	---------	---------

AMBU® RES-CUE MASK

- 100% Latex-free
- Can be used alone or with a manual resuscitator
- Prevents direct contact with patient
- Transparent dome allows visual check of patients airway status
- Complete version includes: mask, valve/filter, gloves, wipes, head strap and oxygen inlet
- Oxygen Inlet facilitates increased oxygen supply
- Basic version includes: mask, valve/filter, gloves and wipes
- Single-use



000252102	with O ₂ Inlet	50/case	\$17.49	\$16.62	\$16.62	\$15.74
000252103	without O ₂ Inlet	50/case	\$14.49	\$13.77	\$13.77	\$13.04

FLEXICARE® VENTICAIRE® POCKET RESUSCITATION MASK

- Protect patients against cross-contamination with the flexicare single patient use resuscitator mask.
- Secure, durable carrying case with easy open release catch
- Easy-to-follow instructions
- Non-return valve prevents reflux contamination while allowing gas escape through side vents
- Integrated bacterial/viral filter protects both patient and resuscitator
- Pre-filled cushion provides maximum seal
- Latex-free strap assists in maintaining mask in secure position
- Built-in oxygen line port with cap allows connection of oxygen if required
- Single-use



010999U	Resuscitation Mask	each	\$12.49	\$11.87	\$11.87	\$11.24
---------	--------------------	------	---------	---------	---------	---------

CPR MICROSHIELD® BARRIER

- Flexible, foldable
 - Clear, soft barrier
 - One-way anti-reflux valve
 - In sealed, orange pouch with instructions
- Regular size, Opens to 25 sq inches



PLS70150H	CPR Barrier	each	\$16.49	\$15.67	\$15.67	\$14.84
-----------	-------------	------	---------	---------	---------	---------



PHYSIO CONTROL LIFEPAK CR® PLUS AED

The LIFEPAK CR PLUS AED is your complete solution for large scale AED program implementations, facilities looking for a device that is simple to use or with long-term product support.



The CR Plus Features:

- A fully automatic model with a two-step operation and no need to push a shock button
- One of the industry's longest warranties at 8 years

The Benefit of Fully Automatic:

The fully automatic LIFEPAK CR Plus AED can be used with just two simple steps: turn on the device and apply the electrode pads to the victim. The CR Plus then takes over, giving ClearVoice™ prompts and delivering the shock without anyone pushing a button. A recent study shows that safety is not compromised when untrained rescuers use a fully automatic AED compared to a semi-automatic AED.¹

¹ Hosmans T, et al. Safety of fully automatic external defibrillation by untrained lay rescuers in the presence of a bystander. Resuscitation. 2008;77-2:216-219.

LIFEPAK CR Plus AED

- Includes: 2 sets of electrodes, 1 CHARGE-PAK™ battery charger, a soft shell carry case and Ambu® Res-Cue Mask® kit

80403000149			\$1,930.99	\$1,834.44	\$1,834.44	\$1,737.89
-------------	--	--	------------	------------	------------	------------

Replacement Kits for CHARGE-PAK Battery Charger

- Includes: 1 or 2 sets of electrodes (see corresponding item numbers below), battery charger, replacement instructions, and discharger for safe disposal of used CHARGE-PAK charger

11403000001	w/2 sets of electrodes	each	\$148.29	\$140.88	\$140.88	\$133.46
11403000002	w/1 set of electrodes	each	\$126.49	\$120.17	\$120.17	\$113.84

Infant/Child Reduced Energy Defibrillation Electrodes

For use on children less than 8 years of age or less than 55 lbs.

- Includes: 1 zippered storage pouch, 1 set of electrodes and instruction card

11101000017	Pediatric starter kit	each	\$204.49	\$194.27	\$194.27	\$184.04
-------------	-----------------------	------	----------	----------	----------	----------

Replacement Infant/Child Reduced Energy Defibrillation Electrodes

11101000016	1 set of electrodes	each	\$135.99	\$129.19	\$129.19	\$122.39
-------------	---------------------	------	----------	----------	----------	----------

Wall Mount Bracket

11210000021		each	\$117.49	\$111.62	\$111.62	\$105.74
-------------	--	------	----------	----------	----------	----------

ANESTHESIA SUPPLIES/EMERGENCY RESPONSE

2



ZOLL® AED PLUS

When a cardiac arrest occurs, the fact is that only half of the victims will need a shock. The other half requires CPR.

The AED Plus features Real CPR Help®, a tool that is able to actually see what you are doing and provide feedback to help you do it well. Audio and visual prompts help you rescue with confidence and clarity unmatched by any other AED.

- Not pushing hard enough? It will tell you when to push harder.
- Pushing hard enough? It will say, "Good compression."
- Not pushing fast enough? A metronome will lead you to the right rate. It will even show you the depth of each compression, in real time.
- Not yet started? The AED Plus will tell you again to get started.
- Compressions stopped? It will tell you to continue.

Ordinary AEDs require that the rescuer use two confusing electrodes. ZOLL's Full-Rescue AED features CPR-D•padz®, a one-piece, pre-connected pair of electrodes with Real CPR Help. By using CPR hand placement as the landmark, it is easy to place the CPR-D•padz on the victim quickly.

Consumer Batteries for Convenience and Savings - The AED Plus is the only AED to use inexpensive, consumer lithium batteries that are easy to replace and last five years in stand-by mode.

- AED Plus is backed by a 5-year warranty

AED Plus

- Includes: 1 Set of CPR-D•padz, 1 Set of Lithium Batteries, and a Soft Carry Case

102011010		\$1,999.99	\$1,899.99	\$1,899.99	\$1,799.99
-----------	--	------------	------------	------------	------------

Replacement CPR-D•padz

8900080001	CPR-D•padz	each	\$235.99	\$224.19	\$224.19	\$212.39
------------	------------	------	----------	----------	----------	----------

Pedi•Padz® II Pediatric Pads

8900081001		each	\$126.29	\$119.98	\$119.98	\$113.66
------------	--	------	----------	----------	----------	----------

Replacement Lithium Battery

800080701		10/box	\$84.99	\$80.74	\$80.74	\$76.49
-----------	--	--------	---------	---------	---------	---------

Wall Mount

800080901		each	\$75.29	\$71.53	\$71.53	\$67.76
-----------	--	------	---------	---------	---------	---------



CARDIAC SCIENCE POWERHEART® G5 AED

Whether you are an experienced rescuer or a first-time responder, you will appreciate how the Powerheart® G5 AED combines ease of use, reliability and advanced technology to help victims survive sudden cardiac arrest.

The Powerheart® G5 AED with Intellisense™ CPR provides corrective feedback to ensure delivery of high-quality compressions. One button dual-language functionality and comprehensive self-tests give rescuers the confidence to act.

- **Intellisense™ CPR Feedback** Real-time voice and text corrective prompts for CPR quality.
- **Rescue Ready®** Conducts regular, automatic self-tests on pads, battery and internal circuitry.
- **RescueCoach™** Guides users through each critical step of the rescue with user paced voice/text prompts.
- **Dual Language** Change from primary to secondary language in the push of a button.
- **90 minutes of rescue data**, accessible by USB cable or USB memory stick
- **Non-polarized Intellisense™ defibrillation pads** offer a 2-year shelf life
- **Intellisense® Lithium battery** offers a 4-year shelf life
- Powerheart® G5 AED is backed by a 8-year warranty

Powerheart® G5 AED

- Includes: Pads, Batteries, Soft Carry Case, AED Manager CD and USB cable

G5A80CS		\$1,995.49	\$1,895.72	\$1,895.72	\$1,795.94
---------	--	------------	------------	------------	------------

Adult Defibrillator Pads with CPR Device

XELAED001B	each	\$73.99	\$70.29	\$70.29	\$66.59
------------	------	---------	---------	---------	---------

Pediatric Defibrillator Pads

For patients 8 years or younger weighing less than 55 lbs (25kg)

XELAED003C	each	\$112.79	\$107.15	\$107.15	\$101.51
------------	------	----------	----------	----------	----------

Intellisense® Replacement Battery

XBTAED001A	each	\$423.29	\$402.13	\$402.13	\$380.96
------------	------	----------	----------	----------	----------

BONE GRAFTING

3

alloOss® ALLOGRAFT PARTICULATE, BLOCK AND PUTTY

alloOss Allograft is offered in particulate, block and putty. alloOss is a safe, predictable and effective solution for providing a scaffolding for cell ingrowth and promoting the growth of new bone.

Natural Allograft Bone

alloOss mineralized Allograft Particulates and Blocks preserve natural collagen and minerals with the bony structure of allograft. The existence of these natural minerals and collagen facilitate remodeling of the host bone suitable for implantation. alloOss Allograft Putty combines demineralized bone with cancellous bone and an organic based carrier. alloOss products are designed to promote new bone formation.

Safe and Effective

The complete line of alloOss products have been processed and sterilized to the highest standards set forth by the American Association of Tissue Banks (AATB®) and the U.S. Food and Drug Administration (FDA) regulations. alloOss allograft is obtained and processed in conjunction with AATB accredited tissue banks.

Functional Quality Comparison: Cancellous and Cortical

Particulate Qualities	Cancellous	Cortical
Osteoconductive Function	★ ★ ★	★ ★
Space Maintenance	★ ★	★ ★ ★
Density	★ ★	★ ★ ★
Volume Enhancement	★ ★	★ ★ ★



Cleansing Process for alloOss

Phase 1: Stringent screening procedures for donor selection.



Phase 2: Initial physical inspection and debridement to eliminate substandard tissue and aseptic placement in canisters, reducing risk of tissue contamination. Tissue is processed via minimal manipulation in cGTP conditions in ISO 5 (Class 100) clean rooms using procedures and screening criteria that meets the requirements of the AATB.



Phase 3: Cleansing reaches deep within the tissues with a process that sustains consistency. The process, which uses mechanical energy and rotation, is customized based on specific tissue type and weight. Microbiological testing is performed before and after processing to ensure safety of the final product.



Phase 4: Low dose low temperature irradiation provides the tissue with a Sterility Assurance Level (SAL) of 10⁻⁶ without compromising the mechanical integrity of the allograft tissue.



Phase 5: Full water immersion extraction, as opposed to batch swab, is taken for increased accuracy of microbial presence reading.

AATB® is a registered service mark of the American Association of Tissue Banks.

alloOss® MINERALIZED PARTICULATE

alloOss Allograft Particulates offer the structural strength required to maintain space and volume during the remodeling phase. Used alone or as a composite graft, alloOss Particulates are an effective alternative to intraoral graft harvesting.



- Provides scaffold for cell ingrowth
- Predictable Guided Bone Regeneration
- Fills voids and extensive bone defects
- Remodels completely
- Stabilizes implants after 4-6 months
- Regeneration of periodontal defects

Mineralized Cancellous Particulate

01107001	250-1,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01107002	250-1,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01107003	250-1,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23
01107101	1,000-2,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01107102	1,000-2,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01107103	1,000-2,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23

Mineralized Cortical Particulate

01160001	250-1,000 mic.	0.5cc	\$77.29	\$77.29	\$61.83	\$61.83
01160002	250-1,000 mic.	1.0cc	\$101.79	\$101.79	\$81.43	\$81.43
01160003	250-1,000 mic.	2.0cc	\$162.29	\$162.29	\$129.83	\$129.83
01160101	1,000-2,000 mic.	0.5cc	\$77.29	\$77.29	\$61.83	\$61.83
01160102	1,000-2,000 mic.	1.0cc	\$101.79	\$101.79	\$81.43	\$81.43
01160103	1,000-2,000 mic.	2.0cc	\$162.29	\$162.29	\$129.83	\$129.83

Mineralized 50-50 Cortico/Cancellous Particulate

01104000	250-1,000 mic.	0.25cc	\$71.29	\$71.29	\$57.03	\$57.03
01104001	250-1,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01104002	250-1,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01104003	250-1,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23
01104101	1,000-2,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01104102	1,000-2,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01104103	1,000-2,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23

**alloOss®
NATURAL BLEND
PARTICULATE**

alloOss Natural Blend is a perfect blend of both cortical and cancellous bone particles, allowing for a more natural mirroring of the composition of human bone; making it a perfect solution for promoting the growth of new bone.



Mineralized Cortico/Cancellous Particulate

01108200	500-1,000 mic.	0.25cc	\$71.29	\$71.29	\$57.03	\$57.03
01108201	500-1,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01108202	500-1,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01108203	500-1,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23

**alloOss®
DEMINEALIZED
PARTICULATE**

alloOss Demineralized Cortical Bone is appropriate for use in any procedure where bone healing is desired. Routine testing for osteoinductivity is performed to ensure the highest quality and efficacy for your regeneration procedures.



Demineralized Cortical Bone

01170001	125-710 mic.	0.25cc	\$30.29	\$30.29	\$24.23	\$24.23
01170002	125-710 mic.	0.5cc	\$51.79	\$51.79	\$41.43	\$41.43
01170003	125-710 mic.	1.0cc	\$77.79	\$77.79	\$62.23	\$62.23
01170004	125-710 mic.	3.0cc	\$200.79	\$200.79	\$160.63	\$160.63



Mineralized Cortico/Cancellous Particulate In A Syringe

01108204	500-1,000 mic.	0.5cc	\$115.49	\$115.49	\$92.39	\$92.39
01108205	500-1,000 mic.	0.75cc	\$134.99	\$134.99	\$107.99	\$107.99

Mineralized/Demineralized Cortico/Cancellous Particulate

01109200	500-1,000 mic.	0.25cc	\$71.29	\$71.29	\$57.03	\$57.03
01109201	500-1,000 mic.	0.5cc	\$96.79	\$96.79	\$77.43	\$77.43
01109202	500-1,000 mic.	1.0cc	\$126.29	\$126.29	\$101.03	\$101.03
01109203	500-1,000 mic.	2.0cc	\$202.79	\$202.79	\$162.23	\$162.23

**alloOss®
BLOCK GRAFT**

A cortico/cancellous block for restoring bone volume in severely resorbed ridges without the harvesting of autogenous block graft from the patient.

- Approximate healing in 5-6 months
- Cancellous component for rapid healing and revascularization of host bone to the new graft
- Maintains strength and volume with predictable remodeling
- Cortical component for strength and rigidity of the graft with volume enhancement



Cortico/Cancellous Block Graft

01241006	18 x 6 x 8mm	each	\$597.79	\$597.79	\$478.23	\$478.23
01241010	18 x 10 x 8mm	each	\$693.49	\$693.49	\$554.79	\$554.79
01241015	18 x 15 x 8mm	each	\$776.99	\$776.99	\$621.59	\$621.59

**alloOss® DEMINEALIZED
ALLOGRAFT PUTTY**

alloOss Putty is a combination of DBM (demineralized bone matrix), cancellous FDBA (freeze dried bone allograft), and a carrier medium of CMC (carboxymethyl-cellulose). This mixture results in a highly viscous dental bone grafting putty with easy handling, inductive potential and excellent osteoconductive scaffolding.

Characteristics

- High osteoconductive ability
- High viscosity
- Excellent biocompatibility
- Forms freely and is completely moldable
- Osteoinductive ability

Demineralized Allograft Putty

01105002	150-850 mic.	0.5cc	\$129.29	\$129.29	\$103.43	\$103.43
01105003	150-850 mic.	1.0cc	\$232.79	\$232.79	\$186.23	\$186.23



BONE GRAFTING

3

alloOss® PLUS ALLOGRAFT BONE MATRIX

Introducing alloOss® Plus

This is a truly exciting addition to the alloOss allograft family. alloOss Plus is a bone viable matrix that is aligned and priced for dental use. It contains a proprietary mix of bone fibers, bone chips and particles that retains its native bone-forming cells¹. alloOss Plus is simple to use, just thaw and go without any decanting. The fibrous structure is what makes this unique product exceptionally pliable and easy to shape and form around defects.

alloOss Plus Contains:¹

- Viable cancellous bone with endogenous bone forming cells, including osteoblast, osteoprogenitor, and mesenchymal stem cells (MSC).
- Cortical-derived demineralized bone (DBM) fibers and particles.
- Smaller bone chips aligned for dental use.
- Cortical-derived DBM fibers provide better scaffold than DBM particles and prevent less cancellous bone chip migration from implant site.

Handling Characteristics:

- AlloOss Plus Bone Viable Matrix (BVM) fibrous material allows for easy forming around defects.



Safety Profile

- Aseptic tissue processing.
- Antimicrobial treatment of tissue.
- Sterility cultures performed on every lot.
- Selective immunodepletion process.
- Non-immunogenic evidenced by the inhibition proliferation of lymphocytes.¹

Donor Testing

- Hepatitis B Surface Antigen
- Hepatitis B Core Antibody
- Hepatitis C Virus Antibody
- Hepatitis C Nucleic Acid Test
- HIV-1 and 2 Antibody
- HIV-1 Nucleic Acid Test
- HBV NAT
- Syphilis (Treponema Pallidum)

Donor Screening

- Medical and social history evaluation.
- Physical examination.
- Medical record evaluation, including autopsy report (if performed).
- Licensed physician review and release of every donor record.

Quality Assurance

- cGTP compliant according to FDA regulations.
- All applicable state licensure / registrations.
- AATB® Accreditation

Packaging and Logistics:

- AlloOss Plus is available in a 1cc or 2cc size. It is provided in a shelf box containing a vial sealed in a sterile pouch, IFU and patient labels.
- AlloOss Plus is stored at -75°C and shipped on dry ice. The product should remain sealed in the shipper until time of use.
- The sealed shippers are validated to maintain appropriate storage temperature for either 6 days (if 3 or 8 units or less are ordered) or for 4 days (if up to 20 units are ordered).
- The product is shipped overnight, please plan on ordering the product to arrive the day before the surgery. Please be mindful of holidays when planning your cases.
- AlloOss Plus can be stored in a cryogenic freezer (-75°C) until the expiration date of the product.



Handling:

- Do not overpack the site. AlloOss Plus can be compressed, but once the pressure is removed, the material will bounce back.
- AlloOss Plus can be thawed in 5-15 minutes at room temperature.
- AlloOss Plus should be used within 2 hours post thaw.
- AlloOss Plus can be used immediately after thawing, no additional wash steps are required.
- AlloOss Plus has no liquid to decant.
- Unlike particulate, alloOss Plus is a cortical fiber-based product, so it cannot easily be mixed with other bone particulate. If multiple bone grafting product are required for a case, we recommend layering the bone grafts.
- May be used with or without a membrane, based on clinical judgement.

Important Ordering Details

- Due to the nature of the product, the order and delivery of alloOss Plus needs to be carefully coordinated. alloOss Plus is not available for purchase online. alloOss Plus requires Priority Shipping and Handling, which is in addition to prices shown. Please call for ordering details.

alloOss® Plus Allograft Bone Matrix

01103001	1cc vial	\$399.99	\$399.99	\$319.99	\$319.99
01103002	2cc vial	\$799.99	\$799.99	\$639.99	\$639.99

alloOss® Plus Bone Viable Matrix is manufactured for ACE Surgical Supply Co., Inc. AATB® is a registered service mark of the American Association of Tissue Banks.

Footnote References:

- ¹ Data on file at manufacturer.
- * See Instructions For Use (IFU) for full prescribing information.

Consider the Benefits of our In-Office

CRYOGENIC FREEZER PROGRAM

- Save on alloOss Plus shipping charges
- Maximize treatment planning options
- Safely store alloOss Plus for up to 2 years (or until the expiration date of the product)
- Fits under counter-tops
- Inquire for complete program details

Freezer Storage Specifications:

Temperature: -50° to -122°F (-45° to -86°C)
 Dimensions: 29"h x 20"w x 22"d
 Capacity: 1 Cubic Foot Storage



OSTEOCEL® PLUS CELLULAR ALLOGRAFT BONE MATRIX

Osteocel® Plus mimics the biological profile of autograft. The proprietary processing technology used to produce Osteocel® Plus results in a viable bone matrix that preserves the native Mesenchymal Stem Cells (MSCs) and Osteoprogenitors found in marrow-rich bone. It is the first product to have the desired beneficial osteoconductive, osteoinductive and osteogenic properties of autograft. Osteocel® Plus allows surgeons to offer a complete solution without the added risk and cost of performing a secondary autograft procedure.

Viable Cell Content

Each lot is tested for cell concentration, cell viability, and cell activity (osteogenic potential). The osteogenic potential arises from the mesenchymal and osteoprogenitor cells in Osteocel® Plus. Following processing of marrow-rich bone, release testing demonstrates osteogenic potential according to the following criteria:

- Rich supply of osteopotent cells: Greater than 250,000 cells/cc
- Viability: Greater than 70% cell viability
- Positive osteogenicity: In vitro cell culture assay

Safety Profile

- Aseptic tissue processing.
- Antimicrobial treatment of tissue.
- Sterility cultures performed on every lot.
- Selective immunodepletion process.
- No evidence of ectopic tissue formation or inflammatory cellular response with experimental high doses of MSCs cells.5

Donor Testing

- Hepatitis B Surface Antigen
- Hepatitis B Core Antibody
- Hepatitis C Virus Antibody
- Hepatitis C Nucleic Acid Test
- HIV-1 and 2 Antibody
- HIV-1 Nucleic Acid Test
- Human T-Lymphotropic Virus Antibody I/II
- Syphilis (Treponema Pallidum)

Donor Screening

- Medical and social history evaluation.
- Physical examination.
- Medical record evaluation, including autopsy report (if performed).
- Licensed physician review and release of every donor record.

Quality Assurance

- GTP compliant according to FDA regulations.
- All applicable state licensure / registrations.
- AATB® Accreditation



Important Ordering Details

- Due to the nature of the product, the order and delivery of Osteocel® Plus needs to be carefully coordinated. Osteocel® Plus is not available for purchase online. Osteocel® Plus requires Priority Shipping and Handling, which is in addition to prices shown. Please call for ordering details.

Osteocel® Plus Allograft Bone Matrix

5013001A	1cc vial	\$751.99	\$751.99	\$639.19	\$639.19
5013005A	5cc vial	\$3,385.79	\$3,385.79	\$2,877.92	\$2,877.92

Osteocel® and Osteocel® Plus are federally registered trademarks of NuVasive, Inc.

Footnote References:

- 1 As of March 2017. Data on file with NuVasive®, Inc.
- 2 Mesenchymal Stem Cells Avoid Allogeneic Rejection. – Ryan, Barry, Murphy, Mahon. – Journal of Inflammation, July 2005
- 3 Histologic Evaluation of a Stem Cell Based Sinus Augmentation Procedure: A Case Series. – McAllister, Haghigat, Gonshor. – Journal of Perio., April 2009
- 4 Data on file with NuVasive, Inc.
- 5 Data on file with Osiris Therapeutics

Ridge Augmentation Clinical Case



Consider the Benefits of our In-Office

CRYOGENIC FREEZER PROGRAM

- Save on Osteocel® Plus shipping charges
- Maximize treatment planning options
- Safely store Osteocel® Plus for up to 5 years
- Fits under counter-tops
- Inquire for complete program details

Freezer Storage Specifications:

Temperature: -50° to -122°F (-45° to -86°C)
 Dimensions: 29"h x 20"w x 22"d
 Capacity: 1 Cubic Foot Storage



BONE GRAFTING

3

NuOss® XENOGRAFT PARTICULATE

A Reliable Bone Substitute Used by Dental Professionals Worldwide

NuOss natural cancellous and cortical bone matrix is produced by the removal of all organic components from bovine bone. Due to its natural porous structure, NuOss is physically and chemically comparable to the mineralized matrix of human bone.

The proprietary manufacturing process of NuOss allows for maintenance of the naturally occurring bone mineral structure, including macroscopic and microscopic structures. Due to this trabecular architecture with interconnecting macro and micro pores, the ingrowth of new bone and associated vascularization at the grafting site is optimized.

NuOss may be used with autogenous bone where there is an insufficient quantity of bone to fulfill the needs of the proposed surgical procedure.

Handling: The structure of NuOss aids in particulate hydration and delivery of the product to the defect.

Particle Distribution: NuOss maximizes the distribution of the preferred particle sizing within desired ranges.

Particle Structure and Porosity: NuOss maintains the high porosity and trabecular structure of natural mineralized bone and maximizes the potential for angiogenesis.¹

NuOss Cancellous Particulate

5099001	0.25-1.0mm	0.25g	0.6cc	\$86.29	\$86.29	\$69.03	\$69.03
5099002	0.25-1.0mm	0.5g	1.2cc	\$138.79	\$138.79	\$111.03	\$111.03
5099003	0.25-1.0mm	1.0g	2.4cc	\$232.29	\$232.29	\$185.83	\$185.83
5099004	0.25-1.0mm	2.0g	4.7cc	\$402.99	\$402.99	\$322.39	\$322.39
5099102	1.0-2.0mm	0.5g	1.7cc	\$138.79	\$138.79	\$111.03	\$111.03
5099104	1.0-2.0mm	2.0g	6.8cc	\$402.99	\$402.99	\$322.39	\$322.39

NuOss Cortical Particulate

5091001	0.50-1.0mm	0.25g	0.4cc	\$86.29	\$86.29	\$69.03	\$69.03
5091002	0.50-1.0mm	0.5g	0.8cc	\$138.79	\$138.79	\$111.03	\$111.03
5091003	0.50-1.0mm	1.0g	1.6cc	\$232.29	\$232.29	\$185.83	\$185.83
5091004	0.50-1.0mm	2.0g	3.2cc	\$402.99	\$402.99	\$322.39	\$322.39



NuOss DS Cancellous granules in a preloaded Delivery Syringe

Particulate: NuOss® DS features our best-selling small particle NuOss® cancellous.

Delivery Syringe: The preloaded NuOss Delivery Syringe makes it easier to hydrate and place the graft material, minimizing preparation time and particulate migration at the surgical site.

NuOss Cancellous Particulate in a Delivery Syringe

5099025	0.25-1.0mm	0.25cc	\$67.99	\$67.99	\$54.39	\$54.39
5099050	0.25-1.0mm	0.50cc	\$98.29	\$98.29	\$78.63	\$78.63

NuOss is manufactured for ACE Surgical.

1. Klenke FM, Liu Y, Yuan H, Hunziker EB, Siebenrock KA, Hofstetter W. Impact of pore size on the vascularization and osseointegration of ceramic bone substitutes in vivo. J Biomed Mater Res A. 2008 Jun 1; 85(3): 777-86.

NuOss® COLLAGEN PRODUCTS

Just like NuOss granules, naturally structured NuOss Collagen combines the mineral properties of cancellous particulate with the space maintenance of collagen. This highly manageable combination provides a convenient approach to filling periodontal and maxillofacial defects. With NuOss Collagen products, there is no mixing or loose handling.



NuOss Plug Anorganic Bovine Bone with Collagen

NuOss Plug contains low heat processed anorganic bone that will remodel and support placement of an implant within, on average, 6-9 months. NuOss Plug is easy to place and prevents particulate migration once placed into the surgical site.

NuOss Plug

5094001	Tapered, 9 x 15mm	5/box	\$275.79	\$275.79	\$220.63	\$220.63
5094002	Straight, 10 x 20mm	5/box	\$274.99	\$274.99	\$219.99	\$219.99


NuOss Collagen Block Anorganic Bovine Bone with Collagen

NuOss Collagen Block: Available in 3 sizes, NuOss Collagen can be easily trimmed to shape the defect. Hydration can be done prior to insertion or once it has been placed.



NuOss Collagen Block

5099100	6 x 7 x 8mm	100mg	\$98.29	\$98.29	\$78.63	\$78.63
5099250	8 x 9 x 9mm	250mg	\$183.49	\$183.49	\$146.79	\$146.79
5099500	10 x 11 x 12mm	500mg	\$305.79	\$305.79	\$244.63	\$244.63



NuOss XC Self-Expanding composite

Expansion: Expanding composite material allows for placement in a compressed form.

Upon hydration/contact with blood, the material will self-expand to fill the entire defect.

NuOss XC Socket: Preshaped and ready to place directly into the extraction cavity.

NuOss XC Sinus: Upon successful elevation of the sinus membrane, NuOss XC Sinus is placed without pre-hydration, into the sinus cavity.

NuOss XC Socket

5092001	6 x 15mm expanded state	\$92.49	\$92.49	\$78.62	\$78.62
5092002	9 x 15mm expanded state	\$147.49	\$147.49	\$125.37	\$125.37

NuOss XC Sinus

5093001	17 x 10mm expanded state	\$147.49	\$147.49	\$125.37	\$125.37
5093002	22 x 12mm expanded state	\$238.49	\$238.49	\$202.72	\$202.72

**synOss™
SYNTHETIC BONE
PARTICULATE**

synOss Synthetic Mineral is an osteoconductive calcium phosphate based bone graft material with a carbonate apatite structure similar to the structure of human bone.

Material Composition

The material composition is calcium phosphate in a carbonate apatite structure. The product is supplied in granular form, and is sterile, non-pyrogenic, and for single use only.

- Particulate size: 0.35-1.0mm



When synOss is hydrated it holds to the delivery spoon (#1000216), even when it's inverted.

5095000	0.25 gram	0.5cc	\$74.29	\$74.29	\$63.15	\$63.15
5095100	0.5 gram	1.0cc	\$110.79	\$110.79	\$94.17	\$94.17
5095200	1.0 gram	2.0cc	\$202.29	\$202.29	\$171.95	\$171.95
5095300	2.0 gram	3.5cc	\$353.99	\$353.99	\$300.89	\$300.89

**CaSO₄
CALCIUM SULFATE
HEMIHYDRATE**

CaSO₄ is one of the first bone materials ever used in a grafting procedure and maintains a safety record dating back more than 100 years. CaSO₄ is considered one of the clinician's first choices among synthetic bone grafting materials. Calcium Sulfate is sterile, highly purified and continues its bioactivity through the entire healing process. CaSO₄ is totally resorbed at the completion of the healing period.

CaSO₄ Features

- Osteoconductive
- Promotes wound healing
- Easy to set-up
- Easy to place
- Combines well with bone grafting particulate
- Inhibits non-osteogenic cells
- Completely resorbable
- Easy to handle
- Cost effective

CaSO₄ Kit Contents

- 0.5 gram or 1.0 gram calcium sulfate powder in a sterile mixing dish
- Regular set liquid

5052002	0.5 gram	each	\$43.49	\$41.32	\$41.32	\$39.14
5052001	1.0 gram	each	\$51.29	\$48.73	\$48.73	\$46.16



BE THE OFFICE HERO

Save Time and Money with...



REPLENISH

Our Replenish program allows you to set up scheduled deliveries for the products your practice uses on a regular basis.

Ensure that your practice never runs out of essential items that keep your office running smoothly.

SHOP from hundreds of eligible Replenish products

CHOOSE your delivery frequency; every week, every month, you're in control

SAVE an additional 3% on all Replenish-eligible products

*See page 155 for complete details or visit acesouthern.com/replenish



BONE GRAFTING

3



alloGuide® GF GROWTH FACTOR, MEMBRANE

AlloGuide GF is a lyophilized placental barrier membrane that preserves its endogenous growth factors associated with the repair of oral wounds. AlloGuide GF is an ideal membrane to support Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR) dental procedures as it retains the high levels of regenerative growth factors naturally present in the placental membrane.

The GF stands for Growth Factor

Placental Tissue is a rich source of extracellular matrix proteins and growth factors that have been shown to be involved in the natural wound repair process¹ and is commonly used in a wide array of therapeutic applications, including wound care and orthopedics.^{2,3} The processing and preservation of the placental tissue allograft will influence the preservation of those growth factors.^{4,5}

A Sampling of the Preserved Growth Factors⁶

Growth Factor	Present	Growth Factor	Present
BDNF	✓	HGF	✓
bFGF	✓	OPG	✓
FGF	✓	PDGF	✓
BMP-4	✓	SCF	✓
BMP-5	✓	TGF	✓
BMP-7	✓	TGFβ-3	✓
EGF	✓	VEGF	✓

Characteristics of AlloGuide GF:

- Preserves protein structures and maintains bioactive capabilities of important wound healing elements of placental tissue
- Preserves extracellular matrix
- Maintains natural tissue thickness
- Excellent handling and ease of application.

AlloGuide GF Safety:

- AATB® standards, FDA regulations and stringent internal processes, policies and procedures
- Processed under aseptic conditions
- Sterility testing performed on every lot

AlloGuide GF Processing:

- AlloGuide GF is processed using a proprietary freeze-drying solution during lyophilization
- Preserves protein structures and maintains bioactive capabilities of all the important wound healing elements in the placental tissue
- Preserves the extracellular matrix
- Maintains the natural tissue thickness
- Provides excellent handling characteristics that allow for ease of application.

AlloGuide GF

1250808	8 x 8mm	each	\$59.29	\$59.29	\$50.40	\$50.40
1251010	10 x 10mm	each	\$80.99	\$80.99	\$68.84	\$68.84
1251520	15 x 20mm	each	\$163.79	\$163.79	\$139.22	\$139.22
1252030	20 x 30mm	each	\$255.49	\$255.49	\$217.17	\$217.17

alloGuide® GF is manufactured for ACE Surgical Supply.

References

- ¹ Chopra and Thomas 2013. J Biomim Biomater and Tissue Eng
- ² Brantley and Verla 2015. Advances in Wound Care
- ³ Heckmann et al 2016. American Journal of Orthopedics
- ⁴ McQuilling et al 2017. Wounds
- ⁵ Paolin et al 2016. Cell Tissue Bank
- ⁶ Data on File



MemGuide® RESORBABLE PORCINE MEMBRANE

MemGuide Porcine Dental Membrane is a non-pyrogenic porcine derived dental barrier membrane. MemGuide is sourced from porcine peritoneum which is a naturally thicker membrane than porcine small intestine submucosa¹. The profile of MemGuide results in outstanding handling characteristics. MemGuide retains Growth Factors VEGF and TGF² VEGF and TGFβ are known to be key contributors to tissue repair.³ The Growth Factors in porcine peritoneum share homology with human growth factors which can induce response in human cells.^{4,6} MemGuide is fully resorbable, eliminating the need for a second surgical procedure and is fully resorbable with a 6-month resorption profile.

Features

- Fully resorbable with a resorption time of 6 months
- Guides healing and regeneration of bone and surrounding tissue
- Place and trim in a hydrated or dry state
- May be sutured when required
- Available in 3 convenient sizes
- Packaged sterile in a pouch in pouch configuration
- Store at room temperature. 3-year shelf life

Advantages:

- Fully resorbable - no second surgery required for removal
- Long resorption time for sustained function during bone repair
- Flexible clinical procedure allows for trimming and placing in hydrated or dry state
- Excellent handling characteristics
- Endogenous extracellular matrix proteins and growth factors are preserved^{3,7}
- Porcine growth factor share homology with human growth factors can induce response in human cells.⁴⁻⁶

MemGuide

7351520	15 x 20mm	each	\$134.49	\$134.49	\$114.32	\$114.32
7352030	20 x 30mm	each	\$165.79	\$165.79	\$140.92	\$140.92
7353040	30 x 40mm	each	\$267.29	\$267.29	\$227.20	\$227.20

MemGuide® is manufactured for ACE Surgical Supply.

References

- ¹ Shi and Ronfard Int J Burns Trauma. 2013; 3(4): 173–179.
- ² Hoganson, et al, Biomaterials. 2010 (27): 6934-40.
- ³ Chopra and Thomas J Biomim Biomater and Tissue Eng. 2013 18(1)
- ⁴ Stroobant and Waterfield. The EMBO Journal 1984 3(12).2963-2967
- ⁵ Taylor Journal of Leukocyte Biology 2009. 85
- ⁶ Holmes and Zachery Genome Biology 2005. 6(2) Article 209
- ⁷ Luo et al Journal of Biomedical Materials Research B: Applied Bio Materials 2015

conFORM RESORBABLE COLLAGEN MEMBRANE

conFORM is a resorbable collagen membrane matrix made from Type I bovine collagen. Its pliable characteristics and predictable resorption time of 3-4 months make it an ideal membrane for guiding the healing of bone and surrounding tissue. When hydrated, conFORM will naturally adapt to surgical defects and maintain a high tensile strength, which makes it the ideal membrane for suturing.

With its in vivo stability, permeability and mechanical strength, conFORM is a reliable choice for various clinical applications. When used in conjunction with a bone grafting particulate, conFORM can help mitigate the migration of the particulate from the surgical site. Its resorbability eliminates the need for surgical removal. conFORM is packaged sterile, non-pyrogenic, and for single use.

Features

- White, non-friable membrane matrix engineered from highly purified Type I Collagen
- Fully resorbable with a predictable resorption time of 3-4 months¹
- Guides healing and regeneration of bone and surrounding tissue
- Unique fiber orientation provides high tensile strength for in situ stability of membrane
- Cell occlusive, designed to retard and/or prevent epithelium down growth

Advantages

- Conforms to surgical site without re-shaping
- Easy to handle, does not adhere to itself
- Sufficient resorption time for sustained function in bone repair
- Can easily be trimmed during surgery for a final fit
- May be placed with either side against the tissue or bone

conFORM

5091520	15 x 20mm	each	\$139.99	\$139.99	\$118.99	\$118.99
5092030	20 x 30mm	each	\$172.49	\$172.49	\$146.62	\$146.62
5093040	30 x 40mm	each	\$277.99	\$277.99	\$236.29	\$236.29

References

1. D.Yuen, J. Ulreich, G. Zuclich, H-B. Lin and S-T Li "Prediction of In Vivo Stability of a Resorbable, Reconstituted Type I Collagen Membrane by In Vitro Methods." 2000 Society for Biomaterials - Sixth World Congress.



RCM6[®] RESORBABLE COLLAGEN MEMBRANE

RCM6[®] Resorbable Collagen Membrane allows you to guide healing and regeneration of bone and surrounding tissue within your time frame. RCM6 is intended to aid in the wound healing process within oral surgical procedures, such as: dental implants, bone defects, periodontal defects and ridge augmentations. When used in conjunction with a bone grafting particulate, RCM6 can help mitigate the migration of the particulate from the surgical site. RCM6 is packaged sterile, non-pyrogenic, and for single use. Its resorbability eliminates the need for surgical removal.


Features

- White, non-friable membrane matrix engineered from highly purified Type I Collagen
- Fully resorbable with a resorption time of up to 6 months¹
- Guides healing and regeneration of bone and surrounding tissue
- Unique fiber orientation provides high tensile strength for in situ stability of membrane
- Cell occlusive, designed to retard and/or prevent epithelium down growth

Advantages

- Long resorption time, for sustained function in bone repair
- High mechanical strength for membrane stabilization
- Each piece can be easily trimmed during surgery for a final fit
- May be placed with either side against the tissue or bone

Resorption Comparison *



* Approximate in vivo resorption time.

RCM6

5096000	15 x 20mm	each	\$113.99	\$113.99	\$96.89	\$96.89
5093000	20 x 30mm	each	\$139.99	\$139.99	\$118.99	\$118.99
5098000	30 x 40mm	each	\$209.49	\$209.49	\$178.07	\$178.07

References

1. D.Yuen, J. Ulreich, G. Zuclich, H-B. Lin and S-T Li "Prediction of In Vivo Stability of a Resorbable, Reconstituted Type I Collagen Membrane by In Vitro Methods." 2000 Society for Biomaterials - Sixth World Congress.
2. Reported from the product's Instructions for Use.
3. RCM6 Data on file

BioMend[®] Extend is a registered trademarks of Zimmer Dental.
Bio-Gide[®] is a registered trademark of ED. GEISTLICH SOHNE AG.



INTEGRA[®] HELIMEND[®] ABSORBABLE COLLAGEN MEMBRANES

Integra absorbable collagen membranes are white, compressed, non-friable matrices fabricated from collagen derived from bovine-deep flexor (Achilles) tendon. These collagen products provide wound stabilization and create space during guided tissue regeneration procedures where an extended resorption profile is desired. These membranes are completely absorbable, eliminating the need for the second surgery that is often required to remove a non-resorbable membrane. The HeliMend membrane incorporates into the surrounding tissue and is generally reabsorbed within 4-8 weeks, and the HeliMend Advanced collagen membrane is generally reabsorbed within 18 weeks.

- Can be cut to any size or shape in the wet or dry state, without tearing or fragmenting
- Textured surface for easy handling, non-slippery when wet
- Wet characteristics greatly facilitate positioning & handling during surgical procedures
- Semi-occlusive, it allows essential nutrients to pass through the membrane
- Have mechanical strength to be sutured in place
- Prevents epithelial cell migration and maintains space for periodontal ligament and enhances clinician's ability to promote bone regeneration
- Functions as a barrier during the critical period of wound healing

HeliMend Absorbable Collagen Membrane


- Absorbed in 4 to 8 weeks

62203	15 x 30mm	each	\$120.49	\$120.49	\$102.42	\$102.42
62204	20 x 30mm	each	\$151.79	\$151.79	\$129.02	\$129.02
62205	30 x 40mm	each	\$241.79	\$241.79	\$205.52	\$205.52

HeliMend Advanced Absorbable Collagen Membrane

- Absorbed in 4 to 8 weeks

62206	15 x 30mm	each	\$137.29	\$137.29	\$116.70	\$116.70
62207	20 x 30mm	each	\$169.49	\$169.49	\$144.07	\$144.07
62208	30 x 40mm	each	\$248.79	\$248.79	\$211.47	\$211.47



BONE GRAFTING

3



RESORBABLE COLLAGEN PLUG, FOAM AND TAPE WOUND DRESSINGS

- Provides an excellent matrix to assist in wound healing
- Beneficial in protecting the oral wound site
- Controls minor bleeding due to the intrinsic hemostatic properties of collagen fibers

Features	Plug	Foam	Tape
Fully resorbable	✓	✓	✓
Ideal for applications in oral wounds	✓	✓	✓
Provides an excellent matrix to assist in wound healing	✓	✓	✓
Protects the oral wound bed	✓	✓	✓
Ideal for extraction sockets	✓		
Controls minor bleeding due to the intrinsic hemostatic properties of collagen fibers	✓	✓	✓

RC Plug Collagen Plug Wound Dressing

- Sterile

5094000	1cm x 2cm	10/box	\$123.29	\$123.29	\$104.80	\$104.80
---------	-----------	--------	----------	----------	----------	----------

RC Foam Collagen Foam Wound Dressing


- Sterile

5092000	2cm x 4cm x 3mm	10/box	\$135.79	\$135.79	\$115.42	\$115.42
---------	-----------------	--------	----------	----------	----------	----------


RC Tape Collagen Tape Wound Dressing

- Sterile


5091000	2.5 x 7.5cm x 1mm	10/box	\$185.49	\$176.22	\$176.22	\$166.94
---------	-------------------	--------	----------	----------	----------	----------



62-202



62-200



62-201

INTEGRA® ABSORBABLE COLLAGEN WOUND DRESSINGS FOR DENTAL SURGERY

Integra absorbable collagen wound dressing products are soft, white, pliable, non-friable, absorbent sponges made from collagen obtained from bovine deep flexor (Achilles) tendon. The tendon is one of the purest sources of collagen that can be readily obtained and processed in commercial quantities.

These dressings are intended for use on moist or bleeding wounds where a stabilized blood clot can help protect the surface of the wound from further injury. Collagen wound dressing products are typically reabsorbed by the body within 2 to 4 weeks after placement.

- Dressings retain their structural integrity even when wet
- Very porous, the dressings are highly absorbent, holding many times their own weight of saline solution
- Controls bleeding usually within two to five minutes
- In situations where frequent aspiration is required to remove accumulated blood and exudate, the Absorbable Collagen Wound Dressings can be used to dress bleeding or oozing wounds and thereby improve patient comfort
- Easily removed without damaging the wound bed, they can be used as temporary dressings prior to surgical closure
- The structure of the dressings provides a three-dimensional matrix for additional strengthening of the blood clot
- Designed to be totally absorbable when left in situ after hemostasis
- Require no preparation time and are easy to handle

HeliPlug Collagen Plug Wound Dressing

- Extraction sites
- Biopsy sites

62202	3/8" x 3/4"	10/box	\$141.99	\$141.99	\$120.69	\$120.69
-------	-------------	--------	----------	----------	----------	----------

HeliCote Collagen Cote Wound Dressing

- Palatal donor sites
- Mucosal flaps

62201	3/4" x 1-1/2"	10/box	\$154.79	\$154.79	\$131.57	\$131.57
-------	---------------	--------	----------	----------	----------	----------

HeliTape Collagen Tape Wound Dressing

- Minor oral wounds
- Closure of grafted sites
- Repair of Schneiderian Membranes

62200	1" x 3"	10/box	\$208.29	\$208.29	\$177.05	\$177.05
-------	---------	--------	----------	----------	----------	----------



SOCKET PRESERVATION KIT

All-In-One: NuOss®, RC Plug, conFORM and Violet Braided PGA Suture

The Socket Preservation Kit is the perfect combination of bone grafting, membrane and closure products for your single tooth extractions. This kit is designed to make socket preservation convenient and successful each and every time.

Contents

- 1 - 0.5 gram NuOss Cancellous Granules
- 1 - RCP Resorbable Collagen Plug
- 1 - conFORM 15 x 20mm Resorbable Collagen Membrane
- 1 - 3-0 Violet Braided PGA Suture, DS18 Needle

0034227	Single-use kit	each	\$260.79	\$260.79	\$221.67	\$221.67
---------	----------------	------	----------	----------	----------	----------


Instruments shown holding the products above are not included in this kit.

BONE GRAFTING

3

NON-TEXTURED

TEXTURED on both sides



CYTOFLEX® TEFGUARD® PTFE BARRIER MEMBRANES

Cytoflex® Tefguard® is a non-resorbable membrane for use as a space-maintaining barrier in guided tissue and guided bone regeneration procedures. Tefguard® is composed of a proprietary micro porous polytetrafluoroethylene (PTFE) thin film.

The construction of the surface and porosity of the material is suitable for adhesion of host cells to the material while preventing passage and integration of bacteria within the interstices of the material. The supple membrane conforms easily to tissue contours, and yet offers sufficient stiffness to maintain space over the bony defect.

Both smooth and micro-texture surfaces are designed to prevent tissue and bacteria ingrowth. Tefguard® Textured design allows for increased flap attachment for additional localized stability. Both Tefguard® options are intended for use as a barrier membrane in guided tissue and guided bone regeneration procedures.

Tefguard® Advantages:

- Excellent handling, rigidity and easily adaptable
- Can be easily trimmed to desired shape
- Can be easily fixated using either bone tacks or bone screws
- Easily retrieved as one piece
- Resists fibroblasts and bacteria penetration
- Increased capability of host tissue attachment with Tefguard® Textured
- Facilitates tissue regeneration despite flap recessions, or incomplete primary closure
- 0.2mm thick

Non-textured A

C010101	25 x 30mm	each	\$69.49	\$69.49	\$59.07	\$59.07
C010105		5/pkg	\$279.29	\$279.29	\$237.40	\$237.40

Non-textured B Single extraction

C010301	12 x 24mm	each	\$51.49	\$51.49	\$43.77	\$43.77
C010305		5/pkg	\$222.99	\$222.99	\$189.54	\$189.54

Textured C

C040201	25 x 30mm	each	\$85.79	\$85.79	\$72.92	\$72.92
C040205		5/pkg	\$340.29	\$340.29	\$289.25	\$289.25

Textured D Single extraction

C040101	12 x 24mm	each	\$66.99	\$66.99	\$56.94	\$56.94
C040105		5/pkg	\$275.49	\$275.49	\$234.17	\$234.17

ETHICON VICRYL™ KNITTED MESH

Ethicon Vicryl™ (polyglactin 910) Knitted Mesh for temporary wound or organ support. Made of an absorbable material, it can provide temporary support during your patients' healing process.

- Undyed
- 6" x 6" (15cm x 15cm)
- Sterile

VKMM	3/box	\$3,424.29	\$3,253.08	\$3,253.08	\$3,081.86
------	-------	------------	------------	------------	------------



Small and large membranes are shown actual size.

CYTOFLEX® TEFGUARD® PTFE TI-ENFORCED MEMBRANES

TITANIUM REINFORCED PTFE BARRIER MEMBRANES FOR GTR AND GBR PROCEDURES

Cytoflex® Tefguard® PTFE Textured Membranes with titanium reinforcements are non-resorbable membranes for use as a space-making barrier in guided tissue and guided bone regeneration procedures.

- Improved handling, ideal for modeling and shaping around defects
- Dual textured surface allows membrane to be placed either side up or down
- Range of sizes and configurations for large and small defects
- Can be easily trimmed to desired shape (It is recommended to avoid cutting titanium struts when trimming reinforced membranes)
- Can be easily fixated using either bone tacks or bone screws
- Sterile

Ti-Reinforced Tefguard PTFE Barrier Membrane

C050301	E 12 x 24mm	each	\$128.79	\$128.79	\$109.47	\$109.47
C050901	F 19 x 26mm	each	\$157.99	\$157.99	\$134.29	\$134.29
C051701	G 25 x 30mm	each	\$215.79	\$215.79	\$183.42	\$183.42
C053901	H 25 x 30 x 41mm	each	\$256.29	\$256.29	\$217.85	\$217.85
C054101	I 32 x 40mm	each	\$262.29	\$262.29	\$222.95	\$222.95



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

BONE GRAFTING

3

TRUFIX FIXATION SYSTEM

The All-In-One Bone Fixation and Membrane Stabilization Solution


The TruFIX System from ACE Surgical Supply provides you with everything you need for both bone fixation and membrane stabilization. This system incorporates all the necessary components to pick up and drive the ACE TruSCREW and TruTACK. The TruFIX system eliminates the need for multiple systems and unnecessary components, making it user-friendly for your practice.

- True centered patented locking connection
- True axial alignment with pickup each time
- TruSCREW self-drilling screw developed for easy insertion with maximum fixation
- TruSCREW Patented Removal Sleeve for disengaging the screw from the driver without damaging screw head
- TruTACK hexagon driver head and barbed tip to pierce without drilling for simple tack insertion and easy screw-like removal

The TruSCREW, with its aggressive cutting flutes, is the ideal bone screw for the fixation of small bone within the oral and maxillofacial environment. These cutting flutes (in most instances) eliminate the need for any pre-drilling. The patented design of the screw ensures an effortless insertion into all types of bone.

The TruTACK makes the stabilization of membranes quick and trouble-free. Our unique tack incorporates a hexagon on its head and threads on its shaft, allowing for easy removal. The truTACK is placed like a tack and removed like a screw, a feature that you will not find in any other system.


The TruTENT is a refinement of the truSCREW. Its raised collar and wider head is designed to support a titanium mesh during augmentation procedures.



The TruFIX driver's true-centered locking connection holds the TruSCREW securely.

The removal sleeve disengages the driver from the screw without compromising the screw or the stability of the graft.


45418015	truFIX Complete Fixation System	1 kit	\$1,184.99	\$1,184.99	\$1,007.24	\$1,007.24
System Includes:						
45419001	A truFIX Driver Handle	each	\$332.29	\$332.29	\$299.06	\$299.06
45417001	B truTACK Driver Tip (comes with Tip Cover)	each	\$100.79	\$100.79	\$90.71	\$90.71
45415001	C truSCREW Driver Tip	each	\$78.49	\$78.49	\$70.64	\$70.64
45417901	C truSCREW Driver Removal Sleeve	each	\$81.49	\$81.49	\$73.34	\$73.34
45415201	D truSCREW Driver Tip, Contra-Angle	each	\$89.49	\$89.49	\$80.54	\$80.54
45440203	E 2-Step Countersink Bur, Ø 0.8 and 1.6mm, Contra-Angle	each	\$64.79	\$64.79	\$58.31	\$58.31
45440202	F Ø 1.1mm Pilot Twist Drill, Contra-Angle, 29.8mm Length	each	\$75.29	\$75.29	\$67.76	\$67.76
502700045	G Ø 0.45mm Pilot Bur, Contra-Angle, 27mm Length	each	\$49.49	\$49.49	\$44.54	\$44.54
20157702	H Hand Driver, 0.88mm External Hex	each	\$77.79	\$77.79	\$70.01	\$70.01
45450201	I truSCREW Packaging Removal Tool	each	\$49.49	\$49.49	\$44.54	\$44.54
45418501	truFIX Tray	each	\$388.79	\$388.79	\$349.91	\$349.91
TruFIX Small Driver Handle (sold separately)						
This smaller, thinner, and lighter handle is often preferred by surgeons with smaller hands.						
45419501	J Small Driver Handle, 89mm (3.5") long	each	\$132.29	\$132.29	\$119.06	\$119.06



Handles shown actual size.

TruTACK, TruSCREW AND TruTENT

Bone Tacks, Self Drilling Bone Screws and Tenting Screws



A truTACK 2.5mm head diameter

- Sterile

see chart 10/cassette **\$149.79** **\$149.79** **\$127.32** **\$127.32**

Total Length	Thread Ø	Number
3.0mm	0.7mm Thread Ø	9600313
5.0mm	0.7mm Thread Ø	9600314

B truSCREW 3mm head diameter

- Sterile

see chart 5/cassette **\$145.49** **\$145.49** **\$123.67** **\$123.67**

Total Length	1.2mm Thread Ø	1.5mm Thread Ø	2.0mm Thread Ø
4.5mm	45427202	—	—
6.0mm	45427203	45427502	45428002
7.5mm	45427204	45427503	45428003
9.0mm	45427205	45427504	45428004
10.5mm	45427206	45427505	45428005
12.0mm	—	45427506	45428006
13.5mm	—	45427507	45428007
15.0mm	—	45427508	45428008

C truTENT 5mm head diameter

- Sterile

see chart 1/package **\$52.99** **\$52.99** **\$45.04** **\$45.04**

Total Length	Collar Height	1.5mm Thread Ø
10.0mm	4mm	454391001
12.0mm	6mm	454391002
14.0mm	8mm	454391003





BONE GRAFTING

ACE TITANIUM MICRO MESH

Customizable Titanium Mesh for Guided Bone Regeneration

- See next page for details
- Mesh images are shown at actual size.
- Circled detail images are enlarged 400%.

3



Recessed holes in mesh D and E are compatible with \varnothing 3mm heads of the ACE TruScrew.

Holes will accept screws up to a 2mm thread diameter.

ACE TITANIUM MICRO MESH

Customizable Titanium Mesh for Guided Bone Regeneration

See previous page for images

- Manufactured from commercially pure titanium (6Al-4V ELI)
- Perforated • Non-sterile
- Can be secured with ACE TruFIX Self-Drilling bone screws
- Trim to size using ACE Wire Cutting Scissors



0039429	A	34 x 25mm	0.1mm thick	each	\$64.79	\$64.79	\$55.07	\$55.07
0039430	B	120 x 60mm	0.1mm thick	each	\$250.99	\$250.99	\$213.34	\$213.34
0039433	C	152 x 66mm	0.2mm thick	each	\$270.49	\$270.49	\$229.92	\$229.92
00394501	D	98.5 x 77.5mm	0.6mm thick	each	\$270.49	\$270.49	\$229.92	\$229.92
00394701	E	98.5 x 77.5mm	0.2mm thick	each	\$270.49	\$270.49	\$229.92	\$229.92


ACE TITANIUM TENTING MESH

Pre-Shaped Titanium Mesh for Guided Bone Regeneration

- **F G H** feature a 12 point center locking design, which supports the connection on top of standard 2.7mm External Hexagon Implants and ACE Tenting Screws
- Manufactured from commercially pure titanium (6Al-4V ELI), perforated, non-sterile
- Can be secured with ACE TruFIX self-drilling bone screws

0039434	F	7 x 14mm	0.2mm thick	each	\$61.49	\$61.49	\$61.49	\$52.27
0039440	G	13 x 33mm	0.2mm thick	each	\$76.79	\$76.79	\$76.79	\$65.27
0039442	H	13 x 33mm custom	0.2mm thick	each	\$88.29	\$88.29	\$88.29	\$75.05
0039444	I	30 x 80mm	0.2mm thick	each	\$115.79	\$115.79	\$115.79	\$98.42

Mesh inset images are enlarged 400%



Tenting Mesh all shown 100% actual size

ACE TITANIUM FOIL

- Non-perforated
- Sterile



0039435	16 x 25mm	each	\$33.29	\$33.29	\$28.30	\$28.30
0039436	76 x 100mm	each	\$280.29	\$280.29	\$238.25	\$238.25

ACE MEMBRANE PUNCH

- Revised punch design minimizes the flared edge when punching titanium mesh membranes.
- The first punch flares a mesh membrane.
- Flip the instrument so the second punch flattens out the flared edge.
- Stainless



9600303	3-1/4" (8.5cm)	each	\$80.29	\$80.29	\$72.26	\$72.26
---------	----------------	------	--------------------	----------------	--------------------	----------------

BONE GRAFTING

3



Regular Bone



Autografted Bone Histology of Compressed Cortico/Cancellous Bone

ACE AUTOGRAFTER™

Simple method for obtaining autogenous bone

The patented design of the ACE Autografter allows the surgeon to filter and retain particles of bone which have been suctioned along with standard irrigation during surgery. The graft may be obtained from almost any bony site; examples include dental implant preparation sites, the mandible, maxilla, calvarium, etc. The collected graft material can be used either in a wet amorphous form, or a compressed osseous particulate. The drier particulate,

which is easier to handle, can be obtained by compressing the collected material in a screen with the titanium plunger. This forces out excess moisture through the screen.

NOTE: The ACE Autografter™ is not a routine suction device; therefore, bleeding should be controlled prior to use. The ilium should not be used as clots will clog this device.

ACE Autografter™

0034215		Autografter™ Start Up Kit, (includes: A, C, D and G)	1 kit	\$448.29	\$448.29	\$403.46	\$403.46
0034212	A	ACE Autografter™, for HVAC Tip, Autoclavable	each	\$218.79	\$218.79	\$196.91	\$196.91
0034210	B	ACE Autografter™	each	\$190.49	\$190.49	\$171.44	\$171.44

ACE Autografter™ Accessories

0900024	C	HVAC Suction Tip, 4mm Tip, 45° Angle, 7" (18cm)	each	\$32.29	\$32.29	\$29.06	\$29.06
0034211	D	Titanium Plunger, 4-1/4", 11cm	each	\$119.49	\$119.49	\$107.54	\$107.54
9600449	E	Square Titanium Bowl, 12mm Deep, 50 x 50mm	each	\$105.99	\$105.99	\$95.39	\$95.39
1000216	F	Implant Curette with Round Cups, Titanium, 5-1/2" (14cm) Cups are 5 x 10mm & 6 x 20mm	each	\$115.49	\$115.49	\$103.94	\$103.94
6462002	G	Screen Assembly, Disposable, Sterile, 12cc Capacity	10/box	\$134.79	\$134.79	\$121.31	\$121.31
0030276		Suction Connection Tubing, Sterile, 1' (3cm)	each	\$4.79	\$4.55	\$4.55	\$4.31
0034238		Suction Connection Tubing, Sterile, 3' (.9m)	20/case	\$60.49	\$57.47	\$57.47	\$54.44
0034239		Suction Connection Tubing, Sterile, 6' (1.8m)	20/case	\$82.29	\$78.18	\$78.18	\$74.06
Y2606	H	Y Connector, 7.5mm I.D., 11.5mm O.D.	each	\$5.79	\$5.50	\$5.50	\$5.21

OSSEOUS TECHNOLOGIES
DISPOSABLE BONE COLLECTOR

- Harvest valuable autogenous bone
- Correct small defects around implant sites and periodontal comprised teeth
- Shut off valve cuts off suction eliminating the chance of drying out bone
- Sterile collectors are sold in packs of 10
- Single patient use only
- Includes: 10 Collectors, each with 1 Filter, per Bag

BC100	bag	\$183.99	\$183.99	\$165.59	\$165.59
-------	-----	----------	----------	----------	----------



BONE GRAFTING

3

ACE REUSABLE BONE SCRAPERS

The tool for simple and fast recovery of allogeneic bone grafts.

- Made of surgical-grade stainless steel, free of biocompatibility hazards
- Easy to clean and sterilize
- One re-usable cutting blade included, replacements are available
- Stainless Steel



Bone Scraper, Straight

39895700	Straight	each	\$276.99	\$276.99	\$249.29	\$249.29
39895701	Replacement blade	each	\$134.29	\$134.29	\$120.86	\$120.86

Bone Scraper, Angulated

39895900	Angulated	each	\$330.79	\$330.79	\$281.17	\$281.17
39895901	Replacement blade	each	\$134.29	\$134.29	\$120.86	\$120.86



OSSEOUS TECHNOLOGIES OSSEOUS GLIDER

Disposable bone scraper

- Always has a fresh blade
- 360° blade is marked in 4 sections
- Expose new cutting surface with the twist of any .050 hex driver
- Recommended for the external oblique ridge
- Hold at a 35° angle on the ridge.

BG205	Osseous Glider	each	\$54.99	\$54.99	\$49.49	\$49.49
-------	----------------	------	---------	---------	---------	---------

INTEGRA® BONE BASIN

- Stainless steel deep dish with weighted bottom to prevent tipping
- Holds 20 grams of wetted material



6271	Bone Basin	each	\$68.29	\$68.29	\$61.46	\$61.46
------	------------	------	---------	---------	---------	---------

ACE BONE CRUSHER

The Bone Crusher introduces an alternative to the Bone Mill. Insert the bone piece into the crusher and use the pestle in order to crush the bone. If necessary, you may additionally use a mallet.

- Stainless steel
- Crushing surface: 20mm



9600050	Bone Crusher	each	\$279.29	\$279.29	\$251.36	\$251.36
---------	--------------	------	----------	----------	----------	----------

ACE GRAFT DELIVERY SYRINGES

- With funnels
- Stainless steel




Syringe openings, available sizes



9800035	3.5mm, curved	each	\$128.79	\$128.79	\$115.91	\$115.91
9800045	4.5mm, curved	each	\$128.79	\$128.79	\$115.91	\$115.91
9600025	5.0mm, straight	each	\$128.79	\$128.79	\$115.91	\$115.91
9800060	6.0mm, curved	each	\$128.79	\$128.79	\$115.91	\$115.91

A TITAN BONE FUNNELS

- Stainless steel
- Produced in the USA
- Thumb rings are color coded to indicate size
- Graft material loads at both ends to prevent clogging



3.5mm - Holds 1.55 cc of Bone Material (Blue)

BFS35MM	Straight	each	\$93.79	\$93.79	\$84.41	\$84.41
BFC35MM	Curved	each	\$93.79	\$93.79	\$84.41	\$84.41

4.5mm - Holds 2.0 cc of Bone Material (Red)

BFS45MM	Straight	each	\$93.79	\$93.79	\$84.41	\$84.41
BFC45MM	Curved	each	\$93.79	\$93.79	\$84.41	\$84.41

6.0mm - Holds 2.6 cc of Bone Material (Black)

BFS60MM	Straight	each	\$93.79	\$93.79	\$84.41	\$84.41
BFC60MM	Curved	each	\$93.79	\$93.79	\$84.41	\$84.41

ACE STERILE DISPOSABLE SYRINGES

- Carries material to surgical site
- Single use
- Can be used for delivery of any material to surgical site



0032019		10/pack	\$43.79	\$43.79	\$39.41	\$39.41
---------	--	---------	---------	---------	---------	---------

EXEL DISPOSABLE BONE GRAFTING SYRINGE

- 1 cc sterile bone delivery syringe used for delivery of material to the surgical site
- Disposable



OSTD51		10/pack	\$24.99	\$23.74	\$23.74	\$22.49
--------	--	---------	---------	---------	---------	---------

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

BONE GRAFTING

3

BONE GRAFT PLUGGERS, COMPACTORS AND CONDENSERS

A-Titan Bone Graft Spoon & Pluggers

To carry and pack bone



SINUS7	each	\$106.29	\$106.29	\$95.66	\$95.66
SINUS7S With Smaller Spoon	each	\$106.29	\$106.29	\$95.66	\$95.66

Hu-Friedy® Black Line Grisdale Bone Condenser

PLGGR1X	each	\$67.99	\$67.99	\$61.19	\$61.19
---------	------	---------	---------	---------	---------



Hu-Friedy® Black Line Labanca Plugger

PLGLABANCX 3mm x 1.5mm	each	\$59.29	\$59.29	\$53.36	\$53.36
------------------------	------	---------	---------	---------	---------



Integra™ Bone Compactors

Features a serrated tip and bent angle to reach difficult areas

6261 6 and 8mm	each	\$49.29	\$49.29	\$44.36	\$44.36
MIL6260 3 and 4mm	each	\$49.29	\$49.29	\$44.36	\$44.36



BONE GRAFT MATERIAL SPOONS

Integra™ Bone Graft Material Spoon

6250	each	\$56.79	\$56.79	\$51.11	\$51.11
------	------	---------	---------	---------	---------



Integra™ Bone Graft Material Spoon

6251 With 3mm serrated packer	each	\$112.29	\$112.29	\$101.06	\$101.06
-------------------------------	------	----------	----------	----------	----------



Integra™ Bone Graft Material Spoon

6252 With 4mm serrated packer	each	\$112.29	\$112.29	\$101.06	\$101.06
-------------------------------	------	----------	----------	----------	----------



MEMBRANE PLACEMENT INSTRUMENTS

ACE Buser-type Elevator

0300139 7"	each	\$39.49	\$39.49	\$35.54	\$35.54
------------	------	---------	---------	---------	---------



Hu-Friedy® Membrane Placement Instrument

PMPI	each	\$59.29	\$59.29	\$53.36	\$53.36
------	------	---------	---------	---------	---------



Hu-Friedy® Black Line Membrane Placement Instrument

PMPIX	each	\$83.29	\$83.29	\$74.96	\$74.96
-------	------	---------	---------	---------	---------



Integra™ Membrane Placement Instrument

6257	each	\$60.79	\$60.79	\$54.71	\$54.71
------	------	---------	---------	---------	---------



BONE GRAFT AND TISSUE FORCEPS

ACE Bone Graft Holding Forceps

9600026 8" (20cm)	each	\$108.29	\$108.29	\$97.46	\$97.46
9600027 6" (15cm)	each	\$108.29	\$108.29	\$97.46	\$97.46



ACE Kuehne Membrane Holding Forceps

• 4-1/4" (10.5cm)					
9600322 Angled, no lock	each	\$36.99	\$36.99	\$33.29	\$33.29
9600324 Straight, with lock	each	\$67.99	\$67.99	\$61.19	\$61.19
9600326 Straight, no lock	each	\$36.99	\$36.99	\$33.29	\$33.29



ACE Flap Holding Plier

9600330 5-3/4" (14.5cm) serrated	each	\$58.79	\$58.79	\$52.91	\$52.91
----------------------------------	------	---------	---------	---------	---------



Integra™ Tissue Graft Forceps

6256 6" (15.4cm)	each	\$60.79	\$60.79	\$54.71	\$54.71
------------------	------	---------	---------	---------	---------



BONE GRAFTING

3

BONE GRAFT CHISELS

Hu-Friedy® 1 Chandler Bone Chisel

CC1	5mm mono-bevel tip	each	\$55.29	\$55.29	\$49.76	\$49.76
-----	--------------------	------	---------	---------	---------	---------



Hu-Friedy® Black Line Periodontal Bone Chisel - 1/2 Buser Modified

PPBUSER12X	3mm/4mm	each	\$89.79	\$89.79	\$80.81	\$80.81
------------	---------	------	---------	---------	---------	---------



Hu-Friedy® Black Line Periodontal Bone Chisel - 1/3 Buser Modified

PPBUSER13X	3mm/4mm	each	\$89.79	\$89.79	\$80.81	\$80.81
------------	---------	------	---------	---------	---------	---------



Hu-Friedy® Black Line Periodontal Bone Chisel - 2/3 Buser Modified

PPBUSER23X	4mm/5mm	each	\$89.79	\$89.79	\$80.81	\$80.81
------------	---------	------	---------	---------	---------	---------



Integra™ Graduated Chisels

- Feature reference graduation marks used to precisely widen ridges.
- With 7-10-13-15-18mm Markings

6220	3.8mm, straight	each	\$59.79	\$59.79	\$53.81	\$53.81
6221	3.8mm, curved	each	\$59.79	\$59.79	\$53.81	\$53.81

Integra™ Freer Chisel - Single Cut Blade

6241	4mm	each	\$57.29	\$57.29	\$51.56	\$51.56
------	-----	------	---------	---------	---------	---------



Integra™ Freer Chisel - Double Cut Blade

MIL6243	2mm	each	\$57.29	\$57.29	\$51.56	\$51.56
---------	-----	------	---------	---------	---------	---------



MIL6244	4mm	each	\$57.29	\$57.29	\$51.56	\$51.56
---------	-----	------	---------	---------	---------	---------



6245	6mm	each	\$57.29	\$57.29	\$51.56	\$51.56
------	-----	------	---------	---------	---------	---------




ACE MISCH RAMUS RETRACTORS

Unique design to retract tissues of the posterior mandible in preparation for obtaining graft from lateral mandibular external oblique ridge.

- Curved Down


0502056	8-½" (21cm) right	each	\$120.49	\$120.49	\$108.44	\$108.44
0502057	8-½" (21cm) left	each	\$120.49	\$120.49	\$108.44	\$108.44



LOOKING FOR
ADDITIONAL
INSTRUMENTS?

Flip over to the
Instruments Section
starting on Page 156

And for the most
comprehensive instrument
listing, please visit
acesouthern.com today



BONE GRAFTING

3



ACE SINUS LIFT OSTEOTOME KIT

- All osteotomes are 6", 15cm
- Concave tip, sharply honed cutting surface
- Concave tip can be used to carry and place bone grafting material into site
- All osteotomes and tray may be purchased separately.

Kit Includes:

- Sinus Lift Osteotomes Step**
- 4545001 1.6 mm/2.0 mm Step
 - 4545101 2.0 mm/2.8 mm Step
 - 4545201 2.8 mm/3.3 mm Step
 - 4546001 3.3 mm/4.0 mm Step
 - 4546101 4.0 mm/4.35 mm Step
 - 4546301 4.35 mm/5.0 mm Step

- Sinus Lift Osteotomes Final**
- 4545701 2.8 mm Final
 - 4545301 3.3 mm Final
 - 4545401 4.0 mm Final
 - 4546201 4.35 mm Final
 - 4546401 5.0 mm Final
 - 7651060 Autoclavable Tray



ACE Sinus Lift Osteotome Kit

4545000	Set of 11 with tray	each	\$805.29	\$805.29	\$724.76	\$724.76
---------	---------------------	------	----------	----------	----------	----------

Also Available, Not Included in the Kit

4546701	5.0 mm/5.4 mm ftep	each	\$73.29	\$73.29	\$65.96	\$65.96
4546801	5.4 mm final	each	\$73.29	\$73.29	\$65.96	\$65.96




ACE SITE DILATOR KIT

- All dilators are 6", 15cm
- Convex tip compresses bone, limiting the amount of bone that must be removed from the site.
- All dilators and tray may be purchased separately.

Kit Includes:

- Site Dilating Osteotomes Step**
- 4543001 1.6 mm/2.0 mm Step
 - 4543101 2.0 mm/2.8 mm Step
 - 4543201 2.8 mm/3.3 mm Step
 - 4544001 3.3 mm/4.0 mm Step
 - 4544101 4.0 mm/4.35 mm Step
 - 4544301 4.35 mm/5.0 mm Step

- Site Dilating Osteotomes Final**
- 4543701 2.8 mm Final
 - 4543301 3.3 mm Final
 - 4543401 4.0 mm Final
 - 4544201 4.35 mm Final
 - 4544401 5.0 mm Final
 - 7651055 Autoclavable Tray



ACE Site Dilator Kit

4543000	Set of 11 with tray	each	\$805.29	\$805.29	\$724.76	\$724.76
---------	---------------------	------	----------	----------	----------	----------

ACE OFFSET OSTEOTOME KIT

- All Osteotomes are 6", 15cm
- Concave Tip Design - See Above Detail
- All osteotomes and tray may be purchased separately.

Kit Includes:

Sinus Lift Osteotomes Step

- 4542001 1.6 mm/2.0 mm Step
- 4542101 2.0 mm/2.8 mm Step
- 4542201 2.8 mm/3.3 mm Step
- 4542301 3.3 mm/4.0 mm Step

Sinus Lift Osteotomes Final

- 4542401 2.8 mm Final
- 4542601 3.3 mm Final
- 4542801 4.0 mm Final
- 7651080 Autoclavable Tray



ACE Offset Osteotome Kit

4542000	Set of 7 with tray	each	\$763.79	\$763.79	\$687.41	\$687.41
---------	--------------------	------	----------	----------	----------	----------

Also Available, Not Included in the Kit

4542501	4.0 mm/4.35 mm Step	each	\$94.99	\$94.99	\$85.49	\$85.49
4547101	4.35 mm/5.0 mm Step	each	\$94.99	\$94.99	\$85.49	\$85.49
4542901	4.35 mm final	each	\$94.99	\$94.99	\$85.49	\$85.49
4547201	5.0 mm final	each	\$94.99	\$94.99	\$85.49	\$85.49

INTEGRA™ CONCAVE OSTEOTOMES WITH STOPS

- Concave osteotomes break the floor of the maxillary sinus through the implant cavity
- Concave tips also allow for introduction of graft material
- Curved, with 6-11-15mm Markings

62152	2mm	each	\$81.79	\$81.79	\$73.61	\$73.61
62153	3mm	each	\$81.79	\$81.79	\$73.61	\$73.61
MIL62154	3.5 mm	each	\$81.79	\$81.79	\$73.61	\$73.61
62155	4mm	each	\$81.79	\$81.79	\$73.61	\$73.61



BONE GRAFTING

3



ACE SINUS LIFT KIT

- Kit Includes:** shown in order top to bottom
- 03-110-88 Shepherd Sinus Lift Retractor
 - 03-110-23 Freer Periosteal Elevator-Strong Curve
 - 03-110-77 Clepper Tatum Antralplasty Elevator
 - 10-002-36 Sinus Lift Curette-Standard Bend, 3mm/5mm
 - 10-002-37 Sinus Lift Curette-Severe Bend, 3mm/5mm
 - 714-0243 Depth Packing Gauge (Titanium), DE
 - 765-1050 Tray

Tray Included. All instruments may be purchased separately

ACE Sinus Lift Kit

0034194	Set of 6 with tray	each	\$401.49	\$401.49	\$361.34	\$361.34
---------	--------------------	------	----------	----------	----------	----------

A-TITAN SINUS LIFT INSTRUMENT-5



90/90/60 degree sinus elevator has a reverse angle and has extremely smooth tips to ensure the sinus wall isn't punctured when separating the lining from the wall. The handle is anodized red. Made in the USA from 100% 440A USA stainless steel.

SINUS5	each	\$101.29	\$101.29	\$91.16	\$91.16
--------	------	----------	----------	---------	---------

HU-FRIEDY® BLACK LINE SINUS LIFT INSTRUMENTS



Hu-Friedy® Black Line Kramer-Nevins Sinus Lift

IMP6578SX	Obtuse	each	\$123.29	\$123.29	\$110.96	\$110.96
IMP6577SCX	Acute	each	\$111.79	\$111.79	\$100.61	\$100.61



Hu-Friedy® Black Line Sinus Curette 1

SINC1X		each	\$103.49	\$103.49	\$93.14	\$93.14
SINC1SX	Short handle	each	\$103.49	\$103.49	\$93.14	\$93.14



Hu-Friedy® Black Line Sinus Curette 2

SINC2X		each	\$103.49	\$103.49	\$93.14	\$93.14
SINC2SX	Short handle	each	\$103.49	\$103.49	\$93.14	\$93.14

INTEGRA™ SINUS ELEVATORS

- Designed to aid in lifting the sinus membrane during an implant site preparation, the sinus lift instruments have a lightweight, wider diameter stainless steel handle
- Feature a more round, smooth tip for lifting the sinus with minimal puncture risk

62180	lift no. 1	each	\$52.79	\$52.79	\$47.51	\$47.51
-------	------------	------	---------	---------	---------	---------

62181	lift no. 2	each	\$52.79	\$52.79	\$47.51	\$47.51
-------	------------	------	---------	---------	---------	---------

62182	lift no. 3	each	\$52.79	\$52.79	\$47.51	\$47.51
-------	------------	------	---------	---------	---------	---------

62183	lift no. b1	each	\$49.99	\$49.99	\$44.99	\$44.99
-------	-------------	------	---------	---------	---------	---------

62184	lift no. b2	each	\$49.99	\$49.99	\$44.99	\$44.99
-------	-------------	------	---------	---------	---------	---------

62185	lift no. b3	each	\$49.99	\$49.99	\$44.99	\$44.99
-------	-------------	------	---------	---------	---------	---------

BONE GRAFTING

3

ACE RIDGE EXPANDING KIT

- Discover a safe, controllable way to widen a thin questionable ridge to one capable of accommodating an implant up to 4.5mm without compromising cortical bone.
- The slowly revolving dilators of increasing diameter gently widen the bone in a risk free procedure to the exact depth and direction of the final implant.
- After each dilator is inserted, allow 30-60 seconds before removing it. This allows time for the elasticity of the bone to respond.
- Laser markings provide accurate guidance for the placement of implants into the newly expanded ridge.



ACE Ridge Expanding Kit

05560000	Set of 4 with surgical kit	each	\$621.79	\$621.79	\$559.61	\$559.61
----------	----------------------------	------	----------	----------	----------	----------

Kit Components Individually

05560201	2.5mm Expander	each	\$54.29	\$54.29	\$46.15	\$46.15
05560301	3.1mm Expander	each	\$54.29	\$54.29	\$46.15	\$46.15
05560401	3.6mm Expander	each	\$54.29	\$54.29	\$46.15	\$46.15
05560501	4.3mm Expander	each	\$54.29	\$54.29	\$48.86	\$48.86
04973401	Connection driver tool (Hex), for CA handpiece	each	\$102.79	\$102.79	\$87.37	\$87.37
9600228	Universal ratchet	each	\$154.99	\$154.99	\$131.74	\$131.74
20155811	Driver extender (Hex)	each	\$102.79	\$102.79	\$92.51	\$92.51
9600449	Square titanium bowl	each	\$105.99	\$105.99	\$95.39	\$95.39
20195012	Autoclavable tray	each	\$133.99	\$133.99	\$120.59	\$120.59

ACE RIDGE EXPANDING CHISEL KIT

- All Ridge Expanding Chisels are 6", 15cm
- All chisels and tray may be purchased separately.
- Markings at 8, 10, 13, 15 and 18mm



ACE Ridge Expanding Chisel Kit

4544700	Chisel kit	kit	\$336.79	\$336.79	\$303.11	\$303.11
---------	------------	-----	----------	----------	----------	----------

Kit Components Individually

4544701	3.0 mm Ridge Chisel	each	\$72.29	\$72.29	\$65.06	\$65.06
4544801	4.0 mm Ridge Chisel	each	\$72.29	\$72.29	\$65.06	\$65.06
4544901	5.0 mm Ridge Chisel	each	\$72.29	\$72.29	\$65.06	\$65.06
4544501	9.5 mm Ridge Chisel	each	\$72.29	\$72.29	\$65.06	\$65.06
4541501	Chisel Handle (Uses 1.25mm Hex Driver)	each	\$56.29	\$56.29	\$50.66	\$50.66
7651090	Autoclavable Tray	each	\$91.99	\$91.99	\$82.79	\$82.79



FIBRIJET® RATIO APPLICATOR

FibriJet® Ratio Applicator Assembly 11:1 Ratio

- Designed to facilitate biomaterials application by isolating the two components until they passively mix at the treatment site

SA1001	Single use, sterile	each	\$38.29	\$36.38	\$36.38	\$34.46
--------	---------------------	------	---------	---------	---------	---------

FibriJet® Stream Tip

- Malleable metal shaft can be bent or shaped for better access to treatment site
- Luer lock connections prevent accidental dislodging of applicator tips
- Dispenses in a stream
- Dual Cannula - 20G x 2"

SA3615	Single use, sterile	each	\$38.29	\$36.38	\$36.38	\$34.46
--------	---------------------	------	---------	---------	---------	---------

FibriJet® Spray Tip

- Applies two liquids simultaneously to a treatment site
- Offers maximum versatility by providing multiple options for spraying an assortment of materials to open procedures
- Dual Spray

SA3660	Single use, sterile	each	\$38.29	\$36.38	\$36.38	\$34.46
--------	---------------------	------	---------	---------	---------	---------

BONE GRAFTING

3



Tubes shown are not included with centrifuge

LABNET SPECTRAFUGE™ 6C CENTRIFUGE

The Spectrafuge™ 6C has been designed for quick production of PPP (platelet poor plasma) and PRP (platelet rich plasma) as well as other applications in the clinical or research laboratory.

The removable fixed angle (30°) rotor of the Spectrafuge™ 6C accepts 6x15ml tubes, and with the adapters will accommodate all common blood draw tubes. The variable speed allows for a variety of applications to be performed. A timer controls operation from 30 seconds to 30 minutes. The time can also be set to continuous, by selecting the "On" position, or momentary, by pressing and holding the start knob.

- Supplied with 6 x 10/15ml rotor
- Compact design: 8.25"(W) x 9.5"(L) x 7" (H)
- Variable speed: 300 rpm to 6,500 rpm with digital control, Max. g-force: 4,000 x g
- Capable of centrifuging 15ml conical and all common blood draw tubes: 15ml (17 x 125mm), 10ml (16 x 100mm), 7ml (13 x 100mm) and 5ml (12 x 75mm)
- Warranty: 1 yr. warranty against defects in material and workmanship

C0060	each	\$1,358.79	\$1,290.85	\$1,290.85	\$1,222.91
-------	------	------------	------------	------------	------------

ACE PRP (PLATELET RICH PLASMA) BLOOD COLLECTION KIT

Kit Includes:	(Qty) included in kit	Qty when re-ordering
• 7287	(1) - Vacutainer Safetylok Blood Collection Set, 21G x 3/4", 7" tubing, sterile	50/box
• 368607	(1) - Safety Needle - 21G x 1-1/4"	48/box
• 364815	(1) - Vacutainer® One Use Needle Holder, latex-free	250/bag
• 364606	(6) - Whole Blood Tube with ACD, yellow, 8.5ml, sterile	100/box
• 366430	(6) - Blood Tube For Serum, red, 10ml, sterile	100/box
• 309646	(3) - Syringe Only, 5cc, sterile, no cap	125/box
• 0900269	(2) - Blunt Cannula, 18G x 2-1/2", sterile *	each
• 0900270	(1) - Blunt Cannula, 18G x 3", sterile *	each
• 1089	(1) - Drape, 22" x 25", sterile	2/pkg
• 3322	(1) - Gauze Sponge, 2" x 2", 8 ply, sterile	3,000/case
• 3601	(2) - Adhesive Sheer Strip Bandage, 3/4" x 3"	100/box
• 6818	(2) - Alcohol Prep Pad, 1-1/4" x 1-1/2"	200/box
• 0031153	(1) - Tourniquet, 1" x 18", latex-free	10/pkg

* Item is sterile in this PRP kit, but will be non-sterile when re-ordered individually. See our web site for replacement products images & descriptions.

PRP Collection Kit

• Single use, disposable, non-sterile					
0031702	each	\$42.49	\$40.37	\$40.37	\$38.24



GLOBAL SCIENTIFIC INC. BLOOD COLLECTION TUBE RACK

- 36 16mm-holes
- Autoclavable
- 5" x 5" x 2.7"

456545B	Blue	each	\$30.49	\$28.97	\$28.97	\$27.44
---------	------	------	---------	---------	---------	---------



COSMETIC AND RESTORATIVE DENTISTRY

**BOSWORTH CORA-CAINE®
DENTURE PAIN-RELIEVING
ADHESIVE OINTMENT**

Cora-Caine® is a pain-relieving dental analgesic adhesive ointment for new or intermediate dentures. The active ingredient is Benzocaine 16%. Cora-Caine® provides immediate relief and a secure fit by soothing tender ridges and relieving the discomfort of sore spots. Acting as a gel-like cushion between the denture and the tissue, Cora-Caine is applied directly onto the denture to provide positive retention and enhance stability.



16621	4g Tube	36/box	\$70.79	\$67.25	\$67.25	\$63.71
-------	---------	--------	---------	---------	---------	---------

**UNIPACK CAVEX
TEMPORARY CEMENT KIT**

- Short setting time, 1-1/2 - 2 minutes in mouth
- For cementation of temporary crowns and bridges
- Easy and convenient use
- Yellowish color is easy to detect and remove
- Eugenol-free and will not inhibit the setting of composite cement or composite filling material



Kit Includes: 1 Tube 35g Base, 1 Tube 16g Catalyst and Mixing Pad

CE090	kit	\$21.79	\$20.70	\$20.70	\$19.61
-------	-----	---------	---------	---------	---------

VISTA DENTAL BEST-ETCH™

(37% Phosphoric Acid Etch)

- Applicator tips allow for precise delivery
- Designed to etch and prepare dentin and enamel



Kit Includes: (4) 1.2ml Prefilled Syringes, (25) 22g Pre-Bent Applicator Tips

502115	kit	\$19.79	\$18.80	\$18.80	\$17.81
--------	-----	---------	---------	---------	---------

SEPTODONT SEPTOBOND™

Light Cure Adhesive

Universal bonding agent with an acetone base that offers greater effectiveness than ethanol. With SeptoBond, sensitivity is minimized.



Kit Includes: (4) 1.2ml Prefilled Syringes, (25) 22g Pre-Bent Applicator Tips

01C2030	kit	\$59.29	\$56.33	\$56.33	\$53.36
---------	-----	---------	---------	---------	---------

**HAGER WORLDWIDE QUICK-TIPS™
MICRO-APPLICATOR TIPS**

- For accurate application of minute amounts of liquids and pastes: cavity lining, etching, bonding, fluoride varnishing, hemostat and disclosing solutions
- Lint-free, non-absorbing fiber tips release practically all liquid
- Flexible tapered ends maintain desired angle for application in difficult areas
- Hygienic and convenient dispenser pack



605220	100/pack	\$15.99	\$15.19	\$15.19	\$14.39
--------	----------	---------	---------	---------	---------

**VISTA DENTAL
MICRO-APPLICATOR
BRUSHES**



- Ideal for applying etchants, sealants, bonding agents and hemostatic agents
- Bendable micro-applicator
- Lint-free, non-absorbant brush tip

Regular

320802	Silver	100/box	\$12.29	\$11.68	\$11.68	\$11.06
--------	--------	---------	---------	---------	---------	---------

4

**COLTÈNE WHALEDENT COLTOLUX®
LED CURING LIGHT**

Slim, powerful, cordless

The ColtLux LED combines a slim curing probe with high output in a cordless lightweight pen-style design.



- Cures 2mm in only 10 seconds
- Dome shaped lens minimizes light diffusion
- LED light output matches the absorption spectrum of Camphorquinone (CQ) Photoinitiators
- Audible signal every 10 seconds and automatic power saving sleep mode after 20 seconds of operation
- Hundreds of 10-second cures with a fully charged battery
- Test block provided for determining accurate cure times

Includes: Handpiece, Charger Base, Power Supply, 3 Light Shields, 5 Light Lenses, 100 Barrier Sleeves, Curing Test Fixture

C7970100115	each	\$508.49	\$483.07	\$483.07	\$457.64
-------------	------	----------	----------	----------	----------

Barrier Sleeves

C7957	1000/pack	\$49.99	\$47.49	\$47.49	\$44.99
-------	-----------	---------	---------	---------	---------

Snap-On Amber Light Shields

C7990	25/box	\$64.29	\$61.08	\$61.08	\$57.86
-------	--------	---------	---------	---------	---------

PALMERO HEALTHCARE BONDING LENS SHIELD

This shield easily adapts to almost any pistol or wand type VLCs (visible light cure). Helps to protect eyes from the harmful blue light used in curing and bleaching procedures. Constructed of an unbreakable polycarbonate with protection of 530nm.



125	each	\$17.79	\$16.90	\$16.90	\$16.01
-----	------	---------	---------	---------	---------

**PALMERO HEALTHCARE
OVERSHIELD EYEWEAR**

Attractive over-the-glass style. Fits over most prescription eyewear. High-impact polycarbonate lenses with built-in side shield. Anti-fog coated lens. Certified ANSI Z87.1-2003 (high impact) and CSA Z94.3-2007. Filters over 95% of blue and UV light. Black frame/reddish-orange (bonding) lens.



3602BOND	each	\$13.29	\$12.63	\$12.63	\$11.96
----------	------	---------	---------	---------	---------

DISPOSABLES

CROSSTEX PATIENT TOWEL
PRODUCT INFORMATION

- **Polyback®** - 3-ply of soft facial tissue bonded to a layer of poly film results in a towel with maximum absorbency qualities and a reliable moisture barrier
- **Polygard®** - 3" longer than Polyback - provides 25% more protection
- **Polygard® Deluxe** - Largest 3-ply plus towel, provides the maximum area of protection
- **Polyback® Contour Bib** - Oversized 3-ply plus poly towel with a cut-out that allows it to fit the contour of the patient's neck, providing a more secure fit
- **Econoback®** - 2-ply of unsurpassed facial quality tissue plus a layer of poly film offer quick absorbency with a reliable moisture barrier
- **UltraGard®** - Same as Econoback towels, but 3" longer - provides 25% more protection
- **Advantage & Advantage Plus®** - Heavy basis weight, facial quality tissue; These towels are specially treated to maximize both strength and absorption
- **Proback®** - Constructed of an extra heavy ply of facial tissue bonded to a layer of poly film; Provides an excellent absorbent surface backed by a reliable moisture barrier
- **Professional Regular®** - 3-ply of facial quality tissue with unique properties of wet or dry tensile strength combined with facial softness and super absorbency
- **Professional Thrift®** - 2-ply made from the same facial quality tissue as Professional Regular towels; Economical alternative with multiple uses
- **Sani-Tab Chain-Free®** - Exclusive patented adhesive tabs eliminate the need for chains forever; Simply position the towel anywhere on the patient and apply pressure to the tab; Maximize infection control and patient comfort while saving time and money

3-Ply Tissue with Poly Patient Towels

Polyback® Crosstex Horizontal

19" x 13"	500/case	\$32.79	\$31.15	\$31.15	\$29.51
-----------	----------	---------	---------	---------	---------

WPXAQ	Aqua	WPXBG	Beige	WPXBL	Blue	WPXDR	Dusty rose
WPXGR	Green	WPXLV	Lavender	WPXPE	Peach	WPXSI	Silver grey
WPXWH	White	WPXYE	Yellow				

Polyback® Crosstex Horizontal

19" x 16"	500/case	\$43.29	\$41.13	\$41.13	\$38.96
-----------	----------	---------	---------	---------	---------

WGXBL	Blue	WGXWH	White
-------	------	-------	-------

Polygard® Deluxe Vantage Pebble

18" x 30"	300/case	\$46.99	\$44.64	\$44.64	\$42.29
-----------	----------	---------	---------	---------	---------

W3D1930BL	Blue
-----------	------

Polyback® Contour Vantage Pebble

18" x 22"	400/case	\$45.99	\$43.69	\$43.69	\$41.39
-----------	----------	---------	---------	---------	---------

WERUECBL	Blue	WERUECDR	Dusty rose	WERUECGR	Green
WERUECWH	White				

Advantage Plus® Vantage Vertical

18" x 13"	500/case	\$31.79	\$30.20	\$30.20	\$28.61
-----------	----------	---------	---------	---------	---------

WPXAAQ	Aqua	WPXABG	Beige	WPXABL	Blue	WPXADR	Dusty rose
WPXAGR	Green	WPXALV	Lavender	WPXAYE	Yellow		

Sani-Tab Chain-Free® Crosstex Horizontal

19" x 13"	400/case	\$43.99	\$41.79	\$41.79	\$39.59
-----------	----------	---------	---------	---------	---------

WPXBTL	Blue	WPXDRT	Dusty rose	WPXWHT	White
--------	------	--------	------------	--------	-------

Patient Towel Color Chart



Patient Towel Embossing Patterns



Patient Towel Folds



VANTAGE FOLD



CROSSTEX FOLD

Sani-Tab



2-Ply Tissue with Poly Patient Towels

Econoback® Crosstex Horizontal

19" x 13"	500/case	\$27.79	\$26.40	\$26.40	\$25.01
-----------	----------	---------	---------	---------	---------

WEXAQ	Aqua	WEXBG	Beige	WEXBL	Blue	WEXDR	Dusty rose
WEXGR	Green	WEXLV	Lavender	WEXPE	Peach	WEXWH	White
WEXYE	Yellow						

WEXYE	Pink ribbons	\$27.79	\$26.40	\$26.40	\$25.01
-------	--------------	---------	---------	---------	---------

UltraGard® Crosstex Horizontal

19" x 16"	500/case	\$39.49	\$37.52	\$37.52	\$35.54
-----------	----------	---------	---------	---------	---------

WUXBL	Blue	WUXWH	White	WUXYE	Yellow
-------	------	-------	-------	-------	--------

Advantage Vantage Vertical

18" x 13"	500/case	\$25.79	\$24.50	\$24.50	\$23.21
-----------	----------	---------	---------	---------	---------

WEXABG	Beige	WEXABL	Blue	WEXADR	Dusty rose	WEXAGR	Green
WEXALV	Lavender	WEXASI	Silver grey				

Sani-Tab Chain-Free® Crosstex Horizontal

19" x 13"	400/case	\$35.29	\$33.53	\$33.53	\$31.76
-----------	----------	---------	---------	---------	---------

WEXBLT	Blue	WEXDRT	Dusty rose	WEXLVT	Lavender
WEXWHT	White				

Extra Heavy Tissue with Poly Patient Towels

Proback® Crosstex Horizontal

19" x 13"	500/case	\$23.29	\$22.13	\$22.13	\$20.96
-----------	----------	---------	---------	---------	---------

WEPBAQ	Aqua	WEPBBL	Blue	WEPBDR	Dusty rose
WEPBGR	Green	WEPBLV	Lavender	WEPBPE	Peach
WEPBSI	Silver grey	WEPBWH	White	WEPBYE	Yellow

3-Ply Tissue Patient Towels

Professional® Regular Crosstex Horizontal

19" x 13"	500/case	\$25.79	\$24.50	\$24.50	\$23.21
-----------	----------	---------	---------	---------	---------

WTXBG	Beige	WTXBL	Blue	WTXDR	Dusty rose
WTXGR	Green	WTXWH	White	WTXYE	Yellow

2-Ply Tissue Patient Towels

Professional® Thrift Crosstex Horizontal

19" x 13"	500/case	\$23.49	\$22.32	\$22.32	\$21.14
-----------	----------	---------	---------	---------	---------

WT5BL	Blue	WT5WH	White
-------	------	-------	-------

5

DISPOSABLES

TIDI® POLYBACKED TISSUE PATIENT TOWELS

- Offers strength and convenience with the addition of a leak-proof poly backing
- Waffle embossed for enhanced softness and absorbency
- Also available in yellow, gray, mauve, beige, peach and teal. Please call a customer service rep for details.



2-Ply							
13" x 18"	500/case	\$33.29	\$31.63	\$31.63	\$29.96		
917463	Blue	917462	Green	917459	Lavender	917461	White

3-Ply							
13" x 18"	500/case	\$35.79	\$34.00	\$34.00	\$32.21		
917403	Blue	917402	Green	917400	Lavender	917401	White

TIDI® TISSUE PATIENT TOWELS

3-Ply							
Pro-Towels Embossed							
918101	White, 18" x 13"	500/case	\$28.49	\$27.07	\$27.07	\$25.64	
918103	Blue, 18" x 13"	500/case	\$46.29	\$43.98	\$43.98	\$41.66	

TIDI® POLYBACKED TISSUE ORAL SURGERY BIBS

3-Ply							
18" x 25"	250/case	\$38.79	\$36.85	\$36.85	\$34.91		
917677	Blue	917678	Mauve				

TIDI® PROPHY BIBS TISSUE/POLY/TISSUE

- Creates optimal absorbency and barrier protection
- Contour neck extends to the waist for exceptional patient coverage



29" x 21"	100/case	\$30.79	\$29.25	\$29.25	\$27.71		
917483	Blue	917480	Lavender	917486	Mauve		
29" x 42", Knee Length	50/case	\$30.29	\$28.78	\$28.78	\$27.26		
917493	Blue						

TIDI® CHAINLESS TISSUE/POLY BIBS WITH TIES

18" x 25"	250/case	\$46.79	\$44.45	\$44.45	\$42.11		
917673	Blue	917676	Mauve				

UNIPACK POLYBACKED TISSUE PATIENT BIBS

- 2-ply tissue, 1-ply poly provides excellent absorbency
- High quality poly backing keeps patients dry
- Non-sterile

2-Ply							
13-1/2" x 18"	500/case	\$18.49	\$17.57	\$17.57	\$16.64		
UPB6101	Blue	UPB6102	White	UPB6103	Mauve		

PALMERO 15" CHAIN NAPKIN BIB HOLDER

- Autoclavable, chrome-plated alligator clips with durable spring and sharp tooth jaws
- Beautifully finished with ball chain or clear plastic "chill-free" sleeves
- Available in 6 colors (blue, white, yellow, green, red, and black)



310W	White	each	\$4.29	\$4.29	\$3.86	\$3.86	
Chill-Free							
311W	White	each	\$4.29	\$4.29	\$3.86	\$3.86	

BUFFALO DENTAL NAPKIN HOLDERS

'Tite Grip,' with plastic sleeves



51591X72	14" - Blue	each	\$6.29	\$6.29	\$5.66	\$5.66	
51621	14" - White	each	\$6.29	\$6.29	\$5.66	\$5.66	

BUFFALO DENTAL TITE GRIP

- Alligator-type clips securely attached to the bead chain
- All parts finished in polished chrome
- Can be sterilized in the usual manner



51451X72	10 1/4" - Chrome	each	\$4.49	\$4.49	\$4.04	\$4.04	
----------	------------------	------	---------------	---------------	---------------	---------------	--

ACE DENTAL PLASTIC NAPKIN HOLDER

9700044	Blue	each	\$14.99	\$14.99	\$13.49	\$13.49	
9700045	Yellow	each	\$14.99	\$14.99	\$13.49	\$13.49	



CERTOL INTERNATIONAL® PROBARRIER DISPOSABLE CHAIN™ DISPOSABLE BIB CHAIN

- Disposable - helps reduce the risk of cross-contamination. Use once and discard with each patient. Eliminates the need to disinfect metal bib chains.
- Economical - (300) 3/4" x 16-1/2" chains per box. 15% more strips than the competition.
- Patient Comfort - Soft, lightweight non-latex fiber fits comfortably on patient and eliminates cold chains that snag hair. Adhesive strips securely hold chain to bib.



DC0115S	300/box	\$14.79	\$14.05	\$14.05	\$13.31		
---------	---------	----------------	----------------	----------------	----------------	--	--

DISPOSABLES

PRECISION DYNAMICS PILLOW PROTECTOR - ZIPPER TYPE

- Plastibak Nylon

1503511MGD	21" x 27"	each	\$6.29	\$5.98	\$5.98	\$5.66
------------	-----------	------	--------	--------	--------	--------

CROSSTEX PILLOWCASES

- Facial tissue combined with poly for patient comfort and protection
- Soft, strong and absorbent
- Fluid resistant barrier between patient and pillow
- Latex-free
- White, 21" x 30"



PC2130	100/case	\$36.49	\$34.67	\$34.67	\$32.84
--------	----------	---------	---------	---------	---------

TIDI® PILLOWCASES

- Latex-free
- White, 21" x 30"



Spunbond Polypropylene

- Nonwoven, resists fluid

980915	100/case	\$91.49	\$86.92	\$86.92	\$82.34
--------	----------	---------	---------	---------	---------

Embossed Poly

- Provides a strong, leak-proof barrier for dependable infection control

919350	100/case	\$38.79	\$36.85	\$36.85	\$34.91
--------	----------	---------	---------	---------	---------

Fabricel

- Made of Fabricel tissue fused with a layer of polyethylene film on the back
- Very soft, super strong and tear-resistant
- Leak-proof for reliable barrier protection

919355	100/case	\$51.79	\$49.20	\$49.20	\$46.61
--------	----------	---------	---------	---------	---------

Tissue/Poly

- Provides a strong, leak-proof barrier for dependable infection control

919365	100/case	\$39.49	\$37.52	\$37.52	\$35.54
--------	----------	---------	---------	---------	---------

TIDI® EXAMINATION TABLE PAPER

- Wide variety - to fit every examination table
- All rolls electronically metered for exact length
- Individually poly-wrapped
- Bright, smooth white finish
- Strong, economical, and durable
- One side machine glazed for moisture repellence



980914	21" x 225'	12 rolls/case	\$67.49	\$64.12	\$64.12	\$60.74
--------	------------	---------------	---------	---------	---------	---------

DYNAREX PAPER TRAY COVERS

- Sized right for dental tray
- Minimizes noise
- Provides infection control barrier
- Rounded corners
- Easy dispenser box



DY2259	8¼" x 12¼" White	1000/box	\$19.79	\$18.80	\$18.80	\$17.81
2263	8¼" x 12¼" Blue	1000/box	\$19.79	\$18.80	\$18.80	\$17.81

CROSSTEX BRACKET TRAY COVERS

- Manufactured from a unique, heavy weight paper stock
- Impedes flow of moisture
- Minimizes seep-through



Size A - 9-1/2" x 12-1/4" - Weber Chayes

FABL	Blue	1000/box	\$38.29	\$36.38	\$36.38	\$34.46
FAGR	Green	1000/box	\$38.29	\$36.38	\$36.38	\$34.46
FAWH	White	1000/box	\$38.29	\$36.38	\$36.38	\$34.46

Size B - 8-1/2" x 12-1/4" - Ritter

FBAQ	Aqua	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBBG	Beige	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBBL	Blue	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBGR	Green	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBLV	Lavender	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBPE	Peach	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBSI	Silver	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
FBWH	White	1000/box	\$37.49	\$35.62	\$35.62	\$33.74

Size C - 11" x 17-1/2" - Weber Hi-Lo

FCBL	Blue	1000/box	\$60.29	\$57.28	\$57.28	\$54.26
FCGR	Green	1000/box	\$60.29	\$57.28	\$57.28	\$54.26
FCWH	White	1000/box	\$60.29	\$57.28	\$57.28	\$54.26

Size D - 10-1/4" x 15-3/4" - SS White #20

FDBL	Blue	1000/box	\$66.29	\$62.98	\$62.98	\$59.66
FDGR	Green	1000/box	\$66.29	\$62.98	\$62.98	\$59.66
FDWH	White	1000/box	\$66.29	\$62.98	\$62.98	\$59.66

Size E - 9" x 13-1/2" - Midwest

FEWH	White	1000/box	\$36.29	\$34.48	\$34.48	\$32.66
------	-------	----------	---------	---------	---------	---------

Size M - 5" x 8" - Mini

FMINI	White	1000/box	\$34.79	\$33.05	\$33.05	\$31.31
-------	-------	----------	---------	---------	---------	---------

Size X - 11" x 11" - Cox

FCOX	White	1000/box	\$37.49	\$35.62	\$35.62	\$33.74
------	-------	----------	---------	---------	---------	---------

CROSSTEX LIQUID-PROOF BRACKET TRAY COVERS

- Polyethylene-backed paper prevents seep-through
- Dual purpose: poly side down for moisture barrier, poly side up for quick mixing



Size B - 8-1/2" x 12-1/4" - Ritter

FBC	White	500/box	\$29.99	\$28.49	\$28.49	\$26.99
-----	-------	---------	---------	---------	---------	---------

Size F - 8-1/2" x 11" - ADEC

FF	White	500/box	\$28.49	\$27.07	\$27.07	\$25.64
----	-------	---------	---------	---------	---------	---------

DISPOSABLES

CROSSTEX POLY COATED MIXING PADS

- Extra heavy with high poly finish
- Non-skid foam base
- Use for mixing cements, composites, rubber base material and silicates



3" x 3"					
KPAD33	100 sheets/pad	\$3.29	\$3.13	\$3.13	\$2.96
3" x 6"					
KPAD36	100 sheets/pad	\$3.99	\$3.79	\$3.79	\$3.59
6" x 6"					
KPAD66	100 sheets/pad	\$6.99	\$6.64	\$6.64	\$6.29

HT BERRY 1/2 OZ SOUFFLE CUPS



FSWPSC050	1/2 oz	5000/case	\$92.79	\$88.15	\$88.15	\$83.51
-----------	--------	-----------	---------	---------	---------	---------

DYNAREX 1 OZ MEDICINE CUPS

- Manufactured of strong flexible polypropylene
- Easy to read graduations



4258	Economy	100/sleeve	\$1.49	\$1.42	\$1.42	\$1.34
4252	Standard	100/sleeve	\$6.49	\$6.17	\$6.17	\$5.84

DYNAREX PAPER CUPS

- Contemporary design to fit any environment
- Quality construction
- Packaged in sleeves of 100



4335	3 oz	100/sleeve	\$3.79	\$3.60	\$3.60	\$3.41
4336	5 oz	100/sleeve	\$4.49	\$4.27	\$4.27	\$4.04

CROSSTEX 5 OZ PLASTIC CUPS

- Double rolled lip for easy separation
- Unique embossing for greater strength and easy dispensing
- Double coated for lustrous, high gloss finish



CXBG	Beige	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXBL	Blue	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXDR	Dusty rose	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXGR	Green	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXLV	Lavender	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXPE	Peach	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXSI	Silver grey	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXWH	White	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXYE	Yellow	1000/case	\$36.49	\$34.67	\$34.67	\$32.84
CXCL	Clear	2500/case	\$96.79	\$91.95	\$91.95	\$87.11

DYNAREX 5 OZ. DRINKING CUPS

High value, high quality drinking cups made from high impact polypropylene that is both sturdy and shatterproof, yet flexible.



DY4236	White	1000/case	\$25.99	\$24.69	\$24.69	\$23.39
4237	Blue	1000/case	\$25.99	\$24.69	\$24.69	\$23.39

UNIPACK 5 OZ PLASTIC CUPS



UBC6201	Blue	1000/case	\$28.79	\$27.35	\$27.35	\$25.91
UBC6202	White	1000/case	\$28.79	\$27.35	\$27.35	\$25.91
UBC6203	Mauve	1000/case	\$28.79	\$27.35	\$27.35	\$25.91
UBC6204	Lavender	1000/case	\$28.79	\$27.35	\$27.35	\$25.91

KENDALL CURITY™ COTTON BALLS

- 100% cotton balls designed for a variety of applications
- Non-sterile



2600	Medium	500/bag	\$3.99	\$3.79	\$3.79	\$3.59
2601	Large	200/bag	\$4.49	\$4.27	\$4.27	\$4.04

DYNAREX COTTON-TIPPED APPLICATORS

- Cotton tips are highly absorbent
- Sterile applicators are packaged in convenient peel-down pouches
- Bag is autoclavable with a colored autoclaving process marker turning reddish pink to brown



Non-Sterile						
4301	3"	1000/box	\$4.49	\$4.27	\$4.27	\$4.04
4302	6"	1000/box	\$8.29	\$7.88	\$7.88	\$7.46
Sterile						
4303	3"	2/pk, 100 pks/box	\$4.49	\$4.27	\$4.27	\$4.04
4304	6"	1/pk, 100 pks/box	\$5.79	\$5.50	\$5.50	\$5.21

DYNAREX TONGUE DEPRESSORS

- Precision cut, polished smooth edges



Non-Sterile						
DY4311	Junior, 5-1/2"	500/box	\$7.29	\$6.93	\$6.93	\$6.56
DY4312	Regular, 6"	500/box	\$8.29	\$7.88	\$7.88	\$7.46
Sterile						
DY4313	Junior, 5-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
4314	Regular, 6"	100/box	\$6.29	\$5.98	\$5.98	\$5.66

DISPOSABLES

KIMBERLY CLARK PROFESSIONAL® PAPER PRODUCTS

SCOTT® C-Fold Towels

- White
- 10.125" x 13.15"
- 12 boxes of 200/case



01510	2400/case	\$45.79	\$43.50	\$43.50	\$41.21
-------	-----------	---------	---------	---------	---------

SCOTT® Multi-Fold Towels

- White
- 9.2" x 9.4"
- 6 boxes of 250/case



01804	4000/case	\$43.79	\$41.60	\$41.60	\$39.41
-------	-----------	---------	---------	---------	---------

KLEENEX® C-Fold Towels

- Dependable, folded towels that you can look to for softness, absorbency and performance
- FSC Certified Product
- 10.125" x 13.15"
- 16 boxes of 150/case



01500	2400/case	\$54.99	\$52.24	\$52.24	\$49.49
-------	-----------	---------	---------	---------	---------

KLEENEX® Multi-Fold Towels

- Soft and absorbent
- Meets EPA standards with a minimum of 40% post consumer waste
- FSC Certified Product
- White
- 9.2" x 9.4"
- 16 boxes of 150/case



01890	2400/case	\$53.79	\$51.10	\$51.10	\$48.41
-------	-----------	---------	---------	---------	---------

Scott® Rolls - Perforated Paper Towels

- White



41482	20 roll/case	\$49.29	\$46.83	\$46.83	\$44.36
-------	--------------	---------	---------	---------	---------

Scott® Facial Tissues - Surpass

- White
- 30 boxes of 100/sheets

21340	3000/case	\$54.29	\$51.58	\$51.58	\$48.86
-------	-----------	---------	---------	---------	---------

KLEENEX® Facial Tissue

- White
- 8.4" x 8.2"
- 36 boxes of 100/case



21400	3600/case	\$61.79	\$58.70	\$58.70	\$55.61
-------	-----------	---------	---------	---------	---------

- White
- 8.4" x 8"
- 48 boxes of 125/case

21606	6000/case	\$94.79	\$90.05	\$90.05	\$85.31
-------	-----------	---------	---------	---------	---------

KLEENEX COTTONELLE® - Bathroom Tissue

- FSC Certified Product
- White



17713	60 rolls/case	\$72.99	\$69.34	\$69.34	\$65.69
-------	---------------	---------	---------	---------	---------

Scott® Toilet Tissue

- 2-ply Standard Bathroom Tissue



04460	80 rolls/case	\$72.79	\$69.15	\$69.15	\$65.51
-------	---------------	---------	---------	---------	---------

MEDLINE LOW DENSITY WASTE CAN LINERS

- Assorted sizes for step-on cans, linen hampers or garbage drums
- Superior stretch resistance, ideal for heavier or sharper waste



Black, 40" x 46", 45 Gal, 1.8 mil

NONTGG46XXH	100/case	\$94.49	\$89.77	\$89.77	\$85.04
-------------	----------	---------	---------	---------	---------

Black, 40" x 55", 55 Gal, 1.8 mil

NONTGG55XXH	100/case	\$111.79	\$106.20	\$106.20	\$100.61
-------------	----------	----------	----------	----------	----------

Clear, 24" x 23", 7-10 Gal, 0.35 mil

NON022207	500/case	\$27.99	\$26.59	\$26.59	\$25.19
-----------	----------	---------	---------	---------	---------

White, 30" x 36", 20-30 Gal, 0.8 mil

NONXPW36X	200/case	\$53.49	\$50.82	\$50.82	\$48.14
-----------	----------	---------	---------	---------	---------

White, 33" x 39", 33 Gal, 0.75 mil

NONCRW39X	150/case	\$41.49	\$39.42	\$39.42	\$37.34
-----------	----------	---------	---------	---------	---------

HT BERRY WASTE CAN LINERS

- Black

30" x 36", 20-30 Gallon

ACT3036H	250/case	\$59.79	\$56.80	\$56.80	\$53.81
----------	----------	---------	---------	---------	---------

23" x 39", 31-33 Gallon

ACT3339H	200/case	\$59.79	\$56.80	\$56.80	\$53.81
----------	----------	---------	---------	---------	---------

38" x 58", 31-33 Gallon

ACT3858XH	100/case	\$49.29	\$46.83	\$46.83	\$44.36
-----------	----------	---------	---------	---------	---------

CLINTON INDUSTRIES PREMIUM RED WASTE RECEPTACLES

- Heavy gauge steel construction
- Ultra-durable, easy-clean, non-textured, smooth, powder-coated finish
- Full width lid hinge
- Heavy duty lid operating system
- Self-closing lid with bumper strips
- Interior bag holder
- Easy-to-locate, extra-wide, foot pedal
- Rubber base strip protects flooring



24 Quart (6-Gallons), 12" x 13.25" x 16.5"

TP24R	each	\$134.29	\$127.58	\$127.58	\$120.86
-------	------	----------	----------	----------	----------


32 Quart (8-Gallons), 12" x 13.25" x 21.5"

TP32R	each	\$138.79	\$131.85	\$131.85	\$124.91
-------	------	----------	----------	----------	----------

ULINE STORAGE BINS

- Polypropylene construction
- Autoclavable

Size Comparison



S12414BLU	Small, 7.5" x 4" x 3"	each	\$5.79	\$5.50	\$5.50	\$5.21
-----------	-----------------------	------	--------	--------	--------	--------

S12415BLU	Medium, 11" x 5.5" x 5"	each	\$10.49	\$9.97	\$9.97	\$9.44
-----------	-------------------------	------	---------	--------	--------	--------

S12416BLU	Large, 11" x 8" x 7"	each	\$15.99	\$15.19	\$15.19	\$14.39
-----------	----------------------	------	---------	---------	---------	---------

5

DISPOSABLES

CROSSTEX SPARKLE® PROPHY ANGLES

Features

- Smooth-running
- Vibration-free
- No noise
- No over-heating
- Non-locking

Improved Cup

- Dark colored soft cups - better visibility and contrast with teeth
- Perfect cup flare - improves subgingival access
- Reduced splatter - helps minimize cross contamination
- 100% latex-free

A More Ergonomic Design

- Non-slip contour - reduces hand fatigue and micro-traumas
- Light colored shell - for improved visibility in the mouth



Soft Cup, Purple

TPASSC1	100/box	\$41.99	\$39.89	\$39.89	\$37.79
TPASSC5	500/box	\$202.49	\$192.37	\$192.37	\$182.24

Firm Cup, White

TPASFC1	100/box	\$42.49	\$40.37	\$40.37	\$38.24
---------	---------	---------	---------	---------	---------

PREVENTECH PIVOT® DISPOSABLE PROPHY ANGLES

- Designed to fit all slow-speed, doriot (straight nose cone) style handpieces
- Works smoothly and quietly with less vibration, no noisy chatter
- Shorter profile cup/head allows improved access to hard-to-reach surfaces
- Non-slip grip



Non-Latex, Soft Cup

1100043	144/box	\$81.79	\$81.79	\$73.61	\$73.61
1100047	1200/box	\$584.79	\$584.79	\$526.31	\$526.31

DYNAREX ORAL SWABSTICKS

- Gently cleans, moistens and refreshes mouth, teeth and gums
- Raised foam ridges remove particles and mucous for enhanced oral hygiene
- Individually wrapped



1217	Mint flavored	250/box	\$27.29	\$25.93	\$25.93	\$24.56
------	---------------	---------	---------	---------	---------	---------

HAGER WORLDWIDE HAPPY MORNING XYLITOL

Single-Use Toothbrush with Xylitol Toothpaste in Bristles

- Toothbrushes pre-pasted with mint flavored Xylitol toothpaste
- Individually wrapped for ultimate hygiene and to help reduce risk of spreading infection
- Just wet bristles to activate the toothpaste



605496		50/box	\$21.49	\$20.42	\$20.42	\$19.34
--------	--	--------	---------	---------	---------	---------

TESS® PERIO BRUSHES

- Black bristles mask blood and reduce patient alarm
- 0.005" ultra-soft bristles for sensitive teeth and gums
- 38 Tuft



Handle Colors: Black, White, Purple, Orange, Pink, Red, Green, Yellow, Blue

3810	Assorted	72/box	\$53.79	\$51.10	\$51.10	\$48.41
------	----------	--------	---------	---------	---------	---------

TESS® ADULT BRUSHES

Accent Brushes

- Compact head is suited for smaller mouths and hard-to-reach areas
- 0.006" bristles surround 0.007" bristles for soft, comfortable cleaning
- Contoured clear handle features a rubber thumb rest for easy gripping and control
- 36 Tuft



Thumb Rest Colors: Blue, Lavender, Green, Pink, Red

3600	Assorted	72/box	\$68.79	\$65.35	\$65.35	\$61.91
------	----------	--------	---------	---------	---------	---------

Concept Curve Brushes

- Handle features an oval head, flexible neck and indented thumb rest for comfort and ease of access
- 0.007" bristles for gentle but effective cleaning
- 39 Tuft



Pearlescent Handle Colors: Mint, Sage, Blue, Lavender, Rose

3920	Assorted	72/box	\$48.49	\$46.07	\$46.07	\$43.64
------	----------	--------	---------	---------	---------	---------

TESS® CHILD BRUSHES

Junior Brush

- 6" long for ages 4-8
- Combination of 0.006" and 0.007" bristles for soft and effective cleaning
- Outer row bristles are extra soft for comfort along the gumline and between teeth
- 28 Tuft



Handle Colors: Purple, Orange, Pink, Red, Green, Yellow, Blue, White

2800	Assorted	72/box	\$45.29	\$43.03	\$43.03	\$40.76
------	----------	--------	---------	---------	---------	---------

SHERMAN SPECIALTY TAKE HOME BAG

- Super patient pattern



98ZX	Teeth 7"x10"	100/box	\$31.29	\$29.73	\$29.73	\$28.16
------	--------------	---------	---------	---------	---------	---------



infinity

Dental Implant Systems



infinity
OCTAGON BLT

infinity
OCTAGON BL

infinity
OCTAGON TL

infinity
INTERNAL HEX

infinity
TRI-CAM

infinity
EXTERNAL HEX



QUALITY, COMPATIBLE IMPLANTS. UNBEATABLE VALUE.

infinity Dental Implant Systems are built to respect the relationship between your practice and your referring dentists. Our system is backed by industry leading engineering, high quality manufacturing, and our superior customer service. Supported by a lifetime warranty, we are confident that infinity Dental Implant Systems can seamlessly fit into your practice.

To learn more, visit www.infinitydentalimplant.com or call to speak with one of our implant specialists at **800.441.3100**.



www.infinitydentalimplant.com



HANDPIECES AND BURS

IMPLANTMED PLUS SURGERY UNIT

Accurate. Powerful. Uncomplicated.

Your reliable partner in implantology and oral surgery

Implantmed Plus has been specially designed for use in dental implantology.

Implantmed Plus gives you the tools for oral surgery procedures in the fields of implantology and maxillofacial surgery with maximum precision.

Safe - New implant placing & stability determining feature

Simple - Intuitive & clear operation

Flexible - Upgrade anytime

Advantages

- Intuitive interface with glass color touch screen
- Ideal for group practices with 6 user settings
- Documentation function to USB drive
- Automatic torque control
- Thread cutter function
- New pump design simplifies irrigation tube loading
- 24 month Warranty

SAVE HUNDREDS* ON PACKAGE BUNDLES!
choosing is simple as

1 2 3

*compared to buying the base Package 1 and upgrading later

1 CHOOSE: EM-19 motor with/without the LC (Light Control)

Both the EM-19 and the EM-19 LC motor's feature:

- 10mm shorter and over 21g lighter than previous motors
- 200 - 40,000 RPM
- Built only for the Implantmed SI-1015 Surgery Unit
- Compatible with LED G and LG which generates light by motors rotation

The EM-19 LC difference:

- LC (Light Contacts) powers the new "L" series LED lights electrically, not by motors rotation, providing full brightness, even at low speeds.

2 CHOOSE: With or without the Osstell SI-SQ Module Implant Stability Quotient (ISQ)

- Check implants primary stability
- Monitor implants osseointegration
- Immediate, early or delayed placement?
- Document implant stability and ISQ
- Peace of mind
- ISQ = Implant Stability Quotient
- See page 104 for additional ISQ details



3 CHOOSE: Wired or wireless foot control option

Foot Control features:

- Wired foot pedal with multifunctional buttons and variable speed control
- Control the functions of the Implantmed Plus SI-1015, SI-915 white unit, Elcomed SA-310 and Piezomed SA-320 with only one foot control

Wireless only features:

- No cables – no mess, freedom of movement
- Surgical team safety
- Long battery life
- Familiar interactive programmability, control locations and pedal sensitivity



Implantmed Plus SI-1015 Packages

		1 EM-19 EM-19 LC	2 Osstell	3 Wired Wireless				
90000207	Package 1	✓ -	-	✓ -	\$4,240.49	\$4,240.49	\$3,816.44	\$3,816.44
90000208	Package 2	✓ -	✓	✓ -	\$7,879.49	\$7,879.49	\$7,091.54	\$7,091.54
90000209	Package 3	✓ -	-	- ✓	\$4,769.49	\$4,769.49	\$4,292.54	\$4,292.54
90000210	Package 4	✓ -	✓	- ✓	\$8,188.29	\$8,188.29	\$7,369.46	\$7,369.46
90000211	Package 5	- ✓	-	- ✓	\$4,847.29	\$4,847.29	\$4,362.56	\$4,362.56
90000212	Package 6	- ✓	✓	✓ -	\$8,265.49	\$8,265.49	\$7,438.94	\$7,438.94
90000213	Package 7	- ✓	-	- ✓	\$5,070.49	\$5,070.49	\$4,563.44	\$4,563.44
90000214	Package 8	- ✓	✓	- ✓	\$8,320.49	\$8,320.49	\$7,488.44	\$7,488.44

Replacement Motor and Cable for Implantmed Plus SI-1015

30185000	EM-19 - 5.9' (1.8m)	\$1,559.99	\$1,559.99	\$1,403.99	\$1,403.99
30281000	EM-19 LC - 5.9' (1.8m)	\$2,167.79	\$2,167.79	\$1,951.01	\$1,951.01
30281001	EM-19 LC - 11' (3.5m)	\$2,540.99	\$2,540.99	\$2,286.89	\$2,286.89

W&H Osstell SI-SQ module for Implantmed Plus

30210000	with USB Documentation	\$3,639.79	\$3,639.79	\$3,275.81	\$3,275.81
----------	------------------------	------------	------------	------------	------------

ELCOMED SURGERY UNIT

Safe, reliable and powerful

The W&H Elcomed has been designed specifically to meet the requirements of maxillofacial surgery. In addition to the expert knowledge of the surgeon, the drive unit is also a pivotal factor in the success of a surgical or implant procedure.

#1 in the surgical operator

Elcomed from W&H is a high-quality innovative solution for maxillofacial surgery, dental surgery and implantology.

Powerful, versatile motor

Shortest and lightest motor with 50,000 rpm and 7 Ncm, can be used with all instruments that have an e-type / ISO 3964 handpiece connection.

One display, many possibilities

Elcomed makes complicated processes as simple as possible. All necessary parameters can be set easily using just one operating stage. In addition to revolutions per minute and torque, the supply of irrigation can also easily be adjusted during use. Six programs for the most frequent recurring operational procedures can be saved at the first start-up.

Advantages


- All information in one user level
- Optional documentation via W&H USB flash drive
- Shortest and lightest motor with 50,000 rpm
- Up to 80 Ncm on the rotary instrument
- 6 programs with individually adjustable parameters
- Automatic thread cutter function
- Motor with cable is thermo washer disinfectable and sterilizable

Each Elcomed unit includes:

- Control Unit with wired S-N1 Foot Control
- Motor and 5.9' (1.8 m) Cable
- Irrigation Stand with Motor/Handpiece Support
- 3 Disposable Irrigation Tubing Sets and 5 Tubing Clips
- Mains Cable, 2 Locking Pins
- 1 Year Warranty



OPTIONAL UPGRADE



SA-310 Elcomed Surgery Unit with Motor and Foot Control

15933100	with USB Documentation	\$8,506.99	\$8,506.99	\$7,656.29	\$7,656.29
15933102	without USB Documentation	\$7,979.99	\$7,979.99	\$7,181.99	\$7,181.99

Replacement Motor and Cable for Elcomed SA-310

05513400	5.9' (1.8m)	\$2,389.29	\$2,389.29	\$2,150.36	\$2,150.36
06205800	11' (3.5m)	\$2,428.79	\$2,428.79	\$2,185.91	\$2,185.91

W&H CONTRA ANGLES AND HANDPEICES can be found on **pages 102-103**

GENUINE W&H CASES & ACCESSORIES can be found on **pages 104-105**

HANDPIECES AND BURS

CONTRA-ANGLE HANDPIECES



Durable, ergonomic and robust.

W&H contra-angle handpieces are specifically designed for fatigue-free and flexible operation in the areas of oral, maxillofacial surgery and implantology. The high quality stainless steel of each handpiece with special scratch-resistant coating makes these handpieces particularly robust.

Providing The Best Treatment

- The high standards required for the achievement of optimum surgical results are met by the W&H range of surgical contra-angles, whose design is the result of close attention to detail, high operational safety and stainless steel construction.
- Able to be fully dismantled for thorough cleaning, which allows W&H handpieces to meet the most stringent hygiene requirements (except WI-75 E/KM)
- Thermo washer disinfectable and sterilizable
- Wide range for all surgical applications
- Coupling system: ISO 3964 (DIN 13.940)
- All W&H surgical and L G handpieces work with all E-type motors, although the L series work specifically with the EM-19 LC motor
- 1 year warranty (2 year on L series), 6 months warranty on W&H factory repairs

L G Series Mini LED+ Handpieces

W&H L G (Light Generator) Mini LED+ handpieces provide daylight quality lighting where you need it most – directly at the surgical sight, producing optimal visibility and natural red tones. Light is delivered through an independent power supply with an integrated generator built into the handpiece. The handpieces RPM creates the light, so the LED light gets brighter at higher RPM speeds and there is no light at low speeds.

L Series Mini LED+ Handpieces

The new L series Mini LED+ handpieces also provide daylight quality lighting and natural red tones, while shining at full brightness, even at slower “implant-placement” speeds. The L series will only work with the new EM-19 motors and new Implantmed SI-1015, but they will require the EM-19 LC (Light Contacts) motor to power L series LED+ light.

WS-56 Surgical Contra-angle, 1:1

- Designed for difficult surgical procedures
- 3 models available: **Standard** (no light), **L G LED+**, and **L LED+** (for EM-19 LC only)
- Push-button coupling system, uses surgical burs and cutters with contra-angled latch-type shank (Ø 2.35 mm surgical length burs), single spray
- Maximum speed: 50,000 RPM (40,000 for WS-56 L G)
- Indications: maxillary and mandibular osteotomy, gernectomy, sequestrectomy



WS-75 Surgical Contra-angle, 20:1

- 3 models available: **Standard** (no light), **L G LED+**, and **L LED+** (for EM-19 LC only)
- Hexagon chucking system uses surgical burs and cutters with contra-angled latch-type shank (Ø 2.35 mm surgical length burs, shank hardness >50 HRC, >520 HV)
- Kirschner/Meyer internal cooling system with left or right side external irrigation green spray clips, single spray
- Maximum speed: 50,000 RPM (40,000 for WS-75 L G)
- Indications: implantation, osteotomy



WI-75 E/KM Implantology Contra-angle, 20:1

- Standard version only (LED not available)
- Hexagon chucking system uses surgical burs and cutters with contra-angled latch-type shank (Ø 2.35 mm surgical length burs, shank hardness >50 HRC, >520 HV)
- Kirschner/Meyer internal cooling system with left side only external irrigation black spray clips, single spray
- Maximum speed: 50,000 RPM
- Indications: implantation, osteotomy
- Cannot be dismantled



WS-92 Surgical Contra-angle, 1:2.7

- 3 Models available: **Standard** (no light), **L G LED+**, and **L LED+** (for EM-19 LC only)
- Push-button chucking system uses surgical burs and cutters with FG (Friction Grip) shank (Ø 1.6 mm surgical length burs)
- A triple spray guarantees sufficient cooling of the rotary instrument and the tooth
- Maximum speed: 50,000 RPM (40,000 for WS-92 L G)
- Indications: hemisection, extraction of wisdom teeth



WS-91 Surgical Contra-angle with 45° Head, 1:2.7

- 45° head offers better visibility of the treatment site and makes the access to the maxillary and mandibular molars much simpler
- 3 models available: **Standard** (no light), **L G LED+**, and **L LED+** (for EM-19 LC only)
- Push-button chucking system uses surgical burs and cutters with FG (Friction Grip) shank (Ø 1.6 mm surgical length burs)
- A triple spray guarantees sufficient cooling of the rotary instrument and the tooth
- Maximum speed: 50,000 RPM (40,000 for WS-91 L G)
- Indications: hemisection, extraction of wisdom teeth



	What motors will this Contra-angle work with? >	Works with EM-19	Works with EM-19 LC	Implantmed older motors	Elcomed motor				
30038000	WS-56, Standard - without light	Yes	Yes	Yes	Yes	\$1,503.49	\$1,503.49	\$1,353.14	\$1,353.14
30039000	WS-56 LG, LED+ Self-generating LED	Yes *	Yes *	Yes	Yes	\$1,927.49	\$1,927.49	\$1,734.74	\$1,734.74
30040000	WS-56 L, LED+ (EM-19*/EM-19 LC motors only)	Yes *	Yes	No	No	\$1,588.79	\$1,588.79	\$1,429.91	\$1,429.91
30065000	WS-75, Standard - without light	Yes	Yes	Yes	Yes	\$1,764.29	\$1,764.29	\$1,587.86	\$1,587.86
30032000	WS-75 LG, LED+ Self-generating LED	Yes *	Yes *	Yes	Yes	\$1,960.49	\$1,960.49	\$1,764.44	\$1,764.44
30033000	WS-75 L, LED+ (EM-19*/EM-19 LC motors only)	Yes *	Yes	No	No	\$1,617.79	\$1,617.79	\$1,456.01	\$1,456.01
10207550	WI-75 E/KM Contra-angle	Yes	Yes	Yes	Yes	\$1,430.79	\$1,430.79	\$1,287.71	\$1,287.71
30041000	WS-92, Standard - without light	Yes	Yes	Yes	Yes	\$1,601.29	\$1,601.29	\$1,441.16	\$1,441.16
30042000	WS-92 LG, LED+ Self-generating LED	Yes *	Yes *	Yes	Yes	\$1,927.49	\$1,927.49	\$1,734.74	\$1,734.74
30043000	WS-92 L, LED+ (EM-19*/EM-19 LC motors only)	Yes *	Yes	No	No	\$1,476.99	\$1,476.99	\$1,329.29	\$1,329.29
30047000	WS-91, Standard - without light	Yes	Yes	Yes	Yes	\$1,601.29	\$1,601.29	\$1,441.16	\$1,441.16
30048000	WS-91 LG, LED+ Self-generating LED	Yes *	Yes *	Yes	Yes	\$1,927.49	\$1,927.49	\$1,734.74	\$1,734.74
30049000	WS-91 L, LED+ (EM-19*/EM-19 LC motors only)	Yes *	Yes	No	No	\$1,588.79	\$1,588.79	\$1,429.91	\$1,429.91

* Will also work with the EM-19 and EM-19 LC motors as the "LG" LED is self-generating.

* Will work with the EM-19 motor – but the LED light will only work with the EM-19 LC motor.

HANDPIECES AND BURS

STRAIGHT AND ANGLED HANDPIECES



Durable, ergonomic and robust.

The surgical straight handpieces have been specially designed for oral and maxillo-facial surgery and implantology. The high-quality stainless steel with special coating makes them particularly robust.

Providing the Best Treatment

W&H has developed a series of straight and angled handpieces for oral and maxillofacial surgery and implantology.

- Fatigue-free working thanks to optimal ergonomics
- Durable and robust due to high-quality stainless steel and scratch-resistant coating
- Wide range for all surgical applications
- Able to dismantle for thorough cleaning
- Thermo washer disinfectable and sterilizable
- Coupling system: ISO 3964 (DIN 13.940)
- All W&H surgical and L G handpieces work with all E-type motors, although the S-11 L works specifically with the EM-19 LC motor
- 1 year warranty (2 years on L series)
- 6-months warranty on W&H factory repairs

S-9 Angled, 1:1

- Standard and L G LED+ models available
- Lever chuck uses surgical burs and cutters Ø 2.35 mm (also Stryker system) with maximum length = 45 mm, external spray
- Maximum speed: 50,000 RPM (40,000 for S-9 LG)
- Indications: apical resection, osteotomy, apical trepanning, bone modeling

S-10 Angled, 1:1

- Standard version only (LED not available)
- Lever chuck, surgical burs and cutters Ø2.35 mm with recommended length of 65-70mm, external spray
- Maximum speed: 50,000 RPM
- Indications: osteotomy, fenestration, osteosynthesis in distal areas



S-11 Straight, 1:1

- 3 Models available: **Standard** (no light), **L G LED+**, and **L LED+** (for EM-19 LC only)
- Lever chuck, surgical burs and cutters Ø 2.35 mm (also Stryker system) with Maximum length = 45 mm, external spray
- Maximum speed: 50,000 RPM (40,000 for S-11 LG)
- Indications: apical resection, osteotomy, apical trepanning, bone modeling



S-12 Angled, 1:2

- Standard version only (LED not available)
- Lever chuck, surgical burs and cutters Ø 2.35 mm with recommended length of 65-70 mm, external spray
- Maximum speed: 40,000 RPM
- Indications: osteotomy, fenestration, osteosynthesis in distal areas



S-16 Angled, 1:2

- Standard version only (LED not available)
- Lever chuck, surgical burs and cutters Ø 2.35 mm with recommended length of 45 mm, external spray
- Maximum speed: 40,000 RPM
- Indications: premolar and molar area

SZ-75 L Zygoma Angled, 20:1

- Standard version only (LED not available)
- Twist chuck, surgical burs and cutters Ø 2.35 mm, external spray
- Maximum speed: 50,000 RPM
- Indications: Implantation of Zygoma implants osteotomy

	What motors will this Contra-angle work with? >	Works with EM-19	Works with EM-19 LC	Implantmed older motors	Elcomed motor				
30053000	S-9 Handpiece	Yes	Yes	Yes	Yes	\$1,731.29	\$1,731.29	\$1,558.16	\$1,558.16
30054000	S-9 L G LED+ Handpiece, self-generating LED	Yes *	Yes *	Yes	Yes	\$2,206.99	\$2,206.99	\$1,986.29	\$1,986.29
30059000	S-10 Handpiece	Yes	Yes	No	No	\$1,713.49	\$1,713.49	\$1,542.14	\$1,542.14
30056000	S-11 Handpiece	Yes	Yes	Yes	Yes	\$1,339.99	\$1,339.99	\$1,205.99	\$1,205.99
30057000	S-11 L G LED+ Handpiece, self-generating LED	Yes *	Yes *	Yes	Yes	\$2,057.49	\$2,057.49	\$1,851.74	\$1,851.74
30058000	S-11 L LED+ (EM-19*/EM-19 LC motor only)	Yes *	Yes	No	No	\$1,589.29	\$1,589.29	\$1,430.36	\$1,430.36
30061000	S-12 Handpiece	Yes	Yes	Yes	Yes	\$1,862.29	\$1,862.29	\$1,676.06	\$1,676.06
30411000	S-16 Handpiece	Yes	Yes	Yes	Yes	\$1,765.99	\$1,765.99	\$1,589.39	\$1,589.39
30095000	SZ-75 L Zygoma Handpiece	Yes	Yes	Yes	Yes	\$2,266.49	\$2,266.49	\$2,039.84	\$2,039.84

* Will also work with the EM-19 and EM-19 LC motors as the "LG" LED is self-generating.

* Will work with the EM-19 motor – but the LED light will only work with the EM-19 LC motor.



WORK IN DAYLIGHT QUALITY

W&H Surgical Handpieces & Contra-Angles with Mini LED+

Excellent lighting conditions facilitate accurate diagnoses and perfect treatment results. The surgical handpieces with Mini LED+ produce LED light independently. Now available in two models:

L G Series Mini LED+ Handpieces

W&H L G (LED Generator) Mini LED+ Handpieces provide daylight quality lighting where you need it most, directly at the surgical sight, producing optimal visibility and natural red tones. Light is delivered through an independent power supply with an integrated generator built into the handpiece. The RPM speed of the handpiece creates the light, so the light gets brighter at higher speeds and there is no light at low speeds.

L Series Mini LED+ Handpieces

The new W&H L series Mini LED+ handpieces also provide daylight quality lighting and natural red tones, while shining at full brightness, even at slower "implant-placement" speeds. The L series will only work with the new EM-19 motors and new Implantmed SI-1015, but they will require the EM-19 LC (Light Contacts) motor to power L series LED+ light.

HANDPIECES AND BURS




OSSTELL ISQ

Measurement result

Implant Stability Quotient (ISQ) is an objective industry standard for measuring implant stability. It is based on Resonance Frequency Analysis (RFA). The result is presented as an ISQ value of 1-100. The higher the ISQ, the more stable the implant.

Thus, osseointegration can be observed during the treatment progress to diagnose the optimum moment to load the implant.



How does it work?

The technology is based on the tuning fork principle. A new, sterile and disposable SmartPeg is attached to the implant and made to vibrate, just like a tuning fork. The purpose is to find the resonance frequency i.e. the frequency with the strongest vibration. The higher the resonance frequency, the higher the ISQ value and the more stable the implant is. This measurement is done quickly; it takes only a few seconds and is non-invasive. The result is presented as an ISQ value between 1-100.

Why measure implant stability?

- Help determine when an implant is ready for loading and avoid premature loading.
- Manage patients with risk factors.
- Avoid unnecessarily long treatment times.
- Achieve more predictable outcomes.
- Identify situations in which it is best to unload a provisional and delay a final restoration.
- Optimize communication and increase trust with colleagues and patients.
- Improve case documentation and quality assurance.

W&H OSSTELL BEACON

Opt for even better treatment

The Osstell Beacon shows you within a matter of seconds when an implant is osseointegrated and can therefore be loaded. Because it operates wirelessly, it gives you greater freedom of movement in your dental practice. The data can be safely stored and compared on the free OsstellConnect web site and or used with the Implantmed Plus SI-1015 unit to show and document each patients ISQ reading.



Key Features:

- Monitoring of osseointegration and evaluation of implant stability
- Unique technology based on over 1,000 scientific published articles
- Optimum user-friendliness and intuitive ISQ measurements
- Results that are easy to interpret using a colour scheme
- Highly suited to a digital working environment
- Included with the Osstell Beacon measuring device:
- USB Key for Bluetooth communication with a computer and or Implantmed Plus SI-1015 unit
- 1 Osstell TestPeg
- 1 SmartPeg mount
- Power / charging cord
- Carrying / storage case
- Instruction manual

08022260	each	\$2,969.79	\$2,969.79	\$2,672.81	\$2,672.81
----------	------	------------	------------	------------	------------

W&H OSSTELL SI-SQ MODULE

Implantmed SI-1015 Attachment

The Osstell ISQ technology shows you within a matter of seconds when an implant is osseointegrated and can therefore be loaded. Implant stability can be measured precisely with the W&H Osstell SI-SQ module, available as an accessory for the Implantmed SI-1015 only.



- The measured Implant Stability Quotient is confirmed by over 1,100 published articles.
- Included with the Osstell SI-SQ Module:
- Sterilizable measuring probe with cable
- 1 Osstell TestPeg
- 1 SmartPeg mount

30210000	SI-1015 only	each	\$3,639.79	\$3,639.79	\$3,275.81	\$3,275.81
----------	--------------	------	------------	------------	------------	------------

This SmartPeg Mount is included with the Beacon and the SI-SQ Module, but is also available separately



W&H OSSTELL SMARTPEGS

Sterile, Single-use

To measure the ISQ of an implant, you first screw a SmartPeg into the placed implant. The SmartPeg, with its magnet on top, works like a small tuning fork. The magnet on the SmartPeg is 'hit' with magnetic pulses, from the probe, which makes the SmartPeg vibrate. Due to the stiffness in the interface between the implant surface and the bone the SmartPeg will vibrate accordingly.

Includes:

- 5 Sterile SmartPegs (select the implant system and the platform)

Visit www.osstell.com to determine which SmartPeg will work with your specific implant – if not listed below



Number	Description	infinity Implant System		Platform
08009830	SmartPeg Type 12	TRI-CAM	3.5	
08009840	SmartPeg Type 13		4.3	
08009840	SmartPeg Type 13		5.0	
08010020	SmartPeg Type 32	INTERNAL HEX	3.5	
08009980	SmartPeg Type 27		4.5	
08010120	SmartPeg Type 4	OCTAGON TL	RP WP	
08010240	SmartPeg Type 53		NP	
08010250	SmartPeg Type 54		RP	

SmartPeg


see chart for order number	5/pack	\$128.49	\$128.49	\$115.64	\$115.64
----------------------------	--------	----------	----------	----------	----------

Replacement SmartPeg Mount

07460300	SmartPeg Mount	each	\$30.99	\$30.99	\$27.89	\$27.89
----------	----------------	------	---------	---------	---------	---------



HANDPIECES AND BURS



DISPOSABLE IRRIGATION TUBING

Each irrigation tube set is packaged sterile for single use.

04363600	A	2.2 m length, Implantmed & Elcomed SA-310	6/box	\$88.99	\$88.99	\$80.09	\$80.09
04364100	B	3.8 m length, Implantmed & Elcomed SA-310	6/box	\$101.79	\$101.79	\$91.61	\$91.61
04365100	C	3.8 m length, Elcomed SA-200 Motor Unit only	6/box	\$109.79	\$109.79	\$98.81	\$98.81

Green Luer Lock Irrigation Tubes

00929300		7' 2" (2.2 m) Implantmed SI-915 & Elcomed SA-310	10/pkg	\$85.79	\$85.79	\$77.21	\$77.21
00929400		12' 4" (3.8 m) Elcomed Units w/3.5 m motors only	10/pkg	\$100.29	\$100.29	\$90.26	\$90.26

FOOT CONTROLS FOR IMPLANTMED SURGERY UNITS

Replacement Foot Control Units are available for most generations of the W&H Surgery Units. Please know your unit model information before calling to order the specific foot control for your W&H unit.

30264000	D	WiFi, S-NW with CAN-Bus, includes 1 dongle, SI-1015 only	each	\$1,389.79	\$1,389.79	\$1,250.81	\$1,250.81
30285000	E	Wired, S-N2, for Implantmed SI-1015 only	each	\$860.99	\$860.99	\$774.89	\$774.89



GENUINE W&H ACCESSORIES FOR IMPLANTMED & ELCOMED UNITS

04013500	E	Sterilization Cassette, 27.5 cm x 17.25 cm x 2.75 cm	each	\$254.29	\$241.58	\$241.58	\$228.86
06338400	F	Irrigation Spike with Roller Clamp Autoclavable, for SI-95, SI-915 and SA-310 only	each	\$91.99	\$91.99	\$82.79	\$82.79
04006800	G	Locking Pin	each	\$13.99	\$13.99	\$12.59	\$12.59
04013900	H	Pump Tube Complete, Autoclavable fits SI-95, SI-915 and SA-310 Only	each	\$44.79	\$44.79	\$40.31	\$40.31
04014000	I	Spare Pump Feed Tubes, 85 mm, 4 mm O.D. fits SI-95, SI-915 and SA-310 only	3/set	\$21.99	\$21.99	\$19.79	\$19.79
04019000	J	Tube Clamps, for Implantmed SI-915 and older series only	5/pkg	\$25.99	\$25.99	\$23.39	\$23.39
06290600		Tube Clamps, for Implantmed SI-1015 and Elcomed SA-310	5/pkg	\$24.79	\$24.79	\$22.31	\$22.31
06016400		Y-Connector with tubing	3/set	\$56.29	\$56.29	\$50.66	\$50.66
04014700	K	1.25 Amp 250V Fuse, for Implantmed SI-95 and SI-915	each	\$7.29	\$7.29	\$6.56	\$6.56
04005900	L	Irrigation Stand	each	\$42.49	\$42.49	\$38.24	\$38.24
02821400	M	Main unit power cord	each	\$41.49	\$41.49	\$37.34	\$37.34

Coolant Tubes

02610500	O	Internal Bur Coolant Tube	each	\$41.49	\$41.49	\$37.34	\$37.34
05369700		Movable, for old S-11 LED G, SI-11 LED G, metal connection	each	\$68.49	\$68.49	\$61.64	\$61.64
06879500		Movable, for new S-9, S-9LG, S-11 & S-11LG, green	each	\$61.79	\$58.70	\$58.70	\$55.61
06903200		Movable, for new S-10 & S-12 handpieces, green connection	each	\$61.99	\$61.99	\$55.79	\$55.79

GENUINE W&H IRRIGATION SPRAY CLIPS


04757100	P	Internal/External, WI-75 E/KM (older 20:1's), black	each	\$32.29	\$32.29	\$29.06	\$29.06
06946300	Q	Internal/External, all three WS-75 models, left, green	3/pkg	\$48.29	\$48.29	\$43.46	\$43.46
06948300	R	Internal/External, all three WS-75 models, right, green	3/pkg	\$48.29	\$48.29	\$43.46	\$43.46
06948400	S	External only, left, green	3/pkg	\$48.29	\$48.29	\$43.46	\$43.46
06949500	T	External only, right, green	3/pkg	\$48.29	\$48.29	\$43.46	\$43.46

The **BLACK** spray clips will not work on handpieces designed for green clips and **GREEN** spray clips will not work on handpieces designed for black spray clips.



HANDPIECE LUBRICANTS

MD400	U	W&H F1 Service Oil Spray (spray cap sold separately)	400ml can	\$84.29	\$80.08	\$80.08	\$75.86
02038200	V	W&H Spray Cap for F1 Service Oil, reusable, for e-type handpieces	each	\$12.29	\$11.68	\$11.68	\$11.06
02036100	V	W&H Spray Cap for F1 Service Oil, reusable, for WS-91 / WS-92	each	\$10.99	\$10.44	\$10.44	\$9.89
1795	W	Nou-Clean Spray, for Nouvag and E-coupling handpieces	500ml can	\$113.29	\$107.63	\$107.63	\$101.96
OILWNOZZLE	X	Osada Silicone Spray Oil, for E-coupling handpieces	6 oz can	\$51.29	\$48.73	\$48.73	\$46.16
1600064006	Y	Bien-Air Lubrifluid® Spray Lubricant (blue can)	500ml can	\$73.99	\$70.29	\$70.29	\$73.89
1600036006	Y	Bien-Air Spraynet® Cleaning Spray (red can)	500ml can	\$73.99	\$70.29	\$70.29	\$66.59
AHP09	Z	Aseptico AseptiSpray, cleaner/lubricant, withstands autoclaving	6 oz can	\$41.79	\$39.70	\$39.70	\$37.61



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

HANDPIECES AND BURS

HALL® MICROFREE® HIGH SPEED BATTERY OPERATED DRILL



When comparing Hall's MicroFree® versus competitors, you will feel the difference and see the savings. MicroFree is lighter and ergonomic, sterilization methods take less time and energy – making it less expensive to process than competitors.

Automated Cleaning/Washing

IPX 6, 8, & 9 rating. Protected against high-temperature and high-pressure water jets and continuous submersion in water over one meter deep.

Sterilization Methods

Steam sterilization is safe, eliminating the need for more expensive STERRAD battery sterilization.

One-Step Sterilization

MicroFree Handpieces and batteries can be processed together in one tray, saving time and energy.

Weight and Ergonomics

MicroFree Drill and Power Unit weighs 380.7 grams; 55.20 grams less than Stryker's F1® drill and was engineered for optimal balance, reducing muscle fatigue.

HALL® MicroFree® Drill

- MicroFree Lithium power unit (battery) sold separately from the handpiece.
- 1 year warranty on Handpiece and Charger, 6 month warranty on battery

PRO8000SB	Battery powered drill	\$5,499.29	\$5,499.29	\$4,949.36	\$4,949.36
L3500	Lithium power unit adapter	\$300.29	\$300.29	\$270.26	\$270.26
L3500SB	Lithium power unit	\$2,000.29	\$2,000.29	\$1,800.26	\$1,800.26
L3000	Lithium battery charger	\$4,499.29	\$4,499.29	\$4,049.36	\$4,049.36

Stryker F1 is a registered trademark of Stryker

HALL® MICROPOWER+ ORAL MAX ELITE HIGH SPEED ELECTRIC DRILL



When you pick up a handpiece, you want the highest degree of control and precision and you want to have the stability and comfort in your hand that gives you exacting control over your instrument. With the ConMed MicroPower+ System, you get all this and more.

HALL® MicroPower+ Oral Max Elite

- 1 year warranty on Handpiece and Drive System

PRO7000DE	Electric drill	\$5,899.29	\$5,899.29	\$5,309.36	\$5,309.36
D4000	drive system	\$6,000.29	\$6,000.29	\$5,400.26	\$5,400.26
D4000A	drive system with irrigation	\$6,249.99	\$6,249.99	\$5,624.99	\$5,624.99
MC5057	universal electric cord	\$789.29	\$789.29	\$710.36	\$710.36
5020053	electric footswitch	\$1,300.29	\$1,300.29	\$1,170.26	\$1,170.26
E9418	irrigation tubing set	\$328.79	\$328.79	\$295.91	\$295.91

ACE SURGI-DRILL I HALL-TYPE AIR DRILL SYSTEM

Proven workhorse in Oral Surgery

- Speed: 100,000 RPM
- Operating pressure: 90-110 PSI
- For both delicate and complex bone surgery
- Compact • Lightweight • Convenient
- Low Maintenance • Safe • Surgically efficient
- Accepts standard 2.35 mm handpiece surgical burs 44 mm or longer
- 90-days Warranty on HANDPIECES: covers manufacturing defects and workmanship
- 60-days Warranty on REPAIRS



ACE Surgi-Drill I

Burs and bur guards shown sold separately.

00320000	Pneumatic air drill	\$3,073.79	\$3,073.79	\$2,766.41	\$2,766.41
----------	---------------------	------------	------------	------------	------------

ACE Surgi-Drill I Air Hose

90-day Warranty on HOSES: covers manufacturing defects and workmanship

00320104	Surgi-Drill I 10'	\$501.29	\$501.29	\$451.16	\$451.16
----------	-------------------	----------	----------	----------	----------

ACE Handpiece Lubricant

00320105	Pneumatic lubricant	\$15.49	\$14.72	\$14.72	\$13.94
----------	---------------------	---------	---------	---------	---------

HALL® SURGAIRTOME TWO® HIGH SPEED PNEUMATIC DRILL

The Surgairtome Two® has nearly a half century of dependability and performance and is among the most reliable oral high speed drills on the market today.



- Lightweight and balanced for optimal control, feel and visibility
- Highly efficient at cutting bone with minimal gas consumption
- Proven reliability with Oral and Maxillofacial Surgeons

HALL® Surgairtome Two®

- 1 year warranty on Handpiece

5058001	Pneumatic air drill	\$4,899.29	\$4,899.29	\$4,409.36	\$4,409.36
5052010	10' pneumatic hose	\$649.99	\$649.99	\$584.99	\$584.99
1375028	Pneumatic footswitch	\$821.49	\$821.49	\$739.34	\$739.34

ACE Surgi-Drill II Air Hose

90-days Warranty on Hoses: covers manufacturing defects and workmanship

00320125	Surgi-Drill II 10'	\$501.29	\$501.29	\$451.16	\$451.16
----------	--------------------	----------	----------	----------	----------

HANDPIECES AND BURS

ACCESSORIES FOR ACE SURGI-DRILL I, II AND HALL® HIGH SPEED DRILLS



ACE Bur Guards - stainless steel, ACE bur guards fit all air driven handpieces

00320101	Medium, 1-3/4"	each	\$199.99	\$199.99	\$169.99	\$169.99
00320102	Long, 2-1/2"	each	\$251.49	\$251.49	\$213.77	\$213.77



Hall® Bur Guards

1375012	Medium	each	\$299.99	\$299.99	\$269.99	\$269.99
1375011	Long	each	\$299.99	\$299.99	\$269.99	\$269.99
1375023	X-long	each	\$366.29	\$366.29	\$329.66	\$329.66

ACE Bur Brush - stainless steel



00200106	5-1/4" (135 mm)	each	\$35.79	\$35.79	\$30.42	\$30.42
----------	-----------------	------	---------	---------	---------	---------



ACE Autoclave Case - stainless steel, fits Hall drill, hose and bur guard

00320103	10" L x 6.5" W x 2.5" H	each	\$338.49	\$338.49	\$304.64	\$304.64
----------	-------------------------	------	----------	----------	----------	----------

ACE HALL DRILL WALL MOUNTS



Handpiece Holder-Receptacle Complete Autoclavable

00320118	Complete receptacle	each	\$123.29	\$123.29	\$110.96	\$110.96
----------	---------------------	------	----------	----------	----------	----------

Replacement Holder-Receptacle Components Autoclavable

0034243	Receptacle, plastic	each	\$8.99	\$8.99	\$8.09	\$8.09
0030003	Mounting plate w/screws	each	\$38.29	\$38.29	\$32.55	\$32.55

Handpiece Holder - Metal

Autoclavable, use as-is or with plastic replacement receptacle (above)



00320120	Metal HP holder	each	\$81.29	\$81.29	\$69.10	\$69.10
----------	-----------------	------	---------	---------	---------	---------

Handpiece Holder - Delrin

Cold sterilize only, comes with 2 screws



HHS	Delrin HP holder	each	\$73.29	\$73.29	\$65.96	\$65.96
-----	------------------	------	---------	---------	---------	---------



PALISADES DENTAL IMPACT AIR 45® HANDPIECE

The IMPACT AIR 45®'s unique 45° angled head allows maximum access and visibility. It is available with or without fiberoptics.



- Cut teeth and bone
- Apicoectomy
- Osteoplasty
- Odontoplasty
- Root Resection
- Compressed air or nitrogen driven (400-500,000 RPM)
- Integral irrigation with water stream only
- Uses standard friction grip burs
- Fits standard 2, 4, or 5 hole fiberoptic dental hose couplings
- 1 year handpiece warranty on manufacturer defects; 3 months warranty on turbine

PDT403	4 hole	each	\$872.79	\$872.79	\$785.51	\$785.51
PDT401	5 hole w/fiberoptics	each	\$1,029.29	\$1,029.29	\$926.36	\$926.36

Impact Air 45® Turbine Replacement

3 month warranty



603	for push button hp's	each	\$272.29	\$272.29	\$245.06	\$245.06
-----	----------------------	------	----------	----------	----------	----------

Glenlube Lubricant for Impact Air 45®



61802LUBE	2 oz bottle	each	\$88.79	\$84.35	\$84.35	\$79.91
-----------	-------------	------	---------	---------	---------	---------

REJUVENATE YOUR INSTRUMENTS

One Visit, and They'll Come Back Sharper and Stronger.



- Sharpening & Re-tipping
- Air Drills & Bur Guards
- Calibration Services
- Electric Surgical Motors
- Electrosurgery Units
- Sterilizers
- and more

It's like a SPA DAY for your INSTRUMENTS!

See page 205

HANDPIECES AND BURS

ACE CUSTOM FOOT PEDAL FOR AIR DRILL CONNECTION

CUSTOM
YOU CHOOSE THE CONNECTIONS & ADAPTERS

Custom foot actuator for nitrogen or compressed air operated handpieces

- You choose the connections & adapters
- Can also be used as a back up high speed surgical unit

Can be used with:

- Impact Air 45 degree handpiece
- ACE Hall Type Air Drill
- Henry Schein 45 degree handpiece
- Adec
- Midland



Custom Foot Pedals - Note: specify wall outlet and handpiece type

0031717	w/gauge, pedal housing	each	\$561.99	\$561.99	\$505.79	\$505.79
0031720	no gauge or pedal housing	each	\$501.29	\$501.29	\$451.16	\$451.16

ALLIED HEALTHCARE NITROGEN REGULATOR

Nitrogen Adapter and Regulator Assembly

- Includes regulator, gauge and quick-disconnect for Dynacon outlet
- Can be converted to any outlet



11200002		each	\$651.99	\$651.99	\$586.79	\$586.79
----------	--	------	----------	----------	----------	----------

ACE SINGLE NITROGEN REGULATOR

with Female Schrader Outlet

- Regulator for application where flow is measured in psi (pounds per square inch) is required.
- Adjustable flow (0-200 psi)
- Requires H Size Tank



00320300		each	\$386.29	\$386.29	\$347.66	\$347.66
----------	--	------	----------	----------	----------	----------

COMPRESSED AIR/NITROGEN ADAPTERS



Bay Corporation Schrader Air Adapters

SCF164	Schrader Air Check with 1/4" NPT male adapter	each	\$121.99	\$115.89	\$115.89	\$109.79
SH1617S	Swivel quick-connect to 1/4" serrated hose stem	each	\$49.99	\$47.49	\$47.49	\$44.99

Allied Healthcare Schrader Air Adapters

12061110	Male quick-connect non swivel to 1/4 NPT male	each	\$125.29	\$119.03	\$119.03	\$112.76
----------	---	------	----------	----------	----------	----------

Imperial Eastman Adapter

DCI0015	1/4" male with 1/4" poly-flow connection	each	\$48.79	\$46.35	\$46.35	\$43.91
---------	--	------	---------	---------	---------	---------

MACAN RADIOSURGE ELECTROSURGERY UNIT



Unit Includes:

- Autoclavable handpiece and cord
- Indifferent plate and an 8' extension cord (003-3306)
- Six Ultraflex electrodes (see next page, **KIT** = included)
- "The Electrosurgery Handbook"
- Three dental waveforms (cutting, cut-coagulate, coagulate)
- Foot pedal
- Units are compact and lightweight
- 110V (220V available)
- Dimensions: 6.5" wide x 2.75" high x 11.5" deep
- 3-year warranty on RadioSurge unit (excludes handpiece, pedal, plate and electrodes)

MAAMC6A	RadioSurge unit	\$2,391.99	\$2,272.39	\$2,272.39	\$2,152.79
---------	-----------------	------------	------------	------------	------------

BIPOLAR JEWELER'S FORCEPS



- Insulated
- Fine Serrated Tip
- Sold with Cord
- Microfine, Straight, Smooth

BPFS1	4-1/2" (11.5 cm)	each	\$356.79	\$338.95	\$338.95	\$321.11
-------	------------------	------	----------	----------	----------	----------

BOVIE® HIGH TEMPERATURE CAUTERIES



Bovie manufactures more cauteries than anyone in the world. Therefore, Bovie has the ability to put the right cautery in the physician's hand regardless of the procedure. High temperature (2200° F) cauteries are popular for many types of procedures including, pinpoint hemostasis, controlling diffuse bleeders, and many others.

- Individually packaged, sterile, disposable
- Latex-free

AA25	Micro-fine tip	10/box	\$149.29	\$141.83	\$141.83	\$134.36
AA01	Fine tip	10/box	\$114.99	\$109.24	\$109.24	\$103.49
AA03	Loop tip	10/box	\$125.29	\$119.03	\$119.03	\$112.76
AA05	Loop tip, 2" Shaft	10/box	\$138.29	\$131.38	\$131.38	\$124.46

MEDLINE HIGH TEMP CAUTERY PENS



- Intended for coagulating tissue or arresting bleeding from small vessels by heat conducted through the wire tip
- Alkaline battery-powered cautery available with various tip styles
- Protective cap covers the switch and top
- Tip heats in approximately 1.0 second
- Tip glows a dull gray at low temperatures, glows a bright red at high temperatures.



ESCT001	Fine tip	each	\$12.99	\$12.34	\$12.34	\$11.69
ESCT003	Loop tip	each	\$11.99	\$11.39	\$11.39	\$10.79

HANDPIECES AND BURS

MACAN FULGURATION ADAPTOR



For use with RadioSurge MC6A and MC-4A low voltage units, produces high voltage needed for fulguration. The adaptor raises voltage for pin-point soft coagulation, treating enucleated cysts or fistulous tracts, for treating tumor excision sites, and for the desiccation of small pre-cancerous lesions. About 1mm arc length is available which suggests the use of an optical loupe during use. The adaptor limits therapeutic current in the event of accidental tissue contact. Use a maximum power setting of "6" and COAG when using the fulguration adaptor.

ME6000	Fulguration adaptor	each	\$153.49	\$145.82	\$145.82	\$138.14
--------	---------------------	------	----------	----------	----------	----------



MACAN RADIOSURGE ACCESSORIES

DPP1	Dispersive "ground" plate, metal, autoclavable, 4" x 4"	each	\$74.79	\$71.05	\$71.05	\$67.31
DPC2	8' extension cord – red	each	\$40.49	\$38.47	\$38.47	\$36.44
HPAC1	Handpiece and 8' extension cord – black	each	\$146.99	\$139.64	\$139.64	\$132.29
NCS1	Nose cap for handpiece	each	\$28.49	\$27.07	\$27.07	\$25.64

MACAN RIGID ELECTRODES

- Use with Macan RadioSurge (MAAMC6A)
- Not shown to scale

Wires	\$27.99	\$26.59	\$26.59	\$25.19
RC55	each			
Blades	\$29.29	\$27.83	\$27.83	\$26.36
R71	each			
Balls	\$27.99	\$26.59	\$26.59	\$25.19
RC51	each			
RC52	each			
Needles	\$27.99	\$26.59	\$26.59	\$25.19
RF10	each			
RF12	each			
RSN62	each			
RSN61	each			
Loops	\$27.99	\$26.59	\$26.59	\$25.19
RL32	each			
RL33	each			
RL34	each			
R45	each			
RL35	each			
RL36	each			
RL39	each			
RL30	each			
RL42	each			
RL22	each			
RL37	each			

MACAN ULTRA-FLEX BENDABLE ELECTRODES

- **KIT** = included with initial Unit
- Use with Macan RadioSurge (MAAMC6A)

Adjusto-Tip	\$36.99	\$35.14	\$36.99	\$36.99
--------------------	---------	---------	---------	---------

- Adjustable length of tip permits optimal control.
- Predetermine depth of cut in any surgical incision.
- Precise depth control is obtained by adjusting the cutting wire in or out of the isolated guide sheath and electrode shaft.
- Especially good for use anteriorly in thin tissue and for working around crowns, bridges and implants.
- Not shown to scale

UFF1	each	KIT		
------	------	------------	--	--

Fine Needles	\$32.29	\$30.68	\$30.68	\$29.06
---------------------	---------	---------	---------	---------

UFF10	each			
UFF12	each			
UFF11	each			
UFF13	each			

Wire	\$32.29	\$30.68	\$30.68	\$29.06
-------------	---------	---------	---------	---------

UF54	each			
Standard Needles	\$32.29	\$30.68	\$30.68	\$29.06
UFSN61	each			
UFSN62	each	KIT		
Balls	\$32.29	\$30.68	\$30.68	\$29.06
UFC51	each	KIT		

Loops	\$32.29	\$30.68	\$30.68	\$29.06
--------------	---------	---------	---------	---------

UFL22	each	KIT		
UFL32	each			
UFL33	each			
UFL34	each	KIT		
UFL45	each			
UFL35	each	KIT		
UFL36	each			
UFL38	each			
UFL31	each			
UFL37	each			

HANDPIECES AND BURS



Shanks above and drill images below have been enlarged to show detail. See descriptions for actual sizes.

SURGICAL BURS

6

CARBIDE		shank	head material	head diameter	bur number	total length	cutting length	features	qty	regular	member	premier	premier+
ROUND BURS	21441960	HP	Carbide	1.8 mm	6	44 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21501960	HP	Carbide	1.8 mm	6	50 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21651960	HP	Carbide	1.8 mm	6	65 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21441980	HP	Carbide	2.3 mm	8	44 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21441980S	SN	Carbide	2.3 mm	8	44 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21501980	HP	Carbide	2.3 mm	8	50 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
	21651980	HP	Carbide	2.3 mm	8	65 mm	-	Round	10	\$40.49	\$38.47	\$38.47	\$36.44
702 TAPER BURS	61441902	HP	Carbide	1.6 mm	702	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	61501902	HP	Carbide	1.6 mm	702	50 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	61651902	HP	Carbide	1.6 mm	702	65 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	62441902	HP	Carbide	1.6 mm	1702	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	62501902	HP	Carbide	1.6 mm	1702	50 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	62651902	HP	Carbide	1.6 mm	1702	65 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
1702 TAPER DOME BURS	62441902ST	ST	Carbide	1.6 mm	1702	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	62591902ST	ST	Carbide	1.6 mm	1702	59 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	61441909	HP	Carbide	2.1 mm	703	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	61501909	HP	Carbide	2.1 mm	703	50 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
703 TAPER BURS	61651909	HP	Carbide	2.1 mm	703	65 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	63441903	HP	Carbide	2.1 mm	1703	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	63501903	HP	Carbide	2.1 mm	1703	50 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
1703L TAPER DOME BURS	63651903	HP	Carbide	2.1 mm	1703	65 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	63441903ST	ST	Carbide	2.1 mm	1703	44 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	63591903ST	ST	Carbide	2.1 mm	1703	59 mm	5 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	64441903	HP	Carbide	2.1 mm	1703L	44 mm	7 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
1703L TAPER DOME BURS	64501903	HP	Carbide	2.1 mm	1703L	50 mm	7 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	64651903	HP	Carbide	2.1 mm	1703L	65 mm	7 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	64441903ST	ST	Carbide	2.1 mm	1703L	44 mm	7 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	64591903ST	ST	Carbide	2.1 mm	1703L	59 mm	7 mm	Tapered	10	\$40.49	\$38.47	\$38.47	\$36.44
	40500016	HP	Carbide	1.5 mm	-	50 mm	7.5 mm	-	10	\$109.79	\$104.30	\$104.30	\$98.81
	CYLINDER BUR	504080501	HP	Carbide	4 mm	-	48 mm	8 mm	Cross Cut	1	\$36.79	\$34.95	\$34.95
31680031		HP	Carbide	4 mm	-	70 mm	8 mm	Cross Cut	1	\$43.49	\$41.32	\$41.32	\$39.14
31440023		HP	Carbide	5.5 mm	-	44 mm	10 mm	Cross Cut	1	\$42.49	\$40.37	\$40.37	\$38.24
31500023		HP	Carbide	5.5 mm	-	50 mm	10 mm	Cross Cut	1	\$42.49	\$40.37	\$40.37	\$38.24
504080505		HP	Carbide	5.5 mm	-	70 mm	10 mm	Cross Cut	1	\$43.49	\$41.32	\$41.32	\$39.14
31440020S		SN	Carbide	4 mm	-	44 mm	8 mm	Cross Cut	1	\$38.29	\$36.38	\$36.38	\$34.46
504080501S		SN	Carbide	4 mm	-	48 mm	8 mm	Cross Cut	1	\$38.29	\$36.38	\$36.38	\$34.46
504502ST		ST	Carbide	4 mm	-	54 mm	8 mm	Cross Cut	1	\$38.29	\$36.38	\$36.38	\$34.46
31440023ST		ST	Carbide	5.5 mm	-	44 mm	10 mm	Cross Cut	1	\$44.79	\$42.55	\$42.55	\$40.31
31500023ST		ST	Carbide	5.5 mm	-	50 mm	10 mm	Cross Cut	1	\$44.79	\$42.55	\$42.55	\$40.31
31310023		CA	Carbide	4 mm	-	31 mm	8 mm	Cross Cut	1	\$48.79	\$46.35	\$46.35	\$43.91
31310024		CA	Carbide	5.5 mm	-	31 mm	10 mm	Cross Cut	1	\$48.79	\$46.35	\$46.35	\$43.91
SINUS LIFT EGG BURS		00350028	CA	Carbide	2.8 mm	-	35 mm	5 mm	Cross Cut	1	\$36.99	\$35.14	\$35.14
	00450028	HP	Carbide	2.8 mm	-	45 mm	5 mm	Cross Cut	1	\$36.99	\$35.14	\$35.14	\$33.29
	00450028ST	ST	Carbide	2.8 mm	-	45 mm	5 mm	Cross Cut	1	\$39.49	\$37.52	\$37.52	\$35.54
	00450028S	SN	Carbide	2.8 mm	-	45 mm	5 mm	Cross Cut	1	\$39.49	\$37.52	\$37.52	\$35.54

continued on following page



HANDPIECES AND BURS

													SPECIALTY
	shank	head material	head diameter	bur number	total length	cutting length	features	qty	regular	member	premier	premier+	
10440166	HP	SS	2 mm	-	44 mm	10 mm	-	1	\$30.49	\$28.97	\$28.97	\$27.44	LINDEMAN BURS 
10600168	HP	SS	2 mm	-	60 mm	24 mm	-	1	\$41.49	\$39.42	\$39.42	\$37.34	
10650168	HP	SS	2.3 mm	-	65 mm	24 mm	-	1	\$41.49	\$39.42	\$39.42	\$37.34	
10800169	HP	SS	2.3 mm	-	80 mm	35 mm	-	1	\$50.99	\$48.44	\$48.44	\$45.89	
10350166	CA	SS	2 mm	-	35 mm	10 mm	-	1	\$27.49	\$26.12	\$26.12	\$24.74	
10440166S	SN	SS	2 mm	-	44 mm	10 mm	-	1	\$33.29	\$31.63	\$31.63	\$29.96	
10650168S	SN	SS	2.3 mm	-	65 mm	24 mm	-	1	\$38.29	\$36.38	\$36.38	\$34.46	
10650168ST	ST	SS	2.3 mm	-	65 mm	24 mm	-	1	\$38.29	\$36.38	\$36.38	\$34.46	
19260162	FG	SS	1.6 mm	-	26 mm	9 mm	-	1	\$27.49	\$26.12	\$26.12	\$24.74	
88330009	CA	SS	1.6-2.35 mm	-	34 mm	15 mm	Internal Irr.	1	\$58.29	\$55.38	\$55.38	\$52.46	
88330010	CA	SS	3.15-4 mm	-	34 mm	15 mm	Internal Irr.	1	\$61.49	\$58.42	\$58.42	\$55.34	
00500013	HP	SS	6 mm	-	50 mm	12 mm	Contouring	1	\$36.99	\$35.14	\$35.14	\$33.29	PEAR BUR 
00500014	HP	SS	6 mm	-	50 mm	12 mm	Contouring	1	\$36.99	\$35.14	\$35.14	\$33.29	BARREL BUR 
00500015	HP	SS	4 mm	-	50 mm	14 mm	Contouring	1	\$36.99	\$35.14	\$35.14	\$33.29	CONE BURS 
45444501	HP	SS	0.81-1.6-2.1mm	-	48 mm	-	-	1	\$64.79	\$61.55	\$61.55	\$58.31	3-STEP PILOT DRILLS  ideal for use in lag screw placement
45444001	CA	SS	0.81-1.6-2.1mm	-	30 mm	-	-	1	\$64.79	\$61.55	\$61.55	\$58.31	
454281001	CA	SS	0.81-1.6 mm	-	30 mm	-	-	1	\$64.79	\$61.55	\$61.55	\$58.31	2-STEP PILOT DRILLS 
454282001	HP	SS	0.81-1.6 mm	-	47 mm	-	-	1	\$64.79	\$61.55	\$61.55	\$58.31	
45440201	CA	SS	1.1 mm	-	23.6 mm	-	-	1	\$75.29	\$71.53	\$71.53	\$67.76	PILOT TWIST DRILLS 
45440501	CA	SS	0.8 mm	-	22 mm	-	-	1	\$75.29	\$71.53	\$71.53	\$67.76	
11260357	CA	SS	3.75-7 mm	-	26 mm	7 mm	-	1	\$120.49	\$114.47	\$114.47	\$108.44	CRESTOTOME BUR 
00500035	HP	Titanium	3.5 mm	-	50 mm	14 mm	Contouring	1	\$36.99	\$35.14	\$35.14	\$33.29	TITANIUM CUTTING BURS to form Alveolar ridge 
00500060	HP	Titanium	6 mm	-	50 mm	14 mm	Contouring	1	\$36.99	\$35.14	\$35.14	\$33.29	
													DIAMOND
	shank	head material	head diameter	bur number	total length	cutting length	features	qty	regular	member	premier	premier+	
00500080	HP	Diamond	5 mm	-	44 mm	-	0.6 mm thick	1	\$36.99	\$35.14	\$35.14	\$33.29	fine to medium grit DIAMOND CIRCULAR SAWS 
00230080	CA	Diamond	5 mm	-	22 mm	-	0.6 mm thick	1	\$36.99	\$35.14	\$35.14	\$33.29	
00500080ST	ST	Diamond	5 mm	-	44 mm	-	0.6 mm thick	1	\$36.99	\$35.14	\$35.14	\$33.29	
00500080S	SN	Diamond	5 mm	-	44 mm	-	0.6 mm thick	1	\$36.99	\$35.14	\$35.14	\$33.29	
90250014	CA	Diamond	1.4 mm	4	22 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	ROUND DIAMOND BURS 
90250018	CA	Diamond	1.8 mm	6	22 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	
90250023	CA	Diamond	2.3 mm	8	22 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	
92250014	FG	Diamond	1.4 mm	4	25 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	
92250018	FG	Diamond	1.8 mm	6	25 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	
92250023	FG	Diamond	2.3 mm	8	25 mm	-	Round	10	\$107.29	\$101.93	\$101.93	\$96.56	
90440018	HP	Diamond	1.8 mm	6	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	ROUND DIAMOND BURS 
90440023	HP	Diamond	2.3 mm	8	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	
90500018	HP	Diamond	1.8 mm	6	50 mm	-	Round	1	\$51.99	\$49.39	\$49.39	\$46.79	
90500023	HP	Diamond	2.3 mm	8	50 mm	-	Round	1	\$51.99	\$49.39	\$49.39	\$46.79	
90500100	HP	Diamond	2.5 mm	-	50 mm	-	Round	1	\$51.99	\$49.39	\$49.39	\$46.79	
90600200	HP	Diamond	5 mm	-	60 mm	-	Round	1	\$60.49	\$57.47	\$57.47	\$54.44	
90700045	HP	Diamond	1.5 mm	-	70 mm	-	Round	1	\$60.49	\$57.47	\$57.47	\$54.44	
90700100	HP	Diamond	2.5 mm	-	70 mm	-	Round	1	\$60.49	\$57.47	\$57.47	\$54.44	
90440018ST	DT	Diamond	1.8 mm	6	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	
90440023ST	ST	Diamond	2.3 mm	8	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	
90440018S	SN	Diamond	1.8 mm	6	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	
90440023S	SN	Diamond	2.3 mm	8	44 mm	-	Round	1	\$35.99	\$34.19	\$34.19	\$32.39	
90500023S	DN	Diamond	2.3 mm	8	50 mm	-	Round	1	\$51.99	\$49.39	\$49.39	\$46.79	

6

Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

HANDPIECES AND BURS



ORAL SURGERY BURS

Stryker® and Hall® Compatible



Round	Shape	Head Ø	SHANK 1†	SHANK 1††	SHANK 2††	SHANK 3††	SHANK 4††	SHANK 5††	SHANK 6††	SLRA†	SLFG†
	1	1.0	14823								14001
	2	1.2	14826							14112	14002
	4	1.4	14832							14114	14004
	6	1.8	14838	34838	30020	30035	30050	30060	30000	14116	14006
	8	2.3	14844	34844	30021	30036	30051	30061	30001	14118	14008
	10	2.7	14846	34846							

Cross Cut Fissure

	557	1.0	14871							14157	14057
	558	1.2	14874							14158	14058
	560	1.6	14880				30052		30002		
	1557	1.0	14898							14159	14059
	1558	1.2	14901								14060
	700	1.0	14883								14070
	701	1.2	14886	34916	30022		30053	30062			14071
	701L	1.2	14910								
	702	1.6	14889	34889	30023	30037	30054	30063	30003	14172	14072
	702L	1.6	14890		30024		30055	30064	30004		
	1702	1.6	14916		30029		30048	30069		14174	14074
	1702L	1.6	14917		30030		30049	30070			
	703	2.1	14892	34892	30025	30038	30056	30065	30005	14173	14073
	703L	2.1	14893		30026	30039	30057	30066	30006		
	1703	2.1	14914	34914	30027	30040	30058	30067	30007	14175	14075
	1703L	2.1	14915	34915	30028		30059	30068	30008		

† Non-Sterile, †† Sterile

SHANK SELECTION Reference the matrix above for item numbers

SHANK 1	Non-sterile - 44.5mm, fits Hall* & Most Standard Handpieces for Medium Bur Guard	10/box	\$23.99	\$22.79	\$22.79	\$20.39
SHANK 1	Sterile - 44.5mm, fits Hall* & Most Standard Handpieces for Medium Bur Guard	10/box	\$39.29	\$37.33	\$37.33	\$35.36
SHANK 2	Sterile - 51mm, fits Hall* and MicroAire*** Handpieces with Medium Bur Guard	10/box	\$57.29	\$54.43	\$54.43	\$51.56
SHANK 3	Sterile - 65mm, fits Hall* and MicroAire*** Handpieces with Long Bur Guard	10/box	\$57.29	\$54.43	\$54.43	\$51.56
SHANK 4	Sterile - 44.5mm, fits Stryker** 100K & Impaction Handpieces	10/box	\$57.29	\$54.43	\$54.43	\$51.56
SHANK 5	Sterile - 59mm, fits Stryker** 100K & Impaction Handpieces	10/box	\$57.29	\$54.43	\$54.43	\$51.56
SHANK 6	Sterile - 44.5mm, fits Stryker** 50K & Osteotome Drill J-Notch	10/box	\$57.29	\$54.43	\$54.43	\$51.56
SLRA	Non-sterile - surgical length right angle - 1.023" diameter .0925"	5/pkg	\$19.49	\$18.52	\$18.52	\$17.54
SLFG	Non-sterile - surgical length friction grip - 1.004" diameter .0630"	5/pkg	\$19.49	\$18.52	\$18.52	\$17.54


*Hall® is a registered trademark of Zimmer, Inc. **Stryker® is a registered trademark of Stryker Instruments. ***MicroAire® is a registered trademark of MicroAire Surgical Instruments, Inc.

HANDPIECES AND BURS

SS WHITE GREAT WHITE® BURS

Tackle your toughest material cutting challenges and save time with the multi-functional performance of Great White® carbide burs. Great White burs cut all hard materials smoothly, making the treatment procedure an easier experience for your patient and reducing rotary bur inventory in your office.

- All SLFG - surgical length friction grip


15062	#2 (GW-2)	10/pkg	\$39.99	\$37.99	\$37.99	\$35.99	
18212	#245 (GW245)	10/pkg	\$39.99	\$37.99	\$37.99	\$35.99	
18218	#557 (GW557)	10/pkg	\$39.99	\$37.99	\$37.99	\$35.99	
13091	#557 (GW557)	100/pkg	\$395.79	\$376.00	\$376.00	\$356.21	
18220-5	#557SL (GWSL557) long	5/pkg	\$25.49	\$24.22	\$24.22	\$22.94	
18222	#558 (GW558)	10/pkg	\$39.99	\$37.99	\$37.99	\$35.99	
18221	#1557 (GW1557)	10/pkg	\$39.99	\$37.99	\$37.99	\$35.99	

ADVANCE PLUS® STERILE ORAL SURGERY BURS

- Designed and manufactured for optimum strength, performance and durability
- Packaged in individual gamma sterile pouches
- Extensive carbide offering
- Wide variety of surgical carbides for bone cutting, oral surgery and periodontal procedures




					
Us No.	8	702	703	703L	1703L
Figure No.	H1	H33	H33	H33L	H33LR
Head Size (1/10 mm)	023	016	021	021	021
Cutting Length (mm)	-	4.4	4.8	7.5	7.5
Shank #2	•	•	•	-	-
Shank #5	-	-	•	•	•



Shank #	Cross Reference		See chart above	Head	Length	Shank	Quantity	Pricing					
	SS WHITE	Brasseler						Regular Price	Member Price	Premier Price	Premier+ Price		
Shank 2	ADV8210S	30021	000114S0	H1.50.023	#8 Round	51mm	2	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV82100S	—	000114S1	H1.50.023	#8 Round	51mm	2	100/box	\$446.49	\$424.17	\$424.17	\$401.84	
Shank 2	ADV702210S	30023	000465S0	H33.50.016	#702	51mm*	2	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV7022100S	—	000465S1	H33.50.016	#702	51mm*	2	100/box	\$446.49	\$424.17	\$424.17	\$401.84	
Shank 2	ADV703210S	30025	000466S0	H33.50.021	#703	51mm*	2	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV7032100S	—	000466S1	H33.50.021	#703	51mm*	2	100/box	\$446.49	\$424.17	\$424.17	\$401.84	
Shank 5	ADV703L510S	30066	5023881S0	H33L.95.021	#703L	59mm*	5	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV703L5100S	—	5023881S1	H33L.95.021	#703L	59mm*	5	100/box	\$446.49	\$424.17	\$424.17	\$401.84	
Shank 5	ADV1703L510S	30068	5023877S0	H33LR.95.021	#1703L	59mm†	5	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV1703L5100S	—	5023877S1	H33LR.95.021	#1703L	59mm†	5	100/box	\$446.49	\$424.17	\$424.17	\$401.84	
Shank 5	ADV703510S	30065	5023875S0	H33.95.021	#703	59mm*	5	10/bag	\$46.99	\$44.64	\$44.64	\$42.29	
	ADV7035100S	—	5023875S1	H33.95.021	#703	59mm*	5	100/box	\$446.49	\$424.17	\$424.17	\$401.84	

HANDPIECES AND BURS



ACE TREPHINE BURS - CONTRA-ANGLE

- For harvesting bone in autogenous bone graft placement.
- May be used for retrieval of implants.
- Laser marked score lines.
- Manufactured from Titanium Alloy (Ti-6Al-4V ELI).
- Internally irrigated to minimize the risk of overheating the bone during drilling.
- Can be used with the Drill Guide Paralleling System to eliminate drill chatter during graft or implant harvesting.

		I.D. mm	O.D. mm	Length mm	Cutting length	Markings						
04948001	A	2	2.8	35	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948002	B	2.75	3.5	35	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948201	C	3.75	4.6	35	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948202	D	4.25	5.10	35	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948401	E	5.25	6.10	32	15.5 mm	8,10,13,15 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948501	F	6	6.9	32	15.7 mm	8,10,13,15 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948502	G	7	7.9	32	15.8 mm	8,10,13,15 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948701	H	8	8.9	32	14 mm	8,10,13 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04948702	I	10	10.8	32	14 mm	8,10,13 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	

ACE TREPHINE BURS – HANDPIECE

- For harvesting bone in autogenous bone graft placement.
- Laser marked score lines.
- Burs may be notched for the Stryker 100K Impaction Drill – Flat or Stepped – or J-Notched for an additional charge.
- Stryker notched burs may be used in straight Hall type handpieces.
- Manufactured from Titanium Alloy (Ti-6Al-4V ELI).



		I.D. mm	O.D. mm	Length mm	Cutting length	Markings						
04949201	J	3.75	4.6	52	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04949202	K	4.25	5.1	52	18.5 mm	8,10,13,15,18 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	
04949501	L	6	6.9	49	15.7 mm	8,10,13,15 mm	each	\$85.79	\$85.79	\$77.21	\$77.21	

SUGGESTED TREPHINES TO USE FOR RETRIEVE FAILED IMPLANTS For Implants up to.....

3.3 mm use 04948201

3.5–4.0 mm use 04948202


4.5–5.0 mm use 04948401

5.0–6.0 mm use 04948501

ACE TREPHINE BURS SURGICAL KIT

- Includes the 9 Contra Angle Trephines listed above
- Includes the 3 Handpiece Trephines listed above
- Includes an ACE 2" Square Titanium Bowl and Delrin Case

0494512	kit	\$873.79	\$873.79	\$786.41	\$786.41
---------	-----	----------	----------	----------	----------



ACE DR. SMILER'S TREPHINE BURS – POSTERIOR

- For harvesting autogenous bone in the posterior maxillary region.
- Cutting length is 6.75 mm, total length is 23 mm.
- Laser marked score line at 6 mm depth.
- Manufactured from Titanium Alloy (Ti-6Al-4V ELI).




		I.D. mm	O.D. mm									
04948601	M	6	6.9	each	\$85.79	\$85.79	\$77.21	\$77.21				
04948602	N	7	7.9	each	\$85.79	\$85.79	\$77.21	\$77.21				

ACE LATCH TO SQUARE DRIVER CONVERTER

- Secures and locks contra shaft for hand or ratchet drive

20155411	14 mm long	each	\$49.49	\$49.49	\$44.54	\$44.54
----------	------------	------	---------	---------	---------	---------



* = Flat-End Taper Cross-Cut Fissure
† = Round-End Taper Cross-Cut Fissure

HANDPIECES AND BURS

ACE EASYRETRIEVE TREPHINE SYSTEM

The ACE EasyRetrieve Trephine System makes using trephines more practical and predictable.

- The Drill Guide centering pin, when used after the included 2mm pilot drill, precisely centers the trephine while controlling and virtually eliminating the walking and chatter most commonly associated with trephine bur use.
- The EasyRetrieve Trephines comes apart after harvesting to allow easy access to the retrieved bone core.

ACE EasyRetrieve Basic Kit

0494510	Two Piece Trephine Kit w/Drill Guide	kit	\$907.99	\$907.99	\$817.19	\$817.19
---------	--------------------------------------	-----	----------	----------	----------	----------

Additional EasyRetrieve Positive Depth Stop Pins

04941001	Contra Angle Driver Shaft	each	\$42.49	\$42.49	\$38.24	\$38.24
04953251	2.0mm Twist Drill	each	\$72.79	\$72.79	\$65.51	\$65.51
20151001	Installation & Removal Tool, 1.25 mm Hex	each	\$37.79	\$37.79	\$32.12	\$32.12
20155411	CA Ratchet Extension/Hand Driver Tool	each	\$49.49	\$49.49	\$44.54	\$44.54
04941301	Two Piece Trephine Knurled Holder	each	\$56.29	\$56.29	\$50.66	\$50.66
04941202	Two Piece Bone Removal Tool	each	\$37.79	\$37.79	\$34.01	\$34.01
04941801	Bone Core Separator Tool, 4.5" length	each	\$91.49	\$91.49	\$82.34	\$82.34
9600449	Square Titanium Bowl	each	\$105.99	\$105.99	\$95.39	\$95.39
04941108	22mm Drill Guide (for 4.25 I.D. & smaller)	each	\$32.29	\$32.29	\$27.45	\$27.45
04941109	23mm Drill Guide (for 3.5 I.D. & larger)	each	\$32.29	\$32.29	\$27.45	\$27.45
20195013	Delrin Tray	each	\$270.49	\$270.49	\$243.44	\$243.44
see below	Choose your first EasyRetrieve Trephine	each	included	included	included	included
see below	Choose your second EasyRetrieve Trephine	each	included	included	included	included

Patent Pending. I.D. = Internal Diameter O.D. = Outer Diameter

ACE EASYRETRIEVE TREPHINES

2 EasyRetrive trephines included with purchase of the ACE EasyRetrieve Trephine Kit (above)

- Contra angle driver shaft (04941001) sold separately



	I.D. mm	O.D. mm	Markings				
04942001	A 2.0	2.8	8,10,13,15 mm	each	\$89.49	\$89.49	\$80.54 \$80.54
04942002	B 2.5	3.2	8,10,13,15 mm	each	\$89.49	\$89.49	\$80.54 \$80.54
04942003	C 2.75	3.5	8,10,13,15 mm	each	\$89.49	\$89.49	\$80.54 \$80.54
04942004	D 3.2	4.0	8,10,13,15 mm	each	\$89.49	\$89.49	\$80.54 \$80.54
04942201	E 4.25	5.0	8,10,13 mm	each	\$89.49	\$89.49	\$80.54 \$80.54
04942202	F 5.25	6.0	8,10,13 mm	each	\$89.49	\$89.49	\$80.54 \$80.54

I.D. = Internal Diameter O.D. = Outer Diameter

ACE EASYRETRIEVE POSITIVE DEPTH STOP KIT


- Inserts securely into driver shaft before trephine is attached, stops trephine when the Depth Stop hits bone.
- If there was a 2mm pilot hole drilled the trephine will advance 2mm deeper, please plan accordingly.

04941100	Depth Stop Kit, includes: 4, 6, 8 and 10mm	4/set	\$100.29	\$100.29	\$90.26	\$90.26
----------	--	-------	----------	----------	---------	---------

EasyRetrieve Trephine	Depth Drilling Guide				Trephines cut deeper when using pilot holes.
	with 2mm pilot		without pilot hole		
	A - D	E - F	A - D	E - F	Trephines cut deeper with shorter Depth Stops.
These guides are available as part of the kit only	4 mm	17 mm	15 mm	13 mm	
	6 mm	15 mm	13 mm	11 mm	
	8 mm	13 mm	11 mm	9 mm	
	10 mm	11 mm	9 mm	7 mm	

Additional EasyRetrieve Positive Depth Stop Pins	each	\$32.29	\$32.29	\$29.06	\$29.06
--	------	---------	---------	---------	---------

04941103	13 mm	8 mm	6 mm	4 mm	
04941104	15 mm	6 mm	4 mm	2 mm	
04941105	17 mm	4 mm	2 mm	0 mm	



The key component of the kit is the two-piece Removal Tool (04-9413-01 and 04-9412-02), which replaces the CA Driver Shaft.

Drill to desired depth using the Positive Depth Stop Pins, or drill to a Laser Marking. The Harvested Bone Core will usually stay intact making it ideal for grafting and clinical case evaluations.

Slowly thread the Bone Removal Tool through the Knurled Holder. Bone core is slowly and controllably pushed out of the trephine.

ACE EASYRETRIEVE CASE SUGGESTIONS

EasyRetrieve 2 mm Trephine: Specially Designed for Clinical Case Evaluation and Research Studies

- Used to retrieve quick and easy cores for clinical research and studies
- See EasyRetrieve Trephine A above

04942001	each	\$89.49	\$89.49	\$80.54	\$80.54
----------	------	---------	---------	---------	---------

EasyRetrieve Dr. Langer Sinus Lift Procedure: used for Implant Placement

- Use the EasyRetrieve Basic Kit with appropriate implant bone taps:

ACE 3.75 mm Bone Tap

04956001	each	\$73.49	\$73.49	\$62.47	\$62.47
----------	------	---------	---------	---------	---------

ACE 4.75 mm Bone Tap

04955001	each	\$73.49	\$73.49	\$62.47	\$62.47
----------	------	---------	---------	---------	---------

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

HANDPIECES AND BURS



ACE DRILL GUIDE PARALLELING SYSTEM

Kit Includes:

- 20973201 2.8mm Threaded Parallelizing Pin
- 20973002 2.0/2.8mm Straight Parallelizing Pins (4)
- 20977001 Parallelizing Drill Guide, 8mm
- 20977501 Parallelizing Drill Guide, 20mm
- 20977801 Parallelizing Drill Guide, 28mm
- 04953251 2mm Twist Drill
- 7651030 Autoclavable Tray



20973000	Complete Kit	each	\$526.29	\$526.29	\$473.66	\$473.66
----------	--------------	------	----------	----------	----------	----------

6



ACE DRILL, BUR AND TREPHINE TRAY

- Ideal autoclavable carry case for drill, tissue punches and trephine burs
- Holds 6 instruments up to 33 mm total length

09991006	3-1/2" (8.5cm)	each	\$95.49	\$95.49	\$85.94	\$85.94
----------	----------------	------	---------	---------	---------	---------

SS WHITE HUNTER SCRATCH BUR BRUSH



SW26040		each	\$11.29	\$11.29	\$10.16	\$10.16
---------	--	------	---------	---------	---------	---------

SS WHITE UNIVERSAL ALUMINUM BLOCK

- Autoclavable



16391		each	\$31.29	\$31.29	\$28.16	\$28.16
-------	--	------	---------	---------	---------	---------

INTEGRA™ STAINLESS MESH BUR AND DIAMOND HOLDER



6528	2-1/2" inside diameter	each	\$16.49	\$16.49	\$14.84	\$14.84
------	------------------------	------	---------	---------	---------	---------

HU-FRIEDY® BUR CUSHION

- Hollow, rectangular, silicone base with a stainless steel lid
- Designed to hold burs securely in place during cleaning and sterilization
- Hollow design allows drainage and proper cleaning
- For use in ultrasonic cleaners



IMS1372T	Tall: holds 12 burs 3.3" x 1.22" x 0.9"	each	\$34.49	\$34.49	\$29.32	\$29.32
IMS1372TH	Tall-half: holds 6 burs 3.3" x 2.44" x 0.9"	each	\$26.49	\$26.49	\$22.52	\$22.52

ENCOUNTER BOOKS

- Step by Step Consultation
- Quality For Your Office
- PC Ready Forms for Your Practice- Including Exam Forms, IV Sedation Forms, and Consent Forms and so much more...
- Proper Procedure Coding



9910002	Dental Implants	50/pkg	\$179.99	\$170.99	\$170.99	\$161.99
9910001	Wisdom Teeth	100/pkg	\$179.99	\$170.99	\$170.99	\$161.99

A-TITAN EASY X-TRAC SYSTEM

Includes:

- (1) Sterilization Cassette
- (1) Easy X-Tractor
- (1) Set of 3 Drills
- (3) 28mm Extraction Screws
- (3) 33mm Extraction Screws
- (1) Ratchet
- (1) Protector Plate-Left
- (1) Protector Plate-Right
- (1) Bur Block
- (1) Instructional CD-Rom



K8000	kit	\$1,847.29	\$1,847.29	\$1,662.56	\$1,662.56
-------	-----	------------	------------	------------	------------

Also Sold Individually



XT120	X-TRAC Ratchet	3/box	\$215.79	\$215.79	\$194.21	\$194.21
-------	----------------	-------	----------	----------	----------	----------



A-Titan X-TRAC Drills

- set of 3 drills: 1.1mm, 1.3mm, 1.5mm

XT110	Includes 3 drills	kit	\$84.49	\$84.49	\$76.04	\$76.04
-------	-------------------	-----	---------	---------	---------	---------

A-Titan X-TRAC Screws

XT128	28mm length	each	\$84.49	\$84.49	\$76.04	\$76.04
XT130	30mm length	each	\$84.49	\$84.49	\$76.04	\$76.04
XT133	33mm length	each	\$84.49	\$84.49	\$76.04	\$76.04
XT135	35mm length, 2mm diameter for enlarged root canals	each	\$84.49	\$84.49	\$76.04	\$76.04

A-TITAN IMPLANTATION TREPHINE DRILL KIT

Designed to assist in the placement of CAMLOG screwline implants

- Implantation trephine drills automatically create a hole in the bone with the precise bottom (apical) diameter for implant placement
- The special 'form drills' build the side walls up to diameter of the implant
- All drills have depth markings of 9.0mm, 11.0mm and 13.0mm and are 14mm long



K363	Includes the 6 items below	kit	\$498.99	\$498.99	\$449.09	\$449.09
------	----------------------------	-----	----------	----------	----------	----------

Also Sold Individually

36302	2.4mm for 3.3mm implant	each	\$86.29	\$86.29	\$77.66	\$77.66
36303	3.3mm for 3.8mm implant	each	\$86.29	\$86.29	\$77.66	\$77.66
36304	3.6mm for 4.3mm implant	each	\$86.29	\$86.29	\$77.66	\$77.66
36305	4.3mm for 5.0mm implant	each	\$86.29	\$86.29	\$77.66	\$77.66
36306	5.2mm for 6.0mm implant	each	\$86.29	\$86.29	\$77.66	\$77.66

Sterilization Block (not sold individually)

RIGHT ANGLE DRIVER KIT

The manually driven right-angled driver makes it easy to place posterior screws in the mandible or maxilla.

- Kit Includes: RAD, ATA and BITSKIT (listed below)

RADKIT	kit	\$1,099.99	\$1,099.99	\$989.99	\$989.99
--------	-----	------------	------------	----------	----------

Also Sold Individually

Right Angle Driver

- Maximize access and comfort with the ability to rotate the head to the exact angle needed
- Pre-set torque adapters allow you to achieve a specific torque value
- Torque adapter is sold separately (below)



RAD	each	\$421.99	\$421.99	\$379.79	\$379.79
-----	------	----------	----------	----------	----------

Torque Adapter

- 7 torque values: 10, 15, 20, 25, 30, 32, 35 Ncm
- Easily connects to Right Angled Driver (RAD)
- Autoclavable



ATA	each	\$539.29	\$539.29	\$485.36	\$485.36
-----	------	----------	----------	----------	----------

Bits Kit

- Compatible with all major implant systems (RAD or any standard latch-type contra-angle system)



BITSKIT	set of 7	\$335.29	\$335.29	\$301.76	\$301.76
---------	----------	----------	----------	----------	----------

IMPRESSION MATERIAL



BOSWORTH DENTUSIL™

Silicone Soft Reline Material Kit

Cartridge dispensed silicone soft reline material provides denture patients with a comfortable, reliable reline in just one office visit. Flexible and resilient, DentuSil's silicone base absorbs less water so it's more resistant to stains and odors. Translucent pink color makes for a more natural, healthy appearance.

- Adhesive holds material in place, securing the reline in just minutes
- Contains 10 applications
- DentuSil™ lasts in the mouth for approximately 1 year
- Dispenser not included in this kit

Kit Includes: 50 ml Cartridge, Mixing Tips and 0.25 oz Liquid

0921276	kit	\$171.99	\$163.39	\$163.39	\$154.79
---------	-----	----------	----------	----------	----------



GC AMERICA COE-COMFORT™

- Extremely soft temporary denture lining for use beneath ill-fitting dentures
- Permits the patient to wear dentures while irritated mouth tissues return to a healthy form and condition for new impressions

Kit Includes: 6 oz Powder, 6 oz Liquid

341001	kit	\$190.99	\$181.44	\$181.44	\$171.89
--------	-----	----------	----------	----------	----------



VOCO QUICK UP

- Composite for bonding in attachments such as ball, Locator® and telescopic anchors in composite-based dentures.
- Reattachment of secondary elements such as bar retainers in acrylic-based dentures.
- Very simple chair-side application
- High strength values
- Complete set for immediate use including control and blocking out material for primary structures that are to be protected, such as implants
- Gingiva-colored

Kit Includes: 7.5g QuickMix Syringe, 2g Syringe Quick Up LC, 4ml Adhesive, 5ml QuickMix Syringe Fit Test C & B

1625	kit	\$199.99	\$189.99	\$189.99	\$179.99
------	-----	----------	----------	----------	----------

Refill: 7.5g Quick Up QuickMix Syringe with 10 Mixing Tips

VOC1626	each	\$64.29	\$61.08	\$61.08	\$57.86
---------	------	---------	---------	---------	---------



GC AMERICA COE ALGINATE™

- ADA Certified
- Eliminates can fluffing
- Low flow minimizes gagging
- High compressive and tear strength
- Sufficient powder for 24 full arch impressions

Fast Set – 1-3/4 Minutes, 1lb Canister

GCA120101	Mint	each	\$24.79	\$23.55	\$23.55	\$22.31
121101	Wild cherry	each	\$24.79	\$23.55	\$23.55	\$22.31



SEPTODONT PLASTALGIN

Dust-free alginate for accurate impression taking. Impression for case study and orthodontic models, fabrication of models for mouthguards, provisional restorations, bleaching trays, partial denture frameworks and occlusal splint appliances.

Type II Normal Setting, 1 lb Bag

01C1100	Apple	each	\$13.79	\$13.10	\$13.10	\$12.41
---------	-------	------	---------	---------	---------	---------

Ortho Type I Quick Setting, 1 lb Bag

01C1000	Apricot	each	\$13.79	\$13.10	\$13.10	\$12.41
---------	---------	------	---------	---------	---------	---------

IMPRESSION MATERIAL

ZHERMACK HYDROGUM® 5

Extra fast setting dust-free high stability alginate. Setting time: 1'50".

Clinical Applications

- Temporary crowns & bridges
- Study models
- Removable prosthesis and chomes
- Antagonists in fixed and removable prosthesis

Features

- High dimensional stability alginate
- Very long impression storage - 5 days stability (120 hours)
- Homogenous mass, smooth, compact surfaces
- 5-year shelf life
- Fast working and setting times
- Thixotropic consistency
- Lilac color, mangustan tropical fruit flavor

Intro Kit Includes: 453g (1 lb) Bag, 1 Storage Canister, 1 Measuring Set

C302072	kit	\$22.99	\$21.84	\$21.84	\$20.69
---------	-----	---------	---------	---------	---------

Refill: 453g (1 lb) Single Bag

C302070	each	\$17.29	\$16.43	\$16.43	\$15.56
---------	------	---------	---------	---------	---------



ZHERMACK® TROPICALGIN

Fast setting chromatic dust-free alginate. Setting time: 2'35".

Clinical Applications

- Removable prosthesis and chomes
- Antagonists in fixed and removable prosthesis
- Study models

Features

- Fast working and setting times
- Mango flavor
- Easy mixing
- Long impression storage
- Three-phase chromatic alginate: Red phase: mixing time, Orange phase: working time, Yellow phase: time in mouth
- Thixotropic
- Fast water absorption
- Homogenous mass, smooth, compact surfaces
- Packaged in modified atmosphere

Kit Includes: 2 x 453g (1 lb) Bags, 1 Storage Canister, 1 Measuring Set

C302242	kit	\$36.79	\$34.95	\$34.95	\$33.11
---------	-----	---------	---------	---------	---------

Refil: 453g (1 lb) Single Bag

C302240	each	\$14.49	\$13.77	\$13.77	\$13.04
---------	------	---------	---------	---------	---------



ZHERMACK OCCLUFAS®T ROCK

A-Silicone for bite registration: Maximum intercuspitation bite registration, centric, protrusion or lateral registrations, gnathological evaluations, diagnostic evaluation of craniocervical mandibular disorders, and orthodontic centric occlusion

Features

- Imperceptible consistency
- Fast 60 seconds: minimal time in the mouth ensures greater patient comfort and reduces possible errors due to minor jaw movements
- Extremely rigid after setting (95 Shore A)

Advantages

- Easy to use and no air bubbles
- Minimum thickness, virtually equal to zero
- Control of positioning, does not drip
- Easy to trim or mill to eliminate excess material

Includes: 2 - 50ml Cartridges (base + catalyst) + 12 Green Mixing Tips

C200726	kit	\$47.79	\$45.40	\$45.40	\$43.01
---------	-----	---------	---------	---------	---------



IMPRESSION MATERIAL



PARKELL BLU-MOUSSE®

Precision bite registration is only the beginning

- Exceptional accuracy
- Dependable, no-slump consistency
- From conventional bite registration techniques to modern-day digital impressions
- Plaster-like hardness (85-durometer)
- Choice of setting speeds
- Guns sold separately, not available from ACE SOUTHERN

Kits Include: 2 - 50ml Cartridges and 10 Mixing Tips

S438S	Scent-free classic, in split cartridge	kit	\$63.29	\$60.13	\$60.13	\$56.96
S438	Scent-free classic, in standard cartridge	kit	\$63.29	\$60.13	\$60.13	\$56.96
S457S	Scent-free super-fast, in split cartridge	kit	\$63.29	\$60.13	\$60.13	\$56.96
S457	Scent-free super-fast, in standard cartridge	kit	\$63.29	\$60.13	\$60.13	\$56.96

7

ZHERMACK ELITE® HD+

A- Silicone Impression Material

- Hydrophilic: Immediate hydrophilic behavior after just 2 seconds.
- Flowability: Nanoparticles in the fillers improve flowability in the gingival sulcus.
- Detailed: The synergistic action between immediate hydrophilia and nano-technology ensures a consistent reproduction of almost any kind of detail.
- Versatility: A complete range of viscosities and setting times suitable for any impression technique.
- Practical: The 5:1 system for automatic mixing machines fits straight into the machine. It does not need a cartridge sleeve.
- Complete: A range of products that satisfies all demands.



Fast Setting, Light Body - Includes: 2 - 50ml Cartridges (Base + Catalyst) + 12 Yellow Mixing Tips

C203040		kit	\$42.29	\$40.18	\$40.18	\$38.06
---------	--	-----	---------	---------	---------	---------

Normal Setting, Light Body - Includes: 2 - 50ml Cartridges (Base + Catalyst) + 12 Yellow Mixing Tips

C203030		kit	\$42.29	\$40.18	\$40.18	\$38.06
---------	--	-----	---------	---------	---------	---------

Tray Material, Heavy Body - Includes: 2 - 50ml Cartridges (Base + Catalyst) + 6 Green Mixing Tips

C202032		kit	\$42.29	\$40.18	\$40.18	\$38.06
---------	--	-----	---------	---------	---------	---------



BAUSCH ARTICULATING PAPER

- No-smudge
- Blue

BK09	Extra-thin, 0.0016" (0.04 mm) blue	200/book	\$20.49	\$19.47	\$19.47	\$18.44
------	------------------------------------	----------	---------	---------	---------	---------



GC AMERICA MIXING PADS

159112	COE parchment, mini, 6" x 6"	100 sheets/pad	\$12.29	\$11.68	\$11.68	\$11.06
--------	------------------------------	----------------	---------	---------	---------	---------



J&J INSTRUMENTS PLASTER SPATULAS

07863	3R 3.75" Stiff blade	each	\$8.79	\$8.79	\$7.91	\$7.91
07870	10R 4.25" Extra stiff blade	each	\$8.79	\$8.79	\$7.91	\$7.91

IMPRESSION MATERIAL

BUFFALO DENTAL SPATULAS FOR PLASTER AND STONE

- Rosewood handle

79390	8R Spatula, 4" blade, stiff, stainless steel	each	\$16.49	\$16.49	\$14.84	\$14.84
55570	6R Knife, 1-5/8" blade, plaster	each	\$16.49	\$16.49	\$14.84	\$14.84



PALMERO HEALTHCARE ALGINATE SPATULA

A white plastic, broad blade spatula with medium flex. Autoclavable up to 270 degrees.

1519	7-1/2" Length	each	\$4.79	\$4.79	\$4.31	\$4.31
------	---------------	------	--------	--------	--------	--------



PALMERO HEALTHCARE CANDEEZ FLEXIBLE MIXING BOWLS

Manufactured of a soft yet strong PVC material that is extremely pliable and allows for easier and more efficient mixing of plaster, stone, investments or alginates. The ultra smooth surface of the bowl makes for a bubble and void free mix that eases mixing effort and eliminates waste. The "self-cleaning" bowl cuts clean-up time in half with a unique flexible design that allows the bowl to be folded inside out and easily wiped clean. All of these bowls may also be used with most automatic alginate mixers and are a great addition to any dental operator or laboratory.

- Available in six vibrant colors; red, blue, yellow, purple, green, pink bubble gum scent (medium only) and clear

350 cc, Medium, 4"W x 2-1/2"D

1530R	Red	each	\$7.79	\$7.40	\$7.40	\$7.01
1530B	Blue	each	\$7.79	\$7.40	\$7.40	\$7.01
1530Y	Yellow	each	\$7.79	\$7.40	\$7.40	\$7.01
1530P	Purple	each	\$7.79	\$7.40	\$7.40	\$7.01
1530G	Green	each	\$7.79	\$7.40	\$7.40	\$7.01
1530C	Clear	each	\$7.79	\$7.40	\$7.40	\$7.01
1530BP	Pink w/bubblegum scent	each	\$7.79	\$7.40	\$7.40	\$7.01



600 cc, Large, 4-1/2"W x 3"D

1531R	Red	each	\$9.49	\$9.02	\$9.02	\$8.54
1531B	Blue	each	\$9.49	\$9.02	\$9.02	\$8.54
1531Y	Yellow	each	\$9.49	\$9.02	\$9.02	\$8.54
1531P	Purple	each	\$9.49	\$9.02	\$9.02	\$8.54
1531G	Green	each	\$9.49	\$9.02	\$9.02	\$8.54
1531C	Clear	each	\$9.49	\$9.02	\$9.02	\$8.54

ZHERMACK® RUBBER MIXING BOWL FOR ALGINATES

350 cc, Medium

U113280	each	\$8.49	\$8.07	\$8.07	\$7.64
---------	------	--------	--------	--------	--------



COLTÈNE WHALEDENT HYGENIC FLEXIBOLE

600 cc, Large, 4-1/4" Dia x 3-1/2" H

H00565	each	\$29.79	\$28.30	\$28.30	\$26.81
--------	------	---------	---------	---------	---------



IMPRESSION MATERIAL



ZHERMACK HYDROGUM® 5 MEASURING SET

Includes: Powder Scoop, Water Measure

D105187	each	\$3.29	\$3.13	\$3.13	\$2.96
---------	------	---------------	---------------	---------------	---------------



ZHERMACK® DISPENSING GUN AND TIPS

Intraoral Tips

C202090	Yellow	48/box	\$19.99	\$18.99	\$18.99	\$17.99
---------	--------	--------	----------------	----------------	----------------	----------------

Mixing Tips

C202070	Small, yellow	48/box	\$56.29	\$53.48	\$53.48	\$50.66
C202080	Standard, green	48/box	\$56.29	\$53.48	\$53.48	\$50.66
C205500	Medium, pink	48/box	\$56.29	\$53.48	\$53.48	\$50.66

1:1 Universal 50ml Dispensing Gun


C202100	each	\$111.29	\$105.73	\$105.73	\$100.16
---------	------	-----------------	-----------------	-----------------	-----------------



GC AMERICA CARTRIDGE DISPENSER

Used to extrude impression materials. Designed to work with either 48-50 mL or 75 mL delivery systems. Will not work with Blu-Mousse Classic in standard cartridge.

152101	each	\$157.49	\$149.62	\$149.62	\$141.74
--------	------	-----------------	-----------------	-----------------	-----------------



VISTA DENTAL IMPRESSION TRAYS - CLEAR

Ideal for all types of impression material

- Clear tray made of rigid acrylic providing superior visibility and added rigidity for greater accuracy
- Key-way perforations and side wall slots ensure material retention
- Posterior rimlock feature to minimize secondary flow or gagging
- Occlusal surface space maintainer
- Finger rest ledges for easier patient removal

VD310001	#1, Large upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310002	#2, Large lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310003	#3, Medium upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310004	#4, Medium lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310005	#5, Small upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310006	#6, Small lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310007	#7, Upper right, lower left	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310008	#8, Upper left, lower right	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
310009	#9, Anterior	12/bag	\$11.29	\$10.73	\$10.73	\$10.16


IMPRESSION MATERIAL

VISTA DENTAL IMPRESSION TRAYS - GREEN

Ideal for all types of impression material

- Key-way perforations and side wall slots ensure material retention
- Posterior rimlock feature to minimize secondary flow or gagging
- Occlusal surface space maintainer
- Finger rest ledges for easier patient removal

311001	#1, Large upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311002	#2, Large lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311003	#3, Medium upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311004	#4, Medium lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311005	#5, Small upper	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311006	#6, Small lower	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311007	#7, Upper right, lower left	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311008	#8, Upper left, lower right	12/bag	\$11.29	\$10.73	\$10.73	\$10.16
311009	#9, Anterior	12/bag	\$11.29	\$10.73	\$10.73	\$10.16



BOSWORTH BOTRAYS®

An economical choice for consistent, unsurpassed impressions. Alginate “grippers” and slotted walls keep material in the tray to reduce gag effect. Made from rigid, high-impact FDA approved styrene and organic FDA approved dye. Green in color.

0921364	#1, Perforated large upper	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921365	#2, Perforated large lower	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921366	#3, Perforated medium upper	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921367	#4, Perforated medium lower	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921368	#5, Perforated small upper	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921369	#6, Perforated small lower	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921370	#7, Perforated partial (upper right, lower left)	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921371	#8, Perforated partial (upper left, lower right)	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921372	#9, Perforated lower anterior	12/bag	\$12.29	\$11.68	\$11.68	\$11.06
0921373	#10, Perforated upper anterior	12/bag	\$12.29	\$11.68	\$11.68	\$11.06



GC AMERICA COE NICKEL-PLATED REGULAR PERFORATED TRAYS

- Industry standard
- Long lasting
- Ideal for alginate impressions
- Available in upper, lower and quadrant

Set of 8 Perforated Trays Includes: #1, 3, 4, 5, 7, 20, 21, 22

260018	set	\$229.29	\$217.83	\$217.83	\$206.36
--------	-----	----------	----------	----------	----------

Sizes are also sold individually. Call customer service or visit our website for more information.



GC AMERICA COE TRAY CLEANER

- Cleans: Alginate, Zinc Oxide, Eugenol, Modeling Compound, Waxes
- Powder mixed with water

250010	1 lb	each	\$38.29	\$36.38	\$36.38	\$34.46
--------	------	------	---------	---------	---------	---------



INFECTION CONTROL

EARLOOP MASKS

HALYARD EARLOOP PROCEDURE MASKS

- Protective three-layer construction
- Appropriate for situations where exposure to blood and/or bodily fluids is not a risk
- Pleat style



47080	Blue	50/box	\$17.99	\$17.09	\$17.09	\$16.19
47095	Pink	50/box	\$17.99	\$17.09	\$17.09	\$16.19
47117	Yellow	50/box	\$17.99	\$17.09	\$17.09	\$16.19

RICHMOND DENTAL ANTIFLUID EARLOOP MASKS

- 100% latex-free, including the earloops
- Engineered with 3-ply spunbound fabric that provides comfort, filtration and excellent breathability
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort
- ASTM Level 1



400502	Blue	50/box	\$13.49	\$12.82	\$12.82	\$12.14
RD400501	Pink	50/box	\$13.49	\$12.82	\$12.82	\$12.14
400509	Lavender	50/box	\$13.49	\$12.82	\$12.82	\$12.14

8

ACE PROCEDURE MASK

- Fluid resistant
- Latex-free
- High BFE, over 99%
- ASTM Level 1



0032418	Blue	50/box	\$15.99	\$15.19	\$15.19	\$14.39
---------	------	--------	----------------	---------	---------	---------

HALYARD FLUIDSHIELD® ASTM EAR-LOOP PROCEDURE MASKS

- Meets ASTM F2100-11 Level 1 standards for fluid-resistance, filtration and breathability
- Fog-free
- SO SOFT* Lining
- ASTM Level 1



41802	White	50/box	\$12.79	\$12.15	\$12.15	\$11.51
-------	-------	--------	----------------	---------	---------	---------

CARDINAL HEALTH SURGICAL MASK

- Pleat style with ties
- Easy to breath through
- For users preferring adjustable masks
- ASTM Level 1



AT7505P	Blue	50/box	\$12.79	\$12.15	\$12.15	\$11.51
---------	------	--------	----------------	---------	---------	---------

3M™ EARLOOP PROCEDURE MASK

- Flat, earloop face mask for general patient care applications
- Standard pleats



1820	Blue	50/box	\$12.99	\$12.34	\$12.34	\$11.69
------	------	--------	----------------	---------	---------	---------

DYNAREX EARLOOP ISOLATION MASK

- Offers a flexible nose piece that helps form a closure for better protection against exposure
- 99% filtration level
- Pleated earloop with glass-free filter
- Latex-free
- ASTM Level 1



DY2208	Yellow	50/box	\$24.79	\$23.55	\$23.55	\$22.31
--------	--------	--------	----------------	---------	---------	---------

DYNAREX EARLOOP PROCEDURE FACE MASK

- Offers a flexible nose piece that helps form a closure for better protection against exposure
- 99% filtration level
- Pleated earloop with glass-free filter
- Latex-free
- ASTM Level 1



DY2201	Blue	50/box	\$22.29	\$21.18	\$21.18	\$20.06
--------	------	--------	----------------	---------	---------	---------

UNIPACK EARLOOP MASKS

- PFE exceeds 99%
- With nose and chin straps
- 3-ply pleated
- Extra soft and comfortable
- High breathability
- Fluid resistant
- Latex and fiberglass free
- ASTM Level 1



ULM6181	Blue	50/box	\$5.29	\$5.03	\$5.03	\$4.76
ULM6189	Pink	50/box	\$5.29	\$5.03	\$5.03	\$4.76
ULM6184	Lavender	50/box	\$5.29	\$5.03	\$5.03	\$4.76
ULM6185	Green	50/box	\$5.29	\$5.03	\$5.03	\$4.76

EARLOOP MASKS

CROSSTEX ISOFLUID® EARLOOP MASKS

- Fluid resistant outer layer
- White tissue inner layer
- ASTM Level 1



GCIBL	Blue	50/box	\$10.79	\$10.25	\$10.25	\$9.71
GCIGR	Green	50/box	\$10.79	\$10.25	\$10.25	\$9.71
GCIPK	Pink	50/box	\$10.79	\$10.25	\$10.25	\$9.71
GCILV	Lavender	50/box	\$10.79	\$10.25	\$10.25	\$9.71
GCITE	Teal	50/box	\$10.79	\$10.25	\$10.25	\$9.71

ADVANCE PLUS EARLOOP PROCEDURE MASK

- Bacterial Filtration Efficiency (BFE): Rated at 99.8%
- Particulate Filtration Efficiency (PFE) rated at 99.9%
- Fluid resistant up to 160mm HG
- Cooler and more breathable than comparable Level 3 masks
- Features an aluminum nose wire for a customized fit
- Conical earloops for better comfort behind the ear
- ASTM Level 3



ADV2316	Blue	50/box	\$8.99	\$8.54	\$8.54	\$8.09
---------	------	--------	--------	--------	--------	--------

CROSSTEX SECURE FIT® EARLOOP PROCEDURAL MASKS

- Conforms to any size or shape face
- Forms a tight seal with aluminum nose and chin pieces that significantly reduce gapping at the top, sides and bottom of the mask
- ASTM Level 2



GCPBSF	Blue	50/box	\$13.29	\$12.63	\$12.63	\$11.96
GCPLVSF	Lavender	50/box	\$13.29	\$12.63	\$12.63	\$11.96

CROSSTEX ULTRA® EARLOOP MASK

- Fluid resistant outer layer and white inner layer
- ASTM Level 3



GCFCXU	Blue	40/box	\$13.99	\$13.29	\$13.29	\$12.59
--------	------	--------	---------	---------	---------	---------

CROSSTEX SECURE FIT® ULTRA® EARLOOP FACE MASK

- Conforms to any size or shape face
- Forms a tight seal with aluminum nose and chin pieces that significantly reduce gapping at the top, sides and bottom of the mask
- ASTM Level 3



GCFCXUSF	Blue	50/box	\$19.79	\$18.80	\$18.80	\$17.81
----------	------	--------	---------	---------	---------	---------

HALYARD PROCEDURE MASK

- Pleat style with earloops
- Fluid-resistant
- ASTM Level 2



HLY62115	Blue	50/box	\$17.29	\$16.43	\$16.43	\$15.56
----------	------	--------	---------	---------	---------	---------

TIDI PROCEDURE MASK

- Pleat style with earloops
- Fog-free design with foam band
- ASTM Level 2



9030	Blue	50/box	\$15.99	\$15.19	\$15.19	\$14.39
------	------	--------	---------	---------	---------	---------

RICHMOND DENTAL PATH-O-GUARD™ ANTI-FLUID EARLOOP MASKS

A Level III Earloop Face Mask You Can Trust

- Exceeds current ASTM-2100 standard for high barrier surgical face masks
- 3-ply
- Comfortable and non-irritating
- Latex-free
- Bacterial filtration efficiency of >98%
- Anti-fluid
- ASTM Level 3



400581	Pink	50/box	\$13.29	\$12.63	\$12.63	\$11.96
400582	Blue	50/box	\$13.29	\$12.63	\$12.63	\$11.96
400589	Lavender	50/box	\$13.29	\$12.63	\$12.63	\$11.96

INFECTION CONTROL

ANTI-FLUID, ANTI-FOG EARLOOP MASKS

RICHMOND DENTAL ANTI-FLUID, ANTI-FOG EARLOOP MASKS

- 100% latex-free, including the earloops
- Engineered with 3-ply spunbound fabric that provides comfort, filtration and excellent breathability
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort
- ASTM Level 1



400503	Pink	50/box	\$13.49	\$12.82	\$12.82	\$12.14
400504	Blue	50/box	\$13.49	\$12.82	\$12.82	\$12.14
400510	Lavender	50/box	\$13.49	\$12.82	\$12.82	\$12.14

CROSSEX ISOFLUID FOGFREE® EARLOOP MASK

- Fluid resistant outer layer
- White tissue inner layer
- Unique FogFree (Comfort Cushion) feature on inner layer
- ASTM Level 1



GCICXB	Blue	40/box	\$14.99	\$14.24	\$14.24	\$13.49
--------	------	--------	---------	---------	---------	---------

RICHMOND DENTAL PATH-O-GUARD™ ANTI-FLUID, ANTI-FOG EARLOOP MASKS

A Level III Earloop Face Mask You Can Trust

- Exceeds current ASTM-2100 standard for high barrier surgical face masks
- 3-ply
- Comfortable and non-irritating
- Latex-free
- Bacterial filtration efficiency of >98%
- Anti-fluid/Anti-fog
- ASTM Level 3



400583	Pink	50/box	\$12.29	\$11.68	\$11.68	\$11.06
400584	Blue	50/box	\$12.29	\$11.68	\$11.68	\$11.06

HALYARD FLUIDSHIELD® FOG-FREE PROCEDURE EARLOOP MASK

This Fluidshield Procedure Mask offers the protection of patented LONCET® breathable film, which helps protect healthcare workers against exposure to blood and bodily fluids.

- ASTM Level 3



Pleat Style, Foam Band

47107	Orange	40/box	\$14.49	\$13.77	\$13.77	\$13.04
-------	--------	--------	---------	---------	---------	---------

CROSSTEX ULTRA® FOGFREE™ EARLOOP MASK

- Fluid resistant outer layer and white inner layer
- Unique No-Fog feature on inner and outer layers
- Exclusive Comfort Cushion
- ASTM Level 3



GCFCX	Blue	40/box	\$18.49	\$17.57	\$17.57	\$16.64
-------	------	--------	---------	---------	---------	---------

EARLOOP MASKS WITH SHIELD

DYNAREX EARLOOP PROCEDURE MASK WITH PLASTIC SHIELD

- Offers a flexible nose piece that helps form a closure for better protection against exposure
- 99% filtration level
- Pleated earloop with glass-free filter
- Latex-free
- ASTM Level 1



DY2202	Blue	50/box	\$31.29	\$29.73	\$29.73	\$28.16
--------	------	--------	---------	---------	---------	---------

RICHMOND DENTAL ANTI-FLUID, ANTI-FOG EARLOOP MASKS WITH SHIELD

- 100% latex-free, including the earloops
- Engineered with 3-ply spunbound fabric that provides comfort, filtration and excellent breathability
- Soft aluminum nosepiece maintains shape when formed to face and is fully enclosed for comfort
- ASTM Level 1



400505	Pink	25/box	\$31.29	\$29.73	\$29.73	\$28.16
400506	Blue	25/box	\$31.29	\$29.73	\$29.73	\$28.16
400511	Lavender	25/box	\$31.29	\$29.73	\$29.73	\$28.16

CROSSTEX ISOFLUID FOGFREE® EARLOOP MASK WITH SHIELD

- Fluid resistant outer layer
- White tissue inner layer
- Attached optically clear, distortion free shield
- ASTM Level 1



GCIPWB	Blue	25/box	\$39.29	\$37.33	\$37.33	\$35.36
--------	------	--------	---------	---------	---------	---------

continued on following page

EARLOOP MASKS WITH SHIELD

CROSSTEX ULTRA® FOGFREE™ EARLOOP MASK WITH SHIELD

- Fluid resistant outer layer
- Fluid resistant white inner layer
- Unique No-Fog feature on inner and outer layers
- Attached optically clear, distortion free shield
- ASTM Level 3



GCPW	Blue	25/box	\$48.79	\$46.35	\$46.35	\$43.91
------	------	--------	---------	---------	---------	---------

HALYARD FLUIDSHIELD® FOG-FREE PROCEDURE EARLOOP MASK WITH SPLASHGUARD® VISOR

This Fluidshield Procedure Mask offers the protection of patented LONCET® breathable film, which helps protect healthcare workers against exposure to blood and bodily fluids.

- ASTM Level 3



Pleat Style with Wraparound Splashguard® Visor, Foam Band

47147	Orange	25/box	\$49.99	\$47.49	\$47.49	\$44.99
00146	Blue	25/box	\$49.99	\$47.49	\$47.49	\$44.99

MÖLNLYCKE BARRIER® SOF-LOOP™ SURGICAL MASK WITH VISOR

- BFE≥98% Protection Level- Not Rated
- Convenient ear loop design
- Splash resistant
- Anti-fog mask



42671	Blue	50/box	\$68.99	\$65.54	\$65.54	\$62.09
-------	------	--------	---------	---------	---------	---------

SENSITIVE EARLOOP MASKS

VALUMAX TISSUESOFT SENSITIVE EARLOOP MASK

This procedure mask provides protection from transfer of microorganisms, fluids and particulate materials while remaining lightweight and very comfortable.

- Soft facial tissue inner layer is extra gentle to skin
- Latex-free
- Easy breathability - Delta P<2.5
- High filtration: PFE 99.7% at 0.1 micron; BFE 99.5%



TissueSoft

5230EWH	White	50/box	\$9.29	\$8.83	\$8.83	\$8.36
---------	-------	--------	--------	--------	--------	--------

SENSITIVE EARLOOP MASKS

VALUMAX EXTRA-SAFE SENSITIVE EARLOOP MASKS

Procedure mask provides excellent protection from transfer of microorganisms, fluids and particulate materials while remaining lightweight and very comfortable.

- Inner layer made of cellulose material, both sensitive and breathable
- Latex-free
- Easy breathability - Delta P<2.5
- High filtration: PFE 99.7% at 0.1 micron; BFE 99.5%



Extra-Safe

5430ESB	Sky blue	50/box	\$11.29	\$10.73	\$10.73	\$10.16
5430ETE	Teal	50/box	\$11.29	\$10.73	\$10.73	\$10.16
5430EWH	White	50/box	\$11.29	\$10.73	\$10.73	\$10.16

CROSSTEX ULTRA® SENSITIVE EARLOOP MASK

- Fluid resistant outer layer
- White hypoallergenic cellulose inner layer
- Extra soft - will NOT lint, tear, shred
- ASTM Level 3



GCFCXS	White	50/box	\$17.29	\$16.43	\$16.43	\$15.56
--------	-------	--------	---------	---------	---------	---------

CROSSTEX ULTRA® SENSITIVE EARLOOP MASK W/ SECURE FIT® TECHNOLOGY

- No chemicals, no inks, no dyes
- Hypoallergenic = comfort without allergic reactions
- ASTM Level 3



Fluid Resistant Outer Layer

GCFCXSSF	White	50/box	\$18.79	\$17.85	\$17.85	\$16.91
----------	-------	--------	---------	---------	---------	---------

FogFree™ With FogFree™ Strip on Inner Layer

GCFCXSFSF	White	40/box	\$20.99	\$19.94	\$19.94	\$18.89
-----------	-------	--------	---------	---------	---------	---------

VALUMAX ULTRA 3-IN-1 SENSITIVE EARLOOP MASKS WITH SPLASH VISOR

- Cellulose material inner layer, completely lint free, no fuzzy tickle feeling on nose and face
- Extra wet strength keeps masks away from nose and mouth
- High filtration, BFE 99.5%, PFE 99.7% at 0.1 micron
- Sensitive, no dye, no chemical, gentle to skin, latex-free
- Breathable, Delta P < 2.2, fiberglass free, odorless



5673ESB	Blue	25/box	\$27.79	\$26.40	\$26.40	\$25.01
---------	------	--------	---------	---------	---------	---------

Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INFECTION CONTROL

MOLDED SURGICAL MASKS

DYNAREX MOLDED SURGICAL MASK


- Manufactured of soft, cool, comfortable materials
- Glass-free filter
- Fluid resistant
- Flared soft edges to fit all facial contours
- Adjustable aluminum nose piece
- Dual head straps for tight facial seal
- Latex-free



DY2203	Blue	50/box	\$15.99	\$15.19	\$15.19	\$14.39
--------	------	--------	----------------	----------------	----------------	----------------

ACE CONE MASK

- Protection Level- Not Rated
- Fiberglass free
- Easy to breath through
- Soft contoured edges
- High BFE, over 99%
- Polyester coated latex strap
- Fluid repellent
- Aluminum nosepiece for better fit



0032420	Blue	50/box	\$9.99	\$9.49	\$9.49	\$8.99
---------	------	--------	---------------	---------------	---------------	---------------

3M™ PARTICULATE RESPIRATOR AND SURGICAL MASK

This health care particulate respirator and surgical mask helps provide respiratory protection against certain airborne biological particles. It is disposable and fluid resistant to splash and spatter of blood and other infectious material.

- NIOSH Approved: N95
- FDA cleared for use as surgical mask
- Helps protect against certain airborne biological particles
- >99% BFE (Bacterial Filtration Efficiency) according to ASTM F2101
- Fluid resistant according to ASTM F1862
- Respirator contains no components made from natural rubber latex
- Disposable



1860	Blue	20/box	\$24.49	\$23.27	\$23.27	\$22.04
------	------	--------	----------------	----------------	----------------	----------------

SURGICAL MASKS WITH TIES

HALYARD CLASSIC® SURGICAL MASK WITH TIES

- Protective three-layer construction
- For protection when exposure to blood and/or bodily fluids is not a concern



Pleat Style

48201	Blue	50/box	\$10.99	\$10.44	\$10.44	\$9.89
-------	------	--------	----------------	----------------	----------------	---------------

3M™ TIE-ON SURGICAL MASK

- With horizontal ties, this mask is ideal for surgical and general-purpose use
- Lightweight, three-layer construction



1818	Blue	50/box	\$14.49	\$13.77	\$13.77	\$13.04
------	------	--------	----------------	----------------	----------------	----------------

HALYARD DUCKBILL® SURGICAL MASK WITH TIES

- Protective three-layer construction
- For protection when exposure to blood and/or bodily fluids is not a concern
- Larger breathing chamber



Pouch Style

48220	Blue	50/box	\$10.99	\$10.44	\$10.44	\$9.89
-------	------	--------	----------------	----------------	----------------	---------------

HALYARD FOG-FREE SURGICAL MASK WITH TIES

- Special construction to reduce fogging on glasses and goggles



Pleat Style

49214	Blue	50/box	\$11.29	\$10.73	\$10.73	\$10.16
-------	------	--------	----------------	----------------	----------------	----------------

DYNAREX SURGICAL MASK WITH TIES

- Pleated w/glass-free filter
- Latex-free
- ASTM Level 1



DY2205	Blue	50/box	\$21.49	\$20.42	\$20.42	\$19.34
--------	------	--------	----------------	----------------	----------------	----------------

SURGICAL MASKS WITH TIES

HALYARD FLUIDSHIELD® FOG-FREE SURGICAL MASK WITH TIES AND WRAPAROUND SPLASHGUARD® ANTI-GLARE VISOR

This Fluidshield Surgical Mask offers the protection of patented LONCET® breathable film, which helps protect healthcare workers against exposure to blood and bodily fluids.

- ASTM Level 3



Pleat Style

48247	Orange	25/box	\$34.79	\$33.05	\$33.05	\$31.31
-------	--------	--------	---------	---------	---------	---------

MÖLNLYCKE SURGINE® ANTI-FOG MASK

- BFE≥98% Protection Level- Not Rated
- Prevents fogging
- Filters to 0.6 micron



Pleated, Non-Glass Design

42381	Green	50/box	\$21.79	\$20.70	\$20.70	\$19.61
-------	-------	--------	---------	---------	---------	---------

MÖLNLYCKE SURGINE® II FACE MASK

- BFE≥98% Protection Level- Not Rated
- Hypoallergenic inner facing
- Filters to 0.6 micron



Pleated, Non-Glass Design

42301	Blue	50/box	\$18.29	\$17.38	\$17.38	\$16.46
-------	------	--------	---------	---------	---------	---------

MÖLNLYCKE SURGINE® II ANTI-FOG MASK

- BFE≥98% Protection Level- Not Rated
- Hypoallergenic inner facing
- Tie string
- Foam backing
- Prevents fogging
- Filters to 0.6 micron



Pleated, Non-Glass Design

42311	Blue	50/box	\$24.29	\$23.08	\$23.08	\$21.86
-------	------	--------	---------	---------	---------	---------

SHARN INC. DISPOSABLE SPLASH SHIELDS

- Anti-fog
- Lightweight for comfort with new SOFT STRAP
- 13 inch impenetrable foam barrier for wide coverage



FS4505	Full face	24/pk	\$76.29	\$72.48	\$72.48	\$68.66
--------	-----------	-------	---------	---------	---------	---------

PALMERO OVERSHIELD EYEWEAR

Specially designed to wear over prescription eyewear, OverShield's impact-resistant, one-piece polycarbonate lens with integral side shields wraps around for a wide field of vision and complete protection. Telescoping and ratcheting arms ensure a great fit no matter what size your prescription frames. Lightweight and comfortable. Meets the new ANSI Z87.1-2003 standard for safety eyewear. Blocks up to 95% of UVA and UVB rays. Black frame/clear lens.



PAL3602		each	\$12.49	\$11.87	\$11.87	\$11.24
---------	--	------	---------	---------	---------	---------

SS WHITE® REUSABLE PROTECTIVE FACE SHIELDS

- Full face coverage
- Made in the USA
- May be used over other PPE- such as masks, goggles or glasses
- Works with or without loupes, Anti-fog, Optically clear
- Cleanable (up to 10 times if cleaning guidance is properly followed)
- Extra Width, 360mm



PPEFSXL30	Face shield	each	\$24.79	\$23.55	\$23.55	\$22.31
PPERSXL10	Replacement shield	each	\$10.79	\$10.25	\$10.25	\$9.71

INFECTION CONTROL

PALMERO PROVISION™ SEE-BREEZ EYEWEAR

Ultra-lightweight frames with sculpted dual lenses and telescoping arms ensure that these are the most comfortable fitting eyewear you will ever experience. They are also great for adults with smaller framed faces or for children. Completely fog-free (even with a mask on), scratch-resistant and adjustable. Meets the new ANSI Z87.1-2003 standard for safety eyewear.



Colored Frame - Clear Lens

3560BL	Black	each	\$14.79	\$14.05	\$14.05	\$13.31
3560B	Blue	each	\$14.79	\$14.05	\$14.05	\$13.31
3560PL	Platinum	each	\$14.79	\$14.05	\$14.05	\$13.31

Black Frame - Tinted Lens

3561BL	Grey lens	each	\$15.29	\$14.53	\$14.53	\$13.76
--------	-----------	------	---------	---------	---------	---------

Replacement Lens Set

3560R	Clear	each	\$10.99	\$10.44	\$10.44	\$9.89
3561R	Grey	each	\$11.99	\$11.39	\$11.39	\$10.79

PALMERO PROVISION™ ECONO WRAP EYEWEAR

Feather-weight and low-cost eye protection that offers a universal fit and unsurpassed protection. Features a sleek, compact appearance and superior optical quality. Perfect for patient and staff use. Fog-free and scratch-resistant. Meets the new ANSI Z87.1-2003 standard for safety eyewear.



3601C	Clear lens	each	\$5.99	\$5.69	\$5.69	\$5.39
3601G	Grey lens	each	\$5.99	\$5.69	\$5.69	\$5.39
3601CB	Clear lens	10/box	\$56.49	\$53.67	\$53.67	\$50.84

PALMERO MINI ECONO WRAP EYEWEAR

Same styling as our Econo Wrap Eyewear except with a narrower width for children and people with smaller framed faces. Meets the new ANSI Z87.1-2003 standard for safety eyewear.



3607C	Clear lens	each	\$5.99	\$5.69	\$5.69	\$5.39
3607G	Grey lens	each	\$5.99	\$5.69	\$5.69	\$5.39

PALMERO TECH SPECS™ SAFETY EYEWEAR

This eyewear has exceptional styling in a lightweight design that provides excellent protection and comfort. "Hugs" the face with its contoured design and rubberized (non-latex) attachments on the nosepiece and ear rests. Fog-free and scratch-resistant.



With Grey Frame

3710C	Clear lens	each	\$16.29	\$15.48	\$15.48	\$14.66
3710G	Grey lens	each	\$16.29	\$15.48	\$15.48	\$14.66

PALMERO PROVISION™ CRYSTAL CLEAR OPTICAL CLEANING TOWELETTES

Specially designed towelettes are premoistened with an anti-fog/anti-static lens cleaning and polishing solution for the ultimate in convenience. Quickly and effectively cleans dirt, grease, oils and blood-borne pathogens from eyewear with its fast evaporating formula.



3534	6" x 6.75"	160/canister	\$15.79	\$15.00	\$15.00	\$14.21
------	------------	--------------	---------	---------	---------	---------

DYNAREX NURSE CAPS

- Lightweight, cool, comfortable material
- Soft elastic headband



Large, 21"

2103	Blue	100/box	\$13.49	\$12.82	\$12.82	\$12.14
------	------	---------	---------	---------	---------	---------

X-Large, 24"

DY2113	Blue	100/box	\$14.99	\$14.24	\$14.24	\$13.49
--------	------	---------	---------	---------	---------	---------

MEDLINE BOUFFANT CAPS

- Latex-free



CRI1002	White, 21"	100/box	\$14.99	\$14.24	\$14.24	\$13.49
CRI1003	White, large, 24"	100/box	\$15.99	\$15.19	\$15.19	\$14.39
CRI1004	Blue, large, 24"	100/box	\$15.99	\$15.19	\$15.19	\$14.39

HALYARD BOUFFANT CAP

- Made of breathable one-layer Spunbond material
- Appropriate for use when no potential fluid contact is anticipated



Large

69801	Blue	100/box	\$22.29	\$21.18	\$21.18	\$20.06
-------	------	---------	---------	---------	---------	---------

HALYARD PROTECTIVE SURGICAL CAP

- Made from a protective three-layer fabric
- Best suited for use when there is light fluid contamination from splashing or spraying anticipated
- Universal size
- Tie-on



69520	Blue	100/box	\$29.79	\$28.30	\$28.30	\$26.81
-------	------	---------	---------	---------	---------	---------

INFECTION CONTROL

MÖINLYCKE BARRIER® SURGICAL CAP

- Traditional design
- Cool, comfortable, absorbent
- Wide ties provide secure closure

MOL620500	Green	100/box	\$38.79	\$36.85	\$36.85	\$34.91
-----------	-------	---------	---------	---------	---------	---------

ALLEGIANCE (BAXTER) EASY-TIE™ SURGEONS CAP

- Non-woven cap featuring traditional tie-on styling
- Made of absorbent Kaycel® material with a spunbonded polypropylene top for breathability and comfort
- Universal size

4359	Blue	100/box	\$25.79	\$24.50	\$24.50	\$23.21
------	------	---------	---------	---------	---------	---------



PRECEPT® SURGEONS CAP

- Universal size
- Kaycel material

5200	Blue	125/box	\$31.29	\$29.73	\$29.73	\$28.16
------	------	---------	---------	---------	---------	---------

GRAFCO® LATEX FINGER COTS - NON-MEDICAL

- For industrial use only: easy grip for assembly, preventing finger prints, covering finger wounds.
- ⚠️ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov



3908M	Medium	144/box	\$4.29	\$4.08	\$4.08	\$3.86
3908L	Large	144/box	\$4.29	\$4.08	\$4.08	\$3.86

MEDLINE SPUNBONDED NON-SKID SHOE COVERS



CRI2002	Blue, universal size	100/box	\$20.79	\$19.75	\$19.75	\$18.71
---------	----------------------	---------	---------	---------	---------	---------

HALYARD X-TRA TRACTION® SHOE COVER

- 3-layer fabric, adhesive traction



69252	Blue, universal size	100/box	\$24.79	\$23.55	\$23.55	\$22.31
-------	----------------------	---------	---------	---------	---------	---------

HALYARD BASIC SHOE COVERS

- Traction
- No potential fluid contact



69551	Blue, universal size	300/case	\$68.99	\$65.54	\$65.54	\$62.09
-------	----------------------	----------	---------	---------	---------	---------

DYNAREX SHOE COVERS

- Manufactured of spunbound polypropylene
- Sewn seams, not glued - offers greater strength
- Elastic opening provides a snug comfortable fit
- Non-conductive
- Latex-free
- Universal Size



2131	Blue	300/case	\$43.99	\$41.79	\$41.79	\$39.59
DY2132	Blue, non-skid	300/case	\$45.79	\$43.50	\$43.50	\$41.21

LATEX-FREE EXAM GLOVES

DYNAREX SAFE-TOUCH BLUE NITRILE EXAM GLOVES

- Textured fingertips, and a minimum thickness of 3 mil. at the palm and 4 mil. at the fingertip
- Blue color



Latex-Free, Powder-Free

2511	Small	100/box	\$26.49	\$25.17	\$25.17	\$23.84
2512	Medium	100/box	\$26.49	\$25.17	\$25.17	\$23.84
2513	Large	100/box	\$26.49	\$25.17	\$25.17	\$23.84
2514	X-Large	100/box	\$26.29	\$24.98	\$24.98	\$23.66

DYNAREX SAFE-TOUCH BLACK NITRILE EXAM GLOVES

- Textured fingertips, and a minimum thickness of 4 mil.
- Black color



Latex-Free, Powder-Free

2521	Small	100/box	\$24.99	\$23.74	\$23.74	\$22.49
2522	Medium	100/box	\$24.99	\$23.74	\$23.74	\$22.49
2523	Large	100/box	\$24.99	\$23.74	\$23.74	\$22.49
2524	X-Large	100/box	\$24.99	\$23.74	\$23.74	\$22.49

continued on following page

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

INFECTION CONTROL

LATEX-FREE EXAM GLOVES

DYNAREX DYNAPLUS™ NITRILE EXAM GLOVES

- Superior comfort
- Textured fingertips
- Periwinkle color



Latex-Free, Powder-Free

2515	X-Small	200/box	\$24.99	\$23.74	\$23.74	\$22.49
2516	Small	200/box	\$24.99	\$23.74	\$23.74	\$22.49
2517	Medium	200/box	\$24.99	\$23.74	\$23.74	\$22.49
2518	Large	200/box	\$24.99	\$23.74	\$23.74	\$22.49
2519	X-Large	180/box	\$24.99	\$23.74	\$23.74	\$22.49

HALYARD KC100 LAVENDER NITRILE EXAM GLOVES

- Non-sterile
- Textured fingertips, beaded cuff
- Higher count/box saves space and waste
- Lavender color



Latex-Free, Powder-Free

52816	X-small	250/box	\$30.29	\$28.78	\$28.78	\$27.26
52817	Small	250/box	\$30.29	\$28.78	\$28.78	\$27.26
52818	Medium	250/box	\$30.29	\$28.78	\$28.78	\$27.26
52819	Large	250/box	\$30.29	\$28.78	\$28.78	\$27.26
52820	X-Large	230/box	\$30.29	\$28.78	\$28.78	\$27.26

HALYARD PURPLE NITRILE® EXAM GLOVES

- Textured fingertips, beaded cuff
- Purple color



9.5" Long, Latex-Free, Powder-Free

55081	Small	100/box	\$25.99	\$24.69	\$24.69	\$23.39
55082	Medium	100/box	\$25.99	\$24.69	\$24.69	\$23.39
55083	Large	100/box	\$25.99	\$24.69	\$24.69	\$23.39
55084	X-Large	90/box	\$25.99	\$24.69	\$24.69	\$23.39

9.5" Long, Latex-Free, Powder-Free, Sterile

55091	Small	100/box	\$74.79	\$71.05	\$71.05	\$67.31
55092	Medium	100/box	\$74.79	\$71.05	\$71.05	\$67.31
55093	Large	100/box	\$74.79	\$71.05	\$71.05	\$67.31

Latex-Free, Powder-Free

53430	X-Small	100/box	\$22.79	\$21.65	\$21.65	\$20.51
53431	Small	100/box	\$22.79	\$21.65	\$21.65	\$20.51
53432	Medium	100/box	\$22.79	\$21.65	\$21.65	\$20.51
53433	Large	100/box	\$22.79	\$21.65	\$21.65	\$20.51

DASH® ALASTA™ NITRILE EXAM GLOVES

- Performs like latex
- Extreme comfort with precision fit
- Textured fingertips
- Available in 100/box and 200/box
- Periwinkle color



Latex-Free, Powder-Free

ALS100XS	X-Small	100/box	\$26.79	\$25.45	\$25.45	\$24.11
ALS100S	Small	100/box	\$26.79	\$25.45	\$25.45	\$24.11
ALS100M	Medium	100/box	\$26.79	\$25.45	\$25.45	\$24.11
ALS100L	Large	100/box	\$26.79	\$25.45	\$25.45	\$24.11
ALS100XL	X-Large	100/box	\$26.79	\$25.45	\$25.45	\$24.11
ALS200XS	X-Small	200/box	\$50.79	\$48.25	\$48.25	\$45.71
ALS200S	Small	200/box	\$50.79	\$48.25	\$48.25	\$45.71
ALS200M	Medium	200/box	\$50.79	\$48.25	\$48.25	\$45.71
ALS200L	Large	200/box	\$50.79	\$48.25	\$48.25	\$45.71
ALS200XL	X-Large	200/box	\$50.79	\$48.25	\$48.25	\$45.71

DASH® BLACK MAXX™ NITRILE EXAM GLOVES

- Textured surface
- Black color



Latex-Free, Powder-Free

BMN100S	Small	100/box	\$41.79	\$39.70	\$39.70	\$37.61
BMN100M	Medium	100/box	\$41.79	\$39.70	\$39.70	\$37.61
BMN100L	Large	100/box	\$41.79	\$39.70	\$39.70	\$37.61
BMN100XL	X-Large	100/box	\$41.79	\$39.70	\$39.70	\$37.61

HALYARD PURPLE NITRILE® XTRA STERILE EXAM GLOVES

- Sterile
- Maximum coverage and protection
- Textured fingertips
- For moderate to high exposure to fluids or chemicals
- Natural rubber



Latex-Free, Powder-Free

14260	Small	50 pairs/box	\$79.79	\$75.80	\$75.80	\$71.81
14261	Medium	50 pairs/box	\$79.79	\$75.80	\$75.80	\$71.81
14262	Large	50 pairs/box	\$79.79	\$75.80	\$75.80	\$71.81
14263	X-Large	50 pairs/box	\$79.79	\$75.80	\$75.80	\$71.81

LATEX-FREE EXAM GLOVES

HALYARD PURPLE NITRILE® MAX EXAM GLOVES

The new standard for protection in the sterile processing department. Designed specifically to meet the increased protection needs of Central Sterilization technicians. Durable, easy-to-don PURPLE NITRILE® MAX offers CS techs these critical safety features:

- **Longer Length** - PURPLE NITRILE® MAX cuff extends a full 16 inches, providing the right level of protection from dangerous microorganisms, all the way to the mid-forearm.
- **Heavier Weight** - At 20 grams, PURPLE NITRILE® MAX is more than twice the weight of PURPLE NITRILE-XTRA®, for added confidence when handling instruments and decontamination chemicals.
- **Thicker Fingertips** - Its increased fingertip thickness provides almost double the protection of any other exam glove in the Halyard portfolio.
- **Fully Textured** - PURPLE NITRILE® MAX features fully textured fingertips and palm, designed for excellent wet and dry grip.



HALYARD® AQUASOFT NITRILE EXAM GLOVES

- Ultra-thin with tactile sensitivity comparable to latex
- Textured Palm and Fingertips
- Violet blue color



Latex-Free, Powder-Free

43932	X-Small	300/box	\$33.79	\$32.10	\$32.10	\$30.41
43933	Small	300/box	\$33.79	\$32.10	\$32.10	\$30.41
43934	Medium	300/box	\$33.79	\$32.10	\$32.10	\$30.41
43935	Large	300/box	\$33.79	\$32.10	\$32.10	\$30.41

Latex-Free, Powder-Free

44992	Small	50/box	\$68.99	\$65.54	\$65.54	\$62.09
44993	Medium	50/box	\$68.99	\$65.54	\$65.54	\$62.09
44994	Large	50/box	\$68.99	\$65.54	\$65.54	\$62.09
44995	X-Large	50/box	\$68.99	\$65.54	\$65.54	\$62.09

HALYARD® BLACK-FIRE NITRILE EXAM GLOVES

- NFPA 1999-2013 Certified
- Patented nitrile formulation increases comfort and strength
- Reversible
- Black/orange color




Latex-Free, Powder-Free

44756	Small	150/box	\$34.29	\$32.58	\$32.58	\$30.86
44757	Medium	150/box	\$34.29	\$32.58	\$32.58	\$30.86
44758	Large	150/box	\$34.29	\$32.58	\$32.58	\$30.86

HALYARD STERLING® NITRILE EXAM GLOVES

A new manufacturing process boosts the tensile strength of nitrile, and also allows for a reduction in glove thickness which provides a latex-free glove with the flexibility and sensitivity of natural rubber latex.

- Non-sterile
- Textured fingertips, beaded cuff
- Higher count, saves space
- Gray color



Latex-Free, Powder-Free

50705	X-Small	200/box	\$42.29	\$40.18	\$40.18	\$38.06
50706	Small	200/box	\$42.29	\$40.18	\$40.18	\$38.06
50707	Medium	200/box	\$42.29	\$40.18	\$40.18	\$38.06
50708	Large	200/box	\$42.29	\$40.18	\$40.18	\$38.06
50709	X-Large	170/box	\$42.29	\$40.18	\$40.18	\$38.06

SEMPERMED STARMED® PLUS™ NITRILE EXAM GLOVES

- Textured fingertips
- Beaded cuff
- Violet blue color



Latex-Free, Powder-Free

SMNP301	X-Small	300/box	\$43.29	\$41.13	\$41.13	\$38.96
SMNP302	Small	300/box	\$43.29	\$41.13	\$41.13	\$38.96
SMNP303	Medium	300/box	\$43.29	\$41.13	\$41.13	\$38.96
SMNP304	Large	300/box	\$43.29	\$41.13	\$41.13	\$38.96
SMNP305	X-Large	250/box	\$43.29	\$41.13	\$41.13	\$38.96

continued on following page

INFECTION CONTROL

LATEX-FREE EXAM GLOVES

SEMPERMED STARMED® ROSE NITRILE EXAM GLOVES WITH ALOE AND VITAMIN E

StarMed® Rose Nitrile with Aloe and Vitamin E has what it takes to brighten any environment. An aloe and vitamin E coating soothes hands and helps prevent skin irritation, while StarMed® Rose's formulation provides enhanced protection, comfort and tactile sensitivity.



- Aloe and vitamin E coating
- Beaded cuff
- Bright and cheery rose color

Latex-Free, Powder-Free

SMNR201	X-Small	200/box	\$58.79	\$55.85	\$55.85	\$52.91
SMNR202	Small	200/box	\$58.79	\$55.85	\$55.85	\$52.91
SMNR203	Medium	200/box	\$58.79	\$55.85	\$55.85	\$52.91

(Your purchase of this product helps Sempermed proudly support the Prostate Cancer Foundation, www.pcf.org.)

MICROFLEX ULTRAFORM® NITRILE EXAM GLOVES

- First ambidextrous glove available in half sizes
- Beaded cuff
- Textured fingertips
- Smart packaging
- Cobalt blue color



Latex-Free, Powder-Free

UF524XS	X-Small	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524XSS	X-Small/small	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524S	Small	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524SM	Small/medium	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524M	Medium	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524ML	Medium/large	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524L	Large	300/box	\$94.79	\$90.05	\$90.05	\$85.31
UF524XL	X-Large	250/box	\$94.79	\$90.05	\$90.05	\$85.31

SEMPERMED TENDER TOUCH® NITRILE EXAM GLOVES

The right combination of fit, feel, tactile sensitivity, strength and value in a latex-free exam glove.



- Synthetic nitrile rubber (NBR)
- Textured fingertips
- Beaded cuff
- Periwinkle color

Latex-Free, Powder-Free

TTNF201	X-Small	200/box	\$44.49	\$42.27	\$42.27	\$40.04
TTNF202	Small	200/box	\$44.49	\$42.27	\$42.27	\$40.04
TTNF203	Medium	200/box	\$44.49	\$42.27	\$42.27	\$40.04
TTNF204	Large	200/box	\$44.49	\$42.27	\$42.27	\$40.04

ANSELL MICRO-TOUCH® MICRO-THIN NITRILE EXAM GLOVES

- Thinner to enhance sensitivity
- Textured fingertips
- Beaded cuff
- Blue color



Latex-Free, Powder-Free

6034300	X-Small	200/box	\$37.79	\$35.90	\$35.90	\$34.01
6034301	Small	200/box	\$37.79	\$35.90	\$35.90	\$34.01
6034302	Medium	200/box	\$37.79	\$35.90	\$35.90	\$34.01
6034303	Large	200/box	\$37.79	\$35.90	\$35.90	\$34.01
6034304	X-Large	200/box	\$37.79	\$35.90	\$35.90	\$34.01

MICROFLEX XCEED® NITRILE EXAM GLOVES

- Proven to reduce muscle effort
- Industry leading protection and strength
- Textured fingertips
- Blue color



Latex-Free, Powder-Free

XC310XS	X-Small	250/box	\$48.29	\$45.88	\$45.88	\$43.46
XC310SM	Small	250/box	\$48.29	\$45.88	\$45.88	\$43.46
XC310MED	Medium	250/box	\$48.29	\$45.88	\$45.88	\$43.46
XC310LG	Large	250/box	\$48.29	\$45.88	\$45.88	\$43.46
XC310XL	X-Large	230/box	\$48.29	\$45.88	\$45.88	\$43.46

ANSELL MICRO-TOUCH® AFFINITY™ POLYCHLOROPRENE EXAM GLOVES

- Textured fingertips provide excellent touch sensitivity
- 100% neoprene
- Green color



Latex-Free, Powder-Free

PE3770	X-Small	100/box	\$22.29	\$21.18	\$21.18	\$20.06
3771	Small	100/box	\$22.29	\$21.18	\$21.18	\$20.06
3772	Medium	100/box	\$22.29	\$21.18	\$21.18	\$20.06
3773	Large	100/box	\$22.29	\$21.18	\$21.18	\$20.06
3774	X-Large	100/box	\$22.29	\$21.18	\$21.18	\$20.06

LATEX-FREE EXAM GLOVES

HALYARD FLEXAPRENE® GREEN CHLOROPRENE EXAM GLOVES

The new FLEXAPRENE® GREEN exam glove is an advanced chloroprene technology that gives you the feel and tactile sensitivity of latex without the risk of Type I allergic reactions to natural rubber latex proteins - and with the protection you'd expect from a synthetic.

- Tactile sensitivity compared to latex
- Soft to the touch, yet powerful in grip
- Polymer coated for improved donning
- Not made with natural rubber latex
- Powder-free chloroprene technology
- Green color



Latex-Free, Powder-Free

Product Code	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
44792	X-Small	200/box	\$28.79	\$27.35	\$27.35	\$25.91
44793	Small	200/box	\$28.79	\$27.35	\$27.35	\$25.91
44794	Medium	200/box	\$28.79	\$27.35	\$27.35	\$25.91
44795	Large	200/box	\$28.79	\$27.35	\$27.35	\$25.91
44796	X-Large	180/box	\$28.79	\$27.35	\$27.35	\$25.91

MICROFLEX NEOGARD® CHLOROPRENE EXAM GLOVES

Neogard is made from a unique, non-latex formulation that delivers outstanding grip in wet and dry environments, perfect for dealing with explorers and handpieces. With unparalleled tactile sensitivity and strength, Neogard gloves are ideal for every task, from cleanings to orthodontics.

- Soft and stretchy chloroprene formulation is highly comfortable with exceptional tactile sensitivity
- Polymer coated for easy donning, especially on damp hands
- Non-latex formula
- Green color



Latex-Free, Powder-Free

Product Code	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
C520	X-Small	200/box	\$18.99	\$18.04	\$18.04	\$17.09
C521	Small	200/box	\$18.99	\$18.04	\$18.04	\$17.09
C522	Medium	200/box	\$18.99	\$18.04	\$18.04	\$17.09
C523	Large	200/box	\$18.99	\$18.04	\$18.04	\$17.09

LATEX EXAM GLOVES

DYNAREX SAFE-TOUCH™ LATEX EXAM GLOVES

- Rugged construction provides tough strength for light-medium use
- Comfortable fit with acute sensitivity for better functionality
- Cream color



Latex, Powder-Free

Product Code	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
2336	Small	100/box	\$14.99	\$14.24	\$14.24	\$13.49
2337	Medium	100/box	\$14.99	\$14.24	\$14.24	\$13.49
DY2338	Large	100/box	\$14.99	\$14.24	\$14.24	\$13.49
2339	X-Large	100/box	\$14.99	\$14.24	\$14.24	\$13.49

MICROFLEX DIAMOND GRIP PLUS™ LATEX EXAM GLOVES

- The #1 branded latex in the dental market*
- A perfect multi-purpose latex glove
- Fully textured for superior grip even when wet
- Natural color



Latex, Powder-Free

Product Code	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
DGP350S	Small	100/box	\$16.79	\$15.95	\$15.95	\$15.11
DGP350M	Medium	100/box	\$16.79	\$15.95	\$15.95	\$15.11
DGP350L	Large	100/box	\$16.79	\$15.95	\$15.95	\$15.11
DGP350XL	X-Large	100/box	\$16.79	\$15.95	\$15.95	\$15.11

*Strategic Dental Marketing

SEMPERMED STARMED® LATEX EXAM GLOVES

StarMed powder-free latex exam gloves are designed to provide the perfect combination of protection, tactile sensitivity and value. Our advanced manufacturing technology and online chlorination process creates a glove that is strong, flexible and goes on with ease.

- Textured
- Cream color



Latex, Powder-Free

Product Code	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
SM101	X-Small	100/box	\$15.49	\$14.72	\$14.72	\$13.94
SM102	Small	100/box	\$15.49	\$14.72	\$14.72	\$13.94
SM103	Medium	100/box	\$15.49	\$14.72	\$14.72	\$13.94
SM104	Large	100/box	\$15.49	\$14.72	\$14.72	\$13.94
SM105	X-Large	90/box	\$15.49	\$14.72	\$14.72	\$13.94

continued on following page

INFECTION CONTROL

LATEX EXAM GLOVES

ANSELL MICRO-TOUCH™ PLUS LATEX EXAM GLOVES

- Beaded cuff
- Fully Textured
- Polymer Coated Lining
- Cream color



Latex, Powder-Free

6015300	X-Small	150/box	\$32.29	\$30.68	\$30.68	\$29.06
6015301	Small	150/box	\$32.29	\$30.68	\$30.68	\$29.06
6015302	Medium	150/box	\$32.29	\$30.68	\$30.68	\$29.06
6015303	Large	150/box	\$32.29	\$30.68	\$30.68	\$29.06
6015304	X-Large	150/box	\$32.29	\$30.68	\$30.68	\$29.06

SEMPERMED SEMPERCARE® LATEX EXAM GLOVES

SemperCare® Latex is the second generation of Sempermed's popular co-polymer latex glove, boasting enhanced elasticity, tactile sensitivity and sensation. This glove has a slick inner coating for ease of donning and added comfort.



Features

- 50 micrograms or less of total water-extractable protein per gram
- Polymer coating
- Slick inner coating
- Micro-rough textured surface

Benefits

- Reduced potential for allergic reaction
- Superior elasticity and comfortable fit
- Enhanced strength and durability, quick and easy to don
- Improved tactile sensitivity with increased gripping characteristics
- Cream color

Latex, Powder-Free

SCLT101	X-Small	100/box	\$19.99	\$18.99	\$18.99	\$17.99
SCLT102	Small	100/box	\$19.99	\$18.99	\$18.99	\$17.99
SCLT103	Medium	100/box	\$19.99	\$18.99	\$18.99	\$17.99
SCLT104	Large	100/box	\$19.99	\$18.99	\$18.99	\$17.99
SCLT105	X-Large	100/box	\$19.99	\$18.99	\$18.99	\$17.99

DYNAREX ALOETEX™ LATEX EXAM GLOVES WITH ALOE

- Exam gloves that help soothe and moisturize dry, chapped hands
- Coated with skin-soothing aloe
- Finger textured
- Tapered cuff for a more snug and comfortable fit around the wrist
- Green color



Latex, Powder-Free

6716	X-Small	100/box	\$18.79	\$17.85	\$17.85	\$16.91
DY6717	Small	100/box	\$18.79	\$17.85	\$17.85	\$16.91
DY6718	Medium	100/box	\$18.79	\$17.85	\$17.85	\$16.91
6719	Large	100/box	\$18.79	\$17.85	\$17.85	\$16.91

SEMPERMED BEST TOUCH® LATEX EXAM GLOVES W/ ALOE & VITAMIN E

Best Touch Latex with Aloe and Vitamin E provides superior protection and comfort with the added benefit of an aloe and vitamin E inner coating. This exclusive formula allows for an enhanced inner coating which helps prevent your skin from becoming damaged or irritated.



- Micro-rough textured surface
- Cream color

Latex, Powder-Free

BTLA101	X-Small	100/box	\$15.29	\$14.53	\$14.53	\$13.76
BTLA102	Small	100/box	\$15.29	\$14.53	\$14.53	\$13.76
BTLA103	Medium	100/box	\$15.29	\$14.53	\$14.53	\$13.76
BTLA104	Large	100/box	\$15.29	\$14.53	\$14.53	\$13.76
BTLA105	X-Large	90/box	\$15.29	\$14.53	\$14.53	\$13.76

ANSELL MICRO-TOUCH NEXTSTEP® LATEX EXAM GLOVES W/GLYCEROL

- Powder-free to help minimize allergenic and aeroallergen transmission
- Textured fingertips
- Green color



Latex, Powder-Free

PE3200	X-Small	100/box	\$20.79	\$19.75	\$19.75	\$18.71
3201	Small	100/box	\$20.79	\$19.75	\$19.75	\$18.71
3202	Medium	100/box	\$20.79	\$19.75	\$19.75	\$18.71
3203	Large	100/box	\$20.79	\$19.75	\$19.75	\$18.71
3204	X-Large	80/box	\$20.79	\$19.75	\$19.75	\$18.71

8

LATEX-FREE SURGICAL GLOVES

**ANSELL GAMMEX®
NON-LATEX
SURGICAL GLOVES**

- 100% neoprene; Smooth finish
- Beaded cuff; SureFit Technology™
- Polyurethane inner coating
- Green color



Latex-Free, Powder-Free

SKU	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
8511	Size 5.5	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8512	Size 6.0	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8513	Size 6.5	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8514	Size 7.0	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8515	Size 7.5	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8516	Size 8.0	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8517	Size 8.5	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29
8518	Size 9.0	50pr/box	\$156.99	\$149.14	\$149.14	\$141.29

**ANSELL GAMMEX® NON-LATEX
PI SURGICAL GLOVES**

- 100% polyisoprene; Smooth finish
- Beaded cuff; SureFit Technology™
- Polyurethane inner coating
- Pale green color



Latex-Free, Powder-Free

SKU	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
20685255	Size 5.5	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685260	Size 6.0	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685265	Size 6.5	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685270	Size 7.0	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685275	Size 7.5	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685280	Size 8.0	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685285	Size 8.5	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51
20685290	Size 9.0	50pr/box	\$162.79	\$154.65	\$154.65	\$146.51

**SEMPERMED® SYNTEGRA
CR NON-LATEX
CHLOROPRENE SURGICAL
GLOVES**

For the latex sensitive, Sempermed Syntegra CR offers a similar fit and feel to latex. These quality surgical gloves are manufactured with polychloroprene, also known as neoprene, a trusted material with over 70 years of proven performance. Excellent tactile sensitivity, while the textured surface and beaded cuff enhance gripping and protection.



- Beaded cuff
- Cream color

Latex-Free, Powder-Free

SKU	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
SCR600	Size 6.0	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR650	Size 6.5	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR700	Size 7.0	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR750	Size 7.5	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR800	Size 8.0	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR850	Size 8.5	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16
SCR900	Size 9.0	40pr/box	\$101.29	\$96.23	\$96.23	\$91.16

**MEDLINE
SENSICARE® PI
SLT LATEX-FREE
SURGICAL
GLOVES**

- Standard lightly textured (SLT) for surgeons who prefer a surgical glove with more sensitivity
- Isolex™ polyisoprene material
- Anti-roll-down cuff technology significantly reduces the roll down common with many surgical gloves
- Cream color



Latex-Free, Powder-Free

SKU	Size	Quantity	Regular Price	Member Price	Premier Price	Premier+ Price
MSG1555	Size 5.5	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1560	Size 6.0	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1565	Size 6.5	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1570	Size 7.0	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1575	Size 7.5	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1580	Size 8.0	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1585	Size 8.5	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81
MSG1590	Size 9.0	50pr/box	\$149.79	\$142.30	\$142.30	\$134.81

INFECTION CONTROL

LATEX SURGICAL GLOVES

SEMPERMED® SUPREME TEXTURED LATEX SURGICAL GLOVES

When a glove feels this right, nothing can get in the way of your performance. Sempermed Supreme surgical gloves are anatomically designed for comfort and precision. Powder-free to reduce potential aerosolization of latex proteins, these latex gloves have a special multi-layer copolymer composition to provide strength and protection during surgery. Secured with a beaded cuff.



- Cream color

Latex, Powder-Free

SPFP600	Size 6.0	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP650	Size 6.5	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP700	Size 7.0	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP750	Size 7.5	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP800	Size 8.0	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP850	Size 8.5	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36
SPFP900	Size 9.0	50pr/box	\$69.29	\$65.83	\$65.83	\$62.36

ANSELL ENCORE® SENSI-TOUCH® PF LATEX SURGICAL GLOVES

Provides Superior Strength, Reduced Hand Fatigue, with a Beaded Cuff for Extra Protection



- Natural rubber latex
- Beaded cuff
- Damp hand donnability; Specially designed to don like a powdered glove
- Anti-slip smooth finish for a comfortable, secure grip
- White color

Latex, Powder-Free

7821PF	Size 5.5	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7822PF	Size 6.0	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7823PF	Size 6.5	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7824PF	Size 7.0	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7825PF	Size 7.5	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7826PF	Size 8.0	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7827PF	Size 8.5	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99
7828PF	Size 9.0	50pr/box	\$49.99	\$47.49	\$47.49	\$44.99

ANSELL ENCORE® PERRY® STYLE 42® PF LATEX SURGICAL GLOVES

Exceptional Value Without Sacrificing Comfort

- Natural rubber latex
- Inverted, beaded cuff
- Damp hand donnability; Specially designed to don like a powdered glove
- Anti-slip smooth finish for a comfortable, secure grip
- White color



Latex, Powder-Free

5711100PF	Size 5.5	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711101PF	Size 6.0	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711102PF	Size 6.5	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711103PF	Size 7.0	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711104PF	Size 7.5	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711105PF	Size 8.0	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711106PF	Size 8.5	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94
5711107PF	Size 9.0	50pr/box	\$55.49	\$52.72	\$52.72	\$49.94

DYNAREX LATEX SURGICAL GLOVES

- One pair (one right, one left) per paper wallet, each wallet packed in a peel-open pouch
- Maximum strength with no loss of tactile sensitivity
- Textured surface in palm and finger tips
- Cream color



Latex, Powder-Free

DY2360	Size 6.0	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
2365	Size 6.5	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
2370	Size 7.0	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
2375	Size 7.5	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
DY2380	Size 8.0	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
2385	Size 8.5	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09
2390	Size 9.0	50pr/box	\$38.99	\$37.04	\$37.04	\$35.09

8

DYNAREX ISOLATION GOWN

- Non-sterile
- Offers economical, comfortable and reliable protection
- Full length for maximum coverage
- Elastic cuffs
- Sewn seams provide superior strength
- Extra length on waist ties allows gown to be easily secured in front
- Latex-free
- Yellow



DY2141	Universal	50/case	\$115.49	\$109.72	\$109.72	\$103.94
--------	-----------	---------	----------	----------	----------	----------

MAYTEX ISOLATION GOWN

Made of non-woven polypropylene material coated with water repellent agent for added protection. Material is breathable which prevents irritating hot temperature build up. Waist ties are designed to tie in front for ease of use. For most routine medical and dental procedures.



- One size fits most
- Knit cuff
- Latex-free
- Blue
- AAMI Level 1

CV6010B	One size fits most	50/case	\$96.79	\$91.95	\$91.95	\$87.11
---------	--------------------	---------	---------	---------	---------	---------

VALUMAX® ISOLATION GOWN

- Non-sterile
- Lightweight non-woven fabric
- Splash-resistant
- Tie back at neck and waist
- Regular size fits most
- Knit cuffs
- Economical
- Latex-free
- Blue
- AAMI Level 1

CV3260B		50/case	\$219.99	\$208.99	\$208.99	\$197.99
---------	--	---------	----------	----------	----------	----------

MEDLINE NONSTERILE NONREINFORCED ECLIPSE GOWNS

- High-quality, SMS fabric enhances the Eclipse line with increased breathability, lighter weight, higher hydrostatic head and better wicking/moisture evaporation
- Blue
- AAMI Level 2 protection



16836NS	X-Large	60/case	\$293.79	\$279.10	\$279.10	\$264.41
---------	---------	---------	----------	----------	----------	----------

BE THE OFFICE HERO

Save Time and Money with...



REPLENISH

Our Replenish program allows you to set up scheduled deliveries for the products your practice uses on a regular basis.

Ensure that your practice never runs out of essential items that keep your office running smoothly.

SHOP from hundreds of eligible Replenish products

CHOOSE your delivery frequency; every week, every month, you're in control

SAVE an additional 3% on all Replenish-eligible products

*See page 155 for complete details or visit acesouthern.com/replenish



INFECTION CONTROL

HALYARD TRI-LAYER OVER-THE-HEAD ISOLATION GOWN WITH THUMB HOOKS AND PERFORATED OPEN BACK

- Protective Tri-Layer SMS Fabric
- Cloth-like, yet tough and abrasion-resistant
- Spunbound polypropylene material
- Micro fibers filter out many water-based liquids
- Tape-tab neck closure
- Yellow



29549	Universal size	10/bag	\$14.29	\$13.58	\$13.58	\$12.86
-------	----------------	--------	---------	---------	---------	---------

HALYARD EVOLUTION 4 STANDARD NON-REINFORCED SURGICAL GOWNS, BOOK FOLD

- Sterile
- EVOLUTION 4 fabric meets AORN standards for fluid, flame and lint resistance
- Raglan sleeves
- Adjustable neckline with yellow neck binding - hook and loop
- Book fold
- With towel
- Blue



90012	Large	36/case	\$149.49	\$142.02	\$142.02	\$134.54
90042	X-Large	34/case	\$149.99	\$142.49	\$142.49	\$134.99

HALYARD TRI-LAYER AAMI2 ISOLATION GOWNS

- Non-sterile
- For facility-wide use when risk of exposure to fluid is between low and moderate
- Medium-weight SMS fabric
- Elastic cuffs
- Tape-tab neck closure
- Blue
- Meets AAMI Level 2 requirements

69981	Universal size	100/case	\$179.29	\$170.33	\$170.33	\$161.36
69987	X-Large	100/case	\$237.99	\$226.09	\$226.09	\$214.19

TIDI HIRISK™ ISOLATION GOWN

Impervious gown made from soft embossed polyethylene which provides a level of protection against infectious materials.

- Latex-free
- Closed loop neck for easy donning and tear off
- Open back for increased air flow
- Thumb loop cuffs for wrist protection
- Passes ASTM F1671 and ASTM F1670 standards
- Light blue
- Impervious material passes AAMI Level 4 protection



8576	One size fits most	75/case	\$123.79	\$117.60	\$117.60	\$111.41
------	--------------------	---------	----------	----------	----------	----------

HALYARD SCRUB SHIRTS

- Protective 3-layer SMS fabric
- Non-sterile
- Single Use
- Blue



69701	Medium	each	\$6.29	\$5.98	\$5.98	\$5.66
69702	Large	each	\$6.29	\$5.98	\$5.98	\$5.66
69703	X-Large	each	\$6.29	\$5.98	\$5.98	\$5.66

MEDLINE STERILE SURGICAL GOWNS

- Sterile, Nonabsorbent
- Individually packaged with a hand towel
- Long sleeves with cuff
- Velcro closure on back collar
- Tie-on waist closure
- Blue



DYNJP2005S	Small	each	\$6.29	\$5.98	\$5.98	\$5.66
DYNJP2001S	Large	each	\$6.29	\$5.98	\$5.98	\$5.66
DYNJP2002S	X-Large	each	\$6.79	\$6.45	\$6.45	\$6.11

HALYARD EVOLUTION 4 SET-IN-SLEEVE NON-REINFORCED SURGICAL GOWN

- Sterile
- Evolution 4 fabric meets AORN standards for fluid, flame and lint resistance
- Set-in-sleeves
- Adjustable neckline with yellow neck binding - hook and loop
- Book fold
- With towel
- Blue



90018	Large	36/case	\$151.99	\$144.39	\$144.39	\$136.79
90048	X-Large	32/case	\$166.29	\$157.98	\$157.98	\$149.66

8

INFECTION CONTROL

HALYARD AERO BLUE STERILE PERFORMANCE SURGICAL GOWN

- Performance Gown with breathable back panel
- Meets ASTM D4966 for abrasion resistance
- Hook-and-Loop Adjustable Neckline
- AERO BLUE meets both the AAMI Level 3 Liquid Barrier Standard and ASTM F1670 standards for resistance of materials used in protective clothing



41732	Small	34/case	\$181.99	\$172.89	\$172.89	\$163.79
41733	Large	32/case	\$183.79	\$174.60	\$174.60	\$165.41
41734	X-Large	30/case	\$185.29	\$176.03	\$176.03	\$166.76

HALYARD ULTRA® NON-REINFORCED SURGICAL GOWNS

- Sterile
- Fluid resistance: barrier protection to prevent blood strikethrough and fluid contamination
- Flame resistance: meets NFPA industry standard for low ignition
- Lint and abrasion resistance: reduces the risk of lint contamination in the wound and related post-operative complications
- Non-reinforced, for low to moderate fluid contact
- Raglan sleeves
- With towel
- Blue
- AAMI Level 3

95101	Small	34/case	\$150.29	\$142.78	\$142.78	\$135.26
95111	Large	32/case	\$209.29	\$198.83	\$198.83	\$188.36
95121	X-Large	30/case	\$215.99	\$205.19	\$205.19	\$194.39
95131	XX-Large	28/case	\$233.29	\$221.63	\$221.63	\$209.96

ALLEGIANCE (BAXTER) OPTIMA® UNREINFORCED BACK GOWN

Non-woven, fluid repellent Optima fabric provides repellency to help prevent strike-through, yet retains the desirable characteristics of linen in the areas of softness, comfort and drapability. Sizes are generous for greater freedom of movement, and polyester cuffs fit comfortably to the wrist.

- Sterile
- With hand towel
- Blue
- AAMI Level 3

9545	X-Large	20/case	\$176.29	\$167.48	\$167.48	\$158.66
------	---------	---------	----------	----------	----------	----------

TIDI® PATIENT GOWN

Tidi gowns are constructed to provide patients modesty and comfort while ensuring convenient access to your patients and hygienic coverage.

- Embossed
- Poly waist ties; stretch to desired length
- Prevents leak through



Tissue/Poly/Tissue, Front or Back Opening

910520	30" x 42", Blue	50/box	\$45.29	\$43.03	\$43.03	\$40.76
910420	30" x 42", White	50/box	\$45.29	\$43.03	\$43.03	\$40.76

CROSSTEX PATIENT BARRIER APRONS

Impervious to liquids to protect the patient. Features two alligator clips attached to the front of the neck band to hold a disposable patient towel/bib and a touch-tab Velcro® neck closure for fast and easy fastening.



- Barrier aprons may be wiped clean quickly and easily with any antibacterial product or warm soap and water.

Plastic/Knee Length

NP4BL	27" X 44", Blue	each	\$14.79	\$14.05	\$14.05	\$13.31
NP4WH	27" X 44", White	each	\$14.79	\$14.05	\$14.05	\$13.31

Plastic/Waist Length

NP3BL	27" X 30", Blue	each	\$14.79	\$14.05	\$14.05	\$13.31
NP3WH	27" X 30", White	each	\$14.79	\$14.05	\$14.05	\$13.31

PETRA MANUFACTURING REUSABLE PLASTIC PATIENT DRAPES

- Practical protection for patient's clothing during surgical procedures, high speed drill operations, prophylaxis, impression techniques, irrigation, extraction
- Neatly tailored neck bond is adjustable for comfort
- Extra features:
 - Tabs perfectly placed on each garment make it convenient to slip on towel in just-right position before patient is in chair.
 - Brushed silk embossed plastic with touch-on (Velcro) neck adjustment.
- Soft, supple, comfortable
- Stains swish away – surfaces stay cleaner

Non-Sterile Ankle-Length Drapes

850BLUE	36" x 54", Blue	each	\$9.49	\$9.02	\$9.02	\$8.54
---------	-----------------	------	--------	--------	--------	--------

Non-Sterile Knee-Length Drapes

PEM851CLEAR	27" x 44", Clear	each	\$9.49	\$9.02	\$9.02	\$8.54
851BLUE	27" x 44", Blue	each	\$9.49	\$9.02	\$9.02	\$8.54

Non-sterile Waist-Length Drapes

852CLEAR	27" x 30", Clear	each	\$9.49	\$9.02	\$9.02	\$8.54
852BLUE	27" x 30", Blue	each	\$9.49	\$9.02	\$9.02	\$8.54

INFECTION CONTROL

STRADIS HEALTHCARE SURGICAL IMPLANT PACKS

- Latex-free

Basic

1 - 24" x24" CSR Wrap	1 - Suture Bag	4 - Sani-Sleeves
2 - Large Gowns	1 - Headrest Cover	1 - 6" Saliva Ejector
2 - Blue Masks	4 - 4" x 6" Barriers	1 - Aspirator Tip
3 - Bouffant Caps	25 - 2" x 2" Sponges	1 - HVE Tip
2 - ABS Towels	8 - 4" x 4" Sponges	2 - 6" CTAs*
1 - 45" x 45" CSR Wrap	1 - OR Towel	3 - Utility Drapes
1 - Fenestrated Drape	1 - X-Ray Cassette Cover	1 - Mayo Cover

SM40561BIP	4/case	\$221.79	\$210.70	\$210.70	\$199.61
------------	--------	----------	----------	----------	----------

Split Drape

1 - 24" x 24" CSR Wrap	1 - Suture Bag	4 - Sani-Sleeves
2 - Large Gowns	1 - Headrest Cover	1 - 6" Saliva Ejector
2 - Blue Masks	4 - 4" x 6" Barriers	1 - Surg-o-Vac II Tip
3 - Bouffant Caps	25 - 2" x 2" Sponges	1 - HVE Tip
2 - ABS Towels	8 - 4" x 4" Sponges	2 - 6" CTAs*
1 - 45" x 45" CSR Wrap	1 - OR Towel	4 - Utility Drapes
1 - 55" x 77" Split Drape w/22" adh strip	1 - X-Ray Cassette Cover	1 - Mayo Cover

40809SIP	4/case	\$193.99	\$184.29	\$184.29	\$174.59
----------	--------	----------	----------	----------	----------

*Cotton tipped applicators

HALYARD SURGICAL PACKS

Come Standard With:

- Back Table Cover Zone-Reinforced, 44" x 90"
- Mayo Stand Cover, Reinforced, 23" x 54"

Also Include Components as Listed Below:

Basic Pack V

ULTRA® Surgical Gown with Towel (in Overwrap), Large • 1 ULTRA® Surgical Gown, X-Large • 4 Utility Drapes, 26" x 15", with Tape • Large Drape, 60" x 76"

88151	8/case	\$154.79	\$147.05	\$147.05	\$139.31
-------	--------	----------	----------	----------	----------

Basic Pack VI

ULTRA® Surgical Gown, with Towel (in Overwrap), Large • 1 ULTRA® Surgical Gown, X-Large • 4 Utility Drapes, 26" x 15", with Tape • Large Drape, 60" x 76" • Half Drape, 60" x 44"

88161	7/case	\$157.49	\$149.62	\$149.62	\$141.74
-------	--------	----------	----------	----------	----------

EENT Pack I

ULTRA® Surgical Gown, with Towel (in Overwrap), Large • Absorbent Towel, 15" x 36" • Head Drape, 50" x 36" • EENT Split Drape, 76" x 124", Suture Bag

88711	5/case	\$126.79	\$120.45	\$120.45	\$114.11
-------	--------	----------	----------	----------	----------

EENT Pack IV

Surgical Gown with Set-In Sleeve with Towel (in Overwrap), Large • Surgical Gown with Set-In Sleeve with Towel, Large • Absorbent Towel, 15" x 36" • Head Drape, 41" x 32" • EENT Split Drape, 76" x 124", Suture Bag

88719	5/case	\$147.29	\$139.93	\$139.93	\$132.56
-------	--------	----------	----------	----------	----------



ACE BASIC DENTAL IMPLANT PACK

Includes:

- 1 - Split Sheet, 78" x 110" w/Adhesive Split
- 2 - Handpiece & Power Cord Drapes, 3" x 96"
- 1 - Head/Bar Drape w/Tape, 27" x 40"
- 1 - ACE Power Unit Drape, 24" x 24"
- 4 - Towels, 18" x 33", 100% Cotton, Slate
- 1 - ACE Light Handle Cover, 7" x 4"
- 1 - Small Gown w/Hand Towel
- 1 - Large Gown w/Hand Towel
- 1 - Half Sheets, 40" x 58"
- 1 - Absorbent Towel, 15" x 23"
- 10 - ACE All-Gauze Sponges, 2" x 2", 12-ply
- 1 - Foam Tray
- 10 - ACE All-Gauze Sponges, 4" x 4", 12-ply
- 1 - Mayo Stand Cover, 23" x 55", Zone Grid w/Absorbent Pad, 18" x 24"

0033144	each	\$77.29	\$73.43	\$73.43	\$69.56
---------	------	---------	---------	---------	---------



ACE BASIC DENTAL IMPLANT PACK I

Includes:

- 1 - Split Sheet, 78" x 110" with Adhesive Split
- 1 - ACE Power Unit Drape, 24" x 24"
- 1 - Head/Bar Drape with Tape, 27" x 40"
- 2 - ACE Light Handle Covers, 7" x 4"
- 4 - Towels, 17" x 24", 100% Cotton, Blue
- 1 - Half Sheets, 40" x 58"
- 10 - ACE All-Gauze Sponges, 2" x 2", 12-ply
- 1 - Absorbent Towel, 15" x 23"
- 10 - ACE All-Gauze Sponges, 4" x 4", 12-ply
- 1 - Foam Tray
- 2 - Handpiece & Power Cord Drapes, 3" x 96"
- 1 - Mayo Stand Cover, 23" x 55", Zone Grid with Absorbent Pad, 18" x 24"

0033143	each	\$64.29	\$61.08	\$61.08	\$57.86
---------	------	---------	---------	---------	---------



ACE DENTAL IMPLANT SURGICAL PACK III

Includes:

- 1 - Foam Tray, 12" x 16"
- 1 - Split Sheet, 72" x 115"
- 4 - Handpiece and Power Cord Drape, 3" x 96"
- 3 - Bouffant Caps
- 2 - Light Handle Covers, 7" x 4"
- 1 - Power Unit Drapes, 24" x 24"
- 1 - Large Gown
- 1 - X-Large Gown
- 1 - Head/Bar Drape
- 1 - Header Bag, 16" x 22"

0033891	each	\$51.79	\$49.20	\$49.20	\$46.61
---------	------	---------	---------	---------	---------



8

INFECTION CONTROL

CARDINAL HEALTH™ CONVERTORS® SURGICAL PACKS

Come Standard With:

- Suture Bag
- Outer Wrap/Reinforced Poly Table Cover, 44" x 90", with 24" x 90" Reinforcement
- 23"W Reinforced Poly Mayo Stand Cover

Also Include Components as Listed Below:

Basic Pack

4 Drape Towels, Adhesive, Non-Absorbent • Fan-Folded Drape Sheet, 57" x 76"

9112	10/case	\$149.49	\$142.02	\$142.02	\$134.54
------	---------	----------	----------	----------	----------

Head and Neck Pack w/Reinforced Drape

Split Sheet Drape, Reinforced, 25" x 2-1/2" Adhesive; Instrument Pad; Cord-Holding Tabs • Head/Bar Drape with Adhesive and Turban • Scrub Nurse Gown with Towel in Outer Wrap, Unreinforced, Large • Absorbent Towel on Gown • Gown, Sterile Back, Unreinforced, Large • 2 Drape Towels, Adhesive

9194	6/case	\$169.49	\$161.02	\$161.02	\$152.54
------	--------	----------	----------	----------	----------

MEDLINE STERILE ABSORBENT TOWELS

- White Embossed

DYNJP2220	15" x 23"	2/pk	\$3.49	\$3.32	\$3.32	\$3.14
-----------	-----------	------	--------	--------	--------	--------

DYNAREX STERILE TOWEL DRAPES

Poly-interlined, three-ply construction - tissue/poly/tissue - provides absorbency while preventing possible moisture or bacteria migration



Plain, White

DY4410	18" x 26"	50/box	\$13.99	\$13.29	\$13.29	\$12.59
--------	-----------	--------	---------	---------	---------	---------

Fenestrated, White

4409	18" x 26"	50/box	\$13.99	\$13.29	\$13.29	\$12.59
------	-----------	--------	---------	---------	---------	---------

3M STERI-DRAPES

Steri-Drape™ Utility Sheet

1089	22" x 25"	2/pk	\$6.29	\$5.98	\$5.98	\$5.66
------	-----------	------	--------	--------	--------	--------


Clear Steri-Drape™ with Adhesive Strips

1000	17" x 11"	10/box	\$30.49	\$28.97	\$28.97	\$27.44
------	-----------	--------	---------	---------	---------	---------

3M1010	17" x 23"	10/box	\$42.99	\$40.84	\$40.84	\$38.69
--------	-----------	--------	---------	---------	---------	---------



ACE STERILE IMPLANT DRAPES WITH OPENINGS



0033134	40" x 52"	1/pk	\$6.29	\$5.98	\$5.98	\$5.66
---------	-----------	------	--------	--------	--------	--------

CARDINAL HEALTH™ ASTOUND® STERILE DRAPE

Half Sheet

9358	40" x 57"	20/case	\$55.79	\$53.00	\$53.00	\$50.21
------	-----------	---------	---------	---------	---------	---------

3M STERILE DRAPES WITH OPENINGS


- 2-1/2" round opening
- Adhesive area: 3-7/8" x 4-7/8"

3M Small Aperture Clear Drape

1020	15" x 15"	10/box	\$43.99	\$41.79	\$41.79	\$39.59
------	-----------	--------	---------	---------	---------	---------

HALYARD STERILE DRAPES WITH OPENINGS

- 3-1/8" circular opening
- Control® Plus Fabric Reinforcement



89471	31" x 44"	50/case	\$183.29	\$174.13	\$174.13	\$164.96
-------	-----------	---------	----------	----------	----------	----------

HALYARD UTILITY DRAPES WITH TAPE

Sterile, Tape Across Long Side

8973135	15" x 26"	2/pk	\$3.79	\$3.60	\$3.60	\$3.41
---------	-----------	------	--------	--------	--------	--------

89331	60" x 70"	each	\$8.49	\$8.07	\$8.07	\$7.64
-------	-----------	------	--------	--------	--------	--------

ALLEGIANCE (BAXTER) STERILE TIBURON® DRAPE SHEET

Square Folded

9461	38-1/2" x 38-1/2"	30/case	\$73.29	\$69.63	\$69.63	\$65.96
------	-------------------	---------	---------	---------	---------	---------

INFECTION CONTROL

MEDLINE STERILE HEAD DRAPES

- Sterile
- Attached utility drape with tape



DYNJP7009	27" x 40"	each	\$7.29	\$6.93	\$6.93	\$6.56
-----------	-----------	------	--------	--------	--------	--------

ALLEGIANCE (BAXTER) CONVERTORS STERILE TIBURON® ORTHOPEDIC SPLIT DRAPE

- IsoBac® Absorbent Reinforcement, 28" x 38"
- Split 37" x 2-1/2" wide

Square Folded

29431	77" x 108"	10/case	\$119.99	\$113.99	\$113.99	\$107.99
-------	------------	---------	----------	----------	----------	----------

HALYARD STERILE SPLIT SHEETS

- No Adhesive on Split
- Reinforced
- Absorbent Pad and Tube Holder

89301	76" x 120"	17/case	\$202.29	\$192.18	\$192.18	\$182.06
-------	------------	---------	----------	----------	----------	----------

MEDLINE STERILE SPLIT SHEETS

- 3-1/2" x 24" Adhesive Split
- Absorbent Pad and Tube Holder

DYNJP7001	76" x 110"	each	\$12.79	\$12.15	\$12.15	\$11.51
-----------	------------	------	---------	---------	---------	---------

HALYARD STERILE MAYO STAND COVERS

Standard Reinforced Flexible Fabric



89601	23" x 54"	54/case	\$189.29	\$179.83	\$179.83	\$170.36
89602	30" x 57"	44/case	\$130.79	\$124.25	\$124.25	\$117.71

TIDI® NON-STERILE EVERYDAY EQUIPMENT DRAPES

Economical single-use products that provide functional everyday protection.

- Tissue/poly material
- Offer strength and fluid holdout
- Embossed for softness
- Blue

980929	40" x 90"	50/case	\$43.99	\$41.79	\$41.79	\$39.59
980926	40" x 60"	100/case	\$57.29	\$54.43	\$54.43	\$51.56

ALLEGIANCE (BAXTER) CONVERTORS® STERILE TABLE COVER

- Reinforced poly

8376	44" x 75"	22/case	\$93.29	\$88.63	\$88.63	\$83.96
8377	44" x 90"	22/case	\$68.29	\$64.88	\$64.88	\$61.46

MEDLINE STERILE TABLE COVERS

44" x 76" with 23" x 55" Absorbent Pad

DYNJP2309		each	\$4.99	\$4.74	\$4.74	\$4.49
-----------	--	------	--------	--------	--------	--------

HALYARD STERILE TABLE COVERS

Back Table Cover, Reinforced, 44" x 90"

42217		each	\$4.99	\$4.74	\$4.74	\$4.49
-------	--	------	--------	--------	--------	--------

ACE HUCK TOWELS

- Non-sterile
- High quality huck towels
- 100% cotton



0032996	Blue, 17" x 24"	12/pk	\$30.49	\$28.97	\$28.97	\$27.44
0032997	White, 17" x 31"	12/pk	\$42.99	\$40.84	\$40.84	\$38.69

MEDLINE STERILE BLUE O.R. TOWELS

- Prewashed and delinted to provide cleanliness and a low particulate count
- Made of highly absorbent 100% cotton



MDT2168282Z		2/pack	\$3.29	\$3.13	\$3.13	\$2.96
MDT2168284Z		4/pack	\$4.99	\$4.74	\$4.74	\$4.49
MDT2168286		6/pack	\$6.79	\$6.45	\$6.45	\$6.11

INFECTION CONTROL

**MEDICAL ACTION ACTI-SORB™
STERILE BLUE O.R. TOWELS**

- 100% cotton fabric for maximum absorbency
- Vacuum delinted
- Convenient packaging contains correct O.R. fold
- Easy open pouch-pack for aseptic delivery



17" x 26"					
701	1/pack	\$1.79	\$1.70	\$1.70	\$1.61
702	2/pack	\$3.29	\$3.13	\$3.13	\$2.96
704	4/pack	\$6.29	\$5.98	\$5.98	\$5.66
706	6/pack	\$8.49	\$8.07	\$8.07	\$7.64

COVIDIEN UNDERPADS

- With polypropylene back
- Fluff fill

1093	23" x 36"	150/case	\$59.29	\$56.33	\$56.33	\$53.36
------	-----------	----------	---------	---------	---------	---------

TIDI UNDERPADS

- Soft nonwoven top ply provides non-irritating skin protection from wetness
- Highly absorbent tissue-filled center
- Waterproof poly backing
- 3-ply tissue-filled
- Blue



16650	17" x 24"	50/bag	\$8.49	\$8.07	\$8.07	\$7.64
16652	23" x 24"	50/bag	\$15.99	\$15.19	\$15.19	\$14.39
16651	24" x 36"	50/bag	\$19.29	\$18.33	\$18.33	\$17.36

**TIDI® EQUIPMENT DRAPE
SHEET/STRETCHER SHEETS**

Tissue/Poly, White

918211	30" x 48"	100/case	\$37.79	\$35.90	\$35.90	\$34.01
--------	-----------	----------	---------	---------	---------	---------

Tissue/Poly, Blue

918213	30" x 48"	100/case	\$37.79	\$35.90	\$35.90	\$34.01
980924	40" x 48"	100/case	\$42.29	\$40.18	\$40.18	\$38.06



TIDI® POLY PATIENT DRAPE

Provides substantial coverage and poly barriers to protect during endodontic and oral surgery procedures.



White						
917670	36" x 40"	100/case	\$35.49	\$33.72	\$33.72	\$31.94

**TIDI® NON-STERILE PATIENT
DRAPE SHEETS**

Tissue/Poly/Tissue, Blue

919374	40" x 48"	50/case	\$27.79	\$26.40	\$26.40	\$25.01
919375	40" x 60"	50/case	\$31.99	\$30.39	\$30.39	\$28.79



ACE STERILE DRAPES

- Individually Wrapped

Half Length						
0038900	40" x 32"	6/bag	\$20.79	\$19.75	\$19.75	\$18.71

Full Length						
0038901	40" x 48"	20/case	\$66.49	\$63.17	\$63.17	\$59.84



**ACE NON-STERILE
PATIENT DRAPES**

- Disposable
- Tie-on strings on the back collar
- White 3-ply tissue/poly backing

Half Length Patient Drape						
0034101	40" x 32"	12/pkg	\$20.79	\$19.75	\$19.75	\$18.71

Full Length Patient Drape						
0034100	40" x 52"	12/pkg	\$20.79	\$19.75	\$19.75	\$18.71



**MEDLINE STERILE FLEXIBLE
LIGHT HANDLE COVERS**

DYNJLHS1	1/pk	80 pks/case	\$58.79	\$55.85	\$55.85	\$52.91
DYNJLHS2	2/pk	160 pks/case	\$106.79	\$101.45	\$101.45	\$96.11



**ASPEN SURGICAL FLEXIBLE LIGHT
HANDLE COVERS**

- Provides a sterile interface between the non-sterile operating room surgical lamps and the sterile surgeon or nurse
- Eliminates costly reprocessing of reusable light handles
- Sterile
- Not made of natural rubber latex

LTF01B	1/pouch, 100 pouches/box	\$101.79	\$96.70	\$96.70	\$91.61
LTF02	2/pouch, 80 pouches/box	\$121.79	\$115.70	\$115.70	\$109.61



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INFECTION CONTROL

PINNACLE CHAIR SLEEVES™

- Tapered design conforms neatly to your chair
- Disposable, so no chair clean up is required
- Eliminates the need for a headrest cover
- Reduces usage of disinfectant and extends the life of your chair



Extra-Long w/Special Gripping Plastic

3850XL	29"W x 80"L	125/box	\$54.49	\$51.77	\$51.77	\$49.04
--------	-------------	---------	---------	---------	---------	---------

Fits A-dec's Cascade and Other Similar Shapes

3825	32"W x 32"L	200/box	\$56.29	\$53.48	\$53.48	\$50.66
------	-------------	---------	---------	---------	---------	---------

Universal

PIN3800	27-1/2"W x 24"L	225/box	\$48.29	\$45.88	\$45.88	\$43.46
---------	-----------------	---------	---------	---------	---------	---------

No-Slip

This unique two-sided dental chair cover is manufactured from two distinctly different plastics laminated together to produce a slippery top surface and a gripping bottom surface.

SNG150	Clear, 48"x 56"	150/box	\$59.99	\$56.99	\$56.99	\$53.99
--------	-----------------	---------	---------	---------	---------	---------

CROSSTEX CHAIR SLEEVE

- Universal clear plastic sleeve for patient chairs
- Covers chair back and chair controls
- Easily slides on and off
- Clear



BCCS	27.5" x 24"	225/dispenser roll	\$35.29	\$33.53	\$33.53	\$31.76
------	-------------	--------------------	---------	---------	---------	---------

DYNAREX CHAIR SLEEVES

- Designed to easily fit over most chairs
- Highly durable
- Easy to use



Half Chair

DY2151	28" x 24"	225/box	\$35.29	\$33.53	\$33.53	\$31.76
--------	-----------	---------	---------	---------	---------	---------

Full Chair

DY2152	29" x 80"	125/box	\$48.79	\$46.35	\$46.35	\$43.91
--------	-----------	---------	---------	---------	---------	---------

UNIPACK HALF CHAIR COVER

27.5" x 24"

UBC8023	Clear	225/box	\$23.99	\$22.79	\$22.79	\$21.59
---------	-------	---------	---------	---------	---------	---------

UNIPACK HEADREST COVERS

Plastic - Clear

UBC8024	11" x 9-1/2" x 2"	250/box	\$7.29	\$6.93	\$6.93	\$6.56
UBC8025	14" x 9-1/2" x 2"	250/box	\$8.79	\$8.35	\$8.35	\$7.91

Tissue/Poly - White

UBC80262	10" x 10"	500/box	\$23.99	\$22.79	\$22.79	\$21.59
UBC80272	10" x 13"	500/box	\$27.29	\$25.93	\$25.93	\$24.56

CROSSTEX HEADREST COVERS

Non-Woven Headrest Covers

- Soft, breathable non-woven fabric
- Quiet - does not rattle
- Provides maximum patient comfort



Standard - 10" x 10"

LOC	White	500/case	\$48.29	\$45.88	\$45.88	\$43.46
-----	-------	----------	---------	---------	---------	---------

Extra Large - 10" x 13"

L3C	White	500/case	\$54.79	\$52.05	\$52.05	\$49.31
-----	-------	----------	---------	---------	---------	---------

Plastic Headrest Covers

- Clear, soft 1 mil thick polyethylene
- 3" gusset - lays flat on headrest
- Impervious barrier between patient and headrest

Regular - 9.5" x 11"

LOCP	Clear	250/disp. box	\$11.79	\$11.20	\$11.20	\$10.61
------	-------	---------------	---------	---------	---------	---------

Jumbo - 9.5" x 14"

L3CP	Clear	250/disp. box	\$13.79	\$13.10	\$13.10	\$12.41
L3CPW	White	250/disp. box	\$16.79	\$15.95	\$15.95	\$15.11

Polycoated Headrest Covers

- Soft facial tissue combined with poly
- Fluid resistant barrier
- Strong and absorbent

Standard - 10" x 10"

LOBL	Blue	500/case	\$42.99	\$40.84	\$40.84	\$38.69
LODR	Dusty rose	500/case	\$42.99	\$40.84	\$40.84	\$38.69
LOLV	Lavender	500/case	\$42.99	\$40.84	\$40.84	\$38.69
LOWH	White	500/case	\$42.99	\$40.84	\$40.84	\$38.69

Extra Large - 10" x 13"

L3BL	Blue	500/case	\$49.79	\$47.30	\$47.30	\$44.81
L3DR	Dusty rose	500/case	\$49.79	\$47.30	\$47.30	\$44.81
L3GR	Green	500/case	\$49.79	\$47.30	\$47.30	\$44.81
L3LV	Lavender	500/case	\$49.79	\$47.30	\$47.30	\$44.81
L3WH	White	500/case	\$49.79	\$47.30	\$47.30	\$44.81

Jumbo - 13" x 13"

L13BL	Blue	500/case	\$67.29	\$63.93	\$63.93	\$60.56
L13WH	White	500/case	\$67.29	\$63.93	\$63.93	\$60.56

Extra Jumbo - 11.5" x 15"

L15WH	White	500/case	\$72.29	\$68.68	\$68.68	\$65.06
-------	-------	----------	---------	---------	---------	---------

INFECTION CONTROL

TIDI® HEADREST COVERS



Fabricel/Poly - White

919511	10" x 10"	500/case	\$56.79	\$53.95	\$53.95	\$51.11
919512	13" x 10"	500/case	\$63.29	\$60.13	\$60.13	\$56.96
919513	13" x 13"	500/case	\$81.99	\$77.89	\$77.89	\$73.79

Tissue/Poly - White

919711	10" x 10"	500/case	\$45.29	\$43.03	\$43.03	\$40.76
919811	13" x 10"	500/case	\$50.29	\$47.78	\$47.78	\$45.26

All-Poly - Clear

919612	10" x 10"	500/case	\$30.79	\$29.25	\$29.25	\$27.71
919622	13" x 10"	500/case	\$42.49	\$40.37	\$40.37	\$38.24

BUFFALO HEADREST COVERS

- Regular with gussets

White

31830	12.5" x 8.5"	500/box	\$42.99	\$40.84	\$40.84	\$38.69
-------	--------------	---------	--------------------	--------------------	--------------------	--------------------

Clear

31820	12.5" x 8.5"	500/box	\$42.99	\$40.84	\$40.84	\$38.69
31850	15.5" x 8.5"	400/box	\$43.99	\$41.79	\$41.79	\$39.59

PINNACLE TRAY SLEEVE™

Instruments are placed directly on top of the disposable plastic sleeve, which is made to replace the present paper tray cover used on instrument trays. Unique fold lock end keeps the tray safely locked inside the sleeve. Minimizes staff time spent on cleaning the tray. Fits mayo trays.



Clear

3500M	14" x 19.5"	500/box	\$50.49	\$47.97	\$47.97	\$45.44
-------	-------------	---------	--------------------	--------------------	--------------------	--------------------

DYNAREX BARRIER FILM

- Apply to almost any surface without worrying about leaving sticky residues behind.
- Provides all-purpose cleanliness by protecting difficult-to-clean areas.
- Roll of 4"x6" perforated sheets - tear off single or multiple sheets as needed.



4488	Clear, 4" x 6"	1200 sheets/roll	\$15.99	\$15.19	\$15.19	\$14.39
DY4489	Blue, 4" x 6"	1200 sheets/roll	\$15.99	\$15.19	\$15.19	\$14.39


UNIPACK EASY-PEEL BARRIER FILM



4" x 6"

UBC8048	Blue	1200/roll	\$12.99	\$12.34	\$12.34	\$12.99
UBC8049	Clear	1200/roll	\$12.99	\$12.34	\$12.34	\$12.99

CROSSTEX SURFACE BARRIERS



- Surface barriers reduce the risk of cross-contamination caused by airborne aerosols, splash and spatter
- Barrier Film with Finger Lift Edge
- Disposable film for areas/surfaces that are hard to reach, clean or disinfect
- Light tack adhesive easily attaches to surfaces without leaving residue
- Unique 3/16" non-adhesive Finger Lift Edge allows quick, easy placement/removal

4" x 6"

BFBL	Blue	1200/roll	\$22.29	\$21.18	\$21.18	\$20.06
BFCL	Clear	1200/roll	\$22.29	\$21.18	\$21.18	\$20.06

UNIPACK TUBE SLEEVE

- Clear poly material to protect instruments from cross-contamination
- Barriers effectively cover equipment with near custom fit
- Pen and syringe sleeves have pre-cut openings for easy application

UBC80394	4" x 1200"	4/case	\$48.49	\$46.07	\$46.07	\$43.64
----------	------------	--------	--------------------	--------------------	--------------------	--------------------

SANI-SLEEVE™

A disposable drape made of tough, thin polyethylene in a unique, U.S. Patented, FDA-accepted design for medical, surgical and dental procedures.



- Sani-Sleeve in PEEL POUCH
- Diameter 2-3/8" x 66"
- Latex-free

SS9335	Sterile	25/box	\$76.29	\$72.48	\$72.48	\$68.66
1001NS	Non-sterile	25/box	\$48.49	\$46.07	\$46.07	\$43.64

ACE STERILE SANI-SLEEVES AND DRAPES



- Ideal for cords and tubing

0033139	Handpiece / Cord 3" x 96"	5/pk	\$71.79	\$71.79	\$64.61	\$64.61
6411000	Sani-Sleeve 2" x 60"	10/pk	\$45.29	\$43.03	\$43.03	\$40.76

Ideal for motors and monitors, etc.

0033138	Power Unit Drape 24" x 24"	5/pk	\$37.29	\$35.43	\$35.43	\$33.56
---------	-------------------------------	------	--------------------	--------------------	--------------------	--------------------

INFECTION CONTROL

UNIPACK UNIVERSAL X-RAY COVER

15" x 26"

UBC8033	Clear	250/box	\$18.49	\$17.57	\$17.57	\$16.64
---------	-------	---------	----------------	----------------	----------------	----------------

CROSSTEX X-RAY SLEEVE

- Universal clear plastic sleeve fits most x-ray machines
- Will not adversely affect x-ray performance
- May also be used as a keyboard cover



15" x 26"

BCXR		250/dispenser box	\$30.49	\$28.97	\$28.97	\$27.44
------	--	-------------------	----------------	----------------	----------------	----------------

UNIPACK SYRINGE SLEEVE WITH OPENING

2.5" x 10"

UBC8028	Clear	500/box	\$8.29	\$7.88	\$7.88	\$7.46
---------	-------	---------	---------------	---------------	---------------	---------------



CROSSTEX SYRINGE SLEEVES

- Clear plastic sleeve easily fits over three-way syringe, high or low speed handpieces, and/or air abrasion handpieces
- Without Opening: Simply place over syringe/handpiece and push through
- With Opening: Convenient, precut opening for easy insertion
- With Opening & Clear Adhesive Strip: Convenient, precut opening for easy insertion with clear adhesive strip that attaches to tubing to prevent sliding of sleeve creating a closed barrier



Without Opening, 2.5" x 10"

BCSS	Clear	500/dispenser box	\$15.29	\$14.53	\$14.53	\$13.76
------	-------	-------------------	----------------	----------------	----------------	----------------

With Opening, 2.5" x 10"

BCSSO	Clear	500/dispenser box	\$15.79	\$15.00	\$15.00	\$14.21
-------	-------	-------------------	----------------	----------------	----------------	----------------

BCSSBO	Blue	500/dispenser box	\$15.79	\$15.00	\$15.00	\$14.21
--------	------	-------------------	----------------	----------------	----------------	----------------

With Opening and Clear Adhesive Strip, 2.5" x 11.5"

BCSSA	Clear	500/dispenser box	\$18.49	\$17.57	\$17.57	\$16.64
-------	-------	-------------------	----------------	----------------	----------------	----------------

CROSSTEX LIGHT SLEEVE

- Patented clear plastic sleeve for T-Style handles
- White indicator arrow for quick location of die-cut
- Will not fall off during use



4" x 5.75"

BCLST	Clear	500/dispenser box	\$20.79	\$19.75	\$19.75	\$18.71
-------	-------	-------------------	----------------	----------------	----------------	----------------

ACE LIGHT HANDLE COVERS



Sterile

0034048	7" x 4"	20/bag	\$24.29	\$23.08	\$23.08	\$21.86
---------	---------	--------	----------------	----------------	----------------	----------------

Non-sterile

0034049	7" x 4", Bulk	100/case	\$39.99	\$37.99	\$37.99	\$35.99
---------	---------------	----------	----------------	----------------	----------------	----------------

DEFEND® LIGHT HANDLE COVERS

- Light Handle Sleeves made of clear plastic. Easy to apply.
- Made to fit most light handles.

BF5000	4" x 5-3/4"	500/box	\$13.99	\$13.29	\$13.29	\$12.59
--------	-------------	---------	----------------	----------------	----------------	----------------

UNIPACK T-HANDLE COVER

4" x 5.75"

UBC8031	Clear	500/box	\$10.49	\$9.97	\$9.97	\$9.44
---------	-------	---------	----------------	---------------	---------------	---------------



CAREFUSION EXIDINE® 4% SCRUB

- Generic Hibiclens®



29900401	*Gallon	each	\$84.49	\$80.27	\$80.27	\$76.04
----------	---------	------	----------------	----------------	----------------	----------------

MOLNLYCKE HEALTHCARE HIBICLENS® SURGICAL SOAP



057516	16 oz	each	\$15.79	\$15.00	\$15.00	\$14.21
057532	32 oz, no pump	each	\$18.29	\$17.38	\$17.38	\$16.46
057591	*Gallon, no pump	each	\$63.29	\$60.13	\$60.13	\$56.96

*A freight charge of \$5.50 will be added to the price of this product.

INFECTION CONTROL

DYNAREX POVIDONE IODINE PREP AND SCRUB SOLUTIONS

- Antiseptic treatment of skin and mucosa
- Surgical site preparation
- Venipuncture
- Care of abrasions, lacerations and minor burns



Prep Scrub

10% Prep Solution

DY1414	8 oz	each	\$1.99	\$1.89	\$1.89	\$1.79
DY1415	16 oz	each	\$3.79	\$3.60	\$3.60	\$3.41
DY1416	Gallon*	each	\$25.29	\$24.03	\$24.03	\$22.76

7.5% Scrub Solution

1425	16 oz	each	\$3.79	\$3.60	\$3.60	\$3.41
------	-------	------	--------	--------	--------	--------

DYNAREX POVIDONE IODINE SWABSTICKS

- Each swabstick is saturated with a 10% povidone iodine solution
- Ideal for venipuncture, IV starts, and pre-op prepping



DY1201	1 Swabstick/pk	50pk/box	\$8.29	\$7.88	\$7.88	\$7.46
DY1202	3 Swabsticks/pk	25pk/box	\$8.29	\$7.88	\$7.88	\$7.46

PURDUE BETADINE® SOLUTION

- Aqueous, film-forming solution of povidone-iodine developed as a unique topical microbicide
- For disinfection of wounds, lacerations and abrasions
- For post-operative applications to incisions to protect against infection



015017	16 oz	each	\$8.49	\$8.07	\$8.07	\$7.64
PU015001	Gallon*	each	\$56.79	\$53.95	\$53.95	\$51.11

DYNAREX POVIDONE IODINE PREP PADS

- Each pad is saturated with 10% povidone iodine solution



DY1108	Medium	100/box	\$7.79	\$7.40	\$7.40	\$7.01
--------	--------	---------	--------	--------	--------	--------

PURDUE BETADINE® SURGICAL SCRUB

- An antiseptic, sudsing skin cleanser
- Produces a rich, golden lather



015117	16 oz	each	\$15.99	\$15.19	\$15.19	\$14.39
015132	32 oz	each	\$21.79	\$20.70	\$20.70	\$19.61
PU015101	Gallon*	each	\$54.79	\$52.05	\$52.05	\$49.31

DYNAREX ALCOHOL PREP PADS

- Sterile alcohol prep pads saturated with 70% isopropyl alcohol



DY1113	Medium	200/box	\$2.99	\$2.84	\$2.84	\$2.69
1116	Large	100/box	\$2.29	\$2.18	\$2.18	\$2.06

COVIDIEN WEBCOL™ ALCOHOL PREP PADS

- 2 ply
- Non-woven pad saturated with 70% isopropyl alcohol

6818	Medium	200/box	\$4.29	\$4.08	\$4.08	\$3.86
5110	Large	200/box	\$5.29	\$5.03	\$5.03	\$4.76

*A freight charge of \$5.50 will be added to the price of this product.

INFECTION CONTROL

3M AVAGARD™

(Chlorhexidine Gluconate 1% Solution and Ethyl Alcohol 61% w/w)

Surgical and Healthcare Personnel Hand Antiseptic with Moisturizers. Fast, effective antimicrobial protection without water or scrubbing.

Persistence and cumulative activity meets and exceeds all FDA criteria at all measured time points for surgical hand scrubs. Fast-acting and effective, Avagard surgical hand prep provides persistent antimicrobial protection for up to six hours without water or scrubbing.

Because Avagard surgical hand prep requires no water or scrub brush, it can be applied in a fraction of the time required to perform a traditional hand scrub*. You spend less time - without compromising efficacy. Application is fast and easy, eliminating the need for water, scrub brushes and sterile towels. Three pumps from a no-touch, simple dispensing system is all it takes.

**Data on file, 3M Healthcare, Instructions for Use.*

Dispenser Bottle - 16.9 oz (500ml)

9200	each	\$75.29	\$71.53	\$71.53	\$75.29
------	------	---------	---------	---------	---------

Wall Bracket and Foot Pump (use with 9200)

9201	each	\$67.99	\$64.59	\$64.59	\$67.99
------	------	---------	---------	---------	---------

Wall Bracket and Hand Pump (use with 9200)

9202	each	\$46.49	\$44.17	\$44.17	\$46.49
------	------	---------	---------	---------	---------



BD® E-Z SCRUB™ PREOPERATIVE SURGICAL SCRUB BRUSHES

- Soft sponge provides gentle cleaning on irritated areas or when firm mechanical action is not required
- Ergonomic brush design allows greater user comfort and effectiveness



With Povidone Iodine

372053	30/box	\$29.79	\$28.30	\$28.30	\$26.81
--------	--------	---------	---------	---------	---------

With 3% Chloroxylenol (PCMX)

371163	30/box	\$32.29	\$30.68	\$30.68	\$29.06
--------	--------	---------	---------	---------	---------

With 4% Chlorhexidine Gluconate

371073	30/box	\$49.99	\$47.49	\$47.49	\$44.99
--------	--------	---------	---------	---------	---------

Dry Surgical Scrub Brush

371603	30/box	\$28.49	\$27.07	\$27.07	\$25.64
--------	--------	---------	---------	---------	---------

GREEN SOAP SOLUTION

- Effective Detergent/Cleaner for Removal of Dried Blood and Protein Soils from skin and scalp
- Effective for Cleaning Surgical Instruments and Surgical Apparatus before Sterilization



1029999	Gallon*	each	\$18.29	\$17.38	\$17.38	\$16.46
---------	---------	------	---------	---------	---------	---------

**A freight charge of \$5.50 will be added to the price of this product. **A freight charge of \$6.00 will be added to the price of this product.*

METREX VIONEXUS® FOAMING SOAP W/ VITAMIN E

A refreshing plumeria-apple fragrance that is pH balanced and hypoallergenic. This soap is non-irritating, leaving hands feeling smooth and soft. The product is dispensed as an economical foam and is recommended for repeated use.



** 1 Liter Pump Bottle

102001	6/case	\$111.49	\$105.92	\$105.92	\$100.34
--------	--------	----------	----------	----------	----------

METREX VIONEX® ANTIMICROBIAL LIQUID SOAP

Intended for handwashing between routine patient contact, procedures, after exposure to potential pathogens and after removal of personal protective equipment. It is formulated to be gentle and non-irritating for hand washing 30-50 times per day.



- Effective germ reduction in 30 seconds
- Lotion based formula for extra softness; Mild and non-drying to skin, contains conditioners, emollients and softeners
- Active ingredient: chloroxylenol (PCMX)

1518	18 oz Pump Bottle	each	\$13.49	\$12.82	\$12.82	\$12.14
------	-------------------	------	---------	---------	---------	---------

1500	* Gallon	each	\$43.29	\$41.13	\$41.13	\$38.96
------	----------	------	---------	---------	---------	---------

CLOROX HEALTHCARE ALOEGUARD® ANTIMICROBIAL SOAP

- Kills pathogenic microorganisms on hands
- Reduces risk of cross-contamination
- Moisturizes and softens skin with aloe vera
- Gentle, soothing, non-irritating
- Pleasant fragrance



7760	18 oz Pump bottle	each	\$12.29	\$11.68	\$11.68	\$11.06
------	-------------------	------	---------	---------	---------	---------

7720	800ml Bag-in-box	each	\$12.99	\$12.34	\$12.34	\$11.69
------	------------------	------	---------	---------	---------	---------

7740	* Gallon	each	\$40.49	\$38.47	\$38.47	\$36.44
------	----------	------	---------	---------	---------	---------

3M AVAGARD™ D INSTANT HAND ANTISEPTIC

- Sensitive to hands; Leaves hands feeling smooth, not sticky or greasy; Fragrance free
- 3M's patented Liquid Crystalline Technology enables greater skin moisturization, counteracting the drying effects of alcohol
- Broad-spectrum activity kills more than 99% of harmful bacteria within 15 seconds†
- Chlorhexidine Gluconate (CHG) compatible
- Effective against MRSA and VRE antibiotic resistant bacteria†

†See www.3m.com/healthcare for additional information



9222	16.9 oz Pump bottle	each	\$11.29	\$10.73	\$10.73	\$10.16
------	---------------------	------	---------	---------	---------	---------

INFECTION CONTROL

CLOROX HEALTHCARE GBG ALOEGEL® HAND SANITIZER

- Kills 99.9% of harmful microorganisms in 15 seconds, including MRSA, VRE, E. coli and Streptococcus faecalis
- Contains moisturizers and skin conditioners including Vitamin E and Aloe Vera
- No sticky residue



7776	18 oz Pump bottle	each	\$12.79	\$12.15	\$12.15	\$11.51
7777	800ml Refill	each	\$12.79	\$12.15	\$12.15	\$11.51

PURELL® TOUCH-FREE DISPENSER AND REFILLS

- Kills 99.99% of the most common germs that cause illness
- Automatically dispenses PURELL, nothing for the user to touch
- Extra-large window to see when it's time to reload
- Up to 30,000 uses from three standard C-size batteries
- Fully ADA compliant, UL/CE registered
- Uses 1,200 ml. refills (sold separately)
- 11¼" H x 6½" W x 4½" D



5908301	Dispenser only	each	\$33.29	\$31.63	\$31.63	\$29.96
539202	1200ml Refill, foam	each	\$46.29	\$43.98	\$43.98	\$41.66
5909129	1200ml Refill, gel	each	\$21.79	\$20.70	\$20.70	\$19.61

CLOROX HEALTHCARE™ HYDROGEN PEROXIDE CLEANER AND DISINFECTANT WIPES

- Proprietary and patented activated hydrogen peroxide technology
- 30-60 second contact time for most bacteria and viruses (see product label for specific kill claims)
- Kill claims for 32 different microorganisms, including germs of high concern such as nine antibiotic resistant bacterial strains (including MRSA, VRE), Norovirus, Rhinovirus, Rotavirus, RSV, Influenza and E. coli
- High-quality tear resistant wipes for greater durability
- No harsh chemical odor or fumes
- Non-corrosive, non-flammable



Wipes						
30824	6.75" x 9"	95/tub	\$7.79	\$7.40	\$7.40	\$7.01
Cleaner						
30828	32 oz Spray	each	\$7.99	\$7.59	\$7.59	\$7.19

METREX CAVIWIPES® AND CAVICIDE® PRODUCTS

- Effective against TB and HBV
- Disinfectant, cleaner and decontaminant, all in one
- Compatible with most medical device material
- Safe for use on non-porous surfaces
- Towels stay fully saturated and won't dry out during use
- Sturdy and durable, non-woven towelette will not bunch up during use



CaviWipes Disinfecting Towelettes

131100	6" x 6.75"	160/tub	\$12.79	\$12.15	\$12.15	\$11.51
131150	9" x 12"	65/tub	\$14.99	\$14.24	\$14.24	\$13.49

CaviCide Surface Disinfectant/Decontaminant Cleaners

131024	24 oz Spray bottle	each	\$11.49	\$10.92	\$10.92	\$10.34
131000	*Gallon	each	\$29.49	\$28.02	\$28.02	\$26.54
131025	**2.5 Gallon	each	\$74.29	\$70.58	\$70.58	\$66.86

METREX CAVIWIPES1® & CAVICIDE1® PRODUCTS

CaviCide1 kills TB and 25 other microorganisms including Norovirus and A. aumannii in just 1 minute and with only 1 step. It is a multi-purpose, non-bleach, low alcohol cleaner and disinfectant designed for everyday use.



- 1 minute kill claim - Increased compliance and decreased risks of cross contamination
- 1-step cleaning and disinfecting - Faster room turnover
- Ready-to-use - No dilution required
- Fragrance-free - Ideal for respiratory sensitive patients and healthcare workers
- Low alcohol formula - Offers material compatibility

CaviWipes1 Disinfecting Towelettes

135100	6" x 6.75"	160/tub	\$13.29	\$12.63	\$12.63	\$11.96
135150	9" x 12"	65/tub	\$15.49	\$14.72	\$14.72	\$13.94

CaviCide1 Surface Disinfectant/Decontaminant Cleaners

135024	24 oz Spray bottle	each	\$11.99	\$11.39	\$11.39	\$10.79
135000	*Gallon	each	\$30.49	\$28.97	\$28.97	\$27.44
135025	**2.5 Gallon	each	\$72.79	\$69.15	\$69.15	\$65.51

*A freight charge of \$5.50 will be added to the price of this product. | **A freight charge of \$6.00 will be added to the price of this product.

INFECTION CONTROL

METREX CAVIWIPES™ 2.0

- **Prepared for Today's and Tomorrow's Pathogens**
Qualifies for the EPA's Emerging Viral Pathogens claim for all virus types (enveloped, large and small non-enveloped viruses) to meet current and potentially future infection prevention needs.
- **2-minute Contact Time**
2-minute universal contact time for bacteria, viruses, and pathogenic fungi.
- **Enhanced Efficacy Claims**
Effective against over 40 pathogens including Norovirus, SARS-CoV-2 and Candida auris for added confidence in preventing healthcare associated infections.
- **Excellent Material Compatibility**
Compatible with hard, non-porous surfaces commonly used in clinical settings.



141100	6" x 6.75"	160/tub	\$13.29	\$12.63	\$12.63	\$11.96
141150	9" x 12"	65/tub	\$15.29	\$14.53	\$14.53	\$13.76

MEDLINE MICRO-KILL+™

Effective Against Listed Organisms in 2 Minutes

Hard-surface, germicidal wipes; disinfecting, deodorizing, cleaning wipes with alcohol. For use on hard, non-porous surfaces only, such as stainless steel, Formica, glass tables, carts, baskets, counters, cabinets and telephones. Effective against: Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella enterica, MRSA, Escherichia coli, VRE, Mycobacterium tuberculosis, Herpes simplex virus 2, Influenza A, HIV-1, Human Hepatitis C Virus, Human Hepatitis B Virus.



MSC351200	6" x 6.75"	160/tub	\$8.49	\$8.07	\$8.07	\$7.64
MSC351210	10" x 10"	65/tub	\$10.49	\$9.97	\$9.97	\$9.44

PDI SANICLOTH® AF3 GERMICIDAL DISPOSABLE WIPES

- Fragrance, alcohol, bleach, acid, phenol, acetone, and ammonia free
- Effective against 44 microorganisms in 3 minutes including the following MDROs, bloodborne pathogens and viruses¹: Klebsiella pneumoniae - NDM-1, Klebsiella pneumoniae - Carbapenem Resistant, ESBL producing E. coli, CA-MRSA USA 300 and 400, RSV, HBV, HCV and HIV
- Fragrance free formulation; no harsh odors
- Perfect for equipment and other hard non-porous surfaces sensitive to alcohol
- Meets CDC, OSHA and CMS Tag F441 guidelines
- Bactericidal, tuberculocidal, virucidal



P13872	6" x 6.75"	160/tub	\$12.99	\$12.34	\$12.34	\$11.69
P63884	7.5" x 15"	65/tub	\$14.29	\$13.58	\$13.58	\$12.86

¹See Technical Data Bulletin for a complete list of microorganisms

PDI SANI-CLOTH® HB GERMICIDAL DISPOSABLE WIPES

- Alcohol free formula (Technical data available from Professional Disposables International at www.pdihc.com)
- An EPA registered hospital-grade disinfectant to accommodate the many situations and hard non-porous environmental surfaces found in healthcare settings



Q08472	6" x 6.75"	160/tub	\$12.79	\$12.15	\$12.15	\$11.51
Q85484	7.5" x 15"	65/tub	\$17.29	\$16.43	\$16.43	\$15.56

PDI SANI-CLOTH® PLUS GERMICIDAL DISPOSABLE WIPES

- Low Alcohol (14.85%) formula
- An EPA registered hospital-grade disinfectant to accommodate the many situations and hard non-porous environmental surfaces found in healthcare settings



Q89072	6" x 6.75"	160/tub	\$12.29	\$11.68	\$11.68	\$11.06
Q85084	7.5" x 15"	65/tub	\$15.29	\$14.53	\$14.53	\$13.76

PDI SUPER SANI-CLOTH® GERMICIDAL DISPOSABLE WIPES

- High alcohol (55%) formula
- An EPA registered hospital-grade disinfectant to accommodate the many situations and hard non-porous environmental surfaces found in healthcare settings



Q55172	6" x 6.75"	160/tub	\$11.29	\$10.73	\$10.73	\$10.16
Q86984	7.5" x 15"	6 tubs of 65/case	\$15.99	\$15.19	\$15.19	\$14.39

PALMERO DISCIDE® ULTRA

DisCide® ULTRA is a hospital-level, one-step, ready-to-use quaternary ammonium, high-level alcohol-based disinfectant that is laboratory-proven to kill deadly pathogens[†] in one minute. It has one of the fastest kill time on the market. DisCide is non-corrosive and non-staining and leaves behind a pleasant herbal scent with no unsightly residue. DisCide is EPA-registered and meets the disinfection requirements of OSHA's Bloodborne Pathogen Standard.



DisCide ULTRA Towlettes

60DIS	6" x 6.75"	160/tub	\$12.29	\$11.68	\$11.68	\$11.06
-------	------------	---------	---------	---------	---------	---------

DisCide ULTRA Disinfectant

3565Q	**Quart sprayer	12/case	\$164.49	\$156.27	\$156.27	\$148.04
3565G	*Gallon	each	\$34.79	\$33.05	\$33.05	\$31.31

[†]TB, MRSA, HIV, (H1N1, HBV, HBC, VRE, RSV, H3N2, influenza A, adenovirus, HSV-2, coronavirus, Pseudomonas aeruginosa, salmonella, staph, E. coli, athlete's foot fungus, among others.

*A freight charge of \$5.50 will be added to the price of this product. | **A freight charge of \$12.00 will be added to the price of this product.

INFECTION CONTROL

**LEGEND BRANDS
BOTANICLEAN® DISINFECTANT**

Use BotaniClean anywhere the hazards of contamination create risk. Featuring a unique, nature-inspired thymol blend, the product is EPA-registered as a Category IV low-toxicity disinfectant, requiring no warning labels and no PPE. This broad-spectrum antimicrobial is economical, ready-to-use and lets you clean, disinfect and deodorize – all in one easy step. And the product doesn't separate, helping to ensure consistent application every time.



BotaniClean has demonstrated effectiveness against viruses similar to SARS-CoV-2 (as the surrogate for COVID-19) under the current EPA Emerging Pathogens Standard, so is included on List N: Disinfectants for Use Against SARS-CoV-2 (listed under EPA no. 34810-25).

- Cleans and disinfects high-touch devices and hard surfaces - won't streak or stain stainless steel, chrome, glass, Plexiglas®, glazed porcelain, painted or varnished surfaces, or plastics
- Deep cleans and deodorizes carpets, upholstery and fabrics
- Mild pleasant fragrance – no harsh fumes
- Ready-to-use AND economical

224006000	*3 Liter bottle	each	\$22.49	\$21.37	\$21.37	\$20.24
-----------	-----------------	------	---------	---------	---------	---------

**METREX ENVIROCID® SURFACE
DISINFECTANT/CLEANER**

A multi-purpose disinfectant intended for use in cleaning, decontaminating and disinfecting hard non-porous, inanimate surfaces and non-critical instruments in hospitals, laboratories, and other critical care areas where environmental control of cross contamination between treated surfaces is important.



133324	*24 oz Spray bottle	each	\$14.29	\$13.58	\$13.58	\$12.86
133300	*Gallon	each	\$41.49	\$39.42	\$39.42	\$37.34

**CERTOL® INTERNATIONAL PROSPRAY™
READY-TO-USE SURFACE
DISINFECTANT/CLEANER**

- Light lemon scent
- Effective against MRSA-CA and HA
- Ready to use format is safe and convenient-no dilution required



Intended Use:

ProSpray is a broad spectrum antimicrobial disinfectant and cleaner. It is tuberculocidal and meets CDC guidelines for intermediate-level disinfectant use when blood and body soils are likely to be present. It may be used on surfaces, as a cleaning solution for instruments and as a disinfectant for dental impressions and appliances.

Hard on Germs-Not on Bare Hands!

So safe it is now EPA approved for use without gloves (in absence of biohazards)

Safe for Surfaces:

ProSpray is a water based disinfectant with non-alcohol actives. Water based disinfectant formulas are preferred by upholstery manufacturers. And, unlike alcohol based formulas, ProSpray will not soften or haze vinyls and plastic. ProSpray will not corrode metals or stainless steel counters as may occur with bleach and hydrogen peroxide based products. It will not stain or bleach surfaces even with repeated use.

PSC240	24 oz Spray bottle	each	\$8.49	\$8.07	\$8.07	\$7.64
PSC128	*Gallon Refill	each	\$23.79	\$22.60	\$22.60	\$21.41

**CETYLITE CETYLICIDE-II® HARD
SURFACE DISINFECTANT**

- EPA-registered, broad-spectrum hospital disinfectant
- Germicidal, fungicidal and virucidal
- Effective against 140+ disease-causing organisms including: HIV, MRSA, VRE, HBV and Influenza A/Brazil Virus (H1N1).**
- 1/2 oz bottle makes 1 quart of ready-to-use solution
- 2 oz of Cetylicide-II Concentrate makes 1 gallon of use solution
- 1 quart of Cetylicide-II Concentrate makes 16 gallons of solution
- No mess dispensing with easy tip-and-pour bottle
- Broad spectrum surface disinfectant/deodorizer
- Complies with OSHA Bloodborne Pathogen Rule
- Pleasant lemon scent
- Economical, up to 80% less expensive than ready-to-use formulas
- Saves inventory space
- No waste - mix only what you need
- Dilutes completely in ordinary tap water



0153	1/2 oz Concentrate	each	\$11.79	\$11.20	\$11.20	\$10.61
0152	One quart concentrate	each	\$128.49	\$122.07	\$122.07	\$115.64

**Influenza A/Brazil Virus (H1N1) - EPA believes, "based on available scientific information, that the currently registered Influenza A virus products will be effective against the 2009-H1N1 flu strain and other Influenza A virus strains on hard, non-porous surfaces" <http://www.epa.gov/oppad001/influenza-disinfectants.html>

**CLOROX CITRACE® HOSPITAL
DISINFECTANT AND DEODORIZER**

- Eliminates odors
- Kills 99.9% of bacteria on hard surfaces in 30 seconds
- Disinfects hard surfaces in 5 minutes or less
- Inhibits the growth of mold and mildew on fabric



HEA49100	14 oz Can	each	\$8.49	\$8.07	\$8.07	\$7.64
----------	-----------	------	--------	--------	--------	--------

LYSOL® I.C. DISINFECTANT SPRAY

- Highly effective hospital disinfectant



58345029	19 oz Can	each	\$15.99	\$15.19	\$15.19	\$14.39
----------	-----------	------	---------	---------	---------	---------

LYSOL® DISINFECTANT FOAM

95524	24 oz Can	each	\$11.79	\$11.20	\$11.20	\$10.61
-------	-----------	------	---------	---------	---------	---------



*A freight charge of \$5.50 will be added to the price of this product.

INFECTION CONTROL

**ASEPSIS INTERNATIONAL ASEPTIFOAM
GERMICIDAL CLEANER AND DEODORANT**



AI201218	18 oz Can	each	\$13.99	\$13.29	\$13.29	\$12.59
----------	-----------	------	----------------	--------------------	--------------------	----------------

MEDLINE BIOHAZARD LINERS

- Color coding
- Biohazard symbols and bilingual markings to help keep materials safety segregated



1.2 mil

NON122424	10 Gallon, 24" x 24"	500/case	\$113.29	\$107.63	\$107.63	\$101.96
-----------	----------------------	----------	-----------------	---------------------	---------------------	-----------------

1.5 mil

NON151717	4 Gallon, 17" x 17"	500/case	\$86.99	\$82.64	\$82.64	\$78.29
-----------	---------------------	----------	----------------	--------------------	--------------------	----------------

NON152426	10 Gallon, 24" x 26"	200/case	\$69.29	\$65.83	\$65.83	\$62.36
-----------	----------------------	----------	----------------	--------------------	--------------------	----------------

**CERTOL PRO TECTOR®
INFECTIOUS WASTE BAGS**

- Bright red waste bags, heavy-duty polyethylene with universal biohazard symbol
- Meets State and OSHA Standards
- Imprinted with required biohazard symbol.
- Meets ASTM D-1709, ASTM D-1922, NY and FL "Red Bag" standards




PW15011	10 Gallon, 24" x 23"	100/box	\$34.29	\$32.58	\$32.58	\$30.86
---------	----------------------	---------	----------------	--------------------	--------------------	----------------

PW15051	12 Gallon, 24" x 30"	100/box	\$43.79	\$41.60	\$41.60	\$39.41
---------	----------------------	---------	----------------	--------------------	--------------------	----------------

**ULINE
STORAGE BINS**

- Polypropylene construction
- Autoclavable




S12414BLU	Small, 7.5" x 4" x 3"	each	\$5.79	\$5.50	\$5.50	\$5.21
-----------	-----------------------	------	---------------	-------------------	-------------------	---------------

S12415BLU	Medium, 11" x 5.5" x 5"	each	\$10.49	\$9.97	\$9.97	\$9.44
-----------	-------------------------	------	----------------	-------------------	-------------------	---------------

S12416BLU	Large, 11" x 8" x 7"	each	\$15.99	\$15.19	\$15.19	\$14.39
-----------	----------------------	------	----------------	--------------------	--------------------	----------------

BECAUSE WE JUST CAN'T FIT EVERYTHING INTO THIS CATALOG!

SHOP ONLINE AT ACESOUTHERN.COM



BE THE OFFICE HERO

Save Time and Money with...



REPLENISH

Our **REPLENISH** program allows you to set up scheduled deliveries for the products your practice uses on a regular basis. Ensure that your practice never runs out of essential items that keep your office running smoothly.

The ACE SOUTHERN Replenish Program, easy as 1-2-3

- 1 SHOP** from hundreds of eligible **Replenish** products
- 2 CHOOSE** your delivery frequency; every week, every month, **you're in control**
- 3 SAVE** an additional 3% on all **Replenish**-eligible products

Easy to Manage

Review, revise or skip orders at your convenience through your online account page. There, you can manage your Replenish items, change the quantities, change the frequency or update products. Prior to each delivery, we will send you an email showing the products that will be shipped. Put the Replenish program to work for you today and see how easy it is to save time and money.



For more details, scan the qr-code, or go to acesouthern.com/replenish



ACE SOUTHERN
Surgical Solutions

800.624.5926 | 800.441.3100

info@acesouthern.com


[ACESOUTHERN.COM](https://www.acesouthern.com) | [in](#) [f](#) [t](#)

INSTRUMENTS

EXTRACTION FORCEPS





MAXILLARY

For more detailed information and images, please visit www.acesouthern.com



MAXILLARY									
6300301	A	#1 Standard, upper anterior	7"	ACE	\$92.99	\$92.99	\$83.69	\$83.69	
1PLUS	B	#1 Upper anterior, PLUS series		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
1PPLUS		#1 Upper anterior, pedo, PLUS series		A-Titan	\$182.79	\$182.79	\$164.51	\$164.51	
F1XS	C	#1 Standard, upper incisors, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
FAF1XS	D	#1 Apical, upper incisors, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
DEF1STD	E	#1 Upper incisors and cuspids		Integra™	\$195.49	\$195.49	\$175.94	\$175.94	
DEF10S	F	#10 Upper molars, universal		Integra™	\$145.79	\$145.79	\$131.21	\$131.21	
V9010S	G	#10 Upper molar, universal - Vantage® econo		Integra™	\$54.49	\$54.49	\$49.04	\$49.04	
6300343	H	#16 Upper molar	7"	ACE	\$92.99	\$92.99	\$83.69	\$83.69	
ATT1700	I	#17 Upper left molar, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
ATT1800	J	#18 Upper right molar, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
6300335	K	#32 Upper molar	7-½"	ACE	\$58.79	\$58.79	\$52.91	\$52.91	
ATT3400	L	#34 Upper anterior, notched beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
ATT3500	M	#35 Upper universal, notched beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
3508	N	#35 Upper universal, beveled beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
F35XS	O	#35 Upper canines & premolars, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
FAF35XS	P	#35 Apical upper canines & premolars, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
ATT4415	Q	#44 Universal root, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
ATT5115	R	#51 Upper root, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
5197	S	#51 Upper root fragment, narrow, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
51PLUS	T	#51 Upper root, PLUS series		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
51PPLUS		#51 Upper root, pedo, PLUS series		A-Titan	\$182.79	\$182.79	\$164.51	\$164.51	
FX51XS	U	#51 Upper roots, serrated, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
6300307	V	#65 Upper root	7-¼"	ACE	\$92.99	\$92.99	\$83.69	\$83.69	
6717	W	#67 Upper 3rd left molar, 3 prong, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
6718	X	#67 Upper 3rd right molar, 3 prong, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
6701	Y	#67 Upper molar, universal, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
F67XS	Z	#67 Apical upper molars, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	

continued on following page



MAXILLARY - continued							
DEF88L	A	#88L 1st & 2nd upper molars, left	Integra™	\$195.49	\$195.49	\$175.94	\$175.94
V9088L		#88L 1st & 2nd upper molars, left, Vantage® econo	Integra™	\$54.49	\$54.49	\$49.04	\$49.04
DEF88R	B	#88R 1st & 2nd upper molars, left	Integra™	\$195.49	\$195.49	\$175.94	\$175.94
V9088R		#88R 1st & 2nd upper molars, left, Vantage® econo	Integra™	\$54.49	\$54.49	\$49.04	\$49.04
6300453	C	#89 Upper Cowhorn	7" ACE	\$115.49	\$115.49	\$103.94	\$103.94
6300452	D	#90 Upper Cowhorn	7" ACE	\$115.49	\$115.49	\$103.94	\$103.94
6300444	E	#150 Upper molar, smooth	7" ACE	\$85.79	\$85.79	\$77.21	\$77.21
6315309	F	#150 Upper molar, serrated	7" ACE	\$115.49	\$115.49	\$103.94	\$103.94
150PLUS	G	#150 Upper universal, PLUS series	A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
150PLUS		#150 Upper universal, pedo, PLUS series	A-Titan	\$182.79	\$182.79	\$164.51	\$164.51
FAF150	H	#150 Apical standard, universal, Presidential®	Hu-Friedy	\$235.29	\$235.29	\$211.76	\$211.76
DEF150	I	#150 Upper incisors, cuspids, bicuspid *and roots; universal	Integra™	\$145.79	\$145.79	\$131.21	\$131.21
V90150		#150 Upper incisors, cuspids, bicuspid and roots; universal, Vantage® econo	Integra™	\$54.49	\$54.49	\$49.04	\$49.04
6300314	J	#150A Upper molar	7" ACE	\$92.99	\$92.99	\$83.69	\$83.69
F150A	K	#150A Cryer parallel beaks, upper incisors, canines & premolars, universal, Presidential®	Hu-Friedy	\$208.79	\$208.79	\$187.91	\$187.91
6315442	L	#150AB Upper molar, split beak, serrated	7" ACE	\$72.79	\$72.79	\$65.51	\$65.51
DEF150AS	M	#150AS Upper anteriors, split beaks, serrated	Integra™	\$195.49	\$195.49	\$175.94	\$175.94
6300468	N	#150AS Upper molar, pedo, serrated	6" ACE	\$72.79	\$72.79	\$65.51	\$65.51
6316473	O	#150S Upper	6" ACE	\$85.79	\$85.79	\$77.21	\$77.21
210PLUS	P	#210 Upper molar, PLUS series	A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
210PLUS		#210 Upper molar, pedo, PLUS series	A-Titan	\$182.79	\$182.79	\$164.51	\$164.51
6300333	Q	#210S Upper 3rd molar	7" ACE	\$92.99	\$92.99	\$83.69	\$83.69
F210S	R	#210S Upper 3rd molar, Presidential® series	Hu-Friedy	\$208.79	\$208.79	\$187.91	\$187.91
0700	S	Upper universal, 1st molar, X-TRAC	A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
2001170	T	Steiglitz root, straight	5-½" ACE	\$68.99	\$68.99	\$62.09	\$62.09
2002171		Steiglitz root, 45 degree angle	5-¼" ACE	\$72.79	\$72.79	\$65.51	\$65.51
2002172	U	Steiglitz root, 90 degree angle	5" ACE	\$72.79	\$72.79	\$65.51	\$65.51

Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INSTRUMENTS

EXTRACTION FORCEPS

MANDIBULAR


For more detailed information and images, please visit www.acesouthern.com



MANDIBULAR									
6300407	A	#13 Lower root forceps	6"	ACE	\$102.29	\$102.29	\$92.06	\$92.06	
F13XS	B	#13 Lower premolar, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
6300344	C	#17 Lower molar forceps	6-3/4"	ACE	\$92.99	\$92.99	\$83.69	\$83.69	
6300345	D	#23 Lower molar forceps	6-1/2"	ACE	\$83.29	\$83.29	\$74.96	\$74.96	
ATT2300	E	#23 Cowhorn lower, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
F23XS	F	#23 Cowhorn lower molars, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
F23	G	#23 Cowhorn, 1st & 2nd lower molar	6-1/2"	Hu-Friedy	\$208.79	\$208.79	\$187.91	\$187.91	
DEF23	H	#23 1st & 2nd lower molar; Cowhorn; universal		Integra™	\$195.49	\$195.49	\$175.94	\$175.94	
V9023	I	#23 1st & 2nd lower molar; Cowhorn; universal, Vantage® economy grade		Integra™	\$54.49	\$54.49	\$49.04	\$49.04	
6300413	J	#33 Lower root forceps	5-1/2"	ACE	\$102.29	\$102.29	\$92.06	\$92.06	
ATT3600	K	#36 Lower anterior, notched beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
3608	L	#36 Lower anterior, beveled beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
F36XS	M	#36 Lower incisors, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
FAF36XS	N	#36 Apical lower premolars/incisors, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
ATT4415	O	#44 Universal root, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
ATT4515	P	#45 Lower root, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
4507	Q	#45 Lower root fragment, narrow, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
FX45XS	R	#45 Lower roots, serrated, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
7301	S	#73 Lower molar, Ash pattern, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
FX73XS	T	#73 Lower molar, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	
FX74N	U	#74N European-style roots, narrow beaks, serrated		Hu-Friedy	\$227.99	\$227.99	\$205.19	\$205.19	
6300417	V	#74N Lower incisor forceps	6"	ACE	\$91.99	\$91.99	\$82.79	\$82.79	
7900	W	#79 Lower molar, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81	
FX79XS	X	#79 Lower molar, Atraumair		Hu-Friedy	\$265.79	\$265.79	\$239.21	\$239.21	

continued on following page

INSTRUMENTS



MANDIBULAR - continued


6300320	A	#151 Lower molar forceps	6-½"	ACE	\$95.49	\$95.49	\$85.94	\$85.94
151PLUS	B	#151 Lower universal, PLUS series		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
151PPLUS	C	#151 Lower universal, pedo, PLUS series		A-Titan	\$182.79	\$182.79	\$164.51	\$164.51
F151	D	#151 Cryer universal, lower incisors, canines, premolars & roots, Presidential®		Hu-Friedy	\$208.79	\$208.79	\$187.91	\$187.91
DEF151	E	#151 Lower incisors, cuspids, bicuspid and roots; universal		Integra™	\$145.79	\$145.79	\$131.21	\$131.21
V90151	F	151 Lower incisors, cuspids, bicuspid and roots; universal, Vantage® economy		Integra™	\$54.49	\$54.49	\$49.04	\$49.04
6300325	G	#151A Lower molar forceps	6-½"	ACE	\$87.99	\$87.99	\$79.19	\$79.19
6300326	H	#151AB Lower anterior forceps, serrated	7"	ACE	\$92.99	\$92.99	\$83.69	\$83.69
6316472	I	#151S Lower molar forceps	6"	ACE	\$85.79	\$85.79	\$77.21	\$77.21
6300349	J	#217 Lower molar forceps	7"	ACE	\$116.49	\$116.49	\$104.84	\$104.84
6300350	K	#222 Lower molar forceps	7"	ACE	\$92.99	\$92.99	\$83.69	\$83.69
7990	L	Lower universal, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
ATT2190	M	Lower universal, notched beaks, X-TRAC		A-Titan	\$229.79	\$229.79	\$206.81	\$206.81
6300427	N	#MD3 Mead lower anterior forceps	5-½"	ACE	\$115.49	\$115.49	\$103.94	\$103.94

X-TRAC® Root Tip Kit

Includes: Upper Root Tip Forceps ATT5115, Lower Root Tip Forceps ATT4515, and Bone Preservation Elevator 63

K4000	O	Lower universal, notched beaks, X-TRAC		A-Titan	\$459.49	\$459.49	\$413.54	\$413.54
-------	---	--	--	---------	----------	----------	----------	----------

INSTRUMENTS



SURGICAL SCISSORS

For more detailed information and images, please visit www.acesouthern.com

DEAN								
4000001	A	#9, smooth	6-1/2"	ACE	\$72.79	\$72.79	\$65.51	\$65.51
4015002		#9, serrated	6-1/2"	ACE	\$66.49	\$66.49	\$59.84	\$59.84
861		Serrated blades	7"	A-TITAN	\$105.99	\$105.99	\$95.39	\$95.39
S9		1 serrated blade, Standard series	6"	Hu-Friedy	\$124.29	\$124.29	\$111.86	\$111.86
5264		1 serrated blade	6-3/4"	Integra™	\$145.79	\$145.79	\$131.21	\$131.21
12100	B	1 serrated blade, angled	7"	J&J Instru.	\$16.99	\$16.99	\$15.29	\$15.29
4005201		#9, serrated, Tungsten carbide	6-3/4"	ACE	\$167.99	\$167.99	\$151.19	\$151.19

GOLDMAN-FOX								
4001028		#16, straight, 1 blade serrated	5-1/2"	ACE	\$44.79	\$44.79	\$40.31	\$40.31
S20SCX	C	Straight, super-cut, Black Line	5-1/4"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
S16SCX		Curved, super-cut, Black Line	5-1/4"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
5320TC	D	1 serrated and 1 CARB-N-SERT blade	5"	Integra™	\$232.79	\$232.79	\$209.51	\$209.51
4002029		#16, curved, 1 blade serrated	5-1/2"	ACE	\$44.79	\$44.79	\$40.31	\$40.31
5SC320	E	Gum scissors, curved	5"	Integra™	\$202.29	\$202.29	\$182.06	\$182.06

IRIS								
4001021		Straight, very delicate	4"	ACE	\$48.79	\$48.79	\$43.91	\$43.91
4001019	F	#17, straight, delicate,	4-1/2"	ACE	\$46.29	\$46.29	\$41.66	\$41.66
S17SCX	G	Straight, delicate, super-cut, Black Line	4-3/4"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
S17		Straight, delicate, Standard series	4-1/2"	Hu-Friedy	\$80.29	\$80.29	\$72.26	\$72.26
869		Straight	4-1/2"	A-TITAN	\$65.79	\$65.79	\$59.21	\$59.21
5300		Straight	4"	Integra™	\$71.29	\$71.29	\$64.16	\$64.16
5304		Straight	4-1/2"	Integra™	\$60.99	\$60.99	\$54.89	\$54.89
V95304	H	Straight, Vantage® economy	4-1/2"	Integra™	\$13.29	\$13.29	\$11.96	\$11.96
221500	I	Straight	4-1/2"	J&J Instru.	\$7.79	\$7.79	\$7.01	\$7.01
4005210		Straight, Tungsten carbide	4-1/2"	ACE	\$123.79	\$123.79	\$111.41	\$111.41
4002022		Curved, very delicate	4"	ACE	\$48.79	\$48.79	\$43.91	\$43.91
4002020		#18, curved, delicate	4-1/2"	ACE	\$46.29	\$46.29	\$41.66	\$41.66
S18SCX		Curved, delicate, super-cut, Black Line	6-3/4"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
ATT840		Curved, extra fine	4-1/2"	A-TITAN	\$79.79	\$79.79	\$71.81	\$71.81
840SC		Curved, extra fine, super cut	4-1/2"	A-TITAN	\$90.99	\$90.99	\$81.89	\$81.89
865SC	J	Curved, super cut	4-1/2"	A-TITAN	\$90.99	\$90.99	\$81.89	\$81.89
865	K	Curved	4-1/2"	A-TITAN	\$65.79	\$65.79	\$59.21	\$59.21
5302		Curved	4"	Integra™	\$74.29	\$74.29	\$66.86	\$66.86
5306		Curved	4-1/2"	Integra™	\$60.99	\$60.99	\$54.89	\$54.89
V95306		Curved, Vantage® economy	4-1/2"	Integra™	\$14.49	\$14.49	\$13.04	\$13.04
221510		Curved	4-1/2"	J&J Instru.	\$8.49	\$8.49	\$7.64	\$7.64
4005211		Curved, Tungsten carbide	4-1/2"	ACE	\$123.79	\$123.79	\$111.41	\$111.41
8111C	L	Curved, carbide tip	4-1/2"	A-TITAN	\$220.79	\$220.79	\$198.71	\$198.71

continued on following page

9



INSTRUMENTS

SURGICAL SCISSORS - continued


For more detailed information and images, please visit www.acesouthern.com

KELLY								
4001003		Straight, smooth	6-¼"	ACE	\$67.99	\$67.99	\$61.19	\$61.19
5250		Straight, sharp/sharp points	6-¼"	Integra™	\$112.29	\$112.29	\$101.06	\$101.06
12215	A	Straight	6-¼"	J&J Instru.	\$10.29	\$10.29	\$9.26	\$9.26
12275		Straight, 1 serrated blade	6-¼"	J&J Instru.	\$10.29	\$10.29	\$9.26	\$9.26
4002005	B	#1 curved, smooth	6-¼"	ACE	\$69.99	\$69.99	\$62.99	\$62.99
4002006		#1 curved, serrated	6-¼"	ACE	\$69.99	\$69.99	\$62.99	\$62.99
859	C	Curved, serrated	6-¼"	A-TITAN	\$90.99	\$90.99	\$81.89	\$81.89
S1		Curved, pointed, 1 serrated, Std. series	6-¼"	Hu-Friedy	\$100.79	\$100.79	\$90.71	\$90.71
5252	D	Curved, sharp/sharp points	6-¼"	Integra™	\$115.99	\$115.99	\$104.39	\$104.39
12235	E	Curved	6-¼"	J&J Instru.	\$10.29	\$10.29	\$9.26	\$9.26
12295		Curved, 1 serrated blade	6-¼"	J&J Instru.	\$10.49	\$10.49	\$9.44	\$9.44
4005203		Curved, Tungsten carbide	6-¼"	ACE	\$167.99	\$167.99	\$151.19	\$151.19
4002007		Angled, smooth	6-¼"	ACE	\$67.99	\$67.99	\$61.19	\$61.19
4002008		Angled, serrated	6-¼"	ACE	\$67.99	\$67.99	\$61.19	\$61.19
LaGRANGE								
4002015		#14, curved, serrated	4-¼"	ACE	\$61.99	\$61.99	\$55.79	\$55.79
877SC	F	Curved, with super cut	4-½"	A-TITAN	\$119.99	\$119.99	\$107.99	\$107.99
877		Curved	4-½"	A-TITAN	\$107.79	\$107.79	\$97.01	\$97.01
221540	G	Curved	4-½"	J&J Instru.	\$9.79	\$9.79	\$8.81	\$8.81
S14SCX		Double-curved, super-cut, Black Line	4-½"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
LAHEY METZENBAUM, METZENBAUM								
4002063		Lahey Metzenbaum, dissecting, curved	5-½"	ACE	\$59.29	\$59.29	\$53.36	\$53.36
S19SCX	H	Metzenbaum super-cut, Black Line	5-½"	Hu-Friedy	\$226.49	\$226.49	\$203.84	\$203.84
LOCKLIN								
4000011	I	#12	6-¼"	ACE	\$67.99	\$67.99	\$61.19	\$61.19
882	J	Open handle, serrated blades	6-¼"	A-TITAN	\$105.99	\$105.99	\$95.39	\$95.39
880	K	Closed handle, serrated blades	6-¼"	A-TITAN	\$105.99	\$105.99	\$95.39	\$95.39

continued on following page



INSTRUMENTS



SURGICAL SCISSORS - continued

For more detailed information and images, please visit www.acesouthern.com

OPERATING									
V952	A	Sharp/sharp points - Vantage® econo	4-½"	Integra™	\$8.79	\$8.79	\$7.91	\$7.91	
V956		Sharp/sharp points - Vantage® econo	5-½"	Integra™	\$9.99	\$9.99	\$8.99	\$8.99	
V958		Sharp/sharp points - Vantage® econo	6-½"	Integra™	\$14.79	\$14.79	\$13.31	\$13.31	
SULLIVAN									
4000009	B	#15 Sullivan, smooth	5-½"	ACE	\$67.99	\$67.99	\$61.19	\$61.19	
SUTURE SCISSORS									
899	C	Suture, straight, with loop	6"	A-TITAN	\$89.29	\$89.29	\$80.36	\$80.36	
809	D	Suture, curved, with loop	4-½"	A-TITAN	\$81.79	\$81.79	\$73.61	\$73.61	
STITCH SCISSORS									
4000030	E	#13 Littauer	6"	ACE	\$67.99	\$67.99	\$61.19	\$61.19	
4000032		#13S Littauer	5"	ACE	\$35.99	\$35.99	\$32.39	\$32.39	
9102	F	Littauer, junior	4-½"	Integra™	\$85.79	\$85.79	\$77.21	\$77.21	
9104		Littauer, standard	5-½"	Integra™	\$92.49	\$92.49	\$83.24	\$83.24	
V99102		Littauer, junior, Vantage® economy	4-½"	Integra™	\$11.29	\$11.29	\$10.16	\$10.16	
V99104		Littauer, standard, Vantage® economy	5-½"	Integra™	\$12.29	\$12.29	\$11.06	\$11.06	
4000033		Sickle stitch	5"	ACE	\$48.79	\$48.79	\$43.91	\$43.91	
MIL9100	G	Spencer, delicate	3-½"	Integra™	\$80.99	\$80.99	\$72.89	\$72.89	
V99100		Spencer, delicate, Vantage® economy	3-½"	Integra™	\$15.79	\$15.79	\$14.21	\$14.21	

MICROSURGICAL SCISSORS

CASTROVIEJO									
707	H	Tissue, curved	4-¾"	A-Titan	\$271.99	\$271.99	\$244.79	\$244.79	
SPVX	I	with super cut, Black Line series	6-½"	Hu-Friedy	\$510.79	\$510.79	\$459.71	\$459.71	
WESTCOTT									
4102042	J	Tenotomy, curved, sharp	5-½"	ACE	\$102.29	\$102.29	\$92.06	\$92.06	
4106032		Straight, sharp	4-¾"	ACE	\$102.29	\$102.29	\$92.06	\$92.06	
4106033		Stitch, curved, sharp	4-¾"	ACE	\$102.29	\$102.29	\$92.06	\$92.06	
4102003		Micro Swedish scissors, curved, sharp	4"	ACE	\$54.49	\$54.49	\$49.04	\$49.04	



MISC. SCISSORS

For more detailed information and images, please visit www.acesouthern.com



BANDAGE AND UTILITY SCISSORS

American Diagnostic Co. Lister Bandage Scissors

301	A	Economy, regular blade	5-1/2"	Amer. Diag	\$2.79	\$2.79	\$2.51	\$2.51
320BK	B	Medicut™ economy, black, razor sharp, hardened	7-1/4"	Amer. Diag	\$5.29	\$5.29	\$4.76	\$4.76
320NP		Medicut™ economy, pink, razor sharp, hardened	7-1/4"	Amer. Diag	\$5.29	\$5.29	\$4.76	\$4.76
320R		Medicut™ economy, red, razor sharp, hardened	7-1/4"	Amer. Diag	\$5.29	\$5.29	\$4.76	\$4.76
320V		Medicut™ economy, purple, razor sharp, hardened	7-1/4"	Amer. Diag	\$5.29	\$5.29	\$4.76	\$4.76



Integra™ Bandage and Utility Scissors

5998	C	Black, needle destroyer, serrated, autoclavable	6-1/2"	Integra™	\$25.49	\$25.49	\$22.94	\$22.94
51000		Black, needle destroyer, serrated, autoclavable	7-1/2"	Integra™	\$25.49	\$25.49	\$22.94	\$22.94

WIRE CUTTERS, SCISSORS AND TWISTERS

Wire Cutters

0805192	D	WCPU Wire cutter, Tungsten carbide	5"	ACE	\$135.79	\$135.79	\$122.21	\$122.21
---------	----------	------------------------------------	----	-----	----------	----------	----------	----------



Wire Cutting scissors

0800211	E	Wire cutting scissors	4-3/4"	ACE	\$53.79	\$53.79	\$48.41	\$48.41
0805212		Tungsten carbide	4-3/4"	ACE	\$144.49	\$144.49	\$130.04	\$130.04
6614052		Economy	4-3/4"	ACE	\$8.99	\$8.99	\$8.09	\$8.09
0800213		Smith, 1 blade serrated	6"	ACE	\$71.79	\$71.79	\$64.61	\$64.61
0805260		Smith, 1 blade serrated, Tungsten carbide	7"	ACE	\$155.99	\$155.99	\$140.39	\$140.39
9124	F	1 serrated blade, angled	4-3/4"	Integra™	\$78.79	\$78.79	\$70.91	\$70.91
V99124		1 serrated blade, angled, Vantage® economy	4-3/4"	Integra™	\$18.29	\$18.29	\$16.46	\$16.46
9D133	G	Collar and Crown, straight	4-1/4"	Integra™	\$71.79	\$71.79	\$64.61	\$64.61
9D135		Collar and Crown, curved	4-1/4"	Integra™	\$73.49	\$73.49	\$66.14	\$66.14



Wire Twisters


0808236	H	Blunt	6"	ACE	\$65.49	\$65.49	\$58.94	\$58.94
0805237		Blunt, Tungsten carbide	6"	ACE	\$117.99	\$117.99	\$106.19	\$106.19

INSTRUMENTS



MICROSURGICAL NEEDLEHOLDERS

For more detailed information and images, please visit www.acesouthern.com



CASTROVIEJO								
3105206	A	Serrated, straight	5-½"	ACE	\$99.29	\$99.29	\$89.36	\$89.36
181828	B	Serrated, straight, with lock	5-½"	Integra™	\$515.29	\$515.29	\$463.76	\$463.76
V9181828		Serrated, straight, Vantage® economy	5-½"	Integra™	\$90.79	\$90.79	\$81.71	\$81.71
3105205	C	Serrated, straight, fine tip, Tungsten carbide	5-½"	ACE	\$221.99	\$221.99	\$199.79	\$199.79
3105110		Serrated, straight, narrow, Tungsten carbide	5-½"	ACE	\$201.79	\$201.79	\$181.61	\$181.61
3105010		Serrated, straight, broad, Tungsten carbide	5-½"	ACE	\$201.79	\$201.79	\$181.61	\$181.61
3105016		Serrated, straight, broad, Tungsten carbide	7"	ACE	\$242.29	\$242.29	\$218.06	\$218.06
3105018		Serrated, curved, broad, Tungsten carbide	7"	ACE	\$242.29	\$242.29	\$218.06	\$218.06
3105112		Serrated, curved, narrow, Tungsten carbide	5-½"	ACE	\$221.99	\$221.99	\$199.79	\$199.79
187TCC	D	Curved round handle, Tungsten carbide	5-½"	A-TITAN	\$393.79	\$393.79	\$354.41	\$354.41
188TC	E	Flat handle, Tungsten carbide	5-½"	A-TITAN	\$393.79	\$393.79	\$354.41	\$354.41

NEEDLEHOLDERS

For more detailed information and images, please visit www.acesouthern.com




BAUMGARTNER								
840		Straight, serrated	5-½"	Integra™	\$103.99	\$103.99	\$93.59	\$93.59
V9840	F	Straight, serrated, Vantage® economy	5-½"	Integra™	\$19.79	\$19.79	\$17.81	\$17.81
B5		Straight, serrated, Tungsten carbide	5-½"	A-TITAN	\$159.49	\$159.49	\$143.54	\$143.54
840TC	G	Straight, serrated, Tungsten carb-n-sert	5-½"	Integra™	\$183.49	\$183.49	\$165.14	\$165.14

CRILE WOOD								
3001012	H	Straight, serrated	6"	ACE	\$56.29	\$56.29	\$50.66	\$50.66
256420		Straight, serrated	6"	J&J Instru.	\$12.49	\$12.49	\$11.24	\$11.24
MIL850		Straight, serrated	6"	Integra™	\$87.79	\$87.79	\$79.01	\$79.01
V9850		Straight, serrated, Vantage® economy	6"	Integra™	\$18.79	\$18.79	\$16.91	\$16.91
3005117		Straight, Tungsten carbide	6"	ACE	\$99.29	\$99.29	\$89.36	\$89.36
C6C		Straight, serrated, Tungsten carbide	6"	A-TITAN	\$159.49	\$159.49	\$143.54	\$143.54
850TC	I	Straight, serrated, Tungsten carbide	6"	Integra™	\$186.99	\$186.99	\$168.29	\$168.29
3005118		Straight, serrated, Tungsten carbide	7"	ACE	\$99.29	\$99.29	\$89.36	\$89.36
3002016		Curved, serrated	5-½"	ACE	\$87.99	\$87.99	\$79.19	\$79.19

DEBAKEY								
3215010	J	Straight, serrated, Tungsten carbide	7"	ACE	\$115.49	\$115.49	\$103.94	\$103.94

continued on following page




INSTRUMENTS

NEEDLEHOLDERS - continued

For more detailed information and images, please visit www.acesouthern.com


HALSEY								
3015001	A	Straight, serrated	5-1/4"	ACE	\$49.99	\$49.99	\$44.99	\$44.99
3005100		Straight, serrated, Tungsten carbide	5"	ACE	\$99.79	\$99.79	\$89.81	\$89.81
MAYO HEGAR								
842	B	Straight, serrated	5"	Integra™	\$77.49	\$77.49	\$69.74	\$69.74
256290		Straight, serrated	5"	J&J Intsr.	\$10.99	\$10.99	\$9.89	\$9.89
3001018		Straight, serrated	6"	ACE	\$49.99	\$49.99	\$44.99	\$44.99
844		Straight, serrated	6"	Integra™	\$81.79	\$81.79	\$73.61	\$73.61
V9844	C	Straight, serrated, Vantage® economy	6"	Integra™	\$19.79	\$19.79	\$17.81	\$17.81
256360		Straight, serrated	6"	J&J Intsr.	\$12.49	\$12.49	\$11.24	\$11.24
3001019		Straight, serrated	7"	ACE	\$56.29	\$56.29	\$50.66	\$50.66
846		Straight, serrated	7"	Integra™	\$91.29	\$91.29	\$82.16	\$82.16
V9846		Straight, serrated, Vantage® economy	7"	Integra™	\$22.79	\$22.79	\$20.51	\$20.51
256370	D	Straight, serrated	7"	J&J Intsr.	\$15.49	\$15.49	\$13.94	\$13.94
3001020		Straight, serrated	8"	ACE	\$56.29	\$56.29	\$50.66	\$50.66
842TC		Straight, serrated, Tungsten carbide	5"	Integra™	\$182.99	\$182.99	\$164.69	\$164.69
M6	E	Straight, serrated, Tungsten carbide	6"	A-TITAN	\$159.49	\$159.49	\$143.54	\$143.54
844TC		Straight, serrated, Tungsten carbide	6"	Integra™	\$185.29	\$185.29	\$166.76	\$166.76
3005109	F	Straight, serrated, Tungsten carbide	6 1/4"	ACE	\$99.29	\$99.29	\$89.36	\$89.36
3005110		Straight, serrated, Tungsten carbide	7"	ACE	\$128.79	\$128.79	\$115.91	\$115.91
M7	K	Straight, serrated, Tungsten carbide	7"	A-TITAN	\$164.29	\$164.29	\$147.86	\$147.86
846TC		Straight, serrated, Tungsten carbide	7"	Integra™	\$202.79	\$202.79	\$182.51	\$182.51
MICRO RYDER (SMALLER TIP)								
M6C	G	Straight, serrated, Tungsten carbide	6"	A-TITAN	\$199.79	\$199.79	\$179.81	\$179.81
OLSEN HEGAR (BUILT-IN SCISSORS)								
O5C	H	Straight, serrated, Tungsten carbide	5-1/2"	A-TITAN	\$189.79	\$189.79	\$170.81	\$170.81
SWEDISH MICRO (SMALLER TIP)								
S6	I	Straight, serrated, Tungsten carbideslightly thinner beaks than a Crilewood	6"	A-TITAN	\$68.49	\$68.49	\$61.64	\$61.64
WEBSTER								
3005103		Serrated, Tungsten carbide	4-3/4"	ACE	\$95.49	\$95.49	\$85.94	\$85.94
86	J	Smooth jaws, standard pattern, delicate	5"	Integra™	\$97.79	\$97.79	\$88.01	\$88.01



INSTRUMENTS

PLIERS AND FORCEPS

For more detailed information and images, please visit www.acesouthern.com



CORN SUTURE PLIERS

2000107	A	Angled	6"	ACE	\$58.29	\$58.29	\$52.46	\$52.46
106		Angled	6"	A-Titan	\$83.29	\$83.29	\$74.96	\$74.96
70153	B	Angled	6"	Integra™	\$57.49	\$57.49	\$51.74	\$51.74

DRESSING FORCEPS/PLIERS

66317		Cotton, serrated 14-14.5mm jaw		Integra™	\$31.29	\$31.29	\$28.16	\$28.16
66317S		Cotton, smooth, 14-14.5mm jaw		Integra™	\$30.29	\$30.29	\$27.26	\$27.26
66317A		Cotton, serrated, extra-fine, 11.5-12mm jaw		Integra™	\$32.29	\$32.29	\$29.06	\$29.06
08200		College, serrated	6"	J&J Instru.	\$5.79	\$5.79	\$5.21	\$5.21
08205	C	College, serrated, narrow tips	6"	J&J Instru.	\$5.99	\$5.99	\$5.39	\$5.39
08210	D	College, serrated, self-locking	6"	J&J Instru.	\$6.79	\$6.79	\$6.11	\$6.11
2002145	E	#17 Dressing	6-¼"	ACE	\$20.79	\$20.79	\$18.71	\$18.71
8817034		#17 Dressing, angular, economy	6-¼"	ACE	\$7.79	\$7.79	\$7.01	\$7.01
DPU17		U17 utility pick-up dressing	6-¼"	Integra™	\$18.99	\$18.99	\$17.09	\$17.09
2002146	F	#18 Dressing	6-½"	ACE	\$21.99	\$21.99	\$19.79	\$19.79
8815005		Dressing, economy	6"	ACE	\$5.79	\$5.79	\$5.21	\$5.21
66	G	Dressing, standard pattern, serrated	5"	Integra™	\$25.99	\$25.99	\$23.39	\$23.39
692		Dressing, delicate pattern, 1 x 2 teeth	5"	Integra™	\$37.79	\$37.79	\$34.01	\$34.01
696	H	Dressing, delicate pattern, 1 x 2 teeth	6"	Integra™	\$45.49	\$45.49	\$40.94	\$40.94

TISSUE FORCEPS

6118	I	Adson, delicate, serrated	4-¾"	Integra™	\$47.99	\$47.99	\$43.19	\$43.19
V96118		Adson, delicate, serrated, Vantage® economy	4-¾"	Integra™	\$12.79	\$12.79	\$11.51	\$11.51
2015115	J	Adson #41, serrated	4-¾"	ACE	\$21.99	\$21.99	\$19.79	\$19.79
485		Adson, serrated	4-¾"	A-Titan	\$56.29	\$56.29	\$50.66	\$50.66
6121		Adson, cross serrated	4-¾"	Integra™	\$59.79	\$59.79	\$53.81	\$53.81
V96120		Adson, cross serrated, Vantage® economy	4-¾"	Integra™	\$12.79	\$12.79	\$11.51	\$11.51
2000118		Adson #42, 1x2 teeth	4-¾"	ACE	\$28.49	\$28.49	\$25.64	\$25.64
486		Adson, 1x2 teeth, serrated	4-¾"	A-Titan	\$56.29	\$56.29	\$50.66	\$50.66
6120A		Adson, 1x2 teeth, delicate pattern, angled	4-¾"	Integra™	\$65.29	\$65.29	\$58.76	\$58.76
2000119		Adson, 1x2 teeth	6"	ACE	\$30.49	\$30.49	\$27.44	\$27.44
2000127		Adson, 1x2 teeth, thumb grips	4-¾"	ACE	\$44.79	\$44.79	\$40.31	\$40.31
355020		Adson, 1x2 teeth		J&J Instru.	\$6.29	\$6.29	\$5.66	\$5.66
642		Adson, 1x2 teeth	5"	Integra™	\$27.79	\$27.79	\$25.01	\$25.01
V9646		Adson, 1x2 teeth, Vantage® economy	6"	Integra™	\$10.79	\$10.79	\$9.71	\$9.71
2100038		Micro Adson, 1x2 teeth .8mm	6"	ACE	\$39.49	\$39.49	\$35.54	\$35.54
2115037		Micro Adson, 1x2 teeth, straight	4-¾"	ACE	\$39.49	\$39.49	\$35.54	\$35.54
2115035		Micro Adson, straight, serrated	4-¾"	ACE	\$39.49	\$39.49	\$35.54	\$35.54

continued on following page

INSTRUMENTS

PLIERS AND FORCEPS - continued

For more detailed information and images, please visit www.acesouthern.com

TISSUE FORCEPS (continued)								
168	A	Allis tissue forceps, 5x6 teeth	6"	Integra™	\$85.29	\$85.29	\$76.76	\$76.76
243050		Allis tissue forceps, 5x6 teeth	6"	J&J Instru.	\$16.49	\$16.49	\$14.84	\$14.84
6124	B	Brown-Adson, 7 x 7, Side Grasping Teeth		Integra™	\$70.99	\$70.99	\$63.89	\$63.89
V96124		Brown-Adson, 7 x 7, Side Grasping Teeth, Vantage® economy		Integra™	\$16.79	\$16.79	\$15.11	\$15.11
355090		Brown-Adson, 7 x 7 teeth		J&J Instru.	\$8.79	\$8.79	\$7.91	\$7.91
2000048	C	Debakey, 2mm wide	6"	ACE	\$85.79	\$85.79	\$77.21	\$77.21
2000049		Debakey, 2mm wide	7-¾"	ACE	\$85.79	\$85.79	\$77.21	\$77.21
24570XL	D	Debakey, atraumatic, fenestrated handles	6"	Integra™	\$163.79	\$163.79	\$147.41	\$147.41
2000104		Gerald, 1x2 teeth, straight	7"	ACE	\$41.49	\$41.49	\$37.34	\$37.34
6182		Gerald, 1x2 teeth, straight	7"	Integra™	\$82.99	\$82.99	\$74.69	\$74.69
2002105	E	Gerald, 1x2 teeth, curved	7"	ACE	\$41.49	\$41.49	\$37.34	\$37.34
2102022	F	Iris, delicate, half curve, serrated	4"	ACE	\$29.49	\$29.49	\$26.54	\$26.54
6104	G	Semkin, serrated	5"	Integra™	\$36.99	\$36.99	\$33.29	\$33.29
V96104		Semkin, serrated, Vantage® economy	5"	Integra™	\$14.49	\$14.49	\$13.04	\$13.04
6109		Semkin, serrated	6"	Integra™	\$45.29	\$45.29	\$40.76	\$40.76
6109A	H	Semkin, 1x2 teeth	6"	Integra™	\$51.99	\$51.99	\$46.79	\$46.79

TISSUE FORCEPS - RUSSIAN								
2001108	I	Russian tissue forceps	6"	ACE	\$40.49	\$40.49	\$36.44	\$36.44
6142		Russian tissue forceps	6"	Integra™	\$57.49	\$57.49	\$51.74	\$51.74
2001109		Russian tissue forceps	8"	ACE	\$41.49	\$41.49	\$37.34	\$37.34
6144		Russian tissue forceps	8"	Integra™	\$67.79	\$67.79	\$61.01	\$61.01

MOSQUITO FORCEPS								
2001201	J	#4 Halsted, straight	4-¾"	ACE	\$40.49	\$40.49	\$36.44	\$36.44
2002202		#3 Halsted, curved	4-¾"	ACE	\$42.49	\$42.49	\$38.24	\$38.24
74		Halsted, curved	5"	Integra™	\$60.99	\$60.99	\$54.89	\$54.89
V974		Halsted, curved, Vantage®, economy	5"	Integra™	\$14.49	\$14.49	\$13.04	\$13.04
726	K	Hartman, curved	3-½"	Integra™	\$72.99	\$72.99	\$65.69	\$65.69
V9726		Hartman, curved, economy	3-½"	Integra™	\$14.29	\$14.29	\$12.86	\$12.86
738	L	Kelly forceps, curved	5-½"	Integra™	\$61.29	\$61.29	\$55.16	\$55.16
V9738		Kelly forceps, curved, Vantage®, economy	5-½"	Integra™	\$15.49	\$15.49	\$13.94	\$13.94
2001220	M	Providence Kelly, straight	5-½"	ACE	\$40.49	\$40.49	\$36.44	\$36.44
2002221		Providence Kelly, curved	5-½"	ACE	\$39.49	\$39.49	\$35.54	\$35.54
2002225	N	Crile, curved	5-½"	ACE	\$28.49	\$28.49	\$25.64	\$25.64
241245		Crile, curved	5-½"	J&J Instru.	\$10.29	\$10.29	\$9.26	\$9.26

continued on following page



9

INSTRUMENTS

PLIERS AND FORCEPS - continued

For more detailed information and images, please visit www.acesouthern.com



9

WYNNMAN CROWN GRIPPER

6037210		Curved, Tungsten carb-n-sert	6-½"	ACE	\$126.29	\$126.29	\$113.66	\$113.66
7210TC	A	Curved, Tungsten carb-n-sert	5-½"	Integra™	\$248.99	\$248.99	\$224.09	\$224.09

TOWEL FORCEPS/CLAMPS

2000305	B	Backhaus towel forceps	3-½"	ACE	\$20.79	\$20.79	\$18.71	\$18.71
2000306		Backhaus towel forceps	5"	ACE	\$21.99	\$21.99	\$19.79	\$19.79
4413031		Backhaus Roeder w/ball stops, economy	5-¼"	ACE	\$9.49	\$9.49	\$8.54	\$8.54
7504		Backhaus towel clamps	3-½"	Integra™	\$54.29	\$54.29	\$48.86	\$48.86
7506		Backhaus towel clamps	5-½"	Integra™	\$54.99	\$54.99	\$49.49	\$49.49
V97504		Backhaus towel clamps, Vantage®, economy	3-½"	Integra™	\$14.79	\$14.79	\$13.31	\$13.31
V97506		Backhaus towel clamps, Vantage®, economy	5-½"	Integra™	\$15.29	\$15.29	\$13.76	\$13.76
2000309	C	Jones Towel Forceps	3-½"	ACE	\$25.99	\$25.99	\$23.39	\$23.39

ARTICULATING PAPER FORCEPS

763	D	Holds strips straight, economy grade	6"	Integra™	\$13.79	\$13.79	\$12.41	\$12.41
-----	---	--------------------------------------	----	----------	---------	---------	---------	---------

MAGILL FORCEPS

4419051	E	Magill forceps, child	7-½"	ACE	\$14.49	\$14.49	\$13.04	\$13.04
315	F	Magill forceps, child	8"	Amer. Diag.	\$9.79	\$9.79	\$8.81	\$8.81
4423050		Magill forceps, adult	9-¾"	ACE	\$19.79	\$19.79	\$17.81	\$17.81
316	G	Magill forceps, adult	9-¾"	Amer. Diag.	\$9.79	\$9.79	\$8.81	\$8.81

STERILIZING FORCEPS

2001320	H	Universal sterilizing forceps, straight	8"	ACE	\$29.49	\$29.49	\$26.54	\$26.54
2000301		3-prong sterilizing forceps	7-¾"	ACE	\$44.79	\$44.79	\$40.31	\$40.31

YOUNG TONGUE SEIZING FORCEPS

0100151	I	with rubber jaws	6-¼"	ACE	\$66.49	\$66.49	\$59.84	\$59.84
0100152		Replacement tips for 0100151 (latex-free)	2/set	ACE	\$12.29	\$12.29	\$11.06	\$11.06
317		with rubber jaws	6-¼"	Amer. Diag.	\$9.29	\$9.29	\$8.36	\$8.36

INSTRUMENTS

RONGEURS

For more detailed information and images, please visit www.acesouthern.com

#1 RONGEUR								
1000001	A	#1 Rongeur	6-½"	ACE	\$147.99	\$147.99	\$133.19	\$133.19

BEYER RONGEUR								
1000043	B	Beyer, double action	7"	ACE	\$262.79	\$262.79	\$236.51	\$236.51

BLUMENTHAL RONGEUR								
1000009	C	Standard Blumenthal, 30° Angle	6"	ACE	\$160.49	\$160.49	\$144.44	\$144.44
ATT300		Blumenthal	5-½"	A-Titan	\$266.29	\$266.29	\$239.66	\$239.66
22480		Blumenthal, 30°	6"	Integra™	\$355.79	\$355.79	\$320.21	\$320.21
11030		Blumenthal, 30°	6"	J&J Instru.	\$43.29	\$43.29	\$38.96	\$38.96
11045	D	Blumenthal, 45°	6"	J&J Instru.	\$43.29	\$43.29	\$38.96	\$38.96
22D422		Mini Blumenthal, 25°	5-¼"	Integra™	\$327.79	\$327.79	\$295.01	\$295.01
11020		Mini Blumenthal, 30°	4-½"	J&J Instru.	\$43.29	\$43.29	\$38.96	\$38.96


CLEVELAND RONGEUR								
1000021	E	#4 Cleveland	5-½"	ACE	\$147.99	\$147.99	\$133.19	\$133.19

FRIEDMAN RONGEUR								
1000034	F	Friedman	5-½"	ACE	\$160.49	\$160.49	\$144.44	\$144.44
ATT310		Friedman	5-¼"	A-Titan	\$266.29	\$266.29	\$239.66	\$239.66
22494		Friedman, very delicate	5-½"	Integra™	\$370.29	\$370.29	\$333.26	\$333.26
11110		Friedman, 30°	5-½"	J&J Instru.	\$43.29	\$43.29	\$38.96	\$38.96
1000035	G	Mini Friedman	5"	ACE	\$165.49	\$165.49	\$148.94	\$148.94
11115		Micro Friedman	5-½"	J&J Instru.	\$43.29	\$43.29	\$38.96	\$38.96
1010036		Micro Friedman	5-½"	ACE	\$179.99	\$179.99	\$161.99	\$161.99
310M		Micro Friedman	6"	A-Titan	\$266.29	\$266.29	\$239.66	\$239.66
174801		Micro Friedman, curved jaws, 15°	5-½"	Integra™	\$411.49	\$411.49	\$370.34	\$370.34
17D4802	H	Micro Friedman, 90°	5-½"	Integra™	\$395.79	\$395.79	\$356.21	\$356.21

KERRISON RONGEUR								
1000075	I	Kerrison Rongeur, 3.5mm Bite, 4" Shaft	5-¾"	ACE	\$533.79	\$533.79	\$480.41	\$480.41

LUER RONGEUR								
1000018		Luer Rongeur, 4mm Cup	6"	ACE	\$147.99	\$147.99	\$133.19	\$133.19

UNIVERSAL SIDE CUTTING RONGEUR								
1000025	J	Universal Side Cutting Rongeur	6"	ACE	\$151.99	\$151.99	\$136.79	\$136.79



INSTRUMENTS

SURGICAL MALLETS

ACE Surgical Mallet

Stainless steel

6400312S	6-½"	each	\$81.29	\$81.29	\$73.16	\$73.16
----------	------	------	---------	---------	---------	---------



ACE Mead Mallet

with 2 interchangeable plastic facings

6400307	8"	each	\$102.29	\$102.29	\$92.06	\$92.06
---------	----	------	----------	----------	---------	---------



6400305	Replacement facing	each	\$15.49	\$15.49	\$13.94	\$13.94
---------	--------------------	------	---------	---------	---------	---------

ACE Mead Mallet

with 2 interchangeable plastic discs

6400308	7-½"	each	\$80.29	\$80.29	\$72.26	\$72.26
---------	------	------	---------	---------	---------	---------



6400310	Replacement facing	2/set	\$32.79	\$32.79	\$29.51	\$29.51
---------	--------------------	-------	---------	---------	---------	---------

Integra™ Mead Mallet

with 2 interchangeable nylon caps

19796	7-½"	each	\$177.49	\$177.49	\$159.74	\$159.74
-------	------	------	----------	----------	----------	----------



PERIODONTAL BONE CHISELS

Hu-Friedy® Black Line - Buser Modified

PPBUSER12X	1/2 Buser, 3mm/4mm	each	\$89.79	\$89.79	\$80.81	\$80.81
PPBUSER23X	2/3 Buser, 4mm/5mm	each	\$89.79	\$89.79	\$80.81	\$80.81
PPBUSER13X	1/3 Buser, 3mm/5mm	each	\$89.79	\$89.79	\$80.81	\$80.81

CHISELS

ACE Cottle Nasal Chisel

0502211	4mm wide, 6-¼"	each	\$67.99	\$67.99	\$61.19	\$61.19
---------	----------------	------	---------	---------	---------	---------

0502017	6mm wide, 6-¼"	each	\$67.99	\$67.99	\$61.19	\$61.19
---------	----------------	------	---------	---------	---------	---------



ACE Flat Chisel

0501220	45° Angle, 6"	each	\$67.99	\$67.99	\$61.19	\$61.19
---------	---------------	------	---------	---------	---------	---------



BONE FILES

ACE #12 Howard Bone File

6400004	7", 7.6 and 8mm wide	each	\$39.49	\$39.49	\$35.54	\$35.54
---------	----------------------	------	---------	---------	---------	---------



Integra™ #12 Bone File

401381	#12	each	\$83.49	\$83.49	\$75.14	\$75.14
--------	-----	------	---------	---------	---------	---------

401382	#12A		\$83.49	\$83.49	\$75.14	\$75.14
--------	------	--	---------	---------	---------	---------



J&J Instru. #12 Howard Bone File

02120	#12	each	\$14.29	\$14.29	\$12.86	\$12.86
-------	-----	------	---------	---------	---------	---------

02100	#12A, cross serrations	each	\$14.29	\$14.29	\$12.86	\$12.86
-------	------------------------	------	---------	---------	---------	---------



ACE. #21 Miller Bone File

6400009	7-¼"	each	\$39.49	\$39.49	\$35.54	\$35.54
---------	------	------	---------	---------	---------	---------



Hu-Friedy® #21 Miller Bone File

Straight-Cut Serrations

BF21	#522 handle	each	\$82.29	\$82.29	\$74.06	\$74.06
------	-------------	------	---------	---------	---------	---------




J&J Instru. #21 Miller Bone File

02210		each	\$14.29	\$14.29	\$12.86	\$12.86
-------	--	------	---------	---------	---------	---------



ACE #45 Miller Bone File

6400003	7"	each	\$39.49	\$39.49	\$35.54	\$35.54
---------	----	------	---------	---------	---------	---------



ACE #52 Miller Bone File


6400002	7"	each	\$39.49	\$39.49	\$35.54	\$35.54
---------	----	------	---------	---------	---------	---------



INSTRUMENTS

ACE BONE SPREADER OSTEOTOME KIT

- Single end handle can be tapped with a mallet.
- Markings at 13mm, 16mm and 19mm.
- D-contoured tips to avoid fracturing weaker boney wall.
- Instruments and tray are Stainless Steel.



4542106	A	Kit with mallet and tray	kit	\$463.29	\$463.29	\$416.96	\$416.96
---------	----------	--------------------------	-----	----------	----------	----------	----------

Bone Spreader Osteotomes and Accessories available separately

4542102	B	4.3mm 27mm scoop	each	\$102.79	\$102.79	\$92.51	\$92.51
4542103	C	5.45mm 26mm Scoop	each	\$102.79	\$102.79	\$92.51	\$92.51
4542104	D	7.5mm 26mm Scoop	each	\$102.79	\$102.79	\$92.51	\$92.51
4542105	D	10.0mm 27mm Scoop	each	\$102.79	\$102.79	\$92.51	\$92.51
6400312S		Surgical mallet, 6-½" (16.5cm)	each	\$81.29	\$81.29	\$73.16	\$73.16
9400036		Instrument tray, 175 x 140 x 24mm	each	\$135.79	\$135.79	\$122.21	\$122.21



Especially designed for taking Chin and Ramus monocortical block grafts.

In Use: Here the Bone Spreader Osteotome is in place, separating the block graft from the ramus on the mandible. The site was prepared and the crest burnished. The crest is then perforated to outline the graft using a .75mm round bur. The outline should be 4-5mm medial from the buccal plate. The length of the graft is determined using a caliper.

A 702L bur is used to connect the perforations to the length of the graft. Vertical anterior and posterior cuts are made with the 702L bur 10-12mm in length. At the vertical depth, the buccal plate is scored using a #8 round bur. The block graft is out-fractured using the bone spreader osteotome.



OSTEOTOMES

ACE Spatula Osteotome

0500006	5 x 10mm, 5-½"	each	\$35.99	\$35.99	\$32.39	\$32.39
0500010	4 x 30mm, 6"	each	\$44.79	\$44.79	\$40.31	\$40.31
0500011	5 x 35mm, 6"	each	\$44.79	\$44.79	\$40.31	\$40.31

ACE Maxillary Osteotome

0501026	6 x 1.2mm, Straight, 7"	each	\$80.29	\$80.29	\$72.26	\$72.26
---------	-------------------------	------	---------	---------	---------	---------



**ACE ALLOGRAFTS.
ACE XENOGRAFTS.
ACE ALLOPLASTS.**

BONE REGEN BEGINS ON PAGE 70

OR PURCHASE ANYTIME AT ACESOUTHERN.COM



Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price


INSTRUMENTS

CALIPERS AND GAUGES

All instruments Stainless Steel unless noted.

Mitutoyo Digital Caliper

- Cannot be autoclaved
- With thumb wheel
- 0-6", 0-150mm



8647A42	9" (23cm)	each	\$221.99	\$221.99	\$199.79	\$199.79
---------	-----------	------	----------	----------	----------	----------

Boley Gauges

- 0-100mm



0800010	ACE 5-1/2" (14cm)	each	\$102.29	\$102.29	\$92.06	\$92.06
68694	Integra™ 5-1/2" (14cm)	each	\$57.79	\$57.79	\$52.01	\$52.01

ACE Sliding Caliper

- 0-95mm



0869301	5" (13cm)	each	\$245.99	\$245.99	\$221.39	\$221.39
---------	-----------	------	----------	----------	----------	----------

ACE Jameson Caliper

- 4" (10cm)



0800007	4" (10cm)	each	\$138.29	\$138.29	\$124.46	\$124.46
---------	-----------	------	----------	----------	----------	----------

ACE Castroviejo Calipers - short beak

0800015	3-1/4" 0-20mm, straight	each	\$201.79	\$201.79	\$181.61	\$181.61
0800016	3-1/4" 0-20mm, 45° angle	each	\$234.99	\$234.99	\$211.49	\$211.49
0800017	3-1/4" 0-40mm, straight	each	\$234.99	\$234.99	\$211.49	\$211.49
0800018	3-1/2" 0-40mm, 45° angle	each	\$261.79	\$261.79	\$235.61	\$235.61

ACE Castroviejo Calipers - long beak

0800019	6-1/2" 0-40mm, straight	each	\$261.79	\$261.79	\$235.61	\$235.61
0800020	6-1/2" 0-40mm, 45° angle	each	\$281.29	\$281.29	\$253.16	\$253.16

Integra™ Castroviejo Calipers - short beak

18650	3-1/4" 0-22mm, straight	each	\$376.79	\$376.79	\$339.11	\$339.11
-------	-------------------------	------	----------	----------	----------	----------

ACE Bone Implant Caliper

- 0-20mm

0800050	3-1/2" (9cm)	each	\$188.79	\$188.79	\$169.91	\$169.91
---------	--------------	------	----------	----------	----------	----------

Integra™ Crown Gauge

- w/ 2 Sided Scale

68693		each	\$106.29	\$106.29	\$95.66	\$95.66
-------	--	------	----------	----------	---------	---------

CALIPERS AND GAUGES - continued

All instruments Stainless Steel unless noted.

ACE Ridge Mapping Caliper

- 0-25mm



0800059	6" (15cm)	each	\$248.29	\$248.29	\$223.46	\$223.46
---------	-----------	------	----------	----------	----------	----------

ACE Depth Packing Gauge

- For use with 3.3mm, 3.75mm and 4.0mm self-tapping implants
- Contoured ends allow use as packing instrument and in sinus elevation procedures
- 8, 10, 13, 15, 18mm markings • Titanium (6AI-4V-ELI)



7140243	6-1/2" (16.5cm)	each	\$111.79	\$111.79	\$100.61	\$100.61
---------	-----------------	------	----------	----------	----------	----------

Implant Depth Gauge with Ball End

- Ball end is useful for examination of bone preparation and in cases where the floor of the bone osteotomy site has been perforated
- Use at the 2.0mm drilling stage
- 7, 8, 10, 13, 15, 18, 20mm markings
- Titanium (6AI-4V-ELI)



9600639	ACE, 4" (10cm)	each	\$74.79	\$74.79	\$67.31	\$67.31
62170	Integra™, 4" (10cm)	each	\$135.29	\$135.29	\$121.76	\$121.76

ACE LACRIMAL PROBES

- Sterling
- 5" (13cm)

Actual Size	•	•	•	•	•	•	•	•	•	•	•	•
Ø mm	0.6	0.7	0.8	0.9	1.0	1.1	1.2	1.3	1.4	1.5	1.6	1.7
Ø Inches	.023	.027	.031	.035	.039	.043	.047	.051	.055	.059	.063	.066
Probe Size	0000	000	00	0	1	2	3	4	5	6	7	8

Bowman – Standard (no bulb)

6122540	Size 0000-000	each	\$40.49	\$40.49	\$36.44	\$36.44
6122520	Size 00-0	each	\$40.49	\$40.49	\$36.44	\$36.44
6122512	Size 1-2	each	\$40.49	\$40.49	\$36.44	\$36.44
6122534	Size 3-4	each	\$40.49	\$40.49	\$36.44	\$36.44
6122556	Size 5-6	each	\$40.49	\$40.49	\$36.44	\$36.44
6122578	Size 7-8	each	\$40.49	\$40.49	\$36.44	\$36.44

Williams (with bulb)

6125040	Size 0000-000	each	\$40.49	\$40.49	\$36.44	\$36.44
6125020	Size 00-0	each	\$40.49	\$40.49	\$36.44	\$36.44
6125012	Size 1-2	each	\$40.49	\$40.49	\$36.44	\$36.44
6125034	Size 3-4	each	\$40.49	\$40.49	\$36.44	\$36.44
6125056	Size 5-6	each	\$40.49	\$40.49	\$36.44	\$36.44
6125078	Size 7-8	each	\$40.49	\$40.49	\$36.44	\$36.44

INSTRUMENTS

PERIODONTAL PROBES

DE = double ended

ACE Probe

- For measuring bony osteotomies in hard-to-reach areas
- 8, 10, 13, 15, 18mm Markings
- Titanium (6Al-4V-ELI)

6301105	6" (15cm)	each	\$90.29	\$90.29	\$81.26	\$81.26
---------	-----------	------	---------	---------	---------	---------

American Eagle Instruments Probes

- Precisely calibrated probes provide accuracy for intraoral measurements
- 2 Handle Options: Stainless Steel (SS), standard or EagleLite™ Resin (X) (see sample images to the left)
- Handles are available as single ended or double ended (DE)



3/8" - 9.5mm EagleLite™ Resin (X)



7/32" - 5.6mm Stainless Steel (SS), standard

Marquis

AEP12BX	12, resin X, DE	each	\$25.49	\$25.49	\$22.94	\$22.94
AEP12B	12, SS, DE	each	\$27.29	\$27.29	\$24.56	\$24.56



Marquis

AEP12YX	12 Yellow, resin X	each	\$29.79	\$29.79	\$26.81	\$26.81
AEP12Y	12 Yellow, SS	each	\$28.79	\$28.79	\$25.91	\$25.91



South Dakota

AEPSD12YX	12 Yellow, resin X	each	\$29.79	\$29.79	\$26.81	\$26.81
AEPSD12Y	12 Yellow, SS	each	\$28.79	\$28.79	\$25.91	\$25.91



Yellow/Green

AEP12Y/GX	12, resin X	each	\$29.79	\$29.79	\$26.81	\$26.81
AEP12Y/G	12, SS	each	\$28.79	\$28.79	\$25.91	\$25.91



Hu-Friedy® Goldman-Fox/Williams Probe

PGFW6	DE, satin steel™ handle	each	\$32.49	\$32.49	\$29.24	\$29.24
-------	-------------------------	------	---------	---------	---------	---------



PERIODONTAL PROBES - continued

DE = double ended

Integra™ Periodontal Probes

696212	DE	each	\$33.99	\$33.99	\$30.59	\$30.59
696218	X2 probe	each	\$33.99	\$33.99	\$30.59	\$30.59
696216	C probe	each	\$33.99	\$33.99	\$30.59	\$30.59

J&J Instruments Periodontal Probes



Williams Round

16540	1, 2, 3, 5, 7, 8, 9, 10mm	each	\$10.29	\$10.29	\$9.26	\$9.26
-------	---------------------------	------	---------	---------	--------	--------

#12 Probe

16540	3, 6, 9, 12mm	each	\$10.29	\$10.29	\$9.26	\$9.26
-------	---------------	------	---------	---------	--------	--------

PERIODONTAL EXPLORERS

DE = double ended

American Eagle Instruments Explorers

- Providing the best tactile feel, American Eagle explorers are extremely sharp and have excellent memory
- 2 Handle Options: Stainless Steel (SS), standard or EagleLite™ Resin (X) (see sample images to the left)
- Handles are available as single ended or double ended (DE)

Pigtail 2

AEEXP2X	11-12, resin X, DE	each	\$22.79	\$22.79	\$20.51	\$20.51
AEEXP2	11-12, SS, DE	each	\$21.29	\$21.29	\$19.16	\$19.16



DE 5

AEEXP5X	23-17, resin X, DE	each	\$22.79	\$22.79	\$20.51	\$20.51
AEEXP5	23-17, SS, DE	each	\$21.29	\$21.29	\$19.16	\$19.16



Shepard's Hook

AEEXP23X	Single end, resin X	each	\$16.29	\$16.29	\$14.66	\$14.66
AEEXP23	Single end, SS	each	\$14.99	\$14.99	\$13.49	\$13.49



J&J Instruments Explorers

04850	#5, DE	each	\$7.49	\$7.49	\$6.74	\$6.74
04820	#2, DE	each	\$7.49	\$7.49	\$6.74	\$6.74




Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INSTRUMENTS

ACE CEPHALOMETRIC MEASURING INSTRUMENTS

- Cold sterilization only



ACE Protractor

0032936	11-1/2" (29cm)	each	\$17.49	\$17.49	\$15.74	\$15.74
---------	----------------	------	---------	---------	---------	---------




ACE Incisor/Molar Tracing Template

0032938		each	\$11.29	\$11.29	\$10.16	\$10.16
---------	--	------	---------	---------	---------	---------

ACE PLASTIC INTERMAXILLARY RULER

- Specifically designed to confirm intermaxillary opening
- Millimeter increments
- Cold sterilization only
- 4-3/4" x 3-3/4"



0500201	0-60mm markings	each	\$15.49	\$15.49	\$13.94	\$13.94
---------	-----------------	------	---------	---------	---------	---------



ACE STAINLESS STEEL RULERS


0800003	1/2" x 4"	each	\$13.49	\$13.49	\$12.14	\$12.14
0800011	1/2" x 6"	each	\$13.49	\$13.49	\$12.14	\$12.14



INTEGRA™ FLEXIBLE STAINLESS RULER

- Graduated in 1/32" fractions and mm

18660	6" (15 cm)	each	\$10.99	\$10.99	\$9.89	\$9.89
-------	------------	------	---------	---------	--------	--------



ACE PLASTIC RULER

- Standard and metric increments
- Cold sterilization only


0032937	7/8" x 6"	10/pkg	\$18.79	\$18.79	\$16.91	\$16.91
---------	-----------	--------	---------	---------	---------	---------

ACE PATIENT OPENING GAUGE

- Allows patient to monitor progress at home
- Cold sterilization only




0500204		each	\$8.29	\$8.29	\$7.46	\$7.46
---------	--	------	--------	--------	--------	--------



CHECK OUT THE WEB STORE FOR A MORE COMPLETE SELECTION OF SURGICAL INSTRUMENTS!

VISIT ACESOUTHERN.COM TODAY

INSTRUMENTS



ACE X-RAY MAGNIFICATION CALCULATOR™

- Specially designed for dental implantology and orthodontics.
- Eliminates the complex calculations for determining panoramic x-ray distortion.
- Measurements can easily be obtained to determine the precise magnification of the x-ray and to accurately determine the amount of bone available for implant placement.
- Can be used with the Delrin Bite Block with Stainless Steel Pin (01-004-48), (01-004-56) or Exacta Mouth Prop® with 6mm Marking Ball (01-004-50).
- Actual size: 3-3/4" x 9-1/2", 9.5cm x 24cm.
- Cold sterilization only.

0031967	each	\$30.49	\$30.49	\$27.44	\$27.44
---------	------	---------	---------	---------	---------

ACE EXACTA® MOUTH PROPS

For use with #0031967 ACE X-Ray Magnification Calculator (above)

Exacta® Mouth Prop with Calculator Pin

- Exacta mouth prop has a 20mm pin



0100456	20mm calculator pin	each	\$54.79	\$54.79	\$49.31	\$49.31
---------	---------------------	------	---------	---------	---------	---------

Exacta® Mouth Prop with 6mm Marking Ball

- Exacta mouth prop has a 6mm spherical marking ball



0100450	6mm marking ball	each	\$54.79	\$54.79	\$49.31	\$49.31
---------	------------------	------	---------	---------	---------	---------

Exacta® Mouth Prop with 5mm and 10mm Pin

- The size of the pin on a panoramic x-ray is comparable to the original pin length. The amount of elongation is determined and the exact measurement extrapolated. (If a 10mm pin increases to 12mm, a 20% elongation is realized – the x-ray bone height is reduced by 20%.)



0100350	5mm and 10mm pins	each	\$54.79	\$54.79	\$49.31	\$54.79
---------	-------------------	------	---------	---------	---------	---------

Exacta® Mouth Prop with 10mm Pin

- Allows simple and exact determination of vertical bone height in dental implantology
- Autoclavable
- Contains a 10mm pin
- For single tooth or narrow edentulous area



0100448	10mm pin	each	\$28.49	\$28.49	\$25.64	\$25.64
---------	----------	------	---------	---------	---------	---------



MCKESSON TYPE SILICONE MOUTH PROPS

Silicone mouth props support the mandible and relieve muscle tension by providing a soft, pacifying surface to bite on. Reduces patient gagging and fatigue. Ideal for children and adults. Finger channel for insertion. Deep corrugation for secure stabilization. Virtually indestructible silicone will provide years of reusable service.

- Autoclavable
- Not shown actual size
- Silicone — Latex-free

Small Child

0116035	A Black	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0116039	Blue	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0116036	Tan	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
PAL101	#01Tan	Palmero	pair	\$18.79	\$18.79	\$16.91	\$16.91

Child

0116037	B Black	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0116040	C Blue	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0116038	D Tan	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
101A	Tan	Palmero	pair	\$18.79	\$18.79	\$16.91	\$16.91

Adult

0117041	E Black	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0117043	Blue	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
0117042	Tan	ACE	pair	\$25.99	\$25.99	\$23.39	\$23.39
101B	Tan	Palmero	pair	\$18.79	\$18.79	\$16.91	\$16.91

DCI RUBBER MOUTH PROPS

- Rubber over brass construction
- No latex used
- Gray, sterilizable
- Autoclavable



2436	E Large adult	DCI	pair	\$76.49	\$76.49	\$68.84	\$68.84
2437	F Edentulous	DCI	pair	\$133.29	\$133.29	\$119.96	\$119.96

HU-FRIEDY® MARKEL MOUTH PROPS

- Latex-free



MP52	Large adult	Hu-Friedy®	pair	\$54.29	\$54.29	\$48.86	\$48.86
MP54	Edentulous	Hu-Friedy®	pair	\$54.29	\$54.29	\$48.86	\$48.86

INSTRUMENTS



INTEGRA™ VIENNA LIGHTWEIGHT NASAL SPECULA

206	Adult, 5-3/4"	each	\$156.49	\$156.49	\$140.84	\$140.84
207	Child, 5-3/4"	each	\$156.49	\$156.49	\$140.84	\$140.84

ACE PLASTIC LIP/CHEEK RETRACTORS


- Autoclavable
- Can be steam sterilized to 270°F (130°C)
- Dry heat not recommended

Double-Ended

0100030	E Universal, 4" (40mm)	2/pkg	\$28.49	\$27.07	\$27.07	\$25.64
---------	-------------------------------	-------	---------	---------	---------	---------

Single Ended

0100038	F Adult, 4-1/4" (45mm)	2/pkg	\$29.49	\$28.02	\$28.02	\$26.54
---------	-------------------------------	-------	---------	---------	---------	---------



MOUTH GAGS

Autoclavable

MOLT MOUTH GAGS

0117016	Adult 5-1/2"	ACE	each	\$111.79	\$111.79	\$100.61	\$100.61
0116018	A Child 4-1/2"	ACE	each	\$115.49	\$115.49	\$103.94	\$103.94
6310	Ped./small child	Integra™	each	\$355.29	\$355.29	\$355.29	\$319.76
6315	Child	Integra™	each	\$355.29	\$355.29	\$355.29	\$319.76
6320	B Adult	Integra™	each	\$355.29	\$355.29	\$355.29	\$319.76
07900	Child	J&J Instru.	each	\$58.99	\$58.99	\$58.99	\$53.09
07910	C Adult	J&J Instru.	each	\$55.49	\$55.49	\$55.49	\$49.94

Replacement Covers for Molt Mouth Gags

0117023	Child/adult, contains latex	ACE	2/pack	\$8.29	\$8.29	\$8.29	\$7.46
PTRTC	Child, silicone, latex-free	Hu-Friedy®	2/pack	\$11.99	\$11.99	\$11.99	\$10.79
PTRTA	Adult, silicone, latex-free	Hu-Friedy®	2/pack	\$11.99	\$11.99	\$11.99	\$10.79
6310A	Ped./small child, latex-free	Integra™	2/pack	\$7.99	\$7.99	\$7.99	\$7.19
6315A	Child, silicone, latex-free	Integra™	2/pack	\$7.99	\$7.99	\$7.99	\$7.19
6320A	Adult, silicone, latex-free	Integra™	2/pack	\$7.99	\$7.99	\$7.99	\$7.19
07905	Child/adult, silicone latex-free	J&J Instru.	2/pack	\$5.79	\$5.79	\$5.79	\$5.21

DENHARDT MOUTH GAGS

0100025	D 5"	ACE	each	\$128.79	\$128.79	\$115.91	\$115.91
---------	-------------	-----	------	----------	----------	----------	----------

SELF-SPANNING PLASTIC LIP/CHEEK RETRACTORS

- Autoclavable
- Latex-free
- Can be steam sterilized to 250°F (121°C) for 3 to 3-1/2 minutes
- Dry heat not recommended

ACE Self-Spanning

0116032	G Child, 35mm	each	\$24.99	\$24.99	\$22.49	\$22.49
0117031	Adult, 50mm	each	\$21.99	\$21.99	\$19.79	\$19.79


ACE Self-Spanning with Holders

0116025	Child, 35mm	each	\$28.29	\$28.29	\$25.46	\$25.46
0117025	H Adult, 50mm	each	\$27.49	\$27.49	\$24.74	\$24.74

Hager Worldwide Spandex Lip and Cheek Retractors

- Two-sided
- Allows clear view of work area
- Non-reflective

605455	Child	2/pack	\$40.79	\$40.79	\$36.71	\$36.71
605454	Adult	2/pack	\$42.49	\$42.49	\$38.24	\$38.24



INSTRUMENTS

**STAINLESS STEEL
LIP/CHEEK/TONGUE RETRACTORS**

- Autoclavable

Apfel

0100163	A	6-¼" (16cm)	ACE	each	\$77.79	\$77.79	\$70.01	\$70.01
---------	----------	-------------	-----	------	---------	---------	---------	---------

Austin, #2

0100158	B	3" (7.5cm)	ACE	each	\$52.49	\$52.49	\$47.24	\$47.24
07500	C	#2	J&J Instru.	each	\$12.99	\$12.99	\$11.69	\$11.69

Austin Retractor, University of Minnesota

0100160	D	3" (7.5cm)	ACE	each	\$48.79	\$48.79	\$43.91	\$43.91
---------	----------	------------	-----	------	---------	---------	---------	---------

Bishop, Tongue and Cheek

0100077	E	6" (15cm)	ACE	each	\$41.49	\$41.49	\$37.34	\$37.34
---------	----------	-----------	-----	------	---------	---------	---------	---------

Branemark® Type

6100390	F	6-¼" (16cm)	ACE	each	\$111.79	\$111.79	\$100.61	\$100.61
---------	----------	-------------	-----	------	----------	----------	----------	----------

Catania Cheek

Double Ended

0100107	G	5-½" (14cm)	ACE	each	\$50.99	\$50.99	\$45.89	\$45.89
---------	----------	-------------	-----	------	---------	---------	---------	---------

Columbia

Double Ended

0100391	H	5" (13cm)	ACE	each	\$34.99	\$34.99	\$31.49	\$31.49
---------	----------	-----------	-----	------	---------	---------	---------	---------

Feldman

Double Ended

0100394	I	5-½" (14cm)	ACE	each	\$35.99	\$35.99	\$32.39	\$32.39
---------	----------	-------------	-----	------	---------	---------	---------	---------

Ronnau

0100162	J	6" (15cm)	ACE	each	\$77.79	\$77.79	\$70.01	\$70.01
---------	----------	-----------	-----	------	---------	---------	---------	---------

Simplex, Lip/Cheek

0100404	K	4-½" (11.5cm)	ACE	each	\$48.79	\$48.79	\$43.91	\$43.91
---------	----------	---------------	-----	------	---------	---------	---------	---------


University of Minnesota, Tongue/Cheek

Double Ended


0100088	L	5-½" (14cm)	ACE	each	\$28.79	\$28.79	\$25.91	\$25.91
0100092		6" (15cm) Bishop	ACE	each	\$28.49	\$28.49	\$25.64	\$25.64
437		5-½" (14cm) Cawood	A-Titan	each	\$36.79	\$36.79	\$33.11	\$33.11
22D209		5-½" (14cm)	Integra™	each	\$31.79	\$31.79	\$28.61	\$28.61
07505	M		J&J Instru.	each	\$13.99	\$13.99	\$12.59	\$12.59

Weider

0116095	N	5" (13cm) Child	ACE	each	\$57.29	\$57.29	\$51.56	\$51.56
0117094		5" (13cm) Adult	ACE	each	\$57.29	\$57.29	\$51.56	\$51.56
WEIDER2		Serrated	A-Titan	each	\$62.79	\$62.79	\$56.51	\$56.51
222	O	28mm, sm blade	Integra™	each	\$67.49	\$67.49	\$60.74	\$60.74
224		36mm, lg blade	Integra™	each	\$67.49	\$67.49	\$60.74	\$60.74



INSTRUMENTS



ACE INTRAORAL OCCLUSAL PHOTO MIRRORS WITH MALLEABLE HANDLE

- With malleable handle
- Single sided only



5961002	3 3/4" (9.5cm)	each	\$101.29	\$101.29	\$91.16	\$91.16
---------	----------------	------	----------	----------	---------	---------

ACE INTRAORAL PHOTO MIRRORS

Chromium Coated, Front Surface - Not shown actual size

- Material: float glass, thickness 1/8" (3mm) ground edges; absolutely flat surfaces even on edges (unlike metallic mirrors)
- Coating: chromium coating on both sides by sputtering vacuum metallization
- Reflection: reflection capacity = 65%; sharp images without any deformation and without double reflection (front surface)
- Sterilization: autoclavable up to 428°F (220°C)
- Can be used on both sides

Lingual Intraoral

0032969	A Adult, 6"	each	\$61.99	\$61.99	\$55.79	\$55.79
0033002	B Adult, narrow, 6"	each	\$61.99	\$61.99	\$55.79	\$55.79

Buccal Intraoral

0032970	C Adult, 6"	each	\$61.99	\$61.99	\$55.79	\$55.79
0032989	D Adult, 6"	each	\$61.99	\$61.99	\$55.79	\$55.79

Occlusal Intraoral

0032971	E Adult, 4"	each	\$61.99	\$61.99	\$55.79	\$55.79
0032972	F Pediatric, 3 1/4"	each	\$61.99	\$61.99	\$55.79	\$55.79
5961008	G Oversized, 5"	each	\$85.79	\$85.79	\$77.21	\$77.21
5961009	H Oversized, 5 1/4"	each	\$85.79	\$85.79	\$77.21	\$77.21

INTEGRA™ INTRA ORAL MIRRORS

Lingual Intraoral

01726807	I Adult, #2	each	\$126.99	\$126.99	\$114.29	\$114.29
----------	--------------------	------	----------	----------	----------	----------


Buccal Intraoral

01726806	J Adult, #1	each	\$126.99	\$126.99	\$114.29	\$114.29
----------	--------------------	------	----------	----------	----------	----------

Occlusal Intraoral

01726805	K Adult, #3	each	\$126.99	\$126.99	\$114.29	\$114.29
01726808	L Child, #3	each	\$126.99	\$126.99	\$114.29	\$114.29

9



DISPOSABLE DENTAL MIRRORS

Tess® Mouth Mirrors

Mouth mirrors enable you to see at home what the dentist sees in the examination chair. The ergonomic handle is made of a strong but lightweight plastic. The mirror is glass for a clearer image and resists scratching unlike the reflective stickers you see on many.

- Assorted Colors: Blue, Green, Pink, Purple, Yellow

46072ASST	72/box	\$34.79	\$33.79	\$33.05	\$31.31
-----------	--------	---------	---------	---------	---------

Palmero Health Care Mouth Mirrors

Clear, distortion-free glass mirrors with an ergonomic handle and cheek-retracting strength. Ultra-thin glass for enhanced imaging. Ideal for high-volume offices where interim sterilization is not practical. 7/8" diameter.



93	72/box	\$43.49	\$43.49	\$39.14	\$39.14
----	--------	---------	---------	---------	---------

INSTRUMENTS

Plane mirrors may produce a "double image" as light reflects off both the glass surface and the silver coating behind the glass. Front surface mirrors have rhodium surfaces, which prevent this double reflection, providing clear, distortion-free images.



ACE DENTAL MIRRORS

Simple Stem		Cone Socket				
9700010	Size 3, plane	each	\$4.79	\$4.79	\$4.31	\$4.31
9700012	Size 4, plane	each	\$4.29	\$4.29	\$3.86	\$3.86
9700014	Size 5, plane	each	\$4.29	\$4.29	\$3.86	\$3.86

Cone Socket						
9700011	Size 3, plane	each	\$4.79	\$4.79	\$4.31	\$4.31
9700013	Size 4, plane	each	\$4.29	\$4.29	\$3.86	\$3.86
9700015	Size 5, plane	each	\$4.29	\$4.29	\$3.86	\$3.86

Simple Stem Mirror Handle						
9700030		each	\$8.99	\$8.99	\$8.09	\$8.09

Cone Socket Mirror Handle						
9700031	(mirror shown sold separately)	each	\$8.99	\$8.99	\$8.09	\$8.09

AMERICAN EAGLE INSTRUMENTS MOUTH MIRRORS

Cone Socket Dental Mirrors						
AEM4	#4	each	\$8.49	\$8.49	\$7.64	\$7.64
AEM5	#5	each	\$11.29	\$11.29	\$10.16	\$10.16

Mirror Handle EagleLite®						
AEMHRX	3/8" resin, round	each	\$17.99	\$17.99	\$16.19	\$16.19
AEMHRZ	3/8" stainless, round	each	\$17.99	\$17.99	\$16.19	\$16.19
AEMHR	1/4" stainless, round	each	\$16.29	\$16.29	\$14.66	\$14.66

INTEGRA™ MILTEX DENTAL MIRRORS

676984	size 4	dozen	\$27.49	\$27.49	\$24.74	\$24.74
--------	--------	-------	---------	---------	---------	---------

HU-FRIEDY® DENTAL MIRRORS

Cone Socket, Front Surface						
MIR4	size 4	each	\$7.29	\$7.29	\$6.56	\$6.56
MIR5	size 5	each	\$7.29	\$7.29	\$6.56	\$6.56

Double-Sided Cone Socket, Front Surface

- Two back to back, front surface mirrors
- Due to contra-angle design, these mirrors also act as a cheek, tongue and lip retractor
- Provide better visibility
- Autoclavable



MIR5DS3	size 5	3/box	\$25.29	\$25.29	\$22.76	\$22.76
---------	--------	-------	---------	---------	---------	---------

Cone Socket Mirror Handles

MH1	#1 5" (12.5cm)	each	\$18.79	\$18.79	\$16.91	\$16.91
MH24	#24 5" (12.5cm)	each	\$18.79	\$18.79	\$16.91	\$16.91

J&J INSTRUMENTS MOUTH MIRRORS

Cone Socket Mouth Mirrors

07210	#4	each	\$18.29	\$18.29	\$16.46	\$16.46
07230	#5	each	\$18.29	\$18.29	\$16.46	\$16.46

Cone Socket Mirror Handles

07405	Octagonal	each	\$2.49	\$2.49	\$2.24	\$2.24
07420	Round	each	\$2.49	\$2.49	\$2.24	\$2.24



A INSTRUMENTS

ACE ARCH BARS

- Meter roll (approximately 36")
- Not MRI compatible



0800110	each	\$37.79	\$37.79	\$34.01	\$34.01
---------	------	---------	---------	---------	---------

BAKER® WINTER SPLINT KIT

Baker Winter Fracture Splint Kit: Rhodium covered.

Comes complete with easily adjustable upper and lower arch bars, tie wires, waxed silk ligatures and rubber bands. Directions included.



130801	each	\$244.29	\$244.29	\$219.86	\$219.86
--------	------	----------	----------	----------	----------

Baker Winter Arch Bar



HDS130901	each	\$113.99	\$113.99	\$102.59	\$102.59
-----------	------	----------	----------	----------	----------

ACE UTILITY ARCH BAR WAX

0800160	2" Ropes	10 boxes/pkg	\$48.79	\$48.79	\$43.91	\$43.91
---------	----------	--------------	---------	---------	---------	---------



ACE SPOOL WIRE ON AUTOCLAVABLE DELRIN SPOOLS


- 8% nickel in tube and spool line
- 4 oz. Spool
- 200 feet
- Stainless steel



0800087	28G, .012"	each	\$44.79	\$44.79	\$40.31	\$40.31
0800088	26G, .016"	each	\$44.79	\$44.79	\$40.31	\$40.31
0800089	25G, .018"	each	\$44.79	\$44.79	\$40.31	\$40.31
0800090	24G, .020"	each	\$44.79	\$44.79	\$40.31	\$40.31
0800091	22G, .025"	each	\$44.79	\$44.79	\$40.31	\$40.31
0800092	20G, .032"	each	\$44.79	\$44.79	\$40.31	\$40.31

ACE TUBE WIRE

- 8% nickel in tube wire
- Stainless steel
- 18" strands



0800182	28G, .012"	50/tube	\$35.99	\$35.99	\$32.39	\$32.39
0800100	26G, .016"	50/tube	\$35.99	\$35.99	\$32.39	\$32.39
0800101	25G, .018"	50/tube	\$35.99	\$35.99	\$32.39	\$32.39
0800102	24G, .020"	50/tube	\$35.99	\$35.99	\$32.39	\$32.39

WALL MEDICAL MAXILLOFACIAL FRACTURE KITS

Both Wall and Erich Splints are available in dozen packs, as are wire coils and elastic bands.

- The unique design of the Unisplint Kit contributes substantially to the ease in treating fractures of the mandible and maxilla
- Its simplicity lends itself to emergency treatment
- Splints are preformed and have the correct malleability for easy custom contouring to the patient's dental arch
- All surfaces have a very high polish; The lugs are smooth and rounded, reducing the incidence of abrasion




Unisplint Maxillofacial Kit (Wall Universal Splint)

Includes: two 5" Wall universal splints, 3 yards of 25g dead soft temper wire and approximately fifty 1/8" I.D. latex rubber bands

U100	1 kit	\$64.99	\$64.99	\$58.49	\$58.49
------	-------	---------	---------	---------	---------

Wall Universal Maxillofacial Fracture Splints

The Wall Splint has a universal design that prevents upside-down placement. Its smooth, trapezoidal lugs are placed so that either elastic bands or wires may be attached with ease.



U101	each	\$92.99	\$92.99	\$83.69	\$83.69
------	------	---------	---------	---------	---------

Erich Maxillofacial Kit

Includes: Two 5" Erich splints, 3 yards of 25g dead soft temper stainless steel wire and approximately fifty 1/8" I.D. latex rubber bands

U400	1 kit	\$64.99	\$64.99	\$58.49	\$58.49
------	-------	---------	---------	---------	---------

Erich Maxillofacial Fracture Splint

U401	5"	12/pack	\$127.29	\$127.29	\$114.56	\$114.56
------	----	---------	----------	----------	----------	----------



Erich Arch Bar Coil

U402	1 meter coil	each	\$55.29	\$55.29	\$49.76	\$49.76
------	--------------	------	---------	---------	---------	---------

1/8" I.D. Latex Rubber Bands

approx 50/pack

U102	12 packs/box	\$21.49	\$21.49	\$19.34	\$19.34
------	--------------	---------	---------	---------	---------

25g Dead Soft Temper Stainless Steel Wire

U103	3 yard coils	12 coils/box	\$51.49	\$51.49	\$46.34	\$46.34
------	--------------	--------------	---------	---------	---------	---------



ACE ELASTICS

- Contains Latex

0800128	1/8" Light	100/pack	\$4.79	\$4.55	\$4.55	\$4.31
0800133	3/16" Medium	100/pack	\$4.79	\$4.55	\$4.55	\$4.31
0800134	3/16" Heavy	100/pack	\$4.79	\$4.55	\$4.55	\$4.31
0800135	3/16" Extra Heavy	100/pack	\$4.79	\$4.79	\$4.31	\$4.31
0800136	1/4" Light	100/pack	\$4.79	\$4.55	\$4.55	\$4.31
0800137	1/4" Medium	100/pack	\$4.79	\$4.79	\$4.31	\$4.31
0800138	1/4" Heavy	100/pack	\$4.79	\$4.79	\$4.31	\$4.31
0800142	5/16" Heavy	100/pack	\$4.79	\$4.79	\$4.31	\$4.31
0800146	3/8" Heavy	100/pack	\$4.79	\$4.79	\$4.31	\$4.31

INSTRUMENTS

ORTHO TRACTION PADS

Button Pad Oval Curved

Oval pad, gently curved for cuspids.
Chain swivels 360 degrees.

BPOC	4mm x 3.1mm	each	\$26.29	\$26.29	\$23.66	\$23.66
------	-------------	------	--------------------	---------	---------	---------



Button Pad Round

Small round pad with very slight concave shape.
Chain swivels 360 degrees.

BPRO	3.4mm diameter	10/pack	\$250.29	\$250.29	\$225.26	\$225.26
------	----------------	---------	---------------------	----------	----------	----------



Large Pad Low Profile

Large rectangular pad, nicely curved for cuspids.
Low profile so nothing interferes with tissue

LPLP	4.3mm x 3.5mm	5/pack	\$167.99	\$167.99	\$151.19	\$151.19
------	---------------	--------	---------------------	----------	----------	----------



Mesh Pad Eyelet Profile

Mesh screen for great bonding. Your adhesive flows into the mesh. Eyelet profiles so it's easy to pick up.

MPEP	3.9mm diameter	5/pack	\$167.99	\$167.99	\$151.19	\$151.19
------	----------------	--------	---------------------	----------	----------	----------



Mesh Pad Low Profile

Mesh screen is for great bonding. Low profile allows no interference with tissue. Trim and bend to work anywhere.

MPLP	3.9mm diameter	5/pack	\$167.99	\$167.99	\$151.19	\$151.19
------	----------------	--------	---------------------	----------	----------	----------



Round Pad Eyelet Profile

Small round pad with very slight concave shape. Eyelet profile so it's easy to pick up. Great in "holes."

RPEP	3.4mm diameter	5/pack	\$167.99	\$167.99	\$151.19	\$151.19
------	----------------	--------	---------------------	----------	----------	----------



Small Pad Eyelet Profile

Small rectangular pad, less curvey for non-cuspids. Eyelets profile so it's easy to pick up.

SPEP	3.5mm x 3mm	5/pack	\$167.99	\$167.99	\$151.19	\$151.19
------	-------------	--------	---------------------	----------	----------	----------



ASEPTICO® TECHNIQUE REFILL CHAIN & GRID

The Technique Refill traction system, or chain ligation system, is made for the treatment of altered eruption and impaction of teeth. It is available in three styles: Cuspid, Bicuspid, and Incisor exposing brackets, or orthodontic buttons.

The ligation chain is 7/8" (22m) in length, and its 14k gold construction reduces irritation for the patient during the treatment of erupted and impacted teeth.



Cuspid Type Refill

AMA09A	10/box	\$387.79	\$387.79	\$349.01	\$349.01
--------	--------	---------------------	----------	----------	----------

Bicuspid Type Refill

AMA09A2	10/box	\$387.79	\$387.79	\$349.01	\$349.01
---------	--------	---------------------	----------	----------	----------

Incisor Type Refill

AMA09A3	10/box	\$387.79	\$387.79	\$349.01	\$349.01
---------	--------	---------------------	----------	----------	----------

HU-FRIEDY® DIRECT BOND BRACKET HOLDER TWEEZERS

for orthodontic chains

678212	each	\$71.49	\$71.49	\$64.34	\$64.34
--------	------	--------------------	---------	---------	---------



INTEGRA™ DIRECT BOND BRACKET HOLDER


74212	5"	each	\$38.29	\$38.29	\$34.46	\$34.46
-------	----	------	--------------------	---------	---------	---------



**GET INSIDER ACCESS
ORDER. TRACK. FAVORITES.
MANAGE. AND MORE.**

SIGN UP FOR WEB ACCESS TODAY

[ACESOUTHERN.COM/CUSTOMER/ACCOUNT/CREATE/](https://acesouthern.com/customer/account/create/)



INSTRUMENTS

**PREMIER UNI-PUNCH®
DISPOSABLE BIOPSY
PUNCHES**

- Sterile, seamless, razor sharp blades



9033509	1.5mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033502	2.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033511	2.5mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033503	3.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033510	3.5mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033504	4.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033505	5.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033506	6.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36
9033508	8.0mm	25/box	\$49.29	\$49.29	\$44.36	\$44.36

Assortment

Includes: 5 each of 2.0mm, 3.0mm, 3.5mm, 4.0mm, & 5.0mm

9033525	25/box	\$49.29	\$49.29	\$44.36	\$44.36
---------	--------	---------	---------	---------	---------

**INTEGRA™
STERILE DISPOSABLE
BIOPSY PUNCHES**

Feature a seamless stainless steel razor sharp cutting edge and a ribbed handle for maximum control and comfort. Each punch is sealed in a color-coded, peel-apart, see-through blister pack.



3331AA	1.0mm	50/box	\$118.49	\$118.49	\$106.64	\$106.64
3331A	1.5mm	50/box	\$118.49	\$118.49	\$106.64	\$106.64
3331	2.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3332	3.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3333	3.5mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3334	4.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3335	5.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3336	6.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06
3337	8.0mm	50/box	\$132.29	\$132.29	\$119.06	\$119.06

Assortment

Includes: 10 each of 2.0mm, 3.0mm, 4.0mm, 5.0mm, & 6.0mm

3338	50/box	\$128.79	\$128.79	\$115.91	\$115.91
------	--------	----------	----------	----------	----------

**INTEGRA™
SHORT HANDLED STERILE
DISPOSABLE BIOPSY PUNCHES**

Designed to enable dentists to take tissue samples in all areas of the mouth. Featuring the same razor-sharp blade as the original Integra biopsy punches, the shorter handle was modified to allow dental professionals access in even the hardest to reach, posterior areas of the mouth.



3333SH	3.5mm	50/box	\$116.49	\$116.49	\$104.84	\$104.84
3334SH	4.0mm	50/box	\$116.49	\$116.49	\$104.84	\$104.84
3335SH	5.0mm	50/box	\$116.49	\$116.49	\$104.84	\$104.84
3336SH	6.0mm	50/box	\$116.49	\$116.49	\$104.84	\$104.84

TISSUE PUNCH SIZE CHART

1.5mm	2mm	2.5mm	3mm	3.5mm	4mm	5mm	6mm	8mm
-------	-----	-------	-----	-------	-----	-----	-----	-----



ACE REUSABLE TISSUE PUNCHES

- Provides a precise easy method of soft tissue removal to insure positive seating of healing abutments
- Diameter of instrument is designed to extend over the implant for simplified tissue removal
- Can be used for tissue biopsies
- Honed cutting edge
- Hand-held
- Surgical grade stainless steel



Long Handle Straight Tissue Punches

- 4" (10cm)

9000003	3.0mm	each	\$35.99	\$35.99	\$32.39	\$32.39
9000004	4.0mm	each	\$35.99	\$35.99	\$32.39	\$32.39
9000005	5.0mm	each	\$35.99	\$35.99	\$32.39	\$32.39

Long Handle Angled Tissue Punches

- 4" (10cm)

9000060	3.0mm	each	\$56.29	\$56.29	\$50.66	\$50.66
9000080	4.0mm	each	\$56.29	\$56.29	\$50.66	\$50.66

Short Handle Half-Cut Tissue Punches

- 1-1/4" (3cm)

20178501	3.5mm	each	\$41.49	\$41.49	\$37.34	\$37.34
20178502	4.0mm	each	\$41.49	\$41.49	\$37.34	\$37.34
20178503	5.0mm	each	\$41.49	\$41.49	\$37.34	\$37.34
20178504	6.0mm	each	\$41.49	\$41.49	\$37.34	\$37.34

**ACE REUSABLE
CONTRA ANGLE
MUCOSA PUNCHES**

- 1" (25mm)



9200002	2.0mm	each	\$71.29	\$71.29	\$60.60	\$60.60
9200003	3.0mm	each	\$71.29	\$71.29	\$64.16	\$64.16
9200001	3.5mm	each	\$71.29	\$71.29	\$64.16	\$64.16
9200004	4.0mm	each	\$71.29	\$71.29	\$64.16	\$64.16
9200005	5.0mm	each	\$71.29	\$71.29	\$64.16	\$64.16
9200006	6.0mm	each	\$71.29	\$71.29	\$64.16	\$64.16

Tissue Punch Driver

20155001	Thumb knob	each	\$35.49	\$35.49	\$30.17	\$30.17
20155411	Ratchet extension	each	\$49.49	\$49.49	\$44.54	\$44.54

**Complete Set of
Mucosa Tissue Punches**

- Includes 6 Punches (above)
- Plus ratchet extension, thumb knob (see inset) and delrin tray




9200020	set of 6	\$537.99	\$537.99	\$484.19	\$484.19
---------	----------	----------	----------	----------	----------

INSTRUMENTS

A-TITAN CYLINDRICAL TISSUE PUNCH KIT

Tissue Punches are made from specially formulated BioRez, 100% stainless steel, to stay sharper for an extended period of time and are specifically designed for slow speed handpieces.

Includes: Set of 4 Stainless Steel Tissue Punches: 3mm, 4mm, 5mm, 6mm, and a Custom Designed Autoclavable Block and Lid to House all 4 Tissue Punches




K3800C	set of 4	\$288.79	\$288.79	\$259.91	\$259.91
--------	----------	----------	----------	----------	----------

UNIGRAFT KNIFE™ CONNECTIVE TISSUE HARVESTER

For taking precise keratinized epithelial and connective tissue grafts.

Allows for the harvesting of a predictable, uniform, keratinized epithelial graft from the palate to build and ensure adequate marginal or crestal gingiva.



The harvested piece of connective tissue has a consistent 1.5mm thickness and a 9mm width. The length is controlled by the ACE CT Graft Harvester advance and is 10mm for this example.



SINGLE SHOE

DOUBLE SHOE

Unigraft Knife™

The same handle is used for both the Single and Double shoe blades

6500014	Handle only	each	\$232.79	\$232.79	\$209.51	\$209.51
---------	-------------	------	----------	----------	----------	----------


Unigraft Knife™ Single Shoe Blade

6500005	1.5mm x 6mm	each	\$64.29	\$64.29	\$57.86	\$57.86
6500007	1.25mm x 9mm	each	\$64.29	\$64.29	\$57.86	\$57.86
6500008	1.5mm x 9mm	each	\$64.29	\$64.29	\$57.86	\$57.86

Unigraft Knife™ Double Shoe Blade for Connective Tissue Harvesting

6500009	1.5mm x 9mm	each	\$85.79	\$85.79	\$77.21	\$77.21
---------	-------------	------	---------	---------	---------	---------

* Replacement Single and Double Shoes are Interchangeable.



Did you know 99% of ACE SOUTHERN customers receive their order within 1-2 business days of shipment? ACE SOUTHERN has four order fulfillment facilities throughout the United States to ensure your order is **TO YOU IN 2***.

*See page 261 for complete details or visit acesouthern.com/shipping-information



ASPEN SURGICAL SKIN MARKER WITH RULER

- Sterile
- Taper Tip
- Not made with natural rubber latex

100200PDG	each	\$4.29	\$4.08	\$4.08	\$3.86
-----------	------	--------	--------	--------	--------

9

INSTRUMENTS

BLADE HANDLES

#3 Flat Blade Handles

- For surgical blades 10 through 15

0030992	A	5" with metric markings	ACE	each	\$14.49	\$13.77	\$13.77	\$13.04
0030993		5"	ACE	each	\$13.49	\$12.82	\$12.82	\$12.14
371030	B	5" with metric markings	Bard Parker	each	\$25.99	\$24.69	\$24.69	\$23.39
1013003		5" with metric markings	Hu-Friedy®	each	\$34.49	\$32.77	\$32.77	\$31.04
47	C		Integra™	each	\$21.29	\$21.29	\$19.16	\$19.16
47R		Economy, rustless, nickel-plated	Integra™	each	\$9.99	\$9.99	\$8.99	\$8.99
193000	D		J&J Instru.	each	\$4.29	\$4.29	\$3.86	\$3.86

#4 Flat Blade Handles

- For surgical blades 20 through 25

48	E	Extra-fine	Integra™	each	\$21.29	\$21.29	\$19.16	\$19.16
48R		Economy, rustless, nickel-plated	Integra™	each	\$9.99	\$9.99	\$8.99	\$8.99
194000	F		J&J Instru.	each	\$4.29	\$4.29	\$3.86	\$3.86

#5 Round-Style Blade Handles

- For surgical blades 10 through 15

0030995	G	6-1/4" straight	ACE	each	\$35.49	\$35.49	\$31.94	\$31.94
1013005		5-3/4" straight	Hu-Friedy®	each	\$87.99	\$87.99	\$79.19	\$79.19
0039356	H	5-1/2" angled 45°	ACE	each	\$44.79	\$42.55	\$42.55	\$40.31
0039358	I	5-3/4" angled 45° universal <i>Blade may be reversed to approach the surgical site from either direction</i>	ACE	each	\$44.79	\$44.79	\$40.31	\$40.31
0039360		5-1/2" angled 45° reversed <i>For blade placement on the side</i>	ACE	each	\$44.79	\$42.55	\$42.55	\$40.31

#7 Flat Blade Handles

- For surgical blades 10 through 15

0030997	J	6-1/2" straight	ACE	each	\$22.49	\$22.49	\$20.24	\$20.24
49	K	Extra-fine	Integra™	each	\$32.79	\$32.79	\$29.51	\$29.51
197000	L		J&J Instru.	each	\$4.29	\$4.29	\$3.86	\$3.86

MINI BLADE HANDLES

- For miniature surgical blades only

ACE Mini Blade Handles

7861330		4" round knurled		each	\$23.49	\$23.49	\$21.14	\$21.14
7861350	M	5" round knurled		each	\$23.49	\$23.49	\$21.14	\$21.14
7861320		6" round knurled		each	\$23.49	\$23.49	\$21.14	\$21.14

Beaver Mini Blade Handles

374310		4" #3K, round knurled		each	\$64.29	\$61.08	\$61.08	\$57.86
371360		5" #3D, round knurled		each	\$78.49	\$74.57	\$74.57	\$70.64
371350		5" #3KD, round knurled		each	\$71.79	\$68.20	\$68.20	\$64.61
374340		6" #3HL, hexed handle		each	\$104.99	\$99.74	\$99.74	\$94.49
371305		Replacement chuck collet for 371360		each	\$16.49	\$16.49	\$14.84	\$14.84


Integra™ Mini Blade Handles

4400	N	3-3/4" round knurled, with chuck		each	\$56.79	\$56.79	\$51.11	\$51.11
4401	O	3-3/4" round knurled, with self-locking chuck		each	\$57.49	\$57.49	\$51.74	\$51.74



#3 Handles

#4 Handles




#5 Handles

#7 Handles



Mini Handles

ACE DOUBLE BLADE HANDLE FOR SUBEPITHELIAL TISSUE GRAFTS



#3 Flat Handle with Metric Markings

0039362	5", cuts 1.5mm wide	each	\$102.29	\$102.29	\$92.06	\$92.06
---------	---------------------	------	----------	----------	---------	---------

#5 Round Handle

0039372	5-3/4", cuts 1mm wide	each	\$102.29	\$97.18	\$97.18	\$92.06
---------	-----------------------	------	----------	---------	---------	---------

BE THE OFFICE HERO

Save Time and Money with...



REPLENISH

*See page 155 for complete details or visit acesouthern.com/replenish



DISPOSABLE STERILE STAINLESS SCALPELS

Cincinnati Surgical Disposable Scalpels

CS0511	#11 blade	10/box	\$14.29	\$13.58	\$13.58	\$12.86
CS0512	#12 blade	10/box	\$14.29	\$13.58	\$13.58	\$12.86
CS0515	#15 blade	10/box	\$14.29	\$13.58	\$13.58	\$12.86



Integra™ Disposable Safety Scalpels

- Complies with OSHA directive for the Needlestick Safety and Prevention Act
- Securely locks blade into both extended and retracted position
- Simple, retractable mechanism provides precision and protection
- Reliable sharpness of the Miltex stainless steel blade
- Metric ruler assists in determining length of incision

4510	#10 blade	10/box	\$34.49	\$32.77	\$32.77	\$31.04
4511	#11 blade	10/box	\$34.49	\$32.77	\$32.77	\$31.04
4515	#15 blade	10/box	\$34.49	\$32.77	\$32.77	\$31.04
4515C	#15C blade	10/box	\$34.49	\$32.77	\$32.77	\$31.04



Mycro Reli®-Cut Plus Disposable Scalpels

- Top activation allows for ambidextrous use
- Audible "clicks" help confirm proper blade position
- Fully extended blade enables clear line of sight
- Non-slip, ergonomic handle with integrated metric scale
- Scalpel blades produced from high-quality Swiss stainless steel

6008TR10	#10 blade	10/box	\$12.79	\$12.15	\$12.15	\$11.51
6008TR11	#11 blade	10/box	\$12.79	\$12.15	\$12.15	\$11.51
6008TR12	#12 blade	10/box	\$12.79	\$12.15	\$12.15	\$11.51
6008TR15	#15 blade	10/box	\$12.79	\$12.15	\$12.15	\$11.51
6008TR20	#20 blade	10/box	\$12.79	\$12.15	\$12.15	\$11.51



Swann-Morton Disposable Scalpels

SM0504	#12 blade	10/box	\$15.29	\$14.53	\$14.53	\$13.76
SM0518	#12D blade	10/box	\$15.29	\$14.53	\$14.53	\$13.76
SM0505	#15 blade	10/box	\$15.29	\$14.53	\$14.53	\$13.76
SM0592	#15T blade	10/box	\$15.29	\$14.53	\$14.53	\$13.76



INSTRUMENTS

**STERILE SURGICAL BLADES
CARBON STEEL**

ACE Carbon Steel

0036112	#12	100/box	\$28.49	\$27.07	\$27.07	\$25.64
0036312	#12D	100/box	\$38.29	\$36.38	\$36.38	\$34.46
0036115	#15*	100/box	\$24.79	\$23.55	\$23.55	\$22.31
0036315	#15C	100/box	\$27.49	\$26.12	\$26.12	\$24.74

* the tip of the 15* blade is slightly different than a standard 15

Bard Parker Carbon Steel

BP1110	#10	50/box	\$39.49	\$37.52	\$37.52	\$35.54
BP1111	#11	50/box	\$39.49	\$37.52	\$37.52	\$35.54
BP1112	#12	50/box	\$127.79	\$121.40	\$121.40	\$115.01
BP1115	#15	50/box	\$39.49	\$37.52	\$37.52	\$35.54

Cincinnati Surgical Carbon Steel

CS0012	#12	100/box	\$43.29	\$41.13	\$41.13	\$38.96
CS0015	#15	100/box	\$43.29	\$41.13	\$41.13	\$38.96

Hu-Friedy® Carbon Steel

1025515	#15	100/box	\$144.29	\$137.08	\$137.08	\$129.86
1025615	#15C	100/box	\$215.29	\$204.53	\$204.53	\$193.76

Miltex Carbon Steel

4112B	#12B	100/box	\$101.29	\$96.23	\$96.23	\$91.16
4115C	#15C	100/box	\$103.29	\$98.13	\$98.13	\$92.96

MYCO GLASSVAN® Carbon Steel


2001T11	#11	100/box	\$20.29	\$19.28	\$19.28	\$18.26
2001T12	#12	100/box	\$20.29	\$19.28	\$19.28	\$18.26
2001T12B	#12B	100/box	\$20.29	\$19.28	\$19.28	\$18.26
2001T15	#15	100/box	\$20.29	\$19.28	\$19.28	\$18.26
2001T15C	#15C	100/box	\$20.29	\$19.28	\$19.28	\$18.26

Propper Carbon Steel

122012	#12 Carbon Steel	150/box	\$78.49	\$74.57	\$74.57	\$70.64
122015	#15 Carbon Steel	150/box	\$75.79	\$72.00	\$72.00	\$68.21

Swann-Morton Carbon Steel

00SM10	#10	100/box	\$45.29	\$43.03	\$43.03	\$40.76
00SM11	#11	100/box	\$45.29	\$43.03	\$43.03	\$40.76
00SM12	#12	100/box	\$45.29	\$43.03	\$43.03	\$40.76
00SM12D	#12D	100/box	\$53.29	\$50.63	\$50.63	\$47.96
00SM15	#15	100/box	\$45.29	\$43.03	\$43.03	\$40.76
00SM15C	#15C	100/box	\$53.29	\$50.63	\$50.63	\$47.96
00SM15T	#15T	100/box	\$45.29	\$43.03	\$43.03	\$40.76



**STERILE SURGICAL BLADES
STAINLESS STEEL**

ACE Stainless Steel

0037215	#15*	100/box	\$27.49	\$26.12	\$26.12	\$24.74
0037315	#15C	100/box	\$27.49	\$26.12	\$26.12	\$24.74

* the tip of the 15* blade is slightly different than a standard 15

Bard Parker Stainless Steel

BP1210	#10	50/box	\$43.49	\$41.32	\$41.32	\$39.14
1211	#11	50/box	\$43.49	\$41.32	\$41.32	\$39.14
BP1212	#12	50/box	\$154.99	\$147.24	\$147.24	\$139.49
BP1712	#12B	50/box	\$283.99	\$269.79	\$269.79	\$255.59
BP1215	#15	50/box	\$37.49	\$35.62	\$35.62	\$33.74

Cincinnati Surgical Stainless Steel

CS0111	#11	100/box	\$43.29	\$41.13	\$41.13	\$38.96
CS0112	#12	100/box	\$43.29	\$41.13	\$41.13	\$38.96
CS0112D	#12D	100/box	\$52.29	\$49.68	\$49.68	\$47.06
CS0115	#15	100/box	\$43.29	\$41.13	\$41.13	\$38.96
CS0115C	#15C	100/box	\$52.29	\$49.68	\$49.68	\$47.06

Miltex Stainless Steel

4312B	#12B	100/box	\$101.29	\$96.23	\$96.23	\$91.16
4315	#15	100/box	\$52.79	\$50.15	\$50.15	\$47.51

MYCO GLASSVAN® Stainless Steel

3001T11	#11	100/box	\$20.29	\$19.28	\$19.28	\$18.26
3001T12	#12	100/box	\$20.29	\$19.28	\$19.28	\$18.26
3001T12B	#12B	100/box	\$20.29	\$19.28	\$19.28	\$18.26
3001T15	#15	100/box	\$20.29	\$19.28	\$19.28	\$18.26

Swann-Morton Stainless Steel

01SM10	#10	100/box	\$45.29	\$43.03	\$43.03	\$40.76
01SM11	#11	100/box	\$45.29	\$43.03	\$43.03	\$40.76
01SM12	#12	100/box	\$45.29	\$43.03	\$43.03	\$40.76
01SM12D	#12D	100/box	\$53.29	\$50.63	\$50.63	\$47.96
01SM15	#15	100/box	\$45.29	\$43.03	\$43.03	\$40.76
01SM15C	#15C	100/box	\$53.29	\$50.63	\$50.63	\$47.96
01SM15T	#15T	100/box	\$45.29	\$43.03	\$43.03	\$40.76

INSTRUMENTS



Beaver Threaded Mini Blades



Beaver threaded Mini Blades are screwed onto a mini blade handle (with the end clamp removed) while still in the protective tray.

Continue to thread handle until the thread is fully engaged and the Blade is released from the tray clips.

STERILE MINI BLADES
STAINLESS STEEL

ACE Mini Stainless Steel

7866200	#62	10/pkg	\$25.99	\$24.69	\$24.69	\$23.39
7866400	#64	10/pkg	\$25.99	\$24.69	\$24.69	\$23.39
7866500	#65	10/pkg	\$25.99	\$24.69	\$24.69	\$23.39
7866700	#67	10/pkg	\$25.99	\$24.69	\$24.69	\$23.39
7866900	#69	10/pkg	\$56.29	\$53.48	\$53.48	\$50.66

Beaver Mini Stainless Steel

376300	#63	6/pkg	\$146.49	\$139.17	\$139.17	\$131.84
BEAVER6400	#64	20/pkg	\$82.49	\$78.37	\$78.37	\$74.24
BEAVER6700	#67	20/pkg	\$83.79	\$79.60	\$79.60	\$75.41
BEAVER6900	#69	20/pkg	\$272.99	\$259.34	\$259.34	\$245.69
BEAVER6910	#6910 Cleft Palate 60° Angle, 180° Cutting Surface	6/pkg	\$152.49	\$144.87	\$144.87	\$137.24

Swann-Morton Mini Stainless Steel

01SM63	#63	25/box	\$97.29	\$92.43	\$92.43	\$87.56
01SM64	#64	25/box	\$60.49	\$57.47	\$57.47	\$54.44
01SM67	#67	25/box	\$60.49	\$57.47	\$57.47	\$54.44
01SM69	#69	25/box	\$97.29	\$92.43	\$92.43	\$87.56



HU-FRIEDY® SCALPEL
BLADE REMOVER

- For use with all blade shapes
- Reusable
- Allows simple, safe removal of blade from scalpel handle

1019900	2-1/4" (5.5cm)	each	\$31.29	\$31.29	\$28.16	\$28.16
---------	----------------	------	---------	---------	---------	---------



CINCINNATI SURGICAL
AUNEX BLADE REMOVER

- The only blade remover designed to remove ALL standard sized surgical blades, and is the first blade remover for size 60 and 70 blades. Remove used blades with a single-handed technique. No fingers, forceps or re-sheathing and NO INJURIES.

48CS	each	\$26.29	\$26.29	\$23.66	\$23.66
------	------	---------	---------	---------	---------



MYCO
QUICKSMART®
BLADEFLASK
SURGICAL BLADE
REMOVAL SYSTEM

Adheres to OSHA Bloodborne Pathogen Safety Standard "1910.1030"

- NO OVER FILLING – Quicksmart has a counter and shut-off mechanism when the container reaches its 100-blade capacity
- SAFE – Quicksmart allows single-handed safe scalpel blade removal and disposal, avoiding the risk and cost of injury and infection
- CONVENIENT –The secure sharps container can be mounted to a wall, trolley, or bench via a universal bracket
- EASY –The scalpel is inserted into the slot in the container until the click is heard; the handle is withdrawn, and the blade remains sealed inside for later disposal
- COMPATIBLE –Quicksmart has been tested and is suitable for use with blades using the following most commonly used handles: no. 3, no. 4, no. 5, no. 7 and no. 9

QFRGEN	Blade Remover	each	\$39.99	\$39.99	\$35.99	\$35.99
QBRGEN	Wall mounting bracket	each	\$11.99	\$11.99	\$10.79	\$10.79

9

Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INSTRUMENTS



ELEVATORS

ACE

6300004	A	#301 Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300008	B	#303 Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300018	D	#8 Crane Pick	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300011		#46R Heidbrink Elevator	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300059		#71 Miller Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300060		#72 Miller Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300061		#73 Miller Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300062		#74 Miller Apexo	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300014		#34 Seldin	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6300140	C	#34S Seldin	6" each	\$46.29	\$46.29	\$41.66	\$41.66
6315015		#34S Seldin, Serrated	6" each	\$46.29	\$46.29	\$41.66	\$41.66

A-Titan

K5000		Bone Preservation Elevator Kit , includes the next 4 items	each	\$291.79	\$291.79	\$262.61	\$262.61
63	E	Spade Tip, Straight	each	\$79.79	\$79.79	\$71.81	\$71.81
60B	F	Spade Tip, Angled	each	\$79.79	\$79.79	\$71.81	\$71.81
60R	G	Spade Tip, Right Molar	each	\$79.79	\$79.79	\$71.81	\$71.81
ATT60L	H	Spade Tip, Left Molar	each	\$79.79	\$79.79	\$71.81	\$71.81

K3000		Molar Elevator Kit , includes the next 2 items	each	\$149.29	\$149.29	\$134.36	\$134.36
60R	G	Spade Tip, Right Molar	each	\$79.79	\$79.79	\$71.81	\$71.81
ATT60L	H	Spade Tip, Left Molar	each	\$79.79	\$79.79	\$71.81	\$71.81

60U	I	Spade Tip, Back Action	each	\$79.79	\$79.79	\$71.81	\$71.81
ATT301	J	Small, Straight - 2.5mm	each	\$79.79	\$79.79	\$71.81	\$71.81
304	K	Medium, Straight - 3.5mm	each	\$79.79	\$79.79	\$71.81	\$71.81
34S	L	Large, Straight - 3.75mm	each	\$79.79	\$79.79	\$71.81	\$71.81
34ST	M	Large, Straight, Thin - 4.5mm	each	\$79.79	\$79.79	\$71.81	\$71.81
77R	N	Angled, Molar Region	each	\$79.79	\$79.79	\$71.81	\$71.81
46R	O	Angled, Serrated Molar Region	each	\$79.79	\$79.79	\$71.81	\$71.81

Hu-Friedy®

E77R	P	#77R, 3mm, serrated	each	\$94.29	\$94.29	\$84.86	\$84.86
E190	Q	#190, 3mm	each	\$94.29	\$94.29	\$84.86	\$84.86
E191	R	#191, 3mm	each	\$94.29	\$94.29	\$84.86	\$84.86
EA	S	A, Cogswell	each	\$94.29	\$94.29	\$84.86	\$84.86
EB	T	B, Cogswell	each	\$94.29	\$94.29	\$84.86	\$84.86
EA30	U	A30, Cryer	each	\$94.29	\$94.29	\$84.86	\$84.86
EA31	V	A31, Cryer	each	\$94.29	\$94.29	\$84.86	\$84.86

Integra™

DEL301	W	#301 Apical	each	\$79.79	\$79.79	\$71.81	\$71.81
DEL34	X	#34 Apical	each	\$79.79	\$79.79	\$71.81	\$71.81
DEL46	Y	#46 Apical	each	\$79.79	\$79.79	\$71.81	\$71.81
DEL77R	Z	#77R, serrated	each	\$76.49	\$76.49	\$68.84	\$68.84
DEL44	1	#44 Cryer	each	\$83.79	\$83.79	\$75.41	\$75.41
DEL45	2	#45 Cryer	each	\$83.79	\$83.79	\$75.41	\$75.41

J&J Instruments

04520	3	#1L Seldin	each	\$19.29	\$19.29	\$17.36	\$17.36
04540	4	#1R Seldin	each	\$19.29	\$19.29	\$17.36	\$17.36
04990	5	#190 General	each	\$19.29	\$19.29	\$17.36	\$17.36
04991	6	#191 General	each	\$19.29	\$19.29	\$17.36	\$17.36

LUXATING ELEVATORS

A-Titan D-Lux

K5500A		Curved Kit , includes the 4 curved X-Otomes below	each	\$384.49	\$384.49	\$346.04	\$346.04
LUX2A	A	2mm curved A-T 1	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX3A	B	3mm curved A-T 2	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX4A	C	4mm curved A-T 3	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX5A	D	5mm curved A-T 4	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX2	E	2mm straight A-T 5	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX3	F	3mm straight A-T 6	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX4	G	4mm straight A-T 7	each	\$115.29	\$115.29	\$103.76	\$103.76
LUX5	H	5mm straight A-T 8	each	\$115.29	\$115.29	\$103.76	\$103.76
K7800	I	Spade Kit , includes the next 5 X-Otomes below and a sterilization case (T10)	each	\$609.49	\$609.49	\$548.54	\$548.54
XO77		Spade, beveled, straight	each	\$126.79	\$126.79	\$114.11	\$114.11
XO77M		Spade, beveled, angled mesial	each	\$126.79	\$126.79	\$114.11	\$114.11
XO77D		Spade, beveled, angled distal	each	\$126.79	\$126.79	\$114.11	\$114.11
XO77AL		Spade, beveled, angled left	each	\$126.79	\$126.79	\$114.11	\$114.11
XO77AR		Spade, beveled, angled right	each	\$126.79	\$126.79	\$114.11	\$114.11
ATTT10		Sterilization case	each	\$129.29	\$129.29	\$116.36	\$116.36

Hu-Friedy® Black Line

EL2SX	J	2mm, straight	each	\$113.79	\$113.79	\$102.41	\$102.41
EL2CX		2mm, curved	each	\$113.79	\$113.79	\$102.41	\$102.41
EL3SX		3mm, straight	each	\$113.79	\$113.79	\$102.41	\$102.41
EL3CX		3mm, curved	each	\$113.79	\$113.79	\$102.41	\$102.41
EL4SX		4mm, straight	each	\$113.79	\$113.79	\$102.41	\$102.41
EL4CX		4mm, curved	each	\$113.79	\$113.79	\$102.41	\$102.41
EL5SX		5mm, straight	each	\$113.79	\$113.79	\$102.41	\$102.41
EL5CX		5mm, curved	each	\$113.79	\$113.79	\$102.41	\$102.41



Hu-Friedy®

EL2C		L2C, 2mm, #510 handle	each	\$94.29	\$94.29	\$84.86	\$84.86
EL3C	K	L3C, 3mm, #510 handle	each	\$94.29	\$94.29	\$84.86	\$84.86

Integra™

DELLUXS2	L	2mm, straight tip	each	\$69.79	\$69.79	\$62.81	\$62.81
DELLUXC3	M	3mm, curved tip	each	\$69.79	\$69.79	\$62.81	\$62.81
DELLUXS3	N	3mm, straight tip	each	\$69.79	\$69.79	\$62.81	\$62.81
DELLUXC4	O	4mm, curved tip	each	\$69.79	\$69.79	\$62.81	\$62.81
DELLUXS4	P	4mm, straight tip	each	\$69.79	\$69.79	\$62.81	\$62.81

INSTRUMENTS



9

Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price

INSTRUMENTS



ROOT TIP ELEVATORS / PICKS

ACE

6300113	#1 Heidbrink	6-½"	each	\$29.49	\$29.49	\$26.54	\$26.54
6300114	A #2 Heidbrink	6-½"	each	\$29.49	\$29.49	\$26.54	\$26.54
6300115	B #3 Heidbrink	6-½"	each	\$29.49	\$29.49	\$26.54	\$26.54
6300118	#13/14 Heidbrink	7-½"	each	\$25.99	\$25.99	\$23.39	\$23.39
6315130	C Double-ended, serrated	7-¾"	each	\$28.49	\$28.49	\$25.64	\$25.64

Integra™

DELH1	D #1		each	\$62.29	\$62.29	\$56.06	\$56.06
DELH2	E #2		each	\$62.29	\$62.29	\$56.06	\$56.06
DELH3	F #3		each	\$62.29	\$62.29	\$56.06	\$56.06
DELH23	G #2-3		each	\$62.29	\$62.29	\$56.06	\$56.06
DEL78	H #78		each	\$81.99	\$81.99	\$73.79	\$73.79
DEL79	I #79		each	\$81.99	\$81.99	\$73.79	\$73.79
DEL80	J #80		each	\$81.99	\$81.99	\$73.79	\$73.79

T-BAR HANDLE ELEVATORS

ACE

6300063	#6X Potts, T-bar handle	6"	each	\$46.29	\$46.29	\$41.66	\$41.66
6300064	#7X Potts, T-bar handle	6"	each	\$46.29	\$46.29	\$41.66	\$41.66

Hu-Friedy®

E11LX	L 11L Winter cross-bar handle		each	\$94.29	\$94.29	\$84.86	\$84.86
E11RX	K 11R Winter cross-bar handle		each	\$94.29	\$94.29	\$84.86	\$84.86

Integra™

DEL6X	M #6X Potts, T-bar handle		each	\$65.29	\$65.29	\$58.76	\$58.76
DEL7X	N #7X Potts, T-bar handle		each	\$65.29	\$65.29	\$58.76	\$58.76

Integra™

04591	O #1 cross-bar handle		each	\$20.49	\$20.49	\$18.44	\$18.44
04592	#2 cross-bar handle		each	\$20.49	\$20.49	\$18.44	\$18.44

INSTRUMENTS

PERIOSTEAL ELEVATORS

ACE

0300080	A	#9A Allen	6-½" each	\$38.29	\$38.29	\$34.46	\$34.46
0300139	B	Buser	7" each	\$39.49	\$39.49	\$35.54	\$35.54
0311023	C	Freer, strong curve	7-½" each	\$49.99	\$49.99	\$44.99	\$44.99
0300022	D	#15, Freer, curved ends	8" each	\$46.29	\$46.29	\$41.66	\$41.66
0300023	E	#16, Freer, strong curved ends	8" each	\$46.29	\$46.29	\$41.66	\$41.66
0300073	F	Molt #9, double ended	7" each	\$38.29	\$38.29	\$34.46	\$34.46
0314042		Molt, double ended, SB	7-½" each	\$53.79	\$53.79	\$48.41	\$48.41
0300072	G	Molt #9, double ended, titanium	7-½" each	\$134.79	\$134.79	\$121.31	\$121.31
0300001		#23 Seldin, flat handle, double ended	7-¾" each	\$53.79	\$53.79	\$48.41	\$48.41
0300004	H	#23 Seldin, octagonal handle, double ended	7-½" each	\$54.79	\$54.79	\$49.31	\$49.31
0300065	I	Pritchard, double ended	7-¾" each	\$53.79	\$53.79	\$48.41	\$48.41
0300016	J	#24G Raspatory straight, double ended	6-¾" each	\$58.29	\$58.29	\$52.46	\$52.46

A-Titan

G24	K	#24 cutting edge, A.Titan handle	each	\$75.99	\$75.99	\$68.39	\$68.39
BUSER	L	Buser, double ended	each	\$90.79	\$90.79	\$81.71	\$81.71

Hu-Friedy® Black Line

PPAELX	M	Allen	each	\$124.29	\$124.29	\$111.86	\$111.86
PPAELSX		Allen, short handle	each	\$124.29	\$124.29	\$111.86	\$111.86
PPAELAX		Allen, anterior	each	\$124.29	\$124.29	\$111.86	\$111.86
PPAELASX		Allen, anterior, short handle	each	\$124.29	\$124.29	\$111.86	\$111.86
PPAELPX		Allen, posterior	each	\$127.29	\$127.29	\$114.56	\$114.56
PPBUSERX	N	Buser	each	\$95.29	\$95.29	\$85.76	\$85.76
PPBUSERSX		Buser, short handle	each	\$95.29	\$95.29	\$85.76	\$85.76
PKN152X	O	#152 Kramer-Nevis	each	\$95.29	\$95.29	\$85.76	\$85.76
PKN152SX		#152 Kramer-Nevis, short handle	each	\$95.29	\$95.29	\$85.76	\$85.76
P9X	P	#9 Molt	each	\$95.29	\$95.29	\$85.76	\$85.76
PPR3X	Q	#3 Prichard	each	\$108.99	\$108.99	\$98.09	\$98.09

Hu-Friedy®

P23	R	#23 Seldin	each	\$105.79	\$105.79	\$95.21	\$95.21
PPAELP	S	Allen, posterior, #6 handle	each	\$101.79	\$101.79	\$91.61	\$91.61

Integra™

DELM9	T	#9 Molt	each	\$79.79	\$79.79	\$71.81	\$71.81
DELS23	U	#23 Seldin	each	\$88.79	\$88.79	\$79.91	\$79.91
DELPR3	V	#3 Pritchard	each	\$79.29	\$79.29	\$71.36	\$71.36

J&J Instruments

04120	W	#9 Molt	each	\$16.99	\$16.99	\$15.29	\$15.29
04160	X	#23 Seldin	each	\$16.99	\$16.99	\$15.29	\$15.29




9

INSTRUMENTS



INTERCHANGEABLE TIP PERIOTOMES

ACE PERIOTOME HANDLE & BLADES KIT



- Periotome is used to remove ligament attachment from the root of the tooth and the alveolar housing.
- Single end handle can be tapped with a mallet or used with hand pressure.
- Includes 6-Blade Storage Holder: stands up for easy access during procedures and folds flat for sterilization and storage.
- Components can be purchased separately.


3441770	A	Periotome Kit, includes the next 9 items	kit	\$379.29	\$379.29	\$341.36	\$341.36
16901001	B	Single-end handle only, 6" (15cm)	each	\$135.79	\$135.79	\$122.21	\$122.21
3441651	C	P3, 1.5mm, angled, 1.75" (4.5cm)	each	\$51.99	\$51.99	\$46.79	\$46.79
3441631	D	P1, 2mm, straight, 1.5" (4cm)	each	\$51.99	\$51.99	\$46.79	\$46.79
3441641	F	P2, 1.5mm, straight, 1.5" (4cm)	each	\$51.99	\$51.99	\$46.79	\$46.79
3441691	G	2.5mm, straight, semi rigid, 1.5" (4cm)	each	\$51.99	\$51.99	\$46.79	\$46.79
3441661	H	Arrow tip	each	\$51.99	\$51.99	\$46.79	\$46.79
6400312S		Surgical mallet, 6-½" (16.5cm)	each	\$81.29	\$81.29	\$73.16	\$73.16
3441665	I	Blade storage holder, 6" (15.2cm)	each	\$62.49	\$62.49	\$56.24	\$56.24
9400036		Instrument tray, 175 x 140 x 24mm	each	\$135.79	\$135.79	\$122.21	\$122.21

A-TITAN INTERCHANGEABLE PERIOTOME KIT

K6100		Interchangeable Periotome Kit, includes the next 4 items	kit	\$271.99	\$271.99	\$244.79	\$244.79
P0	J	Single-end periotome handle	each	\$152.79	\$152.79	\$137.51	\$137.51
P1	K	Periotome Tip - Wide	each	\$57.49	\$57.49	\$51.74	\$51.74
P2	L	Periotome Tip - Narrow	each	\$57.49	\$57.49	\$51.74	\$51.74
P3	M	Periotome Tip - Angled	each	\$67.49	\$67.49	\$60.74	\$60.74


HU-FRIEDY® INTERCHANGEABLE PERIOTOME KIT

PTIKIT		Interchangeable Periotome Kit, includes the next 6 items	kit	\$325.29	\$325.29	\$292.76	\$292.76
PTIH	N	Interchangeable Periotome Handle	each	\$190.79	\$190.79	\$171.71	\$171.71
PTI7		Tip, Straight Wide	each	\$67.49	\$67.49	\$60.74	\$60.74
PTI8		Tip, Straight Narrow	each	\$67.49	\$67.49	\$60.74	\$60.74
PTI9	O	Tip, Angulated	each	\$67.49	\$67.49	\$60.74	\$60.74
PTI10		Tip, Mesial/Distal Angulated	each	\$67.49	\$67.49	\$60.74	\$60.74
IMS1372TH		IMS Bur Cushion w/lid, holds 6	each	\$32.29	\$32.29	\$29.06	\$29.06



INSTRUMENTS

FIXED TIP PERIOTOMES



ACE PERIODONTAL LIGAMENT/LUXATING ELEVATOR KIT

- Single end handle can be tapped with a mallet.
- Instruments and tray are stainless steel.

4547111	A	Periodontal Ligament/Luxator Elevator Kit , includes the next 7 items	kit	\$699.99	\$699.99	\$629.99	\$629.99
4547701	B	Elevator, brown, angled "OUT"	each	\$114.29	\$114.29	\$102.86	\$102.86
4547901	C	Elevator, green, curved "OUT"	each	\$114.29	\$114.29	\$102.86	\$102.86
4547401	D	Elevator, blue, straight 2.5mm	each	\$114.29	\$114.29	\$102.86	\$102.86
4547501	E	Elevator, red, arrow tip	each	\$114.29	\$114.29	\$102.86	\$102.86
4547601	F	Elevator, gold, angled "IN"	each	\$114.29	\$114.29	\$102.86	\$102.86
6400312S		Surgical mallet, 6-1/2" (16.5cm)	each	\$81.29	\$81.29	\$73.16	\$73.16
9400036		Instrument tray, 175 x 140 x 24mm	each	\$135.79	\$135.79	\$122.21	\$122.21

A-TITAN X-OTOMES KIT

K6600		X-Otomes Kit Periotope Kit includes the next 7 items	kit	\$848.49	\$848.49	\$763.64	\$763.64
X010	G	Straight small	each	\$126.79	\$126.79	\$114.11	\$114.11
X020	H	Curved small	each	\$126.79	\$126.79	\$114.11	\$114.11
X030	I	Curved angled mesial	each	\$126.79	\$126.79	\$114.11	\$114.11
X040	J	Angled distal extra fine	each	\$126.79	\$126.79	\$114.11	\$114.11
X050	K	Angled mesial	each	\$126.79	\$126.79	\$114.11	\$114.11
X060	L	Angled distal	each	\$126.79	\$126.79	\$114.11	\$114.11
X077	M	Straight beveled spade	each	\$126.79	\$126.79	\$114.11	\$114.11

Additional X-Otomes available separately

X070	N	Straight universal	each	\$126.79	\$126.79	\$114.11	\$114.11
X077AR	O	Straight beveled angled right	each	\$126.79	\$126.79	\$114.11	\$114.11
X077AL	P	Straight beveled angled left	each	\$126.79	\$126.79	\$114.11	\$114.11
X077D	Q	Angled distal w/beveled spade tip	each	\$126.79	\$126.79	\$114.11	\$114.11
X077M	R	Angled mesial w/beveled spade tip	each	\$126.79	\$126.79	\$114.11	\$114.11
X080	S	Straight extra fine	each	\$126.79	\$126.79	\$114.11	\$114.11
X090	T	Angled mesial large	each	\$126.79	\$126.79	\$114.11	\$114.11
X0100	U	Angled distal large	each	\$126.79	\$126.79	\$114.11	\$114.11
X0110	V	Angled 2mm	each	\$126.79	\$126.79	\$114.11	\$114.11

A-TITAN FIXED-TIP PERIOTOME KIT

- Hybrid instrumentation, a cross between an elevator & periotope
- Custom sure-grip handle for power transfer during extraction
- Honed narrow tips to sever periodontal fibers
- Preserve tissue structures; eliminate bone fracture
- 100% stainless steel; made in the U.S.A.

K6000		Fixed-Tip Periotope Kit , includes the next 7 items	kit	\$215.79	\$215.79	\$194.21	\$194.21
P100	W	Wide blade - universal	each	\$89.29	\$89.29	\$80.36	\$80.36
P200	X	Narrow blade - anterior	each	\$89.29	\$89.29	\$80.36	\$80.36
P300	Y	Angled blade - posterior	each	\$89.29	\$89.29	\$80.36	\$80.36

Additional X-Otomes available separately

P500	Z	Angled blade - interproximal	each	\$89.29	\$89.29	\$80.36	\$80.36
------	----------	------------------------------	------	---------	---------	---------	---------



9

Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price

INSTRUMENTS



FIXED TIP PERIOTOMES - continued

HU-FRIEDY® LUXATING HYBRID KIT

Designed with thin, precisely angled tips to aid in preserving the buccal plate during atraumatic extractions.

- Improved apical and interproximal access with thin working tips
- Easily severs the periodontal ligament with finely honed, sharp cutting edges
- Enhanced oral adaptation with 7 precisely manufactured tip options
- Longevity and corrosion resistance afforded by handcrafting from the finest grade surgical stainless steel
- With #6 handle

EPTKIT	A	Kit with the next 9 items	kit	\$1,047.29	\$1,047.29	\$942.56	\$942.56
EPTSSP	B	straight spade	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTSMS	C	small straight	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTSMC	D	small curved	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTSMD	E	small distal	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTSMF	F	small mesial	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTLGD	G	large distal	each	\$142.79	\$142.79	\$128.51	\$128.51
EPTLGM	H	large mesial	each	\$142.79	\$142.79	\$128.51	\$128.51
SSLUX		Ceramic sharpening stone, fine grit	each	\$35.79	\$35.79	\$32.21	\$32.21
IM5080		Signature series cassette, 8-instrument, small, red	each	\$145.79	\$145.79	\$131.21	\$131.21

HU-FRIEDY® PERIOTOMES

PT1	I	Posterior, double-ended	each	\$87.29	\$87.29	\$78.56	\$78.56
PT2		Anterior, double-ended	each	\$87.29	\$87.29	\$78.56	\$78.56
PT3		Posterior, right, single-ended	each	\$87.29	\$87.29	\$78.56	\$78.56
PT4		Posterior, left, single-ended	each	\$87.29	\$87.29	\$78.56	\$78.56
PT5		Anterior, angled, single-ended	each	\$87.29	\$87.29	\$78.56	\$78.56
PT6		Anterior, straight, single-ended	each	\$87.29	\$87.29	\$78.56	\$78.56

HU-FRIEDY® BLACK LINE PERIOTOMES

PT1X	J	Posterior, Double-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PT2X		Anterior, Double-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PT3X		Right, Single-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PT4X		Left, Single-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PT5X		Anterior Angled, Single-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PT6X		Anterior Straight, Single-Ended	each	\$108.99	\$108.99	\$98.09	\$98.09
PEREURX		Wise	each	\$108.99	\$108.99	\$98.09	\$98.09

PERIODONTAL KNIVES

HU-FRIEDY® BLACK LINE PERIOTOMES

TKN1X	L	Tunneling 1	each	\$124.29	\$124.29	\$111.86	\$111.86
TKN1SX		Tunneling 1, Short Handle	each	\$124.29	\$124.29	\$111.86	\$111.86
TKN2X	M	Tunneling 2	each	\$124.29	\$124.29	\$111.86	\$111.86
TKN2SX		Tunneling 2, Short Handle	each	\$124.29	\$124.29	\$111.86	\$111.86
K012KP03AX	N	1/2 Allen Orban	each	\$102.29	\$102.29	\$92.06	\$92.06
K012KP03RX		1/2 Allen Orban, Round	each	\$102.29	\$102.29	\$92.06	\$92.06
KK1516X		15/16 Kirkland	each	\$102.29	\$102.29	\$92.06	\$92.06
KK1516SX		15/16 Kirkland, short handle	each	\$102.29	\$102.29	\$92.06	\$92.06
KPAX	O	Allen End-Cutting Intrasulcular	each	\$102.29	\$102.29	\$92.06	\$92.06
K012X		1/2 Orban	each	\$102.29	\$102.29	\$92.06	\$92.06
KGf11X	P	11 Goldman-Fox	each	\$102.29	\$102.29	\$92.06	\$92.06

INSTRUMENTS

PERIO & ORAL SURGERY CURETTES

DE = double ended

COLUMBIA

SC4R4L9E2	A	#4R/4L Everedge® 2.0, #9 handle, DE	Hu-Friedy®	\$46.79	\$46.79	\$42.11	\$42.11
SC13149E2	B	#13/14, Everedge® 2.0, #9 handle, DE	Hu-Friedy®	\$46.79	\$46.79	\$42.11	\$42.11

LUCAS

612	C	#75, straight, DE	Integra™	\$73.99	\$73.99	\$66.59	\$66.59
1000209		#85	7" ACE	\$29.49	\$29.49	\$26.54	\$26.54
614		#85, angular, small, DE	7" Integra™	\$73.99	\$73.99	\$66.59	\$66.59
1000210	D	#86	7" ACE	\$29.49	\$29.49	\$26.54	\$26.54
616	E	#86, angular, medium, DE	Integra™	\$73.99	\$73.99	\$66.59	\$66.59
CL86X	F	#86, Black Line, DE	Hu-Friedy®	\$92.29	\$92.29	\$83.06	\$83.06
1000211		#87	7" ACE	\$29.49	\$29.49	\$26.54	\$26.54
CL87X		#87, Black Line, DE	Hu-Friedy®	\$92.29	\$92.29	\$83.06	\$83.06
618		#87, angular, large, DE	Integra™	\$73.99	\$73.99	\$66.59	\$66.59
1000212		#88	7" ACE	\$29.49	\$29.49	\$26.54	\$26.54
1000224	G	Angled, 3.5mm, serrated 5-teeth	6-½" ACE	\$53.79	\$53.79	\$48.41	\$48.41
1000226	H	Angled, 4.5mm, serrated 7-teeth	6-½" ACE	\$53.79	\$53.79	\$48.41	\$48.41

MCCALL

SM1314S9E2	9	#13S/14S McCall Curette		\$60.29	\$60.29	\$54.26	\$54.26
------------	----------	-------------------------	--	---------	---------	---------	---------

MILLER

1000203	J	#9	7-¼" ACE	\$28.49	\$28.49	\$25.64	\$25.64
1000204	K	#10	7-¼" ACE	\$28.49	\$28.49	\$25.64	\$25.64
1000205	L	#11	7" ACE	\$28.49	\$28.49	\$25.64	\$25.64
1000206	M	#12	7" ACE	\$28.49	\$28.49	\$25.64	\$25.64

MOLT

1000218		#1, single ended	6" ACE	\$35.99	\$35.99	\$32.39	\$32.39
1000219	N	#2, single ended	6" ACE	\$35.99	\$35.99	\$32.39	\$32.39
MIL6112	O	#2	Integra™	\$94.29	\$94.29	\$84.86	\$84.86
02292	P	#2, straight, medium	J&J Instru.	\$18.79	\$18.79	\$16.91	\$16.91
1000221	Q	#2/4, DE	7" ACE	\$53.79	\$53.79	\$48.41	\$48.41
CM24X	R	#2/4, Black Line, DE	Hu-Friedy®	\$118.79	\$118.79	\$106.91	\$106.91
1000220		#4, single ended	6" ACE	\$35.99	\$35.99	\$32.39	\$32.39
CM4	S	#4, single ended, #502 handle	Hu-Friedy®	\$74.79	\$74.79	\$67.31	\$67.31
02294	T	#4	J&J Instru.	\$18.79	\$18.79	\$16.91	\$16.91

SHEPHERD SINUS LIFT

1000248	U	Shepherd Sinus Lift, DE		\$49.99	\$49.99	\$44.99	\$44.99
---------	----------	-------------------------	--	---------	---------	---------	---------



9

For a full selection of scalers and curettes visit www.acesouthern.com

Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

INSTRUMENTS



SCALERS - IMPLANTS

A-TITAN TITANIUM IMPLANT SCALER KIT

Implant maintenance instrumentation with color coded handles. Sterilizable titanium tips. Heat treated and fully autoclavable. Made from 6AL-4V titanium alloy, the same material as implants/abutments which eliminates damage during cleaning.

- Guaranteed sharpness
- Made in the U.S.A.

K1800	A Kit with the next 6 items	kit	\$257.79	\$257.79	\$232.01	\$232.01
G1112TI	Gracey 11/12	each	\$51.79	\$51.79	\$46.61	\$46.61
G1314TI	Gracey 13/14	each	\$51.79	\$51.79	\$46.61	\$46.61
COL1314TI	Columbia 13/14	each	\$51.79	\$51.79	\$46.61	\$46.61
204STI	204S	each	\$51.79	\$51.79	\$46.61	\$46.61
H6H7TI	H6-H7	each	\$51.79	\$51.79	\$46.61	\$46.61


Sterilization tray

9

HU-FRIEDY® TITANIUM IMPLANT SCALERS

Perform at your best with Hu-Friedy's state of the art line of Titanium Implant Scalers. Made from the same titanium alloy as implants and abutments, these scalers are expertly designed and manufactured for implant maintenance, debridement, biofilm removal and can be used both supra and subgingival.

- Improved visual acuity with Hu-Friedy's signature teal colored anodized titanium
- Optimized sharpness with cutting edges that are finely honed and sharpened post anodization process
- Enhanced instrument characteristics for implant maintenance, the ability to remove calculus and biofilm
- Increased instrument value with the ability to be sharpened at any time to extend its useful life
- Ergonomic sensitivity featuring Hu-Friedy's ergo handle that offers a large diameter and cross/hatch design for maximum tactile grip and support



IMP204SDT	B 204SD	each	\$71.79	\$71.79	\$64.61	\$64.61
IMPLG12T	C Langer 1/2	each	\$71.79	\$71.79	\$64.61	\$64.61
IMPM12T	D Mini five 1/2	each	\$71.79	\$71.79	\$64.61	\$64.61
IMPM1112T	E Mini five 11/12	each	\$71.79	\$71.79	\$64.61	\$64.61
IMPM1314T	F Mini five 13/14	each	\$71.79	\$71.79	\$64.61	\$64.61

TESS IMPLANT PROPHY+™ KIT WITH SHARPENING STONE

- Resharpenable and autoclavable, sharpening stone included
- Graphite-free, safe on all implant surfaces
- Can be used on root-sensitive areas, composites and veneers
- Used on implant surfaces, with no unsightly lines on veneers or composite surfaces like metal devices; plastic contains no glass or graphite fillers
- Plastic tool and sharpening stone are autoclavable

502	Kit includes: Gracey 5/6, Gracey 11/12, Gracey 13/14, Columbia 13/14U and a Sharpening Stone	kit	\$102.49	\$102.49	\$92.24	\$92.24
-----	---	-----	----------	----------	---------	---------



SCALERS & CURETTES

HU-FRIEDY® THE EVEREDGE® 2.0 COLLECTION

Hu-Friedy's EverEdge 2.0 Technology is unlike anything you've experienced before in a scaler. Hu-Friedy has combined state-of-the-art technology in metallurgy, heat treatment and cryogenics to create a superior stainless steel that keeps your scalers and curettes sharper longer. That means less frequent sharpening, less hand fatigue, and greater comfort throughout the day.

STAYS SHARPER LONGER

Revolutionary stainless steel alloy is super-durable and stays sharper longer.

PROPRIETARY PROCESS

Proprietary processing ensures superior edge retention will last the entire life of the instrument. EverEdge 2.0 is not a superficial coating, so instrument blades can be resharpened.

DIAMOND KNURL

The unique diamond knurl pattern and large diameter handle provides a confident and comfortable grasp. The #9 handle was also designed to be the optimal weight for reduced hand fatigue.

COMFORT ZONE CONE WITH SIGNATURE SERIES GRIP

The cushioned grip provides a smooth transition and increased comfort. With distinctive purple grips, EverEdge 2.0 scalers can be easily identified in a set up.

FINELY-HONED FINISH

EverEdge 2.0 Scalers feature a hand-crafted finish that keeps edges effectively sharp.



EverEdge® 2.0 #9 Handles



EverEdge® 2.0 Resin Handles
8 Colors = Oral Location

- Purple Resin - Anterior Sickle Scalers
- Red Resin - Posterior Sickle Scalers
- Lavender Resin - Universal Curettes
- Gray Resin - Gracey Anterior
- Yellow Resin - Gracey Anterior Scalers
- Green Resin - Gracey Buccal, Lingual Scalers
- Orange Resin - Gracey Mesial Surfaces, Premolars, Molars
- Blue Resin - Gracey Distal Surfaces, Premolars, Molars

Hu-Friedy® EverEdge® 2.0 GRACEY Curettes

SG129E2	A	#1/2 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG12C8E2	B	#1/2 Gracey Curette	Resin Handle <input checked="" type="checkbox"/>	\$55.29	\$55.29	\$49.76	\$49.76
SG349E2	C	#3/4 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG569E2	D	#5/6 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG789E2	E	#7/8 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG9109E2	F	#9/10 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG11129E2	G	#11/12 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG1112R9E2		#11/12R Gracey Rigid Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG11149E2	H	#11/14 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG12139E2	I	#12/13 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG13149E2	J	#13/14 Gracey Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SG1314R9E2	K	#13/14R Gracey Rigid Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26

Hu-Friedy® EverEdge® 2.0 SICKLE Scalers

SH5339E2	L	#H5/33 Scaler	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SH679E2	M	#H6/H7 Scaler	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SH67C8E2		#H6/H7 Scaler	Resin Handle <input checked="" type="checkbox"/>	\$55.29	\$55.29	\$49.76	\$49.76
SN1379E2	N	#137 Sickle/Curette	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
S204S9E2	O	#204S Sickle	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SCNEVI19E2	P	Nevi #1 Anterior Scaler	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
SCNEVI29E2	Q	Nevi #2 Posterior Scaler	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26

Hu-Friedy® EverEdge® 2.0 YOUNGER GOOD Scalers


SYG789E2	R	#7/8 Younger Good Scaler	#9 Handle <input type="checkbox"/>	\$60.29	\$60.29	\$54.26	\$54.26
----------	----------	--------------------------	------------------------------------	---------	---------	---------	---------



INSTRUMENTS

ACE STAINLESS STEEL INSTRUMENT TRAYS

- With 2 Supports and 1 Retainer
- With Silicone Storing Frame



9400035	A	10.75" x 6.75" x 1"	each	\$146.49	\$146.49	\$131.84	\$131.84
9400036	B	7.375" x 5.6" x 1"	each	\$135.79	\$135.79	\$122.21	\$122.21
9400037	C	7" x 3.6" x 1"	each	\$87.99	\$87.99	\$79.19	\$79.19
9400038	D	7.6" x 5.6" x 1.25"	each	\$227.79	\$227.79	\$205.01	\$205.01


ACE POLYMER (PLASTIC) INSTRUMENT TRAYS

- Autoclave up to 270° F



Tray with Finger Mat

7651001		6" x 2-1/2" x 3/4"	each	\$34.99	\$34.99	\$31.49	\$31.49
7651002	E	6-1/2" x 4", 3/4"	each	\$48.79	\$46.35	\$46.35	\$43.91
7651011		10" x 6" x 3/4"	each	\$96.79	\$91.95	\$91.95	\$87.11



2-Tier Sterilization Tray

7651021	F	10" x 6" x 1-1/2"	each	\$201.79	\$201.79	\$181.61	\$181.61
---------	----------	-------------------	------	----------	----------	----------	----------

9

HU-FRIEDY® SIGNATURE SERIES® ORAL SURGERY CASSETTES

- Fits in most 10" chambers or larger sterilizers, most instrument washers and larger ultrasonics.
- IMS® cassettes are made of long-lasting, low-maintenance stainless steel.
- Soft, colored silicone rails securely hold instruments throughout the processing cycle.
- Accommodates 12 instruments, 4 elevators, 3 hinged instruments in the lid and has a separate accessory area.
- Has a built-in needle cap holder and slot for a winged syringe.
- Exterior dimensions: 8" x 14.5" x 1.25".



IMOS0	red	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS1	gray	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS2	white	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS3	orange	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS5	yellow	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS6	lavender	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS8	blue	each	\$280.29	\$280.29	\$252.26	\$252.26
IMOS9	green	each	\$280.29	\$280.29	\$252.26	\$252.26

ACE STERILIZING TRAYS

- Perforated



0033123		10" x 6-1/2" x 2-1/2"	each	\$174.79	\$174.79	\$157.31	\$157.31
---------	--	-----------------------	------	----------	----------	----------	----------

HU-FRIEDY® IMS PARTS BOX

- Designed to keep small items such as rubber dam clamps and matrix band retainers organized during the sterilization process.
- Exterior dimensions: 3.5" x 1.12" x 1"



IMS1273	Large	each	\$57.79	\$57.79	\$52.01	\$52.01
---------	-------	------	---------	---------	---------	---------

INTEGRA™ STAINLESS MESH BUR AND DIAMOND HOLDER



6528		2-1/2" inside diameter	each	\$16.49	\$16.49	\$14.84	\$14.84
------	--	------------------------	------	---------	---------	---------	---------

ACE CLEANING BRUSHES

- Not autoclavable
- 7" (18cm)



6039015	Nylon Bristles	3/pack	\$15.49	\$15.49	\$13.94	\$13.94
6039016	Stainless Steel Bristles	3/pack	\$16.99	\$16.99	\$15.29	\$15.29

INSTRUMENTS

ACE STAINLESS STEEL WARE

- Heavy Gauge, 18/8, 304 stainless steel throughout
- Seamless, crevice-free, sanitary construction
- Uniform flat bottoms
- Contoured for nursing efficiency and patient comfort
- Permanent, brilliant mirror finish



ACE INSTRUMENT TRAYS

Round corners

0531030	A	8-7/8" x 5" x 2"	each	\$25.99	\$25.99	\$23.39	\$23.39
---------	----------	------------------	------	---------	---------	---------	---------

Square corners

0531122	B	12-1/4" x 7-5/8" x 2-1/8"	each	\$32.29	\$32.29	\$29.06	\$29.06
---------	----------	---------------------------	------	---------	---------	---------	---------

ACE INSTRUMENT TRAY WITH STRAP HANDLE COVER



Covered Trays

0531037		8-1/2" x 3" x 1-1/2"	each	\$32.29	\$32.29	\$29.06	\$29.06
0531038		8-7/8" x 5" x 2"	each	\$39.49	\$39.49	\$35.54	\$35.54

ACE INSTRUMENT TRAY WITH RECESSED HANDLE COVER



Covered Trays

0531031		12-1/4" x 7-5/8" x 2-1/8"	each	\$48.79	\$48.79	\$43.91	\$43.91
---------	--	---------------------------	------	---------	---------	---------	---------

MEDLINE INSTRUMENT TRAY WITH RECESSED HANDLE COVER



- Heavy-duty, seamless tray has covered corners for easy sterilization
- Tray and lid sold separately
- Tray dimensions: 8-3/4" x 3-1/4" x 1-1/2"

DYND0583Z	Tray only	each	\$30.99	\$30.99	\$27.89	\$27.89
DYND0583ZZ	Lid only	each	\$23.79	\$23.79	\$21.41	\$21.41

ACE INSTRUMENT JARS



Forceps Jar - Heavy Base

0531027		4-7/16" x 2"	each	\$35.99	\$35.99	\$32.39	\$32.39
0531025	C	7-1/4" x 2"	each	\$40.49	\$40.49	\$36.44	\$36.44

Needle/Ointment Jar

0531008	D	1/4 qt. 3" x 2-1/2"	each	\$17.49	\$17.49	\$15.74	\$15.74
---------	----------	---------------------	------	---------	---------	---------	---------

Thermometer Jar

0531026	E	4" x 1"	each	\$21.99	\$21.99	\$19.79	\$19.79
---------	----------	---------	------	---------	---------	---------	---------

Dressing Jar - Recessed Handle

0531013	F	1/2 qt. 4-1/8" x 3-1/4"	each	\$29.49	\$29.49	\$26.54	\$26.54
0531014		1 qt. 4-1/8" x 5-1/4"	each	\$35.99	\$35.99	\$32.39	\$32.39
0531015		2-1/4 qt. 4-7/8" x 6-7/8"	each	\$48.79	\$48.79	\$43.91	\$43.91
0531016		4-1/2 qt. 7" x 7-3/8"	each	\$56.29	\$56.29	\$50.66	\$50.66

SPONGE BOWLS



ACE

0531009		1 qt. 5-3/4" x 2-1/2"	each	\$29.49	\$29.49	\$26.54	\$26.54
0531010	G	1-3/8" qt. 6-1/2" x 3"	each	\$32.29	\$32.29	\$29.06	\$29.06

Integra™

3903	H	30oz. 6-27/64" x 2-1/4"	each	\$18.79	\$18.79	\$16.91	\$16.91
------	----------	-------------------------	------	---------	---------	---------	---------

ACE EMESIS BASIN



0531033		6" x 3-1/2" x 1-1/2"	each	\$17.49	\$17.49	\$15.74	\$15.74
0531034		8" x 3-3/4" x 1-5/8"	each	\$17.49	\$17.49	\$15.74	\$15.74
0531035		10" x 4-5/8" x 2-5/8"	each	\$17.49	\$17.49	\$15.74	\$15.74
0531036		12" x 5-1/2" x 2-1/4"	each	\$25.99	\$25.99	\$23.39	\$23.39

ACE INSTRUMENT CUPS



Medicine Cup, Graduated

0531005	I	2 oz., 4 tsp., ACE	each	\$7.79	\$7.79	\$7.01	\$7.01
---------	----------	--------------------	------	--------	--------	--------	--------

Iodine or Oil Cup

0531007	J	6 oz., ACE	each	\$8.29	\$8.29	\$7.46	\$7.46
3912	J	6 oz., Integra™	each	\$8.49	\$8.49	\$7.64	\$7.64

Tumbler

0531006	K	7 oz., ACE	each	\$9.49	\$9.49	\$8.54	\$8.54
---------	----------	------------	------	--------	--------	--------	--------

INSTRUMENTS

GRAHAM FIELD COVERED DRESSING JARS

- Made of high-resistance glass
- Unlabeled
- Stainless steel covers
- Can only be cleaned with alcohol free soap and water
- ⚠ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov



GF3460	3" x 3", holds 2"x 2" gauze	each	\$15.99	\$15.19	\$15.19	\$14.39
3461	4" x 4", holds 3"x 3" gauze	each	\$20.79	\$19.75	\$19.75	\$18.71
3462	5" x 5", holds 4"x 4" gauze	each	\$23.49	\$22.32	\$22.32	\$21.14
3463	6" x 6", holds 5"x 5" gauze	each	\$25.99	\$24.69	\$24.69	\$23.39

GRAHAM FIELD LABELED GLASS SUNDRY JARS

- Each jar is 7" high and 4" diameter
- Stainless steel covers
- Can only be cleaned with alcohol free soap and water
- ⚠ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov



3454AC	Cotton jar	each	\$16.49	\$15.67	\$15.67	\$14.84
3454A	Applicator jar	each	\$16.49	\$15.67	\$15.67	\$14.84
3454AG	Gauze jar	each	\$16.49	\$15.67	\$15.67	\$14.84
3454AT	Tongue Depressor jar	each	\$16.49	\$15.67	\$15.67	\$14.84

CLEAR GLASS DAPPEN DISH

Keystone Ind.

- Made of solid heat-resistant glass.
- The thick and heavy base makes dish stable, protects from spills.
- Ideal for mixing acrylics, amalgam, resins, sealants and other materials.



0600610	Clear	each	\$2.79	\$2.79	\$2.51	\$2.51
---------	-------	------	--------	--------	--------	--------

Buffalo Dental

- Cold Sterilization Only.
- 1-1/4" OD, 1 x 1-3/8" Depth.



46520	Clear	each	\$6.79	\$6.79	\$6.11	\$6.11
-------	-------	------	--------	--------	--------	--------

CEMENT SPATULA



Integra™

7376	C-5	each	\$34.29	\$34.29	\$30.86	\$30.86
------	-----	------	---------	---------	---------	---------

J&J Instruments

13124	#24	each	\$6.99	\$6.99	\$6.29	\$6.29
13125	#24A	each	\$6.99	\$6.99	\$6.29	\$6.29

ACE WOODSON PLASTIC FILLING INSTRUMENT

- Double Ended
- Stainless steel

6200097A	#1, 6-1/2"	each	\$28.49	\$28.49	\$25.64	\$25.64
----------	------------	------	---------	---------	---------	---------

KEYSTONE IND. DISPOSABLE DAPPEN DISHES

- Disposable dappen dishes are made of special thermoplastic material which is impervious to all dental materials.
- The dishes are designed for quick and easy mixing of materials such as prophylactic paste, fluoride products, composite materials, acrylics, cements and amalgams.
- Assorted: White, Yellow and Blue



1014410		1,000/box	\$36.49	\$34.67	\$34.67	\$32.84
---------	--	-----------	---------	---------	---------	---------

PULPDENT MINI BOWLS

- Clear, Autoclavable Silicone
- Non-stick coating for easy mixing
- Suction cup bottom for convenient, one-hand mixing



BMS3	Small, 8cc	set of 3	\$21.29	\$21.29	\$19.16	\$19.16
BMM2	Medium, 30cc	set of 2	\$21.29	\$21.29	\$19.16	\$19.16



INTEGRA™ TIP-IT™ INSTRUMENT GUARDS

- Used to protect instrument tips and edges from damage during instrument handling and sterilization
- Steam and gas sterilizable
- Radiopaque and latex-free; assorted colors

32510	50/pack	\$16.99	\$16.99	\$15.29	\$15.29
-------	---------	---------	---------	---------	---------

SHAMROCK INSTRUMENT MARKING TAPE



- 1/4" x 200" rolls

SMT3826	black	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT386	blue	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT388	brown	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT3814	gray	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT383	light green	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT385	orange	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT3813	purple	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT384	red	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT381	white	each	\$16.99	\$16.99	\$15.29	\$15.29
SMT382	yellow	each	\$16.99	\$16.99	\$15.29	\$15.29

HU-FRIEDY® COLOR CODE INSTRUMENT RINGS

- Silicone rubber rings for identifying and organizing instruments
- Will withstand all heat sterilization methods



Regular - For Use on #2 and #4 Handles

IMS1283	orange	50/pack	\$17.49	\$17.49	\$15.74	\$15.74
IMS1289	black	50/pack	\$17.49	\$17.49	\$15.74	\$15.74
IMS1280	assorted colors	90/pack	\$29.49	\$29.49	\$26.54	\$26.54

Regular - For Use on #2 and #4 Handles

IMS1286L	red	50/pack	\$17.49	\$17.49	\$15.74	\$15.74
IMS1288L	blue	50/pack	\$17.49	\$17.49	\$15.74	\$15.74
IMS12810L	pink	50/pack	\$17.49	\$17.49	\$15.74	\$15.74
IMS12811L	lavender	50/pack	\$17.49	\$17.49	\$15.74	\$15.74

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

Introducing the ACE SOUTHERN Alliance Program (ASAP)

When you create an alliance with **ACE SOUTHERN** you have a committed partner in your success. Part of that commitment is to reward you **for your loyalty with significant savings on every dollar you spend.**

As an ASAP Qualified Customer, you will benefit from savings and special promotions across all the products **ACE SOUTHERN** offers to help your practice prosper.

	Premier \$15K/yr spend	Premier Plus* \$25K/year spend
All Rx, Supplies and Equipment	5% saving	10% saving
All Instruments	10% savings	10% savings
Bone Grafting Materials and Implants	15-20% savings	15-20% savings

* AAOMS/AAP members only

Enroll for free today at www.acesouthern.com/asap



INSTRUMENTS



9

INSTRUMENT STANDS

Blickman Mayo Stands **A**

- Maximum utility for the crowded surgical environment
- Foot-operated with 40"-62" height adjustment
- Stainless steel construction
- Low profile permits base to easily slide under surgical table base or beds with low clearance
- Heavy, weighted base
- Cantilever design provides added strength and stability
- Feature two 2" swivel casters
- Drop shipped item, non returnable
- Please Note: Drop ship and shipping fees not included in the pricing.

8867SS	A	12-5/8" x 19-1/8"	each	\$1,021.29	\$1,021.29	\$919.16	\$919.16
8869SS	A	16" x 21"	each	\$1,125.79	\$1,125.79	\$1,013.21	\$1,013.21

Blickman Mayo Instrument Stand

- Chrome steel frames
- Designed for use under low clearance units
- Static free due to contact of metal base with floor
- Height adjustment: 27-1/2" to 46-1/2"
- Offset position extends 20" over table
- Tray size: 12-5/8" x 19-1/8"
- Spread of legs: 18-1/2" x 24"

1515	B		each	\$179.49	\$179.49	\$161.54	\$161.54
------	----------	--	------	----------	----------	----------	----------

Blickman Instrument Stand

- Chrome steel frames
- 5 leg base
- Wall-saver polished aluminum base
- Height adjustment: 27-1/2"-46-1/2"
- Tray size: 12-1/2" x 19" stainless steel
- Offset position extends 15.5" over table
- Spread of legs: 23" diameter

BL11501	C		each	\$195.79	\$195.79	\$176.21	\$176.21
---------	----------	--	------	----------	----------	----------	----------

Brandt Industries Mayo Instrument Stands

- Brandt instrument stands combine stability with mobile convenience
- Chrome steel frames
- 12-1/2" x 19" stainless steel trays
- 2" dual wheel casters

Mayo Instrument Stand

- Height adjustment 33" - 52"

33100	D		each	\$142.49	\$142.49	\$128.24	\$128.24
-------	----------	--	------	----------	----------	----------	----------

Dual Position Instrument Stand with Aluminum Base

- May be used in either offset or center post position
- Supplied with wall-saver polished aluminum base

39101	E		each	\$174.49	\$174.49	\$157.04	\$157.04
-------	----------	--	------	----------	----------	----------	----------

Double Post Instrument Stand

- Adjustable from 34" to 54-1/2"
- For use under low clearance units

BRN33101	F		each	\$151.29	\$151.29	\$136.16	\$136.16
----------	----------	--	------	----------	----------	----------	----------

Replacement Stainless Steel Tray

- For: 33100, 39101, and BRN33101
- Tray size: 12-1/2" x 19" x 3/4" stainless steel

027		Tray Only	each	\$90.29	\$90.29	\$81.26	\$81.26
-----	--	-----------	------	---------	---------	---------	---------

INSTRUMENTS

MAYO TRAYS



ACE Mayo Trays

- Lustrous, heavy-duty instrument trays
- Stainless Steel


0531121	15-1/8" x 10-1/2" x 5/8"	each	\$38.29	\$38.29	\$34.46	\$34.46
0531024	19" x 12" x 5/8"	each	\$46.79	\$46.79	\$42.11	\$42.11

Blickman Mayo Stand Trays

- Heavy-duty

9719126250	12-5/8" x 19-1/8"	each	\$69.49	\$69.49	\$62.54	\$62.54
IT1621	16" x 21"	each	\$154.29	\$154.29	\$138.86	\$138.86

COVIDIEN DEVON®
HANDS-FREE-TRANSFER™
MAGNETIC DRAPES



- Disposable magnetic drapes are specifically designed to attract and hold surgical instruments in the sterile field
- Flexible drape conforms to any shape with non-slip foam backing for secure positioning
- Raised magnets provide strong holding power while allowing instruments to be easily grasped

31140570	10" x 16"	30/case	\$252.79	\$240.15	\$240.15	\$227.51
----------	-----------	---------	----------	----------	----------	----------

MOBILE I.V. INFUSION STANDS

Steelcraft 4-Hook Mobile IV Stand

- Five point all steel base for greater stability
- 2" ball bearing swivel casters
- Full height extension from 60" to 110"
- Shipped completely knocked down in two cartons with simple assembly

399	A	4-hook stand	each	\$178.49	\$169.57	\$169.57	\$160.64
-----	---	--------------	------	----------	----------	----------	----------

ACE I.V. SOLUTION
CEILING SUSPENSION UNIT

- Attaches to track of suspended ceiling
- Can be attached to ceiling timber, using plate and screws provided
- Includes 2' (varying length) chain

0031183		each	\$62.99	\$59.84	\$59.84	\$56.69
---------	--	------	---------	---------	---------	---------

FOOTSTOOLS

- Stainless steel construction
- Flared legs with non-slip conductive vinyl top for firm support without tipping
- Static-electricity conductive rubber feet

Mabis Healthcare Footstool - 13.5" x 18" Base

- 9.5"H x 14"W x 10"D

53919010000	C	each	\$53.49	\$50.82	\$50.82	\$48.14
-------------	---	------	---------	---------	---------	---------

Blickman O.R. Quality, One Step

- 8"H - x 18"W x 12"D

7757MR	D	each	\$285.29	\$271.03	\$271.03	\$256.76
--------	---	------	----------	----------	----------	----------

BREWER ADD-A-TRAY

- Holds instruments or patient items
- 15 lb maximum
- Removable stainless steel tray
- Tray size: 9.75" x 13.5" x .75"
- Clamp sold separately below



43431	Tray only	each	\$138.79	\$131.85	\$131.85	\$124.91
43433	Clamp only	each	\$98.79	\$93.85	\$93.85	\$88.91



INSTRUMENTS



ACE GALETTI ARTICULATOR

- Repairs Available

0500096	each	\$590.79	\$590.79	\$531.71	\$531.71
---------	------	----------	----------	----------	----------



Short-Handled Drivers

Long-Handled Drivers

A-TITAN UNIVERSAL DRIVERS

- Long and short-handled drivers work with all major implant systems (compatible with 95% of all implant systems sold in the United States)
- Long-handled drivers provide outstanding access between teeth and are specifically designed for maxillary cases
- Short drivers are designed for anterior mandible cases
- Drivers are color-coded for easy identification
- Fully autoclavable

Implant Driver Kits Include:


- Straight Screwdriver (Straight Systems Tips)
 - SCS Screwdriver (Straumann Tips)
 - Unigrip Driver (Replace Select Tips)
 - Hex 1.3mm/.050" (Lifecore, Biohorizons, Astra, Zimmer, Camlog Tips)
 - Hex 1.2mm/.048" (3i System Tips)
 - Hex 0.9mm/.035" (Friadent, Ankylos, 3i Cover Screw Tips)
 - Square 1.3mm/.050" (Gold Screw/Nobel Biocare Tips)
- Kit contents also sold individually.

4.25" Long-Handled Drivers Kit

K539	set of 7 drivers/kit	\$314.29	\$314.29	\$282.86	\$282.86
------	----------------------	----------	----------	----------	----------

1.7" Short-Handled Drivers Kit

K540	set of 7 drivers/kit	\$314.29	\$314.29	\$282.86	\$282.86
------	----------------------	----------	----------	----------	----------



In Port

Out Port

Trocar for Superior Chamber

ACE TMJ DYSFUNCTION TREATMENT

- Levaging TM Joint – Shepherd Arthrocentesis Instrument
- A small incision (1/8") is made into the skin; mosquito forceps are used to make a blunt dissection to the lateral ligament; the arthrocentesis instrument may then enter the capsule into the superior chamber.

Reference: Trieger, N, Hoffman, C, Rodriguez, E. The Effect of Arthrocentesis of the Temporomandibular Joint in Patients With Rheumatoid Arthritis. Journal of Oral and Maxillofacial Surgery May 1999; 57:537-541.

Shepherd Arthrocentesis Instrument


0702042	18G x 2.5"	each	\$167.99	\$167.99	\$151.19	\$151.19
---------	------------	------	----------	----------	----------	----------

9

"CAN YOU ORDER SOME MORE OF THESE MASKS?"

ORDER HISTORY, JUST ANOTHER FEATURE ON YOUR DASHBOARD

VISIT ACESOUTHERN.COM TODAY



REJUVENATE YOUR INSTRUMENTS

One Visit, and They'll Come Back Sharper and Stronger.




IN-HOUSE REPAIR DEPARTMENT.

ACE SOUTHERN offers a state-of-the-art repair department, dedicated to reconditioning your surgical instruments. In addition, our repair department will sharpen or retip most of your cutting and drilling instruments. For over 50 years our repair center has maintained the latest in calibration tools, multiple machining centers, torque and tensile testing units, and digital inspection equipment and gauges to keep your tools in shape.

- **Air Drills**
- **Bur Guards**
- **Electric Surgical Motors**
- **Electrosurgery Units**
- **Sterilizers**

IV ADMINISTRATION




EXEL IV ADMINISTRATION SETS

- Sterile
- Latex-free
- Non-toxic
- Non-pyrogenic fluid path
- Vented/non-vented
- Roller clamp
- Occlusion clamp
- Pouch packaging
- Priming volume: 23ml/approx
- Injection site

A 78" IV Set - Luer Slip					
29081	15 drops/ml	50/case	\$54.49	\$51.77	\$49.04
29083	60 drops/ml	50/case	\$61.99	\$58.89	\$55.79
B 78" IV Set - Luer Lock/Slip (Slides to Injection Site)			COMPARE TO V1416		
29086	15 drops/ml	50/case	\$69.29	\$65.83	\$62.36
C 78" IV Set - Luer Lock					
29082	15 drops/ml	50/case	\$62.29	\$59.18	\$56.06
29085	60 drops/ml	50/case	\$63.49	\$60.32	\$57.14
D 105" IV Set - Additional Injection Site And Occlusion Clamp, Priming Volume: 28ml/Approx					
29088	15 drops/ml	50/case	\$66.29	\$62.98	\$59.66


10



B|BRAUN IV ADMINISTRATION SETS

- Universal spike
- Roller clamp
- SPIN-LOCK® connector
- Latex-free


E 73" IV Set - Injection Site 6" Above Distal End, Slide Clamp, Priming Volume: 16ml					
V1416	15 drops/ml	50/case	\$151.29	\$143.73	\$136.16
F 79" IV Set - Slide Clamp, Priming Volume: 14ml					
V1402	15 drops/ml	50/case	\$62.99	\$59.84	\$56.69
G 81" IV Set - Injection Site 6" Above Distal End, Priming Volume: 16ml					
V1390	15 drops/ml	50/case	\$113.29	\$107.63	\$101.96
H 93" IV Set - Injection Site 26" From Distal End, Slide Clamp, Priming Volume: 16ml					
V141515	15 drops/ml	50/case	\$111.49	\$105.92	\$100.34
I 106" IV Set - Injection Sites 82" and 6" Above Distal End, Slide Clamp, Priming Volume: 19ml					
V1423	15 drops/ml	50/case	\$162.79	\$154.65	\$146.51
J 106" IV Set - Injection Sites 28" and 6" Above Distal End, Slide Clamp, Priming Volume: 19ml					
V1426	60 drops/ml	50/case	\$113.29	\$107.63	\$101.96



B|BRAUN IV ADMINISTRATION SET

66" IV Set - Injection Sites 28" Above Distal End, Nonvented Spike, Control Clamp, Slide Clamp, Priming Volume: 12ml

V142701	60 drops/ml	50/case	\$126.29	\$119.98	\$113.66
---------	-------------	---------	----------	----------	----------



HOSPIRA MACROBORE IV ADMINISTRATION SETS

- Convertible piercing pin and sight chamber
- Non-DEHP
- CAIR™ clamp and secure lock male adapter

34" Secondary Set - With Extension Hook

1423028	15 drops/ml	50/case	\$98.79	\$93.85	\$88.91
---------	-------------	---------	---------	---------	---------

100" Primary Set - With Slide Clamp

1422128	15 drops/ml	50/case	\$139.29	\$132.33	\$125.36
---------	-------------	---------	----------	----------	----------

IV ADMINISTRATION

BAXTER IV ADMINISTRATION SETS

Injection site 6" from male luer lock adapter w/ retractable collar

K 103" Solution Set - With DUO-VENT Spike

ACT5435	10 drops/ml	48/case	\$141.99	\$134.89	\$134.89	\$127.79
---------	-------------	---------	----------	----------	----------	----------

L 104" Set - With DUO-VENT Spike, Additional Injection Site 34" From Male Luer Lock Adapter

ACT5439	10 drops/ml	48/case	\$145.29	\$138.03	\$138.03	\$130.76
---------	-------------	---------	----------	----------	----------	----------

M 67" Mini-Drip Set

ACT5432	60 drops/ml	48/case	\$157.49	\$149.62	\$149.62	\$141.74
---------	-------------	---------	----------	----------	----------	----------



AMSINO IV ADMINISTRATION SET

- Universal spike w/ drip chamber
- Slide clamp and roller clamp
- Latex-free injection site
- Rotating male luer lock, vented/non-vented

78" Standard Set

AA3101	15 drops/ml	50/case	\$71.79	\$68.20	\$68.20	\$64.61
--------	-------------	---------	---------	---------	---------	---------



KAWASUMI IV ADMINISTRATION SETS

- Non-DEHP
- Latex-free

N 84" IV Set - Injection Site, Priming Volume: 18ml

IV009Y	10 drops/ml	50/case	\$56.49	\$53.67	\$53.67	\$50.84
--------	-------------	---------	---------	---------	---------	---------

O 100" IV Set - 2 Injection Sites, Priming Volume: 23ml

IV016Y	10 drops/ml	50/case	\$72.99	\$69.34	\$69.34	\$65.69
--------	-------------	---------	---------	---------	---------	---------



UNIVERSAL SYRINGE PUMP SETS

Exel 72"

SA1772		60/case	\$179.99	\$170.99	\$170.99	\$161.99
--------	--	---------	----------	----------	----------	----------

Baxter 72" - Male Luer Lock Adapter, 3.0ml Approximate Fill Volume

2C9218		60/case	\$289.49	\$275.02	\$275.02	\$260.54
--------	--	---------	----------	----------	----------	----------



ACE EXTENSION SETS



Fabricated from high-quality medical tubing, odorless, non-toxic, luer slip. Individually packaged in a sterile see-through pouch.

No Medication Site

0030012	12"	2 cc Priming	100/case	\$78.99	\$75.04	\$75.04	\$71.09
0030020	20"	3 cc Priming	100/case	\$89.29	\$84.83	\$84.83	\$80.36
0030030	30"	4 cc Priming	100/case	\$95.99	\$91.19	\$91.19	\$86.39

IV ADMINISTRATION



EXEL EXTENSION SETS

- Luer lock/slip
- Latex-free

A No Injection Site

SA1712	12", 1.6ml priming volume	100/case	\$55.79	\$53.00	\$53.00	\$50.21
SA1720	20", 2.5ml priming volume	100/case	\$94.49	\$89.77	\$89.77	\$85.04
SA1730	30", 3.6ml priming volume	100/case	\$97.79	\$92.90	\$92.90	\$88.01

B With Injection Site

1720WY	20"	100/case	\$133.49	\$126.82	\$126.82	\$120.14
1730WY	30"	100/case	\$143.79	\$136.60	\$136.60	\$129.41



HOSPIRA MACROBORE EXTENSION SETS

- No injection site
- Female adapter
- Slide clamp
- Secure lock male adapter
- Non-DEHP

20", Priming Capacity: 2.7ml

1265528	50/case	\$116.79	\$110.95	\$110.95	\$105.11
---------	---------	----------	----------	----------	----------

30", Priming Capacity: 4ml

1265628	50/case	\$108.29	\$102.88	\$102.88	\$97.46
---------	---------	----------	----------	----------	---------

10



BAXTER EXTENSION SET

- Ideal for the administration of Propofol
- With luer-lock adapter

61", Priming Capacity: 1.4ml

2C9219	60/case	\$209.29	\$198.83	\$198.83	\$188.36
--------	---------	----------	----------	----------	----------

B|BRAUN EXTENSION SETS



C 60", Microbore - Female and Female Luer Lock Connectors, Priming Volume: 0.5ml

V6222	50/case	\$117.29	\$111.43	\$111.43	\$105.56
-------	---------	----------	----------	----------	----------



D 13", Small Bore - Female Luer Connector, Non-Removable Slide Clamp, Distal SPIN-LOCK®, Priming Volume: 0.4ml

473105	100/case	\$253.99	\$241.29	\$241.29	\$228.59
--------	----------	----------	----------	----------	----------



E 21", Standard Bore - Female Luer Connection, Distal SPIN-LOCK® Connector, Priming Volume: 3ml

V5406	50/case	\$95.49	\$90.72	\$90.72	\$85.94
-------	---------	---------	---------	---------	---------



F 34", Standard Bore - Female Luer Connector, Injection Sites 24" and 6" Above Distal End, Slide Clamp, Distal SPIN-LOCK® Connector, DEHP-Free, Latex-Free, Priming Volume: 5ml


V5409	50/case	\$124.79	\$118.55	\$118.55	\$112.31
-------	---------	----------	----------	----------	----------

NEEDLE-FREE SETS AND ACCESSORIES

BAXTER INTERLINK® PREPIERCED IV ADMINISTRATION SETS

- Male luer lock adapter
- Latex-free

G 76" Needle-Free Solution Set - 1 INTERLINK Injection Site, Lever Lock Cannula (Detached)						
2C6401	10 drops/ml	48/case	\$238.29	\$226.38	\$226.38	\$214.46
H 92" Needle-Free Solution Set - 1 INTERLINK Injection Site, DUO-VENT™ Spike						
2C6419	10 drops/ml	48/case	\$233.29	\$221.63	\$221.63	\$209.96
I 100" Needle-Free Solution Set - 2 INTERLINK Injection Sites						
2C6425	10 drops/ml	48/case	\$281.79	\$267.70	\$267.70	\$253.61
J 76" Minidrip Needle-Free Solution Set - 1 INTERLINK Injection Site, Lever Lock Cannula (Detached)						
2C6402	60 drops/ml	48/case	\$234.29	\$222.58	\$222.58	\$210.86
K 92" Minidrip Needle-Free Solution Set - 2 INTERLINK Injection Sites						
2C6424	60 drops/ml	48/case	\$259.49	\$246.52	\$246.52	\$233.54



BAXTER INTERLINK® INJECTION SITE


Luer Lock						
2N3399		200/case	\$457.99	\$435.09	\$435.09	\$412.19



BD™ INTERLINK® LEVER LOCK CANNULA

Versatile locking cannula for IV connection to injection site or Y-site - The INTERLINK Lever Lock Cannula is specifically designed for use with the INTERLINK injection site, distinctively marked with a colored alert ring. Also compatible with Hospira LifeShield® and B Braun Safeline®.


303370		100/box	\$52.29	\$49.68	\$49.68	\$47.06
--------	--	---------	---------	---------	---------	---------



HOSPIRA MACROBORE PREPIERCED IV ADMINISTRATION SETS

- Non-DEHP

L 78" LifeShield® Prepierced Primary IV Set						
Convertible piercing pin and sight chamber, slide clamp, CAIR™ clamp, prepierced y-site and secure lock male adapter						
1267228	15 drops/ml	48/case	\$121.49	\$115.42	\$115.42	\$109.34
M 80" LifeShield® Prepierced Primary Piggyback IV Set						
Convertible piercing pin and sight chamber, upper backcheck valve, 2 prepierced y-sites, 2 slide clamps, CAIR™ clamp, and secure lock male adapter						
1266528	15 drops/ml	48/case	\$186.29	\$176.98	\$176.98	\$167.66
N 70" Microdrip® Prepierced Primary IV Set						
Convertible piercing pin with sight chamber, slide clamp, CAIR™ clamp, prepierced y-site and secure lock male adapter						
1422228	60 drops/ml	48/case	\$156.79	\$148.95	\$148.95	\$141.11
O 80" LifeShield® Microdrip® Prepierced Primary Piggyback IV Set						
Convertible piercing pin and sight chamber, automatic backcheck valve with prepierced y-site, slide clamp, CAIR™ clamp, prepierced y-site and Option-Lok™ male adapter						
1266928	60 drops/ml	48/case	\$175.29	\$166.53	\$166.53	\$157.76



IV ADMINISTRATION

NEEDLE-FREE SETS AND ACCESSORIES



HOSPIRA LIFESHIELD® PREPIERCED RESEAL MALE ADAPTER PLUG

- For use with Lifeshield blunt cannula, prefilled and emergency syringes
- Attaches to hub of venipuncture device allowing for intermittent IV administration
- Priming volume - 0.1ml
- Latex-free

1130101	120/case	\$371.99	\$353.39	\$353.39	\$334.79
---------	----------	-----------------	-----------------	-----------------	-----------------



HOSPIRA MACROBORE PREPIERCED EXTENSION SET

32" Female Adapter, 2 Slide Clamps, 2 Prepierced Y-Sites and Option-Lok™ Male Adapter

1264828	48/case	\$108.29	\$102.88	\$102.88	\$97.46
---------	---------	-----------------	-----------------	-----------------	----------------

10



AMSINO PREPIERCED/NEEDLE-FREE IV ADMINISTRATION SET

83" Non-Vented, Roller Clamp, Pre-Pierced Y-Site, Sure-Lok® Needle-Free Y-Site, Rotating Male Luer Lock, DEHP-Free, Latex-Free

108306	10 drops/ml	50/case	\$107.99	\$102.59	\$102.59	\$97.19
--------	-------------	---------	-----------------	-----------------	-----------------	----------------



AMSINO PREPIERCED EXTENSION SET

27" Female Luer Lock Adapter, Rotating Male Luer Lock, 2 Injection Sites, 2 Slide Clamps, DEHP-Free, Latex-Free

003201		50/case	\$62.79	\$59.65	\$59.65	\$56.51
--------	--	---------	----------------	----------------	----------------	----------------



EXEL ALUMINUM HUB BLUNT NEEDLES

- Sterile, single-use, individually packaged blunt needle
- Stainless steel laser welded, non-lubricated cannula
- Aluminum hub, luer lock design
- Hub crimp cannula attachment
- Color-coded for easy identification
- 2-piece polypropylene needle cap and sheath

26471	15g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26472	17g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26473	18g x 1"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26474	19g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26475	20g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89

NEEDLE-FREE SETS AND ACCESSORIES

COVIDIEN MONOJECT™ BLUNT CANNULAS

- Stainless steel cannula
- Aluminum luer lock hub
- Smooth rounded blunt tip
- Color-coded needle sheath
- Autoclavable in rigid pack
- Sterile

202314	15g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202322	16g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202330	17g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202348	18g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202355	19g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202363	20g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202371	21g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202389	22g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202397	23g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64



HOSPIRA CLAVE® NEEDLE-FREE MACROBORE IV ADMINISTRATION SETS

A 100" LifeShield® Needle-Free Primary IV Set

Convertible piercing pin with sight chamber, slide clamp. CAIR™ clamp, Clave™ Y-site and Option-Lok™ male adapter, non-DEHP

1266128	15 drops/ml	48/case	\$271.29	\$257.73	\$257.73	\$244.16
---------	-------------	---------	----------	----------	----------	----------

B 100" LifeShield® Needle-Free Primary IV Set


Convertible piercing pin with sight chamber, backcheck valve, 2- Clave™ Y-sites, slide clamp, CAIR™ clamp and Option-Lok™ male adapter, non-DEHP

1266228	15 drops/ml	48/case	\$294.49	\$279.77	\$279.77	\$265.04
---------	-------------	---------	----------	----------	----------	----------

C 100" LifeShield® Needle-Free Microdrip® Primary IV Set

Convertible piercing pin with Microdrip™ sight chamber, slide clamp, CAIR™ clamp and Option-Lok™ male adapter, non-DEHP

2067028	60 drops/ml	48/case	\$315.29	\$299.53	\$299.53	\$283.76
---------	-------------	---------	----------	----------	----------	----------



10

HOSPIRA CLAVE® NEEDLE-FREE CONNECTOR

- Luer slip, lock and rotating luer compatible
- The unique dedicated internal fluid path of the CLAVE prohibits contact of the IV solutions with outside components
- A mechanically closed device which is automatically self-sealing when not accessed to prevent blood or fluid loss
- Requires no additional components

1195601		100/case	\$278.99	\$265.04	\$265.04	\$251.09
---------	--	----------	----------	----------	----------	----------



EXEL NEEDLE-FREE IV ADMINISTRATION SET

- Requires no additional cannulas or adapters for access

78" Needle-Free IV Set

Luer lock/slip, (slides to injection site) vented/non-vented, with needle-free injection site, roller clamp, occlusion clamp, latex-free

27086	15 drops/ml	50/case	\$119.29	\$113.33	\$113.33	\$107.36
-------	-------------	---------	----------	----------	----------	----------



IV ADMINISTRATION

NEEDLE-FREE SETS AND ACCESSORIES


ICU MEDICAL MICROCLAVE® NEEDLE-FREE IV ADMINISTRATION SETS



A 84" Needle-Free IV Set

MicroClave® connector, port location 8" from patient end, rotating luer, non-DEHP tubing

B990055	15 drops/ml	25/case	\$103.79	\$98.60	\$98.60	\$93.41
---------	-------------	---------	-----------------	----------------	----------------	----------------



B 106" Needle-Free IV Set

2 MicroClave® connectors, port locations at 12" and 86" from patient end, rotating luer, backcheck valve, non-DEHP tubing

B9900140	15 drops/ml	25/case	\$137.99	\$131.09	\$131.09	\$124.19
----------	-------------	---------	-----------------	-----------------	-----------------	-----------------



C 100" Needle-Free IV Set

MicroClave® connector, port location 6" from patient end, rotating luer, non-DEHP tubing

B99325	60 drops/ml	50/case	\$219.99	\$208.99	\$208.99	\$197.99
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------



D 106" Needle-Free IV Set

2 MicroClave® connectors, port locations 12" and 74" from patient end, rotating luer, non-DEHP tubing

B9900237	60 drops/ml	50/case	\$239.99	\$227.99	\$227.99	\$215.99
----------	-------------	---------	-----------------	-----------------	-----------------	-----------------

10




ICU MEDICAL MICROCLAVE® NEEDLE-FREE CONNECTOR

- Maintains a closed system whether or not the device has been activated


B99278		100/case	\$143.29	\$136.13	\$136.13	\$128.96
--------	--	----------	-----------------	-----------------	-----------------	-----------------

BAXTER CLEARLINK™ NEEDLE-FREE IV ADMINISTRATION SETS



E 76" Needle-Free IV Set, Priming Volume: 12.5ml

2C8401	10 drops/ml	48/case	\$256.79	\$243.95	\$243.95	\$231.11
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------



F 92" Needle-Free IV Set, Priming Volume: 14.42ml

Duo-Vent Spike

2C8419	10 drops/ml	48/case	\$311.99	\$296.39	\$296.39	\$280.79
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------

B|BRAUN SAFEDAY™ NEEDLE-FREE IV ADMINISTRATION SETS




106" Needle-Free IV Set

Universal spike, slide clamp, backcheck valve, SafeDAY valves 87" and 14" above distal end, roller clamp, sliding SPIN-LOCK connector, non-DEHP, latex-free, priming volume: 18.5ml

352636	15 drops/ml	50/case	\$172.29	\$163.68	\$163.68	\$155.06
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------

B|BRAUN CARESITE® VALVE NEEDLE-FREE IV ADMINISTRATION SETS



G 104" Needle-Free IV Set

Universal spike

354203	15 drops/ml	50/case	\$221.49	\$210.42	\$210.42	\$199.34
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------



H 103" Needle-Free IV Set

354204	60 drops/ml	50/case	\$289.29	\$274.83	\$274.83	\$260.36
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------

- Universal spike
- Roller clamp
- SPIN-LOCK® connector
- Slide clamp
- Backcheck valve
- CARESITE valve injection sites 12" and 87" above distal end
- Latex-free, non-DEHP

NEEDLE-FREE SETS AND ACCESSORIES

B|BRAUN ULTRASITE® NEEDLE-FREE INJECTION SITE

The ULTRASITE valve is the first high-flow, one piece capless valve with the positive displacement design. The patented ULTRASITE valve is a swabbable, capless, 2-way valve for injection/gravity/pump flow or aspiration of fluids. The ULTRASITE valve is a completely closed system. Upon luer removal, the ULTRASITE valve design expels fluid from the system, which helps prevent blood from backing up into the IV system. Latex-free.

- Truly passive needle-free compliance; does not permit needle use
- High flow rates; Flow rates up to 250ml/min



415110		100/case	\$243.99	\$231.79	\$231.79	\$219.59
--------	--	----------	-----------------	-----------------	-----------------	-----------------

KAWASUMI NEEDLE-FREE IV ADMINISTRATION SET

- Non-vented, one injection site
- Latex-free, non-DEHP

84" Needle-Free Standard Set

IV009N	10 drops/ml	50/case	\$112.99	\$107.34	\$107.34	\$101.69
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------



MEDLINE NEEDLE-FREE IV ADMINISTRATION SET

- 3 needle-free ports
- Latex-free, non-DEHP

109" Needle-Free IV Set

DYNDTN1545	15 drops/ml	50/case	\$274.29	\$260.58	\$260.58	\$246.86
------------	-------------	---------	-----------------	-----------------	-----------------	-----------------



B|BRAUN CARESITE® VALVE NEEDLE-FREE IV ADMINISTRATION SET FOR VISTA™ PUMP

111" CARESITE Vista™ Needle-Free Pump Set

Universal spike, integrated free flow protector, roller clamp, SPIN-LOCK® connector, CARESITE valve injection site 6" above distal end

354201	15 drops/ml	50/case	\$202.29	\$192.18	\$192.18	\$182.06
--------	-------------	---------	-----------------	-----------------	-----------------	-----------------



DYNAREX NEEDLE-FREE EXTENSION SET

7" with Needle-Free Luer Lock Connector

7061	15 drops/ml	100/box	\$122.99	\$116.84	\$116.84	\$110.69
------	-------------	---------	-----------------	-----------------	-----------------	-----------------



DYNAREX NEEDLELESS PRN CONNECTOR

- Needleless male adapter
- Not made with natural rubber latex

7054		100/box	\$130.79	\$124.25	\$124.25	\$117.71
------	--	---------	-----------------	-----------------	-----------------	-----------------



IV ADMINISTRATION



ACE 2-WAY AUTOMATIC VALVE

- Plastic
- Sterile
- Allows for transfer of anesthetic
- Disposable
- Completely transparent
- Not autoclavable
- No replacement valves or springs
- Inexpensive, efficient, convenient
- For transferring Brevital from reservoir bottle to patient


0031156	10/pack	\$54.79	\$52.05	\$52.05	\$49.31
---------	---------	---------	---------	---------	---------



ACE 3-WAY DISPOSABLE STOPCOCK

- Plastic
- Sterile
- Disposable
- Economical
- Latex-free

0031161	50/box	\$71.79	\$68.20	\$68.20	\$64.61
---------	--------	---------	---------	---------	---------



MEDFUSION STOPCOCKS

G Small Bore - 3-Way w/ Swivel Male Luer Lock, .035-.065 Priming Volume


5311L	50/case	\$53.79	\$51.10	\$51.10	\$48.41
-------	---------	---------	---------	---------	---------

H Large Bore - 3-Way Hi-Flo™ w/ Male Luer Lock, .126 Priming Volume

4311L	50/case	\$66.49	\$63.17	\$63.17	\$59.84
-------	---------	---------	---------	---------	---------

I 4-Way Hi-Flo™ w/ 33" Standard Bore Extension Tubing w/ Luer Lock

MX92133L	50/case	\$131.29	\$124.73	\$124.73	\$118.16
----------	---------	----------	----------	----------	----------



10



SMITHS MEDICAL JELCO® IV CATHETERS

- FEP polymer construction and "J"-point needle design
- Thin wall catheter provides higher flow rates than standard wall catheters
- Latex-free

4055	18g x 1-1/4"	50/box	\$57.99	\$55.09	\$55.09	\$52.19
4057	20g x 1"	50/box	\$57.99	\$55.09	\$55.09	\$52.19
4056	20g x 1-1/4"	50/box	\$57.99	\$55.09	\$55.09	\$52.19
4050	22g x 1"	50/box	\$57.99	\$55.09	\$55.09	\$52.19
4053	24g x 3/4"	50/box	\$57.99	\$55.09	\$55.09	\$52.19



SMITHS MEDICAL CATHLON® CATHETER

Offers all the features and quality of the JELCO® I.V. Catheter combined with the added feature of a metal hub, allowing for optimal balance, weight and suturing. Only CATHLON® I.V. Catheters, made of clear non-radiopaque, FEP polymer, give you the precise and confident feel of metal.

- Strong metal hub provides optimal balance and weight
- Suture rings for additional securement
- Translucent construction allows for consistent and easy flashback visualization
- Self-venting removable flash plug
- Finely tapered catheter tip eases insertions to minimize tissue trauma
- Thinwall catheter design provides higher flow rates than standard wall catheters, often allowing for the use of a smaller gauge device

CR4426	20g x 1-1/4"	50/box	\$248.29	\$235.88	\$235.88	\$223.46
--------	--------------	--------	----------	----------	----------	----------

IV ADMINISTRATION

TERUMO SURFLO® IV CATHETERS

- Ultra-sharp double-bevel introducer needle plus a medical grade lubricant allows for easier penetration
- Preview chamber gives indication of proper placement
- Flexible thinwall catheter with large inner-diameter
- Latex-free

1632	16g x 1-1/4"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
1832	18g x 1-1/4"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
1851	18g x 2"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
2025	20g x 1"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
2032	20g x 1-1/4"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
2051	20g x 2"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
2225	22g x 1"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
2419	24g x 3/4"	50/box	\$84.79	\$80.55	\$80.55	\$76.31



TERUMO SURFLASH® POLYURETHANE IV CATHETERS

- Made of polyurethane; extremely flexible, kink-resistant and durable
- Improves visualization of difficult venipunctures
- Latex-free

SRFF1451	14g x 2"	50/box	\$89.79	\$85.30	\$85.30	\$80.81
SRFF2025	20g x 1"	50/box	\$89.79	\$85.30	\$85.30	\$80.81
SRFF2032	20g x 1-1/4"	50/box	\$89.79	\$85.30	\$85.30	\$80.81
SRFF2225	22g x 1"	50/box	\$89.79	\$85.30	\$85.30	\$80.81
SRFF2419	24g x 3/4"	50/box	\$89.79	\$85.30	\$85.30	\$80.81

EXEL IV CATHETERS

- Translucent catheter hub and flashback chamber allows for easy detection of blood flashback at vein insertion
- Teflon radiopaque catheter
- "New Point" finished sharp introducer needle
- Latex-free

26735	18g x 1-1/4"	50/box	\$44.29	\$42.08	\$42.08	\$39.86
26741	20g x 1"	50/box	\$44.29	\$42.08	\$42.08	\$39.86
26742	20g x 1-1/4"	50/box	\$44.29	\$42.08	\$42.08	\$39.86
26746	22g x 1"	50/box	\$44.29	\$42.08	\$42.08	\$39.86
26751	24g x 3/4"	50/box	\$44.29	\$42.08	\$42.08	\$39.86



BD™ INSYTE® AUTOGUARD™ SHIELDED IV CATHETERS

The patented BD AutoGuard push-button shielding technology is proven to safeguard healthcare workers from needlestick injuries.

- Radiopaque

381433	20g x 1", Pink	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381434	20g x 1-4/25", Pink	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381534	20g x 1-4/25", Pink winged	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381423	22g x 1", Blue	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381523	22g x 1", Blue winged	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381412	24g x 3/4", Yellow	50/box	\$144.79	\$137.55	\$137.55	\$130.31
381512	24g x 3/4", Yellow winged	50/box	\$144.79	\$137.55	\$137.55	\$130.31



Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

IV ADMINISTRATION



BD SAF-T-INTIMA™ CLOSED IV CATHETER SYSTEM

The BD Saf-T-Intima™ closed IV catheter system incorporates preattached extension tubing to help support reduction in blood exposure.

- Advanced shielding design: A telescoping needle shield passively covers the needle, helping protect clinicians from needlestick injuries
- Removable adapter: A removable PRN adapter allows the use of an ISO - compliant needleless connector
- Closed system: Preattached extension tubing offers a closed system to minimize blood exposure during catheter insertion
- BD Vialon™ biomaterial: The system is constructed with BD Vialon biomaterial, a proprietary catheter material that softens in the vein

383346	18g x 1", Green, Y adapter	25/box	\$130.49	\$123.97	\$123.97	\$117.44
383335	20g x 1", Pink, PRN adapter	25/box	\$130.49	\$123.97	\$123.97	\$117.44
383336	20g x 1", Pink, Y adapter	25/box	\$130.49	\$123.97	\$123.97	\$117.44
383323	22g x 3/4", Blue, Y adapter	25/box	\$130.49	\$123.97	\$123.97	\$117.44
383312	24g x 3/4", Yellow, PRN adapter	25/box	\$130.49	\$123.97	\$123.97	\$117.44

B|BRAUN INTROCAN® SAFETY™ IV CATHETERS

Designed to minimize accidental needlesticks without requiring user activation. A stainless steel clip shields the needle tip. A syringe can be pre-attached to facilitate aspiration and injection during insertion.

- Radiopaque

FEP (Teflon) - Straight

425256002	18g x 1-1/4"	50/box	\$145.79	\$138.50	\$138.50	\$131.21
425254302	20g x 1"	50/box	\$145.79	\$138.50	\$138.50	\$131.21
425253502	20g x 1-1/4"	50/box	\$145.79	\$138.50	\$138.50	\$131.21
425251902	22g x 1"	50/box	\$145.79	\$138.50	\$138.50	\$131.21
425250002	24g x 3/4"	50/box	\$145.79	\$138.50	\$138.50	\$131.21

FEP (Teflon) - Winged

425454602	20g x 1"	50/box	\$149.99	\$142.49	\$142.49	\$134.99
425451102	22g x 1"	50/box	\$149.99	\$142.49	\$142.49	\$134.99
425450302	24g x 3/4"	50/box	\$149.99	\$142.49	\$142.49	\$134.99

Polyurethane - Straight

425165202	20g x 1"	50/box	\$149.99	\$142.49	\$142.49	\$134.99
425164402	20g x 1-1/4"	50/box	\$149.99	\$142.49	\$142.49	\$134.99
425162802	22g x 1"	50/box	\$149.99	\$142.49	\$142.49	\$134.99
425160102	24g x 3/4"	50/box	\$149.99	\$142.49	\$142.49	\$134.99



10

**GET INSIDER ACCESS
ORDER. TRACK. FAVORITES.
MANAGE. AND MORE.**

SIGN UP FOR WEB ACCESS TODAY

ACESOUTHERN.COM/CUSTOMER/ACCOUNT/CREATE/

ACE SOUTHERN



IV ADMINISTRATION

SMITHS MEDICAL PROTECTIV® AND PROTECTIV PLUS SAFETY IV CATHETERS

- Designed to resist accidental activation with you in control
- Safety encases the needle as the catheter is advanced
- Optically clear, yet radiopaque polyurethane material for better flashback visualization
- Unique V-point needle provides ultra-smooth insertions
- Prominent push-off tab allows easy threading
- Needle guard with patented gasket designed to minimize blood contact during insertion



PROTECTIV Safety IV Catheter (FEP Polymer)

3055	18g x 1 1/4"	50/box	\$129.49	\$123.02	\$123.02	\$116.54
3057	20g x 1"	50/box	\$129.49	\$123.02	\$123.02	\$116.54
3056	20g x 1 1/4"	50/box	\$129.49	\$123.02	\$123.02	\$116.54
3050	22g x 1"	50/box	\$129.49	\$123.02	\$123.02	\$116.54
3053	24g x 3/4"	50/box	\$129.49	\$123.02	\$123.02	\$116.54

PROTECTIV PLUS Safety IV Catheter (Polyurethane)

3065	18g x 1 1/4"	50/box	\$115.79	\$110.00	\$110.00	\$104.21
3067	20g x 1"	50/box	\$115.79	\$110.00	\$110.00	\$104.21
3066	20g x 1 1/4"	50/box	\$115.79	\$110.00	\$110.00	\$104.21
3060	22g x 1"	50/box	\$115.79	\$110.00	\$110.00	\$104.21
SMM3063	24g x 3/4"	50/box	\$115.79	\$110.00	\$110.00	\$104.21

SMITHS MEDICAL VIAVALVE™ SAFETY IV CATHETERS WITH BLOOD CONTROL

Designed to protect against blood exposure, contamination and needlestick injury. ViaValve safety IV catheters provide blood control to help reduce the risk of blood exposure and contamination, while maintaining a clean access site.

- A valve inside the catheter hub impedes blood flow during initial venipuncture and needle withdrawal
- The introducer needle retracts into the needle housing and the valve closes to keep the access site clean
- The luer fitting of an IV line compresses and opens the valve allowing unrestricted flow of fluids

Safety Features

- Audible "click" of the needle guard when safety is activated signifies the needle is contained for handling and disposal
- Safety guard surrounds the needle to help prevent unintended contact with the needle

Introducer and Needle Housing

- Allows for one-handed or two-handed insertion, threading and disconnection

3268	14g x 1-1/4"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
SMM3262	16g x 1-1/4"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
3265	18g x 1-1/4"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
SMM3267	* 20g x 1"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
3266	* 20g x 1-1/4"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
3260	* 22g x 1"	50/box	\$133.99	\$127.29	\$127.29	\$120.59
3263	* 24g x 5/8"	50/box	\$133.99	\$127.29	\$127.29	\$120.59

*Features a FLASH-VUE® notched needle for early flash detection.



EXEL INJECTION PLUG

- Latex-free
- With cap
- Caution: Do not access with needles larger than 20 gauge

26539	Luer lock	50/box	\$20.49	\$19.47	\$19.47	\$18.44
-------	-----------	--------	---------	---------	---------	---------



KAWASUMI INTERMITTENT INJECTION SITE

- Allows clinician to perform multiple injections when attached to a winged infusion set

INT01	Male luer connector, latex-free	50/box	\$21.79	\$20.70	\$20.70	\$19.61
-------	---------------------------------	--------	---------	---------	---------	---------



Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

IV ADMINISTRATION



B|BRAUN INTERMITTENT INJECTION CAP

- Opaque yellow body, proximal self-sealing injection port, distal male luer lock connector with recessed luer taper
- Latex-free, priming volume: .16ml, length 1-1/8"

6542374238010	100/box	\$97.79	\$92.90	\$92.90	\$88.01
---------------	---------	----------------	----------------	----------------	----------------

MEDIFIL IV FLUSH SOLUTION

- 0.9% Preservative-free sodium chloride

011130	10ml Syringe	60/box	\$49.99	\$47.49	\$47.49	\$44.99
--------	--------------	--------	----------------	----------------	----------------	----------------

TERUMO SURFLO® WINGED INFUSION SETS

- Flexible interlocking wings readily conform to the body's contour
- Angled wing for ease of needle positioning during venipuncture
- Precisely honed needle assures ease of penetration and continued access with a single insertion

With 12" Tubing

19BLK	19g x 3/4"	50/box	\$28.99	\$27.54	\$27.54	\$26.09
21BLK	21g x 3/4"	50/box	\$28.99	\$27.54	\$27.54	\$26.09
22BLK	22g x 3/4"	50/box	\$28.99	\$27.54	\$27.54	\$26.09
23BLK	23g x 3/4"	50/box	\$28.99	\$27.54	\$27.54	\$26.09
25BLK	25g x 3/4"	50/box	\$28.99	\$27.54	\$27.54	\$26.09

Shorty with 3-1/2" Tubing

21BLS	21g x 3/4"	100/box	\$60.99	\$57.94	\$57.94	\$54.89
23BLS	23g x 3/4"	100/box	\$60.99	\$57.94	\$57.94	\$54.89
25BLS	25g x 3/4"	100/box	\$60.99	\$57.94	\$57.94	\$54.89

EXEL WINGED INFUSION SETS

With 12" Tubing

26702	19g x 3/4"	50/box	\$19.79	\$18.80	\$18.80	\$17.81
26704	21g x 3/4"	50/box	\$19.79	\$18.80	\$18.80	\$17.81
26706	23g x 3/4"	50/box	\$19.79	\$18.80	\$18.80	\$17.81
26708	25g x 3/4"	50/box	\$19.79	\$18.80	\$18.80	\$17.81
26709	27g x 3/4"	50/box	\$19.79	\$18.80	\$18.80	\$17.81

TERUMO SURSHIELD® SAFETY WINGED INFUSION SETS

- Unique needle locking mechanism helps reduce risk of needlestick injuries
- One-handed activation meets key criteria for safety and ease of use
- No change in needle insertion technique required
- Latex-free


With 12" Tubing

S19BL	19g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26
S21BLA	21g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26
S23BL	23g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26
S25BL	25g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26

Shorty with 3-1/2" Tubing

S21BLS	21g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26
S23BLS	23g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26
S25BLS	25g x 3/4"	50/box	\$70.29	\$66.78	\$66.78	\$63.26

10



IV ADMINISTRATION

BD™ VACUTAINER SAFETY-LOK™ BLOOD COLLECTION INFUSION SETS

- Larger wing size permits improved positioning during use
- Blood flashback is seen clearly through translucent shield
- Convenient, integrated multiple-sample luer adapter allows specimen collection from a “closed system”
- After use, needle can be easily and completely retracted into the protective shield and securely locked in place

12" Tubing without Adapter

367296	21g x 3/4"	50/box	\$80.79	\$76.75	\$76.75	\$72.71
367297	23g x 3/4"	50/box	\$80.79	\$76.75	\$76.75	\$72.71

12" Tubing with Luer Adapter

367281	21g x 3/4"	50/box	\$86.79	\$82.45	\$82.45	\$78.11
367283	23g x 3/4"	50/box	\$86.79	\$82.45	\$82.45	\$78.11

7" Tubing with Luer Adapter

7287	21g x 3/4"	50/box	\$90.29	\$85.78	\$85.78	\$81.26
------	------------	--------	---------	---------	---------	---------



BD™ VACUTAINER® BLOOD COLLECTION TUBES

Serum

- Red color
- Conventional closure
- Silicone-coated interior

Plastic Tubes, Increased Silica Act Clot Activator

367820	10.0ml draw volume	100/box	\$29.29	\$27.83	\$27.83	\$26.36
368660	6.0ml draw volume	100/box	\$26.49	\$25.17	\$25.17	\$23.84

Glass Tubes, No Additive

366430	10.0ml draw volume	100/box	\$72.99	\$69.34	\$69.34	\$65.69
--------	--------------------	---------	---------	---------	---------	---------

SST™

- Plastic
- Red-grey color
- Conventional closure
- Silicone-coated interior

With Silica Clot Activator, Polymer Gel

367988	8.5ml draw volume	100/box	\$51.29	\$48.73	\$48.73	\$46.16
--------	-------------------	---------	---------	---------	---------	---------

Acid Citrate Dextrose (ACD)

- Glass
- Yellow
- Conventional closure
- Whole blood tube w/ anticoagulant, ACD Sol A

364606	8.5ml draw volume	100/box	\$109.99	\$104.49	\$104.49	\$98.99
--------	-------------------	---------	----------	----------	----------	---------

BD™ VACUTAINER® BLOOD TRANSFER DEVICE

The use of a needle to transfer venous blood from a syringe to a blood collection tube or blood culture bottle is both a dangerous procedure and an OSHA-prohibited practice. The BD Vacutainer Blood Transfer Device was designed with healthcare workers' safety in mind. This single-use device reduces the risk of transfer-related injuries, while maintaining the specimen integrity required for accurate results.

364880		200/case	\$244.49	\$232.27	\$232.27	\$220.04
--------	--	----------	----------	----------	----------	----------

BD™ VACUTAINER® LUER-LOK™ ACCESS DEVICE

Designed for sterile, secure and safer specimen sampling. This device provides the security of a threaded, locking luer connection - the patented BD Vacutainer® Luer-Lok™ that replaces a luer slip device. The product is also compatible with a female luer connection or needleless IV site designed for luer-lock access, and luer locking Foley Catheter sampling ports.

364902		200/case	\$247.49	\$235.12	\$235.12	\$222.74
--------	--	----------	----------	----------	----------	----------

BD™ VACUTAINER® ONE USE HOLDER

Compatible with the entire BD™ Vacutainer® blood collection system.

364815		250/bag	\$17.99	\$17.09	\$17.09	\$16.19
--------	--	---------	---------	---------	---------	---------



IV ADMINISTRATION

BD™ VACUTAINER® ECLIPSE NEEDLE

Effective way to collect blood and help reduce the possibility of accidental needle stick, allowing for one handed activation technique.



368607	21g x 1 1/4"	48/box	\$29.49	\$28.02	\$28.02	\$26.54
--------	--------------	--------	---------	---------	---------	---------

EXEL LUER LOCK HOLDER WITH PRE-ATTACHED MULTIPLE SAMPLE ADAPTER

- Does not contain natural rubber latex
- Sterile
- Individually packed
- Luer lock connection
- Blue color thread



26532	5-7.5ml Holder	200/case	\$168.99	\$160.54	\$160.54	\$152.09
-------	----------------	----------	----------	----------	----------	----------

EXEL MULTI-SAMPLE LUER ADAPTER

26537		100/box	\$16.79	\$15.95	\$15.95	\$15.11
-------	--	---------	---------	---------	---------	---------



VENOSCOPE® II

The Venoscope II Transilluminator, The Vein Finder™, allows a clinician to easily locate hard to find veins. A clinician can reduce multiple sticks and increase the first stick success rate toward the 100% mark.

The use of the Venoscope II is a non-invasive procedure and is done by placing the Venoscope II on the surface of the skin in such a way that the high intensity light is projected into the subcutaneous tissue of a patient. The light causes the veins to contrast with the surrounding tissue making it easier to locate veins suitable for IV therapy and blood draws.



VT03	*	each	\$280.79	\$266.75	\$266.75	\$252.71
------	---	------	----------	----------	----------	----------

Disposable Covers

DPC050		50/box	\$22.99	\$21.84	\$21.84	\$20.69
--------	--	--------	---------	---------	---------	---------

*Requires 3 AA batteries, not included

DYNAREX IV START KIT

- Latex-free
- Sterile

Contents:

- (1) Alcohol prep pad
- (2) Gauze sponge, 2" x 2"
- (1) PVP prep pad
- (1) Tourniquet, 3/4" x 18"
- (1) Transpore tape, 3/4" x 18"
- (1) ID label
- (1) Transparent dressing



4692		50/case	\$64.99	\$61.74	\$61.74	\$58.49
------	--	---------	---------	---------	---------	---------

ACE LATEX-FREE TOURNIQUET

- Blue
- Flat



0031153	1" x 18"	10/pack	\$6.79	\$6.45	\$6.45	\$6.11
---------	----------	---------	--------	--------	--------	--------

TEMPO NON-LATEX TOURNIQUET

Smooth tourniquet offers high end comfort, ideal for patients who require gentle care due to fragile or at-risk skin.



TS118	1" x 18"	250/box	\$40.99	\$38.94	\$38.94	\$36.89
-------	----------	---------	---------	---------	---------	---------

BD™ VACUTAINER® STRETCH LATEX-FREE TOURNIQUET

Reduce the danger of cross-contamination between patients and safeguard both you and your patients from latex allergens.



367203	1" x 18"	25/box	\$15.99	\$15.19	\$15.19	\$15.99
--------	----------	--------	---------	---------	---------	---------

HEALTH CARE LOGISTICS LATEX-FREE TOURNIQUETS ON A ROLL

Tourniquets are packaged on rolls of 25, and are pre-notched every 18" for easy removal. Just roll out and tear at the notch.



10311B	Neon blue	25/roll	\$35.29	\$33.53	\$33.53	\$31.76
10311G	Neon green	25/roll	\$35.29	\$33.53	\$33.53	\$31.76
10311PK	Neon pink	25/roll	\$35.29	\$33.53	\$33.53	\$31.76
10311PP	Neon purple	25/roll	\$35.29	\$33.53	\$33.53	\$31.76

ADC LATEX-FREE VELCRO® TOURNIQUET

ADC341	Child, 9-1/2" x 1"	each	\$3.29	\$3.13	\$3.13	\$2.96
ADC340	Adult, 14-1/4" x 1"	each	\$3.29	\$3.13	\$3.13	\$2.96
ADC342	Leg, 17" x 1-5/8"	each	\$4.49	\$4.27	\$4.27	\$4.04



PROPPER TOURNIQUET

- Adult size
- Cloth/Rubber
- Seraket
- Adjustable



227001	24"	each	\$17.49	\$16.62	\$16.62	\$15.74
--------	-----	------	---------	---------	---------	---------

IV ADMINISTRATION

ACE CUSTOM SURGICAL ARM BOARD

- Can be customized for any surgical chair, including A-dec, Belmont, X Caliber, Decade Plus and Matrix
- Please send in template: instructions below
- Non-refundable

Features

- A** Padded armrest is comfortable to patient and cleans easily
- B** 2" wide Velcro straps open and close easily while securely immobilizing patient arm
- C** Spring loaded control for elevation, can be elevated to 60°
- D** Control for advancing and rotating arm
- E** All aluminum and stainless steel construction
- F** High density, non-slip cork sheet assures positive positioning

0031158

Call for Quote




INSTRUCTIONS FOR TEMPLATE MAKING

- Step 1:** Take the padded arm off of the dental chair; leaving just metal frame.
 - Step 2:** Trace the outer shape of the arm onto a piece of cardboard, include all screw holes.
 - Step 3:** Measure from the back of the armrest to the backrest (In its upright position).
 - Step 4:** Include model and make of the chair.
- If you have any questions or need assistance please call customer service.


ACE SIDE ARM ATTACHMENT WITH SYRINGE GUARD FOR 10 CC AND 20 CC SYRINGES

0039089 each **\$258.79** **\$245.85** **\$245.85** **\$232.91**



ACE STAINLESS STEEL ARM BOARD BRACKETS

G 0030002	Female bracket	each	\$30.49	\$28.97	\$28.97	\$30.49
H 0030003	Female bracket	each	\$38.29	\$38.29	\$32.55	\$38.29
I 0030004	Male plate	each	\$11.29	\$10.73	\$10.73	\$11.29



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

IV ADMINISTRATION

PROFEX ARMBOARD PAD

- With Velcro straps for attachment to chair
- Flexible



0194	6" x 24" x 21"	each	\$81.29	\$77.23	\$77.23	\$73.16
------	----------------	------	---------	---------	---------	---------

PRECISION DYNAMICS DISPOSABLE IV BOARDS

- Cardboard and foam with plastic covering
- Straps not included



4500311MPF	Child-2" x 9"	each	\$3.79	\$3.60	\$3.60	\$3.41
4500411MPG	Wrist, 3" x 9"	each	\$4.29	\$4.08	\$4.08	\$3.86
4500211MPG	Adult, 3" x 17 1/4"	each	\$3.49	\$3.32	\$3.32	\$3.14

MEDLINE DISPOSABLE IV ARMBOARD

Made of vinyl, these boards have a soft, spongy core for patient comfort and a firm covering for stability. Latex-free.



NON24293	3" x 9"	each	\$1.99	\$1.89	\$1.89	\$1.79
----------	---------	------	--------	--------	--------	--------

MEDLINE FOAM ARMBOARD STRAP

- Disposable



NON081540	1 1/2" x 32"	48/case	\$103.79	\$98.60	\$98.60	\$93.41
-----------	--------------	---------	----------	---------	---------	---------


WALL MEDICAL ADULT ELBOW IMMOBILIZER



WEB100		each	\$138.49	\$131.57	\$131.57	\$124.64
--------	--	------	----------	----------	----------	----------

ACE ELBOW IMMOBILIZERS

- Maintains arm in rigid position during administration of IV solutions



8372161	Child, 7"	each	\$24.99	\$23.74	\$23.74	\$22.49
8372160	Adult, 11"	each	\$24.99	\$23.74	\$23.74	\$22.49

ACE KORBOARD WITH VELCRO STRAPS

- Velfoam pad for patient comfort
- X-ray transparent
- Latex-free



4842003	Child, 7 1/4" x 2 1/2"	each	\$34.99	\$33.24	\$33.24	\$31.49
4842002	Adult, 9" x 3 1/2"	each	\$39.49	\$37.52	\$37.52	\$35.54

MEDLINE GEL DONUT HEAD POSITIONERS

Gel head positioners help protect the head, face, neck, and ears and can be used during a variety of procedures.

Pediatric - 5-1/2"x2-1/4"x1-1/4"

PCD082004		each	\$69.49	\$66.02	\$66.02	\$62.54
-----------	--	------	---------	---------	---------	---------

Adult - 8"x3"x1-1/4"

PCD082000		each	\$70.49	\$66.97	\$66.97	\$63.44
-----------	--	------	---------	---------	---------	---------

DYNAREX IV CATHETER DRESSING

Secures and Stabilizes IV Catheters on Patients

- Offers ability to write on foam border
- Sterile, single use only
- Comfortable, and has a low profile
- Waterproof



4609		100/box	\$47.79	\$45.40	\$45.40	\$43.01
------	--	---------	---------	---------	---------	---------

CONMED VENI-GARD® STERILE IV STABILIZATION DRESSING

- Stabilizing dressing for peripheral lines
- Enhances patient comfort
- Adhesive on foam border allows patients to maintain ADL's
- High MVTR (Moisture Vapor Transmission Rate)



7054431	Veni-gard	100/box	\$60.79	\$57.75	\$57.75	\$54.71
---------	-----------	---------	---------	---------	---------	---------

7304432	Veni-gard Junior	100/box	\$75.79	\$72.00	\$72.00	\$68.21
---------	------------------	---------	---------	---------	---------	---------

3M TEGADERM™ STERILE TRANSPARENT DRESSING

Sterile, waterproof, latex-free, transparent film dressing used to cover and secure IV sites, protect at-risk skin, or as a primary or secondary dressing over acute wounds, surgical sites, and chronic wounds.



Easy Application Frame Style with Label

1624W	2-3/8" x 2-3/4"	100/box	\$51.79	\$49.20	\$49.20	\$46.61
-------	-----------------	---------	---------	---------	---------	---------

1626W	4" x 4-3/4"	50/box	\$84.79	\$80.55	\$80.55	\$76.31
-------	-------------	--------	---------	---------	---------	---------

Original Frame Style

1626	4" x 4-3/4"	50/box	\$86.99	\$82.64	\$82.64	\$78.29
------	-------------	--------	---------	---------	---------	---------

IV ADMINISTRATION

ACE INFUSER BAG

- Use with either a 500ml or 1000ml flexible bag
- For faster fluid flow
- Made of heavy nylon
- Includes bag, gauge, bladder, and bulb



0035100	each	\$183.99	\$174.79	\$174.79	\$165.59
---------	------	----------	----------	----------	----------

DYNAREX 2-HOOK IV STAND

Economic with Removable Hook


- 4 rubber wheels provide for a smooth transport and maneuverability
- Comes with a removable 2-hook hanger
- Locking collar allows for easy height adjustment



10554	each	\$64.79	\$61.55	\$61.55	\$58.31
-------	------	---------	---------	---------	---------

DYNAREX 2-HOOK IV STAND WITH ALUMINUM BASE

- Chrome-plated steel pole provides strength and durability
- Composite 5-leg base and 5 swivel ball bearing casters provide easy, smooth maneuvering on all surfaces
- Comes with a removable 2-hook hanger
- Locking collar allows for easy height adjustment



DY10555	each	\$113.29	\$107.63	\$107.63	\$101.96
---------	------	----------	----------	----------	----------

STEELCRAFT 4-HOOK IV STAND

- Five point all steel base for greater stability
- 2" ball bearing swivel casters
- Full height extension from 60" to 110"
- Shipped completely knocked down in two cartons with simple assembly



399	each	\$178.49	\$169.57	\$169.57	\$160.64
-----	------	----------	----------	----------	----------

BREWER ADD-A-TRAY

- Holds instruments or patient items/15 lb maximum
- Removable stainless steel tray
- Tray size: 9-3/4" x 13-1/2" x 3/4"
- Clamp sold separately below



43431	each	\$138.79	\$131.85	\$131.85	\$124.91
-------	------	----------	----------	----------	----------

Easy Application Frame Style with Label

43433	each	\$98.79	\$93.85	\$93.85	\$88.91
-------	------	---------	---------	---------	---------



INTRAVENOUS FLUIDS

0.9% Sodium Chloride*

50ml Bag	80/case	\$225.79	\$214.50	\$214.50	\$203.21
100ml Bag	80/case	\$249.29	\$236.83	\$236.83	\$224.36
150ml Bag	32/case	\$279.99	\$265.99	\$265.99	\$251.99
250ml Bag	24/case	\$197.99	\$188.09	\$188.09	\$178.19
500ml Bag	24/case	\$171.79	\$163.20	\$163.20	\$154.61
1000ml Bag	12/case	\$142.99	\$135.84	\$135.84	\$128.69

0.45% Sodium Chloride*

500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Bag	12/case	\$121.29	\$115.23	\$115.23	\$109.16

Lactated Ringer's*

250ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Bag	12/case	\$109.79	\$104.30	\$104.30	\$98.81

Sterile Water*

250ml Bag	24/case	\$182.49	\$173.37	\$173.37	\$164.24
500ml Bag	24/case	\$182.49	\$173.37	\$173.37	\$164.24
1000ml Bag	14/case	\$148.79	\$141.35	\$141.35	\$133.91

5% Dextrose*

50ml Bag	80/case	\$284.29	\$270.08	\$270.08	\$255.86
100ml Bag	4/pack	\$14.49	\$13.77	\$13.77	\$13.04
100ml Bag	80/case	\$258.49	\$245.57	\$245.57	\$232.64
150ml Bag	32/case	\$255.29	\$242.53	\$242.53	\$229.76
250ml Bag	36/case	\$225.49	\$214.22	\$214.22	\$202.94
500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Bag	12/case	\$138.49	\$131.57	\$131.57	\$124.64

Dextrose 5% Sodium Chloride 0.9%*

500ml Bag	24/case	\$165.99	\$157.69	\$157.69	\$149.39
1000ml Bag	12/case	\$118.29	\$112.38	\$112.38	\$106.46

Dextrose 5% Sodium Chloride 0.45%*

250ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Bag	12/case	\$105.29	\$100.03	\$100.03	\$94.76

Dextrose 5% Sodium Chloride 0.3%*

500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
-----------	---------	----------	----------	----------	----------

Dextrose 5% Sodium Chloride 0.225%*

250ml Bag	24/case	\$196.99	\$187.14	\$187.14	\$177.29
500ml Bag	24/case	\$196.99	\$187.14	\$187.14	\$177.29

Dextrose 5% Lactated Ringer's*

500ml Bag	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Bag	12/case	\$121.79	\$115.70	\$115.70	\$109.61

*A freight charge of \$8.00 will be added to the price of this product.

PATIENT DIAGNOSTIC

EDAN X SERIES PATIENT MONITORS

Based on its expandability of monitoring parameters, and stunning ultra-slim design, the X series monitors fulfill your primary clinical requirements in various scenarios. They cover all patient types, realizing from neonatal to adult monitoring.

EDAN G2 CO₂ (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB[®] algorithm with intelligent CO₂ pseudo wave identification technology
- Automatic zero calibration within seconds
- Accessories for all patient types

ECG

- Customizable filter for better performance
- Pacemaker detection and Defibrillator/ESU protection
- Industry leading iSEAP[™] algorithm for accurate arrhythmia detection
- SEMIP[®] algorithm with 208 ECG findings over age/gender diversities

NIBP

- Dual dust filter design makes no blockage inside and provides accurate NIBP readings
- Unique cleaning mode for routine maintenance
- iCUFS[™] algorithm with smart deflation technology

SpO₂

- iMAT[™] algorithm with outstanding motion resistance and low perfusion resistance performance
- Reference reading of Perfusion Index (PI) from 0 to 10 according to perfusion changes
- Simultaneous measurements of SpO₂ and NIBP of the same limb

Warranty

- 3-year warranty on unit
- 1-year warranty on accessories

Includes:

- Reusable SpO₂ finger sensor (8')
- Large adult cuff (33-43cm)
- Regular adult cuff (27-35cm)
- Small adult/child cuff (20.5-28cm)
- NIBP hose
- 3-lead ECG cable
- ECG electrodes (3)
- Skin temperature probe
- User manual
- CO₂ sampling cannula (1)
- CO₂ and O₂ divided cannula (1)
- Straight sampling line (1)
- CO₂ water trap (1)
- Power supply
- Rechargeable lithium-ion battery
- Ground cable
- Small screwdriver
- Printer paper roll

Patient Monitor 10.1" Touchscreen, ECG, SpO₂, NIBP, CO₂, RR, Printer

X10	each	\$4,595.29	\$4,365.53	\$4,365.53	\$4,135.76
-----	------	------------	------------	------------	------------

Patient Monitor 12.1" Touchscreen, ECG, SpO₂, NIBP, CO₂, RR, Printer

X12	each	\$4,989.29	\$4,739.83	\$4,739.83	\$4,490.36
-----	------	------------	------------	------------	------------



iM60GT pictured

EDAN iM60 PATIENT MONITORS

- 10.4" color TFT-LCD touchscreen with maximum 11 waveforms
- Touch & Configure: Intuitive operation by clicking on specific parameter or waveform to configure
- Connectivity with: OMS Vision, all DSN products (including Oral Surgery Exec and PerioExec)
- Pacemaker detection
- Electrosurgical interference proof
- Defibrillation protection and defibrillation synchronization
- Pitch tone (pulse-tone modulation)
- Unique iSEAP algorithm which is specially optimized for arrhythmia patients
- Dual-mode anti-interference pulse oximetry
- Clinically validated NIBP
- VGA output
- HL7 support via XML files
- Suitable for adult, pediatric and neonatal patients
- Rechargeable lithium-ion battery (up to 7 hours)
- 2-year warranty; 1-year warranty on finger sensor

Features Respironics CO₂ Sidestream Module w/Accessories for Intubated/Non-Intubated Patients

- Proven ruggedness-Solid state design
- No need to calibrate on a regular basis
- Sidestream sampling rate of 50 ml/min
- The iM60 does not require dehumidification tubing and accepts standard male luer connections for EtCO₂ monitoring.

Both iM60's include:

- Adult reusable SpO₂ sensor (8')
- Large adult cuff (34-43cm)
- Adult cuff (27-35cm)
- Small adult/child cuff (20.5-28cm)
- NIBP tube (3m) w/connector
- 3-lead ECG cable
- Rechargeable Lithium-Ion battery
- Power supply
- Printer paper roll
- User manual
- ECG electrodes (3)

iM60GT also includes:

- Skin temperature probe
- CO₂ water trap (1)
- CO₂ sampling cannula (1)
- CO₂ and O₂ divided cannula (1)
- Straight sampling line (1)

EDAN iM60 Patient Monitors

w/ 3-lead ECG, EDAN SpO₂, NIBP, PR, Touchscreen, Thermal Recorder

IM60PT	each	\$2,829.29	\$2,687.83	\$2,687.83	\$2,546.36
--------	------	------------	------------	------------	------------

w/ 3-lead ECG, RESP, EDAN SpO₂, NIBP, PR, 2-TEMP, Respironics CO₂, Touchscreen and Thermal Recorder

IM60GT	each	\$4,430.29	\$4,208.78	\$4,208.78	\$3,987.26
--------	------	------------	------------	------------	------------

Data Management Kit for M3P and iM60GT

- This software allows transfer of data from the patient monitor to a computer
 - Does not direct connect, manual download of information via USB needed
- Includes: Data Management Software, User Manual, USB Sentinel

0204240050	each	\$154.49	\$146.77	\$146.77	\$139.04
------------	------	----------	----------	----------	----------

Data Management Kit for M3P, iM60GT, X10 and X12

- This software allows transfer of data from the patient monitor to a computer
 - Direct connect from monitor to computer
- Includes: CD with Software, User Manual, USB Dongle

8360261775	each	\$446.29	\$423.98	\$423.98	\$401.66
------------	------	----------	----------	----------	----------



PATIENT DIAGNOSTIC

EDAN M3 VITAL SIGNS MONITOR



- 5.6" high resolution color TFT-LCD screen
- Lightweight & portable design
- User-friendly interface
- Auto/Manual/Continuous/Average BP mode, Average NIBP feature: taking average NIBP reading out of 3 or 5 measurements
- Real-time measurements and trend display
- Built-in rechargeable lithium-ion battery for 8 hours of continuous work
- USB flash disk support
- Suitable for adult, pediatric and neonatal patients
- 2-year warranty; 1-year warranty on accessories

Includes:

- Adult reusable SpO₂ sensor
- Large adult cuff (33-43cm)
- Adult cuff (27-35cm)
- Small adult/child cuff (20.5-28cm)
- NIBP hose
- Rechargeable Lithium-Ion battery
- Power cord
- Printer paper roll
- User manual

With NIBP, SpO₂ and printer

M3P	each	\$2,165.29	\$2,057.03	\$2,057.03	\$1,948.76
-----	------	------------	------------	------------	------------

EDAN iM3 VITAL SIGNS MONITOR



With Touchscreen, SpO₂, NIBP, PR and Printer

Combining together a stunning ultra-slim design with a large capacitive touch screen, and a rich set of leading functions, the iM3 takes the lead and presents itself at the beginning of a new generation of vital signs monitors. iM3 serves for both long term monitoring and spot check applications.

Redefinition of Adaptive Work Modes

Round Mode: Offers a special work mode that is designed to make daily use more efficient. The data can be easily uploaded through a USB connection or downloaded from the EMR system. The color-coded list comes with a square windows interface design so you can view patients status and information at a quick glance.

Monitoring Mode: For continuous patient monitoring with real-time data, real-time alarms, and history trend.

Spot Mode: For spot-check applications, such as out-patient, emergency, etc. Real-time vitals are acquired and displayed, no alarm will be seen or heard.

Warranty

- 2-year warranty on unit
- 1-year warranty on finger sensor

With NIBP, SpO₂, PR and printer

iM3P	each	\$2,265.29	\$2,152.03	\$2,152.03	\$2,038.76
------	------	------------	------------	------------	------------

EDAN PATIENT MONITOR ACCESSORIES

	X10	X12	iM60PT	iM60GT	M3P/iM3	ITEM NO.	UOM				
8' Reusable Adult Finger Sensor						0201210119	each	\$109.29	\$103.83	\$103.83	\$98.36
8' Reusable Adult Finger Sensor (7-Pin)						0257225029	each	\$114.79	\$109.05	\$109.05	\$103.31
Water Trap						0201210520	each	\$30.79	\$29.25	\$29.25	\$27.71
Roll Stand w/ Basket and Mounting Plate						MT20760	each	\$429.79	\$408.30	\$408.30	\$386.81
Roll Stand w/ Basket and Mounting Plate						MT20750	each	\$429.79	\$408.30	\$408.30	\$386.81
Printer Paper w/ Grid						015778035	12/box	\$48.99	\$46.54	\$46.54	\$44.09
3m NIBP Hose						0157471005012	each	\$57.29	\$54.43	\$54.43	\$51.56
10' NIBP Hose						VNIBPTUBE	each	\$61.79	\$58.70	\$58.70	\$55.61
Blood Pressure Cuff - Infant 16-21.5cm						0157040212	each	\$34.79	\$33.05	\$33.05	\$31.31
Blood Pressure Cuff - Child/Small Adult 20.5 - 28cm						0157040211	each	\$34.79	\$33.05	\$33.05	\$31.31
Blood Pressure Cuff - Adult 27 - 35cm						0157040205	each	\$34.79	\$33.05	\$33.05	\$31.31
Blood Pressure Cuff - Large Adult 33 - 43cm						0157040210	each	\$34.79	\$33.05	\$33.05	\$31.31
Blood Pressure Cuff - Thigh 42 - 54cm						CM1305	each	\$47.29	\$44.93	\$44.93	\$42.56

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

PATIENT DIAGNOSTIC

**CRITICARE
NGENUITY™ SERIES
PATIENT MONITORS**



Basic vital signs monitoring in a simple-to-use, lightweight, and portable design. Focus on your patients, not your monitor. The nGenuity patient monitor is the perfect choice for busy practices demanding accurate, reliable performance. Compact and lightweight, the nGenuity patient monitor can monitor adults, pediatrics and neonates.

- Innovative user interface features simplified menus and dedicated function keys for timely set-ups
- Big bright display with large numerics can easily be seen, even from across the room
- Easy-to-read 10.4" active color TFT display
- Lightweight design and versatile mounting options
- Proprietary ComfortCuff™ NIBP measures on inflation for faster, more comfortable readings
- 2-year warranty, 1-year warranty on finger sensor

8100E1 and 8100EP1 Includes:

- 474 - Child/Small Adult Blood Pressure Cuff (18-26cm)
- 475 - Adult Blood Pressure Cuff (25-35cm)
- 476 - Large Adult Blood Pressure Cuff (33-47cm)
- CRI706 - NIBP Hose (10 ft. Straight)
- 93410DN - Adult Finger (Clip) Reusable Sensor (10 ft.)
- 527 - ECG Electrodes
- 1073S - ECG Cable Assembly
- 938FNC - WaterChek™ Plus (2 ea.)
- 625N - Sample Line (5 ea.)
- 989 - Power Cord (Down Angle North American)
- Operator's Manual (CD Version)

w/ ECG, SpO₂, NIBP, sidestream EtCO₂, temperature* and respiratory rate

8100E1	\$6,254.99	\$5,942.24	\$5,942.24	\$5,629.49
---------------	-------------------	-------------------	-------------------	-------------------

w/ ECG, SpO₂, NIBP, sidestream EtCO₂, temperature*, respiratory rate and printer

8100EP1	\$6,386.99	\$6,067.64	\$6,067.64	\$5,748.29
----------------	-------------------	-------------------	-------------------	-------------------

* Does not include temperature probe

8100E and 8100EP Includes:

- 474 - Child/Small Adult Blood Pressure Cuff (18-26cm)
- 475 - Adult Blood Pressure Cuff (25-35cm)
- 476 - Large Adult Blood Pressure Cuff (33-47cm)
- CRI706 - NIBP Hose (10 ft. Straight)
- 93410DN - Adult Finger (Clip) Reusable Sensor (10 ft.)
- 527 - ECG Electrodes
- 1073S - ECG Cable Assembly
- 989 - Power Cord (Down Angle North American)
- Operator's Manual (CD Version)

w/ ECG, SpO₂, NIBP, temperature* and respiratory rate

8100E	\$3,674.99	\$3,491.24	\$3,491.24	\$3,307.49
--------------	-------------------	-------------------	-------------------	-------------------

w/ ECG, SpO₂, NIBP, temperature*, respiratory rate and printer

8100EP	\$4,391.29	\$4,171.73	\$4,171.73	\$3,952.16
---------------	-------------------	-------------------	-------------------	-------------------

**CRITICARE
VITALCARE™
506N3 SERIES
PATIENT MONITORS**



Proprietary ComfortCuff™ NIBP, pulse oximetry, and predictive temperature in one compact, easily transported monitor. Spot-check or continuously monitor three critical parameters on neonates to adults with a single monitor. Designed to travel everywhere you do, the VitalCare 506N3 series patient monitor with integrated carry-handle attaches quickly and easily to an IV pole mount or rolling stand.

- Smart battery management allows up to 100 NIBP measurements on a single charge
- Proprietary ComfortCuff NIBP measures on inflation for faster, more comfortable readings
- Audible and visual alarms are adjustable, providing instant notification
- Optional internal printer for on-demand, automatic interval trends or multi-patient reports
- 2-year warranty

506NP3 Includes:

- 474 - Child/Small Adult Blood Pressure Cuff (18-26cm)
- 475 - Adult Blood Pressure Cuff (25-35cm)
- 476 - Large Adult Blood Pressure Cuff (33-47cm)
- CRI706 - NIBP Hose (10 ft. Straight)
- 989 - Power Cord (Down Angle North American)
- Operator's Manual (CD Version), Help Card w/ NIBP and printer

506NP3	\$2,071.49	\$1,967.92	\$1,967.92	\$1,864.34
---------------	-------------------	-------------------	-------------------	-------------------

506DNP3 Includes:

- 474 - Child/Small Adult Blood Pressure Cuff (18-26cm)
- 475 - Adult Blood Pressure Cuff (25-35cm)
- 476 - Large Adult Blood Pressure Cuff (33-47cm)
- CRI706 - NIBP Hose (10 ft. Straight)
- 93410DN - Adult Finger (Clip) Reusable Sensor (10 ft.)
- 989 - Power Cord (Down Angle North American)
- Operator's Manual (CD Version), Help Card w/ NIBP, CSI SpO₂ and printer

506DNP3	\$2,426.99	\$2,305.64	\$2,305.64	\$2,184.29
----------------	-------------------	-------------------	-------------------	-------------------

506DNVP3 Includes:

- 474 - Child/Small Adult Blood Pressure Cuff (18-26cm)
- 475 - Adult Blood Pressure Cuff (25-35cm)
- 476 - Large Adult Blood Pressure Cuff (33-47cm)
- CRI706 - NIBP Hose (10 ft. Straight)
- 93410DN - Adult Finger (Clip) Reusable Sensor 10 ft.
- 1544 - TurboTemp Oral Temperature Probe
- 446 - Probe Covers
- 202089 Isolation Well for Oral Probe
- 989 - Power Cord (Down Angle North American)
- Operator's Manual (CD Version), Help Card w/ NIBP, CSI SpO₂, turbo temperature and printer

506DNVP3	\$3,265.29	\$3,102.03	\$3,102.03	\$2,938.76
-----------------	-------------------	-------------------	-------------------	-------------------

accessories on next page ...

**CRITICARE SYSTEMS, INC.
MONITOR ACCESSORIES**

		8100 nGenuity™ Series	506N3 Series					
Recorder Paper	5 rolls/box			CSI553	\$63.29	\$60.13	\$60.13	\$56.96
USB Serial Port Adapter, Used with CRI1088 for EMR Interface	each			1550	\$125.29	\$119.03	\$119.03	\$112.76
6' Modem Cable Used w/ 1550 for EMR Interface	each			CRI1088	\$131.29	\$124.73	\$124.73	\$118.16
6' Modem Cable Used w/ 1550 for EMR Interface	each			87016B002	\$76.99	\$73.14	\$73.14	\$69.29
Roll Stand w/ Basket	each			CRI1091	\$881.99	\$837.89	\$837.89	\$793.79
Roll Stand w/ Basket	each			1043	\$801.99	\$761.89	\$761.89	\$721.79
Carrying Case	each			CRI1092	\$275.79	\$262.00	\$262.00	\$248.21
Battery	each			PN80519B001	\$79.79	\$75.80	\$75.80	\$71.81
Battery	each			80518B001	\$196.49	\$186.67	\$186.67	\$176.84
Temperature Probe, Oral - Blue	each			CRI1544	\$267.49	\$254.12	\$254.12	\$240.74
Temperature Probe, Rectal - Red	each			1545	\$267.49	\$254.12	\$254.12	\$240.74
Temperature Probe Covers	200/box			447	\$64.79	\$61.55	\$61.55	\$58.31
Temperature Probe, Rectal/Esophageal - Disposable	each			438	\$16.99	\$16.14	\$16.14	\$15.29
Temperature Probe, Rectal/Esophageal - Reusable	each			420	\$146.99	\$139.64	\$139.64	\$132.29
Temperature Interconnect Cable	each			439	\$96.29	\$91.48	\$91.48	\$86.66
Blood Pressure Cuff - Infant (10-19cm)	each			473	\$56.79	\$53.95	\$53.95	\$51.11
Blood Pressure Cuff - Child/Small Adult (18-26cm)	each			474	\$61.49	\$58.42	\$58.42	\$55.34
Blood Pressure Cuff - Adult (25-35cm)	each			475	\$68.49	\$65.07	\$65.07	\$61.64
Blood Pressure Cuff - Large Adult (33-47cm)	each			476	\$91.99	\$87.39	\$87.39	\$82.79
Blood Pressure Cuff - Thigh (46-66cm)	each			477	\$130.49	\$123.97	\$123.97	\$117.44
Blood Pressure Hose - 10' Straight	each			CRI706	\$57.79	\$54.90	\$54.90	\$52.01
Blood Pressure Hose - 10' Coiled	each			CRI707	\$76.99	\$73.14	\$73.14	\$69.29
Adult Pulse Oximeter Finger Sensor - 10' Reusable w/ DB9 Connection	each			93410DN	\$225.49	\$214.22	\$214.22	\$202.94
Extension Cable for Pulse Oximeter Sensor - 10' w/ DB9 Connection	each			518DD	\$207.79	\$197.40	\$197.40	\$187.01
ECG Wrist/Limp Clips	pair			CRI1166	\$125.99	\$119.69	\$119.69	\$113.39
ECG Cable Assembly - 3 Lead Snap	each			1073S	\$168.29	\$159.88	\$159.88	\$151.46
WaterChek™ Plus Trap/Filter †	each			938FNC	\$16.49	\$15.67	\$15.67	\$14.84
Nasal Cannula CO ₂ /O ₂ - 7' Female Luer	10/box			CRI628	\$181.79	\$172.70	\$172.70	\$163.61
Sample Line - 8' Male/Male Luer Lock	25/box			625N	\$139.79	\$132.80	\$132.80	\$125.81

Note: More Criticare monitor accessories are available.
Contact Customer Service for availability and pricing.

† Recommended to be changed weekly or as needed

Looking for additional accessories for your Capnography Monitor?
See the Anesthesia Section for: Sampling Lines • Dual Nasal Hoods • CO₂ Sampling Cannulas!

Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price

PATIENT DIAGNOSTIC



11

ZOE NIGHTINGALE™ PPM3 VITAL SIGNS MONITOR

The Nightingale PPM3 provides non-invasive vital signs monitoring, including 3 or 5-lead ECG, pulse oximetry, blood pressure, respiration (EtCO2 or impedance), and temperature. An optional configuration includes two channels of invasive pressure monitoring.

The Nightingale PPM3 patient monitor meets the needs of the oral surgery and recovery centers with reliable monitoring of physiological parameters in a clean, compact design. The bright, high-contrast display, together with the simple screen and menu layouts make the device quick to learn and easy to operate.

Product Highlights:

- Fast, accurate, continuous monitoring of key physiological parameters
- Monitors ECG, Respiration, CO₂, SpO₂, NIBP, and Temperature
- Heart Rate (HR) measured from ECG, SpO₂, or NIBP
- ECG beat detect tones with variable pitch (based on SpO₂)
- Connectivity with major practice management systems (ie. Care Stream WinOMS, DSN, etc.)
- Future safe — software upgradable via USB stick
- Printer firmly attached to back of the monitor and prints clear waveform and vitals information
- Portable, Long-Life Lithium Ion battery
- Manufactured in the U.S.A.
- Warranty - 3 years, including the battery

Easy to Use:

- Large, bright display with wide viewing angle
- Simple menus and layout for easy operations
- No CO₂ water trap, no calibration required

Nightingale™ PPM3 Monitor

w/ ECG, SpO₂, NIBP, temperature, respiratory rate and printer

1791000	\$3,942.79	\$3,745.65	\$3,745.65	\$3,548.51
---------	------------	------------	------------	------------

Included with the 1791000 Monitor Purchase:

- Reusable SpO₂ Sensor
- SpO₂ Extender Cable (8')
- NBP Hose (10')
- NBP Cuff, Child (13-20cm)
- NBP Cuff, Small Adult (18-26cm)
- NBP Cuff, Adult (26-35cm)
- ECG 3-lead cable (10')
- Printer 2"/50mm
- Table Stand with Recorder
- Power Cord, US (10 ft)

Nightingale™ PPM3 Monitor

w/ ECG, SpO₂, NIBP, Masimo EtCO₂, temperature, respiratory rate, printer and two Masimo CO₂ Nomoline adapters

1791001	\$5,270.79	\$5,007.25	\$5,007.25	\$4,743.71
---------	------------	------------	------------	------------

Included with the 1791001 Monitor Purchase:

- Masimo ISA EtCO₂ Sidestream Analyzer
- Masimo Nomoline Adapter (2 pcs)
- Adult Nasal Cannula CO₂ Sampling (5 pcs)
- Reusable Envitec SpO₂ Sensor
- SpO₂ Extender Cable (8')
- NBP Hose (10')
- NBP Cuff, Small Adult (18-26cm)
- NBP Cuff, Adult (26-35cm)
- NBP Cuff, Large Adult (32-42cm)
- ECG 3-lead cable (10')
- Table Stand with Recorder
- 50mm Printer for PPM3 Monitor Cable, PPM3-Phasein Interface
- Mounting Bracket Kit for Masimo ISA Analyzer
- NBP Hose 3 meter
- Power Cord Us (10 ft)

Nightingale™ Wirelss PPM3 Monitor

w/ ECG, SpO₂, NIBP, Masimo EtCO₂, temperature, respiratory rate, and printer for Carestream WinOMS & DSN OMS Exec

17910011807600	\$5,646.79	\$5,364.45	\$5,364.45	\$5,082.11
----------------	------------	------------	------------	------------

Included with the 17910011807600 Monitor Purchase:

- Masimo CO₂ Nomoline Adapters (2 pcs)
- Masimo ISA EtCO₂ Sidestream Analyzer
- Masimo Nomoline Adapter (2 pcs)
- Adult Nasal Cannula CO₂ Sampling (5 pcs)
- Reusable Envitec SpO₂ Sensor
- SpO₂ Extender Cable (8')
- NBP Hose (10')
- NBP Cuff, Small Adult (18-26cm)
- NBP Cuff, Adult (26-35cm)
- NBP Cuff, Large Adult (32-42cm)
- ECG 3-lead cable (10')
- Table Stand with Recorder
- 50mm Printer for PPM3 Monitor

accessories on next page ...

PATIENT DIAGNOSTIC

ZOE NIGHTINGALE™ MONITOR ACCESSORIES

ECG, SpO₂ & NIBP Accessories

4210032	ECG 5 Lead Set, Snap-On, Cables 10' Length	each	\$166.99	\$158.64	\$158.64	\$150.29
4210046	ECG 3 Lead Set, Snap-On, Cables 10' Length	each	\$125.79	\$119.50	\$119.50	\$113.21
1803007	SpO ₂ Extender Cable 8' Length, Nellcor	each	\$120.49	\$114.47	\$114.47	\$108.44
1803014	Reusable SpO ₂ Sensor, Envitec (Nellcor Compatible)	each	\$84.79	\$80.55	\$80.55	\$76.31
01030255	NBP Hose, 6' Length	each	\$62.29	\$59.18	\$59.18	\$56.06
1801014	NBP Hose, 10' Length	each	\$66.49	\$63.17	\$63.17	\$59.84
1801021	NBP Cuff, Small Adult (18-26cm")	each	\$36.79	\$34.95	\$34.95	\$33.11
1801022	NBP Cuff, Adult (26-35cm)	each	\$37.79	\$35.90	\$35.90	\$34.01
1801023	NBP Cuff, Large Adult (32-42cm)	each	\$37.79	\$35.90	\$35.90	\$34.01
1804009	Oral Temperature Probe Kit, Long	each	\$272.49	\$258.87	\$258.87	\$245.24

Mounting & Printer Solutions

01020172G	Roll Stand w/Basket only, for PPM3 Must also purchase * items sold separately	each	\$591.79	\$562.20	\$562.20	\$532.61
1808070	* 3" Mount Adapter Kit for PPM3 (attach to roll stand)	each	\$41.49	\$39.42	\$39.42	\$37.34
01020820	* Printer Mounting Kit for Roll Stand	each	\$63.79	\$60.60	\$60.60	\$57.41
1800106	2" Printer Thermal Paper (BPA Free), Each Roll	each	\$7.29	\$6.93	\$6.93	\$6.56

Miscellaneous Accessories

4210052	PPM2/PPM3 to DB9 Serial Cable, 9' (3m)	each	\$47.29	\$44.93	\$44.93	\$42.56
1807600	Wireless Kit, for PPM3 Monitor	each	\$525.79	\$499.50	\$499.50	\$473.21



01020172G
Roll Stand with basket only.
Other products shown
available separately.

Power Supply Holster
shown is not included on
stand, ask for details.



Masimo ISA™ EtCO₂ Accessories

01020488	Masimo CO ₂ Nomoline Adapter	each	\$38.79	\$36.85	\$36.85	\$34.91
1808053	Mounting Bracket Kit for Masimo ISA Analyzer	each	\$25.79	\$24.50	\$24.50	\$23.21

PATIENT DIAGNOSTIC

CONMED® ELECTRODES

- Closed cell foam keeps fluids from penetrating
- Pull-tab makes it easy to apply and remove
- Latex-free

ConMed offers several different foam electrodes that feature their proprietary gel formulation, GEL 2000X. This conductive adhesive gel formulation is designed to minimize impedance compared to earlier ConMed designs, maintain signal quality and adhere well to skin while avoiding skin irritation.

**Foam Electrodes with Conductive Adhesive Gel
Omnitrace® General Purpose Electrodes**

For an economical solution for facility wide utilization, the rectangular shape and size of the Omnitrace electrodes provides optimum adhesion and functionality for use up to three days. Unlike all other adult ECG electrodes that have a shelf life limitation of 30 days once the sealed pouch is opened, the Omnitrace electrodes offer a shelf life of 45 days.



3600003	3/pack, 30/box	\$8.99	\$8.54	\$8.54	\$8.09
3600005	5/pack, 50/box	\$15.49	\$14.72	\$14.72	\$13.94
3600050	50/pack	\$14.49	\$13.77	\$13.77	\$13.04

Positrac® Suretrace® RTL General Purpose Electrodes

The small, tear-drop shape of the Suretrace RTL electrodes provides flexibility for electrode placement and adequate adhesion for general purpose monitoring applications.



1800C003	3/pack, 30/box	\$10.99	\$10.44	\$10.44	\$9.89
1800C005	5/pack, 50/box	\$12.29	\$11.68	\$11.68	\$11.06
1800C050	50/pack	\$11.99	\$11.39	\$11.39	\$10.79

Positrac® RTL Diaphoretic Electrodes

Designed with a larger tear-drop shape and a more robust foam adhesive, the Positrac electrodes provide a more aggressive adhesive for more challenging or diaphoretic applications.



1870C003	3/pack, 30/box	\$10.99	\$10.44	\$10.44	\$9.89
----------	----------------	---------	---------	---------	--------

**Tape Electrodes with Conductive Adhesive Gel
Cleartrace™ and Cleartrace™ RTL
General Purpose Electrodes**

The combination of the low profile perforated tape, GEL 2000X, and the round shape and size, make the Cleartrace RTL electrodes an ideal solution for general purpose monitoring.



Cleartrace™					
1700005	5/pack, 50/box	\$14.49	\$13.77	\$13.77	\$13.04

Cleartrace™ RTL					
1700C003	3/pack, 30/box	\$12.79	\$12.15	\$12.15	\$11.51
1700C050	50/pack	\$17.79	\$16.90	\$16.90	\$16.01

3M RED DOT™ ELECTRODES

**Foam ECG Monitoring Electrodes
Ventilated Tape Backing With Hypoallergenic Adhesive**

- Adhesive foam-tape backing and proprietary sticky gel designed to ensure excellent skin contact.
- Clinically proven to stay adhered to skin that becomes diaphoretic
- Has built-in abradant



2570	Adult	50/bag	\$16.49	\$15.67	\$15.67	\$14.84
------	-------	--------	---------	---------	---------	---------

**Soft Cloth ECG Monitoring Electrodes
Ventilated Tape Backing With Hypoallergenic Adhesive**

- Soft cloth backing is conformable for increased patient comfort
- Patented solid gel conforms to the skin quickly allowing high quality, uninterrupted traces
- Can be used on adults with fragile skin

3M2238	Adult	50/bag	\$11.99	\$11.39	\$11.39	\$10.79
--------	-------	--------	---------	---------	---------	---------

**MEDLINE MEDGEL™
GENERAL MONITORING FOAM
ELECTRODE WITH SOLID GEL**



- Clear tracing every time with high performance, specially formulated gel with extended surface area
- MedGel™ electrodes stay securely on the skin, yet remove cleanly with top-quality adhesive
- Flexible lead wire attachment with perforated strip backing
- Easy removal with full-edge lift tab
- Flexible, fluid-resistant foam material

MDSM611550H		50/pack	\$11.29	\$10.73	\$10.73	\$10.16
-------------	--	---------	---------	---------	---------	---------

**MEDLINE MEDGEL™
WET GEL STRESS DIAPHORETIC
FOAM ELECTRODE**

- Made of comfortable, thin, fluid-impervious, medical-grade foam
- Flexible lead wire attachment with perforated strip backing
- Easy removal with full-edge lift tab
- Instant, high quality trace pickup with liquid gel



MDSM618305		100/box	\$16.79	\$15.95	\$15.95	\$15.11
------------	--	---------	---------	---------	---------	---------

**VERMONT MEDICAL
ELECTRODES**



A100063	Adult, Foam, 2-1/4"	3/pack, 75/box	\$39.49	\$37.52	\$37.52	\$35.54
A10001	Pediatric, Tape, 1-5/8	3/pack, 75/box	\$40.49	\$38.47	\$38.47	\$36.44

PATIENT DIAGNOSTIC

ACE MONITORING ELECTRODES

- 3 electrodes per bag • Pregelled
- Disposable • Economical • Hypoallergenic




Silver/silver chloride monitoring foam adhesive with low-chloride gel

- 3 electrodes per bag - 10 bag per box

ME403	Foam, Round, 2"	30/pack	\$9.29	\$8.83	\$8.83	\$8.36
-------	-----------------	---------	--------	--------	--------	--------

**Similar to 3M Red Dot Electrodes
Silver/silver chloride sensing system**

- 3 electrodes per strip - 10 strips per pouch



ME330	Round Cloth, 1-7/8"	30/pack	\$7.79	\$7.40	\$7.40	\$7.01
-------	---------------------	---------	--------	--------	--------	--------

ME430	Round Foam, 2"	30/pack	\$6.79	\$6.45	\$6.45	\$6.11
-------	----------------	---------	--------	--------	--------	--------

**ACE REUSEABLE
WRIST LIMB LEADS**

Velcro™ Limb Leads

- 2 Wrist, 1 Upper Arm Bicep
- 3 Lead Silver/Silver Eyelet Snaps



0032327		3/pack	\$72.79	\$69.15	\$69.15	\$65.51
---------	--	--------	---------	---------	---------	---------

ACE Wrist Limb Leads

- 3 electrode leads, silver/silver eyelet
- Snap on flush mount



LIMBKITFS		pair	\$93.99	\$89.29	\$89.29	\$84.59
-----------	--	------	---------	---------	---------	---------

**PAPER FOR PRINTERS
AND MONITORS**

**IBS Corporation
For Automated Blood Pressure**



7026	4" wide roll	each	\$8.99	\$8.54	\$8.54	\$8.09
------	--------------	------	--------	--------	--------	--------

Precision Chart Standard Chart Recorder Paper, ECG Single Channel

- Nonglare • Permanent Grid

10050GBC	40mm wide, 100'	3/box	\$14.49	\$13.77	\$13.77	\$13.04
----------	-----------------	-------	---------	---------	---------	---------

Precision Chart Chart Recorder Paper, for use with:

- Cardiotron PC-1, PC-2, PC-3
- Lumiscribe RS-106, RS-100D, RS-100DLT
- Miniscribe American Optical

10049	100' roll	2/box	\$9.79	\$9.30	\$9.30	\$8.81
-------	-----------	-------	--------	--------	--------	--------

PARKER SIGNASPRAY® ELECTRODE SOLUTION AND SKIN PREP

- Ideal electrolyte for use with muscle stimulation equipment
- Superior non-gritty skin prep
- Bacteriostatic
- Hypoallergenic
- Non-staining
- No residue build-up



1825	8.5 fl oz spray	each	\$8.99	\$8.54	\$8.54	\$8.09
------	-----------------	------	--------	--------	--------	--------

**PARKER SIGNAGEL®
ELECTRODE GEL**

- Highly conductive
- Hypoallergenic, bacteriostatic
- Water soluble, non-staining, non-gritty



PK1525	8.5 oz tube	each	\$5.49	\$5.22	\$5.22	\$4.94
--------	-------------	------	--------	--------	--------	--------

PARKER REDUX® GEL

- Soluble formula for high-volume routine ECG recording
- Cleans with tap water
- Recommended for use on transistorized ECG machines



6504	4 oz tube	each	\$8.99	\$8.54	\$8.54	\$8.09
------	-----------	------	--------	--------	--------	--------

**PARKER SPECTRA 360
ELECTRODE GEL**

1208	8 oz tube	each	\$2.49	\$2.37	\$2.37	\$2.24
------	-----------	------	--------	--------	--------	--------



TO YOU
IN 2



Did you know 99% of ACE SOUTHERN customers receive their order within 1-2 business days of shipment? ACE SOUTHERN has four order fulfillment facilities throughout the United States to ensure your order is TO YOU IN 2*.

***See page 261 for complete details or visit acesouthern.com/shipping-information**



PATIENT DIAGNOSTIC



LS1R9RRESP



10937002

NONIN® RESPSENSE™ AND LIFESENSE® CAPNOGRAPHS

The WIDESCREEN™ touch-panel display monitor offers MedAir™ sidestream EtCO₂ technology for a wide range of clinical settings and for both intubated and non-intubated patients.

- WIDESCREEN touch panel display
- Numerical EtCO₂ and respiratory rate display
- MedAir EtCO₂ technology
- 1.5 hours trending of EtCO₂ and respiration rate for RespSense
- 4 hours trending of EtCO₂ and respiration rate for LifeSense
- Backlit LCD display
- Audible and visual alarms
- 3-year warranty, 1-year warranty on finger sensor

RespSense™ Capnograph Includes:

- 1 Moisture Trap
- 3 Adult CO₂ Nasal Cannulas
- 1 Straight T Connector
- Operators Manual on CD-ROM
- 3 Filters
- 1 CO₂ Sample Line
- Power Supply

LS1R9RRESP	each	\$2,200.29	\$2,090.28	\$2,090.28	\$1,980.26
------------	------	------------	------------	------------	------------

LifeSense® Capnograph with Pulse Oximetry Includes:

- 1 Moisture Trap
- 3 Adult CO₂ Nasal Cannulas
- 1 Straight T Connector
- 1 8000AA2 Adult Finger Clip Sensor (2 meter)
- Operators Manual on CD-ROM
- 3 Filters
- 1 CO₂ Sample Line
- Power Supply

LS19RLIFE	each	\$2,725.29	\$2,589.03	\$2,589.03	\$2,452.76
-----------	------	------------	------------	------------	------------

NONIN® RESPSENSE™ II AND LIFESENSE® II CAPNOGRAPHS

WIDESCREEN™ capnographs offer proven PureSAT® SpO₂ and Nonin Medical's sidestream EtCO₂ technologies for a wide range of clinical settings for both intubated and non-intubated patients. Easy to use, accurate and cost effective.

- Enhanced configuration options
- Easy-to-use, large (17cm), color touchscreen display
- 36-hour internal memory with USB port
- Adjustable alarm volume
- Configurable selections for language, units of measurement and alarms
- 1.5 hours trending EtCO₂ and respiration rate for RespSense™
- Industry leading 3-year warranty, 1-year warranty on finger sensor

RespSense™ II Capnograph Includes:

- 1 Moisture Trap
- 3 Adult CO₂ Nasal Cannulas
- 1 Straight T Connector
- Operators Manual on CD-ROM
- 3 Filters
- 1 CO₂ Sample Line
- Power Supply

10938001	each	\$3,558.29	\$3,380.38	\$3,380.38	\$3,380.38
----------	------	------------	------------	------------	------------

LifeSense® II Capnograph with Pulse Oximeter

- Nonin PureSAT® SpO₂ technology
- Four-hour trending of pulse rate, SpO₂, respiration and EtCO₂

Includes:

- 1 Moisture Trap
- 3 Adult CO₂ Nasal Cannulas
- 1 Straight T Connector
- 1 8000AA2 Adult Finger Clip Sensor (2 meter)
- Operators Manual on CD-ROM
- 3 Filters
- 1 CO₂ Sample Line
- Power Supply

10937002	each	\$3,599.99	\$3,419.99	\$3,419.99	\$3,239.99
----------	------	------------	------------	------------	------------

NONIN® CAPNOGRAPH ACCESSORIES

NON5787000	TrendSense™ Memory Module Download Cable and Software*	each	\$313.99	\$298.29	\$298.29	\$282.59
760111	Nonin SpO ₂ Adult Finger Clip Sensor, 9 ft. Cable	each	\$218.29	\$207.38	\$207.38	\$196.46
3211000	nVision Data Software. (Not for use with RespSense™ II and LifeSense® II, they utilize a USB for data transfer)	each	\$364.29	\$346.08	\$346.08	\$327.86
TF110	Single-use Disposable Moisture Trap and Filter	10/box	\$81.29	\$77.23	\$77.23	\$73.16
8517001	Moisture Trap Replacement Filters	25/box	\$42.29	\$40.18	\$40.18	\$38.06
7775001	Mounting Bracket (also need 9588001)	each	\$22.29	\$21.18	\$21.18	\$20.06
9588001	Mounting Clamp (also need 7775001)	each	\$184.79	\$175.55	\$175.55	\$166.31
SENSPSNA8860	Power Supply	each	\$104.29	\$99.08	\$99.08	\$93.86

PATIENT DIAGNOSTIC

WELCH ALLYN® PROPAQ CS, VITAL SIGNS AND ATLAS PATIENT MONITORING ACCESSORIES

Welch Allyn® Blood Pressure Accessories

Small child 12-16 cm, durable 1-pc cuff, w/connectors					
REUSE-08-1SC	each	\$29.49	\$28.02	\$28.02	\$26.54

Child 15-21 cm, durable 1-pc cuff, w/connectors					
REUSE-09-1SC	each	\$30.49	\$28.97	\$28.97	\$27.44

Small adult 20-26 cm, durable 1-pc cuff, w/connectors					
REUSE-10-1SC	each	\$32.29	\$30.68	\$30.68	\$29.06

Adult 24-34 cm, durable 1-pc cuff, w/connectors					
REUSE-11-1SC	each	\$34.99	\$33.24	\$33.24	\$31.49

Large adult 32-43 cm, durable 1-pc cuff, w/connectors					
REUSE-12-1SC	each	\$34.49	\$32.77	\$32.77	\$31.04

Straight pressure hose, 8 ft.					
520012	each	\$35.99	\$34.19	\$34.19	\$32.39

Welch Allyn® Miscellaneous Accessories

ECG, 3 lead cable, 10' cable w/3 attached lead wires, snap					
008031500		\$242.79	\$230.65	\$230.65	\$218.51

Welch Allyn® Atlas Monitor Replacement Parts

Printer paper for Atlas monitor					
620040	5 rolls /box	\$24.99	\$23.74	\$23.74	\$22.49



- This Product Monitors:
 - ✓ BPM - Beats/Minute
 - ✓ Pulse Oximetry

NONIN® AVANT® 9600 TABLETOP/PORTABLE DIGITAL PULSE OXIMETER

- The stylish Avant 9600 is well-suited for a variety of applications due to its simple functionality, quick recharge and extraordinary 115 hours of memory.
- Easy-to-use - Bright LED displays, large display screen allows clear readability from across the room
- 3-year warranty

9600	each	\$1,402.49	\$1,332.37	\$1,332.37	\$1,262.24
------	------	------------	------------	------------	------------

NONIN® BLOOD PRESSURE ACCESSORIES

6191001	BP cuff, small 18-26 cm	each	\$24.99	\$23.74	\$23.74	\$22.49
6191002	BP cuff, standard, 26-35 cm	each	\$24.99	\$23.74	\$23.74	\$22.49
6191003	BP cuff, large, 32-42 cm	each	\$24.99	\$23.74	\$23.74	\$22.49

Nonin SpO₂ Adult Finger Clip Sensor


620040	9 ft. cable	each	\$24.99	\$23.74	\$23.74	\$22.49
--------	-------------	------	---------	---------	---------	---------

11

LOOKING FOR ADDITIONAL ACCESSORIES FOR YOUR CAPNOGRAPHY MONITOR?

See the Anesthesia Section for:
Sampling Lines
Dual Nasal Hoods
CO₂ Sampling Cannulas!

For the most comprehensive listing, please visit acesouthern.com today



Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price

PATIENT DIAGNOSTIC

ADC ADVIEW® 2 MODULAR DIAGNOSTIC STATION

All the features you've loved, now even better.

- Clinical-grade Advantage™ BP technology from SunTech Medical
- Contemporary, compact design with integrated handle and easy-to-read LCD display
- Manual blood pressure mode for real-time verification
- Customizable average mode options
- Next-generation EMR ready with HL7 messaging for electronic data transfer
- Network enabled with built-in Ethernet or optional Wi-Fi
- Lithium-ion power for longer-lasting battery life and a lighter weight
- Built-in memory with up to 99 charting sessions
- New Pediatric Mode to ensure accurate measurements for younger patients
- For spot check, not continuous monitoring

The ADView 2 modular diagnostic station grows with your needs. Start with blood pressure and choose the temperature or pulse oximetry options at the time of purchase. Or add new modules and accessories in the field, anytime. Our user-expandable, connect-and-go construction makes upgrading effortless.

- 3 year warranty on Main Unit, SpO2 Modules, Temp Module
- 2 year warranty on Blood Pressure Cuffs
- 2 year warranty on Desktop Caddy, Wall Mount, Mobile Stand
- 90 day warranty on accessories (SpO2 Sensors, Temp Probes, etc.)



PURCHASE BASE UNIT 1
and upgrade later, or
CUSTOMIZE and SAVE
by adding modules
2 or 3 or BOTH today!

ADC ADVIEW® 2 ACCESSORIES

ADView® 2 Mounting Platforms

Desktop Caddy: optional Caddy keeps cuffs and cords organized

9005D	each	\$46.49	\$44.17	\$44.17	\$41.84
9005W	each	\$103.29	\$98.13	\$98.13	\$92.96
9006M	each	\$324.99	\$308.74	\$308.74	\$292.49

ADView® 2 Cuff Options

All cuffs equipped with male bayonet connector.

90059CGR1MB	each	\$13.49	\$12.82	\$12.82	\$12.14
900510SARB1MB	each	\$15.29	\$14.53	\$14.53	\$13.76
900511AN1MB	each	\$15.49	\$14.72	\$14.72	\$13.94
900511APN1MB	each	\$17.79	\$16.90	\$16.90	\$16.01
900512XBD1MB	each	\$19.79	\$18.80	\$18.80	\$17.81
900513TBR1MB	each	\$22.49	\$21.37	\$21.37	\$20.24
900512XPBD1MB	each	\$22.49	\$21.37	\$21.37	\$20.24

ADView® 2 SpO2 ChipOx Finger Sensor

9005SP	each	\$97.79	\$92.90	\$92.90	\$88.01
9005SPDYS	20/box	\$334.99	\$318.24	\$318.24	\$301.49

ADView® 2 Accessories

9000TP	25 boxes/case	\$39.99	\$37.99	\$37.99	\$35.99
9005BAT	each	\$88.29	\$83.88	\$83.88	\$79.46
9005PC	each	\$18.79	\$17.85	\$17.85	\$16.91

1 ADC ADVIEW® 2 Modular Diagnostic Station Blood Pressure Base Unit

The ADView 2 easily captures automated blood pressure measurements at the touch of a button. The non-invasive blood pressure module uses the oscillometric method, and displays in about 30 seconds. With seven available cuff sizes, ranging from child to thigh, the ADView 2 can be used for patients ages three and older. Adult Plus cuffs come standard.

Includes:

- Rechargeable lithium-ion battery
- U.S. power cord
- 10-foot length of tubing
- Adult Plus and Large Adult Plus cuffs
- Mini USB to USB Cable
- Quick Reference Guide

9005BP	Base unit only	each	\$774.29	\$735.58	\$735.58	\$696.86
--------	----------------	------	----------	----------	----------	----------

1 ADC ADVIEW® 2 Blood Pressure Base Unit with SpO2 Pulse Oximetry Module

Measure oxygen saturation of arterial blood flow in about 10 seconds. The ChipOx Nellcor™-compatible SpO2 module displays pulse-rate signal strength during measurement. Comes with a reusable finger clip sensor and 8-foot cable.

Includes:

- Everything from the Base Unit
- the SpO2 ChipOx Module
- Adult finger sensor and 8-foot extension cable

9005BPS	Base and SpO2	each	\$1,394.79	\$1,325.05	\$1,325.05	\$1,255.31
---------	---------------	------	------------	------------	------------	------------

1 ADC ADVIEW® 2 Blood Pressure Base Unit with Temperature Module

Add a Covidien Filac™ 3000 Thermometry module to your ADView 2 for oral or axillary temperature readings in as little as 3-6 seconds. Features Fahrenheit and Celsius scales, plus a cold mode for lower ranges. Includes a built-in, storage well for the probe probe covers. An optional color coded rectal chamber/probe is available.

Includes:

- Everything from the Base Unit
- the Temperature Module
- Oral/axillary probe and well (blue)
- box of 20 probe covers

9005BPTO	Base and Temp	each	\$1,455.79	\$1,383.00	\$1,383.00	\$1,310.21
----------	---------------	------	------------	------------	------------	------------

1 Blood Pressure Base Unit with Temperature Module and SpO2 Module

Includes:

- Everything from the Base Unit
- The SpO2 Module and accessories that come with 9005BPS
- The Temperature Module and accessories that come with 9005BPTO

9005BPSTO	Base, temp, SpO2	each	\$2,510.29	\$2,384.78	\$2,384.78	\$2,259.26
-----------	------------------	------	------------	------------	------------	------------

PATIENT DIAGNOSTIC

NONIN ONYX® VANTAGE 9590 FINGER PULSE OXIMETER

The only finger pulse oximeter with scientifically proven accuracy and is tested for use on fingers, thumbs and toes.

- Proven Accuracy - Scientifically proven accuracy in the widest range of patients.
- Reliable Results- Accurate readings on your challenging cases due to low perfusion and darker skin tones. In addition, Onyx Vantage 9590 has been tested in motion.
- Versatile - One unit works on pediatric to larger adult patients and can be used on fingers, thumbs or toes.
- Clinician-Facing, Easy-to-Read, Bright LED - SpO2 and pulse rate can be read from any angle, in all lighting conditions, day or night.
- Three-Color Pulse Quality Indicator LED - Green, yellow or red light provides a quick assessment of patient pulse quality.
- 4-year warranty and up to 6,000 spot checks on 2 AAA batteries



9590	each	\$221.29	\$210.23	\$210.23	\$199.16
------	------	----------	----------	----------	----------

SMITHS MEDICAL BCI® 3301 HANDHELD PULSE OXIMETER

- Monitor
- 3044 adult finger sensor
- Three "C" disposable batteries
- Operation manual

Key Features:

- Display SpO2, pulse rate, pulse strength
- Rugged, ergonomic design
- Bright LED display
- Two-button operation
- Storage data for up to 99 patients
- 2-year warranty



3301A1	each	\$554.99	\$527.24	\$527.24	\$499.49
--------	------	----------	----------	----------	----------

Protective boot (yellow) is not included, can be purchased separately

EDAN H100B PULSE OXIMETER

The H100B Pulse Oximeter provides a simple universal tool for patient monitoring, from spot checks to continuous monitoring. All monitoring data can be transferred to PC for storage, management, review and printing.

- SpO2 and pulse rate measurement
- LCD display shows numeric and waveform of SpO2 simultaneously
- Powerful data storage capacity, up to 100 IDs and 300 hours data storage capability for continuous monitoring
- The latest 10 minutes trend graph and table of SpO2 and pulse rate can be reviewed in the screen
- Audible and visible alarm capability
- Long battery life – up to 48 hours on 4 AA, or 36 hours on Ni-MH rechargeable batteries (optional upgrade)
- Data can be transferred to PC for storage, review and printing (optional upgrade)
- Unit comes with carrying case and standard reusable DB9 adult sensor.
- Disposable and reusable pediatric sensors or neonate wraps also available
- 2-year warranty



H100B	Pulse Oximeter with 3' adult finger sensor	each	\$429.29	\$407.83	\$407.83	\$386.36
SH1D	3' adult finger sensor	each	\$92.79	\$88.15	\$88.15	\$83.51
2121064164	Replacement battery	each	\$19.99	\$18.99	\$18.99	\$17.99
0316113589	Battery charger stand	each	\$198.79	\$188.85	\$188.85	\$178.91

ADC DIAGNOSTIX™ FINGERTIP PULSE OXIMETER 2100

- Reads SpO2 and Pulse Rate
- Accuracy: SpO2 ± 2%/HR ± 5 BPM
- 6 changeable display modes
- 10 display brightness settings
- Fits fingers with a thickness of .3" to .8"
- Dual color OLED display
- Complete with carry lanyard and black rubberized bumper
- 2 AAA batteries
- 2-year limited warranty



ADC2100	each	\$129.79	\$123.30	\$123.30	\$116.81
---------	------	----------	----------	----------	----------

NONIN 2500A PALMSAT® DIGITAL HANDHELD PULSE OXIMETER

- Simple - Easy two-button operation
- Alarm Functions - Audible, visible, pulse rate high/low, SpO2 high/low and low perfusion
- Compact Size - 7.5 oz unit
- Memory for 72 hours of data storage
- 60 hours of battery life; 40 hours - rechargeable batteries
- 3-year warranty



Pulse Oximeter Includes:

8000AA - Adult Reusable Finger Clip Sensor, 4 AA alkaline batteries and a CD rom w/instructions

2500A	2500A PalmSAT®	each	\$786.79	\$747.45	\$747.45	\$708.11
2500B	Battery Pack	each	\$60.29	\$57.28	\$57.28	\$54.26
2500B	Universal Charger Stand w/NIMH Battery Pack	each	\$60.29	\$57.28	\$57.28	\$54.26

- This Product Monitors:
- ✓ BPM - Beats/Minute
- ✓ Pulse Oximetry

SMITHS MEDICAL BCI SPECTRO2™ 10 PULSE OXIMETER

Versatility with Multiple Sensors

- The SPECTRO2™ Pulse Oximetry Portfolio draws on proprietary BCI® technology to deliver accurate, trustworthy readings across a full range of patient care settings – even during motion or in low perfusion conditions.
- One of the largest displays on the hand-held market, 50% larger LEDs than comparable hand-held devices
- Simple, one-handed operation
- Intuitive, ergonomic design
- 3-year warranty



WW1000EN	Spectro2™ 10 Pulse Oximeter	each	\$385.29	\$366.03	\$366.03	\$346.76
----------	-----------------------------	------	----------	----------	----------	----------

Spectro2™ Rechargeable Battery Pack

WW1090	fits 10, 20 and 30	each	\$98.49	\$93.57	\$93.57	\$88.64
--------	--------------------	------	---------	---------	---------	---------

PATIENT DIAGNOSTIC

WELCH ALLYN SPHYGMOMANOMETERS

Platinum Series DS58 Hand Aneroids

- Jeweled movement for long life
- Laser-engraved dial for accuracy
- Premium inflation system for comfort
- Premium aesthetics
- Super shock-resistant—can withstand up to 60" drop and remain accurate
- Lifetime calibration warranty



Welch Allyn Tycos DS58 Classic

Hand gauge, thumbscrew release valve

The Welch Allyn Tycos DS58 Classic Premium Aneroid features gear-free DuraShock™ technology. You can count on the DS58 aneroid to remain in calibration longer than traditional aneroids, providing more accurate readings over time and a lower lifetime cost, helping to improve practice efficiencies and your bottom line. The premium inflation system comes with an ergonomic inflation bulb designed to enhance comfort while inflating a cuff.

- Aneroid w/Adult durable one-piece cuff and zipper case

WELDS5811	each	\$255.29	\$242.53	\$242.53	\$229.76
-----------	------	----------	----------	----------	----------

ADC PROSPHYG™ 760 SERIES



- Contemporary black enamel 300mmHg no-pin stop manometer provides day-in, day-out dependability
- Adcuff™ nylon cuff with ADC's proprietary Size Guide™ marking system prevents mis-cuffing
- Dip molded latex inflation bladder withstands heavy use
- Chrome-plated brass air release valve
- Leatherette carrying case with durable nylon zipper
- Assembled, inspected and packaged in the USA from foreign components
- Latex-free

Adult Blood Pressure Cuff

76011A	specify color *	each	\$30.29	\$28.78	\$28.78	\$27.26
--------	-----------------	------	---------	---------	---------	---------

Large Adult Blood Pressure Cuff

76012X	specify color *	each	\$37.79	\$35.90	\$35.90	\$34.01
--------	-----------------	------	---------	---------	---------	---------

Thigh Blood Pressure Cuff

76013T	specify color **	each	\$41.79	\$39.70	\$39.70	\$37.61
--------	------------------	------	---------	---------	---------	---------

* Color selections: black, burgundy, navy, orange

** Color selections: black, brown, navy, orange

11

WELCH ALLYN ECONOMY POCKET ANEROID

- Accuracy at an economy price
- Distinctive shape
- Indexed adult cuff
- One piece latex-free Velcro Cuff
- One tube
- 5-year warranty



DS4509	Child	each	\$95.49	\$90.72	\$85.94	\$85.94
DS4511	Adult	each	\$88.99	\$84.54	\$80.09	\$80.09
DS4512	Large adult	each	\$99.29	\$94.33	\$89.36	\$89.36

ADC MULTIKUF™

General Practice Multicuff Kit

- Manometer exceeds ANSI/AAMI/ISO 81060-1 specifications for accuracy
- Backed by a lifetime calibration warranty
- Ambidextrous, trigger-release air valve
- Quick-release luer connectors for rapid selection of appropriate cuff
- Four-section, smoke-gray plexiglass tray features a small 14" x 3-1/2" footprint
- Tray may be wall mounted with included hardware



General Practice Multicuff Kit Includes:

Small adult, adult and large adult calibrated nylon Adcuffs™

705GPKBK	kit	\$153.29	\$145.63	\$145.63	\$137.96
----------	-----	----------	----------	----------	----------

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

See page 2 for Details

DYNAREX® SPHYGMOMANOMETERS

- Accurate blood pressure readings with state of the art aneroid sphygmomanometer
- Chrome-plated, precision crafted, easy to read 300mmHG with air release valve



7106	Child, 5-9 years old	each	\$15.99	\$15.19	\$15.19	\$14.39
7107	Adult	each	\$15.99	\$15.19	\$15.19	\$14.39
7108	Large adult	each	\$19.99	\$18.99	\$18.99	\$17.99
7109	Thigh, x-large	each	\$18.99	\$18.04	\$18.04	\$17.09

PLEASE NOTE: THIS PAGE IS FOR "MANUAL" BLOOD PRESSURE MONITORING REPLACEMENT PARTS ONLY.

WELCH ALLYN BAUM LATEX BULB

06611890	w/Air Flo® control	each	\$15.49	\$14.72	\$14.72	\$13.94
06611891	Latex bulb	each	\$6.79	\$6.45	\$6.45	\$6.11

WELCH ALLYN BAUM ACCESSORIES

06612910	Extendex® coiled tubing, 8', with connectors	each	\$15.49	\$14.72	\$14.72	\$13.94
06612920	Tubing connectors, 1 male, 1 female	box	\$3.79	\$3.60	\$3.60	\$3.41
06612921	Male tubing connector	each	\$3.79	\$3.60	\$3.60	\$3.41
06612922	Female tubing connector	each	\$3.79	\$3.60	\$3.60	\$3.41
2420	Inflation system wall basket, 6-½"W x 7"H x 3-¼"D	each	\$29.49	\$28.02	\$28.02	\$26.54



WELCH ALLYN BAUM CUFFS

06611870	Adult	each	\$18.79	\$17.85	\$17.85	\$16.91
06611868	Large arm	each	\$24.99	\$23.74	\$23.74	\$22.49



WELCH ALLYN BAUM VELCRO™ CUFF & BAG COMBO

06611881	Child	each	\$24.99	\$23.74	\$23.74	\$22.49
06611880	Adult	each	\$27.49	\$26.12	\$26.12	\$24.74
06611869	Large arm	each	\$36.99	\$35.14	\$35.14	\$33.29

WELCH ALLYN BAUM INFLATION BAG

1841SS	Child, 1 tube	each	\$23.49	\$22.32	\$22.32	\$21.14
1840SS	Adult, 1 tube	each	\$24.99	\$23.74	\$23.74	\$22.49

PATIENT DIAGNOSTIC



WELCH ALLYN INFLATION SYSTEMS

Velcro™ Cuff

508221	Child, 2 tube	each	\$68.99	\$65.54	\$65.54	\$62.09
508223	Large adult, 2 tube	each	\$76.29	\$72.48	\$72.48	\$68.66

WELCH ALLYN VELCRO™ CUFFS

Velcro™ Cuff, Nonconductive

508202	Child	each	\$25.99	\$24.69	\$24.69	\$23.39
508201	Adult	each	\$25.99	\$24.69	\$24.69	\$23.39
508211	Thigh	each	\$46.29	\$43.98	\$43.98	\$41.66
WA508216	Large adult	each	\$43.49	\$41.32	\$41.32	\$39.14



WELCH ALLYN INFLATION BAGS

508904	Child, 1 tube	each	\$38.29	\$36.38	\$36.38	\$34.46
508903	Adult, 1 tube	each	\$40.49	\$38.47	\$38.47	\$36.44
508901	Adult, 2 tubes	each	\$23.49	\$22.32	\$22.32	\$21.14
508918	Thigh, 1 tube	each	\$84.79	\$80.55	\$80.55	\$76.31
508921	Large adult, 1 tube	each	\$70.29	\$66.78	\$66.78	\$63.26

WELCH ALLYN BULBS & VALVES

508701	Valve only	each	\$23.49	\$22.32	\$22.32	\$21.14
508601	Latex inflation bulb only	each	\$17.49	\$16.62	\$16.62	\$15.74



WELCH ALLYN CUFF & BAG

508224	Child, 2 tubes	each	\$65.49	\$62.22	\$62.22	\$58.94
508225	Adult, 2 tubes	each	\$43.39	\$41.22	\$41.22	\$39.05



PATIENT DIAGNOSTIC

3M™ LITTMANN® LIGHTWEIGHT II S.E. STETHOSCOPE

A general purpose stethoscope with an improved design that provides reliable acoustics for basic blood-pressure and limited physical assessment.

- Chestpiece shape designed for better orientation around blood-pressure cuffs and body contours
- Traditional bell combined with a tunable diaphragm that conveniently alternates between low- and high-frequency sounds without turning over the chestpiece
- Comfortably angled, anatomically correct headset
- Single-lumen tubing in a variety of colors
- 2-year warranty



2450	28" tubing, black *	each	\$59.29	\$56.33	\$56.33	\$53.36
------	---------------------	------	---------	---------	---------	---------

* If you prefer another color please specify from the colors below:
Burgundy, Ceil Blue, Lilac, Seafoam Green, Pearl Pink

3M™ LITTMANN® MASTER CLASSIC II STETHOSCOPE

Excellent acoustic sensitivity of a tunable diaphragm is the hallmark of the low-profile Littmann Master Classic II. Award-winning design combines high performance and portable convenience in this stethoscope with a single-sided chestpiece.

- Tunable diaphragm conveniently alternates between the low- and high-frequency sounds without the need to turn over the chestpiece
- Portable convenience of a single-sided chestpiece
- High-performance single-lumen tubing in multiple colors
- Excellent acoustic seal and comfortable fit with patented 3M Littmann Snaught Tight Soft-Sealing Eartips
- Comfortably angled, anatomically correct headset
- Patient-friendly, non-chill rim and diaphragm
- 3-year warranty



2144L	27" tubing, black *	each	\$118.99	\$113.04	\$113.04	\$107.09
-------	---------------------	------	----------	----------	----------	----------

* If you prefer another color please specify from the colors below:
Burgundy, Caribbean Blue, Ceil Blue, Hunter Green, Navy Blue, Pine Green, Purple

3M™ LITTMANN® CLASSIC III™ STETHOSCOPE

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world.

- High acoustic sensitivity for performing general physical assessments
- Tunable diaphragms on both the adult and pediatric sides of the chestpiece
- Signal-piece diaphragm that's easy to attach and easy to clean because the surface is smooth without crevices
- Next-generation tubing with improved resistance to skin oils and alcohol for longer tubing life, less likely to pick up stains



5620	27" tubing, black *	each	\$109.79	\$104.30	\$104.30	\$98.81
------	---------------------	------	----------	----------	----------	---------

* If you prefer another color please specify from the colors below:
Burgundy, Navy Blue, Plum, Chocolate, Gray, Ceil Blue, Raspberry, Caribbean Blue, Pearl Pink, Dark Olive Green, Lavender, Turquoise, Lemon-Lime, Rose Pink

3M™ LITTMANN® CARDIOLOGY IV DIAGNOSTIC STETHOSCOPES

The 3M Littmann® Cardiology IV Stethoscope is designed for clinicians who require outstanding acoustic performance combined with exceptional versatility. Its innovative design provides a single-piece tunable diaphragm on each side of the chestpiece.

- Outstanding acoustics with better audibility of high-frequency sounds (adult side, firm pressure)
- Tunable diaphragm: Hear high or low frequency sounds by slightly adjusting pressure on the chestpiece
- Dual-lumen tubing: Two sound paths in one tube eliminates the rubbing noise of traditional double tubes
- Next-generation tubing lasts longer and resists stains; Not made with natural rubber latex or phthalate plasticizers
- Headset is easily adjusted for individual fit and comfort; Angled ear tubes align with ear canals
- 7-year warranty



6151	22" tubing, black only	each	\$208.99	\$198.54	\$198.54	\$188.09
------	------------------------	------	----------	----------	----------	----------

6152	27" tubing, black *	each	\$208.99	\$198.54	\$198.54	\$188.09
------	---------------------	------	----------	----------	----------	----------

* If you prefer a different color please specify from the colors below:
Hunter Green, Navy Blue, Plum, Raspberry

3M™ LITTMANN® STETHOSCOPE TUBING COLOR SELECTIONS

Stethoscopes may not be available in all colors shown. Colors are representative.

Black	Chocolate	Raspberry	Lilac	Hunter Green	Purple	Peach	Lavender	Rose Pink
Burgundy	Gray	Caribbean Blue	Royal Blue	Plum	Seafoam Green	Pearl Pink	Turquoise	
Navy Blue	Pine Green	Orange	Ocean Blue	Ceil Blue	Red	Dark Olive Green	Lemon-Lime	

PATIENT DIAGNOSTIC

**ADC PROSCOPE™
660 AND 670
STETHOSCOPES**

ADC Proscopes combine lightweight comfort, sensitive acoustics and handsome styling, with truly affordable pricing.

- Lightweight aluminum chestpiece construction
- Proprietary diaphragm enhances acoustic response by 50% over traditional designs
- Adjustable chrome plated brass binaurals with PVC eartips
- Spare pair of black mushroom eartips included
- Inspected and packaged in the USA
- 1-year warranty

Proscope 660

- Non-chill diaphragm only nurse style chestpiece
- Ideally suited for blood pressure measurement
- Weighs 3 oz • Overall length 32"

660	specify color *	each	\$5.79	\$5.50	\$5.50	\$5.21
-----	-----------------	------	--------	--------	--------	--------

Proscope 670

- Combination diaphragm/bell chestpiece for wide frequency response
- Non-chill bell ring and diaphragm
- Weighs 3.5 oz • Overall length 32"

670	specify color *	each	\$7.29	\$6.93	\$6.93	\$6.56
-----	-----------------	------	--------	--------	--------	--------

* Color selections: Black, Burgundy, Gray, Gold, Green, Light Blue, Lavender, Magenta, Navy, Neon Orange, Neon Pink, Neon Yellow, Pink, Purple, Red, Royal Blue, Teal, Yellow



**ADC SPRAGUE
RAPPAPORT TYPE
STETHOSCOPES**

Exclusive design features an improved valve mechanism that virtually eliminates acoustic leakage and a chestpiece key to facilitate after market service.

- Latex-free double tube configuration
- Adjustable chrome plated brass binaurals
- Complete with accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips, and one pair of silicon ear tips
- Threaded chrome plated zinc chestpiece drum permits the use of 5 interchangeable chestpiece fittings - adult & pediatric diaphragm & 3 bells (1-1/4", 1", 3/4")
- Inspected and packaged in the USA

Proscope™ 640

- Traditional 16" tubing
- Weighs 7 oz • Overall length 24"
- 1-year warranty

640BK	Black only	each	\$14.29	\$13.58	\$13.58	\$12.86
-------	------------	------	---------	---------	---------	---------

Adscope™ 641

- Convenient 22" tubing length
- Weighs 8 oz • Overall length 30"
- 5-year warranty

641	specify color *	each	\$17.99	\$17.09	\$17.09	\$16.19
-----	-----------------	------	---------	---------	---------	---------

* Color selections: Black, Boysenberry, Burgundy, Caribbean Blue, Dark Green, Frosted Kiwi, Frosted Lilac, Frosted Magenta, Frosted Peacock, Frosted Plum, Frosted Seafoam, Green, Lavender, Light Blue, Magenta, Navy, Neon Blue, Neon Green, Neon Orange, Neon Pink, Neon Yellow, Pink, Purple, Red, Royal Blue, Teal



**ADC ADSCOPE® 608
CONVERTIBLE CLINICIAN
STETHOSCOPE**

Combining AFD (adjustable frequency) technology with a convertible design makes the new Adscope® 608 the most versatile of ADC's Clinician scopes.

- Convertible chestpiece provides adult diaphragm/bell or adult/pediatric diaphragm configurations
- ADC's AFD technology provides enhanced acoustic response: light pressure for low, bell-type frequency, firm pressure for high, diaphragm-type frequency
- Extra-large bell (1 3/8" diameter) for unsurpassed low-frequency response
- Ultra-sensitive diaphragm (1 3/4" diameter) for greater amplification and crisper high-frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Tubing is soft, flexible, and made without phthalates
- Includes accessory kit with two pairs of spare Adsoft Plus eartips
- Scope ID Tag included
- Weighs 5.8 oz • Overall length 30"
- Limited lifetime warranty covers all parts

608BK	Black	each	\$66.79	\$63.45	\$63.45	\$60.11
-------	-------	------	---------	---------	---------	---------



**ADC ADSCOPE® 603
CLINICIAN STETHOSCOPE**

- Combination diaphragm/bell adult chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability
- Extra large bell for unsurpassed low frequency response
- Ultra-sensitive diaphragm for greater amplification and crisper high frequency response
- Non-chill bell and diaphragm retaining rims for patient comfort
- Overall length 31"
- Lifetime warranty

603	Specify color *	each	\$50.99	\$48.44	\$48.44	\$45.89
-----	-----------------	------	---------	---------	---------	---------

* Color selections: Amethyst, Black, Blue Diamond, Burgundy, Dark Green, Gray, Indigo, Lavender, Light Blue, Metallic Caribbean, Metallic Ceil Blue, Metallic Gray, Metallic Raspberry, Metallic Orchid Haze, Navy, Orange, Pink, Purple, Red, Rose Quartz, Royal Blue, Sapphire Ice, Sea Glass, Teal, Turquoise



**ADC ADSCOPE® 602
TRADITIONAL CARDIOLOGY
STETHOSCOPE**

- Combination chestpiece precisely machined to exacting tolerance from surgical stainless steel
- Extra deep bell for unsurpassed low frequency response
- Ultra sensitive diaphragm for greater amplification and crisper high frequency transmission
- Non-chill bell and diaphragm retaining rim for patient comfort
- Large internal bore stainless steel binaurals are fixed at a 15° angle for greater comfort
- Custom molded bi-lumen 16" PVC tubing stays soft and flexible
- Double leaf internal spring construction reduces bulk and prevents kinking when folded
- Complete with one pair of PVC eartips and Adsoft deluxe PVC eartips in small and large size; Spare diaphragm and Scope ID Tag included
- Weighs 7.5 oz • Overall length 25-1/2"
- Limited lifetime warranty on metal parts

602	Specify color *	each	\$107.29	\$101.93	\$101.93	\$96.56
-----	-----------------	------	----------	----------	----------	---------

* Color selections: Black, Burgundy, Dark Green, Navy



PATIENT DIAGNOSTIC

WELCH ALLYN LIGHTWEIGHT DOUBLE-HEAD STETHOSCOPES



- Excellent for general auscultation and blood pressure
- 4 oz. lightweight chestpieces with nonchill rims
- Chrome-plated, double-head (flat diaphragm) chestpiece
- Rotatable binaurals • Flexible single-channel tubing • Double-Head, 21"
- 2-year warranty

507973	Velvet black	each	\$44.79	\$42.55	\$42.55	\$40.31
507976	Misty blue	each	\$44.79	\$42.55	\$42.55	\$40.31

AMPLIFIED PRECORDIAL STETHOSCOPE

Are you listening? High quality, strong volume auscultation during sedation with freedom of movement.

Want more freedom? Listen without wires using Bluetooth® Wireless Technology.

Sedation Stethoscope-Piezo

- Piezo chest piece with built-in microphone
- Diagnostic quality sound
- Clear heart tones
- Less ambient noise
- Warranty: 1 year on unit, 90 days on headset

Wireless with Bluetooth® and Back-up Earbuds includes: Rechargeable Amplifier with Bluetooth Integrated Technology, Bluetooth Headset, 72" Precordial Pick-Up Cable w/ Piezo Chest Piece, A/C Charger, Batteries, Manual and Amplifier Case

2110001S	Precordial stethoscope	each	\$792.49	\$752.87	\$752.87	\$713.24
----------	------------------------	------	----------	----------	----------	----------

Bluetooth® Headset Upgrade Kit (Piezo)

- Allows use of additional headset
- Includes: Piezo Wireless Amplifier, Splitter and One Plantronics M90 Headset

2100002	Headset upgrade kit	each	\$392.99	\$373.34	\$373.34	\$353.69
---------	---------------------	------	----------	----------	----------	----------

Bluetooth® Speaker Upgrade Kit includes:

JBL Bluetooth Speaker, Rechargeable Wireless Piezo Green Amplifier, Splitter and Instructions for Use

2100003	Speaker upgrade kit	each	\$529.79	\$503.30	\$503.30	\$476.81
21AAAR	Rechargeable Batteries	4/pack	\$13.99	\$13.29	\$13.29	\$12.59



SUNMED COMFORT FIT SOFTEAR™ MONASCOPE

Connects to esophageal and brachial stethoscope and Wenger chest piece

The foam is compressed, inserted into the ear canal and slowly expands to make a perfect, low pressure seal with the canal wall. This seal is critical due to the low volume of sounds transmitted. The low pressure against the skin in the ear canal makes the earpiece more comfortable which encourages more vigilant use. Tubing length is 4 feet.



6103025	each	\$3.29	\$3.13	\$3.13	\$2.96
---------	------	--------	--------	--------	--------

SUNMED WENGER CHEST PIECES

These are solid chest pieces made from stainless steel. They have relatively deep sound chambers and are supplied with luer connectors for attachment to ear pieces. They may be placed on the sternal notch where they may or may not be taped.

- All chest pieces have a luer slip male connector
- Stainless steel - no chrome plating to flake or peel



6002500	Newborn	each	\$20.79	\$19.75	\$19.75	\$18.71
6002501	Infant	each	\$17.49	\$16.62	\$16.62	\$15.74
SUN002502	Child	each	\$15.29	\$14.53	\$14.53	\$13.76
002503	Adult	each	\$16.29	\$15.48	\$15.48	\$14.66

SUNMED DOUBLE ADHESIVE DISKS

- Latex-free
- Single-use
- Disposable
- 25 strips of 4/roll



6002525	100/roll	\$22.79	\$21.65	\$21.65	\$20.51
---------	----------	---------	---------	---------	---------

3M™ DOUBLE STICK DISCS

- Pre-cut, hypoallergenic discs that are ideal for securing 3M Littmann stethoscope chestpieces
- The discs have an 11/32" aperture and a 1-1/4" diameter
- Discs are tabbed for easy handling and removal
- Excellent adhesion for secure attachment to the skin
- Hypoallergenic adhesive reduces the risk of allergic reaction



2181	102/pack	\$21.79	\$20.70	\$20.70	\$19.61
------	----------	---------	---------	---------	---------

PATIENT DIAGNOSTIC

AMSINO® NO-TOUCH INFRARED DIGITAL THERMOMETER

- One-button operation
- One second response time
- Measured distance: 0-2"
- Requires 2 AAA batteries, included
- 1-year warranty



YT1	each	\$77.29	\$73.43	\$73.43	\$69.56
-----	------	---------	---------	---------	---------

MEDLINE NO-TOUCH THERMOMETER

- Able to take human and surface temperatures in either Celcius or Fahrenheit units
- Completely eliminates patient contact
- Easily interpret measurement results using the thermometer's result indicators
- Approved for use in both clinical and home settings
- Compact and innovative design is exclusive to Medline
- Requires one CR2032 3-volt coin cell battery (included)
- 1-year warranty



MDSNOTOUCH	each	\$82.29	\$78.18	\$78.18	\$74.06
------------	------	---------	---------	---------	---------

ADC ADTEMP™ 429 NON-CONTACT THERMOMETER

Easy to use IR thermometer measures temperature instantly at the forehead. Convenient, non-contact measurement makes the Adtemp™ 429 the ideal instrument for use on all patients.

Features:

- Reads in 1 second
- Large 6mm illuminated display
- Measuring range (human) 50.0° F - 122.0° F +/- .5° F (10°C - 50°C +/- .3°C)
- 30 reading memory
- Night mode function mutes audibles for patient care
- Audible beep when measurement complete
- Auto off conserves battery strength
- Compact, ergonomic design - measures 5" (12.7cm) long
- Weighs 3 oz. (87g) including battery
- Powered by 2 AAA alkaline batteries (included)
- Optional wall bracket available
- 2-year warranty



429	each	\$64.29	\$61.08	\$61.08	\$57.86
-----	------	---------	---------	---------	---------

SHARN TEMP-A-STRIP TEMPERATURE INDICATOR



SHARN Forehead Temperature Strips have been used on more than 65 million surgical patients. Easy to use, non-invasive and inexpensive, SHARN temperature strips are a good alternative to electronic probes. Leave them on during your patient's stay in PACU for easy monitoring. These temperature strips are adjusted to display the equivalent of core temperature. Crystalline® forehead strips are ideal for use with laryngeal mask airways, and mac, or regional anesthesia when you don't have to intubate. Stays with patients from pre-op to PACU. The reliability of an electronic probe without the electronics or the probe.

5101 II	25/pack	\$67.99	\$64.59	\$64.59	\$61.19
---------	---------	---------	---------	---------	---------

WELCH ALLYN SURETEMP® PLUS 690 ELECTRONIC THERMOMETER

- Easy-to-read LCD display with flashing icons that clearly communicate in any language, intuitive user interface and ergonomic shape
- Patient safety is improved thanks to an interchangeable probe well that reduces the risk of cross-contamination
- Fast and Accurate - the SureTemp Plus 690 is one of the fastest thermometers on the market, providing 4-second oral, 10-second pediatric axillary, 15-second adult axillary and 10-second rectal measurements
- Use the SureTemp Plus 690 with all your patients, from infants to adults
- Optional wall mount has storage space for 50 additional probe covers
- Mode selection button for selecting oral, pediatric axillary, or adult axillary
- Requires 3 AA batteries (included)
- 2-year warranty on thermometer; 1-year warranty on probe



690200	SureTemp PLUS 690	each	\$326.99	\$310.64	\$310.64	\$294.29
21326	Wall mount	each	\$51.79	\$49.20	\$49.20	\$46.61
05031	Probe covers	250/sleeve	\$18.29	\$17.38	\$17.38	\$16.46

ADC ADTEMP™ 413 60-SECOND DIGITAL THERMOMETER

- Compact size - under 5" long
- Scabbard carry case keeps instrument sanitary



413B	AdTemp™ 413	each	\$5.49	\$5.22	\$5.22	\$4.94
416100	Disposable sheaths	100/box	\$3.79	\$3.60	\$3.60	\$3.41

MEDLINE ORAL DIGITAL STICK THERMOMETER

- Fast and accurate results
- With the press of a button, you can get an easy-to-read temperature in about 30 seconds
- One-piece snap case for greater convenience



MDS9950		each	\$6.29	\$5.98	\$5.98	\$5.66
MDS9607	Disposable covers	100/box	\$6.29	\$5.98	\$5.98	\$5.66

FEVERSCAN® DISPOSABLE FOREHEAD THERMOMETER



The FeverScan Disposable Thermometer with a non-sensitizing adhesive back is ideal for continual temperature monitoring.

430		100/box	\$60.29	\$57.28	\$57.28	\$54.26
-----	--	---------	---------	---------	---------	---------

11

PATIENT DIAGNOSTIC

**WELCH ALLYN®
3.5V HALOGEN
DIAGNOSTIC SET**

Welch Allyn's finest ophthalmoscope and otoscopes

- Portable
- Rechargeable battery



97150M	each	\$783.29	\$744.13	\$744.13	\$704.96
--------	------	----------	----------	----------	----------

**WELCH ALLYN®
3.5V HALOGEN
OP MACROVIEW
OTOSCOPE
DIAGNOSTIC SET**



- 28 lenses and 6 standard apertures
- 30% more magnification and twice the field of view of traditional otoscopes
- 3.5V direct plug in (Accepts all Welch Allyn 3.5V headss)
- Includes replacement NiCad Rechargeable Battery
- Hard Case for 3.5V Set

97150MC	each	\$749.99	\$712.49	\$712.49	\$674.99
---------	------	----------	----------	----------	----------

**WELCH ALLYN®
UNIVERSAL CHARGER
AND DESK SETS**

- Compact desk charger for convenient drop-in recharging and flexible mobility
- Colored LED's indicate power to the units, and that the handles are seated and charging
- Advanced circuitry for constant, automatic, safe charging, and optimum lamp performance and life
- Compact for convenience and mobility



Universal Charger Complete Desk Set

Includes: Universal charger base, two rechargeable NiCad handles, coaxial ophthalmoscope, and MacroView™ Otoscope with specula

71641M	each	\$1,089.29	\$1,034.83	\$1,034.83	\$1,034.83
--------	------	------------	------------	------------	------------

Universal Charger with two Handles

Includes: Universal charger base and two rechargeable NiCad handles

71630	each	\$639.49	\$607.52	\$607.52	\$607.52
-------	------	----------	----------	----------	----------

**WELCH ALLYN®
3.5V STANDARD
OPHTHALMOSCOPE**

- 28 focusing lenses with a range of -25 to +40
- Rubber brow rest prevents scratching of eyeglasses
- 6 apertures for generalist use: micro, small, and large spot sizes, fixation target, slit aperture, and red-free filter
- Halogen light for true tissue color and consistent, long-lasting illumination
- Ophthalmoscope - Head Only



11710	Head only	each	\$321.29	\$305.23	\$305.23	\$289.16
-------	-----------	------	----------	----------	----------	----------

**WELCH ALLYN®
OTOSCOPE 3.5V
WITH SPECULA**

- Halogen HPX™ lamp provides 30% more light output for true tissue color and consistent, long-lasting illumination
- Fiber optics project cool light to the distal tip with no reflections and no obstructions
- Wide-angle viewing lens allows for instrumentation under magnification
- Sealed system allows for pneumatic otoscopy



25020	Head only	each	\$206.79	\$196.45	\$196.45	\$186.11
-------	-----------	------	----------	----------	----------	----------

**WELCH ALLYN® 3.5V
DIRECT PLUG-IN
RECHARGEABLE HANDLE**

- Direct plug-in will fit all Welch Allyn 3.5V instrument heads
- Built-in AC adapter with prongs enables convenient recharging wherever your practice takes you
- Adjustable rheostat controls light intensity
- Long-lasting rechargeable nickel-cadmium battery included
- Smooth, non-abrasive knurled finish ensures a no-slip grip
- Chrome-plated brass construction provides durable and reliable long-life



71000A	each	\$222.99	\$211.84	\$211.84	\$200.69
--------	------	----------	----------	----------	----------

**WELCH ALLYN® 3.5V
RECHARGEABLE NICAD HANDLE
FOR DESK/WELL CHARGERS**

- Accepts all Welch Allyn 3.5V instrument heads
- Recharges in the Welch Allyn Universal Charger and numerous stand consoles and well chargers
- Adjustable rheostat controls light intensity
- Long-lasting rechargeable nickel-cadmium battery included
- Smooth, non-abrasive knurled finish ensures a no-slip grip
- Chrome-plated brass construction provides durable and reliable long-life



71670	each	\$188.79	\$179.35	\$179.35	\$169.91
-------	------	----------	----------	----------	----------

PATIENT DIAGNOSTIC

WELCH ALLYN® 3.5V CONVERTIBLE RECHARGEABLE HANDLE

- A versatile handle for use with both rechargeable and standard C-cell batteries
- Accepts all Welch Allyn 3.5V instrument heads
- Built-in AC adapter with prongs enables convenient recharging wherever your practice takes you
- Converts to C-cell operation quickly and easily, a stop gap if you ever forget or cannot recharge your handle
- Adjustable rheostat controls light intensity
- Long-lasting rechargeable nickel-cadmium battery included
- Chrome-plated brass construction provides durable and reliable long-life construction



71000C	each	\$250.49	\$237.97	\$237.97	\$225.44
--------	------	----------	----------	----------	----------

WELCH ALLYN® REPLACEMENT LAMPS

- For: Otoscope 20000, 25020, 21700, 20200, 20202
- 3.5V Handle Adapter 73500
- Tongue Blade Holder 28100
- Illuminator, Transilluminator 27000, 27050, 26530, 41100, 43300



03100U	each	\$24.49	\$23.27	\$23.27	\$22.04
--------	------	---------	---------	---------	---------

HEALTHLINK® REPLACEMENT BATTERIES

- 3.5V NiCad
- Orange Label, Compare to Welch Allyn® 72300

HEA1230	each	\$43.49	\$41.32	\$41.32	\$39.14
---------	------	---------	---------	---------	---------

WELCH ALLYN® REPLACEMENT BATTERIES

For Handle 70700, 71500

72000	2.5V	each	\$77.99	\$74.09	\$74.09	\$70.19
-------	------	------	---------	---------	---------	---------



For Handle 71000A, 71020A, 71000C, 71020C, 71055C

72300	3.5V	each	\$49.99	\$47.49	\$47.49	\$44.99
-------	------	------	---------	---------	---------	---------




For Handle 71670

72200	3.5V	each	\$49.99	\$47.49	\$47.49	\$44.99
-------	------	------	---------	---------	---------	---------



Actual Size



WELCH ALLYN® REUSABLE POLYPROPYLENE DIAGNOSTIC OTOSCOPE SPECULA

- Reusable specula are lightweight yet durable, and may be autoclaved, boiled, or sterilized in solution
- May be used with 23810, 23820, 20000, 25020, 22800/22820 and 24031



24400U	2.5, 3, 4, and 5mm	4/box	\$14.49	\$13.77	\$13.77	\$13.04
WEL24303U	3mm tip only	each	\$4.79	\$4.55	\$4.55	\$4.31
24304U	4mm tip only	each	\$4.79	\$4.55	\$4.55	\$4.31
WEL24305U	5mm tip only	each	\$4.79	\$4.55	\$4.55	\$4.31

WELCH ALLYN® REUSABLE OTOSCOPE SPECULA

- for Pneumatic, and Consulting Oscopes
- Reusable polypropylene — lightweight but durable
- May be autoclaved, boiled or cleaned with standard germicides



WA22100	2, 3, 4, 5 and 9mm	5/box	\$23.49	\$22.32	\$22.32	\$21.14
22003	3mm tip only	each	\$6.79	\$6.45	\$6.45	\$6.11
22004	4mm tip only	each	\$6.79	\$6.45	\$6.45	\$6.11


WELCH ALLYN® KLEENSPEC® DISPOSABLE SPECULA

- Compatible with any black or beige (diagnostic) Welch Allyn otoscope, model series 200, 250, and PocketScope™

52432	2.75mm	850/box	\$36.99	\$35.14	\$35.14	\$33.29
52434	4.25mm	850/box	\$56.49	\$53.67	\$53.67	\$50.84

BECAUSE WE JUST CAN'T FIT EVERYTHING INTO THIS CATALOG!

VISIT ACESOUTHERN.COM TODAY



Tier Pricing Details see page 2

Regular Price | Member Price | Premier Price | Premier+ Price

PATIENT DIAGNOSTIC

BREWER TASK LIGHTING

LED Eco-Series

Light details

- Standard 1-year warranty
- Eco-friendly long-life LED bulb for lower cost of ownership
- Bright, best-in-class white light drives clinical efficiency
- Cool to the touch and easy to maneuver
- Head rotation: 170° left/right
- High intensity, low voltage, 12-watt LED bulb included
- Brightness: 880 lumens; 5" beam dia. at 18" distance
- Color temperature: 4000° Kelvin
- Bulb life: Approx. 50,000 hrs

Arm details

- Smooth, easy-to-clean flexible arm and head
- Flexible, "no-drift," 23" (58.4cm), gooseneck
- Maximum height: 66" (167.6cm)

Base details

- Power cord is detachable and exits the weighted base at the floor level
- Optional five-caster base available (see below)
- Base diameter: 13" (33cm)
- Cord length: 10' (304.8cm)



BRE19100	Biege	each	\$388.99	\$369.54	\$369.54	\$350.09
----------	-------	------	----------	----------	----------	----------

Accessories

200165	Five Caster Base, beige	each	\$57.79	\$54.90	\$54.90	\$52.01
--------	-------------------------	------	---------	---------	---------	---------

Replacement Bulb

2102755	50,000 hr bulb life	each	\$75.49	\$71.72	\$71.72	\$67.94
---------	---------------------	------	---------	---------	---------	---------

WELCH ALLYN® PROFESSIONAL PENLIGHT

High-quality, halogen illumination in a portable hand-held light

- Halogen light produces true tissue color rendition
- Highly durable solid brass construction
- Operates with AAA batteries
- Welch Allyn will repair or replace, free of charge, any parts of its own manufacture proven to be defective through causes other than misuse, neglect, damage in shipment, or normal wear.



76600		each	\$58.29	\$55.38	\$55.38	\$52.46
-------	--	------	---------	---------	---------	---------

WELCH ALLYN® DIRECT FOCUSING HALOGEN LAMP

02500U	6v vacuum lamp	each	\$25.79	\$24.50	\$24.50	\$23.21
--------	----------------	------	---------	---------	---------	---------

WELCH ALLYN® 460 GREEN SERIES™ PORTABLE HEADLIGHT



Welch Allyn's 460 Green Series Portable Headlight illuminates precisely where it is needed, with cool, bright light that ensures comfort for both patient and physician.

- Cord-free, LED illumination
- Weighing only 4 oz - it's one of the most comfortable lightweight headlights available
- Rechargeable batteries housed in the headband make this a completely portable, compact unit, run time of 60 minutes
- With a color temperature of 5300 degrees Kelvin this light supplies bright (40 lumens), white, shadow-free light
- Its operating life expectancy is more than 10,000 hours

Green Series Portable Headlight

Includes: Terrycloth Headband and Direct Power Supply/Charger

46070	each	\$374.29	\$355.58	\$355.58	\$336.86
-------	------	----------	----------	----------	----------

Rechargeable Battery

72270	each	\$86.29	\$81.98	\$81.98	\$77.66
-------	------	---------	---------	---------	---------

WELCH ALLYN® 490 GREEN SERIES™ MINOR PROCEDURE HEADLIGHT



Welch Allyn's 490 Green Series Minor Procedure Headlight provides cool, bright, white light for procedures and examinations. When you need small, well-defined illumination without the heat of traditional headlights, look to the 490 Green Series Minor Headlight to provide superior quality and reliability.

- Lightweight - Weighing only 6 oz - it's one of the most comfortable lightweight headlights available
- True tissue color - With a color temperature of 5700 degrees Kelvin, this light supplies bright (120 lumens), white, shadow-free light allowing the end user to see tissue characteristics without distortion
- Long operating life - light source has an estimated life up to 50,000 hours
- Direct power or portable - Enjoy the flexibility of operating the headlight using the direct power source or the rechargeable "belt-clip" portable power source, run time of 60 minutes

Green Series Minor Procedure Headlight

Includes: Direct Power Source

49020	each	\$887.79	\$843.40	\$843.40	\$799.01
-------	------	----------	----------	----------	----------

Additional Direct Power Source

75270	each	\$167.49	\$159.12	\$159.12	\$150.74
-------	------	----------	----------	----------	----------

Battery Pack with Charger

75260	each	\$470.99	\$447.44	\$447.44	\$423.89
-------	------	----------	----------	----------	----------

Replacement Cells for Battery Charger

72250	each	\$116.49	\$110.67	\$110.67	\$104.84
-------	------	----------	----------	----------	----------

Portable Power Source

75200	each	\$364.79	\$346.55	\$346.55	\$328.31
-------	------	----------	----------	----------	----------

PATIENT DIAGNOSTIC

MEDLINE EVENCARE® G2 BLOOD GLUCOSE MONITORING SYSTEM WITH VOICE GUIDANCE



Advanced technology makes diabetic care easier.

- Easy-to-use: Voice assistance guides you through the testing process and provides results in either English or Spanish
- No coding required: Minimizes errors and inaccurate results
- Fast test time: Results in 6 seconds
- Small blood sample size: 1µL
- Alternate site testing: Option to test on the forearm, palm or fingertips
- Memory option: Stores 300 test times and can provide 7-, 14- and 30-day averaging
- Alarm option: Set reminders for up to 3 different times
- 5-year warranty



EvenCare G2 Glucose Meter

Includes: Meter, Battery, Carrying Case, Log Book, User's Guide

MPH1540	G2 Meter	each	\$33.29	\$31.63	\$31.63	\$29.96
MPH1550Z	G2 Strips	50/box	\$19.49	\$18.52	\$18.52	\$17.54
MPH1560Z	G2 Hi Lo Control	each	\$15.29	\$14.53	\$14.53	\$13.76

MEDLINE EVENCARE® G3 PROFESSIONAL, MULTI-PATIENT USE BLOOD GLUCOSE MONITORING SYSTEM

The EvenCare G3 Professional Blood Glucose Monitoring System from Medline is designed to meet the unique needs of healthcare facilities and their patients. When used as directed, it's cleared for multi-patient use.



- Quality control lock mode: Ensures quality control testing is performed
- Test strip ejector: Helps reduce cross contamination by eliminating contact with used strip
- No coding required: Minimizes errors and inaccurate results
- Fast test times: Results in 6 seconds
- Small blood sample size: 0.8µL
- Alternate site testing: Option to test on the forearm, palm or fingertips
- 2-year warranty

EvenCare G3 Glucose Meter

Includes: Meter, Battery, Carrying Case, User's Guide

MPH3540	G3 Meter	each	\$31.99	\$30.39	\$30.39	\$28.79
MPH3550Z	G3 Strips	50/box	\$24.99	\$23.74	\$23.74	\$22.49
MPH3560	G3 Hi Lo Control	each	\$18.99	\$18.04	\$18.04	\$17.09

MEDLINE SAFETY LANCETS WITH PRESSURE ACTIVATION

- Simple and safe for both the healthcare professional and the patient
- Needle retracts automatically and locks out for extra safety, protecting against needlestick injuries, and it also helps prevent cross-contamination



MPHPRESS21Z	21g	100/box	\$22.99	\$21.84	\$21.84	\$20.69
MPHPRESS26Z	26g	100/box	\$22.99	\$21.84	\$21.84	\$20.69

LIFESCAN, INC. ONETOUCH® ULTRA BLOOD GLUCOSE MONITORING



OneTouch® Ultra®2 Blood Glucose Monitoring System

A simple way to see how food effects your blood sugar.

Ultra 2 system comes with Delica lancing device,

10 One-Touch Delica 30G sterile lancets and carrying case (batteries included).

- See the effects of food choices over time. Before and after meal results (7, 14 & 30 day averages).
- Fast and easy results in just 5 seconds
- Requires only a speck of blood

5670390	OneTouch Ultra 2	each	\$115.49	\$109.72	\$109.72	\$103.94
---------	------------------	------	----------	----------	----------	----------

OneTouch Delica Lancets

5670081	33g (for 5670390)	100/box	\$16.49	\$15.67	\$15.67	\$14.84
---------	-------------------	---------	---------	---------	---------	---------

OneTouch Ultra 2 Test Strips

5676191	Ultra 2 Test Strips	25/box	\$50.99	\$50.99	\$50.99	\$50.99
5670073	Ultra 2 Test Strips	50/box	\$102.29	\$97.18	\$97.18	\$92.06
5670092	Ultra 2 Test Strips	100/box	\$199.99	\$189.99	\$189.99	\$179.99

QUIDEL QUICKVUE ONE-STEP HCG URINE TEST

Rapid Diagnosis

- Results in 3 minutes
- Positive may become visible sooner
- Simple, one-step procedure uses only 3 drops of sample
- Two-color endpoint for easy-to-read and interpret results

20109		25/box	\$68.29	\$64.88	\$64.88	\$61.46
-------	--	--------	---------	---------	---------	---------

QUIDEL QUICKVUE IN-LINE® STREP A TEST

- Results in 5 minutes or less
- One easy step via in-line extraction
- Two-color endpoint for easy to read and interpret results
- Single reagent method allows for less chance of operator error
- All components necessary to perform a test are contained in the kit, including swabs
- Built-in Positive and Negative controls confirm the test has performed properly
- Positive and Negative external controls are provided in the kit
- Individually packaged for convenient single-use testing and disposal

00343		25/box	\$175.99	\$167.19	\$167.19	\$158.39
-------	--	--------	----------	----------	----------	----------



eELECTRONIC FORM 222
ONLINE CONTROLLED SUBSTANCE
ORDERING SYSTEM

ONLINE ORDERING SYSTEM

CSOS:

The Federal DEA has created the Controlled Substance Ordering System (CSOS) to allow registrants to place orders for Class II narcotics electronically without the use of manual paperwork.

Important Contacts:

FEDERAL DEA CSOS CUSTOMER SUPPORT
877.332.3266

ACE SOUTHERN DEA COMPLIANCE TEAM
888.222.3722 (opt. 5)

**EXPRESS222 WEBSITE
SUPPORT TEAM (LEGISYM)**
254.933.4452 (opt. 2)

AFTER YOUR REGISTRATION IS COMPLETE

- No more paper 222 forms returned/voided due to errors
- Faster transaction times - no mail delays
- Reduced amount of paperwork to complete and retain on file

Four easy steps to get started

Step 1. Enroll

- Visit www.deaecom.gov to enroll.
- Follow the instructions provided. For any questions during this period or to check the progress of your paperwork, contact the Federal DEA at 877.332.3266.

Step 2. Install

- You will receive an Access Code by email and an Access Code Password by postal mail from the Federal DEA.
- Use these to retrieve and install the digital signature onto your computer. For any questions or problems during this portion, please contact the Express222 support team (Legisym) at 254.933.4452 (opt. 2).
Note: this must be done within 60 days of the date printed on your postal mailed notice.

Step 3. Contact us to gain access to the e222 ordering area

- After you have retrieved and installed the digital signature, call a member of the ACE SOUTHERN DEA Compliance Team at 888.222.3722 (opt. 5) to create an ACE SOUTHERN user profile on the Express222 ordering site.

Step 4. Test and Sign In

- Perform a Digital Certificate Test and Email Test on the ordering site: www.express222.com.
- If test was unsuccessful, contact a member of the Express222 support team at 254.933.4452 (opt. 2).
- If the test was successful, proceed to Customer area and place your order.

**FOR MORE INFORMATION VISIT:
www.deaecom.gov**

STERILIZATION

ACE CLEANING BRUSHES

- Not autoclavable
- 7" (18cm)



6039015	Nylon Bristles	3/pack	\$15.49	\$15.49	\$13.94	\$13.94
6039016	Stainless Steel Bristles	3/pack	\$16.99	\$16.99	\$15.29	\$15.29

INTEGRA™ INSTRUMENT CLEANING BRUSHES

7-¼" Nylon Bristle Cleaning Brushes

- Designed to clean instruments without scratching or harming delicate surfaces
- Autoclavable



MIL31000		3/box	\$18.99	\$18.99	\$17.09	\$17.09
----------	--	-------	---------	---------	---------	---------

7-¼" Stainless Steel Bristle Cleaning Brushes

- Designed to clean instrument serrations, rasps, files and burs
- Autoclavable



MIL31001		3/box	\$19.79	\$19.79	\$17.81	\$17.81
----------	--	-------	---------	---------	---------	---------

HU-FRIEDY® LILAC UTILITY GLOVES

Three pack of flocklined nitrile utility gloves for use when cleaning and processing instruments.



40060	small (size 7)	3 pair/bag	\$20.29	\$19.28	\$19.28	\$18.26
40062	medium (size 8)	3 pair/bag	\$20.29	\$19.28	\$19.28	\$18.26
40064	large (size 9)	3 pair/bag	\$20.29	\$19.28	\$19.28	\$18.26
40066	x-large (size 10)	3 pair/bag	\$20.29	\$19.28	\$19.28	\$18.26

CROSSTEX DENTAPURE™ INDEPENDENT WATER BOTTLE CARTRIDGE

Installs onto your unit's existing water bottle pickup tube via included quick-connect fitting.

- 1 installation provides 365 days (or 240L of water if usage records are kept) of safe, compliant dental unit water
- Can use either tap or distilled water
- Once installed, no monitoring/shocking protocol is required
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water



DP365B10	each	\$292.49	\$277.87	\$277.87	\$263.24
----------	------	----------	----------	----------	----------

CROSSTEX DENTAPURE™ MUNICIPAL WATER CARTRIDGE

Easily adapted for installation onto incoming waterlines plumbed directly to city water via included poly-flo fittings.

- 1 installation provides 365 days (or 240L of water if usage records are kept) of safe, compliant dental unit water
- Once installed, no monitoring/shocking protocol is required for the life of the cartridge
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water



DP365M10	each	\$292.49	\$277.87	\$277.87	\$263.24
----------	------	----------	----------	----------	----------

DUKAL ALCOHOL DISPENSER

- Spill-proof
- Polypropylene dispensing unit
- Clear bottle with white top



DUK4024	9 oz. bottle	each	\$12.29	\$11.68	\$11.68	\$11.06
---------	--------------	------	---------	---------	---------	---------

ALCONOX® POWDERED PRECISION CLEANER

- Concentrated to save you money (dilute 1:100, yields up to 52 gallons)
- Biodegradable and readily disposable
- Penetrating wetting power to save you time
- Replaces corrosive acids and hazardous solvents

Used To Remove: Soil, grit, grime, buffing compound, slime, grease, oils, blood, tissue, salts, deposits, particulates, solvents, chemicals, radioisotopes, radioactive contaminations, silicon oils, and mold release agents

Cleaning Method: Soak, brush, sponge, cloth, ultrasonic, flow through clean-in-place; Will foam-not for spray or machine use



21835032	4 lb carton	each	\$34.99	\$33.24	\$33.24	\$31.49
----------	-------------	------	---------	---------	---------	---------

STERILIZATION

INTEGRA™ SURGICAL INSTRUMENT CLEANER & LUBRICANT

- Chemically engineered for use as a presoak or for manual cleaning, washer/decontaminators, ultrasonic cleaners and cart washers
- A special penetrating lubricant is blended into the cleaning formula, protecting the instruments and extending their life
- Detergents and emulsifiers quickly break up and remove organic soil from box locks and other hard to reach places restoring the look and feel of a new instrument



3710	* gallon	each	\$38.29	\$36.38	\$36.38	\$34.46
------	----------	------	--------------------	--------------------	--------------------	--------------------

MEDLINE SURGICAL INSTRUMENT LOW SUDS DETERGENT

- Recommended for use in automatic washers and ultrasonic machines when instruments have been pre-cleaned Hemolytic detergent removes blood and soil, then rinses spot-free, leaving no residue
- Contains lubricating agents and rust inhibitors to help prolong instrument life



MDS88000B1	* gallon	each	\$27.29	\$25.93	\$25.93	\$24.56
------------	----------	------	--------------------	--------------------	--------------------	--------------------

INTEGRA™ SURGICAL INSTRUMENT CLEANER

- Neutral pH and phosphate-free liquid concentrate formula for manual and ultrasonic cleaning or soaking of surgical/dental instruments, accessories and glassware
- Biodegradable and rinses free
- Yields 32 gallons of solution



3720	8 oz bottle	each	\$12.49	\$11.87	\$11.87	\$11.24
------	-------------	------	--------------------	--------------------	--------------------	--------------------

MEDLINE INSTRUMENT LUBRICANT

- Protect your investment and prolong instrument life with lubricants that prevent rust and corrosion
- Non-silicone, patient-safe lubricant keeps box locks working freely and maintains effectiveness of grasping and cutting surfaces
- Protects all types of metal instruments, including carbon steel
- Will not interfere with steam or gas sterilization



MDS88000T2	* gallon	each	\$44.79	\$42.55	\$42.55	\$40.31
------------	----------	------	--------------------	--------------------	--------------------	--------------------

INTEGRA™ SPRAY LUBE

- Pre-mixed, non-silicone, anti-corrosive formula in a ready-to-use 8 oz spray pump bottle
- Steam permeable and water soluble
- Helps prevent spotting, staining and rusting
- Lubricates instrument box locks, screw locks, scissors blades, and any other metal to metal moving parts



3700	8 oz bottle	each	\$12.49	\$11.87	\$11.7	\$11.24
------	-------------	------	--------------------	--------------------	-------------------	--------------------

HU-FRIEDY® INSTRUMENT LUBRICANT SPRAY

- Does not contain silicone.
- Comes in a convenient non-aerosol pump bottle.



ILS	8 oz bottle	each	\$20.79	\$19.75	\$19.75	\$18.71
-----	-------------	------	--------------------	--------------------	--------------------	--------------------

BIEN-AIR SPRAYNET® CLEANING SPRAY

1600036006	500ml	each	\$73.99	\$70.29	\$70.29	\$66.59
------------	-------	------	--------------------	--------------------	--------------------	--------------------



BIEN-AIR LUBRIFLUID® SPRAY LUBRICANT

1600064006	500ml	each	\$73.99	\$70.29	\$70.29	\$65.59
------------	-------	------	--------------------	--------------------	--------------------	--------------------



METREX METRILUBE® INSTRUMENT LUBRICANT CONCENTRATE

When diluted, this oil-based instrument "milk" effectively coats metal instruments with an oil film that inhibits corrosion and keeps hinges working freely. MetriLube helps prevent spotting, staining, rusting and freeze up of instruments.

- Safe for automated washers
- Non-greasy
- Anti-corrosive



103425	24 oz spray	each	\$11.79	\$11.20	\$11.20	\$10.61
103400	* gallon	each	\$39.49	\$37.52	\$37.52	\$35.54

* A freight charge of \$5.50 will be added to the price of this product.

STERILIZATION

**ASEPTICO
ASEPTISPRAY**

- Pressurized handpiece cleaner and lubricant.
- Withstands autoclaving.
- A must for proper cleaning and maintenance of your handpieces.



AHP09	8.8 oz can	each	\$41.79	\$41.79	\$37.61	\$37.61
-------	------------	------	---------	---------	---------	---------

**METREX EMPOWER® FOAM
FOAMING DUAL ENZYMATIC SPRAY**

Ready-to-use foaming dual enzymatic spray that starts the precleaning process by breaking down blood, tissue, mucous and other protein-rich body fluids from instruments.

- Covers instruments, reducing splashes and spills
- Ready to use in a convenient sprayer
- Breaks down bioburden on instruments
- Speeds up the cleaning process



10-4224	24 oz spray	each	\$18.29	\$17.38	\$17.38	\$16.46
---------	-------------	------	---------	---------	---------	---------

**MEDLINE DUAL-ENZYMATIC
SURGICAL INSTRUMENT
DETERGENT AND PRESOAK**

- Recommended for automatic washers, ultrasonic machines or manual cleaning when instruments have not been precleaned
- Highly effective, dual-enzyme formula disintegrates protein deposits, starch and carbohydrates, even in hard-to-reach areas



MDS88000B9	* gallon	each	\$37.79	\$35.90	\$35.90	\$34.01
------------	----------	------	---------	---------	---------	---------

**INTEGRA™ EZ-ZYME®
ALL PURPOSE
ENZYME CLEANER**

Neutral pH multi-enzyme liquid concentrate ideal for instrument presoak, ultrasonic cleaning, laundry presoak, and evacuation system cleaner. Ideal for cleaning hard to remove proteinaceous (mucous & tissue) and lipid (fatty) materials from the surfaces of surgical/dental instruments, linen, water submersible equipment, or where mechanical scrubbing is impractical.



- ¼ oz (each packet yields 1 gallon of solution)

3750	32 packets/box		\$77.79	\$73.90	\$73.90	\$70.01
------	----------------	--	---------	---------	---------	---------

**HU-FRIEDY® ENZYMAX® LIQUID
CONCENTRATED ULTRASONIC
DETERGENT AND PRESOAK**

Enzymax Dual Enzyme Ultrasonic Detergent is specifically formulated for efficient ultrasonic cleaning or pre-soaking of cassettes and instruments. Contains Steelgard agents to protect instruments from harmful minerals naturally found in water that contribute to buildup, spotting and corrosion. Can also be used as an evacuation system cleaner.



- 1 gallon bottle makes 128 gallons of solution

IMS1226	* gallon	each	\$74.79	\$71.05	\$71.05	\$67.31
---------	----------	------	---------	---------	---------	---------

**JOHNSON & JOHNSON ENZOL®
ENZYMATIC DETERGENT**

Allows instruments to be thoroughly cleaned before they are high level disinfected or sterilized.



- Results begin in just 1 minute
- Effective at room temperature, eliminating the need for hot water
- Separate, premeasured pump allows easy dispensing
- Gently removes organic matter from a wide range of instruments
- Easy to rinse detergent does not leave film on instruments or clog scopes
- Low-foaming feature facilitates use in multiple systems including automatic endoscope reprocessors, ultrasonic cleaners, and washer decontaminators according to manufacturer's recommendations
- Mild pH reduces risk of instrument damage

2252	* gallon	each	\$83.79	\$79.60	\$79.60	\$75.41
------	----------	------	---------	---------	---------	---------

**HU-FRIEDY® ENZYMAX® EARTH
POWDER PACKETS - ULTRASONIC CLEANER**

IMS1333	96 packets	each	\$132.79	\$126.15	\$126.15	\$119.51
---------	------------	------	----------	----------	----------	----------

**METREX EMPOWER® DUAL
ENZYMATIC DETERGENT**

A dual-enzymatic detergent that works fast to remove blood, tissue, mucous and other protein-rich body fluids from instruments.


- Contains two proteolytic enzymes for effective cleaning
- Safe for use on instruments
- Less corrosive than distilled water
- Low-foaming formula is ideal for use in automated machines
- Fresh scent



104100	* gallon	each	\$40.29	\$38.28	\$38.28	\$36.26
--------	----------	------	---------	---------	---------	---------

**MEDLINE SINGLE-ENZYMATIC
SURGICAL INSTRUMENT
DETERGENT & PRESOAK**

- Recommended for automatic washers, ultrasonic machines or manual cleaning when instruments have not been precleaned
- Single-enzyme formula removes protein only



MDS88000B91	* gallon	each	\$32.49	\$30.87	\$30.87	\$29.24
-------------	----------	------	---------	---------	---------	---------

* A freight charge of \$5.50 will be added to the price of this product.

STERILIZATION

GREAT LAKES SURGICAL WAVICIDE-01® DISINFECTANT

Wavicide-01 comes pre-mixed in a gallon bottle. There is no chance of a mixing error. The remaining solution in the bottle is good up to the expiration date that is stamped on the gallon bottle. Wavicide-01 is effective at room temperature and can be reused up to 30 days.



0104	* gallon	each	\$55.49	\$52.72	\$52.72	\$49.94
------	----------	------	---------	---------	---------	---------

SULTAN HEALTHCARE SPOROX II STERILIZING AND DISINFECTING SOLUTION

High-level disinfection without irritating glutaraldehydes. Sporox II Solution is formulated with powerful hydrogen peroxide to disinfect and sterilize heat-sensitive instruments that cannot be autoclaved.



- Ready to use - no mixing, heating or activation needed
- Safe, with no noxious odors, does not need a special ventilation system installed
- No special disposal required
- Won't bond infectious proteins to instruments, unlike glutaraldehydes
- Disinfects in 30 minutes, sterilizes in 6 hours
- Re-use the solution for up to 21 days

75156	* gallon	each	\$47.99	\$45.59	\$45.59	\$43.19
-------	----------	------	---------	---------	---------	---------

JOHNSON & JOHNSON CIDEX® OPA SOLUTION

Nothing is better for expensive endoscopes and other medical devices than CIDEX OPA Solution (0.55% ortho-phthalaldehyde), the non-glutaraldehyde solution for disinfection of medical devices.

20390	* gallon	each	\$83.99	\$79.79	\$79.79	\$75.59
-------	----------	------	---------	---------	---------	---------

METREX METRICIDE® OPA PLUS SOLUTION

ortho-Phthalaldehyde (OPA) High-Level Disinfectant

- Ideal for both automated and manual instrument reprocessing
 - When used manually, it disinfects in 12 minutes at 20°C and can be reused for up to 30 days (when monitored according to label instructions) - the longest lasting high-level disinfectant on the market today
 - When used in an AER, it disinfects in 5 minutes at 25°C in an AER and can be reused for up to 14 days (when monitored according to label instructions)
 - MetriCide OPA Plus Solution Test Strips are fast - see results in only 60 seconds
- A convenient, ready-to-use, high-level disinfectant that provides a broad spectrum of kill including Tuberculosis, Hepatitis A (AER), Hepatitis B, HIV-1 and Polio Virus Type 1. It can disinfect a wide range of endoscopes and other heat sensitive semi-critical medical devices.



106000	* gallon	each	\$51.49	\$48.92	\$48.92	\$46.34
--------	----------	------	---------	---------	---------	---------

Test Strips

10602	100/bottle		\$89.49	\$85.02	\$85.02	\$80.54
-------	------------	--	---------	---------	---------	---------

METREX METRICIDE™

MetriCide™ Plus 30® 28 Day 3.4% Glutaraldehyde

MetriCide Plus 30 is a long life 3.4% buffered glutaraldehyde. It is the extra-strength solution for the sterilization and disinfection of immersible medical and dental equipment and instruments.



- No dilution required
- Contains a surfactant
- Recommended for heat-sensitive equipment and instruments

103200	* gallon	each	\$31.49	\$29.92	\$29.92	\$28.34
--------	----------	------	---------	---------	---------	---------

MetriCide™ 28 Day 2.5% Glutaraldehyde

MetriCide 28 is intended for use in a tray system with a variety of noncritical, semi-critical, and critical devices - including anesthesia equipment, respiratory therapy equipment, metallic equipment or instruments, rubber objects, plastic objects, and thermometers.



MetriCide 28 may be reused up to 28 days, assuming the minimum concentration of glutaraldehyde, as measured by a chemical indicator, remains within acceptable parameters and other conditions of use are met.

102800	* gallon	each	\$25.79	\$24.50	\$24.50	\$23.21
--------	----------	------	---------	---------	---------	---------

MetriCide™ 14 Day 3.4% Glutaraldehyde

MetriCide is intended for use in a tray system with a variety of noncritical, semi-critical, and critical devices - including lensed instruments, anesthesia equipment, respiratory therapy equipment, rubber objects, plastic objects, sharp instruments, thermometers, and flexible fiberoptic endoscopes.



MetriCide 14 may be reused up to 14 days, assuming the minimum concentration of glutaraldehyde, as measured by a chemical indicator, remains within acceptable parameters and other conditions of use are met.

101400	* gallon	each	\$22.29	\$21.18	\$21.18	\$20.06
--------	----------	------	---------	---------	---------	---------

MetriTest™ 1.8% Test Strips

- Easy to use monitoring system for users of high-level disinfecting and sterilizing solutions
- Quickly determines if the minimum effective concentration (MEC) of 1.8% glutaraldehyde is present during the 28 day use life of MetriCide Plus 30 28 Day and MetriCide 28 Day



10304	1 bottle	60 strips/bottle	\$79.29	\$75.33	\$75.33	\$71.36
-------	----------	------------------	---------	---------	---------	---------

MetriCide™ 1.5% Test Strips

MetriTest 1.5% Test Strips are glutaraldehyde concentration monitor for use in glutaraldehyde-containing germicide solutions with a minimum effective concentration (MEC) of 1.5% glutaraldehyde. It is recommended to test the glutaraldehyde solution prior to each use to help prevent healthcare-acquired infections (HAIs). Use with MetriCide 14 Day.



10303	60 strips/bottle	2 bottles/case	\$164.29	\$156.08	\$156.08	\$147.86
-------	------------------	----------------	----------	----------	----------	----------

Metrex Dispensing Pump

- For use with Metrex one gallon containers

461000		each	\$5.79	\$5.50	\$5.50	\$5.21
--------	--	------	--------	--------	--------	--------

* A freight charge of \$5.50 will be added to the price of this product.

STERILIZATION

**PASCAL DENTAL BANICIDE®
DISINFECTING SOLUTION**

2.65% Glutaraldehyde

- Ready to Use, No Activation Required
- Eliminates Mixing Mistakes
- Sterilization in 10 hours @ 22c
- High Level Disinfection in 45 Minutes @ 22c
- Does Not Contain Ortho-phenyl Phenol
- Includes Corrosion Inhibitor
- Effective against HIV
- Pleasant Lemon Scent



15200	* gallon	each	\$49.99	\$47.49	\$47.49	\$44.99
-------	----------	------	----------------	----------------	----------------	----------------

**ASEPTICO PROCIDE-D™ PLUS 28-DAY
GLUTARALDEHYDE STERILIZING
AND DISINFECTING SOLUTION**

- Procide Plus is a surfactant-containing glutaraldehyde disinfectant.
- Procide Plus is glutaraldehyde in a concentration of 3.4%.
- The activated solution contains anti-corrosion and antioxidant ingredients.



103260	* gallon	each	\$27.79	\$26.40	\$26.40	\$25.01
--------	----------	------	----------------	----------------	----------------	----------------

GC AMERICA COECIDE™ XL

2.5% Alkaline Glutaraldehyde

- High-level sterilant and disinfectant



551128	* gallon	each	\$38.29	\$36.38	\$36.38	\$34.46
--------	----------	------	----------------	----------------	----------------	----------------

**IMS INSTRUMENT
DAILY CLEAN CONCENTRATE**

IMS1220	1.1 lb. (makes 35 gal.)	each	\$26.49	\$25.17	\$25.17	\$23.84
IMS1218	5 lb. (makes 160 gal.)	each	\$75.29	\$71.53	\$71.53	\$67.76



**CETYLITE® POWER CLEANER
GENERAL PURPOSE ULTRASONIC
CLEANER/HOLDING SOLUTION**

Cetylite Power Cleaner is a highly-concentrated general purpose ultrasonic cleaner that removes blood, body fluids and debris from items prior to disinfection or sterilization. The tip-and-pour bottle measures 1 ounce of concentrate per each gallon of use-solution for a total yield of 32 gallons. Contains powerful rust inhibitor.



0130	1 quart bottle	each	\$81.49	\$77.42	\$77.42	\$73.34
------	----------------	------	----------------	----------------	----------------	----------------

**L&R GENERAL PURPOSE
CONCENTRATE**

- Non-Ammoniated

00228	* gallon	each	\$40.49	\$38.47	\$38.47	\$36.44
-------	----------	------	----------------	----------------	----------------	----------------

**L&R INSTRUMENT-MILK®
LUBRICANT**

- Keeps box-locks free, providing vital lubrication with regular use
- Inhibits instrument rust, discoloration and electrolytic mineral deposit build-up
- Reduces repair and replacement; eliminates instrument freeze-up
- Keeps instruments clean; lubricant film prevents build-up of baked-on protein and mineral deposits for more effective detergent cleaning
- Eliminates sterilization lubricant interference; steam permeable



00076	* gallon	each	\$36.99	\$35.14	\$35.14	\$33.29
00139	8 oz. ready-to-use spray	each	\$19.79	\$18.80	\$18.80	\$17.81

**CETYLITE® CETYLCIDE G
DISINFECTION CONCENTRATE**



0122	32 oz. bottle	each	\$88.99	\$84.54	\$84.54	\$80.09
------	---------------	------	----------------	----------------	----------------	----------------

* A freight charge of \$5.50 will be added to the price of this product.

STERILIZATION

**COLTÉNE WHALEDENT
BIOSONIC® ULTRASONIC
SOLUTION CONCENTRATES**

The solution for your cleaning needs

For use with all ultrasonic cleaning units

- Super concentrated
- Super effective
- Super convenient

Exclusive MeterDose® Bottle

- Dispenses exactly the amount needed; no waste, no spills, no mess
- Simply squeeze, tilt and pour
- No bulky gallon jugs to store; save up to 85% in storage space
- Fast acting; no need to wait for tablets or powder to dissolve



**Super Concentrate General Purpose
Ultrasonic Cleaning Solution**

- Dissolves and removes blood, saliva and microscopic particles
- Perfect for cleaning all dental instruments-even plastic and rubber
- 16oz bottle yields 8 gallons of solution

UC30	16oz	each	\$42.79	\$40.65	\$40.65	\$38.51
------	------	------	---------	---------	---------	---------

Enzymatic Ultrasonic Cleaning Concentrate

- Minimizes hand scrubbing by removing protein soils (dried blood, tissue, saliva)
- Effectively breaks down and loosens hardened debris
- 8oz bottle yields 8 gallons of solution



UC32	8oz	each	\$40.99	\$38.94	\$38.94	\$36.89
------	-----	------	---------	---------	---------	---------

Germicidal Ultrasonic Cleaning Concentrate

- Kills HIV virus and tuberculosis in the bath while it cleans instruments
- EPA registered
- Excellent as a pre-soak
- 16oz bottle yields 16 gallons of solution

UC42	16oz	each	\$86.99	\$82.64	\$82.64	\$78.29
------	------	------	---------	---------	---------	---------

**SANAX PROTECTIVE PRODUCTS
DEFEND® ULTRASONIC
ENZYMATIC TABLETS**

- Easy to use tablets dissolve in a gal. of water



UT1000	64 tablets/box		\$45.79	\$43.50	\$43.50	\$41.21
--------	----------------	--	---------	---------	---------	---------

**TUTTNAUER CLEAN & SIMPLE™
ULTRASONIC CLEANERS**

These machines provide 60kHz of quiet cleaning power which is less aggressive to instruments while still removing all surface contamination. The 1 gallon holds a standard selection of instruments or 2 half size cassettes. The 3 gallon holds 3 full size cassettes.

- High efficiency 60kHz of cavitation for thorough and quieter cleaning
- Durable vinyl clad steel housing is easy to clean and resists scratching
- Stainless steel tank with cover
- 2.5-year warranty

**Clean & Simple
One Gallon Ultrasonic Cleaner**

- Tank Dimensions: 9-3/8"W x 4"H x 5-3/8"D
- External Dimensions: 12-5/8"W x 8-3/4"H x 9-7/8"D
- Tank Capacity: 0.8 gallons



CSU1	Ultrasonic Cleaner	each	\$596.49	\$566.67	\$566.67	\$536.84
CSU1H	with heater	each	\$684.99	\$650.74	\$650.74	\$616.49
CSU1B	1 gallon basket	each	\$107.49	\$102.12	\$102.12	\$96.74

**Clean & Simple
Three Gallon Ultrasonic Cleaner**

- Tank Dimensions: 15-1/2"W x 6"H x 9"D
- External Dimensions: 19"W x 13-1/2"H x 16"D
- Tank Capacity: 3 gallons



CSU3	Ultrasonic Cleaner	each	\$1,425.99	\$1,354.69	\$1,354.69	\$1,283.39
CSU3B	3 gallon basket	each	\$230.79	\$219.25	\$219.25	\$207.71

**TUTTNAUER
CLEAN & SIMPLE™
ULTRASONIC
CLEANING TABLETS**

- Eliminate handling of heavy 1 gallon containers
- No pre-mixing required
- Enzymatic solution effectively breaks down protein



CS0064	64 tablets/box		\$61.49	\$58.42	\$58.42	\$55.34
--------	----------------	--	---------	---------	---------	---------

STERILIZATION

**L&R ULTRASONICS
QUANTREX
ULTRASONIC
CLEANING SYSTEMS**



L&R's Quantrex® Ultrasonic Cleaning Systems represent the most advanced technology in ultrasonic cleaning. The Quantrex machines offer increased power for greater cavitation, more standard features and the most comprehensive warranty in the industry today.

Quantrex systems are manufactured with vinyl-clad steel housings, stainless steel tanks, and only the highest quality components available. Each tank and body is carefully inspected for strength and durability. The contemporary design and quiet operation complement any environment.

The Quantrex Systems' standard features include a stainless steel cover, automatic tuning and a pilot light. The Quantrex 140 and larger machines are equipped with a stainless steel, multi-position drain for convenient, easy solution disposal. A 30-minute/hold timer allows you to preset the cleaning cycle duration, or select continuous cleaning. Further, all Quantrex machines have earned worldwide approvals and meet FCC requirements.

Each Quantrex System is backed by L&R's exclusive "Double Warranty". The entire unit is covered for defects in materials and workmanship for a period of 2-½ years. That's 30 months of continuous use. Additionally, the housing is warranted for 5 years and the tank-to-transducer bonding is covered for 10 years.

- "On-Off" Switch, pilot light and drain standard.
- Timer and heater optional.
- Tank internal dimensions: 9-3/8" x 5-3/8" x 4"
- Overall dimensions: 10-3/8" x 6-3/8" x 8"
- Peak output: 270 watts; 43 kHz freq.
- Tank capacity: 3-3/8 quarts
- Input: 117 VAC, 50/60 Hz Avg. Watts 0.8 amp., 95 watts, with heat 150 watts
- Basket: sold separately
- Weight: 9 pounds

Model Quantrex 140 Ultrasonic Cleaning System

00311	each	\$747.79	\$710.40	\$710.40	\$673.01
-------	------	----------	----------	----------	----------

**L&R ULTRASONICS
QUANTREX ULTRASONIC
STAINLESS STEEL WIRE BASKET**



Stainless Steel Basket for Q310

13518	14-15/16" x 8" x 4-1/2"	each	\$289.79	\$275.30	\$275.30	\$260.81
-------	-------------------------	------	----------	----------	----------	----------

Stainless Steel Wire Basket for Q140

11088	8-3/4" x 4-3/4" x 3"	each	\$129.79	\$123.30	\$123.30	\$116.81
-------	----------------------	------	----------	----------	----------	----------

**L&R QUANTREX
ULTRASONIC BEAKERS**



- For use with specialized cleaning solutions.
- Cover with openings allows beakers to be placed (suspended) in tank for specialized procedures.
- Stainless steel covers contain odors and minimize evaporation.

40006	600ml	each	\$30.49	\$28.97	\$28.97	\$27.44
-------	-------	------	---------	---------	---------	---------

**CETYLITE LABS
SOAKING TUB**



C-Tub Instrument Receptacle

- Holds up to 1 gallon of use-solution
- Requires limited counter space
- Guaranteed against breakage in normal use
- Transparent lid reduces technician exposure to chemicals
- Embossed measuring markers for in receptacle mixing
- Dimensions: Tub = 13.5" Length x 7.5" Width x 6.25" Depth
Inner Basket = 13.5" Length x 5" Width x 4" Depth
Mini Basket = 5" Length x 1.94" Width x 1.44" Depth

CET00106	C-Tub	each	\$71.79	\$68.20	\$68.20	\$64.61
----------	-------	------	---------	---------	---------	---------

**BECAUSE WE JUST
CAN'T FIT EVERYTHING
INTO THIS CATALOG!**

VISIT ACESOUTHERN.COM TODAY



STERILIZATION

TUTTNAUER 2340M AND 2540M MANUAL AUTOCLAVES

Tuttnauer manual autoclaves are designed to compliment any healthcare facility. The bright, easy to use panel was designed with the operator in mind.

Manual



Features:

- Automatic shut-off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- 1-year parts and labor warranty

2340M	each	\$4,642.79	\$4,410.65	\$4,410.65	\$4,178.51
2540M	each	\$5,499.99	\$5,224.99	\$5,224.99	\$4,949.99

	Dimensions			Volume / Capacity			Cycle Time	Power Requirments	
	Chamber D x L	Overall D x W x H	Tray D x W x H	Internal Volume	Standard Cassettes	Number of Trays	Standard Unwrapped	Voltage (V) Frequency (Hz)	Power (W) Current (A)
2340M	9" x 18"	21.5" x 20" x 14.4"	16.3" x 6.7" x 0.8"	5 gal / 19L	2 full 2 half	3	Cold: 27 mins Hot: 13 mins	120V 50/60Hz	1400W 12A
2540M	10" x 19"	21.5" x 20" x 14.4"	16.3" x 6.7" x 0.8"	6 gal / 23L	3 full 3 half	4	Cold: 30 mins Hot: 14 mins	120V 50/60Hz	1400W 12A

TUTTNAUER EZ10™ FULLY AUTOMATIC AUTOCLAVE WITH INTERNAL PRINTER

With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer Fully Automatic Autoclaves fill, sterilize, exhaust and dry at the touch of a button. Forget it once you set it.

Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying
- Hepa filter 0.2µm air filter
- Pre-programmed cycles for: unwrapped instruments; wrapped/packs; glassware and additional drying cycle
- All program parameters can be customized and stored
- Automatic shut off at the end of both the sterilization and dry cycles
- Audible and visual cycle interruption alert
- Double safety locking device prevents door from opening while chamber is pressurized
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir
- Constant monitoring of temperature and pressure for sterilization
- Power outage recovery system
- Door design prevents steam from coming in contact with control panel
- 510k clearance for dental handpieces
- 2-year parts and labor warranty
- Printer includes date, time, temperature and pressure

Automatic



EZ10	each	\$6,199.99	\$5,889.99	\$5,889.99	\$5,579.99
------	------	------------	------------	------------	------------

	Dimensions			Volume / Capacity			Cycle Time	Power Requirments	
	Chamber D x L	Overall D x W x H	Tray D x W x H	Internal Volume	Standard Cassettes	Number of Trays	Standard Unwrapped	Voltage (V) Frequency (Hz)	Power (W) Current (A)
EZ10™	10" x 19"	21.5" x 20" x 14.4"	16.3" x 6.7" x 0.8"	6 gal / 23L	3 full 3 half	4	Cold: 30 mins Hot: 14 mins	120V 50/60Hz	1400W 12A

TUTTNAUER EZ10K™ FULLY AUTOMATIC KWIKLAVE WITH INTERNAL PRINTER

The Tuttnauer EZ10k Kwiklave offers all the proven, dependable features of Tuttnauer's EZ10, with the added benefit of:

- Reducing the total sterilization cycle by about 50%.* The perfect choice for quick turn-around or a high volume office.
- Speed and reliability are accomplished without sacrificing load size.
- Requires a voltage between 220V and 235V for proper operation.
- When applicable, a buck/boost transformer (0.5KVA) may be required.

Automatic




EZ10K	each	\$6,599.99	\$6,269.99	\$6,269.99	\$5,939.99
-------	------	------------	------------	------------	------------

	Dimensions			Volume / Capacity			Cycle Time	Power Requirments	
	Chamber D x L	Overall D x W x H	Tray D x W x H	Internal Volume	Standard Cassettes	Number of Trays	Standard Unwrapped	Voltage (V) Frequency (Hz)	Power (W) Current (A)
EZ10K™	10" x 19"	21.5" x 20" x 14.4"	16.3" x 6.7" x 0.8"	6 gal / 23L	3 full 3 half	4	Cold: 21 mins Hot: 11 mins	230V 50/60Hz	2200W 10A

*The overall cycle is reduced, sterilization and exposure time is never reduced.

STERILIZATION

TUTTNAUER AUTOCLAVE POUCH RACKS



All pouched instruments should be placed in a rack for proper sterilization and drying. This stainless steel pouch rack is convenient for the loading and unloading of bagged instruments.

- 15" long x 6-1/2" wide x 3" high
- For 9" and 10" Chambers

AR910	Pouch rack	2/box	\$102.99	\$97.84	\$97.84	\$92.69
-------	------------	-------	----------	---------	---------	---------

TUTTNAUER AUTOCLAVE REPLACEMENT DOOR GASKETS

02610118	for model #2340	each	\$57.29	\$54.43	\$54.43	\$51.56
02610023	for model #2540	each	\$60.49	\$57.47	\$57.47	\$54.44

TUTTNAUER CHAMBER BRITE™ POWDERED AUTOCLAVE CLEANER



One of the world's most effective and easy to use autoclave cleaners. Chamber Brite™ effectively cleans and descales the chamber, reservoir water tank, and copper tubing of any steam autoclave. Regular cleaning extends the life of autoclaves.

- Pre-measured packs, no pre-mixing required
- Regular usage promotes efficiency
- Odor free

CB0010	10 packs/box		\$61.49	\$58.42	\$58.42	\$55.34
--------	--------------	--	---------	---------	---------	---------

HU-FRIEDY® SPORECHECK™

24 Hour In-Office Test

In office biological monitoring allows you to take control of your steam sterilizer testing and provides you with the peace of mind that your sterilizer is functioning properly with results in 24 hours.



- The SporeCheck™ system is easy to use; one biological indicator (BI) vial works for all types of steam sterilizers, including gravity, pre-vacuum and flash sterilization.
- Test results are easy to interpret with a visible color change for positive results, leaving no room for error.
- SporeCheck biological indicator vials are a self-contained, easy to use system which greatly reduces any chance for errors.
- The CDC guidelines recommend that steam sterilizers should be monitored at least weekly, preferably daily and with every load that contains implantable devices.
- SporeCheck™ Starter Kit
- Contains:
 - Self-Contained Steam Biological Indicator Box (25 Vials)
 - Dry Block Incubator
 - Record Notebook
 - Vial Crusher

IMS1376	SporeCheck™ Starter Kit	each	\$285.49	\$271.22	\$271.22	\$256.94
---------	-------------------------	------	----------	----------	----------	----------

CERTOL INTERNATIONAL PROCHEK® S STEAM STERILIZATION INTEGRATORS



- Integrators are more resistant to color change
- Monitors effective penetration of steam and heat over time into thick pads
- Dramatic color change from yellow to purple when full parameters are met


TS14401	100/box		\$44.79	\$42.55	\$42.55	\$40.31
---------	---------	--	---------	---------	---------	---------

3M ATTEST™ BIOLOGICAL INCUBATOR

- Convenient system provides the products you need to monitor the steam sterilization cycle.
- Self-contained biological indicators significantly reduce the possibility of contamination, minimizing false positives and ensuring more accurate results.
- Incubator automatically crushes vial to release medium, making it safer and easier to use than other systems.
- 14-vial capacity to meet the needs of both small and large offices.

3M116	each		\$203.99	\$193.79	\$193.79	\$183.59
-------	------	--	----------	----------	----------	----------

DEFEND® BY MYDENT INTERNATIONAL IMPULSE HEAT SEALER



- Seals items up to 8" wide x 1/8"
- Unit size: 12" x 3" x 5"
- 110V
- If using paper/plastic sealing material – place the plastic on the bottom, nearest the heating element.

HS3000	each		\$147.99	\$140.59	\$140.59	\$133.19
--------	------	--	----------	----------	----------	----------

FORMING AN ALLIANCE
Rewarding Loyalty



ALLIANCE PROGRAM

See page 2 for Details

Tier Pricing Details see page 2	Regular Price	Member Price	Premier Price	Premier+ Price
---------------------------------	---------------	--------------	---------------	----------------

STERILIZATION

**HEALTHLINK PROSURE®
MAIL-IN STERILIZER
MONITORING SERVICE**

ProSure is a mail-in biological indicator program for third party verification on your steam, dry heat, vapor or EO sterilizers.



3910	12 tests/box	\$145.79	\$138.50	\$138.50	\$131.21
------	--------------	----------	----------	----------	----------

**CROSSTEX CONFIRM®
MAIL-IN STERILIZER
MONITORING SERVICE**

The most comprehensive mail-in service in healthcare, with results in just 24 hours! Everything you need to verify sterilizer performance, even to help you stay on your testing schedule.



- Test results in 24 hours - Steam cycle test results in 24 hours. Your choice of immediate test failure notification by telephone, fax or email. Complete sterilization test results are faxed or emailed quarterly.
- Check results online with our secure, confidential website.
- Free reminder program helps you keep on your regular testing schedule.
- Third party documentation of test results provides a tangible record, as well as a legal document, when proof of sterility must be provided.
- One of the most cost-effective methods of biological monitoring.
- Compatible with steam, dry heat, chemical and EtO.

Mail-In Test

Includes: 3 Bacterial Spore Test Strips (2 Tests/1 Control), Postage Prepaid

CST120	Mail-In Test	12 tests/box	\$115.29	\$109.53	\$109.53	\$103.76
--------	--------------	--------------	----------	----------	----------	----------

Mail-In Value Test

Includes: 2 Bacterial Spore Test Strips (1 Test/1 Control), Postage Not Included

CVT520	Mail-In Value Test	52 tests/box	\$222.29	\$211.18	\$211.18	\$200.06
--------	--------------------	--------------	----------	----------	----------	----------

**HU-FRIEDY® SPORECHECK®
24-HOUR IN-OFFICE BIOLOGICAL
MONITOR TEST INDICATOR VIALS**

SporeCheck is a steam sterilizer test system that provides reliable, complete results in 24 hours. Spore growth can be detected in as few as 2.5 hours. The SporeCheck system can be used with all types of steam sterilizers, including gravity, pre-vacuum and Statim or rapid steam sterilizers.



IMS1373		25/box	\$123.29	\$117.13	\$117.13	\$110.96
---------	--	--------	----------	----------	----------	----------

**CROSSTEX CONFIRM® 10
IN-OFFICE BIOLOGICAL
MONITORING SYSTEM**

A breakthrough in biological monitoring...

final results in just 10 hours! Now you can meet biological monitoring guidelines for steam sterilizers in less than half the time with ConForm 10, the first in-office system that delivers final results in just 10 hours.



- Easy to read - Simply process indicator in a normal autoclave cycle followed by a 10 hour incubation cycle
- Easy to read results - Following the incubation cycle, if the media turns yellow, this indicates a failed test
- Fast - Quickest final read on the market

Starter Kit

C10SK	each	\$481.79	\$457.70	\$457.70	\$433.61
-------	------	----------	----------	----------	----------

Steam Biological Indicators

C10BI25	25/box	\$128.19	\$121.78	\$121.78	\$115.37
---------	--------	----------	----------	----------	----------

ConForm Record Keeper Booklet

CRK02	each	\$6.49	\$6.17	\$6.17	\$5.84
-------	------	--------	--------	--------	--------

**3M ATTEST™ BIOLOGICAL
INDICATORS - 48 HOUR READOUT**

- Biological indicator vials are easy to use and interpret with visual color readout in 48 hours
- Self-contained biological indicators significantly reduce the possibility of contamination, minimizing false positives and assuring more accurate results
- For 250°F/121°C gravity and 270°F/132°C vacuum assisted steam sterilizers



Brown Cap

1262P	25/box	\$78.29	\$74.38	\$74.38	\$70.46
1262	100/box	\$207.49	\$197.12	\$197.12	\$186.74

**3M ATTEST™ BIOLOGICAL
INDICATORS - 24 HOUR READOUT**

- Biological indicator vials are easy to use and interpret with visual color readout in 24 hours
- Self-contained biological indicators significantly reduce the possibility of contamination, minimizing false positives and assuring more accurate results
- For 270°F/132°C gravity displacement steam sterilizers



Blue Cap

1261P	25/box	\$72.49	\$68.87	\$68.87	\$65.24
1262	100/box	\$207.49	\$197.12	\$197.12	\$186.74

3M ATTEST™ LOG BOOKS

- Complete record keeping system that maintains your sterilization process information for steam sterilization
- 50 Record Charts per book

Log Book with 50 Record Charts

- Use in recording results for either 3M Attest™ 1261 or 1262 biological indicators

1266A	50 Charts	each	\$15.49	\$14.72	\$14.72	\$13.94
-------	-----------	------	---------	---------	---------	---------

Log Book with 50 Record Sheets

- Use in recording results for either 3M Attest™ Rapid Readout 1291 or 1292 biological indicators

1266	50 Sheets	each	\$18.49	\$17.57	\$17.57	\$16.64
------	-----------	------	---------	---------	---------	---------

STERILIZATION

**CROSSTEX
STEAMPLUS™
TYPE 5 CHEMICAL
INTEGRATORS**



STEAMPlus™ is an immediate readout integrator for use in gravity, prevacuum, and immediate use cycles.

- This product is tested to ST-60 and ISO 11140-1:2005 performance claims of Type 5 integrating integrators
- 4" strip has distinct pass/fail criteria and allows for easy removal from packs
- The STEAMPlus™ is FDA cleared as equivalent in performance to a biological indicator, and offers the highest level of sterility assurance from a chemical indicator.

The STEAMPlus™ does not replace the routine use of a biological indicator.

SSI100	100/box	\$49.49	\$47.02	\$47.02	\$44.54
--------	---------	---------	---------	---------	---------

**3M COMPLY
(STERIGAGE™)
CHEMICAL INTEGRATORS**



- Integrates all critical parameters of steam sterilization to provide a convenient way to monitor proper exposure conditions
- Distinct "accept" or "reject" readout provides confidence that the results were read correctly
- Immediately identifies problems in specific packs
- Convenient small size to fit all packs
- Non-toxic chemicals are contained in an enclosed system, so there's no need to worry about contamination of sterile items
- Moving-front technology

1243B	2" x ¾"	100/pack	\$16.79	\$15.95	\$15.95	\$15.11
1243A	2" x ¾"	500/pack	\$75.49	\$71.72	\$71.72	\$67.94

**3M COMPLY™ STEAM
CHEMICAL INDICATOR STRIPS**



- Designed to monitor exposure to steam conditions in a pack.
- Provides a visual color change to validate that the sterilant has penetrated to the point of placement in the pack and confirm that sufficient exposure conditions have occurred.
- Steam chemical indicator changes color from white to dark brown/black or as dark or darker than the color standard when exposed to steam sterilization conditions.
- Perforated.

1250	240/box	\$39.29	\$37.33	\$37.33	\$35.36
------	---------	---------	---------	---------	---------

**OK® AUTOCLAVE
STEAM INDICATOR STRIPS**



OK Strips are indicators that monitor all critical sterilization parameters: steam, time, and temperature. In sterilization cycles, OK strips guard against incomplete steam penetration, presence of non-condensing gases or excessive air, as well as from inadequate temperature or insufficient time exposure.

264101	250/box	\$25.79	\$24.50	\$24.50	\$23.21
--------	---------	---------	---------	---------	---------

**3M COMPLY™
LEAD FREE STEAM
INDICATOR AUTOCLAVE TAPE**



Designed to seal packs and provide visual evidence that packs have been exposed to the steam sterilization process.

132212MM	0.47" wide x 60 yds	roll	\$5.79	\$5.50	\$5.50	\$5.21
132218MM	0.70" wide x 60 yds	roll	\$6.79	\$6.45	\$6.45	\$6.11
132224MM	0.94" wide x 60 yds	roll	\$8.29	\$7.88	\$7.88	\$7.46
132248MM	1.89" wide x 60 yds	roll	\$14.79	\$14.05	\$14.05	\$13.31

CROSSTEX LEAD-FREE, STEAM INDICATOR TAPE

- Made with adhesive and ink that is free of natural rubber latex, lead or solvents
- Validated for use with all common dynamic-air-removal (pre-vacuum/SFPP) and gravity steam sterilization cycles, allowing simplified stocking of one tape for all steam sterilization needs
- Tape surface can be written on with indelible non-toxic ink to conveniently document appropriate information (e.g. date, sterilizer, contents, load/cycle)
- Satisfies the requirement to have a process indicator visible on the external sterilization packaging
- Compliant to ISO 11140-1: 2014 Type 1
- Natural color

STLF18MM	0.75" x 60 yds	roll	\$18.29	\$17.38	\$17.38	\$16.46
----------	----------------	------	---------	---------	---------	---------

**3M COMPLY™ M52 TAPE
DISPENSER WITH TABBER**



The Comply M52 Tape Dispenser with tabber has a weighted base and is capable of handling three rolls of 1/2", up to two rolls no wider than 1", or one 2" roll of indicator tape. A tabber can be attached to the tape dispenser, which automatically folds a ¼" length of tape up to a 1" roll width.

M52	each	\$138.49	\$131.57	\$131.57	\$124.64
-----	------	----------	----------	----------	----------

**INSTRUMENT TRAYS CAN
BE FOUND THROUGHOUT
THE INSTRUMENT SECTION**



STERILIZATION

PROPPER CHEX-ALL® II STERILIZATION POUCHES

Instant-sealing sterilization pouches with indicators both inside and outside for Steam and ETOH Gas.



Self-Sealing

- Eliminates the need for heat sealers or tape sealing
- Provides a complete hermetic seal every time
- Makes each workstation self-sufficient

Easy to Use

- Designed to be used by all personnel
- Virtually mistake-proof – self-guiding, prefolded construction

Positive, Strong Seal

- Will not open up in steam autoclaves
- Unaffected by gas
- Permits the same sterile storage time as all pouches

Tamper-Proof

- Produces visible evidence of tampering with seals

The Only Pouch with Pre-Printed Indicators Inside and Outside

- In steam: Yellow indicators turn grey-black, red indicators don't change
- In gas: Red indicators turn green, yellow indicators don't change
- No need to insert or apply separate monitoring strips

Plastic/paper

024008	3" x 8"	250/box	\$44.79	\$42.55	\$42.55	\$40.31
024010	5" x 10"	250/box	\$72.79	\$69.15	\$69.15	\$65.51
024012	7" x 12"	250/box	\$86.79	\$82.45	\$82.45	\$78.11
024018	12" x 18"	125/box	\$121.99	\$115.89	\$115.89	\$109.79

UNIPACK SEE-THROUGH, SELF-SEAL AUTOCLAVE POUCHES

- **Triple Seal & Heat Processed:** For a seal that is consistent and strong, eliminating the risk of shattering or tearing
- **Wide Self-Seal Adhesive Strip:** Easy-seal flap to ensure the integrity of contents after sterilization process
- **Dual External Indicators:** Easy to read "Color Changing Indicators" printed on each pouch, assuring proper sterilization process by steam and by ethylene oxide
- **Built-In Internal Indicator:** For Steam Process-to ensure the proper temperature is reached inside the pouch
- **Complete Film/Paper Separation:** Clean separation (peeling) of the 4-ply see-through film (from the medical grade paper) maintains the sterility of the contents until used
- **Blue Tinted Transparent Film:** For easier identification of contents and possible punctures by contrast of color

UPP225	2-¼" x 4"	200/box	\$6.79	\$6.45	\$6.45	\$6.11
UPP2710	2-¾" x 9"	200/box	\$7.79	\$7.40	\$7.40	\$7.01
UPP3212	3-¼" x 12"	200/box	\$8.79	\$8.35	\$8.35	\$7.91
UPP3565	3-½" x 5-¼"	200/box	\$7.29	\$6.93	\$6.93	\$6.56
UPP3510	3-½" x 9"	200/box	\$7.79	\$7.40	\$7.40	\$7.01
UPP4212	4-¼" x 12"	200/box	\$11.79	\$11.20	\$11.20	\$10.61
UPP5211	5-½" x 10"	200/box	\$12.49	\$11.87	\$11.87	\$11.24
UPP7513	7-½" x 13"	200/box	\$23.79	\$22.60	\$22.60	\$21.41
UPP1015	10" x 15"	200/box	\$38.49	\$36.57	\$36.57	\$34.64
UPP1218	12" x 18"	200/box	\$57.29	\$54.43	\$54.43	\$51.56

DYNAREX SELF-SEALING STERILIZATION POUCHES

- Self-sealing - no tape or heat sealing needed
- Blue-tinted transparent cover makes it easy to detect contamination
- Internal and external color-changing dual indicators assure proper processing in chemical vapor, steam, or ethylene oxide (EO) gas sterilization
- Not made with natural rubber latex



DY4460	2-¼" x 2.¾"	200/box	\$8.49	\$8.07	\$8.07	\$7.64
4461	3-½" x 5-¼"	200/box	\$7.29	\$6.93	\$6.93	\$6.56
4462	3-½" x 9"	200/box	\$10.29	\$9.78	\$9.78	\$9.26
4464	5-¼" x 10"	200/box	\$14.49	\$13.77	\$13.77	\$13.04
4465	7-½" x 13"	200/box	\$25.79	\$24.50	\$24.50	\$23.21
DY4466	12" x 16"	200/box	\$39.49	\$37.52	\$37.52	\$35.54

CROSSTEX DUO-CHECK® SELF-SEAL STERILIZATION POUCHES

- Built-in internal/external process indicators that change color when exposed to a specific temperature
- Manufactured from medical grade paper with extra-wide chevron to protect against instrument breaching, thumb-notch for easy opening, 6 corner tack seals to prevent curling and other potential contaminants from entering pouch upon opening and a perforated fold to help ensure the adhesive strip forms an airtight seal
- Once sterilized, contents will remain sterile for up to 3 years provided the integrity of the pouch has not been compromised
- For use with steam and EtO sterilization processes only
- Not made with natural rubber latex



SCXX	2-¼" x 4"	200/box	\$10.99	\$10.44	\$10.44	\$9.89
SCX	2-¾" x 9"	200/box	\$13.79	\$13.10	\$13.10	\$12.41
SCXS	3-½" x 5-¼"	200/box	\$14.29	\$13.58	\$13.58	\$12.86
SCS	3-½" x 9"	200/box	\$13.99	\$13.29	\$13.29	\$12.59
SCM	5-¼" x 10"	200/box	\$21.29	\$20.23	\$20.23	\$19.16
SCL	7-½" x 13"	200/box	\$41.29	\$39.23	\$39.23	\$37.16
SCL1015	10" x 15"	100/box	\$53.49	\$50.82	\$50.82	\$48.14
SCL1215	12" x 15"	100/box	\$53.99	\$51.29	\$51.29	\$48.59
SCL1218	12" x 18"	100/box	\$53.99	\$51.29	\$51.29	\$48.59

White/Pink

SCSPP	3-½" x 9"	200/box	\$12.79	\$12.15	\$12.15	\$11.51
-------	-----------	---------	---------	---------	---------	---------

CARDINAL HEALTH DUAL PEEL SELF-SEAL STERILIZING POUCHES

- With Steam and Gas Indicators
- Plastic/Paper



92308	3-½" x 8"	200/box	\$17.49	\$16.62	\$16.62	\$15.74
92510	5-¼" x 10"	200/box	\$18.49	\$17.57	\$17.57	\$16.64
92713	7-½" x 13"	200/box	\$30.29	\$28.78	\$28.78	\$27.26

STERILIZATION

CROSSTEX SURE-CHECK® SELF-SEAL STERILIZATION POUCHES



Built-in multi-parameter internal/external process indicators that change color when all three criteria for sterilization have been met - time, temperature, and pressure (steam).

- Manufactured from medical grade paper with extra-wide chevron to protect against instrument breaching, thumb-notch for easy opening, 6 corner tack seals to prevent curling and other potential contaminants from entering pouch upon opening and a perforated fold to help ensure the adhesive strip forms an airtight seal
- Once sterilized, contents will remain sterile for up to 3 years provided the integrity of the pouch has not been compromised
- For use with steam and EtO sterilization processes only

SCXX2	2-¼" x 4"	200/box	\$11.99	\$11.39	\$11.39	\$10.79
SCX2	2-¾" x 9"	200/box	\$13.79	\$13.10	\$13.10	\$12.41
SCXS2	3-½" x 5-¼"	200/box	\$14.79	\$14.05	\$14.05	\$13.31
SCS2	3-½" x 9"	200/box	\$14.29	\$13.58	\$13.58	\$12.86
SCL3222	3-½" x 22"	100/box	\$20.79	\$19.75	\$19.75	\$18.71
SCW2	5-¼" x 6-½"	200/box	\$14.99	\$14.24	\$14.24	\$13.49
SCM2	5-½" x 10"	200/box	\$22.79	\$21.65	\$21.65	\$20.51
SCL5152	5-¼" x 15"	200/box	\$45.79	\$43.50	\$43.50	\$41.21
SCL2	7-½" x 13"	200/box	\$45.79	\$43.50	\$43.50	\$41.21
SCL8162	8" x 16"	200/box	\$73.49	\$69.82	\$69.82	\$66.14
SCL10152	10" x 15"	100/box	\$56.79	\$53.95	\$53.95	\$51.11
SCL10172	10-½" x 17"	100/box	\$52.29	\$49.68	\$49.68	\$47.06
SCL12152	12" x 15"	100/box	\$54.79	\$52.05	\$52.05	\$49.31
SCL12182	12" x 18"	100/box	\$56.79	\$53.95	\$53.95	\$51.11

FOX STERILIZING BAGS



- With Steam or Gas (EtOH) Indicators
- White mini-pore Kraft Paper both sides
- For sterilizing by steam or gas (no dry heat)
- Tear resistant
- Not self-sealing
- Thermal indicator ink marked for evidence of autoclaving

Small

For cotton balls, eye pads, eye sponges, forceps, instruments, sponges, 2cc and 5cc syringes and needles

150000101	2-½" x 1-½" x 7-½"	1000/case	\$73.79	\$70.10	\$70.10	\$66.41
-----------	--------------------	-----------	---------	---------	---------	---------

Medium

For drinking tubes, instruments, forceps, scissors, 10cc-50cc syringes, surgeons brushes and needles

150000103	2-½" x 1-½" x 10-½"	300/case	\$83.29	\$79.13	\$79.13	\$74.96
-----------	---------------------	----------	---------	---------	---------	---------

Large

For forceps, glove envelopes and instruments

100000104	6-¼" x ¾" x 14-¾"	250/case	\$153.99	\$153.99	\$153.99	\$153.99
-----------	-------------------	----------	----------	----------	----------	----------

CROSSTEX SELF-SEALING AUTOCLAVE BAGS



- Color process indicators - steam, chemical vapor, EtO, dry heat
- Unique hospital double-fold seal with pull-off adhesive tab
- Gusset feature provides 50% more inner space

OCS	2-½" x 1-½" x 8-¾"	1,000/box	\$64.99	\$61.74	\$61.74	\$58.49
OC	2-½" x 1-½" x 10-½"	1,000/box	\$71.29	\$67.73	\$67.73	\$64.16

CROSSTEX SANI-TUBE® STERILIZATION TUBING



- Use in steam, chemical vapor, dry heat (up to 425°F/218°C)
- Manufactured from clear, puncture-resistant, 1.5 mil nylon
- Seal with process-indicator tape or heat sealer

Nylon with process-indicators:

SI2	2" x 100'	each	\$11.49	\$10.92	\$10.92	\$10.34
SI3	3" x 100'	each	\$14.99	\$14.24	\$14.24	\$13.49
SI4	4" x 100'	each	\$19.49	\$18.52	\$18.52	\$17.54
SI6	6" x 100'	each	\$27.29	\$25.93	\$25.93	\$24.56

Nylon without process-indicators:

SN2	2" x 100'	each	\$10.79	\$10.25	\$10.25	\$9.71
SN3	3" x 100'	each	\$13.79	\$13.10	\$13.10	\$12.41
SN4	4" x 100'	each	\$18.79	\$17.85	\$17.85	\$16.91
SN6	6" x 100'	each	\$25.79	\$24.50	\$24.50	\$23.21

UNIPACK SEE-THROUGH AUTOCLAVE TUBING

- One-side clear film for easy identification of contents
- One-side medical Kraft paper, good wet strength excellent for steam sterilization
- Available in flat and gusseted types
- With sterilization indicators - If steam or chemi-vapor sterilized, the indicating color changes to brown; If EtO sterilized, the indicating color changes to yellow

UPT002	2" x 100'	each	\$7.29	\$6.93	\$6.93	\$6.56
UPT004	4" x 100'	each	\$9.79	\$9.30	\$9.30	\$8.81
UPT006	6" x 100'	each	\$12.49	\$11.87	\$11.87	\$11.24

CARDINAL HEALTH™ STERILIZATION TUBES



- Dual Peel Sterilization Tubes
- With steam and gas indicators
- Paper on one side, clear on the other side

90003	3"	100'/roll	\$20.29	\$19.28	\$19.28	\$18.26
-------	----	-----------	---------	---------	---------	---------

STERILIZATION

HU-FRIEDY® IMS STERILIZING WRAP

- Color: Blue
- Single-Use
- Nonwoven Fabric



IMS1211	15" x 15"	1,000/case	\$175.29	\$166.53	\$166.53	\$157.76
IMS1212	20" x 20"	500/case	\$155.49	\$147.72	\$147.72	\$139.94
IMS1210	24" x 24"	500/case	\$199.49	\$189.52	\$189.52	\$179.54

CARDINAL HEALTH BIO-SHIELD® REGULAR STERILIZATION WRAP

- Ideal for lighter specialty trays, pitchers, small basins and huck towels

4012	12" x 12"	1,000/case	\$69.29	\$65.83	\$65.83	\$62.36
4015	15" x 15"	500/case	\$66.79	\$63.45	\$63.45	\$60.11
4018	18" x 18"	500/case	\$109.99	\$104.49	\$104.49	\$98.99
4020	20" x 20"	500/case	\$90.29	\$85.78	\$85.78	\$81.26
4024	24" x 24"	500/case	\$106.49	\$101.17	\$101.17	\$95.84
4036	36" x 36"	250/case	\$201.99	\$191.89	\$191.89	\$181.79
4040	40" x 40"	250/case	\$178.29	\$169.38	\$169.38	\$160.46

FOX STERILIZING WRAP

- Made from heavy creped stock high in wet strength
- Color: white
- Does not disintegrate in autoclave
- Allows full steam penetration
- Offers adequate protection of sterile items during storage

100001206	20"	500/case	\$91.99	\$87.39	\$87.39	\$82.79
-----------	-----	----------	---------	---------	---------	---------

STERIMED STERI-WRAP

- Single-Use
- Nonwoven Fabric
- Color: Blue

1010360B1515	15" x 15"	1,000/case	\$158.79	\$150.85	\$150.85	\$142.91
1010360B2020	20" x 20"	500/case	\$122.79	\$116.65	\$116.65	\$110.51
1010360B2424	24" x 24"	500/case	\$141.79	\$134.70	\$134.70	\$127.61
1010360B3030	30" x 30"	250/case	\$141.79	\$134.70	\$134.70	\$127.61

HALYARD QUICK CHECK STERILIZATION WRAP



- Comes with white and blue colors that are bonded together - inspect trays easier and faster than ever
- Superior microbial barrier protection with patented Powerguard technology
- Cleared for all of the same sterilization systems as HALYARD KIMGUARD ONE-STEP Wrap
- Made from recyclable #5 polypropylene fabric

34170	15" x 15"	480/case	\$98.29	\$93.38	\$93.38	\$88.46
34156	18" x 18"	480/case	\$128.79	\$122.35	\$122.35	\$115.91
34192	20" x 20"	480/case	\$167.99	\$159.59	\$159.59	\$151.19
34178	24" x 24"	240/case	\$103.29	\$98.13	\$98.13	\$92.96
34176	30" x 30"	144/case	\$94.49	\$89.77	\$89.77	\$85.04
34188	45" x 45"	120/case	\$204.29	\$194.08	\$194.08	\$183.86

HALYARD H100 STERILIZATION WRAP

Independent laboratory testing assures that HALYARD® Sterilization Wrap delivers the optimum balance of properties needed to ensure the sterility of your instruments:

- Strength and tear resistance
- Flame resistance
- Microbial barrier
- Resistance to linting and abrasion
- Effective sterilant penetration
- Easy to use
- Cost-effective

37046	12" x 12"	100/pkg	\$6.79	\$6.45	\$6.45	\$6.11
37049	18" x 18"	100/bag	\$18.79	\$17.85	\$17.85	\$16.91

HALYARD KIMGUARD® SEQUENTIAL STERILIZATION WRAP



The KIMGUARD line is made of SMS sterilization fabric. It is light blue, softer and stronger than ever and still delivers the barrier protection you expect. Often used for light linen packs and gowns, towel packs, light basins and other light wrapping needs.

10712	12" x 12"	1000/case	\$52.29	\$49.68	\$49.68	\$47.06
10715	15" x 15"	1000/case	\$114.29	\$108.58	\$108.58	\$102.86
10718	18" x 18"	1000/case	\$137.79	\$130.90	\$130.90	\$124.01
10720	20" x 20"	1000/case	\$186.29	\$176.98	\$176.98	\$167.66
10724	24" x 24"	500/case	\$123.29	\$117.13	\$117.13	\$110.96
10730	30" x 30"	300/case	\$104.99	\$99.74	\$99.74	\$94.49
10736	36" x 36"	300/case	\$176.99	\$168.14	\$168.14	\$159.29
10740	40" x 40"	250/case	\$146.49	\$139.17	\$139.17	\$131.84
10745	45" x 45"	250/case	\$166.79	\$158.45	\$158.45	\$150.11
10754	54" x 54"	100/case	\$110.99	\$105.44	\$105.44	\$99.89

Additional Package Sizes Available:

37047	15" x 15"	100/bag	\$10.29	\$9.78	\$9.78	\$9.26
37050	20" x 20"	100/bag	\$17.79	\$16.90	\$16.90	\$16.01
37051	24" x 24"	100/bag	\$22.29	\$21.18	\$21.18	\$20.06

GET YOUR ORDER FASTER...

To You in 2 Days Shipping



DID YOU KNOW 99% OF ACE SOUTHERN customers receive their order within 1-2 business days of shipment? ACE SOUTHERN has four order fulfillment facilities throughout the United States to ensure your order is **TO YOU IN 2.**

ACE SOUTHERN prides itself in fulfilling your order as quickly as possible. In-stock items that are ordered and processed by 3:00pm Eastern Standard Time will be shipped same-day.

**Please note, UPS and other carriers are advising that there may be some delays in the delivery of packages due to the impact of the coronavirus (COVID-19). Thank you for your patience as every company works to adapt to this developing situation. For more information please call 1.800.624.5926.



Learn more about
ACE SOUTHERN by visiting
[ACESOUTHERN.COM](https://acesouthern.com) today!



ACE SOUTHERN
Surgical Solutions

☎ 800.624.5926 | 800.441.3100

✉ info@acesouthern.com

🌐 [ACESOUTHERN.COM](https://acesouthern.com) |   

SUCTION/IRRIGATION



IRRIGATION SOLUTIONS


0.9% Sodium Chloride*

250ml Plastic pour bottle	24/case	\$191.99	\$182.39	\$182.39	\$172.79
500ml Plastic pour bottle	24/case	\$191.99	\$182.39	\$182.39	\$172.79
1000ml Plastic pour bottle	12/case	\$95.79	\$91.00	\$91.00	\$86.21

Sterile Water*

250ml Plastic pour bottle	24/case	\$111.79	\$106.20	\$106.20	\$100.61
500ml Plastic pour bottle	18/case	\$110.29	\$104.78	\$104.78	\$99.26
1000ml Plastic pour bottle	16/case	\$99.99	\$94.99	\$94.99	\$89.99
1500ml Plastic pour bottle	9/case	\$75.79	\$72.00	\$72.00	\$68.21
1000ml Bag	12/case	\$83.29	\$79.13	\$79.13	\$74.96
3000ml Bag	4/case	\$112.99	\$107.34	\$107.34	\$101.69

*A freight charge of \$8.00 will be added to the price of this product.



ACE VARIJET™ SUCTION/IRRIGATION SYSTEM

- Single-handed suction with controlled irrigation.
- Flush and clean site with one convenient instrument.
- Improved vision to the site by using fewer instruments.

0900090	3mm, 8" (20cm)	each	\$299.79	\$284.80	\$284.80	\$269.81
---------	----------------	------	----------	----------	----------	----------

13



ACE SYRINGE TIPS AND BULBS

- Single-handed suction with controlled irrigation
- Flush and clean site with one convenient instrument
- Improved vision to the site by using fewer instruments

Autoclavable Stainless Steel Syringe Tips

0900102	Water syringe tip - 1.9mm opening, 4-3/4" (12cm)	each	\$17.49	\$17.49	\$15.74	\$15.74
0900103	Air syringe tip - 2.2mm opening, 4-1/2" (11cm)	each	\$17.49	\$17.49	\$15.74	\$15.74

Replacement Bulbs for ACE Syringe Tips

Steam sterilize at 270°F for no longer than 20 minutes

0900108	Bulb only, black, 60 cc	10/box	\$68.99	\$65.54	\$65.54	\$62.09
---------	-------------------------	--------	---------	---------	---------	---------



VISTA DENTAL AIR/WATER SYRINGE TIPS

Quick-change

- Autoclavable
- Durable, stainless steel construction

315375		5/box	\$31.79	\$31.79	\$28.61	\$28.61
--------	--	-------	---------	---------	---------	---------



UNIPACK DISPOSABLE AIR/WATER SYRINGE TIPS

- Rounded edge for maximum patient comfort
- Convenient dispenser bottle keeps work area clear and assists in preventing cross-contamination
- Rounded non-scratch end to minimize wear and tear of O-ring in hand piece
- Every bag contains 2 dispenser bottles with 125 each
- No adapter required

UAT9052	Plastic w/ plastic core	250/bag	\$14.49	\$13.77	\$13.77	\$13.04
---------	-------------------------	---------	---------	---------	---------	---------

SUCTION/IRRIGATION

UNIPACK COLOR-BRITE DISPOSABLE AIR/WATER SYRINGE TIPS

- Beveled edge for maximum patient comfort
- Two grooves, one for locking indication, the other for securing a tight lock
- Convenient re-sealable pouch contains an assorted arrangement of colors
- Quad core air-channel for ultimate performance
- Rounded non-scratch end to minimize wear and tear of O-ring in hand piece
- Insert tube stays in place during procedures
- No adapter required



UAT9053	Assorted color tips	250/bag	\$14.49	\$13.77	\$13.77	\$13.04
---------	---------------------	---------	---------	---------	---------	---------

ACE IRRIGATING SYRINGES

Disposable with Tapered Curved Tip

- Improved design for complete volume control
- Plastic



0034300	12 cc	50/box	\$23.49	\$22.32	\$22.32	\$21.14
0034302	20 cc	25/box	\$23.49	\$23.49	\$21.14	\$21.14

EXEL CURVED TIP SYRINGE

- Ungraduated with plunger stop
- Ideal for pressure irrigation



26267	12 cc Irrigation	50/box	\$20.79	\$19.75	\$19.75	\$18.71
-------	------------------	--------	---------	---------	---------	---------

COVIDIEN MONOJECT™ CURVED TIP SYRINGE

- Ungraduated
- Ideal for pressure irrigation



412012	12 cc Irrigation	50/box	\$26.29	\$24.98	\$24.98	\$23.66
--------	------------------	--------	---------	---------	---------	---------



ACE SELF-FILLING STAINLESS STEEL 10 CC IRRIGATION SYRINGE

- Includes one 2" angled cannula (replacements sold individually below)

0900101	1mm opening	each	\$163.29	\$155.13	\$155.13	\$146.96
---------	-------------	------	----------	----------	----------	----------



ACE STAINLESS STEEL IRRIGATION CANNULAS

- Luer slip fits any luer type syringe
- 45° angle

1-1/4" (3.5cm)


0900159	12 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900162	14 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900165	16 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900168	18 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900171	20 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84

2" (5cm)

0900160	12 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900163	14 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900166	16 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900169	18 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900172	20 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900175	22 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84

3" (7.5cm)

0900161	12 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900164	14 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900167	16 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900170	18 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84
0900173	20 Gauge	each	\$16.49	\$16.49	\$14.84	\$14.84



Irrigation cannula for internally irrigated systems 3" (7.5cm)

- Flushes and cleans the pathway of internally irrigated burs
- For use with syringe (see 0900101 above)
- Straight
- Blunt end

0900175S	22 Gauge	each	\$16.49	\$15.67	\$15.67	\$14.84
----------	----------	------	---------	---------	---------	---------

SUCTION/IRRIGATION



DYNAREX SUCTION CANISTERS

- Durable plastic container system designed for collecting and retaining aspirated materials
- Float shut-off with built-in filter increases flow, eliminating damaging backflow to the suction regulator
- Foam guard and float valve ultrasonically welded in place, preventing premature shutoff
- Port caps attached to lid, which saves time
- Patient port designed with 20° angle (1200 cc model) or straight port (800 cc model) to prevent kinked tubing - no separate adapter required
- Universal design easily fits wall mounts.

DY32200	800 cc Hi-flow with lid	each	\$4.79	\$4.55	\$4.55	\$4.31
32220	1200 cc Hi-flow with lid	each	\$5.79	\$5.50	\$5.50	\$5.21
32221	1200 cc Kit, above canister with 3/16" x18"x6' tubing	kit	\$6.79	\$6.45	\$6.45	\$6.11



ALLIED HEALTHCARE HANDI-VAC II DISPOSABLE COLLECTION CANISTER

- Adaptable to most portable floor suction units and wall suction units
- May be hung on the wall with a canister holder adapter for regulators having a hose stem inlet
- Attaches directly to a pipeline suction regulator having a DISS male vacuum inlet
- May be used with portable electric suction units

1200 cc Hi-Flow with Lid

20080004		each	\$5.79	\$5.50	\$5.50	\$5.21
		48/case	\$235.29	\$223.53	\$223.53	\$211.76



CARDINAL HEALTH™ MEDI-VAC® GUARDIAN™ SUCTION CANISTERS

An automatic overflow shutoff valve is located inside the lid to help prevent cross contamination of regulators and wall vacuum outlets, and 90° adapters allow tubing to connect at right angles to help prevent kinking and impeded fluid flow. The locking lid encourages the proper disposal of collected wastes and enhances worker safety. Accessory port and orthopedic ports are included.

1200 cc

65651212		each	\$4.99	\$4.74	\$4.74	\$4.49
----------	--	------	--------	--------	--------	--------



BEMIS® QUICK-FIT™ LINER SYSTEM

Liner System:

- Provides easy insertion and safer removal of liners
- Handle offers easy removal and disposal
- Preassembled, sealed lid saves time and eliminates confusion
- Innovative ribbed design and superior tear strength allow easy insertion and removal
- Color coded lids and brackets allow easy visual selection
- Liner unit comes compressed for packaging for maximum use of storage space
- 1500 cc

Suction Liner w/Attached Lid

1504		each	\$4.79	\$4.55	\$4.55	\$4.31
		50/case	\$323.29	\$307.13	\$307.13	\$290.96

Reusable Outer Canister w/Attached Stopcock Bracket

1500SC01		each	\$32.79	\$31.15	\$31.15	\$29.51
----------	--	------	---------	---------	---------	---------

SUCTION/IRRIGATION

BEMIS® HI-FLOW™ RIGID SUCTION CANISTERS WITH AEROSTAT® FILTERS AND FLOAT VALVE SHUTOFF

- Designed to meet high volume collection and disposal needs
- Effectively filters aerosolized microorganisms and particulate matter
- Positive shutoff valve closes immediately when fluid reaches the top
- Foam guard keeps float valve securely in place; prevents premature shutoff

1200 cc Canister

484410	each	\$4.79	\$4.55	\$4.55	\$4.31
	48/case	\$193.49	\$183.82	\$183.82	\$174.14

2000 cc Canister

494410	each	\$4.79	\$4.55	\$4.55	\$4.31
	48/case	\$198.29	\$188.38	\$188.38	\$178.46



BEMIS® HYDROPHOBIC RIGID SUCTION CANISTERS

- Provide superior protection and maximum performance
- Hydrophobic filter provides final filtration and overflow protection
- Breathing screen separates filter stages for improved suction
- Prefilter protects hydrophobic filter from erroneous shutoff
- Tapered protective housing guards shutoff filter from splash or foam

800 cc Hydrophobic Canister

424410	100/case	\$351.49	\$333.92	\$333.92	\$316.34
--------	----------	----------	----------	----------	----------

1200 cc Hydrophobic Canister

434410	each	\$4.99	\$4.74	\$4.74	\$4.49
	48/case	\$199.49	\$189.52	\$189.52	\$179.54

2000 cc Hydrophobic Canister

435410	each	\$4.99	\$4.74	\$4.74	\$4.49
	48/case	\$209.99	\$199.49	\$199.49	\$188.99



13


BEMIS® BRACKETS FOR SUCTION CANISTERS

Wall Plate

530510	each	\$2.99	\$2.84	\$2.84	\$2.69
--------	------	--------	--------	--------	--------

Plastic Ring for Wall Plate

532710	For 1200 cc, 2000 cc and 3000 cc canisters	each	\$3.79	\$3.60	\$3.60	\$3.41
--------	--	------	--------	--------	--------	--------



530510



532710

ALLIED HEALTHCARE CHEMETRON SUCTION ASSEMBLY UNITS

- Comes complete with suction assembly and bottle
- May be used with Chemetron bottle holders and all mounts
- Metal cover and float shut-off assembly
- Use with DISS male connection or suction tubing arrangement
- Use with 1/4" suction tubing

22104003	Suction assembly, 800-2800 cc, plastic	each	\$192.99	\$183.34	\$183.34	\$173.69
1903563	One-gallon plastic replacement bottle, 800-2800 cc	each	\$46.79	\$44.45	\$44.45	\$42.11
22100002	Suction assembly, 200-1200 cc, glass	each	\$160.99	\$152.94	\$152.94	\$144.89
22930002	1/2 Gallon glass bottle, 200-1200 cc	each	\$19.79	\$18.80	\$18.80	\$17.81



Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

SUCTION/IRRIGATION



ALLIED HEALTHCARE CHEMETRON BOTTLE TRAP ASSEMBLY PARTS

- For units purchased after 1985

22100008	Bottle cap assembly, use with 1/4" suction tubing	each	\$132.29	\$125.68	\$125.68	\$119.06
22930007	Replacement gasket only	each	\$5.79	\$5.50	\$5.50	\$5.21
22930009	Replacement polypropylene float only, white	each	\$9.79	\$9.30	\$9.30	\$8.81



ALLIED HEALTHCARE CHEMETRON VACUUM BOTTLE HOLDER WITHOUT LOCK

- Slides into bracket 22110006

22110004		each	\$23.49	\$22.32	\$22.32	\$21.14
----------	--	------	---------	---------	---------	---------



ALLIED HEALTHCARE CHEMETRON MOUNTING BRACKET FOR BOTTLE HOLDER

- 4" Square (10.16 cm)
- Anodized aluminum
- Bracket for bottle holder 22110004

22110006		each	\$39.49	\$37.52	\$37.52	\$35.54
----------	--	------	---------	---------	---------	---------



ALLIED HEALTHCARE CHEMETRON ADAPTERS AND VACUUM COUPLERS

- Quick-connect inlet
- Rectangular stryker

11010017	Adapter vacuum bottle holder, two outlets	each	\$90.29	\$85.78	\$85.78	\$81.26
11010016	Vacuum bottle holder, one outlet	each	\$86.79	\$82.45	\$82.45	\$78.11



ALLIED HEALTHCARE PORTABLE SUCTION ACCESSORY ITEMS

Half-Gallon Plastic Bottle (1200 cc)

AH01903564	For use on all units incorporating bacteria filter as standard equipment. Can be steam autoclaved a limited number of times before discarding.	each	\$50.99	\$48.44	\$48.44	\$45.89
------------	--	------	---------	---------	---------	---------

One-Gallon Plastic Bottle (2800 cc)

1903563		each	\$46.79	\$44.45	\$44.45	\$42.11
---------	--	------	---------	---------	---------	---------



LAERDAL V-VAC™ MANUAL SUCTION SYSTEM

One-hand operation - no batteries required. The V-Vac unit works the instant you squeeze it!

- Adjustment for high and low-pressure suctioning
- Efficient: maximum suction exceeds 380mmHg
- Peak air flow over 70 liters per minute
- Versatile catheter connection
- Disposable one-piece cartridge provides a new, clean pump, valve, and suction tube
- Reusable handle is immersible for cleaning

Starter Kit

Includes: handle, 2 cartridges, short suction catheter w/adaptor tip, double male connector and directions for use

985000		kit	\$268.99	\$255.54	\$255.54	\$242.09
--------	--	-----	----------	----------	----------	----------

Replacement Cartridge

985001		each	\$54.79	\$52.05	\$52.05	\$49.31
--------	--	------	---------	---------	---------	---------

SUCTION/IRRIGATION

PRECISION MEDICAL EASYGO VAC ASPIRATOR PORTABLE SUCTION UNIT

- High performance unit delivers a vacuum range of 2-21 in Hg
- Battery allows 1.5 hours of run time
- Uses any standard 800 cc canister
- 2-year warranty on the unit and a 6-month warranty on the battery
- A dual pump protection system features two safeguards, a float shut-off in the canister to protect overflow and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump

Suction Unit

Includes: Suction Unit, 6' Patient Suction Tubing, 800 cc Canister, Carrying Bag, AC Hospital Grade Power Cord, Filter Assembly (Filter, Tubing and Elbow), Battery Charger/Power Supply

PM65HG	each	\$618.99	\$588.04	\$588.04	\$557.09
--------	------	----------	----------	----------	----------



ALLIED HEALTHCARE SCHUCO®-VAC 430 ASPIRATOR PORTABLE SUCTION SYSTEM

The regulator assembly has been replaced with a bacteria filter, eliminating the cap and float assembly. This filter will prevent fluids and other contaminants from damaging the pump. This durable unit continues to meet the need for controlled vacuum regulation from 0" up to 22" Hg with a flow rate of 28 lpm.

Standard features:

- Vibration-free easy-to-read regulator gauge
- Free-hanging basket accommodates canister
- Quick-connect tubing for all connections
- Use 1/4" tubing with unit
- 3 filters included with every pump
- Hydrophobic/bacteria filter eliminates pump flooding, and prevents contaminant damage
- AC operation, 115v/60Hz, 1/10 HP (does not run on a battery)

S430A	Portable suction unit	each	\$323.49	\$307.32	\$307.32	\$291.14
S1160RPL1200CC	Disposable replacement canister	4/case	\$48.79	\$46.35	\$46.35	\$43.91
01903100	Aspirator bacteria filters	3/box	\$24.29	\$23.08	\$23.08	\$21.86



ALLIED HEALTHCARE GOMCO OPTIVAC® G180 AC/DC PORTABLE SUCTION UNIT

- Long-life internal rechargeable battery eliminates the need for frequent replacement
- Sealed water-resistant unit is easy to clean and minimizes potential contamination or damage to the compressor and battery
- Angled control panel for easy viewing
- Impact-resistant and durable
- 3-year warranty on the unit, 1-year on the battery
- 3 hours run time at maximum vacuum on a fully-charged battery

Suction Unit

includes: 1500ml disposable canister, disposable suction tubing and 3 disposable hydrophobic bacteria filters

G180	each	\$797.99	\$758.09	\$758.09	\$718.19
------	------	----------	----------	----------	----------

1500ml Disposable Collection Canister, Stem Inlet

20080016	each	\$8.29	\$7.88	\$7.88	\$7.46
	16/case	\$116.49	\$110.67	\$110.67	\$104.84

Disposable Hydrophobic/Bacteria Filter, 1/4" Hose Barb X 1/4" Hose Barb

01903928	3/box	\$19.79	\$18.80	\$18.80	\$17.81
----------	-------	---------	---------	---------	---------

Suction Tubing Kit

01902000	Includes: 1-72" and 1-18" clear	kit	\$20.79	\$19.75	\$19.75	\$18.71
----------	---------------------------------	-----	---------	---------	---------	---------

Collection Canister/Tubing Kit

SU1092	Includes: 1-canister, 1-tubing set, 3-filters	kit	\$60.49	\$57.47	\$57.47	\$54.44
--------	---	-----	---------	---------	---------	---------

Replacement Battery

SU1060	each	\$91.99	\$87.39	\$87.39	\$82.79
--------	------	---------	---------	---------	---------



SUCTION/IRRIGATION

CROSSTEX EVACUATOR TIPS



Combo-Tip®

- Two tips in one- vented end on one side/non-vented end on the other
- Unique "S" tip design with NO sharp edges for maximum patient comfort

ZET	White	100/bag	\$3.79	\$3.60	\$3.60	\$3.41
-----	-------	---------	--------	--------	--------	--------

SANAX PROTECTIVE PRODUCTS VENTED ORAL EVACUATION TIPS



Convenient, single-use disposable tips contain smooth edges for maximum patient comfort. The straight handle provides increased control and dual use, as one end is vented and the other is non-vented. Fits all standard evacuator valves. Non-sterile.

VE1121	6" Long, 1/8" dia.	100/pack	\$5.79	\$5.50	\$5.50	\$5.21
--------	--------------------	----------	--------	--------	--------	--------

VISTA DENTAL SURGICAL EVACUATION TUBES



- Use for macro suction (green or blue) or micro aspiration (white)
- Manufactured in rigid plastic to provide greater cheek retraction and repeated autoclaving
- 6" long and pre-bent to 30° angle for greater access

1/8" Diameter

403201	White	25/bag	\$17.29	\$16.43	\$16.43	\$15.56
--------	-------	--------	---------	---------	---------	---------

1/4" Diameter

403202	Green	25/bag	\$17.29	\$16.43	\$16.43	\$15.56
--------	-------	--------	---------	---------	---------	---------

3/8" Diameter (Non-Vented)

403211	Blue	25/bag	\$17.29	\$16.43	\$16.43	\$15.56
--------	------	--------	---------	---------	---------	---------

KEYSTONE INDUSTRIES SURGICAL ASPIRATOR TIPS



1/8" Diameter

1201110	White	25/bag	\$4.29	\$4.08	\$4.08	\$3.86
---------	-------	--------	--------	--------	--------	--------

1/4" Diameter

1201010	Green	25/bag	\$4.29	\$4.08	\$4.08	\$3.86
---------	-------	--------	--------	--------	--------	--------

SANAX PROTECTIVE PRODUCTS SURGICAL ASPIRATOR TIPS



Made with smooth edges for maximum patient comfort. The surgical aspirator is molded at a 30° angle and fits 11mm suction holders. This item is also autoclavable up to 135 C° or may be disposed of after single use.

ST1020	Blue, 1/16" dia.	25/pack	\$5.79	\$5.50	\$5.50	\$5.21
--------	------------------	---------	--------	--------	--------	--------

DYNAREX SURGICAL ASPIRATOR TIPS



The Dynarex Surgical Aspirator Tips provide a secure fit in all evacuation systems. The Dental Aspirator Tips' versatility make them an ideal accessory for a wide range of treatments.

2216	White, 1/8" dia.	25/box	\$5.79	\$5.50	\$5.50	\$5.21
2217	Green, 1/4" dia.	25/box	\$5.79	\$5.50	\$5.50	\$5.21

COVIDIEN MONOJECT™ 444 SURGICAL ASPIRATION TIP



- Pre-bent taper tip
- Smooth, rounded tip provides patient comfort
- Opaque white plastic eliminates unpleasant viewing of evacuated debris
- Universal coupling permits use with all high speed evacuation units
- Sterile, individually packaged, disposable
- For high speed evacuation

444015		25/box	\$28.49	\$27.07	\$27.07	\$25.64
--------	--	--------	---------	---------	---------	---------

CROSSTEX SALIVA EJECTORS



- Wire-reinforced tube forms and maintains any desired shape
- Provides optimal suction without aspirating tissue
- Soft, non-removable tip - bonded to the tube for maximum patient safety

Premium Plus®

ZCB	Clear/blue	100/bag	\$5.29	\$5.03	\$5.03	\$4.76
ZWW	White/white	100/bag	\$5.29	\$5.03	\$5.03	\$4.76
ZPP	Pink/pink*	100/bag	\$5.29	\$5.03	\$5.03	\$4.76

Advantage

ZBBI	Blue/blue	100/bag	\$4.79	\$4.55	\$4.55	\$4.31
ZCBI	Clear/blue	100/bag	\$4.79	\$4.55	\$4.55	\$4.31
ZCCI	Clear/clear	100/bag	\$4.79	\$4.55	\$4.55	\$4.31
ZWWI	White/white	100/bag	\$4.79	\$4.55	\$4.55	\$4.31

SUCTION/IRRIGATION

CROSSTEX COMFORT PLUS® PREMIUM SALIVA EJECTORS



The Comfort Plus Premium Saliva Ejector uses a patented design featuring rounded edges, smooth surfaces and perfectly placed suction ports that help to enhance patient comfort while protecting delicate mucosal tissue. Unlike traditional blunt-ended saliva ejectors that can lock onto tissue and traumatize patients, the unique streamlined Comfort Plus suction tip design:

- Reduces patient discomfort
- Optimizes fluid removal
- Minimizes trauma to mucosal tissue
- Prevents tissue from blocking suction
- Allows for extended use
- Provides you with greater efficiency
- Facilitates uninterrupted fluid and light debris removal
- Perfectly bendable, yet firmly holds the shape you create

Unscented

ZWWCP	White/white	100/bag	\$5.49	\$5.22	\$5.22	\$4.94
ZCWCP	Clear/white	100/bag	\$5.49	\$5.22	\$5.22	\$4.94

Bubble Gum Scented

ZWPCPBG	White/pink	100/bag	\$5.79	\$5.50	\$5.50	\$5.21
---------	------------	---------	--------	--------	--------	--------

SANAX PROTECTIVE PRODUCTS STANDARD SALIVA EJECTORS



These pliable, shape-maintaining saliva ejectors offer a comfortable tip that minimizes tissue impingement. Fits all standard saliva ejector hose ends.

SE7002	Assorted colors	100/pack	\$6.79	\$6.45	\$6.45	\$6.11
--------	-----------------	----------	--------	--------	--------	--------

DYNAREX YANKAUER SUCTION HANDLES

- Bulb tip to minimize tissue trauma
- Ribbed connector for secure attachment to suction tubing
- Blue tinted, transparent material with slip-resistant handle
- Smooth inner surface for maximum flow
- Latex-free



Non-Vented

4691	each	\$1.29	\$1.23	\$1.23	4691
	50/case	\$41.49	\$39.42	\$39.42	\$37.34

Vented

4690	each	\$1.29	\$1.23	\$1.23	4690
	50/case	\$41.49	\$39.42	\$39.42	\$37.34

CONMED® YANKAUER SUCTION HANDLES WITH BULB TIP

- Single-use, sterile
- Transparency allows clear visualization of fluids
- Side eyes reduce traumatic tissue grab



Without Control Vent

34870	each	\$3.49	\$3.32	\$3.32	\$3.14
	50/case	\$68.79	\$65.35	\$65.35	\$61.91

Without Control Vent, With 6' Pre-Attached Tubing

0034910	each	\$5.79	\$5.50	\$5.50	\$5.21
	20/case	\$81.99	\$77.89	\$77.89	\$73.79

With Control Vent

34880	each	\$3.49	\$3.32	\$3.32	\$3.14
	50/case	\$68.79	\$65.35	\$65.35	\$61.91

With On/Off Control Switch

0038730	each	\$6.79	\$6.45	\$6.45	\$6.11
	50/case	\$144.79	\$137.55	\$137.55	\$130.31

CONMED® YANKAUER SUCTION HANDLES WITH OPEN TIP

- Open tip permits access to tight spaces
- Single-use, sterile
- Transparency allows clear visualization of fluids
- Side eyes reduce traumatic tissue grab



Without Control Vent

34970	each	\$1.79	\$1.70	\$1.70	\$1.61
	50/case	\$68.79	\$65.35	\$65.35	\$61.91

With Control Vent

34980	each	\$1.79	\$1.70	\$1.70	\$1.61
	50/case	\$68.79	\$65.35	\$65.35	\$61.91

CARDINAL HEALTH™ MEDI-VAC® TAPERED BULBOUS TIP SUCTION HANDLE

Medi-Vac® suction handles deliver performance and maneuverability. These lightweight handles are designed to provide optimal control and visibility while helping to minimize tissue trauma.

- Non-Vented



Without Control Vent

K80	each	\$1.79	\$1.70	\$1.70	\$1.61
	50/case	\$59.79	\$56.80	\$56.80	\$53.81

With Control Vent

K82	each	\$1.79	\$1.70	\$1.70	\$1.61
	50/case	\$59.79	\$56.80	\$56.80	\$53.81

SUCTION/IRRIGATION

COVIDIEN ARGYLE™ FLEXIBLE YANKAUER SUCTION CATHETER TUBE

- Sterile
- Disposable
- Latex-free

504001	10"	50/case	\$76.79	\$72.95	\$72.95	\$69.11
--------	-----	---------	---------	---------	---------	---------



COVIDIEN ARGYLE™ YANKAUER SUCTION INSTRUMENTS

Rigid

- Instruments are transparent with smooth tips and eyes
- Shatter-resistant
- Sterile, single-use



Open Tip, Non-Vented

505115	*	each	\$1.79	\$1.70	\$1.70	\$1.61
		50/case	\$70.99	\$67.44	\$67.44	\$63.89

Flexible


- Non-breakable, malleable vinyl construction with smooth inner lumen
- Sterile, single-use

High Capacity Tip, Without Vent

502005	*	50/case	\$123.99	\$117.79	\$117.79	\$111.59
--------	---	---------	----------	----------	----------	----------

*Note: Tubing sold separately

13



QUALITY ASPIRATORS ENDODONTIC ASPIRATORS

Micro suction for root canal and root tip surgery. Fits standard high volume cut-off value (11mm).

Button Provides Finger Rest for Precise Control of Aspirator, Includes One Cleaning Stylet

A	O1BA	0.5mm Opening	each	\$84.29	\$84.29	\$75.86	\$75.86
B	O2BA	1.0mm Opening	each	\$84.29	\$84.29	\$75.86	\$75.86
C	O3BA	2.0mm Opening	each	\$84.29	\$84.29	\$75.86	\$75.86

Endodontic Aspirator Made of 1/8" Diameter Stainless Steel

D	O4A	1.7mm Opening	each	\$36.49	\$36.49	\$32.84	\$32.84
E	O5A	2.3mm Opening	each	\$36.49	\$36.49	\$32.84	\$32.84

Endodontic Aspirator Made of 1/8" Diameter Stainless Steel, Tapered Tip Helps Prevent Clogging

F	O490A	1.7mm Opening With 90 right angle bend	each	\$36.49	\$36.49	\$32.84	\$32.84
G	E04	1.7mm Opening	each	\$36.49	\$36.49	\$32.84	\$32.84
H	E0490	1.7mm Opening	each	\$36.49	\$36.49	\$32.84	\$32.84

HU-FRIEDY COUPLAND ASPIRATORS

The Coupland surgical aspirator is used for precision suction in the oral cavity to remove blood, saliva, and fragments of tissue and bone. The Coupland aspirator has 2 parts - a large handle (ASPCH) and separate tips of varying diameters. The large handle is comfortable to hold during procedures and the separate tips are easily interchangeable during surgery. Made of surgical grade stainless steel.

Handle

ASPCH	6-1/2"	each	\$144.99	\$144.99	\$130.49	\$130.49
-------	--------	------	----------	----------	----------	----------

*9-1/2" with tip attached

Tips

ASPCT1	1.5mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09
ASPCT2	2.0mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09
ASPCT3	2.5mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09
ASPCT4	3.0mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09

Tips with Holes

ASPCT3W	2.5mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09
ASPCT4W	3.0mm Tip diameter	each	\$28.99	\$28.99	\$26.09	\$26.09

SUCTION/IRRIGATION

HU-FRIEDY COGSWELL 3 ASPIRATOR

The Cogswell surgical aspirator is used for precision suction in the oral cavity to remove blood, saliva, and fragments of tissue and bone. The Cogswell aspirator is designed as a one-piece suction tube with a finger hole. The hole is located at the concave area midway along the tube. The intensity of the suction is increased when a finger covers the hole. The aspirator has a 2.3mm diameter opening and is made of surgical steel.

ASPCG3	9"	each	\$101.79	\$101.79	\$91.61	\$91.61
--------	----	------	----------	----------	---------	---------



QUALITY ASPIRATORS SURGICAL ASPIRATORS

Surgical aspirator made of 3/16" diameter stainless steel. Fits directly into any 11mm high volume cut off valve. Tapered tip helps prevent clogging. Fits HVE cut-off valve.

15P1A	A 1.5mm Opening	each	\$54.29	\$54.29	\$48.86	\$48.86
15P2A	2.5mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76
15P3A	3.0mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76

Surgical aspirator made of 3/16" diameter stainless steel. Fits directly into any 11mm high volume cut off valve. Tissue relief slots help prevent grabbing of the tongue, cheek and soft tissue. Fits HVE cut-off valve.

16RSA	B 4.0mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76
-------	------------------------	------	---------	---------	---------	---------

Surgical aspirator made of 1/4" diameter stainless steel. Larger tube size allows for more volume of suction. Fits directly into any 11mm high volume cut off valve. Tapered tip helps prevent clogging. Fits HVE cut-off valve.

46P1A	1.5mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76
46P2A	C 2.5mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76
46P3A	3.0mm Opening	each	\$35.29	\$35.29	\$31.76	\$31.76



ACE STAINLESS STEEL REUSABLE SUCTION TIPS

Oral Surgery

0900040	A Suction handle	6" (15cm)	each	\$59.29	\$56.33	\$56.33	\$53.36
0900041	B Suction tip, size 0, 1mm	4" (10cm)	each	\$20.79	\$19.75	\$19.75	\$18.71
0900042	C Suction tip, size 1, 2mm	4" (10cm)	each	\$20.79	\$20.79	\$18.71	\$18.71
0900043	D Suction tip, size 2, 3mm	4" (10cm)	each	\$20.79	\$20.79	\$18.71	\$18.71
0900044	E Suction tonsil tip, size 3	5-3/4" (14.5cm)	each	\$20.79	\$20.79	\$18.71	\$18.71

Coupland (Screw Type)

0900030	G Suction handle	6 1/2" (16.5cm)	each	\$61.79	\$58.70	\$58.70	\$55.61
0900032	H Suction tip, size 2, 2mm	4 1/2" (11.5cm)	each	\$31.79	\$31.79	\$28.61	\$28.61
0900033	I Suction tip, size 3, 3mm	4 1/2" (11.5cm)	each	\$31.79	\$31.79	\$28.61	\$28.61
0900038	J Replacement nut	1" (2.5cm)	each	\$21.49	\$20.42	\$20.42	\$19.34

Pelton Crane

6834301	K Suction tip, 3mm, w/ flow control	6-1/4" (16cm)	each	\$39.49	\$39.49	\$35.54	\$35.54
6834302	L Suction tip, 6mm, w/o flow control	5-1/2" (14cm)	each	\$39.49	\$39.49	\$35.54	\$35.54

Frazier with Stylet, 30° Angle

0900060	Suction tip, 6 fr	7" (18cm)	each	\$48.79	\$48.79	\$43.91	\$43.91
0900062	M Suction tip, 8 fr	7" (18cm)	each	\$48.79	\$48.79	\$43.91	\$43.91
0900063	Suction tip, 9 fr	7" (18cm)	each	\$48.79	\$48.79	\$43.91	\$43.91
0900064	Suction tip, 10 fr	7" (18cm)	each	\$48.79	\$48.79	\$43.91	\$43.91
0900065	Suction tip, 11 fr	7" (18cm)	each	\$48.79	\$46.35	\$46.35	\$43.91
0900066	Suction tip, 12 fr	7" (18cm)	each	\$51.29	\$51.29	\$46.16	\$46.16
0900070	N Replacement stylet	7-1/4" (18.5cm)	each	\$14.49	\$14.49	\$13.04	\$13.04
0900071	Replacement stylet with bulb	7-1/4" (18.5cm)	each	\$14.49	\$13.77	\$13.77	\$13.04



13

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

SUCTION/IRRIGATION



INTEGRA MILTEX FRAZIER FERGUSON SUCTION TUBES

With Finger Cut-Off, Angled


19570	2mm/6 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19572	2.3mm/7 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19573	2.6mm/8 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19574	3mm/9 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19575	3.3mm/10 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19576	3.6mm/11 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89
19577	4mm/12 Fr	each	\$70.99	\$70.99	\$63.89	\$63.89



CONMED FRAZIER METAL SUCTION TIPS

- 9-1/2"
- Sterile


0033100	10 Fr Olive tip	50/case	\$137.99	\$131.09	\$131.09	\$124.19
0033120	12 Fr with control vent	10/box	\$35.99	\$35.99	\$32.39	\$32.39



ACE AUTOCLAVABLE CLEANING BRUSHES FOR SUCTION TUBES

- Ideal for cleaning stainless steel cannulae, suction tubes, endoscopy scopes
- For use with or without detergent

0916117	Small	3mm bristle diameter	12" L	12/pkg	\$36.99	\$35.14	\$35.14	\$33.29
0917118	Medium	5.3mm bristle diameter	12" L	12/pkg	\$36.99	\$35.14	\$35.14	\$33.29
0918119	Large	11mm bristle diameter	12" L	12/pkg	\$38.29	\$36.38	\$36.38	\$34.46
0916118	Small	3mm bristle diameter	24" L	12/pkg	\$36.99	\$35.14	\$35.14	\$33.29
0917119	Medium	5.3mm bristle diameter	24" L	12/pkg	\$36.99	\$35.14	\$35.14	\$33.29
0918120	Large	11mm bristle diameter	24" L	12/pkg	\$38.29	\$36.38	\$36.38	\$34.46



QUALITY ASPIRATORS ASPIRATOR BRUSHES

Specially designed for aspirator tips and saliva ejectors, these brushes come in various sizes to accommodate a wide range of tubular diameters. The white nylon bristles clean blood and debris from inside walls. The stainless steel wire prevents rusting. These brushes will extend the full length of aspirator tips and are flexible enough to round 45 and 90 degree curvatures.

XSCB	Extra small	3/32" x 12"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76
SCB	Small	1/8" x 12"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76
MCB	Medium	1/4" x 12"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76
LCB	Large	1/2" x 12"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76
XLCB	Extra large	3/4" x 8"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76
SECB	Surgical	3/8" x 6 1/2"	12/pkg	\$35.29	\$33.53	\$33.53	\$31.76

SUCTION/IRRIGATION


ACE LATEX RUBBER TUBING

- 50' dispensing wheel
- Seamless wall tubing
- Tack-free outer surface
- Talc and powder-free surface
- Available in amber and black
- May be steam autoclaved (recommended number of times may vary)
- May not be one continuous roll

1/16" Inside Diameter (Bore)		I.D. x O.D.					
0032940	A Amber tubing	1/16" x 1/32"	50' roll	\$24.99	\$23.74	\$23.74	\$22.49
1/8" Inside Diameter (Bore)		I.D. x O.D.					
0032942	B Amber tubing	1/8" x 1/32"	50' roll	\$24.99	\$23.74	\$23.74	\$22.49
0032944	C Amber tubing	1/8" x 3/64"	50' roll	\$29.49	\$28.02	\$28.02	\$26.54
5/32" Inside Diameter (Bore)		I.D. x O.D.					
0032948	D Amber tubing	5/32" x 3/64"	50' roll	\$30.49	\$28.97	\$28.97	\$27.44
3/16" Inside Diameter (Bore)		I.D. x O.D.					
0032950	E Amber tubing	3/16" x 1/16"	50' roll	\$38.29	\$36.38	\$36.38	\$34.46
0032951	E Black tubing	3/16" x 1/16"	50' roll	\$38.29	\$36.38	\$36.38	\$34.46
0032952	F Amber tubing	3/16" x 3/32"	50' roll	\$43.49	\$41.32	\$41.32	\$39.14
0032953	F Black tubing	3/16" x 3/32"	50' roll	\$44.79	\$42.55	\$42.55	\$40.31
0032954	G Amber tubing	3/16" x 1/8"	50' roll	\$60.49	\$57.47	\$57.47	\$54.44
0032955	G Black tubing	3/16" x 1/8"	50' roll	\$68.99	\$65.54	\$65.54	\$62.09
1/4" Inside Diameter (Bore)		I.D. x O.D.					
0032958	H Amber tubing	1/4" x 1/16"	50' roll	\$41.49	\$39.42	\$39.42	\$37.34
0032959	H Black tubing	1/4" x 1/16"	50' roll	\$41.29	\$39.23	\$39.23	\$37.16
0032960	I Amber tubing	1/4" x 3/32"	50' roll	\$50.99	\$48.44	\$48.44	\$45.89
0032961	I Black tubing	1/4" x 3/32"	50' roll	\$50.99	\$48.44	\$48.44	\$45.89
0032962	J Amber tubing	1/4" x 1/8"	50' roll	\$67.99	\$64.59	\$64.59	\$61.19
0032963	J Black tubing	1/4" x 1/8"	50' roll	\$67.99	\$64.59	\$64.59	\$61.19
5/16" Inside Diameter (Bore)		I.D. x O.D.					
0032964	K Amber tubing	5/16" x 1/16"	50' roll	\$41.49	\$39.42	\$39.42	\$37.34
0032965	K Black tubing	5/16" x 1/16"	50' roll	\$41.49	\$39.42	\$39.42	\$37.34
0032966	L Amber tubing	5/16" x 3/32"	50' roll	\$58.29	\$55.38	\$55.38	\$52.46



I.D. = Inside diameter • O.D. = Outside diameter



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price

SUCTION/IRRIGATION



ACE SUCTION CONNECTING TUBING

- Special female connectors will fit all suction instruments
- Economical
- Excellent quality
- Single use
- Latex-free
- Non conductive

Sterile - Clear with Blue Tint

0034238	3' long (0.9m)	1/4" I.D. (6.4mm)	20/case	\$60.49	\$57.47	\$57.47	\$54.44
0034239	6' long (1.8m)	1/4" I.D. (6.4mm)	20/case	\$82.29	\$78.18	\$78.18	\$74.06
0034240	10' long (3.0m)	1/4" I.D. (6.4mm)	20/case	\$103.29	\$98.13	\$98.13	\$92.96

Sterile - Opaque Blue with Adapter**

- Hides fluid from patient

0034241	6' long (1.8m)	1/4" I.D. (6.4mm)	20/case	\$82.29	\$78.18	\$78.18	\$74.06
0034242	10' long (3.0m)	1/4" I.D. (6.4mm)	20/case	\$111.79	\$106.20	\$106.20	\$100.61

Non-Sterile - Clear with Blue Tint

0034235	3' long (0.9m)	1/4" I.D. (6.4mm)	50/case	\$75.79	\$72.00	\$72.00	\$68.21
0034236	6' long (1.8m)	1/4" I.D. (6.4mm)	25/case	\$62.99	\$59.84	\$59.84	\$56.69
0034237	10' long (3.0m)	1/4" I.D. (6.4mm)	25/case	\$85.79	\$81.50	\$81.50	\$77.21

Non-Sterile - Opaque Blue with Adapter**

- Hides fluid from patient

0034244	6' long (1.8m)	1/4" I.D. (6.4mm)	25/case	\$73.79	\$70.10	\$70.10	\$66.41
0034245	10' long (3.0m)	1/4" I.D. (6.4mm)	25/case	\$102.29	\$97.18	\$97.18	\$92.06

***Adapter allows connecting of dental saliva evacuation unit to suction tubing*

CONMED® SUCTION CONNECTING TUBING

Pre-Cut Suction Tubing

- Non-conductive, molded female connector
- Single-use
- Sterile
- Latex-free



3628	6' long (1.8m)	1/4" I.D., 3/8" O.D.	50/case	\$84.79	\$80.55	\$80.55	\$76.31
3629	10' long (3.0m)	1/4" I.D., 3/8" O.D.	50/case	\$107.29	\$101.93	\$101.93	\$96.56
3655	12' long (3.0m)	1/4" I.D., 3/8" O.D.	20/case	\$64.29	\$61.08	\$61.08	\$57.86
3648	6' long (1.8m)	3/16" I.D., 5/16" O.D.	50/case	\$71.79	\$68.20	\$68.20	\$64.61
3677	10' long (3.0m)	3/16" I.D., 5/16" O.D.	50/case	\$104.29	\$99.08	\$99.08	\$93.86
3656	12' long (3.0m)	3/16" I.D., 5/16" O.D.	20/case	\$59.79	\$56.80	\$56.80	\$53.81
3658	6' long (1.8m)	3/16" I.D.*, 5/16" O.D.	50/case	\$97.29	\$92.43	\$92.43	\$87.56

**With male and female connectors*

100' Roll Suction Tubing

- Non-conductive
- Economy roll tubing
- Clean, non-sterile
- 100 ft long w/integral blister connectors every 3 feet

0048250		3/16" I.D., 5/16" O.D.	100' roll/case	\$13.99	\$13.29	\$13.29	\$12.59
0048260		1/4" I.D., 3/8" O.D.	100' roll/case	\$22.29	\$21.18	\$21.18	\$20.06
0048270		9/32" I.D., 9/16" O.D.	100' roll/case	\$30.79	\$29.25	\$29.25	\$27.71

SUCTION/IRRIGATION

PRIMELINE AMBER SUCTION TUBING

- Made of pure natural latex rubber; seamless, flexible and smooth
- Excellent tensile strength
- Translucent amber color; Individually boxed 50 foot reels
- Non-sterile

035RA	I.D. Wall - 1/4" x 3/32", O.D. - 7/16"	50' roll/box	\$38.29	\$36.38	\$36.38	\$34.46
036RA	I.D. Wall - 1/4" x 1/8", O.D. - 1/2"	50' roll/box	\$50.99	\$48.44	\$48.44	\$45.89

DYNAREX SUCTION TUBING

A cost-effective solution with the quality and reliability you need to remove fluids from the surgical field.

- Sterile, non-conductive, latex-free, durable and flexible tubing with straw connector

DY32104	3/16" x 6'	50/case	\$54.79	\$52.05	\$52.05	\$49.31
32105	3/16" x 10'	50/case	\$79.99	\$75.99	\$75.99	\$71.99
32107	3/16" x 20'	20/case	\$41.49	\$39.42	\$39.42	\$37.34
32108	1/4" x 6'	50/case	\$67.79	\$64.40	\$64.40	\$61.01
32109	1/4" x 10'	50/case	\$97.79	\$92.90	\$92.90	\$88.01
32110	1/4" x 12'	20/case	\$38.29	\$36.38	\$36.38	\$34.46
32111	1/4" x 20'	20/case	\$57.29	\$54.43	\$54.43	\$51.56



CARDINAL HEALTH™ MEDI-VAC® SUCTION TUBING

- Single-use, sterile
- Maxi-Grip™ connectors
- Rigid male connector
- Latex-free


N512	3/16" x 12'	30/case	\$40.29	\$38.28	\$38.28	\$36.26
------	-------------	---------	----------------	----------------	----------------	----------------



DCI INTERNATIONAL CENTRAL VACUUM LEVER VALVE

- Includes cleaning brush

DCI5150	5-3/4"	each	\$89.79	\$85.30	\$85.30	\$80.81
---------	--------	------	----------------	----------------	----------------	----------------



Tier Pricing Details see page 2

Regular Price Member Price Premier Price Premier+ Price


SUCTION/IRRIGATION



COVIDIEN ARGYLE™ SUCTION CATHETERS WITH CHIMNEY VALVE

- Chimney style valve has a large opening on the control vent that minimizes residual suction
- Catheters are available in two tip styles - the smooth atraumatic DeLee tip (5 Fr - 8 Fr) with staggered eyes and the angled Whistle tip (10 Fr - 18 Fr)
- Latex-free

30500	5 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
30888	8 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31000	10 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31200	12 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31400	14 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31600	16 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16
31800	18 Fr	each	\$1.29	\$1.23	\$1.23	\$1.16



ACE ADAPTERS FOR HVAC SUCTION TO SUCTION TUBING

Brass

6592968	A O.D. = 1/4" x 7/16"	each	\$28.99	\$27.54	\$27.54	\$26.09
0032973	B O.D. = 5/16" x 7/16"	each	\$30.49	\$28.97	\$28.97	\$27.44

Stainless Steel

0033008	C O.D. = 5/16" x 3/8"	each	\$27.49	\$27.49	\$23.37	\$23.37
---------	------------------------------	------	---------	---------	---------	---------

ACE SUCTION TUBING ADAPTER SET

- Allows conversion from high-volume to low volume suction
- Surgical tubing is flexible and lightweight
- Ideal for long periodontal procedures
- 1/4" x 1/8" amber tubing, 5 ft. long

0900210	Suction tubing and adapter set	each	\$56.29	\$53.48	\$53.48	\$50.66
0900211	Tubing adapter, 5mm head	each	\$30.49	\$30.49	\$27.44	\$27.44

HU-FRIEDY ASPIRATOR ADAPTER WITH SILICONE TUBING

To connect Hu-Friedy aspirators with high-volume evacuation systems

PTAA	Aspirator adapter with 2" silicone tubing	each	\$28.99	\$27.54	\$27.54	\$26.09
PTRTC	Silicone replacement tube	2/pack	\$11.99	\$11.99	\$10.79	\$10.79

QUALITY ASPIRATORS TUBING ADAPTER SET

The tubing adapter set includes a stainless steel adapter and five feet of surgical tubing which allows you to connect 1/4" surgical tubing to your existing high volume cut-off valve.

Adapter and 5' Autoclavable Silicone Tubing

TASS	A	set	\$37.79	\$35.90	\$35.90	\$34.01
------	----------	-----	---------	---------	---------	---------

Adapter Only

TA	B	each	\$24.29	\$23.08	\$23.08	\$21.86
----	----------	------	---------	---------	---------	---------

Autoclavable Silicone Surgical Tubing

SST	C 1/4" ID x 3/32" Wall thickness	order by the foot	\$3.79	\$3.60	\$3.60	\$3.41
-----	---	-------------------	--------	--------	--------	--------

SUCTION/IRRIGATION

PASCAL ORATIP ADAPTER

- White lexan adapter
- Molded to fit the standard chuck on the hose end of most oral evacuator systems
- 4-¼" long and tapered
- Bent 15° for maximum utility
- Sterilizable in any manner
- Durable



24080		each	\$11.29	\$10.73	\$10.73	\$10.16
-------	--	------	---------	---------	---------	---------

ACE STRAIGHT SUCTION TUBING ADAPTER

- Plastic
- Non-sterile



2844000	Blue, 5"	50/bag	\$19.79	\$18.80	\$18.80	\$17.81
---------	----------	--------	---------	---------	---------	---------

PINNACLE EVAC-U-TRAP®

The Evac-u-Trap is a totally disposable canister, lid, mesh filter and gasket designed to replace the present canister and filter in your central vacuum pump. To remove this canister, screw on the lid and segregate with your contact amalgam for recycling. The infectious material is sealed in the canister where it cannot spill or contaminate the workplace. Nothing needs to be cleaned and you can be assured that you have minimized all contact with infectious material.

Features:

- Identical performance
- Constructed with plastic and stainless steel
- Saves staff time spent cleaning filters

Fits 1, 2 hp Den-Tal-Ez single models, Apollo 1 hp, 1.5 single: tandem models, Air Techniques-Vacstar models. Ohmeda-Matrix all models. Turbine Industries, Techwest-Whirl Vac models, Vacudent-Workhorse single, plus a number of other units.



PIN2200	2-3/4" x 3-5/8"	12/box	\$129.49	\$123.02	\$123.02	\$116.54
---------	-----------------	--------	----------	----------	----------	----------

Fits Den-Tal-Ez dual models made after 1978. Fits Dentsply made after 1991.

PIN2300	3-1/2" x 4-3/8"	8/box	\$129.49	\$123.02	\$123.02	\$116.54
---------	-----------------	-------	----------	----------	----------	----------

Fits Apollo 2 hp, 3 hp single and tandem models, Air Techniques - 1997 and later dual models, Vacudent - Workhorse dual. New finer mesh screen.

PIN2350FS	3-1/2" x 4-3/8"	8/box	\$129.49	\$123.02	\$123.02	\$116.54
-----------	-----------------	-------	----------	----------	----------	----------

SUCTION/IRRIGATION



PINNACLE DISPOS-A-TRAP®

The Dispos-a-Trap filter is a totally disposable filter designed to replace the present trap filter in your evacuation unit. It improves the performance of your suction system and reduces infectious aerosol. It does all this while eliminating one of the most hazardous and unpleasant tasks in the dental office, cleaning the waste that is present in your filter.

Features:

- Eliminates hazardous, unpleasant task of cleaning infectious waste in filters
- Improves suction performance
- Saves staff time spent cleaning filters

Fits A-dec (with black plastic canisters), Ampco, Belmont, DCI, Forest, Marus, MDT, and Schein

5501	2-1/8"	144/box	\$75.49	\$71.72	\$71.72	\$67.94
------	--------	---------	---------	---------	---------	---------

Fits Ampco, Belmont, Biotec, Chayes, Dabi-Atlante, Dansereau, DCI, Den-Tal-Ez, Dentech, Forest, Knight, Pelton & Crane, Prisma and others

5503	1-7/8"	144/box	\$75.49	\$71.72	\$71.72	\$67.94
------	--------	---------	---------	---------	---------	---------



DENTSPLY SIRONA PUREVAC® SC EVACUATION SYSTEM CLEANER

- Purevac® SC Evacuation System Cleaner has a separator compatible formula with a near neutral pH specifically designed for systems with amalgam separators
- Super cleaning - designed to remove a wide range of dental debris
- Scent of citrus
- Dispenses easily with 1 oz. (30ml) pump (included)
- Non-foaming and biodegradable

0021132FG	2 Liter concentrate	each	\$74.29	\$70.58	\$70.58	\$66.86
0021135FG	5 Liter concentrate *	each	\$156.79	\$148.95	\$148.95	\$141.11

*A freight charge of \$5.50 will be added to the price of this product.

13



ENZYME INDUSTRIES SANI-TREET PLUS LIQUID EVACUATION CLEANER

- Keeps evacuation systems clean and free of organic waste odors
- Non-toxic
- Non-caustic – harmless to metal, rubber, and plastic
- Biodegradable
- Peppermint scent

4198	Gallon *	each	\$62.99	\$59.84	\$59.84	\$56.69
------	----------	------	---------	---------	---------	---------

*A freight charge of \$5.50 will be added to the price of this product.

SUCTION/IRRIGATION

PASCAL CITRIZYME™ CONCENTRATE EVACUATION CLEANER

- 300g canister makes 35 gallons

PCL15250	each	\$50.99	\$48.44	\$48.44	\$45.89
----------	------	----------------	----------------	----------------	----------------



INTEGRA MILTEX EZ-ZYME® ALL PURPOSE ENZYME CLEANER

Neutral pH multi-enzyme liquid concentrate ideal for instrument presoak, ultrasonic cleaning, laundry presoak, and evacuation system cleaner. Ideal for cleaning hard to remove proteinaceous (mucous & tissue) and lipid (fatty) materials from the surfaces of surgical/dental instruments, linen, water submersible equipment, or where mechanical scrubbing is impractical.

3/4 oz (each packet yields 1 gallon of solution)

3750	32 packets/box	\$77.79	\$73.90	\$73.90	\$70.01
------	----------------	----------------	----------------	----------------	----------------



METREX THE SOLIDIFIER™

Suction Canister Fluid Solidification System

- Solidifies liquid medical waste in collection vessels, such as canisters
- Helps comply with DOT regulation 49CFR173.196 (a)(4) for Infectious Substances
- Protects healthcare workers and promotes workplace safety
- Eliminates spilling, splashing and aerosolizing

SD1200	Solidifies up to 1200 cc	64 bottles/case	\$127.79	\$121.40	\$121.40	\$115.01
SD1500	Solidifies up to 1500 cc	64 bottles/case	\$133.99	\$127.29	\$127.29	\$120.59
SD3000	Solidifies up to 3000 cc	64 bottles/case	\$209.99	\$199.49	\$199.49	\$188.99



13

FORMING AN ALLIANCE **Rewarding Loyalty**



Introducing the ACE SOUTHERN Alliance Program (ASAP)

When you create an alliance with **ACE SOUTHERN** you have a committed partner in your success. Part of that commitment is to reward you for your loyalty with significant savings on every dollar you spend.

As an **ASAP Qualified Customer**, you will benefit from savings and special promotions across all the products **ACE SOUTHERN** offers to help your practice prosper.

	Premier \$15K/yr spend	Premier Plus* \$25K/year spend
All Rx, Supplies and Equipment	5% saving	10% saving
All Instruments	10% savings	10% savings
Bone Grafting Materials and Implants	15-20% savings	15-20% savings

Enroll for free today at www.acesouthern.com/asap * AAOMS/AAP members only

SYRINGE AND NEEDLE

EXEL SUPERJECT™ DENTAL NEEDLES

- Hub design lowers the chances of needle breakage by 500%
- RED DOT on the needle hub marks the needle bevel position, even when the needle point is not visible
- Self-tapping ribs placed on the outside of the hub and inside threads make for ease of use



26557	27g Long	100/box	\$11.29	\$10.73	\$10.73	\$10.16
26562	30g Long	100/box	\$11.29	\$10.73	\$10.73	\$10.16
26555	30g Short	100/box	\$11.29	\$10.73	\$10.73	\$10.16

MONOJECT™ DENTAL NEEDLE

- Sterile, disposable
- Rigid autoclavable polypropylene hard pack
- Tamper-evident heat stake
- Ultra sharp, tri-beveled, anti-coring lancet point
- Anti-coring butt end point design



Plastic Hub

400033	25g Long (1-1/4")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
400041	25g Short (3/4")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
400058	27g Long (1-1/4")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
400066	27g Short (3/4")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
400074	30g Short (3/4")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
400173	30g X-Short (1/2")	100/box	\$13.29	\$12.63	\$12.63	\$11.96

Metal Hub

401031	25g Long (1-3/8")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
401056	27g Long (1-3/8")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
401064	27g Short (1")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
401072	30g Short (1")	100/box	\$13.29	\$12.63	\$12.63	\$11.96
401171	30g X-Short (1/2")	100/box	\$13.29	\$12.63	\$12.63	\$11.96

SEPTODONT SEPTOJECT® DISPOSABLE DENTAL NEEDLES

Each Septoject triple-beveled needle is siliconized to help ensure a smooth motion during injection and maximize patient comfort. Unlike metal hub needles, Septoject's plastic hub will not damage syringe hubs over time. Septoject's plastic hubs are pre-threaded for easy, straight attachment onto the syringe.



N1252	Red, 25g long	100/box	\$14.49	\$13.77	\$13.77	\$13.04
N1271	Orange, 27g short	100/box	\$14.49	\$13.77	\$13.77	\$13.04
N1272	Yellow, 27g long	100/box	\$14.49	\$13.77	\$13.77	\$13.04
N1300	Purple, 30g extra short	100/box	\$14.49	\$13.77	\$13.77	\$13.04
N1301	Blue, 30g short	100/box	\$14.49	\$13.77	\$13.77	\$13.04

SEPTODONT SEPTOJECT® XL DENTAL NEEDLES

Eliminate the confusion of choosing a size and gauge with only two choices - Septoject XL Block or Infiltration



- Large Bore Technology - (43% wider than a standard needle) to reduce level of pressure during injection, thus leading to a more comfortable injection
- High-Grade Surgical Stainless Steel Tubing - Reduce risk of needle breakage
- Siliconized Triple Bevel Design - Ensures a smooth, easy insertion into tissue
- Bevel Indicator - Orient the needle quickly with bevel indicator

01N1500	Infiltration	100/box	\$16.79	\$15.95	\$15.95	\$15.11
01N1550	Block	100/box	\$16.79	\$15.95	\$15.95	\$15.11

COVIDIEN MONOJECT™ 471 ENDODONTIC IRRIGATION NEEDLES

- Rigid autoclavable polypropylene hard pack
- Color-coded cap quickly identifies correct needle gauge
- Flow control notched tip permits greater irrigation of tooth canal
- Translucent polypropylene hub fits tightly on syringes
- Sterile



1471232	23g x 1-1/4"	25/box	\$14.49	\$13.77	\$13.77	\$13.04
1471273	27g x 1-1/4"	25/box	\$14.49	\$13.77	\$13.77	\$13.04

BD™ FILTER NEEDLES

- 5 micron thin wall - effective filtration of glass, metal, rubber, fiber and other contaminants larger than 5 microns
- Sterile



With BD Nokor™ Point

305201	18g x 1-1/2"	100/box	\$42.29	\$40.18	\$40.18	\$38.06
305200	19g x 1-1/2"	100/box	\$35.79	\$34.00	\$34.00	\$32.21

COVIDIEN MONOJECT™ BLUNT CANNULAS

- Stainless steel cannula
- Aluminum luer lock hub
- Smooth rounded blunt tip
- Color-coded needle sheath
- Autoclavable in rigid pack
- Sterile



202314	15g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202322	16g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202330	17g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202348	18g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202355	19g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202363	20g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202371	21g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202389	22g x 1-1/2"	25/box	\$28.49	\$27.07	\$27.07	\$25.64
202397	23g x 1"	25/box	\$28.49	\$27.07	\$27.07	\$25.64

SYRINGE AND NEEDLE

EXEL ALUMINUM HUB BLUNT NEEDLES

- Sterile, single-use, individually packaged blunt needle
- Stainless steel laser welded, non-lubricated cannula
- Aluminium hub, luer lock design
- Hub crimp cannula attachment
- Color-coded for easy identification
- 2-piece polypropylene needle cap and sheath



26471	15g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26472	17g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26473	18g x 1"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26474	19g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89
26475	20g x 1-1/2"	25/box	\$20.99	\$19.94	\$19.94	\$18.89

EXEL HYPODERMIC NEEDLES



26419	18g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26438	19g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26420	18g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26417	20g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26418	20g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26414	21g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26416	21g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26411	22g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26412	22g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26408	23g x 1	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26403	25g x 5/8"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26405	25g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26406	25g x 1-1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26400	27g x 1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26427	27g x 1-1/4"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26437	30g x 1/2"	100/box	\$5.29	\$5.03	\$5.03	\$4.76
26439	30g x 1"	100/box	\$5.29	\$5.03	\$5.03	\$4.76

COVIDIEN MONOJECT™ SOFTPACK HYPODERMIC NEEDLES



818112	18g x 1-1/2"	100/box	\$10.49	\$9.97	\$9.97	\$9.44
820100	20g x 1"	100/box	\$10.49	\$9.97	\$9.97	\$9.44
827112	27g x 1-1/2"	100/box	\$10.49	\$9.97	\$9.97	\$9.44

COVIDIEN MONOJECT™ RIGID PACK HYPODERMIC NEEDLES

250040	18g x 1"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250016	18g x 1-1/2"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250081	19g x 1"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250172	21g x 1"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250222	22g x 1"	100/box	\$13.49	\$12.82	\$12.17	\$12.14
250206	22g x 1-1/2"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250255	23g x 1"	100/box	\$13.49	\$12.82	\$12.17	\$12.14
250545	25g x 1-1/2"	100/box	\$13.49	\$12.82	\$12.82	\$12.14

Short Bevel

250024	18g x 1-1/2"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
250230	22g x 1"	100/box	\$13.49	\$12.82	\$12.82	\$12.14

Aluminum Hub

200433	25g x 1-1/4"	100/box	\$16.49	\$15.67	\$15.67	\$14.84
--------	--------------	---------	---------	---------	---------	---------

BD™ HYPODERMIC NEEDLES



305197	16g x 1"	100/box	\$25.49	\$24.22	\$24.22	\$22.94
305198	16g x 1-1/2"	100/box	\$25.49	\$24.22	\$24.22	\$22.94
305195	18g x 1"	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305196	18g x 1-1/2"	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305186	19g x 1" TW	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305187	19g x 1-1/2" TW	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305175	20g x 1"	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305176	20g x 1-1/2"	100/box	\$8.99	\$8.54	\$8.54	\$8.09
305165	21g x 1"	100/box	\$12.29	\$11.68	\$11.68	\$11.06
305167	21g x 1-1/2"	100/box	\$12.29	\$11.68	\$11.68	\$11.06
305129	21g x 2"	100/box	\$29.29	\$27.83	\$27.83	\$26.36
305155	22g x 1"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305156	22g x 1-1/2"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305145	23g x 1"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305122	25g x 5/8"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305124	25g x 7/8"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305125	25g x 1"	100/box	\$11.99	\$11.39	\$11.39	\$10.79
305127	25g x 1-1/2"	100/box	\$12.29	\$11.68	\$11.68	\$11.06
305109	27g x 1/2"	100/box	\$13.49	\$12.82	\$12.82	\$12.14
305136	27g x 1-1/4"	100/box	\$16.29	\$15.48	\$15.48	\$14.66
305106	30g x 1/2"	100/box	\$36.79	\$34.95	\$34.95	\$33.11
305128	30g x 1"	100/box	\$37.99	\$36.09	\$36.09	\$34.19

Short Bevel

305199	18g x 1-1/2"	100/box	\$11.49	\$10.92	\$10.92	\$10.34
305159	22g x 1-1/2"	1000/cs	\$113.29	\$107.63	\$107.63	\$101.96

TW - Thin Wall

SYRINGE AND NEEDLE

BD™ SPINAL NEEDLES

- Sterile
- Single-use
- Latex-free

405184	18g x 3-1/2", Pink	25/box	\$94.49	\$89.77	\$89.77	\$85.04
405181	22g x 3-1/2", Black	25/box	\$94.49	\$89.77	\$89.77	\$85.04

SEPTODONT ASPIRATING SYRINGES

Septodont dental aspirating syringes are lightweight and made of surgical grade aluminium and stainless steel. The syringes are durable and contain no removable parts and accept 1.7/1.8ml dental cartridges.

Septodont Standard is equipped with a larger thumb ring and offers the standard length plunger for improved comfort.

Septodont Petite is designed for smaller hands. Its compact design provides a smaller thumb ring, finger grip and shorter plunger length. This provides better control and added comfort when aspirating.



N2110	Standard, bronze	each	\$68.49	\$68.49	\$61.64	\$61.64
N2100	Petite, blue	each	\$66.49	\$66.49	\$59.84	\$59.84
01N2120	Petite, pink	each	\$66.49	\$66.49	\$59.84	\$59.84
01N2105	Fusion, titanium	each	\$66.49	\$66.49	\$59.84	\$59.84
01N2020	Saddle grip	each	\$139.49	\$139.49	\$125.54	\$125.54
N2030	Thumb grip	each	\$132.99	\$132.99	\$119.69	\$119.69

ACE COOK-WAITE TYPE ASPIRATING SYRINGES

- Accepts 1.8ml standard dental cartridges

6000001	5-1/2" (14cm)	each	\$32.29	\$32.29	\$29.06	\$29.06
---------	---------------	------	---------	---------	---------	---------



INTEGRA MILTEX ASPIRATING SYRINGES

Cook-Waite Type, Standard

7670	each	\$44.79	\$44.79	\$40.31	\$40.31
------	------	---------	---------	---------	---------

Self-Aspirating

76100	each	\$68.79	\$68.79	\$61.91	\$61.91
-------	------	---------	---------	---------	---------

Astra Type

7680	each	\$38.29	\$38.29	\$34.46	\$34.46
------	------	---------	---------	---------	---------



VISTA DENTAL ASPIRATING SYRINGES



Cook-Waite Type, Standard

107002	each	\$27.29	\$27.29	\$24.56	\$24.56
--------	------	---------	---------	---------	---------

Astra Type

107002	each	\$27.29	\$27.29	\$24.56	\$24.56
--------	------	---------	---------	---------	---------

ACE IRRIGATING SYRINGES

Disposable with Tapered Curved Tip

- Improved design for complete volume control
- Plastic



0034300	12 cc	50/box	\$23.49	\$22.32	\$22.32	\$21.14
0034302	20 cc	25/box	\$23.49	\$23.49	\$21.14	\$21.14

EXEL CURVED TIP SYRINGE

- Ungraduated with plunger stop
- Ideal for pressure irrigation



26267	12 cc Irrigation	50/box	\$20.79	\$19.75	\$19.75	\$18.71
-------	------------------	--------	---------	---------	---------	---------

COVIDIEN MONOJECT™ CURVED TIP SYRINGE

- Ungraduated
- Ideal for pressure irrigation



412012	12 cc Irrigation	50/box	\$26.29	\$24.98	\$24.98	\$23.66
--------	------------------	--------	---------	---------	---------	---------

COVIDIEN MONOJECT™ CLEAR ORAL MEDICATION SYRINGE

- With permanent, heat-etched dual graduations



907102	10 cc	100/box	\$51.29	\$48.73	\$48.73	\$46.16
--------	-------	---------	---------	---------	---------	---------

MEDLINE BULB IRRIGATION SYRINGE

- Easy to grip and use with only one hand, the control-bulb is designed to stand on end, making it excellent for use in clinical settings, sterile



DYND20125	each	\$1.79	\$1.70	\$1.70	\$1.61
-----------	------	--------	--------	--------	--------

EXEL SYRINGE WITH REDUCED DEAD SPACE

- Reduced dead space for reduced medication waste
- Sterile



Lock Tip w/Cap

26049	1 cc	100/box	\$24.79	\$23.55	\$23.55	\$22.31
-------	------	---------	---------	---------	---------	---------

EXEL SYRINGE ONLY



Luer Slip Tip

26048	1 cc tb*	100/box	\$13.49	\$12.82	\$12.82	\$12.14
26201	3 cc*	100/box	\$9.29	\$8.83	\$8.83	\$8.36
26231	5-6 cc*	100/box	\$17.99	\$17.09	\$17.09	\$16.19
26266	10-12 cc*	100/box	\$18.49	\$17.57	\$17.57	\$16.64
26281	20-25 cc	50/box	\$21.99	\$20.89	\$20.89	\$19.79
26301	50-60 cc	25/box	\$17.79	\$16.90	\$16.90	\$16.01

Luer Lock Tip

26200	3 cc*	100/box	\$8.49	\$8.07	\$8.07	\$7.64
26230	5-6 cc*	100/box	\$12.99	\$12.34	\$12.34	\$11.69
26265	10-12 cc*	100/box	\$17.29	\$16.43	\$16.43	\$15.56
26280	20-25 cc*	50/box	\$28.49	\$27.07	\$27.07	\$25.64
26290	30-35 cc	50/box	\$25.29	\$24.03	\$24.03	\$22.76
26300	50-60 cc*	25/box	\$17.79	\$16.90	\$16.90	\$16.01

Luer Slip Eccentric Tip

26291	30-35 cc*	50/box	\$27.49	\$26.12	\$26.12	\$24.74
-------	-----------	--------	---------	---------	---------	---------

Catheter Tip

26304	50-60 cc*	25/box	\$17.79	\$16.90	\$16.90	\$16.01
-------	-----------	--------	---------	---------	---------	---------

*With syringe tip cap

TERUMO SYRINGE ONLY



Luer Lock Tip

60L	60 cc	25/box	\$16.79	\$15.95	\$15.95	\$15.11
-----	-------	--------	---------	---------	---------	---------

SYRINGE AND NEEDLE

COVIDIEN SOFTPACK SYRINGES



Luer Lock Tip

300777	3 cc	100/box	\$19.79	\$18.80	\$18.80	\$17.81
600777	6 cc	100/box	\$54.49	\$51.77	\$51.77	\$49.04
200777T	12 cc	80/box	\$28.49	\$27.07	\$27.07	\$25.64
500777	35 cc	40/box	\$35.99	\$34.19	\$34.19	\$32.39

COVIDIEN MONOJECT RIGID PACK SYRINGE



Regular Tip

501400	1 cc tb	100/box	\$22.29	\$21.18	\$21.18	\$20.06
512852	12 cc	80/box	\$28.79	\$27.35	\$27.35	\$25.91
520673	20 cc	50/box	\$44.29	\$42.08	\$42.08	\$39.86
560224	60 cc	20/box	\$43.79	\$41.60	\$41.60	\$39.41

Luer Lock Tip

513934	3 cc	100/box	\$14.99	\$14.24	\$14.24	\$13.49
516937	6 cc	50/box	\$21.99	\$20.89	\$20.89	\$19.79
512878	12 cc	100/box	\$26.49	\$25.17	\$25.17	\$23.84
520657	20 cc	80/box	\$44.29	\$42.08	\$42.08	\$39.86
535762	35 cc	50/box	\$24.99	\$23.74	\$23.74	\$22.49
560125	60 cc	20/box	\$23.99	\$22.79	\$22.79	\$21.59

BD™ SYRINGE ONLY



Slip Tip

309659	1 cc tb	200/box	\$31.29	\$29.73	\$29.73	\$28.16
309656	3 cc	200/box	\$24.49	\$23.27	\$23.27	\$22.04
309647	5 cc	125/box	\$29.29	\$27.83	\$27.83	\$26.36
303134	10 cc	200/box	\$51.49	\$48.92	\$48.92	\$46.34
302831	20 cc	48/box	\$34.49	\$32.77	\$32.77	\$31.04
302833	30 cc	56/box	\$41.49	\$39.42	\$39.42	\$37.34
309654	60 cc	40/box	\$43.29	\$41.13	\$41.13	\$38.96

Luer Lock Tip

309628	1 cc	100/box	\$49.29	\$46.83	\$46.83	\$44.36
309657	3 cc	200/box	\$20.79	\$19.75	\$19.75	\$18.71
302995	10 cc	200/box	\$37.29	\$35.43	\$35.43	\$33.56
302830	20 cc	48/box	\$20.99	\$19.94	\$19.94	\$18.89
302832	30 cc	56/box	\$24.99	\$23.74	\$23.74	\$22.49
309653	50 cc	40/box	\$35.99	\$34.19	\$34.19	\$32.39

Insulin, Luer Slip

329654	1 cc	200/box	\$50.99	\$48.44	\$48.44	\$45.89
--------	------	---------	---------	---------	---------	---------

SYRINGE AND NEEDLE

BD™ STERILE SYRINGE TIP CAP



- Latex-free polypropylene, sterile

308341	10/tray	\$1.49	\$1.42	\$1.42	\$1.34
305819	200/box	\$229.49	\$218.02	\$218.02	\$206.54

EXEL SYRINGE WITH NEEDLE COMBINATIONS



- 100/box

26108	3 cc w/20g x 1"	\$9.29	\$8.83	\$8.83	\$8.36
26109	3 cc w/20g x 1-1/2"	\$9.29	\$8.83	\$8.83	\$8.36
26105	3 cc w/21g x 1"	\$9.29	\$8.83	\$8.83	\$8.36
26107	3 cc w/21g x 1-1/2"	\$9.29	\$8.83	\$8.83	\$8.36
26102	3 cc w/22g x 1"	\$9.29	\$8.83	\$8.83	\$8.36
26104	3 cc w/22g x 1-1/2"	\$9.29	\$8.83	\$8.83	\$8.36
26101	3 cc w/23g x 1"	\$9.29	\$8.83	\$8.83	\$8.36
26100	3 cc w/25g x 5/8"	\$9.29	\$8.83	\$8.83	\$8.36
26110	3 cc w/18g x 1-1/2"	\$9.29	\$8.83	\$8.83	\$8.36
26111	3 cc w/25g x 1"	\$9.29	\$8.83	\$8.83	\$8.36
26112	3 cc w/25g x 1-1/2"	\$9.29	\$8.83	\$8.83	\$8.36
26214	5-6 cc w/20g x 1"	\$22.99	\$21.84	\$21.84	\$20.69
26212	5-6 cc w/21g x 1"	\$22.99	\$21.84	\$21.84	\$20.69
26213	5-6 cc w/21g x 1-1/2"	\$22.99	\$21.84	\$21.84	\$20.69
26210	5-6 cc w/22g x 1"	\$22.99	\$21.84	\$21.84	\$20.69
26211	5-6 cc w/22g x 1-1/2"	\$22.99	\$21.84	\$21.84	\$20.69
26254	10-12 cc w/20g x 1"	\$25.29	\$24.03	\$24.03	\$22.76
26255	10-12 cc w/20g x 1-1/2"	\$25.29	\$24.03	\$24.03	\$22.76
26252	10-12 cc w/21g x 1"	\$25.29	\$24.03	\$24.03	\$22.76
26253	10-12 cc w/21g x 1-1/2"	\$25.29	\$24.03	\$24.03	\$22.76
26250	10-12 cc w/22g x 1"	\$25.29	\$24.03	\$24.03	\$22.76
26251	10-12 cc w/22g x 1-1/2"	\$25.29	\$24.03	\$24.03	\$22.76
TB					
26044	1 cc w/25g x 5/8"	\$12.99	\$12.34	\$12.34	\$11.69
26043	1 cc w/26g x 1/2"	\$12.99	\$12.34	\$12.34	\$11.69
EXL26040	1 cc w/27g x 1/2"	\$12.99	\$12.34	\$12.34	\$11.69
Insulin					
26028	1/2 cc w/29g x 1/2"	\$14.79	\$14.05	\$14.05	\$13.31
26029	1 cc w/29g x 1/2"	\$15.29	\$14.53	\$14.53	\$13.76

TERUMO SYRINGE WITH NEEDLE



- Includes Caps
- Bold, precise scale markings
- Clear syringe barrel
- Translucent needle hub for quick observation of flashback.
- Maximum patient comfort delivered through ultimate sharpness/smoothness
- 100/box

Luer Lock Tip – Latex-free

SS03L2138	3 cc w/22g x 1-1/2"	\$13.49	\$12.82	\$12.82	\$12.14
-----------	---------------------	---------	---------	---------	---------

COVIDIEN MONOJECT™ SOFTPACK SYRINGE WITH NEEDLE COMBINATIONS



- 100/box

320100	3 cc w/20g x 1"	\$17.29	\$16.43	\$16.43	\$15.56
320112	3 cc w/20g x 1-1/2"	\$17.29	\$16.43	\$16.43	\$15.56
321100	3 cc w/21g x 1"	\$17.29	\$16.43	\$16.43	\$15.56
321112	3 cc w/21g x 1-1/2"	\$17.29	\$16.43	\$16.43	\$15.56
322100	3 cc w/22g x 1"	\$17.29	\$16.43	\$16.43	\$15.56
322112	3 cc w/22g x 1-1/2"	\$17.29	\$16.43	\$16.43	\$15.56
323100	3 cc w/23g x 1"	\$17.29	\$16.43	\$16.43	\$15.56
325058	3 cc w/25g x 5/8"	\$17.29	\$16.43	\$16.43	\$15.56

TB (Detachable Needle)

127012	1 cc w/27g x 1/2"	\$37.79	\$35.90	\$35.90	\$34.01
--------	-------------------	---------	---------	---------	---------

COVIDIEN MONOJECT™ RIGID PACK SYRINGE WITH NEEDLE COMBINATIONS



- Sterile
- Disposable
- Single-use

• 100/box					
513058	3 cc w/20g x 1-1/2"	\$17.49	\$16.62	\$16.62	\$15.74
513132	3 cc w/21g x 1"	\$17.49	\$16.62	\$16.62	\$15.74
513157	3 cc w/21g x 1-1/2"	\$17.49	\$16.62	\$16.62	\$15.74
513330	3 cc w/23g x 1"	\$17.49	\$16.62	\$16.62	\$15.74
• 50/box					
516135	6 cc w/21g x 1"	\$38.49	\$36.57	\$36.57	\$34.64

SYRINGE AND NEEDLE

BD™ SYRINGE WITH NEEDLE COMBINATIONS

• 100/box



309580	3 cc w/18g x 1-1/2"	\$15.49	\$14.72	\$14.72	\$13.94
309578	3 cc w/20g x 1"	\$15.49	\$14.72	\$14.72	\$13.94
309579	3 cc w/20g x 1-1/2"	\$15.49	\$14.72	\$14.72	\$13.94
309575	3 cc w/21g x 1"	\$17.49	\$16.62	\$16.62	\$15.74
309577	3 cc w/21g x 1-1/2"	\$18.79	\$17.85	\$17.85	\$16.91
309572	3 cc w/22g x 1"	\$18.79	\$17.85	\$17.85	\$16.91
309574	3 cc w/22g x 1-1/2"	\$18.79	\$17.85	\$17.85	\$16.91
309571	3 cc w/23g x 1"	\$18.79	\$17.85	\$17.85	\$16.91
309570	3 cc w/25g x 5/8"	\$18.79	\$17.85	\$17.85	\$16.91
309581	3 cc w/25g x 1"	\$18.79	\$17.85	\$17.85	\$16.91
309582	3 cc w/25g x 1-1/2"	\$18.79	\$17.85	\$17.85	\$16.91
309587	3 cc w/26g x 5/8" sub-q	\$18.79	\$17.85	\$17.85	\$16.91
309634	5 cc w/20g x 1"	\$27.49	\$26.12	\$26.12	\$24.74
309635	5 cc w/20g x 1-1/2"	\$27.49	\$26.12	\$26.12	\$24.74
309632	5 cc w/21g x 1"	\$32.29	\$30.68	\$30.68	\$29.06
309633	5 cc w/21g x 1-1/2"	\$34.79	\$33.05	\$33.05	\$31.31
309630	5 cc w/22g x 1"	\$34.79	\$33.05	\$33.05	\$31.31
309631	5 cc w/22g x 1-1/2"	\$34.79	\$33.05	\$33.05	\$31.31
309644	10 cc w/20g x 1"	\$36.29	\$34.48	\$34.48	\$32.66
309645	10 cc w/20g x 1-1/2"	\$36.29	\$34.48	\$34.48	\$32.66
309642	10 cc w/21g x 1"	\$38.29	\$36.38	\$36.38	\$34.46
309643	10 cc w/21g x 1-1/2"	\$38.29	\$36.38	\$36.38	\$34.46
309640	10 cc w/22g x 1"	\$38.29	\$36.38	\$36.38	\$34.46

Insulin

328438	3/10 cc w/31g x 5/16"	\$54.29	\$51.58	\$51.58	\$48.86
328466	1/2 cc w/30g x 1/2"	\$54.29	\$51.58	\$51.58	\$48.86
328418	1 cc w/31g x 5/16"	\$51.29	\$48.73	\$48.73	\$46.16
328411	1 cc w/30g x 1/2"	\$51.29	\$48.73	\$48.73	\$46.16

TB

309626	1 cc w/25g x 5/8"	\$23.49	\$22.32	\$22.32	\$21.14
309625	1 cc w/26g x 3/8"	\$23.49	\$22.32	\$22.32	\$21.14
309623	1 cc w/27g x 1/2"	\$23.49	\$22.32	\$22.32	\$21.14

TERUMO SURGUARD®3 SAFETY HYPODERMIC NEEDLE

- 3 modes of activation—thumb, finger, or hard surface
- Separate finger/thumb grips with ridges and nonslip surfaces
- Hub and needle locks for security post activation
- Firm injection feel
- Sharpest needles among major brands tested*

Needle Only

SG32025	20g x 1"	100/box	\$33.29	\$31.63	\$31.63	\$29.96
SG32038	20g x 1-1/2"	100/box	\$33.29	\$31.63	\$31.63	\$29.96

EXEL SECURETOUCH® SAFETY NEEDLE

Safety Hypodermic Needle

- One-handed activation
- No solid surface required
- Simple finger stroke activation
- Protective arm keeps fingers behind the needle
- An audible "click" as well as tactile feel and visual indicator confirm safety mechanism engagement
- Needle remains locked inside activated cover
- Fits any brand syringe
- This device does not contain natural rubber latex



Safety shield can be pushed back comparative to the needle, allowing a clear view of the injection site

Needle Only

27420	18g x 1-1/2"	100/box	\$28.79	\$27.35	\$27.35	\$25.91
27427	21g x 1-1/4"	100/box	\$28.79	\$27.35	\$27.35	\$25.91
27406	25g x 1-1/2"	100/box	\$28.79	\$27.35	\$27.35	\$25.91
27426	27g x 1-1/2"	100/box	\$28.79	\$27.35	\$27.35	\$25.91

BD SAFETYGLIDE™ SHIELDING HYPODERMIC NEEDLES

Designed for efficiency, safety and patient comfort

- Bevel-up needle-orientation for performing low angle injections
- BD PrecisionGlide™ Needle Technology - Designed for patient comfort
- Activation-Assist™ Technology - For fast and easy needle tip shielding
- Detachable Needle - Fits on any size conventional syringe
- Stainless Steel Latch - Securely shields needle tip after activation



Needle Only

305915	21g x 1"	50/box	\$21.99	\$20.89	\$20.89	\$19.79
305902	23g x 1"	50/box	\$21.99	\$20.89	\$20.89	\$19.79

BD ECLIPSE™ SAFETY NEEDLES

BD Eclipse Injection Needle Offers:

- Single-handed activation
- No hard surface required for activation
- Bevel-oriented needle allows for low-angle injections
- BD PrecisionGlide™ needles minimize injection pain and patient discomfort

Needle Only

305763	22g x 1-1/2"	100/box	\$41.79	\$39.70	\$39.70	\$37.61
305762	23g x 1"	100/box	\$41.79	\$39.70	\$39.70	\$37.61
305761	25g x 1"	100/box	\$41.79	\$39.70	\$39.70	\$37.61
305757	30g x 1/2"	100/box	\$41.79	\$39.70	\$39.70	\$37.61



SYRINGE AND NEEDLE


BD™ SYRINGES W/ PREATTACHED BLUNT CANNULAS



For use with Interlink® or Prepierced System

303346	3 cc w/ Blunt cannula	100/box	\$39.99	\$37.99	\$37.99	\$35.99
303347	5 cc w/ Blunt cannula	100/box	\$44.29	\$42.08	\$42.08	\$39.86
303348	10 cc w/ Blunt cannula	100/box	\$49.29	\$46.83	\$46.83	\$44.36

EXEL SYRINGES WITH SECURETOUCH® SAFETY NEEDLES



Safety Needle Features

- One-handed activation
- No solid surface required
- Simple finger stroke activation
- Protective arm keeps fingers behind the needle
- An audible "click" as well as tactile feel and visual indicator confirm safety mechanism engagement
- Needle remains locked inside activated cover
- Luer lock syringe
- This device does not contain natural rubber latex
- 50/box

Safety shield can be pushed back comparative to the needle, allowing a clear view of the injection site

27044	1 cc w/25g x 5/8"		\$13.49	\$12.82	\$12.82	\$12.14
EXL27110	3 cc w/18g x 1-1/2"		\$12.49	\$11.87	\$11.87	\$11.24
27108	3 cc w/20g x 1"		\$12.49	\$11.87	\$11.87	\$11.24
27105	3 cc w/21g x 1"		\$12.49	\$11.87	\$11.87	\$11.24
27107	3 cc w/21g x 1-1/2"		\$12.49	\$11.87	\$11.87	\$11.24
27102	3 cc w/22g x 1"		\$12.49	\$11.87	\$11.87	\$11.24
27104	3 cc w/22g x 1-1/2"		\$12.49	\$11.87	\$11.87	\$11.24
27101	3 cc w/23g x 1"		\$12.49	\$11.87	\$11.87	\$11.24
27100	3 cc w/25g x 5/8"		\$12.49	\$11.87	\$11.87	\$11.24
27111	3 cc w/25g x 1"		\$12.49	\$11.87	\$11.87	\$11.24
27112	3 cc w/25g x 1-1/2"		\$12.49	\$11.87	\$11.87	\$11.24

TERUMO SURGUARD®3 SAFETY HYPODERMIC NEEDLES WITH SYRINGES

- Latex-Free
- 100/box

SG301T2516	1 cc w/25g x 5/8"		\$59.49	\$56.52	\$56.52	\$53.54
SG303L2125	3 cc w/21g x 1"		\$37.79	\$35.90	\$35.90	\$34.01
SG303L2516	3 cc w/25g x 5/8"		\$37.79	\$35.90	\$35.90	\$34.01

BD™ SYRINGES W/ SAFETYGLIDE™ SHIELDING NEEDLE



- A single hand, a single movement
- BD SafetyGlide needle for unequalled sharpness and penetration ease
- Patented Activation-Assist™ technology for fast and easy needletip shielding
- Bevel-up needle orientation for performing low angle injections
- Stainless steel latch securely shields needletip after activation

305909	3 cc w/21g x 1-1/2"	50/box	\$25.99	\$24.69	\$24.69	\$23.39
305906	3 cc w/22g x 1-1/2"	50/box	\$25.99	\$24.69	\$24.69	\$23.39
305905	3 cc w/23g x 1"	50/box	\$25.99	\$24.69	\$24.69	\$23.39
305904	3 cc w/25g x 5/8"	50/box	\$25.99	\$24.69	\$24.69	\$23.39
305908	10 cc w/22g x 1-1/2"	50/box	\$30.49	\$28.97	\$28.97	\$27.44

CERTOL INTERNATIONAL PRO TECTOR® NEEDLE SHEATH PROP

User Protection:

- Reduces the possibility of accidental needle sticks
- Protects during disassembly and disposal
- Discourages leaving the needle unsheathed on the tray

Autoclavable:

- For one-time sterile tray set-up

Holds Cap Securely:

- Specifically designed to secure cap and hold it in the ready position

Meets CDC Guidelines:

- CDC guidelines state "Use either a one-handed scoop technique or a mechanical device designed for holding the needle cap when recapping needles." Use of the Certol Pro Tector Needle Sheath Prop is one way to meet CDC guidelines for preventing exposures to bloodborne pathogens.



PNS100	100/box	\$19.99	\$18.99	\$18.99	\$17.99
--------	---------	--------------------	--------------------	--------------------	--------------------

HU-FRIEDY NEEDLE COVER HOLDER

- Stainless Steel



STKSKLS	Needle cover holder	\$40.49	\$40.49	\$36.44	\$36.44
---------	---------------------	--------------------	--------------------	--------------------	--------------------

SEPTODONT AIM SAFE AUTOCLAVABLE NEEDLE RECAPPING DEVICE



Aim Safe sits on any flat surface holding the needle cap in a position that allows for easy recapping when the injection is completed.

- Helps protect against needle-stick injuries
- Economical and safe

01N0015	5/box	\$31.79	\$30.20	\$30.20	\$28.61
---------	-------	--------------------	--------------------	--------------------	--------------------

SYRINGE AND NEEDLE

BD™ MULTI-USE SHARPS COLLECTORS

Convenient, vertical "point first" drop, these robust sharps collectors are preassembled and offer a patented one-way funnel valve to minimize needlesticks and overfilling.



Regular Funnel Entry

305488	3.3 Quart	each	\$6.79	\$6.45	\$6.45	\$6.11
305489	6.9 Quart	each	\$11.79	\$11.20	\$11.20	\$10.61
305490	8.2 Quart	each	\$12.49	\$11.87	\$11.87	\$11.24
5491	20 Quart	each	\$19.79	\$18.80	\$18.80	\$17.81

COVIDIEN SHARPSAFETY IN-ROOM™ SHARPS CONTAINERS WITH ALWAYS-OPEN LID

- Always-Open lids provide convenient, quick and easy point-of-use sharps disposal
- Accommodates larger syringes
- Economical horizontal drop maximizes container volume



Transparent Red

851301	5 Quart	each	\$5.79	\$5.50	\$5.50	\$5.21
85321R	2 Gallon	each	\$7.29	\$6.93	\$6.93	\$6.56

BD™ CLEAR TOP SHARPS COLLECTOR

Clear top provides a visualization of collector fill-levels. Increasing fill-level identification helps prevent collector overfill, a leading cause of disposal injuries.

- Dual access



305487	1.5 Quart	each	\$7.79	\$7.40	\$7.40	\$7.01
--------	-----------	------	--------	--------	--------	--------

COVIDIEN SHARPS CONTAINER WITH ROTOR OPENING LID

- Adjustable rotor opening lids accommodate a variety of sharps sizes and provide temporary and permanent closure
- Containers lock for final disposal
- Nestable containers save valuable storage space
- Plastics and pigments are safe to autoclave or incinerate



Red

8970	2 Gallon	each	\$7.29	\$6.93	\$6.93	\$6.56
------	----------	------	--------	--------	--------	--------

COVIDIEN SHARPS CONTAINERS

- Polished viewing windows on four sides
- Translucent colors with lead-free pigments containing no heavy metals
- Clear lid with extra-wide butterfly closure
- Needle removal notches
- Can be autoclaved or incinerated
- Chimney-top



Small

676236	4 Quart	each	\$4.29	\$4.08	\$4.08	\$3.86
--------	---------	------	--------	--------	--------	--------

Medium

676285	8 Quart	each	\$5.49	\$5.22	\$5.22	\$4.94
--------	---------	------	--------	--------	--------	--------

Large

676434	14 Quart	each	\$9.29	\$8.83	\$8.83	\$8.36
--------	----------	------	--------	--------	--------	--------

Mounting Bracket for Vertical Surfaces

8881676012		each	\$13.49	\$12.82	\$12.82	\$12.14
------------	--	------	---------	---------	---------	---------

COVIDIEN SHARPS-A-GATOR™ SHARPS CONTAINER, TORTUOUS PATH

- One-hand, one-step, no-touch sharps depositing
- Horizontal loading, non-mechanical sharps opening
- Tortuous-path design impedes hand entry
- Various mounting brackets and glove box holders are available



Transparent Red

31144010	5 Quart	each	\$4.79	\$4.55	\$4.55	\$4.31
----------	---------	------	--------	--------	--------	--------

UNIFORMS



Tunnel belt rear view

FASHION SEAL WOMEN'S CROSSOVER TUNICS WITH CONTRASTING TRIM

- Classic fit
- Set-in sleeves
- Side vents
- Two lower pockets with pencil divider and inset contrasting pockets
- Scissor holder in right pocket
- Drawstring back tie with coordinating webbed belting
- 65% polyester/35% cotton

Sizes XS-XL

7001	Burgundy w/black	each	\$22.49	\$21.37	\$21.37	\$20.24
7002	Ceil blue w/navy	each	\$22.49	\$21.37	\$21.37	\$20.24
7027	Cobalt w/green apple	each	\$22.49	\$21.37	\$21.37	\$20.24
7005	Dill w/navy	each	\$22.49	\$21.37	\$21.37	\$20.24
7026	Navy w/ceil blue	each	\$22.49	\$21.37	\$21.37	\$20.24
7147	Pewter w/black	each	\$22.49	\$21.37	\$21.37	\$20.24
7599	Pewter w/pretty pink	each	\$22.49	\$21.37	\$21.37	\$20.24
7006	Plum w/eggplant	each	\$22.49	\$21.37	\$21.37	\$20.24
7003	Sangria w/navy	each	\$22.49	\$21.37	\$21.37	\$20.24
7007	Tan w/chocolate	each	\$22.49	\$21.37	\$21.37	\$20.24

Sizes 2XL-4XL also available, call for pricing

FASHION SEAL WOMEN'S V-NECK TUNICS-MODERN FIT

- Set-in sleeves
- Side vents
- Reinforced, bartacked V-neck
- Two lower pockets
- Concealed pocket in left pocket
- Scissor holder in right pocket
- 65% polyester/35% cotton

Sizes XS-XL

7172	Black	each	\$20.29	\$19.28	\$19.28	\$18.26
7140	Ceil blue	each	\$20.29	\$19.28	\$19.28	\$18.26
6900	Citrus green	each	\$20.29	\$19.28	\$19.28	\$18.26
6902	Cocoa	each	\$20.29	\$19.28	\$19.28	\$18.26
6908	Fir green	each	\$20.29	\$19.28	\$19.28	\$18.26
6901	Flamingo	each	\$20.29	\$19.28	\$19.28	\$18.26
6905	Moss green	each	\$20.29	\$19.28	\$19.28	\$18.26
7142	Navy	each	\$20.29	\$19.28	\$19.28	\$18.26
6907	Pacific blue	each	\$20.29	\$19.28	\$19.28	\$18.26
7143	Peacock green	each	\$20.29	\$19.28	\$19.28	\$18.26
7144	Pewter	each	\$20.29	\$19.28	\$19.28	\$18.26
6903	Raisin	each	\$20.29	\$19.28	\$19.28	\$18.26
7145	Royal	each	\$20.29	\$19.28	\$19.28	\$18.26
6904	Sandstone	each	\$20.29	\$19.28	\$19.28	\$18.26
6906	True red	each	\$20.29	\$19.28	\$19.28	\$18.26

Sizes 2XL-4XL also available, call for pricing

15

*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.

Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

FASHION SEAL WOMEN'S V-NECK TUNICS-CLASSIC FIT

- Set-in sleeves
- Side vents
- Reinforced, bartacked v-neck
- Two lower pockets
- Scissor holder in right pocket
- 65% polyester/35% cotton

Sizes XS-XL

7012	Black	each	\$20.49	\$19.47	\$19.47	\$18.44
7281	Blueberry	each	\$20.49	\$19.47	\$19.47	\$18.44
7324	Burgundy	each	\$20.49	\$19.47	\$19.47	\$18.44
7347	Ceil blue	each	\$20.49	\$19.47	\$19.47	\$18.44
7284	Cobalt	each	\$20.49	\$19.47	\$19.47	\$18.44
7325	Dark teal	each	\$20.49	\$19.47	\$19.47	\$18.44
7326	Eggplant	each	\$20.49	\$19.47	\$19.47	\$18.44
7327	Fir green	each	\$20.49	\$19.47	\$19.47	\$18.44
7329	Navy	each	\$20.49	\$19.47	\$19.47	\$18.44
7239	Pewter	each	\$20.49	\$19.47	\$19.47	\$18.44
7016	Sangria	each	\$20.49	\$19.47	\$19.47	\$18.44
7236	Tan	each	\$20.49	\$19.47	\$19.47	\$18.44
7288	Teal	each	\$20.49	\$19.47	\$19.47	\$18.44
7343	White	each	\$20.49	\$19.47	\$19.47	\$18.44

Sizes 2XL-4XL also available, call for pricing



FASHION SEAL WOMEN'S SIDE FLEX TUNICS

- Modern fit
- Set-in sleeves
- Concealed front pocket
- Moisture wicking flexible side panels
- 65% polyester/35% cotton

Sizes XS-XL

7511	Burgundy w/black	each	\$20.29	\$19.28	\$19.28	\$18.26
7541	Ceil blue w/navy	each	\$20.29	\$19.28	\$19.28	\$18.26
7174	Dark teal w/black	each	\$20.29	\$19.28	\$19.28	\$18.26

Sizes 2XL-4XL also available, call for pricing



FASHION SEAL WOMEN'S DOUBLE V-NECK TUNICS

- Modern fit
- Set-in sleeves
- Princess seams
- Two lower pockets with contrasting trim
- Contrasting V-neck placket
- 65% polyester/35% cotton

Sizes XS-XL

7575	Ceil blue w/navy	each	\$21.79	\$20.70	\$20.70	\$19.61
7572	Chocolate w/pretty pink	each	\$21.79	\$20.70	\$20.70	\$19.61
7576	Eggplant w/plum	each	\$21.79	\$20.70	\$20.70	\$19.61
7578	Navy w/ceil blue	each	\$21.79	\$20.70	\$20.70	\$19.61
7579	Pewter w/black	each	\$21.79	\$20.70	\$20.70	\$19.61

Sizes 2XL-4XL also available, call for pricing



*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.

Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

UNIFORMS



FASHION SEAL WOMEN'S DRAWSTRING FLARE PANTS

- Modern fit
- Coordinating easy glide webbed drawcord
- Elastic waist back
- Two side inset pockets
- Right hip pocket
- Natural waist
- Flared legs with side slit
- 65% polyester/35% cotton
- Regular: 31" inseam*
- Petite: 28" inseam*

Sizes XS-XL

Color	Regular	Petite					
Black	7077	7086	each	\$31.99	\$30.39	\$30.39	\$28.79
Burgundy	7078	7087	each	\$31.99	\$30.39	\$30.39	\$28.79
Ceil blue	7080	7088	each	\$31.99	\$30.39	\$30.39	\$28.79
Navy	7082	7090	each	\$31.99	\$30.39	\$30.39	\$28.79
Pewter	7084	7091	each	\$31.99	\$30.39	\$30.39	\$28.79

Sizes 2XL-4XL also available, call for pricing

FASHION SEAL UNISEX FASHION SCRUB SHIRTS

- Unisex
- Extra long cap sleeves
- Left breast pocket
- Reinforced, bartacked v-neck
- 65% polyester/35% cotton
- Tall length is approximately 4-1/4" longer than regular length

Regular:	\$21.79	\$20.70	\$20.70	\$19.61
Tall:	\$21.79	\$20.70	\$20.70	\$19.61

Sizes XS-XL

Color	Regular	Tall	
Black	78729	78789	each
Blueberry	78725	6830	each
Burgundy	78741	78790	each
Caribbean blue	78726		each
Ceil blue	78737	78791	each
Chocolate	78004		each
Cobalt	78720	78792	each
Dark teal	78734	6831	each
Dill	78002		each
Eggplant	78730	78793	each
Fir green	78731	78797	each

Sizes 2XL-4XL also available, call for pricing

Color	Regular	Tall	
Jade	78000		each
Navy	78744	78794	each
Pewter	78785	78796	each
Plum	78722		each
Sage	78750	78795	each
Salmon	78005		each
Sangria	78003		each
Tan	78748	78798	each
Teal	78723		each
White	78735	78799	each



*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.

Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

FASHION SEAL UNISEX FASHION SCRUB PANTS

- Unisex
- Drawcord closure with easy glide coordinating webbing
- Right back pocket

- 65% polyester/35% cotton
- Regular: 31" inseam*
- Tall: 36" inseam*

Regular:	\$21.79	\$20.70	\$20.70	\$19.61
Tall:	\$21.79	\$20.70	\$20.70	\$19.61

Sizes XS-XL

Color	Regular	Tall	
Black	78829	78889	each
Blueberry	78825		each
Burgundy	78841	78890	each
Caribbean blue	78826		each
Ceil blue	78837	78891	each
Chocolate	78887		each
Cobalt	78851	78892	each
Dark teal	78834		each
Dill	78885		each
Eggplant	78830	78893	each
Fir green	78831	78897	each

Color	Regular	Tall	
Jade	78882		each
Navy	78844	78894	each
Pewter	78877	78896	each
Plum	78852		each
Sangria	78886		each
Sage	78850	78895	each
Salmon	78888		each
Tan	78848	78898	each
Teal	78853		each
White	78835	78881	each



Sizes 2XL-4XL also available, call for pricing

FASHION SEAL UNISEX WARM-UPS

- Unisex
- Raglan sleeves
- Two lower pockets
- Scissor holder in right pocket

- Snap front closure
- Coordinating knit cuffs
- 65% polyester/35% cotton

Sizes XS-XL

7673	Black	each	\$28.49	\$27.07	\$27.07	\$25.64
7688	Blueberry	each	\$28.49	\$27.07	\$27.07	\$25.64
7691	Burgundy	each	\$28.49	\$27.07	\$27.07	\$25.64
7690	Ceil blue	each	\$28.49	\$27.07	\$27.07	\$25.64
7675	Cobalt	each	\$28.49	\$27.07	\$27.07	\$25.64
7692	Dark teal	each	\$28.49	\$27.07	\$27.07	\$25.64
7043	Dill	each	\$28.49	\$27.07	\$27.07	\$25.64
7693	Eggplant	each	\$28.49	\$27.07	\$27.07	\$25.64
7694	Fir green	each	\$28.49	\$27.07	\$27.07	\$25.64
7671	Jade	each	\$28.49	\$27.07	\$27.07	\$25.64
7696	Navy	each	\$28.49	\$27.07	\$27.07	\$25.64
7607	Pewter	each	\$28.49	\$27.07	\$27.07	\$25.64
7682	Plum	each	\$28.49	\$27.07	\$27.07	\$25.64
7698	Sage	each	\$28.49	\$27.07	\$27.07	\$25.64
7672	Tan	each	\$28.49	\$27.07	\$27.07	\$25.64
7683	Teal	each	\$28.49	\$27.07	\$27.07	\$25.64
7699	White	each	\$28.49	\$27.07	\$27.07	\$25.64

Sizes 2XL-4XL also available, call for pricing

*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.



Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

UNIFORMS

FASHION SEAL MEN'S SHORT LAB COATS

- 34" length (based on size 40)
- Button front
- Left breast pocket with piping
- Two lower outside pockets with piping
- 65% polyester/35% cotton fine line twill



Sizes 32-50

447	White	each	\$27.79	\$26.40	\$26.40	\$25.01
-----	-------	------	---------	---------	---------	---------

Sizes 52-56 also available, call for pricing

FASHION SEAL MEN'S KNEE LENGTH LAB COATS

- 41" length (based on size 40)
- Button front
- Left breast pocket
- Two lower outside pockets
- Side slits for pocket access
- 65% polyester/35% cotton poplin



Sizes 32-50

485	White	each	\$28.49	\$27.07	\$27.07	\$25.64
-----	-------	------	---------	---------	---------	---------

Sizes 52-56 also available, call for pricing

FASHION SEAL UNISEX SHORT SLEEVE LAB COATS

- 34" length (based on medium size)
- Button front
- Left breast pocket
- Two lower outside pockets
- 65% polyester/35% cotton poplin



Sizes XXS-XL

3409	White	each	\$20.29	\$19.28	\$19.28	\$18.26
------	-------	------	---------	---------	---------	---------

Sizes 2XL, 3XL also available, call for pricing

FASHION SEAL UNISEX LAB COATS

- 41" length (based on size medium)
- Button front
- Left breast pocket
- Two lower oversized outside pockets
- Side slits for pocket access
- 80% polyester/20% cotton poplin



Sizes XXS-XL

3492	White	each	\$26.29	\$24.98	\$24.98	\$23.66
------	-------	------	---------	---------	---------	---------

Sizes 2XL, 3XL also available, call for pricing

**These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.*

VALUMAX® EXTRA-SAFE™ JACKETS & LAB COATS

Fashionable hip-length jackets and knee-length lab coats available in 32 attractive colors. 3-layer "SMS" fabrics meet OSHA fluid-resistant requirements, yet are lightweight and breathable for comfort. Professional appearance with snap-front closure, 3 pockets, knit cuffs and collar.

Latex-free, autoclavable and washable up to 5 times. Static-free: no clinging.

Please specify size and color when ordering. (See color chart below)

Hip-Length Jackets (Sizes XS-XL)

3630	10/pack	\$36.29	\$34.48	\$34.48	\$32.66
------	---------	---------	---------	---------	---------

Knee-Length Lab Coats (Sizes XS-XL)

3660	10/pack	\$40.79	\$38.75	\$38.75	\$36.71
------	---------	---------	---------	---------	---------

Sizes 2XL-5XL also available, call for pricing

Aqua-AQ	Blueberry-BB	Black-BK	Ceil Blue-CB
Deep Purple-DP	Steel Blue-EB	Eggplant-EG	Hot Pink-HP
Navy Blue-NB	Orange-OR	Olive-OV	Multicolor-PGO
Red-RD	Royal Blue-RY	Sky Blue-SB	True Blue-TB

MAYTEX LAB JACKETS

Professional lab jacket with comfortable, latex-free knit cuffs and collar. 3 pockets and 5 snaps in front. Higher neck provides better protection. Seams and pockets are sewn for garment-quality finish. Made of SMS material (Spunbond, Meltblown and Spunbond fused together for better protection while still allowing excellent breathability.) Disposable. Reusable when autoclaved.

#6200 LAB JACKET (SMS Material)								
	XS	S	M	L	XL	2XL	3XL	4XL
Blue		●	●	●	●	●		
Pink	●	●	●	●				
Teal	●	●	●	●	●	●		
Purple	●	●	●	●	●			
White	●	●	●	●	●	●		
Multi	●	●	●	●	●			



Please specify size and color when ordering

6200	Sizes XS-XL	30/case	\$134.49	\$127.77	\$127.77	\$121.04
------	-------------	---------	----------	----------	----------	----------

Sizes 2XL also available, call for pricing

Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

15

UNIFORMS

FASHION SEAL UNISEX PFAS-FREE TEXTURE SHIELD SNAP FRONT LAB COAT WITH THROAT-GUARD

- Front and sleeves made of Texture Shield single-ply fluid-resistant barrier fabric
- Back made of polyester/cotton blend
- Concealed placket front with snap closure
- Throat-Guard collar with adjustable snap closure for improved neck protection
- Two lower pockets
- Left breast pocket
- Side slits for pocket access
- Back vent
- Knit cuffs



Sizes XS-XL

6418	White	each	\$45.79	\$43.50	\$43.50	\$41.21
6419	Navy	each	\$45.79	\$43.50	\$43.50	\$41.21
6420	Ceill blue	each	\$45.79	\$43.50	\$43.50	\$41.21

Sizes 2XL-6XL also available, call for pricing

Style 6418 also available in XXS, 7XL, and 8XL, call for pricing

FASHION SEAL UNISEX PROTECTIVE COATS WITH "HOOK & LOOP" CLOSURE AT NECK

- Front and sleeves made of single-ply fluid-resistant barrier fabric
- Concealed placket front with snap closure
- Throat-Guard collar with "hook and loop" closure for improved neck protection
- Two lower pockets
- Back vent
- Knit cuffs
- Back constructed of polyester/cotton blend
- Fashion Shield® front



Sizes XS-XL

457	Ceill blue	each	\$41.79	\$39.70	\$39.70	\$37.61
450	White	each	\$41.79	\$39.70	\$39.70	\$37.61

Sizes 2XL, 3XL also available, call for pricing

FASHION SEAL UNISEX PFAS-FREE TEXTURE SHIELD D® APRON GOWN

Entire garment constructed of single-ply fluid-resistant barrier fabric. Backless design adds to the comfort of the worker in non-climate controlled areas (i.e. soiled linen and decontamination) Two sets of tie closures, one at the back of the neck, the other at the waist. With ciel blue knit cuffs.



Sizes S-XL

526	Ceill blue	each	\$30.79	\$29.25	\$29.25	\$27.71
-----	------------	------	---------	---------	---------	---------

Sizes 2XL-4XL also available, call for pricing

FASHION SEAL UNISEX PFAS-FREE FASHION SHIELD® APRON GOWN

- Entire garment constructed of single-ply fluid-resistant barrier fabric of Fashion Shield®, a high count, densely woven fabric employing 100% continuous filament man-made polyester yarns, offering excellent fluid resistance properties
- Backless design adds to the comfort of the worker in non-climate controlled areas (i.e. soiled linen and decontamination)
- Two sets of tie closures; one set at back of neck, one set at waist
- White knit cuffs
- Permanent grid for process cycle control



Sizes S-XL

529	Ceill blue	each	\$25.49	\$24.22	\$24.22	\$22.94
-----	------------	------	---------	---------	---------	---------

Sizes 2XL-4XL also available, call for pricing

FASHION SEAL ANGLE BACK FULL SIZE PATIENT GOWN

- Poly/cotton fashion sheeting meets California Fire Marshall's Title 19
- Two ties at waist for adjustability
- Extra large telemeter pocket
- Snap closure on sleeves
- 51" length
- 6" insleeve
- 76" sweep
- Color-coded placket on right sleeve
- 3-1/2" poly bound fenstration



727	Fits up to 3XL, blue optics print	each	\$28.79	\$27.35	\$27.35	\$25.91
-----	-----------------------------------	------	---------	---------	---------	---------

*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.

Please refer to the Fashion Seal size and color charts on page 294 to confirm your size and color selection before ordering. Fashion Seal scrubs and lab coats are available with embroidery. Please see page 295 for info.

UNIFORMS

FASHION SEAL PRODUCTS

YOUR STYLE, YOUR COLOR, YOUR FIT, YOUR WAY

Step 1: Choose your color from the chart below (actual scrub colors may vary from color illustration)



STEP 2: CHOOSE YOUR SIZE BY REFERRING TO THE CHARTS BELOW FOR SCRUB AND LAB COAT SIZING

Ladies' Sizes (Applies to ladies with figures of average proportions and height of 5'3" - 5'7" without shoes)

Women's Size Chart

Order Size	XXS	XS	S	M	L	XL	2XL**	3XL**	4XL**
Numeric Sizes	2	4-6	8-10	12-14	16-18	20-42	44-46	48-50	52-54
Chest	31-32	33-34	35-36	37-39	40-43	44-47	48-51	52-55	56-59
Waist	21-23	24-25	26-27	28-30	31-34	35-38	39-42	43-46	47-50
Hips	33-34	35-36	37-38	39-41	42-45	46-49	50-53	54-57	58-61
Petite Inseam*	28	28	28	28	28	28	28	28	28
Regular Inseam*	31	31	31	31	31	31	31	31	31
Tall Inseam*	34	34	34	34	34	34	34	34	34

Unisex Size Chart


Order Size	XXS	XS	S	M	L	XL	2XL**	3XL**	4XL**	5XL**
Women's Numeric Sizes	2	4-6	8-10	12-14	16-18	20-42	44-46	48-50	52-54	
Chest	29-31	32-24	35-37	38-40	42-44	46-48	50-52	54-56	58-60	62-64
Waist	21-23	24-26	27-29	30-32	34-36	38-40	42-44	46-48	50-52	54-56
Short Inseam*	29	29	29	29	29	29	29	29	29	29
Regular Inseam*	31	31	31	31	31	31	31	31	31	31
Long Inseam*	36	36	36	36	36	36	36	36	36	36

**Where available, at extra charge

*These are general guidelines, some styles' inseams lengthen at predetermined sizes. Contact customer service for more details.

WOUND CLOSURE/SUTURE

STAINLESS STEEL NEEDLE TYPES

Taper Point	
Taper Cut	
Blunt Point	
Conventional Cutting	
Precision Point	
Reverse Cutting	
Micro-Point Reverse Cutting	

NEEDLE COMPARISON CHART

A Selection of Popular Needles

Needle Description	ACE	Advance Plus	A. Titan	Ethicon	Hu-Friedy	Look	Myco Reli Pro	Riverpoint
Reverse Cutting								
3/8 Circle - 11mm				C-2		C16		
3/8 Circle - 12mm	DS12			C-2	C17	C17	YC-2	
3/8 Circle - 12mm				C-3				
3/8 Circle - 16.2mm	DS16			FS-3	C22	C22	YFS-3	FS-3
3/8 Circle - 19mm	DS18	C6	C-6	FS-2	C6	C6	YFS-2	FS-2
3/8 Circle - 24mm	DS24	C7	C-7	FS-1	C7	C7	YFS-1	
3/8 Circle - 26mm				FS		C8	YFS	
3/8 Circle - 30mm				FSL		C10		
1/2 Circle - 15mm	HS16			M-1	C43	C26		
1/2 Circle - 20mm				J-1		C20		
1/2 Circle - 24mm	HS24	C31	C-9	X-1	C9	C31	YX-1	
1/2 Circle - 26mm				CP-2		C21	YCP-2	
1/2 Circle - 37mm				CP-1		C13	YCP-1	
1/2 Circle - 40mm				CP		C14	YCP	
Precision Reverse Cutting								
3/8 Circle - 11mm	P1		C-1	P-1	C1	C1	YP-1	
3/8 Circle - 13mm	DSM13	C3	C-3	P-3	C3	C3	YP-3	P-3
3/8 Circle - 16mm				PS-3		PC34		
3/8 Circle - 19mm	DSM18		C-31	PS-2		PC31	YPS-2	PS-2
3/8 Circle - 24.3mm				PS-1		PC33		
1/2 Circle - 18mm				PS-4		PC30		
Conventional Cutting Needles								
3/8 Circle - 11mm				PC-1		CP1		
Taper Point								
3/8 Circle - 12mm				C-1		T12		
1/2 Circle - 13mm				TF		T22		TF
1/2 Circle - 18mm	HR17			RB-1		T28		
1/2 Circle - 22mm	HR22			SH-1		T4		
1/2 Circle - 26mm	HR26			SH		T		
1/2 Circle - 26.5mm				CT-2		T5		
1/2 Circle - 36mm				MH		T9		
1/2 Circle - 36.6mm				CT-1		T8		
1/2 Circle - 40mm				CT		T10		
Taper Cut								
3/8 Circle - 17mm				RB-1		T28		
1/2 Circle - 16mm	HRT17			V-5	D14	Y5		
1/2 Circle - 17mm				V-5		D14		



REVERSE CUTTING NEEDLES

	3/8 Circle FS3
	FS2
	C2
	P1
	1/2 Circle M1
	C19

ROUND DIAMOND NEEDLES

	1/2 Circle TF
	RB1

ACE PTFE SUTURES

- Undyed white monofilament
- Non-absorbable
- 12/Box

Size:	ACE	Advance Plus	A. Titan	Ethicon	Hu-Friedy	Look	Myco Reli Pro	Riverpoint	
Size: 2-0				Comparable to ETHICON					
430124020	DS18	18"		FS2		\$119.49	\$113.52	\$113.52	\$107.54
Size: 3-0									
430124030	DS18	18"		FS2		\$119.49	\$113.52	\$113.52	\$107.54
430124031	DS16	18"		FS3		\$119.49	\$113.52	\$113.52	\$107.54
430124032	HS16	18"		M1		\$119.49	\$113.52	\$113.52	\$107.54
Size: 4-0									
430124040	DS16	18"		FS3		\$119.49	\$113.52	\$113.52	\$107.54
430124041	DS12	18"		C2		\$119.49	\$119.49	\$107.54	\$107.54
430124042	TF	18"		TF		\$127.79	\$121.40	\$121.40	\$115.01
430124043	DS18	18"		FS2		\$119.49	\$113.52	\$113.52	\$107.54
Size: 5-0									
430124050	C19	18"		C19		\$119.49	\$119.49	\$107.54	\$107.54
430124051	DS12	18"		C2		\$119.49	\$113.52	\$113.52	\$107.54
430124052	TF	18"		RB1		\$127.79	\$121.40	\$121.40	\$115.01
Size: 6-0									
430124062	DS12	18"		C2		\$119.49	\$113.52	\$113.52	\$107.54

ACE CHROMIC GUT

- Absorbable • 12/box

			Comparable to ETHICON	Comparable to LOOK					
Size: 3-0									
0032484	DS18	18"	636G	-	\$30.29	\$28.78	\$28.78	\$27.26	
0030043	DS18	27"	636G	560B	\$31.29	\$29.73	\$29.73	\$28.16	
0032503	DS24	30"	-	-	\$42.49	\$40.37	\$40.37	\$38.24	
0032505	HR26	30"	G122H	-	\$48.79	\$46.35	\$46.35	\$43.91	
Size: 4-0									
0032485	DS18	18"	-	558B	\$31.29	\$29.73	\$29.73	\$28.16	
0032482	DS18	27"	635G	-	\$31.29	\$29.73	\$29.73	\$28.16	
0032494	DSM13	18"	1654G	595B	\$62.49	\$59.37	\$59.37	\$56.24	
0032498	T4	14"	-	589B	\$32.29	\$30.68	\$30.68	\$29.06	
Size: 5-0									
0032499	T28	14"	-	590B	\$32.29	\$30.68	\$30.68	\$29.06	
0032502	DSM13	18"	-	1248B	\$62.49	\$59.37	\$59.37	\$56.24	

ACE PLAIN GUT

- Absorbable • 12/box

			Comparable to ETHICON	Comparable to LOOK					
Size: 3-0									
0032480	DS18	18"	H822G	592B	\$31.29	\$29.73	\$29.73	\$28.16	
0032504	DS18	27"	-	-	\$30.29	\$28.78	\$28.78	\$27.26	
0032481	HS24	18"	612G	594B	\$29.49	\$28.02	\$28.02	\$26.54	
0032493	DS24	30"	H810H	-	\$27.49	\$26.12	\$26.12	\$24.74	
Size: 4-0									
0032483	DS18	18"	7771	591B	\$31.29	\$29.73	\$29.73	\$28.16	
0032506	DS18	27"	H821G	551B	\$29.49	\$28.02	\$28.02	\$26.54	
Size: 5-0									
0032490	DSM13	18"	686G	596B	\$60.49	\$57.47	\$57.47	\$54.44	

ACE PGA

- Equivalent to Vicryl®** • Absorption is essentially complete within 50-90 days • PGA Coated • 12/box

Size: 4-0									
0034220	Undyed braided		DS18	27"	\$47.49	\$45.12	\$45.12	\$42.74	
0034960	Undyed braided		DSM18	18"	\$83.29	\$79.13	\$79.13	\$74.96	
0030494	Undyed braided		DSM13	18"	\$83.29	\$79.13	\$79.13	\$74.96	
Size: 5-0									
0034930	Undyed braided		DSM13	18"	\$94.99	\$90.24	\$90.24	\$85.49	
Size: 6-0									
0030492	Undyed braided		DSM13	18"	\$99.29	\$94.33	\$94.33	\$89.36	

ACE PGA

- Equivalent to Vicryl®** • Absorption is essentially complete within 50-90 days • PGA Coated • 12/box

Size: 3-0									
0033930	Violet braided		DS18	18"	\$44.79	\$42.55	\$42.55	\$40.31	
Size: 4-0									
0033920	Violet braided		DS18	18"	\$46.29	\$43.98	\$43.98	\$41.66	
0034640	Violet braided		DSM13	18"	\$90.79	\$86.25	\$86.25	\$81.71	
Size: 5-0									
0034630	Violet braided		DSM13	18"	\$94.99	\$90.24	\$90.24	\$85.49	

** Vicryl® is a Registered Trademarks of Ethicon, Inc.

WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

- 3/8 Circle
DS18(FS2)
- DS24(FS1)
- DSM13(P3)
- DSM18(PS2)
- DS12(C2)
- 1/2 Circle
HS24(X1)
- HS16(X8)



TAPER POINT NEEDLES

- 1/2 Circle
HRT17(V5)
- T4(SH1)
- T28(RB1)



WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle
DS18(FS2)

DS24(FS1)

DSM13(P3)

DSM18(PS2)

DS12(C2)

1/2 Circle
HS24(X1)

HS16(X8)



TAPER POINT NEEDLES

1/2 Circle
HRT17(V5) Polyester

T4(SH1)

T28(RB1)



ACE PGA-RAPID ABSORBING*

- Equivalent to Vicryl Rapide[®]**
- Absorption is essentially complete within 49-63 days
- 12/box

Size: 3-0

0039907	Undyed braided	DSM18	18"	\$110.79	\$105.25	\$105.25	\$99.71
---------	----------------	-------	-----	----------	----------	----------	---------

Size: 4-0

0039906	Undyed braided	DSM18	18"	\$110.79	\$105.25	\$105.25	\$99.71
0039900	Undyed braided	HRT17	27"	\$113.29	\$107.63	\$107.63	\$101.96

Size: 5-0

0039904	Undyed braided	DSM13	18"	\$110.79	\$105.25	\$105.25	\$99.71
---------	----------------	-------	-----	----------	----------	----------	---------

* Maintains 50% tensile strength at 5 days; 0% at 14 days (data on file). | ** Vicryl[®] and Vicryl Rapide[®] are Registered Trademarks of Ethicon, Inc.

ACE SILK

- Black braided, non-absorbable
- 12/box

Size: 3-0

			Comparable to ETHICON	Comparable to LOOK				
0032487	DS18	18"	2772	784B	\$25.79	\$24.50	\$24.50	\$23.21
0032491	DS24	18"	684G	786B	\$24.99	\$23.74	\$23.74	\$22.49
0032488	HS24	18"	632G	785B	\$23.49	\$22.32	\$22.32	\$21.14
0032486	HS16	18"	7762	-	\$29.49	\$28.02	\$28.02	\$26.54

Size: 4-0

0032489	DS18	18"	683G	-	\$25.79	\$24.50	\$24.50	\$23.21
0032492	DSM13	18"	641G	-	\$60.49	\$57.47	\$57.47	\$54.44

Size: 5-0

0032479	DSM13	18"	-	-	\$80.29	\$76.28	\$76.28	\$72.26
---------	-------	-----	---	---	---------	---------	---------	---------

ACE NYLON

- Black monofilament, non-absorbable
- 12/box

Size: 5-0

			Comparable to ETHICON	Comparable to LOOK				
0032460	DS12	10"	-	913B	\$35.99	\$34.19	\$34.19	\$32.39

Size: 6-0

0032461	DS12	10"	-	911B	\$35.99	\$34.19	\$34.19	\$32.39
---------	------	-----	---	------	---------	---------	---------	---------

ACE POLYPROPYLENE

- Blue monofilament, non-absorbable
- 12/box

Size: 5-0

			Comparable to ETHICON	Comparable to LOOK				
0032501	DSM13	18"	8698G	-	\$77.29	\$73.43	\$73.43	\$69.56

ACE POLYESTER

- Green braided, non-absorbable
- 12/box

Size: 5-0

			Comparable to ETHICON	Comparable to LOOK				
0032500	HRT17	18"	7773	797B	\$50.99	\$48.44	\$48.44	\$45.89

WOUND CLOSURE/SUTURE

ADVANCE+ PLUS

ADVANCE PLUS CHROMIC GUT

- Absorbable
- 12/box

Size: 2-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV561B	C-7	27"	561B	\$22.49	\$21.37	\$21.37	\$20.24
Size: 3-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV525B	C-6	18"	525B	\$22.49	\$21.37	\$21.37	\$20.24
ADV560B	C-6	27"	560B	\$22.49	\$21.37	\$21.37	\$20.24
Size: 4-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV595B	C-3	18"	595B	\$52.99	\$50.34	\$50.34	\$47.69
ADV558B	C-6	18"	558B	\$22.49	\$21.37	\$21.37	\$20.24
ADV559B	C-6	27"	559B	\$22.49	\$21.37	\$21.37	\$20.24
ADV589B	T-4	14"	589B	\$28.99	\$27.54	\$27.54	\$26.09
Size: 5-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV1248B	C-3	18"	1248B	\$52.99	\$50.34	\$50.34	\$47.69

ADVANCE PLUS PLAIN GUT

- Absorbable
- 12/box

Size: 3-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV592B	C-6	18"	592B	\$22.49	\$21.37	\$21.37	\$20.24
ADV554B	C-6	27"	554B	\$22.49	\$21.37	\$21.37	\$20.24
ADV553B	C-7	27"	553B	\$22.49	\$21.37	\$21.37	\$20.24
ADV594B	C-31	18"	594B	\$22.49	\$21.37	\$21.37	\$20.24
Size: 4-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV591B	C-6	18"	591B	\$22.49	\$21.37	\$21.37	\$20.24
ADV551B	C-6	27"	551B	\$22.49	\$21.37	\$21.37	\$20.24

ADVANCE PLUS SILK

- Black braided, non-absorbable
- 12/box

Size: 3-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV784B	C-6	18"	784B	\$20.99	\$19.94	\$19.94	\$18.89
Size: 4-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV781B	C-6	18"	781B	\$21.99	\$20.89	\$20.89	\$19.79

ADVANCE PLUS POLYGLYCOLIC ACID

- Undyed braided, absorbable
- 12/box

Size: 4-0			Comparable to LOOK	Regular Price	Member Price	Premier Price	Premier+ Price
ADV422B	C-6	27"	422B	\$47.99	\$45.59	\$45.59	\$43.19

REVERSE CUTTING NEEDLES

3/8 Circle

C6

C7

1/2 Circle

C31



REVERSE REVERSE CUTTING NEEDLES

3/8 Circle

C3



TAPER POINT NEEDLES

1/2 Circle

T-4



WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle

C16

C17

C6

PC31

C1

PC34

1/2 Circle

C26

C19

C20

C31

C21

C2



REVERSE CUTTING HAND-HONED NEEDLES

3/8 Circle

C3

1/2 Circle

Y5



CONVENTIONAL CUTTING NEEDLES

3/8 Circle

CP1



LOOK CHROMIC GUT

• 12/box

Size: 2-0

561	Absorbable	C-7	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

Size: 3-0

525	Absorbable	C-6	18"	\$32.49	\$30.87	\$30.87	\$29.24
-----	------------	-----	-----	---------	---------	---------	---------

560	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

Size: 4-0

595	Absorbable	C-3	18"	\$66.99	\$63.64	\$63.64	\$60.29
-----	------------	-----	-----	---------	---------	---------	---------

558	Absorbable	C-6	18"	\$32.49	\$30.87	\$30.87	\$29.24
-----	------------	-----	-----	---------	---------	---------	---------

559	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

546B	Absorbable	C-17	18"	\$28.99	\$27.54	\$27.54	\$26.09
------	------------	------	-----	---------	---------	---------	---------

547	Absorbable	C-20	18"	\$30.79	\$29.25	\$29.25	\$27.71
-----	------------	------	-----	---------	---------	---------	---------

544	Absorbable	C-26	18"	\$30.79	\$29.25	\$29.25	\$27.71
-----	------------	------	-----	---------	---------	---------	---------

589	Absorbable	T-4	14"	\$30.49	\$28.97	\$28.97	\$27.44
-----	------------	-----	-----	---------	---------	---------	---------

Size: 5-0

1248	Absorbable	C-3	18"	\$63.49	\$60.32	\$60.32	\$57.14
------	------------	-----	-----	---------	---------	---------	---------

557	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

543	Absorbable	C-16	18"	\$31.49	\$29.92	\$29.92	\$28.34
-----	------------	------	-----	---------	---------	---------	---------

563	Absorbable	C-17	18"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	------	-----	---------	---------	---------	---------

598	Absorbable	C-26	18"	\$35.99	\$34.19	\$34.19	\$32.39
-----	------------	------	-----	---------	---------	---------	---------

590	Absorbable	T-28	14"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	------	-----	---------	---------	---------	---------

Size: 6-0

1242	Absorbable	C-1	12"	\$142.99	\$135.84	\$135.84	\$128.69
------	------------	-----	-----	----------	----------	----------	----------

LOOK PLAIN GUT

• 12/box

Size: 2-0

555	Absorbable	C-7	27"	\$32.49	\$30.87	\$30.87	\$29.24
-----	------------	-----	-----	---------	---------	---------	---------

Size: 3-0

592	Absorbable	C-6	18"	\$32.49	\$30.87	\$30.87	\$29.24
-----	------------	-----	-----	---------	---------	---------	---------

554	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

553	Absorbable	C-7	27"	\$33.49	\$31.82	\$31.82	\$30.14
-----	------------	-----	-----	---------	---------	---------	---------

593	Absorbable	C-21	18"	\$30.29	\$28.78	\$28.78	\$27.26
-----	------------	------	-----	---------	---------	---------	---------

594	Absorbable	C-31	18"	\$30.49	\$28.97	\$28.97	\$27.44
-----	------------	------	-----	---------	---------	---------	---------

Size: 4-0

1241	Absorbable	C-3	18"	\$61.79	\$58.70	\$58.70	\$55.61
------	------------	-----	-----	---------	---------	---------	---------

591	Absorbable	C-6	18"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

551	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

599	Absorbable	C-31	18"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	------	-----	---------	---------	---------	---------

Size: 5-0

596	Absorbable	C-3	18"	\$66.79	\$63.45	\$63.45	\$60.11
-----	------------	-----	-----	---------	---------	---------	---------

550	Absorbable	C-6	27"	\$30.99	\$29.44	\$29.44	\$27.89
-----	------------	-----	-----	---------	---------	---------	---------

535	Absorbable	C-17	18"	\$34.99	\$33.24	\$33.24	\$31.49
-----	------------	------	-----	---------	---------	---------	---------

Size: 6-0

1415	Absorbable	CP-1	18"	\$138.79	\$131.85	\$131.85	\$124.91
------	------------	------	-----	----------	----------	----------	----------

LOOK POLYGLYCOLIC ACID

• 12/box

Size: 2-0

443	Undyed braided, absorbable	C-7	27"	\$51.29	\$48.73	\$48.73	\$46.16
-----	----------------------------	-----	-----	---------	---------	---------	---------

Size: 3-0

423	Undyed braided, absorbable	C-6	27"	\$51.29	\$48.73	\$48.73	\$46.16
442	Undyed braided, absorbable	C-7	27"	\$53.79	\$51.10	\$51.10	\$48.41
497	Undyed braided, absorbable	PC-31	18"	\$96.29	\$91.48	\$91.48	\$86.66

Size: 4-0

494	Undyed braided, absorbable	C-3	18"	\$96.79	\$91.95	\$91.95	\$87.11
422	Undyed braided, absorbable	C-6	27"	\$55.29	\$52.53	\$52.53	\$49.76
386	Undyed braided, absorbable	C-17	10"	\$53.79	\$51.10	\$51.10	\$48.41
496	Undyed braided, absorbable	PC-31	18"	\$97.79	\$92.90	\$92.90	\$88.01

Size: 5-0

493	Undyed braided, absorbable	C-3	18"	\$104.29	\$99.08	\$99.08	\$93.86
421	Undyed braided, absorbable	C-6	27"	\$51.49	\$48.92	\$48.92	\$46.34
385	Undyed braided, absorbable	C-17	10"	\$54.79	\$52.05	\$52.05	\$49.31

Size: 6-0

492	Undyed braided, absorbable	C-3	18"	\$99.79	\$94.80	\$94.80	\$89.81
-----	----------------------------	-----	-----	---------	---------	---------	---------

LOOK SILK

• 12/box

Size: 2-0

770	Black braided, non-absorbable	C-8	18"	\$24.79	\$23.55	\$23.55	\$22.31
-----	-------------------------------	-----	-----	---------	---------	---------	---------

Size: 3-0

784	Black braided, non-absorbable	C-6	18"	\$27.99	\$26.59	\$26.59	\$25.19
768	Black braided, non-absorbable	C-6	30"	\$31.79	\$30.20	\$30.20	\$28.61
786	Black braided, non-absorbable	C-7	18"	\$28.29	\$26.88	\$26.88	\$25.46
795	Black braided, non-absorbable	C-16	18"	\$29.49	\$28.02	\$28.02	\$26.54
782	Black braided, non-absorbable	C-21	18"	\$24.79	\$23.55	\$23.55	\$22.31
783	Black braided, non-absorbable	C-19	18"	\$27.99	\$26.59	\$26.59	\$25.19
785	Black braided, non-absorbable	C-31	18"	\$25.29	\$24.03	\$24.03	\$22.76

Size: 4-0

792	Black braided, non-absorbable	C-3	18"	\$54.99	\$52.24	\$52.24	\$49.49
781	Black braided, non-absorbable	C-6	18"	\$27.99	\$26.59	\$26.59	\$25.19
776	Black braided, non-absorbable	C-7	18"	\$26.29	\$24.98	\$24.98	\$23.66
780	Black braided, non-absorbable	C-20	18"	\$26.29	\$24.98	\$24.98	\$23.66
787	Black braided, non-absorbable	C-17	18"	\$29.49	\$28.02	\$28.02	\$26.54

Size: 5-0

1265	Black braided, non-absorbable	C-3	18"	\$92.29	\$87.68	\$87.68	\$83.06
774	Black braided, non-absorbable	C-17	18"	\$25.79	\$24.50	\$24.50	\$23.21

Size: 6-0

752	Black braided, non-absorbable	C-17	18"	\$29.49	\$28.02	\$28.02	\$26.54
-----	-------------------------------	------	-----	---------	---------	---------	---------

WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle

C16 _____

C17 _____

C6 _____

PC31 _____

C1 _____

PC34 _____

1/2 Circle

C26 _____

C19 _____

C20 _____

C31 _____

C21 _____

C2 _____



REVERSE CUTTING HAND-HONED NEEDLES

3/8 Circle

C3 _____

1/2 Circle

Y5 _____



CONVENTIONAL CUTTING NEEDLES

3/8 Circle

CP1 _____



WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle

C16

C17

C6

PC31

C1

PC34

1/2 Circle

C26

C19

C20

C31

C21

C2



REVERSE CUTTING HAND-HONED NEEDLES

3/8 Circle

C3

1/2 Circle

Y5



CONVENTIONAL CUTTING NEEDLES

3/8 Circle

CP1



LOOK NYLON

- 12/box

Size: 4-0

1291	Black monofilament, non-absorbable	C-3	18"	\$83.49	\$79.32	\$79.32	\$75.14
922	Black monofilament, non-absorbable	C-6	18"	\$37.49	\$35.62	\$35.62	\$33.74
931	Black monofilament, non-absorbable	DS-12	10"	\$38.29	\$36.38	\$36.38	\$34.46

Size: 5-0

1285	Black monofilament, non-absorbable	C-3	18"	\$88.79	\$84.35	\$84.35	\$79.91
918	Black monofilament, non-absorbable	C-6	18"	\$38.79	\$36.85	\$36.85	\$34.91
917	Black monofilament, non-absorbable	C-17	18"	\$39.79	\$37.80	\$37.80	\$35.81

Size: 6-0

1279	Black monofilament, non-absorbable	C-1	18"	\$97.79	\$92.90	\$92.90	\$88.01
1282	Black monofilament, non-absorbable	C-3	18"	\$93.99	\$89.29	\$89.29	\$84.59
916	Black monofilament, non-absorbable	C-17	18"	\$39.49	\$37.52	\$37.52	\$35.54
915	Black monofilament, non-absorbable	C-22	18"	\$44.29	\$42.08	\$42.08	\$39.86

LOOK POLYPROPYLENE

- 12/box

Size: 4-0

1022	Blue monofilament, non-absorbable	C-6	18"	\$58.79	\$55.85	\$55.85	\$52.91
------	-----------------------------------	-----	-----	---------	---------	---------	---------

Size: 5-0

1085	Blue monofilament, non-absorbable	C-3	18"	\$120.79	\$114.75	\$114.75	\$108.71
1018	Blue monofilament, non-absorbable	C-6	18"	\$61.99	\$58.89	\$58.89	\$55.79
1013	Blue monofilament, non-absorbable	C-17	10"	\$49.99	\$47.49	\$47.49	\$44.99

Size: 6-0

1079	Blue monofilament, non-absorbable	C-1	18"	\$116.49	\$110.67	\$110.67	\$104.84
1082	Blue monofilament, non-absorbable	C-3	18"	\$110.29	\$104.78	\$104.78	\$99.26
LK1011	Blue monofilament, non-absorbable	C-17	10"	\$47.79	\$45.40	\$45.40	\$43.01

LOOK POLYVIOLENE

- 12/box

Size: 4-0

1180	Green braided, non-absorbable	C-6	18"	\$33.29	\$31.63	\$31.63	\$29.96
------	-------------------------------	-----	-----	---------	---------	---------	---------

LOOK™ PTFE DENTAL SUTURE

A trusted brand. Consistent quality and value. The LOOK™ line of dental sutures is now even better with PTFE.

- Soft and supple material
- Excellent handling and knot security
- * Minimal package memory
- White

Size: 2-0

815B	Monofilament, non-absorbable	C-6	18"	12/box	\$89.29	\$84.83	\$84.83	\$80.36
------	------------------------------	-----	-----	--------	---------	---------	---------	---------

Size: 3-0

816B	Monofilament, non-absorbable	C-6	18"	12/box	\$101.79	\$96.70	\$96.70	\$91.61
818B	Monofilament, non-absorbable	C-22	18"	12/box	\$89.29	\$84.83	\$84.83	\$80.36

Size: 4-0

820B	Monofilament, non-absorbable	C-3	18"	12/box	\$97.79	\$92.90	\$92.90	\$88.01
817B	Monofilament, non-absorbable	C-22	18"	12/box	\$89.29	\$84.83	\$84.83	\$80.36
821B	Monofilament, non-absorbable	PC-31	18"	12/box	\$97.79	\$92.90	\$92.90	\$88.01
819B	Monofilament, non-absorbable	T-22	18"	12/box	\$89.29	\$84.83	\$84.83	\$80.36

Size: 5-0

823B	Monofilament, non-absorbable	C-3	18"	12/box	\$111.79	\$106.20	\$106.20	\$100.61
------	------------------------------	-----	-----	--------	----------	----------	----------	----------

Size: 6-0

822B	Monofilament, non-absorbable	C-3	18"	12/box	\$111.79	\$106.20	\$106.20	\$100.61
------	------------------------------	-----	-----	--------	----------	----------	----------	----------

WOUND CLOSURE/SUTURE



RELI PRO REDIGUT® CHROMIC PRO

- Absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
GC638M	YFS-1	30"	638G	\$30.49	\$28.97	\$28.97	\$27.44
GC6364M	YFS-2	18"		\$27.99	\$26.59	\$26.59	\$25.19
GC636M	YFS-2	30"	636H	\$33.49	\$31.82	\$31.82	\$30.14
Size: 4-0				\$27.49	\$26.12	\$26.12	\$24.74
GC6354M	YFS-2	18"		\$27.49	\$26.12	\$26.12	\$24.74
GC635M	YFS-2	30"	635H	\$33.49	\$31.82	\$31.82	\$30.14
GC1654M	YP-3	18"	1654G	\$39.79	\$37.80	\$37.80	\$35.81
Size: 5-0				\$35.79	\$34.00	\$34.00	\$32.21
GC634M	YFS-2	30"	634G	\$35.79	\$34.00	\$34.00	\$32.21
GC687M	YP-3	18"	687G	\$39.79	\$37.80	\$37.80	\$35.81
Size: 6-0				\$39.79	\$37.80	\$37.80	\$35.81
GC621M	YC-2	18"		\$39.79	\$37.80	\$37.80	\$35.81

RELI PRO REDIGUT® PLAIN PRO

- Absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
GP3247M	YFS-1	30"	H810H	\$27.29	\$25.93	\$25.93	\$24.56
GP8881M	YFS-2	18"		\$27.79	\$26.40	\$26.40	\$25.01
GP822M	YFS-2	30"	H822H	\$24.29	\$23.08	\$23.08	\$21.86
GP612M	YX-1	18"	612G	\$26.79	\$25.45	\$25.45	\$24.11
Size: 4-0				\$27.79	\$26.40	\$26.40	\$25.01
GP7771M	YFS-2	18"		\$27.79	\$26.40	\$26.40	\$25.01
GP821M	YFS-2	30"	H821H	\$41.49	\$39.42	\$39.42	\$37.34
Size: 5-0				\$37.49	\$35.62	\$35.62	\$33.74
GP686M	YP-3	18"	686G	\$37.49	\$35.62	\$35.62	\$33.74

RELI PRO REDISORB® PRO POLYGLYCOLIC LACTIDE ACID

- Undyed braided, absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
VU423M	YFS-2	30"	J423H	\$39.29	\$37.33	\$37.33	\$35.36
Size: 4-0				\$74.49	\$70.77	\$70.77	\$67.04
VU494M	YP-3	18"	J494H	\$74.49	\$70.77	\$70.77	\$67.04
VU422M	YFS-2	30"	J422H	\$37.79	\$35.90	\$35.90	\$34.01
VU496M	YFS-2	18"	J496H	\$60.99	\$57.94	\$57.94	\$54.89
Size: 5-0				\$67.29	\$63.93	\$63.93	\$60.56
VU421M	YFS-2	30"	J421H	\$67.29	\$63.93	\$63.93	\$60.56
VU493M	YP-3	18"	J493H	\$66.29	\$62.98	\$62.98	\$59.66
Size: 6-0				\$90.29	\$85.78	\$85.78	\$81.26
VU492M	YP-3	18"	J492H	\$90.29	\$85.78	\$85.78	\$81.26

REVERSE CUTTING NEEDLES

3/8 Circle
YFS-2 (C6)

YFS-1 (C7)

1/2 Circle
YX-1 (C31)



REVERSE REVERSE CUTTING NEEDLES

3/8 Circle
YP-3 (C3)



WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle
YFS-2 (C6)

YFS-1 (C7)

1/2 Circle
YX-1 (C31)



REVERSE REVERSE CUTTING NEEDLES

3/8 Circle
YP-3 (C3)



RELI PRO REDISORB® PRO POLYGLYCOLIC LACTIDE ACID

- Violet braided, absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
V442M	YFS-1	18"		\$33.49	\$31.82	\$31.82	\$30.14
Size: 4-0							
V392M	YFS-2	18"	J392H	\$38.99	\$37.04	\$37.04	\$35.09
Size: 5-0							
V463M	YP-3	18"	J463G	\$66.29	\$62.98	\$62.98	\$59.66

RELI PRO REDISILK® PRO

- Black braided, non-absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
SK684M	YFS-1	18"	684G	\$22.99	\$21.84	\$21.84	\$20.69
SK7772M	YFS-2	18"		\$24.29	\$23.08	\$23.08	\$21.86
SK632M	YX-1	18"	632G	\$23.29	\$22.13	\$22.13	\$20.96
Size: 4-0							
SK683M	YFS-2	18"	683G	\$24.79	\$23.55	\$23.55	\$22.31
SK641M	YP-3	18"	641G	\$35.79	\$34.00	\$34.00	\$32.21

ETHICON

ETHICON CHROMIC GUT

Size: 3-0								
636H	Monofilament, absorbable	FS-2	27"	36/box	\$468.79	\$445.35	\$445.35	\$421.91
1638H	Monofilament, absorbable	PS-2	27"	36/box	\$887.29	\$842.93	\$842.93	\$798.56
U204H	Monofilament, absorbable	RB-1	27"	36/box	\$518.79	\$492.85	\$492.85	\$466.91
G122H	Monofilament, absorbable	SH	27"	36/box	\$471.99	\$448.39	\$448.39	\$424.79
G182H	Monofilament, absorbable	SH-1	27"	36/box	\$513.29	\$487.63	\$487.63	\$461.96
Size: 4-0								
635H	Monofilament, absorbable	FS-2	27"	36/box	\$468.79	\$445.35	\$445.35	\$421.91
724G	Monofilament, absorbable	J-1	18"	12/box	\$149.29	\$141.83	\$141.83	\$134.36
1654G	Monofilament, absorbable	P-3	18"	12/box	\$267.29	\$253.93	\$253.93	\$240.56
1637G	Monofilament, absorbable	PS-2	18"	12/box	\$341.79	\$324.70	\$324.70	\$307.61
1643G	Monofilament, absorbable	PS-4	18"	12/box	\$281.29	\$267.23	\$267.23	\$253.16
U203H	Monofilament, absorbable	RB-1	27"	36/box	\$519.79	\$493.80	\$493.80	\$467.81
G121H	Monofilament, absorbable	SH	27"	36/box	\$513.79	\$488.10	\$488.10	\$462.41
G181H	Monofilament, absorbable	SH-1	27"	36/box	\$513.29	\$487.63	\$487.63	\$461.96
Size: 5-0								
634G	Monofilament, absorbable	FS-2	27"	12/box	\$153.99	\$146.29	\$146.29	\$138.59
1658G	Monofilament, absorbable	P-2	18"	12/box	\$324.79	\$308.55	\$308.55	\$292.31
687G	Monofilament, absorbable	P-3	18"	12/box	\$273.79	\$260.10	\$260.10	\$246.41
1636G	Monofilament, absorbable	PS-3	18"	12/box	\$283.99	\$269.79	\$269.79	\$255.59
1642G	Monofilament, absorbable	PS-4	18"	12/box	\$152.49	\$144.87	\$144.87	\$137.24
1652G	Monofilament, absorbable	PS-5	18"	12/box	\$315.49	\$299.72	\$299.72	\$283.94
U202H	Monofilament, absorbable	RB-1	27"	36/box	\$490.29	\$465.78	\$465.78	\$441.26

Continued on the next page

16



ETHICON CHROMIC GUT - continued

Size: 6-0

796G	Monofilament, absorbable	G-1	18"	12/box	\$404.79	\$384.55	\$384.55	\$364.31
1635G	Monofilament, absorbable	PS-3	18"	12/box	\$296.99	\$282.14	\$282.14	\$267.29
1641G	Monofilament, absorbable	PS-4	18"	12/box	\$288.79	\$274.35	\$274.35	\$259.91
1816G	Monofilament, absorbable	PS-6	18"	12/box	\$285.99	\$271.69	\$271.69	\$257.39

ETHICON PLAIN GUT

Size: 2-0

H823H	Monofilament, absorbable	FS-1	27"	36/box	\$476.29	\$452.48	\$452.48	\$428.66
-------	--------------------------	------	-----	--------	----------	----------	----------	----------

Size: 3-0

H810H	Monofilament, absorbable	FS-1	27"	36/box	\$524.29	\$498.08	\$498.08	\$471.86
H822H	Monofilament, absorbable	FS-2	27"	36/box	\$427.49	\$406.12	\$406.12	\$384.74
G322H	Monofilament, absorbable	SH	27"	36/box	\$495.29	\$470.53	\$470.53	\$445.76
612G	Monofilament, absorbable	X-1	18"	12/box	\$154.29	\$146.58	\$146.58	\$138.86

Size: 4-0

H821H	Monofilament, absorbable	FS-2	27"	36/box	\$489.79	\$465.30	\$465.30	\$440.81
1644G	Monofilament, absorbable	P-3	18"	12/box	\$311.99	\$296.39	\$296.39	\$280.79
1633G	Monofilament, absorbable	PS-4	18"	12/box	\$281.29	\$267.23	\$267.23	\$253.16

Size: 5-0

H820G	Monofilament, absorbable	FS-2	27"	12/box	\$153.99	\$146.29	\$146.29	\$138.59
658G	Monofilament, absorbable	P-2	18"	12/box	\$325.99	\$309.69	\$309.69	\$293.39
686G	Monofilament, absorbable	P-3	18"	12/box	\$289.79	\$275.30	\$275.30	\$260.81
1915G	Monofilament, fast absorbing	PC-1	18"	12/box	\$362.79	\$344.65	\$344.65	\$326.51

ETHICON VICRYL®

- Coated

Size: 2-0

J443H	Undyed braided, absorbable	FS-1	27"	36/box	\$332.29	\$315.68	\$315.68	\$299.06
-------	----------------------------	------	-----	--------	----------	----------	----------	----------

Size: 3-0

J442H	Undyed braided, absorbable	FS-1	27"	36/box	\$332.29	\$315.68	\$315.68	\$299.06
J423H	Undyed braided, absorbable	FS-2	27"	36/box	\$347.29	\$329.93	\$329.93	\$312.56
J683G	Undyed braided, absorbable	PS-1	18"	12/box	\$190.99	\$181.44	\$181.44	\$171.89
J497G	Undyed braided, absorbable	PS-2	18"	12/box	\$173.49	\$164.82	\$164.82	\$156.14
J215H	Undyed braided, absorbable	RB-1	27"	36/box	\$272.49	\$258.87	\$258.87	\$245.24
J416H	Undyed braided, absorbable	SH	27"	36/box	\$222.99	\$211.84	\$211.84	\$200.69
J458H	Undyed braided, absorbable	X-1	27"	36/box	\$301.79	\$286.70	\$286.70	\$271.61

Size: 4-0

J422H	Undyed braided, absorbable	FS-2	27"	36/box	\$340.29	\$323.28	\$323.28	\$306.26
J455H	Undyed braided, absorbable	J-1	27"	36/box	\$279.79	\$265.80	\$265.80	\$251.81
J504G	Undyed braided, absorbable	P-2	18"	12/box	\$226.79	\$215.45	\$215.45	\$204.11
J494G	Undyed braided, absorbable	P-3	18"	12/box	\$191.99	\$182.39	\$182.39	\$172.79
J835G	Undyed braided, absorbable	PC-1	18"	12/box	\$187.29	\$177.93	\$177.93	\$168.56
J845G	Undyed braided, absorbable	PC-3	18"	12/box	\$186.29	\$176.98	\$176.98	\$167.66
J823G	Undyed braided, absorbable	PC-5	18"	12/box	\$187.29	\$177.93	\$177.93	\$168.56
J496G	Undyed braided, absorbable	PS-2	18"	12/box	\$213.29	\$202.63	\$202.63	\$191.96
J507G	Undyed braided, absorbable	PS-4	18"	12/box	\$170.99	\$162.44	\$162.44	\$153.89
J214H	Undyed braided, absorbable	RB-1	27"	36/box	\$239.49	\$227.52	\$227.52	\$215.54
J218H	Undyed braided, absorbable	SH-1	27"	36/box	\$222.99	\$211.84	\$211.84	\$200.69

Continued on the next page

WOUND CLOSURE/SUTURE

ETHICON

REVERSE CUTTING NEEDLES

3/8 Circle
FS

FS-1

FS-2

1/2 Circle
X-1



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle
P-1

P-3

PS-2

PS-3

1/2 Circle
P-2

PS-4

PS-5



CONVENTIONAL CUTTING NEEDLES

3/8 Circle
PC-1

PC-3




TAPER POINT NEEDLES

1/2 Circle
RB-1

SH-1

SH



WOUND CLOSURE/SUTURE

ETHICON

REVERSE CUTTING NEEDLES

3/8 Circle
FS



FS-1



FS-2



1/2 Circle
X-1



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle
P-1



P-3



PS-2



PS-3



1/2 Circle
P-2



PS-4



PS-5



CONVENTIONAL CUTTING NEEDLES

3/8 Circle
PC-1



PC-3



TAPER POINT NEEDLES

1/2 Circle
RB-1



SH-1



SH



ETHICON VICRYL® - continued

• Coated

Size: 5-0

J421H	Undyed braided, absorbable	FS-2	27"	36/box	\$306.79	\$291.45	\$291.45	\$276.11
J490G	Undyed braided, absorbable	P-1	18"	12/box	\$217.29	\$206.43	\$206.43	\$195.56
J503G	Undyed braided, absorbable	P-2	18"	12/box	\$228.29	\$216.88	\$216.88	\$205.46
J493G	Undyed braided, absorbable	P-3	18"	12/box	\$236.29	\$224.48	\$224.48	\$212.66
J844G	Undyed braided, absorbable	PC-3	18"	12/box	\$236.29	\$224.48	\$224.48	\$212.66
J495G	Undyed braided, absorbable	PS-2	18"	12/box	\$176.79	\$167.95	\$167.95	\$159.11
J500G	Undyed braided, absorbable	PS-3	18"	12/box	\$215.29	\$204.53	\$204.53	\$193.76
J213H	Undyed braided, absorbable	RB-1	27"	36/box	\$248.29	\$235.88	\$235.88	\$223.46
J433H	Undyed braided, absorbable	TF	27"	36/box	\$145.99	\$138.69	\$138.69	\$131.39

Size: 6-0

J489G	Undyed braided, absorbable	P-1	18"	12/box	\$195.99	\$186.19	\$186.19	\$176.39
J492G	Undyed braided, absorbable	P-3	18"	12/box	\$196.49	\$186.67	\$186.67	\$176.84
J833G	Undyed braided, absorbable	PC-1	18"	12/box	\$239.29	\$227.33	\$227.33	\$215.36
J499G	Undyed braided, absorbable	PS-3	18"	12/box	\$218.29	\$207.38	\$207.38	\$196.46

Size: 7-0

J488G	Undyed braided, absorbable	P-1	18"	12/box	\$244.79	\$232.55	\$232.55	\$220.31
-------	----------------------------	-----	-----	--------	----------	----------	----------	----------

ETHICON VICRYL®

• Coated

Size: 2-0

J461H	Violet braided, absorbable	X-1	27"	36/box	\$301.79	\$286.70	\$286.70	\$271.61
-------	----------------------------	-----	-----	--------	----------	----------	----------	----------

Size: 3-0

J393H	Violet braided, absorbable	FS-2	18"	36/box	\$306.79	\$291.45	\$291.45	\$276.11
J316H	Violet braided, absorbable	SH	27"	36/box	\$242.29	\$230.18	\$230.18	\$218.06
J460H	Violet braided, absorbable	X-1	27"	36/box	\$301.79	\$286.70	\$286.70	\$271.61

Size: 4-0

J386H	Violet braided, absorbable	C-3	18"	36/box	\$302.99	\$287.84	\$287.84	\$272.69
J451H	Violet braided, absorbable	FS-1	27"	36/box	\$318.79	\$302.85	\$302.85	\$286.91
J392H	Violet braided, absorbable	FS-2	18"	36/box	\$296.99	\$282.14	\$282.14	\$267.29
J464G	Violet braided, absorbable	P-3	18"	12/box	\$191.99	\$182.39	\$182.39	\$172.79
J304H	Violet braided, absorbable	RB-1	27"	36/box	\$227.79	\$216.40	\$216.40	\$205.01
J315H	Violet braided, absorbable	SH	27"	36/box	\$231.49	\$219.92	\$219.92	\$208.34
J310H	Violet braided, absorbable	SH-1	27"	36/box	\$241.29	\$229.23	\$229.23	\$217.16

Size: 5-0

J385H	Violet braided, absorbable	C-3	18"	36/box	\$310.99	\$295.44	\$295.44	\$279.89
J391H	Violet braided, absorbable	FS-2	18"	36/box	\$306.79	\$291.45	\$291.45	\$276.11
J463G	Violet braided, absorbable	P-3	18"	12/box	\$193.79	\$184.10	\$184.10	\$174.41
J593G	Violet braided, absorbable	PS-5	18"	12/box	\$257.29	\$244.43	\$244.43	\$231.56

Size: 6-0

J489G	Violet braided, absorbable	P-1	18"	12/box	\$195.99	\$186.19	\$186.19	\$176.39
J510G	Violet braided, absorbable	PS-6	18"	12/box	\$239.49	\$227.52	\$227.52	\$215.54

Size: 7-0

J546G	Violet braided, absorbable	TG-140-8	18"	12/box	\$343.49	\$326.32	\$326.32	\$309.14
-------	----------------------------	----------	-----	--------	----------	----------	----------	----------

ETHICON VICRYL RAPIDE®

- Coated

Size: 3-0

VR497	Undyed braided, fast absorbing	PS-2	18"	12/box	\$344.79	\$327.55	\$327.55	\$310.31
-------	--------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 4-0

VR494	Undyed braided, fast absorbing	P-3	18"	12/box	\$302.29	\$287.18	\$287.18	\$272.06
VR496	Undyed braided, fast absorbing	PS-2	18"	12/box	\$293.79	\$279.10	\$279.10	\$264.41
VR426	Undyed braided, fast absorbing	PS-2	27"	12/box	\$340.29	\$323.28	\$323.28	\$306.26

Size: 5-0

VR493	Undyed braided, fast absorbing	P-3	18"	12/box	\$363.49	\$345.32	\$345.32	\$327.14
VR834	Undyed braided, fast absorbing	PC-1	18"	12/box	\$375.29	\$356.53	\$356.53	\$337.76

ETHICON PDS II®

- Clear

Size: 3-0

Z497G	Monofilament, absorbable	PS-2	18"	12/box	\$271.99	\$258.39	\$258.39	\$244.79
-------	--------------------------	------	-----	--------	----------	----------	----------	----------

Size: 4-0

Z496G	Monofilament, absorbable	PS-2	18"	12/box	\$259.79	\$246.80	\$246.80	\$233.81
-------	--------------------------	------	-----	--------	----------	----------	----------	----------

Size: 5-0

VR490	Monofilament, absorbable	P-1	18"	12/box	\$371.29	\$352.73	\$352.73	\$334.16
Z493G	Monofilament, absorbable	P-3	18"	12/box	\$306.79	\$291.45	\$291.45	\$276.11

ETHICON PDS II®

- Violet

Size: 6-0

Z487G	Monofilament, absorbable	P-1	18"	12/box	\$287.29	\$272.93	\$272.93	\$258.56
-------	--------------------------	-----	-----	--------	----------	----------	----------	----------

ETHICON MONOCRYL®

- Undyed

Size: 3-0

Y497G	Monofilament, absorbable	PS-2	18"	12/box	\$248.29	\$235.88	\$235.88	\$223.46
-------	--------------------------	------	-----	--------	----------	----------	----------	----------

Size: 4-0

Y494G	Monofilament, absorbable	P-3	18"	12/box	\$326.99	\$310.64	\$310.64	\$294.29
Y496G	Monofilament, absorbable	PS-2	18"	12/box	\$250.99	\$238.44	\$238.44	\$225.89

Size: 5-0

Y493G	Monofilament, absorbable	P-3	18"	12/box	\$271.99	\$258.39	\$258.39	\$244.79
Y844G	Monofilament, absorbable	PC-3	18"	12/box	\$271.99	\$258.39	\$258.39	\$244.79
Y495G	Monofilament, absorbable	PS-2	18"	12/box	\$255.29	\$242.53	\$242.53	\$229.76

Size: 6-0

Y492G	Monofilament, absorbable	P-3	18"	12/box	\$322.29	\$306.18	\$306.18	\$290.06
Y833G	Monofilament, absorbable	PC-1	18"	12/box	\$277.49	\$263.62	\$263.62	\$249.74

ETHICON MONOCRYL®

- Violet

Size: 5-0

Y463G	Monofilament, absorbable	P-3	18"	12/box	\$271.99	\$258.39	\$258.39	\$244.79
-------	--------------------------	-----	-----	--------	----------	----------	----------	----------

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

WOUND CLOSURE/SUTURE

ETHICON

REVERSE CUTTING NEEDLES

3/8 Circle

FS _____

FS-1 _____

FS-2 _____

1/2 Circle

X-1 _____



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle

P-1 _____

P-3 _____

PS-2 _____

PS-3 _____

1/2 Circle

P-2 _____

PS-4 _____

PS-5 _____



CONVENTIONAL CUTTING NEEDLES

3/8 Circle

PC-1 _____

PC-3 _____




TAPER POINT NEEDLES

1/2 Circle

RB-1 _____

SH-1 _____

SH _____



WOUND CLOSURE/SUTURE

ETHICON

REVERSE CUTTING NEEDLES

3/8 Circle

FS



FS-1



FS-2



1/2 Circle

X-1



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle

P-1



P-3



PS-2



PS-3



1/2 Circle

P-2



PS-4



PS-5



CONVENTIONAL CUTTING NEEDLES

3/8 Circle

PC-1



PC-3



TAPER POINT NEEDLES

1/2 Circle

RB-1



SH-1



SH



ETHICON ETHILON® NYLON

• Black

Size: 2-0

664G	Monofilament, non-absorbable	FS	18"	12/box	\$83.29	\$79.13	\$79.13	\$74.96
------	------------------------------	----	-----	--------	---------	---------	---------	---------

Size: 3-0

669H	Monofilament, non-absorbable	FS-1	30"	36/box	\$296.29	\$281.48	\$281.48	\$266.66
------	------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 4-0

662G	Monofilament, non-absorbable	FS-2	18"	12/box	\$99.29	\$94.33	\$94.33	\$89.36
------	------------------------------	------	-----	--------	---------	---------	---------	---------

699G	Monofilament, non-absorbable	P-3	18"	12/box	\$217.79	\$206.90	\$206.90	\$196.01
------	------------------------------	-----	-----	--------	----------	----------	----------	----------

1667G	Monofilament, non-absorbable	PS-2	18"	12/box	\$214.29	\$203.58	\$203.58	\$192.86
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 5-0

661G	Monofilament, non-absorbable	FS-2	18"	12/box	\$101.79	\$96.70	\$96.70	\$91.61
------	------------------------------	------	-----	--------	----------	---------	---------	---------

698G	Monofilament, non-absorbable	P-3	18"	12/box	\$163.79	\$155.60	\$155.60	\$147.41
------	------------------------------	-----	-----	--------	----------	----------	----------	----------

1855G	Monofilament, non-absorbable	PC-1	18"	12/box	\$186.29	\$176.98	\$176.98	\$167.66
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

1865G	Monofilament, non-absorbable	PC-3	18"	12/box	\$191.99	\$182.39	\$182.39	\$172.79
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

1666G	Monofilament, non-absorbable	PS-2	18"	12/box	\$216.49	\$205.67	\$205.67	\$194.84
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

1668G	Monofilament, non-absorbable	PS-3	18"	12/box	\$176.79	\$167.95	\$167.95	\$159.11
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

1661G	Monofilament, non-absorbable	PS-4	18"	12/box	\$207.29	\$196.93	\$196.93	\$186.56
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 6-0

697G	Monofilament, non-absorbable	P-1	18"	12/box	\$210.79	\$200.25	\$200.25	\$189.71
------	------------------------------	-----	-----	--------	----------	----------	----------	----------

1698G	Monofilament, non-absorbable	P-3	18"	12/box	\$203.29	\$193.13	\$193.13	\$182.96
-------	------------------------------	-----	-----	--------	----------	----------	----------	----------

1856G	Monofilament, non-absorbable	PC-1	18"	12/box	\$230.29	\$218.78	\$218.78	\$207.26
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

1665G	Monofilament, non-absorbable	PS-3	18"	12/box	\$168.49	\$160.07	\$160.07	\$151.64
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 7-0

1696G	Monofilament, non-absorbable	P-1	18"	12/box	\$224.79	\$213.55	\$213.55	\$202.31
-------	------------------------------	-----	-----	--------	----------	----------	----------	----------

ETHICON PERMA-HAND® SILK

• Black

Size: 2-0

685G	Braided, non-absorbable	FS	18"	12/box	\$86.99	\$82.64	\$82.64	\$78.29
------	-------------------------	----	-----	--------	---------	---------	---------	---------

737G	Braided, non-absorbable	X-1	18"	12/box	\$81.29	\$77.23	\$77.23	\$73.16
------	-------------------------	-----	-----	--------	---------	---------	---------	---------

Size: 3-0

684G	Braided, non-absorbable	FS-2	18"	12/box	\$86.29	\$81.98	\$81.98	\$77.66
------	-------------------------	------	-----	--------	---------	---------	---------	---------

632G	Braided, non-absorbable	X-1	18"	12/box	\$72.79	\$69.15	\$69.15	\$65.51
------	-------------------------	-----	-----	--------	---------	---------	---------	---------

Size: 4-0

736G	Braided, non-absorbable	C-3	18"	12/box	\$92.99	\$88.34	\$88.34	\$83.69
------	-------------------------	-----	-----	--------	---------	---------	---------	---------

629G	Braided, non-absorbable	FS-1	18"	12/box	\$83.29	\$79.13	\$79.13	\$74.96
------	-------------------------	------	-----	--------	---------	---------	---------	---------

683G	Braided, non-absorbable	FS-2	18"	12/box	\$75.79	\$72.00	\$72.00	\$68.21
------	-------------------------	------	-----	--------	---------	---------	---------	---------

734G	Braided, non-absorbable	J-1	18"	12/box	\$99.29	\$94.33	\$94.33	\$89.36
------	-------------------------	-----	-----	--------	---------	---------	---------	---------

641G	Braided, non-absorbable	P-3	18"	12/box	\$172.99	\$164.34	\$164.34	\$155.69
------	-------------------------	-----	-----	--------	----------	----------	----------	----------

1677G	Braided, non-absorbable	PS-2	18"	12/box	\$178.49	\$169.57	\$169.57	\$160.64
-------	-------------------------	------	-----	--------	----------	----------	----------	----------

K871H	Braided, non-absorbable	RB-1	30"	36/box	\$235.99	\$224.19	\$224.19	\$212.39
-------	-------------------------	------	-----	--------	----------	----------	----------	----------

Size: 5-0

682G	Braided, non-absorbable	FS-2	18"	12/box	\$83.29	\$79.13	\$79.13	\$74.96
------	-------------------------	------	-----	--------	---------	---------	---------	---------

640G	Braided, non-absorbable	P-3	18"	12/box	\$173.99	\$165.29	\$165.29	\$156.59
------	-------------------------	-----	-----	--------	----------	----------	----------	----------

Size: 6-0

639G	Braided, non-absorbable	P-1	18"	12/box	\$175.79	\$167.00	\$167.00	\$158.21
------	-------------------------	-----	-----	--------	----------	----------	----------	----------

ETHICON PROLENE®

- Blue

Size: 3-0

8665G	Monofilament, non-absorbable	FS-2	18"	12/box	\$151.99	\$144.39	\$144.39	\$136.79
-------	------------------------------	------	-----	--------	----------	----------	----------	----------

Size: 4-0

8699G	Monofilament, non-absorbable	P-3	18"	12/box	\$253.29	\$240.63	\$240.63	\$227.96
8682G	Monofilament, non-absorbable	PS-2	18"	12/box	\$250.99	\$238.44	\$238.44	\$225.89

Size: 5-0

8698G	Monofilament, non-absorbable	P-3	18"	12/box	\$219.79	\$208.80	\$208.80	\$197.81
8681G	Monofilament, non-absorbable	PS-3	18"	12/box	\$220.79	\$209.75	\$209.75	\$198.71
8661G	Monofilament, non-absorbable	FS-2	18"	12/box	\$126.29	\$119.98	\$119.98	\$113.66

Size: 6-0

8697G	Monofilament, non-absorbable	P-1	18"	12/box	\$319.79	\$303.80	\$303.80	\$287.81
8636G	Monofilament, non-absorbable	PC-3	18"	12/box	\$277.49	\$263.62	\$263.62	\$249.74
8695G	Monofilament, non-absorbable	P-3	18"	12/box	\$272.99	\$259.34	\$259.34	\$245.69
8680G	Monofilament, non-absorbable	PS-3	18"	12/box	\$303.79	\$288.60	\$288.60	\$273.41

ETHICON MERSILENE POLYESTER FIBER

- White

Size: 4-0

R633H	Braided, non-absorbable	FS-2	18"	36/box	\$212.99	\$202.34	\$202.34	\$191.69
-------	-------------------------	------	-----	--------	----------	----------	----------	----------



WOUND CLOSURE/SUTURE

ETHICON

REVERSE CUTTING NEEDLES

- 3/8 Circle
FS
- FS-1
- FS-2
- 1/2 Circle
X-1

PRECISION REVERSE CUTTING NEEDLES

- 3/8 Circle
P-1
- P-3
- PS-2
- PS-3
- 1/2 Circle
P-2
- PS-4
- PS-5

CONVENTIONAL CUTTING NEEDLES

- 3/8 Circle
PC-1
- PC-3

TAPER POINT NEEDLES

- 1/2 Circle
RB-1
- SH-1
- SH

WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle

C-6



C-7

1/2 Circle

C-9




PRECISION REVERSE CUTTING NEEDLES

3/8 Circle

C-1



C-3



A-TITAN CHROMIC GUT

- Absorbable
- 12/box

			Comparable to ETHICON				
Size: 3-0							
6273CG	C-6	27"	636H	\$83.49	\$79.32	\$79.32	\$75.14
Size: 4-0							
3184CG	C-3	18"	1654G	\$182.79	\$173.65	\$173.65	\$164.51
6274CG	C-6	27"	635H	\$83.49	\$79.32	\$79.32	\$75.14
Size: 5-0							
3185CG	C-3	18"	687G	\$182.79	\$173.65	\$173.65	\$164.51

A-TITAN PLAIN GUT

- Absorbable
- 12/box

			Comparable to ETHICON				
Size: 3-0							
6273PG	C-6	27"	H822H	\$83.49	\$79.32	\$79.32	\$75.14
7273PG	C-7	27"	H810H	\$83.49	\$79.32	\$79.32	\$75.14
9183PG	C-9	18"	612G	\$83.49	\$79.32	\$79.32	\$75.14
Size: 4-0							
6184PG	C-6	18"		\$83.49	\$79.32	\$79.32	\$75.14
6274PG	C-6	27"	H821H	\$83.49	\$79.32	\$79.32	\$75.14
Size: 5-0							
3185PG	C-3	18"	686G	\$182.79	\$173.65	\$173.65	\$164.51

A-TITAN PGA

- Undyed braided, absorbable
- 12/box

			Comparable to ETHICON				
Size: 3-0							
6273PGA	C-6	27"	J423H	\$182.79	\$173.65	\$173.65	\$164.51
31183PGA	C-31	18"	J497G	\$162.29	\$154.18	\$154.18	\$146.06
Size: 4-0							
3184PGA	C-3	18"	J494G	\$182.79	\$173.65	\$173.65	\$164.51
6274PGA	C-6	27"	J494G	\$96.79	\$91.95	\$91.95	\$87.11
31184PGA	C-31	18"	J496G	\$162.29	\$154.18	\$154.18	\$146.06
Size: 5-0							
3185PGA	C-3	18"	J493G	\$182.79	\$173.65	\$173.65	\$164.51
31185PGA	C-31	18"	J495G	\$182.79	\$173.65	\$173.65	\$164.51
Size: 6-0							
3186PGA	C-3	18"	J492G	\$182.79	\$173.65	\$173.65	\$164.51

16

WOUND CLOSURE/SUTURE

A-TITAN PGA

- Undyed braided, fast absorbing
- 12/box

Size: 3-0			Comparable to ETHICON				
31183PGAR	C-31	18"	VR497	\$192.29	\$182.68	\$182.68	\$173.06
Size: 4-0							
31184PGAR	C-31	18"	VR496	\$192.29	\$182.68	\$182.68	\$173.06

A-TITAN PGA

- Violet braided, absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
6183PGAV	C-6	18"	J393H	\$96.79	\$91.95	\$91.95	\$87.11
Size: 4-0							
6184PGAV	C-6	18"	J392H	\$96.79	\$91.95	\$91.95	\$87.11
Size: 5-0							
3185PGAV	C-3	18"	J463G	\$182.79	\$173.65	\$173.65	\$164.51
6185PGAV	C-6	18"	J391H	\$182.79	\$173.65	\$173.65	\$164.51

A-TITAN SILK

- Black braided, non-absorbable
- 12/box

Size: 3-0			Comparable to ETHICON				
6183BBS	C-6	18"		\$70.29	\$66.78	\$66.78	\$63.26
9183BBS	C-9	18"	632G	\$70.29	\$66.78	\$66.78	\$63.26
Size: 4-0							
6184BBS	C-6	18"	683G	\$70.29	\$66.78	\$66.78	\$63.26

A-TITAN POLYPROPYLENE

- Non-absorbable
- 12/box

Size: 5-0			Comparable to ETHICON				
3185PPN	C-3	18"		\$182.79	\$173.65	\$173.65	\$164.51
Size: 6-0							
1186PPN	C-1	18"	8697G	\$182.79	\$173.65	\$173.65	\$164.51
3186PPN	C-3	18"		\$182.79	\$173.65	\$173.65	\$164.51
Size: 7-0							
1187PPN	C-1	18"		\$182.79	\$173.65	\$173.65	\$164.51



REVERSE CUTTING NEEDLES

3/8 Circle
C-6

C-7

1/2 Circle
C-9



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle
C-1

C-3



WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES

3/8 Circle

C1 _____

C6 _____

C17 _____

1/2 Circle

C9 _____



HU-FRIEDY PERMA SHARP™ CHROMIC GUT

- Absorbable

Size: 6-0

PSN1816C	C1	18"	12/box	\$172.49	\$163.87	\$163.87	\$155.24
----------	----	-----	--------	----------	----------	----------	----------

HU-FRIEDY PERMA SHARP™ PLAIN GUT

- Absorbable

Size: 3-0

PSNH822A	C6	27"	12/box	\$77.29	\$73.43	\$73.43	\$69.56
----------	----	-----	--------	---------	---------	---------	---------

PSN612A	C9	18"	12/box	\$78.49	\$74.57	\$74.57	\$70.64
---------	----	-----	--------	---------	---------	---------	---------

Size: 4-0

PSNH821A	C6	27"	12/box	\$77.29	\$73.43	\$73.43	\$69.56
----------	----	-----	--------	---------	---------	---------	---------

Size: 6-0

PSN774A	C1	18"	12/box	\$172.49	\$163.87	\$163.87	\$155.24
---------	----	-----	--------	----------	----------	----------	----------

HU-FRIEDY PERMA SHARP™ PGA

- Absorbable

Size: 3-0

PSN393V	Violet braided	C6	18"	12/box	\$88.49	\$84.07	\$84.07	\$79.64
---------	----------------	----	-----	--------	---------	---------	---------	---------

Size: 4-0

PSN392V	Violet braided	C6	18"	12/box	\$88.49	\$84.07	\$84.07	\$79.64
---------	----------------	----	-----	--------	---------	---------	---------	---------

PSN397V	Violet braided	C6	27"	12/box	\$88.49	\$84.07	\$84.07	\$79.64
---------	----------------	----	-----	--------	---------	---------	---------	---------

HU-FRIEDY PERMA SHARP™ SILK

- Non-absorbable

Size: 3-0

PSN772S	Black braided	C6	18"	12/box	\$64.29	\$61.08	\$61.08	\$57.86
---------	---------------	----	-----	--------	---------	---------	---------	---------

PSN632S	Black braided	C9	18"	12/box	\$65.29	\$62.03	\$62.03	\$58.76
---------	---------------	----	-----	--------	---------	---------	---------	---------

Size: 4-0

PSN683S	Black braided	C6	18"	12/box	\$64.29	\$61.08	\$61.08	\$57.86
---------	---------------	----	-----	--------	---------	---------	---------	---------

PSN736S	Black braided	C17	18"	12/box	\$64.29	\$61.08	\$61.08	\$57.86
---------	---------------	-----	-----	--------	---------	---------	---------	---------

Size: 5-0

PSN682S	Black braided	C6	18"	12/box	\$64.29	\$61.08	\$61.08	\$57.86
---------	---------------	----	-----	--------	---------	---------	---------	---------

HU-FRIEDY PERMA SHARP™ POLYPROPYLENE

- Non-absorbable

Size: 4-0

PSN8683P	C6	18"	12/box	\$91.99	\$87.39	\$87.39	\$82.79
----------	----	-----	--------	---------	---------	---------	---------

Size: 5-0

PSN8661P	C6	18"	12/box	\$91.99	\$87.39	\$87.39	\$82.79
----------	----	-----	--------	---------	---------	---------	---------

Size: 6-0

PSN8697P	C1	18"	12/box	\$174.79	\$166.05	\$166.05	\$157.31
----------	----	-----	--------	----------	----------	----------	----------

**RIVERPOINT MEDICAL MONOTEX® DPTFE
POLYTETRAFLUOROETHYLENE SURGICAL SUTURE**

MONOTEX is an ultra-soft, high density surgical suture designed with low surface porosity. The monofilament construction helps reduce bacterial migration along the suture line.

As a soft monofilament, MONOTEX is an excellent suture choice for dental bone grafting and implant procedures. It's clean, comfortable and has little to no package memory, which results in excellent handling and knot security.

- White

Size: 2-0

T8664	Monofilament, Non-Absorbable	FS-2	18"	12/box	\$104.79	\$99.55	\$99.55	\$94.31
-------	------------------------------	------	-----	--------	-----------------	----------------	----------------	----------------

Size: 3-0

T8665	Monofilament, Non-Absorbable	FS-2	18"	12/box	\$104.79	\$99.55	\$99.55	\$94.31
T8000	Monofilament, Non-Absorbable	FS-3	18"	12/box	\$104.79	\$99.55	\$99.55	\$94.31

Size: 4-0

T8001	Monofilament, Non-Absorbable	FS-3	18"	12/box	\$104.79	\$99.55	\$99.55	\$94.31
T8604	Monofilament, Non-Absorbable	P-3	18"	12/box	\$110.79	\$105.25	\$105.25	\$99.71
T8682	Monofilament, Non-Absorbable	PS-2	18"	12/box	\$110.79	\$105.25	\$105.25	\$99.71
T8204	Monofilament, Non-Absorbable	TF	18"	12/box	\$104.79	\$99.55	\$99.55	\$94.31

Size: 5-0

T8698	Monofilament, Non-Absorbable	P-3	18"	12/box	\$110.79	\$105.25	\$105.25	\$99.71
-------	------------------------------	-----	-----	--------	-----------------	-----------------	-----------------	----------------


Size: 6-0

T8695	Monofilament, Non-Absorbable	P-3	18"	12/box	\$110.79	\$105.25	\$105.25	\$99.71
-------	------------------------------	-----	-----	--------	-----------------	-----------------	-----------------	----------------

WOUND CLOSURE/SUTURE



REVERSE CUTTING NEEDLES



ROUND NEEDLES



LOOK™ SILK SPOOLS

- Non-sterile
- 100 Yd spools

Size: 2-0

SP118	Black braided, non-absorbable		each	\$34.99	\$33.24	\$33.24	\$31.49
-------	-------------------------------	--	------	----------------	----------------	----------------	----------------

Size: 3-0

SP117	Black braided, non-absorbable		each	\$28.49	\$27.07	\$27.07	\$25.64
-------	-------------------------------	--	------	----------------	----------------	----------------	----------------

Size: 4-0

SP116	Black braided, non-absorbable		each	\$28.49	\$27.07	\$27.07	\$25.64
-------	-------------------------------	--	------	----------------	----------------	----------------	----------------



LOOK™ UMBILICAL TAPE

- Non-sterile
- Convenient dispensing jars
- Strong, soft and silky braided, mercerized thread

UT75	1/8"		20 yards/jar	\$8.29	\$7.88	\$7.88	\$7.46
------	------	--	--------------	---------------	---------------	---------------	---------------

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

Introducing the ACE SOUTHERN Alliance Program (ASAP)

When you create an alliance with **ACE SOUTHERN** you have a committed partner in your success. Part of that commitment is to reward you **for your loyalty with significant savings on every dollar you spend.**

As an ASAP Qualified Customer, you will benefit from savings and special promotions across all the products **ACE SOUTHERN** offers to help your practice prosper.

	Premier \$15K Annual Spend	Premier Plus \$25K Annual Spend*
5% savings on all Rx, Supplies and Equipment	●	
10% savings on all Instruments	●	●
15-20% savings on Bone Grafting Materials and Implants	●	●
10% savings on all Rx, Supplies and Equipment		●

**AAOMS/AAP members only*

Become an ASAP Member by Enrolling Today • acesouthern.com/ASAP

I'm ready to buy.
Which Price Tier Column Do I Use?

Regular Price

Most customers fall into this category. It's not "suggested retail", but reflects the strong buying power of **ACE SOUTHERN**, which we pass along to you as **everyday low prices.**

Member Price

In partnership with the **AAOMS** and **AAP** societies, their members are given a special 5% discount on general merchandise, including Rx, supplies and equipment.

Premier Price

As shown above, **enrolled members** who have met the required annual spend of \$15,000 are qualified to receive the **Premier Discount.**

Premier+ Price

As a special reward for our **enrolled society members**, once they have met the required annual spend of \$25,000 they will receive the **Premier+ Discount.**

WOUND CLOSURE/SUTURE

3M STERI-STRIP™ SKIN CLOSURES

Reinforced

- Breathable, microporous material
- Hypoallergenic adhesive
- Provides better cosmetic outcome than sutures or staples
- Permits early staple/suture removal
- Reduces the risk of infection
- Non-woven, rayon-backed material reinforced with filaments for strength
- Provides wound support and assists in increasing the tensile strength of the wound

1/4" x 3" Strips

R1541	White	3 strips/pouch, 50 pouches/bx	\$59.49	\$56.52	\$56.52	\$53.54
-------	-------	-------------------------------	---------	---------	---------	---------

1/4" x 4" Strips

R1546	White	10 strips/pouch, 50 pouches/bx	\$92.79	\$88.15	\$88.15	\$83.51
-------	-------	--------------------------------	---------	---------	---------	---------

1/2" x 4" Strips

R1547	White	6 strips/pouch, 50 pouches/bx	\$92.79	\$88.15	\$88.15	\$83.51
-------	-------	-------------------------------	---------	---------	---------	---------



CONMED REFLEX® ONE SKIN STAPLERS

- "Best In Class" ergonomics
- Audible and tactile feedback
- Visible counter
- 35 staples/gun

8035	Regular staples	6 guns/bx	\$63.29	\$60.13	\$60.13	\$56.96
8535	Wide staples	6 guns/bx	\$63.29	\$60.13	\$60.13	\$56.96



3M PRECISE™ MULTI-SHOT DS DISPOSABLE SKIN STAPLER, DS-25

Designed to be a versatile, reliable and cost-effective wound closure tool.

- Quick and easy to use
- Elimination of cross-hatching
- Ease of removal for excellent cosmetic results
- Decreased incidence of infection
- No waste - able to select the number of staples required for each incision
- Reliable performance

3MDS25		25 staples/gun, 12 guns/bx	\$148.99	\$141.54	\$141.54	\$134.09
--------	--	----------------------------	----------	----------	----------	----------



ETHICON DERMABOND ADVANCED™ TOPICAL SKIN ADHESIVE

DERMABOND ADVANCED is intended for topical application only to hold closed easily approximated skin edges of wounds from surgical incisions, including punctures from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations. DERMABOND ADVANCED may be used in conjunction with, but not in place of, deep dermal stitches.

0.7ml

DNX6		6/bx	\$371.29	\$352.73	\$352.73	\$334.16
DNX12		12/bx	\$686.29	\$651.98	\$651.98	\$617.66


Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

WOUND CLOSURE/SUTURE



Regular Viscosity Multiuse Kit



Straight

Angled



High Viscosity Multiuse Kit



Preloaded Applicators

GLUSTITCH PERIACRYL® 90 TISSUE ADHESIVE

Regular Viscosity

- Effective tissue adhesive which is hemostatic and bacteriostatic
- Excellent alternative to conventional sutures in soft-tissue surgery or to conventional sutures
- Used in periodontal surgical procedures including gingivectomy, mucogingival flaps, biopsies, routine extractions and superficial oral ulcerations
- 2ml and 5ml kits include autoclavable trays

PeriAcryl® Violet

PACRYLUV	0.2ml Preloaded applicators	12/box	\$105.99	\$100.69	\$100.69	\$95.39
PACRYL9V	9ml Kit Includes: 3x3 ml bottles, tray, 70 regular pipettes, 20 angled pipettes	kit	\$189.79	\$180.30	\$180.30	\$170.81
PACRYL5V	5ml Multiuse Kit Includes: 5ml MDV, 50 pipettes	kit	\$111.79	\$106.20	\$106.20	\$100.61
PACRYL2V	2ml Mini Kit Includes: 2ml MDV, 20 pipettes	kit	\$73.79	\$70.10	\$70.10	\$66.41

PeriAcryl® Clear

PACRYLUC	0.2ml Preloaded applicators	12/box	\$118.79	\$112.85	\$112.85	\$106.91
PACRYL5C	5ml Multiuse Kit Includes: 5ml MDV, 50 pipettes	kit	\$111.79	\$106.20	\$106.20	\$100.61
PACRYL2	2ml Mini Kit Includes: 2ml MDV, 20 pipettes	kit	\$73.79	\$70.10	\$70.10	\$66.41

GLUSTITCH PERIACRYL® 90 TISSUE ADHESIVE

High Viscosity

- 9 times thicker than regular viscosity PeriAcryl® 90
- Developed using feedback from periodontists
- Increased control for precision application and less waste
- Same trusted set time and strength as PeriAcryl® 90 standard viscosity
- Same user friendly delivery system as PeriAcryl® 90 standard viscosity
- 5ml and 9ml kits include autoclavable trays

PeriAcryl® Violet

PACRYL5VHV	5ml HV Multiuse Kit Includes: 5ml MDV, 50 pipettes	kit	\$118.99	\$113.04	\$113.04	\$107.09
PACRYL9VHV	9ml HV Multiuse Kit Includes: 3x3 ml bottles, tray, 70 regular pipettes, 20 angled pipettes	kit	\$202.79	\$192.65	\$192.65	\$182.51

GLUSTITCH PIPETTES FOR PERIACRYL® MULTIUSE KITS*

PIPETTE	Straight *	100/pack	\$17.99	\$17.09	\$17.09	\$16.19
PIPETTESANGLED	Angled *	50/pack	\$20.79	\$19.75	\$19.75	\$18.71

*Pipettes are empty and do not contain product

16

WOUND MANAGEMENT

JOHNSON & JOHNSON SURGIFOAM™

Absorbable Gelatin Oral Sponge

- Use in oral surgery for hemostasis to fill the dead space created by simple extraction, root amputation and removal of cysts, tumors and impacted teeth
- Rapid hemostasis
- Easy to handle: Compressible, does not require any cutting
- Absorbs up to 40 times its own weight
- Bioresorbable
- Bovine free



1969	Absorbable Gelatin Oral Sponge	24/box	\$155.99	\$148.19	\$148.19	\$140.39
------	--------------------------------	--------	----------	----------	----------	----------

Absorbable Gelatin Sponge

- A sterile, water-insoluble, malleable, porcine gelatin absorbable sponge intended for hemostatic use by applying to a bleeding surface
- Bovine free

1972	Size 12-7, 2cm x 6cm x 7mm	12/box	\$143.29	\$136.13	\$136.13	\$128.96
1973	Size 50, 8cm x 6.25cm x 10mm	4/box	\$157.79	\$149.90	\$149.90	\$142.01
1974	Size 100, 8cm x 12.5cm x 10mm	6/box	\$326.99	\$310.64	\$310.64	\$294.29
1975	Size 100C, 8cm x 12.5cm x 2mm (compressed)	6/box	\$301.79	\$286.70	\$286.70	\$271.61



COREVA HEALTH SCIENCE ACTCEL® HEMOSTATIC GAUZE

- Soluble hemostatic dressing
- Easy-to-use
- Quick clotting time
- Sterile
- Hypoallergenic
- Made from specially treated cellulose

7000228	2" x 2"	20/box	\$119.99	\$113.99	\$113.99	\$107.99
7000440	4" x 4"	20/box	\$206.99	\$196.64	\$196.64	\$186.29

JOHNSON & JOHNSON SURGICEL® ABSORBABLE HEMOSTAT

- Flexible, sheer weave conforms to irregular surfaces
- Excellent for an especially wide variety of procedures

1955	SURGICEL® 1/2" x 2"	12/box	\$569.79	\$541.30	\$541.30	\$512.81
1953	SURGICEL® 2" x 3"	12/box	\$1,245.29	\$1,183.03	\$1,183.03	\$1,120.76
1952	SURGICEL® 4" x 8"	12/box	\$1,886.79	\$1,792.45	\$1,792.45	\$1,698.11
1961	SURGICEL® 1" x 2", Fibrillar	each	\$147.49	\$140.12	\$140.12	\$132.74



INTEGRA™ HELISTAT® ABSORBABLE COLLAGEN HEMOSTAT

Helistat Absorbable Collagen Hemostatic Sponge is a soft, white, pliable, non-friable absorbent sponge. Because of its non-friable coherent sponge structure, the application of Helistat hemostatic sponge to the site where hemostasis is desired is easily controlled. Unwanted dispersal over the operative site is not encountered.

Performance

- Designed to be totally absorbable if left in situ after hemostasis
- Promotes platelet aggregation
- Effectively controls bleeding within 2 to 5 minutes

Better Handling

- Readily placed with dry forceps
- Can be wrapped, packed, sutured or placed on a multitude of bleeding sites

Ease of Use

- Excellent adherence to bleeding site
- Can be used dry or wet and can be used with saline or Thrombin
- Can be cut to size when dry, allowing more efficient application
- May be easily removed without shedding, fragmenting or causing rebleeding at wound site

Compared to Gelfoam®, Helistat showed superior ability to mediate platelet aggregation, deposition and activation, and initiate gross clot formation¹. Helistat outperforms Avitene®, Gelfoam® and Surgicel®².



1690ZZ	Helistat Hemostatic Sponge - 0.5" x 1"	18/box	\$289.99	\$275.49	\$275.49	\$260.99
--------	--	--------	----------	----------	----------	----------

1. Comparative in Vitro Analysis of Topical Hemostatic Agents, Wagner et. al., JOURNAL OF SURGICAL RESEARCH 66, 100-108 (1996) ARTICLE NO. 0379.
2. Evaluation of Hemostatic Agents in Experimental Splenic Lacerations, Coln et al, American Journal of Surgery, Vol 145 Feb 1983.
3. Data on file

WOUND MANAGEMENT



INVOTEC® EQUISPON® ABSORBABLE GELATIN SPONGE USP

EquiSpon® absorbable hemostatic sponges are made of Pharmaceutical Gelatin and are indicated for use to control minimal bleeding by tamponade effect, blood absorption and platelet aggregation, and adhesion prevention following surgery.

- Superior hemostatic capabilities
- Controls minimal bleeding and oozing
- Degradable - no need for painful traumatic removal

EQU602007	Size 12-7, 60mm x 20mm x 7mm	10/box	\$129.99	\$123.49	\$123.49	\$116.99
EQU805010	Size 50, 80mm x 50mm x 10mm	10/box	\$289.99	\$275.49	\$275.49	\$260.99
EQU1258010	Size 100, 125mm x 80mm x 10mm	5/box	\$191.79	\$182.20	\$182.20	\$172.61

LIFESCIENCE PLUS BLOODSTOP® ABSORBENT HEMOSTAT

A natural, biocompatible, non-irritating, animal product-free material that adheres to a wound, absorbs blood, stops bleeding and forms a protective layer to create an environment for wound healing. Its water solubility allows it to be rinsed off without disruption to the clotted surface. Using a proprietary formulation, cotton cellulose is etherized and oxidized to make a highly absorbent, water-soluble hemostatic matrix.

- Acts in less than a minute
- Painless upon application and removal
- Easy to cut, pack or layer to fit the wound
- Single use sterile packaging; several sizes
- Active bleeding control over traditional gauzes
- Adheres to the wound without the discomfort of adhesives
- Safe for all age groups
- Biodegradable

BSIX27	BloodSTOP iX for Surgical Use - 0.5" x 2"	24/box	\$156.29	\$148.48	\$148.48	\$140.66
BSIX14	BloodSTOP iX for Surgical Use - 2" x 2"	12/box	\$269.29	\$255.83	\$255.83	\$242.36
BS10	BloodSTOP for External Use Only - 2" x 2"	20/box	\$125.29	\$119.03	\$119.03	\$112.76
BS12	BloodSTOP for External Use Only - 4" x 4"	10/box	\$148.29	\$140.88	\$140.88	\$133.46
BSMP18	BloodSTOP for External Use Only - 2cm x 2cm	80/box	\$59.29	\$56.33	\$56.33	\$53.36

ORCA ORAL PRODUCTS ORCA™ FOAM

ORCA™ Foam is a sterile, malleable and absorbable gelatin sponge used for rapid hemostatic purposes in oral surgery and dental procedures. It also works well as a carrier medium for treatment liquids such as ORCA Dry Socket or antibiotics.

- Absorbs up to 35 times its weight
- Completely reabsorbs within a week
- Made of water-soluble porcine gelatin
- Small size is ideal for most applications and requires no trimming
- Individually packaged to prevent cross-contamination

09915	ORCA™ Foam - 10mm x 10mm x 10mm	50/box	\$281.49	\$267.42	\$267.42	\$253.34
-------	---------------------------------	--------	-----------------	-----------------	-----------------	-----------------

DAVOL AVITENE™ MICROFIBRILLAR COLLAGEN HEMOSTAT

Avitene™ Microfibrillar Collagen Hemostat is an active collagen hemostat, proven to accelerate clot formation. Collagen accelerates clot formation by enhancing platelet aggregation and the release of proteins to form fibrin, resulting in hemostasis.

- Active absorbable topical hemostatic agent
- 100% collagen: potentiates patient's own clotting mechanism
- Trusted by surgeons for its safety and efficacy for over 40 years
- A collagen hemostat acceptable for all procedures where a topical hemostat is indicated

6101001	Avitene Flour - .5 gm	each	\$156.99	\$149.14	\$149.14	\$141.29
101002	Avitene Flour - 1 gm	each	\$304.79	\$289.55	\$289.55	\$274.31
101059	Avitene Flour - 5 gm	each	\$1,081.99	\$1,027.89	\$1,027.89	\$973.79

ACE COLLATENE™ MICRO FIBRILLAR COLLAGEN

Collatene is an absorbent microfibrillar collagen matrix intended for the management of oral wounds and sores.

- Protects the oral wound bed
- Fully resorbable
- Highly absorbent
- Controls minor bleeding due to intrinsic hemostatic properties of collagen fibers

5091010	Collatene Collagen - 0.1 gram vial	10/box	\$196.99	\$187.14	\$187.14	\$177.29
---------	------------------------------------	--------	-----------------	-----------------	-----------------	-----------------



17

WOUND MANAGEMENT

LOOK™ LUKENS™ BONE WAX

- All-natural softener
- Non-bleached beeswax
- Soft, easy to form and to apply
- Easy-to-handle, safe foil peel packaging

902	Lukens™ Bone Wax - 2.5gm	12/box	\$78.49	\$74.57	\$74.57	\$70.64
-----	--------------------------	--------	---------	---------	---------	---------



ETHICON BONE WAX

- A sterile mixture of beeswax, paraffin, and isopropyl palmitate, a wax-softening agent
- May be used for the control of bleeding from bone surfaces
- Achieves local hemostasis of bone by acting as a mechanical barrier; it does not act biochemically and is non-absorbable
- Should not be used where rapid osseous regeneration and fusion are desired

W31G	Ethicon Bone Wax - 2.5g	12/box	\$281.29	\$267.23	\$267.23	\$253.16
------	-------------------------	--------	----------	----------	----------	----------



CP MEDICAL BONE WAX

Surgical Bone Wax is a simple-to-use local hemostatic agent that controls bleeding on bone surfaces during surgeries and procedures. Bone Wax is a mixture of beeswax, paraffin, and isopropyl palmitate and is opaque in appearance.

CPB31A	CP Medical BONE WAX - 2.5gm	10/box	\$68.99	\$65.54	\$65.54	\$62.09
--------	-----------------------------	--------	---------	---------	---------	---------



ORASOothe® SOCKIT GEL

OraSoothe “Sockit” Gel is an FDA-approved, all-natural, drug-free oral wound hydrogel that comes in an easy-to-use syringe. It is applied in-office during procedures that involve any type of oral wound and is pleasant tasting, optimizes healing, and provides excellent pain relief. After the procedure, the patient takes home the remainder of the gel and applies as needed. Benefits include:

- Optimized healing of oral wounds
- Pain management alternative to prescriptions
- Reduced post-op visits

SI5	10g syringes	5/box	\$81.79	\$77.70	\$77.70	\$73.61
SI25	10g syringes, 1-7 bags	25/bag	\$275.99	\$262.15	\$262.15	\$248.36
	10g syringes, 8 or more bags	25/bag	\$248.99	\$236.54	\$236.54	\$224.09



ORCA ORAL PRODUCTS ORCA™ DRY SOCKET

ORCA™ Dry Socket is an all-natural, unmedicated liquid formulated with essential oils, aloe vera and other proprietary natural ingredients to relieve the symptoms of dry socket. Simply place 8 to 10 drops onto an ORCA™ Foam sponge or similar absorbable sponge, then apply directly to the oral areas experiencing dry socket discomfort.

- Over 95% of patients experience resolution with just one treatment
- Designed for oral surgery and dental applications
- Can also be used with a gauze strip (removal required)

10993	ORCA™ Dry Socket - 10ml bottle	each	\$127.99	\$121.59	\$121.59	\$115.19
-------	--------------------------------	------	----------	----------	----------	----------



ORCA ORAL PRODUCTS ORCA™ PAIN RELIEF PUTTY

- Control post-op pain naturally
- Non-opioid with no harmful chemicals
- Allow your patients to treat dry sockets or discomfort at home
- Convenient packaging with easy application
- Exceptional pricing for peace of mind and improved patient experience
- Same formulation as ORCA Dry Socket but geared for over-the-counter consumers

31099096	ORCA™ Putty	6/box	\$98.29	\$93.38	\$93.38	\$88.46
310990924	ORCA™ Putty	24/box	\$265.79	\$252.50	\$252.50	\$239.21



17

WOUND MANAGEMENT



EUGENOL USP

Eugenol USP is used in medicine as a local antiseptic and anesthetic. Eugenol can be combined with zinc oxide to form a material, known as zinc oxide eugenol, which has restorative and prosthodontic applications in dentistry. The USP grade indicates it is graded suitable for personal care, cosmetic and pharmaceutical applications.

0602077	Eugenol USP - 1 oz. bottle	each	\$8.79	\$8.35	\$8.35	\$7.91
12655	Eugenol USP - 4 oz. bottle	each	\$25.99	\$24.69	\$24.69	\$23.39



GC AMERICA COE-PAK™ PERIODONTAL DRESSING

- ADA Certified
- Does not contain eugenol
- Clean, easy-to-handle
- No burning sensation
- No unpleasant odor or taste
- Contains: Base (90 g) and Accelerator (90 g)

Standard Package - Ready to use in 3 minutes

GCA135001	Standard	each	\$91.49	\$86.92	\$86.92	\$82.34
-----------	----------	------	---------	---------	---------	---------

Hard & Fast Set - Ready to use in 1 minute

GCA135301	Hard & Fast Set	each	\$91.49	\$86.92	\$86.92	\$82.34
-----------	-----------------	------	---------	---------	---------	---------



PULPDENT PERIOCARE

Contains: 13 oz Tube Paste and 13 oz Tube Gel

3840672PC	Perio Care Kit	each	\$44.79	\$42.55	\$42.55	\$40.31
-----------	----------------	------	---------	---------	---------	---------

3M ESPE™ PERIDEX™ CHLORHEXIDINE GLUCONATE 0.12% ORAL RINSE

- Provides antimicrobial activity during oral rinsing
- Reduces certain aerobic and anaerobic bacteria from a range of 54-97% through six months use
- Shows no significant changes in bacterial resistance or adverse changes in the oral microbial ecosystem
- Herbal mint flavor

12134	4oz	each	\$6.99	\$6.64	\$6.64	\$6.29
12133	64oz pump bottle	each	\$41.49	\$39.42	\$39.42	\$37.34

XTTRIUM LABORATORIES, INC. CHLORHEXIDINE ORAL RINSE

For the prevention and treatment of redness, swelling and bleeding of the gums due to gingivitis. Peppermint flavor.

200116	0.12% pint	each	\$6.29	\$5.98	\$5.98	\$5.66
--------	------------	------	--------	--------	--------	--------

WOUND MANAGEMENT

ACE NEW GAUZE STERILE PACKING MATERIAL

- X-ray detectable filament
- Non resorbable
- Used in: antrum, boney cavities, soft tissue, 3rd molar sockets

0034020	Plain, 1" x 5-½ yards	bottle	\$11.29	\$10.73	\$10.73	\$10.16
---------	-----------------------	--------	---------	---------	---------	---------

DERMA SCIENCES PAK-ITS® STERILE GAUZE PACKING STRIPS

Can be used dry to absorb drainage in exudating wounds or soaked with a medication of choice. Pre-cut sizes suit different wound types. A secondary dressing is needed to hold the strips in place and to protect the wound from outside contaminants. The bottle has a special patent-pending strip delivery system that ensures packing material is kept clean and tangle free, with the end always easy to locate at the top of the bottle.

- Latex-free

59220	Plain, 1/2" x 5 yds/bottle	each	\$6.79	\$6.45	\$6.45	\$6.11
-------	----------------------------	------	--------	--------	--------	--------

DYNAREX® STERILE PACKING STRIPS

- Ideal for use as nasal packing, drainage wicks or for packing wounds
- Specially woven gauze will not ravel, lint or fray

DY3401	Plain, 1/4" x 5 yds/bottle	each	\$3.29	\$3.13	\$3.13	\$2.96
DY3402	Plain, 1/2" x 5 yds/bottle	each	\$3.29	\$3.13	\$3.13	\$2.96
3411	Iodoform, 1/4" x 5 yds/bottle	each	\$3.29	\$3.13	\$3.13	\$2.96
3412	Iodoform, 1/2" x 5 yds/bottle	each	\$3.49	\$3.32	\$3.32	\$3.14

CURAD® PACKING MATERIAL

- Packing Strips, Plain
- NOT X-Ray Detectable

NON255145	Plain, ¼" x 5 yards	each	\$4.79	\$4.55	\$4.55	\$4.31
-----------	---------------------	------	--------	--------	--------	--------



Compliance Centered



Our NABP Drug Distribution Accreditation ensures that the pharmaceuticals we ship are 100% authentic, protecting you from counterfeit drugs. And, our turnkey services will help keep you compliant.



WOUND MANAGEMENT

STERILE GAUZE

COVIDIEN™ KENDALL CURITY™ LAPAROTOMY SPONGES, X-RAY DETECTABLE

- Pre-washed 100% cotton, disposable
- Kerlix™ sponges have dual-colored radiopaque handles that remain visible when saturated with blood
- Latex-free
- Sterile



Tray pack

6022	18" x 18"	100/case	\$73.79	\$70.10	\$70.10	\$66.41
------	-----------	----------	---------	---------	---------	---------

Soft Touch Pack

6534	12" x 12"	200/case	\$106.99	\$101.64	\$101.64	\$96.29
------	-----------	----------	----------	----------	----------	---------

COVIDIEN™ KENDALL VISTEC™ X-RAY DETECTABLE SPONGES

- Unique dual-colored element combines blue radiopaque strand with white strand which is visible when sponge is saturated with blood
- Element is latex-free and runs through both sides of sponge to facilitate counting
- Precise, double vacuum delinting process assures low lint
- Packaged sterile in sturdy polystyrene tray with peel-back lid
- Banded 10/tray
- Sterile



Sterile, 12-ply

7318	4" x 8"	800/case	\$160.49	\$152.47	\$152.47	\$144.44
------	---------	----------	----------	----------	----------	----------

16-ply

7317	4" x 4"	1,280/case	\$160.29	\$152.27	\$152.27	\$144.26
------	---------	------------	----------	----------	----------	----------

COVIDIEN™ KENDALL CURITY™ STERILE GAUZE SPONGES

Sterile, 8-ply

6318	4" x 4"	10/tray	\$1.79	\$1.70	\$1.70	\$1.61
		128 trays/case	\$156.29	\$148.48	\$148.48	\$140.66

Sterile, 12-ply

6939	4" x 4"	10/tray	\$1.99	\$1.89	\$1.89	\$1.79
		128 trays/case	\$162.29	\$154.18	\$154.18	\$146.06



DYNAREX® STERILE GAUZE SPONGES

- 100% cotton gauze
- Twin packs
- Packaged in peel-down pouches for convenient, aseptic handling
- Pouches packed in open trays for easy access
- Highly absorbent
- Sterile



Sterile, 12-ply

DY3333	3" x 3"	25 packs/box	\$2.49	\$2.37	\$2.37	\$2.24
DY3343	4" x 4"	25 packs/box	\$3.99	\$3.79	\$3.79	\$3.59

DYNAREX® STERILE GAUZE PADS

Sterile, 12-ply, 1's

DY3353	3" x 3", 12-ply	100/box	\$6.79	\$6.79	\$6.45	\$6.11
DY3354	4" x 4", 12-ply	100/box	\$11.29	\$11.29	\$10.73	\$10.16



DYNAREX® STERILE GAUZE SPONGES

Sterile, 8-ply, 2's

3322	2" x 2"	3,000/case	\$101.29	\$96.23	\$96.23	\$91.16
3342	4" x 4"	50/tray	\$5.79	\$5.50	\$5.50	\$5.21



DYNAREX® NONWOVEN GENERAL USE SPONGES

Sterile, 4-ply, 2's

3362	2" x 2"	3,000/case	\$58.29	\$55.38	\$55.38	\$52.46
3363	3" x 3"	2,400/case	\$66.49	\$63.17	\$63.17	\$59.84
3364	4" x 4"	1,200/case	\$113.79	\$108.79	\$108.79	\$102.41



STERILE GAUZE

WOUND MANAGEMENT

**COVIDIEN™ TELFA™ “OUCHLESS”
NON-ADHERENT STERILE PADS**

- Soft, perforated film allows air to circulate
- Absorbent cotton pad allows fluid collection
- “Ouchless” Telfa dressing won’t disturb healing tissue by sticking to wound
- Perforated carton for easy office dispensing
- Individually wrapped

KD2132	3" x 4"	100/box	\$23.29	\$22.13	\$22.13	\$20.96
--------	---------	---------	---------	---------	---------	---------



COVIDIEN™ KERLIX™ ROLL GAUZE

- Sterile - 6-ply
- Prewashed and fluff dried with permanent crinkle wave
- Soft conformable dressing
- Superior for cleansing wounds, cushioning injured areas against further trauma and prepping
- Available in Curity CSR Pack for convenient aseptic handling
- Length is while stretched

6715	4-½" x 4-½ Yards	1/box	\$3.49	\$3.32	\$3.32	\$3.14
------	------------------	-------	--------	--------	--------	--------



DEROYAL THROAT PACKING

- With radiopaque strip
- 1/pack
- Sterile


31708	2" x 18"	50 packs/case	\$88.29	\$83.88	\$83.88	\$79.46
-------	----------	---------------	---------	---------	---------	---------



**JOHNSON & JOHNSON
CODMAN SURGICAL PATTIES**

- X-ray detectable
- Sterile

801407	1/2" x 3"	10/pkg	\$30.99	\$29.44	\$29.44	\$27.89
--------	-----------	--------	---------	---------	---------	---------



RICHMOND STERIPOCKET™

These sponges are pre-packaged sterile sponges and are ideal to distribute to patients after extractions and oral surgeries.

- Sterile
- 2 per pack

Cotton-Filled

300408	2" x 2"	300 packs/box	\$24.99	\$23.74	\$23.74	\$22.49
--------	---------	---------------	---------	---------	---------	---------

Non-Woven Rayon/Poly

RD303120	2" x 2"	200 packs/box	\$41.29	\$39.23	\$39.23	\$37.16
----------	---------	---------------	---------	---------	---------	---------



WOUND MANAGEMENT

NON-STERILE, NON-WOVEN GAUZE

COVIDIEN™ CURITY™ ALL-PURPOSE SPONGES

- Nonwoven
- Ideal for applying ointments, prepping, wiping needles, cleaning slides and more
- Rayon/polyester blend construction
- Cost-effective
- Highly absorbent



Non-sterile, 4-ply, 2's

9022	2" x 2"	8,000/case	\$99.29	\$94.33	\$94.33	\$89.36
9023	3" x 3"	4,000/case	\$117.99	\$112.09	\$112.09	\$106.19

CROSSTEX ULTRA GAUZE® NON-WOVEN SPONGES

- More absorbent than traditional sponges
- Maintains shape and texture during use
- Easily releases from incision/extraction sites
- Produces little or no lint
- Will not stick to surfaces
- Non-sterile
- 200/bag



4 ply

ENCNWU	2" x 2"	5,000/case	\$92.29	\$87.68	\$83.06	\$83.06
ENCNWUB	2" x 2"	8,000/case	\$143.79	\$136.60	\$129.41	\$129.41
ENC3NWU	3" x 3"	4,000/case	\$154.79	\$147.05	\$139.31	\$139.31
ENC4NWU	4" x 4"	2,000/case	\$123.79	\$117.60	\$111.41	\$111.41

DUKAL CORPORATION CLINISORB NON-WOVEN SPONGES

- Rayon/poly blend construction
- Highly absorbent
- Multi-purpose sponges ideal for general wound care
- Virtually non-linting
- Not made with natural rubber latex
- Non-sterile
- 200/bag



4 ply

DUK2104	4" x 4"	2,000/case	\$45.79	\$43.50	\$43.50	\$41.21
---------	---------	------------	----------------	----------------	----------------	----------------

DYNAREX NON-WOVEN SPONGES

- Ideal for wound dressing, prepping and scrubbing, general purpose cleansing
- Highly absorbent
- Less linting than standard gauze
- Synthetic
- Non-sterile
- 200/bag



4 ply

3252	2" x 2"	4,000/case	\$28.99	\$27.54	\$27.54	\$26.09
3253	3" x 3"	4,000/case	\$56.49	\$53.67	\$53.67	\$50.84
3254	4" x 4"	2,000/case	\$45.29	\$43.03	\$43.03	\$40.76

RICHMOND NON-WOVEN 100% COTTON SPONGES

This sponge features a smooth, textured surface that provides exceptional patient comfort and maximum absorbency. It is also lint-free, so no fibers are left in your patient's mouth.



- Non-sterile
- 200/bag

4 ply

300631	2" x 2"	5,000/case	\$81.79	\$77.70	\$77.70	\$73.61
300632	3" x 3"	4,000/case	\$116.29	\$110.48	\$110.48	\$104.66
300633	4" x 4"	2,000/case	\$91.99	\$87.39	\$87.39	\$82.79

RICHMOND NON-WOVEN RAYON/POLY SPONGES

This sponge is an economical and lint-free alternative. It will not leave fibers in dental work, nor leave lint residue on instruments when used for cleaning.



- Non-sterile
- 200/bag

4 ply

300634	2" x 2"	5,000/case	\$43.29	\$41.13	\$41.13	\$38.96
300635	3" x 3"	4,000/case	\$67.29	\$63.93	\$63.93	\$60.56
300636	4" x 4"	2,000/case	\$59.79	\$56.80	\$56.80	\$53.81

TIDI® OMNISORB NON-WOVEN SPONGES

- Closed weave
- Highly absorbent, low-linting
- For general cleaning, prepping and applying ointments
- Non-sterile
- 200/bag



4 ply

942022	2" x 2"	5,000/case	\$89.99	\$85.49	\$85.49	\$80.99
942033	3" x 3"	4,000/case	\$168.79	\$160.35	\$160.35	\$151.91
942044	4" x 4"	2,000/case	\$113.29	\$107.63	\$107.63	\$101.96

*All gauze and sponges cannot be returned if case is opened

NON-STERILE GAUZE

WOUND MANAGEMENT

ACE GAUZE SPONGES

- All gauze sponges
- For general use
- 200/pkg
- Woven

Non-sterile, 8-ply

0038208	2" x 2"	25 pkgs – 5,000/case	\$43.49	\$41.32	\$41.32	\$39.14
0038212	4" x 4"	20 pkgs – 4,000/case	\$110.79	\$105.25	\$105.25	\$99.71

Non-sterile, 12-ply

0038209	2" x 2"	40 pkgs – 8,000/case	\$72.79	\$69.15	\$69.15	\$65.51
0038211	3" x 3"	20 pkgs – 4,000/case	\$107.29	\$101.93	\$101.93	\$96.56
0038210	4" x 4"	10 pkgs – 2,000/case	\$78.49	\$74.57	\$74.57	\$70.64
0038218	8" x 4"	10 pkgs – 2,000/case	\$166.99	\$158.64	\$158.64	\$150.29

Non-sterile, 16-ply

0038217	4" x 4"	10 pkgs – 2,000/case	\$97.79	\$92.90	\$92.90	\$88.01
---------	---------	----------------------	---------	---------	---------	---------



COVIDIEN™ CURITY™ WOVEN GAUZE SPONGES

- Made from high-grade 100% US grown cotton.
- Exceeds USP Type VII gauze requirements
- Economical dressing in broad range of sizes and plies
- Well suited for variety of applications; wound dressing, general cleaning, minor prepping, wound packing and debriding
- Non-sterile sponges conveniently indexed in 10's for faster, easier handling
- Not cotton-filled

Non-sterile, 12-ply

2346	3" x 3"	4,000/case	\$120.49	\$114.47	\$114.47	\$108.44
2634	4" x 4"	2,000/case	\$103.29	\$98.13	\$98.13	\$92.96

COVIDIEN KENDALL POST-OP SPONGES

- Similar to J&J Topper™ Sponges

Non-sterile

KD1700	3" x 3"	4,000/case	\$144.99	\$137.74	\$137.74	\$130.49
--------	---------	------------	----------	----------	----------	----------

COVIDIEN™ KENDALL DERMACEA™ WOVEN GAUZE SPONGES

- Non-sterile
- 200/bag



8 ply

441217	4" x 4"	4,000/case	\$95.99	\$91.19	\$91.19	\$86.39
--------	---------	------------	---------	---------	---------	---------

12 ply

441219	3" x 3"	4,000/case	\$54.79	\$52.05	\$52.05	\$49.31
--------	---------	------------	---------	---------	---------	---------

COVIDIEN™ DERMACEA™ GAUZE SPONGES

The economy line DERMACEA™ Gauze Sponges offers clinicians a cost effective alternative for non-sterile 100% cotton gauze dressings.



Non-sterile, 12-ply

441215	4" x 4"	2,000/case	\$64.29	\$61.08	\$61.08	\$57.86
--------	---------	------------	---------	---------	---------	---------

DEFEND COTTON DENTAL ROLLS

#2 Medium

- 100% Cotton
- Soft Texture
- Highly Absorbent



CS0200	3/8" x 1-1/2"	2,000/box	\$19.79	\$18.80	\$18.80	\$17.81
--------	---------------	-----------	---------	---------	---------	---------

*All gauze and sponges cannot be returned if case is opened

Tier Pricing Details see page 2

Regular Price	Member Price	Premier Price	Premier+ Price
---------------	--------------	---------------	----------------

WOUND MANAGEMENT

NON-STERILE GAUZE

DYNAREX® GAUZE SPONGES

- 100% cotton gauze
- For a wide range of dressing, cleaning and prepping procedures
- No exposed raw edges - cut edges folded in
- Cotton filled sponges provide maximum absorption
- Non-sterile
- 200/bag



8 ply							
3222	2" x 2"	5,000/case	\$33.49	\$31.82	\$31.82	\$30.14	
3232	3" x 3"	4,000/case	\$59.99	\$56.99	\$56.99	\$53.99	
3242	4" x 4"	4,000/case	\$85.79	\$81.50	\$81.50	\$77.21	
12 ply							
3223	2" x 2"	8,000/case	\$68.49	\$65.07	\$65.07	\$61.64	
3233	3" x 3"	4,000/case	\$69.79	\$66.30	\$66.30	\$59.79	
3243	4" x 4"	2,000/case	\$60.79	\$57.75	\$57.75	\$54.71	
16 ply							
3244	4" x 4"	2,000/case	\$81.29	\$77.23	\$77.23	\$73.16	

DYNAREX® COTTON FILLED GAUZE SPONGES

Non-sterile, 8-ply

3249	2" x 2"	5,000/case	\$47.49	\$45.12	\$45.12	\$42.74	
------	---------	------------	---------	---------	---------	---------	--



RICHMOND ALL-GAUZE SPONGES

All-gauze sponges have a loose weave, giving them a textured, cheesecloth-like appearance and a slightly abrasive surface, which means they provide a good grip. These sponges can be confidently used in a wide variety of applications.

- Non-sterile
- 200/bag

8 ply							
300430	2" x 2"	5,000/case	\$48.79	\$46.35	\$46.35	\$43.91	
300422	4" x 4"	4,000/case	\$122.49	\$116.37	\$116.37	\$110.24	
12 ply							
300420	2" x 2"	5,000/case	\$66.49	\$63.17	\$63.17	\$59.84	
300421	3" x 3"	4,000/case	\$96.99	\$92.14	\$92.14	\$87.29	
300423	4" x 4"	2,000/case	\$91.29	\$86.73	\$86.73	\$82.16	



RICHMOND COTTON FILLED EXODONTIA SPONGES

This sponge looks similar to the all-gauze variety with its cheesecloth-like outer layer, but it contains a layer of pure cotton on the inside.

- Non-sterile
- 200/bag

300410	2" x 2"	5,000/case	\$66.49	\$63.17	\$63.17	\$59.84	
300311	3" x 3"	4,000/case	\$87.79	\$82.79	\$82.79	\$79.01	
300312	4" x 4"	2,000/case	\$104.29	\$99.08	\$99.08	\$93.86	



RICHMOND MULTIPLY™ SPONGES

- Maintains dry surface
- 8-ply layers provide thickness perfect for bite pressure
- Twice as absorbent as gauze, exodontia, or 4 ply non-woven sponges
- Non-sterile
- 100/bag

8 ply							
300637	2" x 2"	3,000/case	\$69.49	\$66.01	\$66.01	\$62.54	
300638	3" x 3"	2,500/case	\$102.99	\$97.84	\$97.84	\$92.69	
300639	4" x 4"	2,000/case	\$151.29	\$143.79	\$143.79	\$136.16	



TIDI POST-OP GAUZE

- Cotton-filled
- Bulk
- Space-Saver Pack

Non-sterile

908224	4" x 4"	2,000/case	\$134.79	\$128.05	\$128.05	\$121.31	
--------	---------	------------	----------	----------	----------	----------	--

UNIPACK™ UNIGAUGE™ SPONGES & GAUZE

8-ply Cotton-Filled Gauze

- Superior absorbency
- No chemical additives
- Maximum patient comfort
- Non-sterile general use
- 200/bag

UCG1022	2" x 2"	5,000/case	\$44.49	\$42.27	\$42.27	\$39.49	
---------	---------	------------	---------	---------	---------	---------	--

4-ply Non-Woven Sponge

- High absorbency
- Soft, smooth, strong and lint-free
- Helps protect and cushion wounds
- Non-sterile general use
- 200/bag

UNW4022	2" x 2"	5,000/case	\$29.29	\$27.83	\$27.83	\$26.36	
UNW3033	3" x 3"	4,000/case	\$38.79	\$36.85	\$36.85	\$34.91	
UNW3044	4" x 4"	2,000/case	\$34.79	\$33.05	\$33.05	\$31.31	



WOUND MANAGEMENT

**BSN MEDICAL COVERLET®
ADHESIVE DRESSING**

- Absorbs liquid faster than competitors**
- Absorbs 10x its own weight in water
- 360-degree adhesive keeps wounds clean
- Elastic fabric



Strips							
00230	¾" x 3"	100/box	\$7.79	\$7.40	\$7.40	\$7.01	

Round Spots							
00301	7/8"	100/box	\$8.99	\$8.54	\$8.54	\$8.09	

**COVIDIEN™ CURITY™
STRETCH BANDAGES**

- Provides softness, low lint, and high absorbency
- Moderate stretch for maximum clinical flexibility
- Cotton/polyester blend



2231	2" x 75"	12 rolls/box	\$8.79	\$8.35	\$8.35	\$7.91
------	----------	--------------	--------	--------	--------	--------

**DYNAREX®
ADHESIVE STRIPS**

- Latex-free
- Sterile



Flexible Fabric							
3612	1" x 3"	100/box	\$5.29	\$5.03	\$5.03	\$4.76	
DY3611	¾" x 3"	100/box	\$5.79	\$5.50	\$5.50	\$5.21	

Sheer Plastic							
3601	¾" x 3"	100/box	\$4.79	\$4.55	\$4.55	\$4.31	

**DYNAREX®
SHEER PLASTIC
ADHESIVE BANDAGE**



3607	7/8" spot	100/box	\$2.79	\$2.65	\$2.65	\$2.51
3602	1" x 3" strip	100/box	\$5.29	\$5.03	\$5.03	\$4.76

**JOHNSON & JOHNSON®
BAND-AIDS**



8 ply							
4434	¾" x 3" strip	100/box	\$7.29	\$6.93	\$6.93	\$6.56	
4444	1" x 3" strip	100/box	\$12.99	\$12.34	\$12.34	\$11.69	

TEMPO® ADHESIVE BANDAGES



- Hypoallergenic
- Latex-free
- Sterile

Plastic							
BAP343	¾" x 3" strip	100/box	\$2.49	\$2.37	\$2.37	\$2.24	
BAP013	1" x 3" strip	100/box	\$2.79	\$2.65	\$2.6	\$2.51	

Fabric							
BAF343	¾" x 3" strip	100/box	\$3.49	\$3.32	\$3.32	\$3.14	
BAF013	1" x 3" strip	100/box	\$3.79	\$3.60	\$3.60	\$3.41	

**3M™ STERI-STRIPS™
SKIN CLOSURES**



R1548	1" x 5", 4/pkg	25 pkg/box	\$109.79	\$104.30	\$104.30	\$98.81
-------	----------------	------------	----------	----------	----------	---------

MASTISOL® LIQUID ADHESIVE

A clear, non-irritating, and non-water-soluble liquid adhesive that secures the most difficult dressings for an extended period of time. Mastisol is 7 to 10 times stickier than Compound Tincture of Benzoin.



052348	2/3 cc vial	48/box	\$207.79	\$197.40	\$197.40	\$187.01
052306	2 oz bottle	each	\$46.29	\$43.98	\$43.98	\$41.66

DETACHOL® ADHESIVE REMOVER

A gentle, non-irritating, adhesive remover recommended for the removal of dressings, wound closure strips, tapes, tubes, and Mastisol® Liquid Adhesive.



051348	2/3 cc vial	48/box	\$121.99	\$115.89	\$115.89	\$109.79
051304	4 oz bottle	each	\$18.29	\$17.38	\$17.38	\$16.46

WOUND MANAGEMENT

TEMPO® CO-EASE COHESIVE BANDAGES

- Lightweight self adhesive bandage
- Extremely comfortable, non-slip construction
- Soft, permeable fabric
- Provides strong and controlled support
- Effortless tear ability
- Water resistant - sweat and water will not cause loosening
- Offers long term compression that won't constrict
- 12 vibrant colors available
- Boxes are an assortment of colors



CE20C	2" x 5 yds	36 rolls/box	\$37.29	\$35.43	\$35.43	\$33.56
CE30C	3" x 5 yds	24 rolls/box	\$33.29	\$31.63	\$31.63	\$29.96

Latex-Free

CE10LFC	1" x 5 yds	60 rolls/box	\$43.79	\$41.60	\$41.60	\$39.41
CE20LFC	2" x 5 yds	36 rolls/box	\$41.99	\$39.89	\$39.89	\$37.79
CE30LFC	3" x 5 yds	24 rolls/box	\$40.99	\$38.94	\$38.94	\$36.89

CURAD ORTHO-POROUS SPORTS ADHESIVE TAPE

- The strong woven cotton backing of this tape provides high tensile strength and support to highly active people
- Powerful, rubber-based adhesive is designed for a long-lasting hold with minimal residue
- Highly porous construction allows the skin to breathe and permits moisture vapor transmission
- Consistently low unwind tension throughout the roll makes tape easy to use



Transparent

NON260301	1"	12/box	\$19.79	\$18.80	\$18.80	\$17.81
NON260302	2"	6/box	\$19.79	\$18.80	\$18.80	\$17.81

KENDALL WET-PRUF® WATERPROOF TAPE

- Waterproof and oil resistant
- High tensile strength



3063	1/2" x 10 yards	24 rolls/box	\$33.29	\$31.63	\$31.63	\$29.96
3142	1" x 10 yards	12 rolls/box	\$33.29	\$31.63	\$31.63	\$29.96

DYNAREX® SURGICAL TAPE

- Latex-free
- 10 yard rolls



Paper						
DY3551	1/2"	24 rolls/box	\$6.29	\$5.98	\$5.98	\$5.66
DY3552	1"	12 rolls/box	\$6.29	\$5.98	\$5.98	\$5.66
DY3553	2"	6 rolls/box	\$6.29	\$5.98	\$5.98	\$5.66

Transparent						
3571	1/2"	24 rolls/box	\$12.79	\$12.15	\$12.15	\$11.51
3572	1"	12 rolls/box	\$12.79	\$12.15	\$12.15	\$11.51
3573	2"	6 rolls/box	\$12.79	\$12.15	\$12.15	\$11.51

Cloth						
3561	1/2"	24 rolls/box	\$13.49	\$12.82	\$12.82	\$12.14
3562	1"	12 rolls/box	\$13.49	\$12.82	\$12.82	\$12.14
3563	2"	6 rolls/box	\$13.49	\$12.82	\$12.82	\$12.14

3M™ MICROPORE™ SURGICAL TAPE

- Gentle to the skin
- Latex-free and hypoallergenic for sensitive patients
- Highly breathable to maintain skin integrity
- Holds well on damp skin for secure placement
- 10 yard rolls



Paper

15300	1/2"	24 rolls/box	\$10.99	\$10.44	\$10.44	\$9.89
15301	1"	12 rolls/box	\$10.99	\$10.44	\$10.44	\$9.89
15302	2"	6 rolls/box	\$10.99	\$10.44	\$10.44	\$9.89

3M™ TRANSPORE™ SURGICAL TAPE

- Easy, straight, bi-directional tear
- Porous for breathability
- Easy to handle with gloves
- Good adhesive to skin and tubing for secure placement
- Transparent for easy monitoring
- Latex-free and hypoallergenic for sensitive patients
- Water-resistant
- 10 yard rolls



Transparent

15270	1/2"	24 rolls/box	\$18.49	\$17.57	\$17.57	\$16.64
15271	1"	12 rolls/box	\$18.49	\$17.57	\$17.57	\$16.64
15272	2"	6 rolls/box	\$18.49	\$17.57	\$17.57	\$16.64
15273	3"	4 rolls/box	\$18.49	\$17.57	\$17.57	\$16.64

3M™ DURAPORE™ SURGICAL TAPE

- Silk-like
- High strength for durability
- Adheres well to dry skin
- Latex-free and hypoallergenic for sensitive patients
- Convenient bi-directional tear for easy tailoring
- Water-resistant
- 10 yard rolls



Cloth

15380	1/2"	24 rolls/box	\$16.29	\$15.48	\$15.48	\$14.66
15381	1"	12 rolls/box	\$16.29	\$15.48	\$15.48	\$14.66
15382	2"	6 rolls/box	\$16.29	\$15.48	\$15.48	\$14.66
15383	3"	4 rolls/box	\$16.29	\$15.48	\$15.48	\$14.66

3M™ MICROFOAM™ TAPE

- Foam hypoallergenic tape
- Latex-free



15281	1" x 5-1/2 yards	12/box	\$29.79	\$28.30	\$28.30	\$26.81
-------	------------------	--------	---------	---------	---------	---------

WOUND MANAGEMENT

**DYNAREX®
COTTON STOCKINETTES**

- Conforms to the shape of the body, preventing wrinkles from forming under the cast
- Manufactured of 100% cotton
- Non-sterile



3692	2" x 25 yd roll	each	\$9.29	\$8.83	\$8.83	\$8.36
3693	3" x 25 yd roll	each	\$12.29	\$11.68	\$11.68	\$11.06
3694	4" x 25 yd roll	each	\$17.79	\$16.90	\$16.90	\$16.01
3696	6" x 25 yd roll	each	\$22.49	\$21.37	\$21.37	\$20.24

**DYNAREX® RUBBER
ELASTIC BANDAGE**

- Washable
- Autoclavable
- Measurement when stretched



3662	2" x 5 yards	10 rolls/box	\$8.99	\$8.54	\$8.54	\$8.09
3663	3" x 5 yards	10 rolls/box	\$10.29	\$9.78	\$9.78	\$9.26
3664	4" x 5 yards	10 rolls/box	\$11.29	\$10.73	\$10.73	\$10.16

COVIDIEN™ ARGYLE™ PENROSE TUBING

- Contain latex
- Sterile
- Radiopaque



12" Flat

513002	7/16"	each	\$1.49	\$1.42	\$1.42	\$1.34
513200	1/2"	each	\$1.49	\$1.42	\$1.42	\$1.34
513408	5/8"	each	\$1.49	\$1.42	\$1.42	\$1.34
513606	3/4"	each	\$1.49	\$1.42	\$1.42	\$1.34
513804	1"	each	\$1.49	\$1.42	\$1.42	\$1.34

18" Flat

514604	7/16"	each	\$1.49	\$1.42	\$1.42	\$1.34
514802	1/2"	each	\$1.49	\$1.42	\$1.42	\$1.34
515007	5/8"	each	\$1.49	\$1.42	\$1.42	\$1.34
515205	3/4"	each	\$1.49	\$1.42	\$1.42	\$1.34
515403	1"	each	\$1.49	\$1.42	\$1.42	\$1.34
515601	1-1/4"	each	\$1.49	\$1.42	\$1.42	\$1.34
515809	1-1/2"	each	\$1.49	\$1.42	\$1.42	\$1.34

**MEDLINE LATEX-FREE SILICONE PENROSE
DRAINS FOR CLOSED WOUND DRAINAGE**

- Made out of medical-grade silicone that is soft, flexible and kink resistant to help decrease potential patient discomfort
- Internal tread design helps pull exudates out of the body
- Smooth exterior surface combined with ribbed inner tubing helps provide ease of insertion and placement and aids gravity fluid drainage
- Sterile
- X-ray detectable

1/4" x 18"

DYND50427	25/case	\$91.29	\$86.73	\$86.73	\$82.16
-----------	---------	---------	---------	---------	---------

1/2" x 18"

DYND50428	25/case	\$123.29	\$117.13	\$117.13	\$110.96
-----------	---------	----------	----------	----------	----------

1" x 18"

DYND50429	25/case	\$160.79	\$152.75	\$152.75	\$144.71
-----------	---------	----------	----------	----------	----------

TOPICAL TREATMENTS

Bacitracin Ointment

006001	14gm tube	each	\$2.79	\$2.65	\$2.65	\$2.51
006003	28gm tube	each	\$3.99	\$3.79	\$3.79	\$3.59
DY11175	4 oz. tube	each	\$19.79	\$18.80	\$18.80	\$17.81
DY11176	1 lb. jar	each	\$62.99	\$59.84	\$59.84	\$56.69
1161	1/32 oz. foil pack	144/box	\$44.79	\$42.55	\$42.55	\$40.31

Neosporin Ointment

073088	.5 oz tube	each	\$6.29	\$5.98	\$5.98	\$5.66
073087	1 oz tube	each	\$8.99	\$8.54	\$8.54	\$8.09

WOUND MANAGEMENT

TEMPO® INSTANT COLD COMPRESS



ICE56M	** 5" x 6"	50/case	\$24.99	\$23.74	\$23.74	\$22.49
--------	------------	---------	---------	---------	---------	---------

CARDINAL HEALTH™ INSTANT COLD PACKS

- Consistent and effective cold therapy
- Easy to activate – just squeeze and shake
- Can be applied directly to skin – no additional cover needed



103	** 5" x 5.5"	50/case	\$44.29	\$42.08	\$42.08	\$39.86
-----	--------------	---------	---------	---------	---------	---------

ACE INSTANT ICE PACK

Care begins in the office and continues at home



Disposable

0031192	** Junior, 5" x 7"	24/case	\$20.79	\$19.75	\$19.75	\$18.71
0031191	** Large, 6" x 9"	24/case	\$23.49	\$22.32	\$22.32	\$21.14

DYNAREX® INSTANT COLD PACKS



DY4511	** 4" x 5"	24/case	\$14.49	\$13.77	\$13.77	\$13.04
DY4512	** 5" x 9"	24/case	\$22.29	\$21.18	\$21.18	\$20.06

COLDSTAR INSTANT HEAT PACK



130104	** 6" x 9"	24/case	\$24.99	\$23.74	\$23.74	\$22.49
--------	------------	---------	---------	---------	---------	---------

DYNAREX® INSTANT HOT PACK



4516	** 5" x 9"	24/case	\$22.99	\$21.84	\$21.84	\$20.69
------	------------	---------	---------	---------	---------	---------

COLDSTAR REUSABLE COLD/HOT PACKS



070204	** 4" x 7"	24/case	\$24.99	\$23.74	\$23.74	\$22.49
070104	** 6" x 9"	24/case	\$29.49	\$28.02	\$28.02	\$26.54

COVIDIEN™ UNI-PATCH™ REUSABLE HOT/COLD GEL PACKS

- Store in freezer for immediate cold therapy
- May be heated in hot water or microwaved in water for hot therapy
- Thick and flexible premium gel
- Formulated with non-toxic and non-caustic ingredients



MH76948	** 4" x 6.5"	48/case	\$49.49	\$47.02	\$47.02	\$44.54
MH73912	* 5" x 10.5"	12/case	\$18.79	\$17.85	\$17.85	\$16.91

MEDICO INTERNATIONAL COOL-JAW® REUSABLE COLD THERAPY GEL PACKS

Round

- Remain pliable even when frozen
- Feature a soft-sided surface to protect the skin



T430CB	** cool blue, 4"	50/case	\$76.49	\$72.67	\$72.67	\$68.84
T430PI	** pink ice, 4"	50/case	\$76.49	\$72.67	\$72.67	\$68.84
T430V	** vanilla, 4"	50/case	\$76.49	\$72.67	\$72.67	\$68.84
T430M	** mango, 4"	50/case	\$76.49	\$72.67	\$72.67	\$68.84
T430K	** kiwi, 4"	50/case	\$76.49	\$72.67	\$72.67	\$68.84
T440CB	* cool blue, 6"	25/case	\$56.99	\$54.14	\$54.14	\$51.29
T440PI	* pink ice, 6"	25/case	\$56.99	\$54.14	\$54.14	\$51.29

Rectangular

- Freeze solid and remain cold for extended periods of time



T410	** 3" x 4"	100/case	\$94.29	\$89.58	\$89.58	\$84.86
T400	* 3" x 5"	50/case	\$50.79	\$48.25	\$48.25	\$45.71

3M™ COVERS ONLY FOR COLD/HOT PACKS



Disposable

1572	4-3/4" x 10-1/2"	100/case	\$84.29	\$80.08	\$80.08	\$75.86
------	------------------	----------	---------	---------	---------	---------

Cool-Jaw Pouches and Sleeves are machine washable. Gel Packs are reusable. Clear gel packs freeze solid. Blue gel packs remain pliable when frozen and can be heated. *A freight charge of \$4.50 will be added to the price of this product. **A freight charge of \$5.50 will be added to the price of this product.

WOUND MANAGEMENT

MEDICO INTERNATIONAL COOL-JAW® FACIAL WRAPS

Cool-Jaw Econo Wrap

- The Econo Wrap is made from a poly blend fabric. It features one continuous sleeve that wraps around the head.
- The T-807 features two rows of Velcro® at the top of the head for an adjustable, universal fit.



T807

Wrap Only

T807S	white	25/case	\$105.29	\$100.03	\$100.03	\$94.76
-------	-------	---------	----------	----------	----------	---------

With 2 Clear Cold Gel Packs

T807	* white	25/case	\$125.79	\$119.50	\$119.50	\$113.21
------	---------	---------	----------	----------	----------	----------

MEDICO INTERNATIONAL COOL-JAW® SOFT STRETCH WRAP

Our T-800 and T-800C Soft Stretch Wraps are made with an elegant, satin spandex fabric that is extremely soft on the patient's skin. Both wraps feature a continuous outer sleeve that extends cheekbone to cheekbone to ensure full coverage over the surgical area. The T-800 and T-800C head wraps both feature two rows of Velcro® closure at the top of the head to provide a universal, adjustable fit. Also available with a built-in chin cup to maintain proper positioning.



T800



T80048

With 2 Clear Cold Gel Packs

T800-2	* white	24/case	\$190.49	\$180.97	\$180.97	\$171.44
--------	---------	---------	----------	----------	----------	----------

With 4 Clear Cold Gel Packs

T800-4	** white	24/case	\$231.99	\$220.39	\$220.39	\$208.79
--------	----------	---------	----------	----------	----------	----------

With 4 Blue Hot/Cold Gel Packs

T8004B	** white	24/case	\$275.29	\$261.53	\$261.53	\$247.76
--------	----------	---------	----------	----------	----------	----------

With Chin Cup and 2 Clear Cold Gel Packs

T800C-2	* white	24/case	\$190.49	\$180.97	\$180.97	\$171.44
---------	---------	---------	----------	----------	----------	----------

With Chin Cup and 4 Clear Cold Gel Packs

T800C-4	** white	24/case	\$213.79	\$203.10	\$203.10	\$192.41
---------	----------	---------	----------	----------	----------	----------

MEDICO INTERNATIONAL COOL-JAW® SOFT STRETCH WRAP-SPLIT STRAP WITH CHIN CUP

The T-900C Soft Stretch Wrap is made with an elegant, satin spandex fabric that is extremely soft on the patient's skin. The wrap features a continuous outer sleeve that extends cheekbone to cheekbone to ensure full coverage over the surgical area.



T900CBLU2B

The T-900C Soft Stretch Wrap features two rows of Velcro® closure at the top of the head to provide a universal, adjustable fit.

Split Wrap w/Chin Cup (Wrap Only)

T900CBLA	black	24/case	\$173.29	\$164.63	\$164.63	\$155.96
----------	-------	---------	----------	----------	----------	----------

T900CBLU	blue	24/case	\$173.29	\$164.63	\$164.63	\$155.96
----------	------	---------	----------	----------	----------	----------

With 2 Clear Cold Gel Packs

T900CBLA2	* black	24/case	\$206.79	\$196.45	\$196.45	\$186.11
-----------	---------	---------	----------	----------	----------	----------

T900CBLU2	* blue	24/case	\$206.79	\$196.45	\$196.45	\$186.11
-----------	--------	---------	----------	----------	----------	----------

With 2 Blue Hot/Cold Gel Packs

T900CBLA2B	* black	24/case	\$246.79	\$234.45	\$234.45	\$222.11
------------	---------	---------	----------	----------	----------	----------

T900CBLU2B	* blue	24/case	\$246.79	\$234.45	\$234.45	\$222.11
------------	--------	---------	----------	----------	----------	----------

With 4 Blue Hot/Cold Gel Packs

T900CBLA4B	** black	24/case	\$278.79	\$264.85	\$264.85	\$250.91
------------	----------	---------	----------	----------	----------	----------

T900CBLU4B	** blue	24/case	\$278.79	\$264.85	\$264.85	\$250.91
------------	---------	---------	----------	----------	----------	----------

Cool-Jaw Pouches and Sleeves are machine washable. Gel Packs are reusable. Clear gel packs freeze solid. Blue gel packs remain pliable when frozen and can be heated.

*A freight charge of \$4.50 will be added to the price of this product. **A freight charge of \$5.50 will be added to the price of this product.

HALYARD BILATERAL FACIAL ICE PACK

Materials: 3-layer STAY-DRY construction with soft outer covering, absorbent middle layer, film inner layer and blue polyester fabric for 5" x 12" strap attachment.



33101	Facial Ice Pack	12/box	\$94.29	\$89.58	\$89.58	\$84.86
-------	-----------------	--------	---------	---------	---------	---------

JOBST FACIAL SLING – SPLINT

- Non-sterile
- Easy to apply with Velcro® fastener
- Pressure dressing to limit edema
- Barton bandage
- Useful for initial treatment of condylar fractures
- Contains Latex



111825	Small	each	\$40.49	\$38.47	\$38.47	\$36.44
111826	Medium	each	\$40.49	\$38.47	\$38.47	\$36.44
111827	Large	each	\$40.49	\$38.47	\$38.47	\$36.44

MEDICO INTERNATIONAL FACIAL SUPPORT

- Made of lightweight spandex
- Extends 3/4" beneath the Adams Apple to apply pressure on the face and lower neck area
- Straps fasten over the head and around the neck



With Ear Openings

T118S	small	each	\$18.79	\$17.85	\$17.85	\$16.91
T118M	medium	each	\$18.79	\$17.85	\$17.85	\$16.91
T118L	large	each	\$18.79	\$17.85	\$17.85	\$16.91

Without Ear Openings

T124S	small	each	\$18.79	\$17.85	\$17.85	\$16.91
T124M	medium	each	\$18.79	\$17.85	\$17.85	\$16.91
T124L	large	each	\$18.79	\$17.85	\$17.85	\$16.91

Allografts

alloss[®]

alloss-Plus

alloguide[®] gf

Osteocel[®]
Plus
The CELLULAR
Advantage

Xenografts

NuOss[®]
PLUG

NuOss[®] DS **NuOss**[®] XC

RCM⁶[®]


CONFORM

RCP **RCF** **RCI**

MemGuide

Alloplasts

CaSO₄ **synoss**



**LAY THE PROPER
GROUNDWORK FOR
A HEALTHY SMILE.**

Allografts – Xenografts – Alloplasts

Solid solutions for maintaining and promoting healthy bone. Backed by years of documented and proven success, our products have been engineered to meet your specific surgical needs.

see pages 70-78 in this catalog

MEDLINK RADIOMAT™ X-RAY FILM

SG Extraoral Film

Full Speed Orthochromatic provides high contrast and excellent detail throughout the entire diagnostic range and rapid processing. This medium speed, high contrast film is compatible with any green-emitting rare-earth intensifying screen. When used with 400-speed green rare-earth screens, images have high clarity, high contrast, and overall image sharpness. When used with 100-speed green rare-earth screens, the film provides outstanding detail in extremity images.

EK771	15cm x 30cm	100 sheets/box	\$64.99	\$61.74	\$61.74	\$58.49
EK8EJ	5" x 12"	100 sheets/box	\$69.99	\$66.49	\$66.49	\$62.99

X-Ray Film Cross Reference

Medlink #	Carestream Dental #
EK771	1869080, 8372716
EK8EJ	1987627

CARESTREAM EXTRAORAL FILM

T-MAT G/RA Film

This dental film provides high-contrast, detailed images of intervening tooth structures while retaining good soft tissue visibility when used with LANEX regular or medium screens.

- Green light sensitive
- 400-speed when used with LANEX regular or EV screens
- The standard for extraoral film in the industry



Panoramic

1987627	TMG-5, 5" x 12"	50 sheets/box	\$145.99	\$138.69	\$138.69	\$131.39
1869080	TMG-15, 5" x 12"	50 sheets/box	\$173.79	\$165.10	\$165.10	\$156.41

EKTAVISION G

This film provides Enhanced Visualization due to its zero crossover control. When paired with either EV screens or LANEX regular screens, it provides a 400 speed system.

- Green light sensitive
- High contrast

Panoramic

1477405	EVG-5, 5" x 12"	50 sheets/box	\$183.99	\$174.79	\$174.79	\$165.59
8372716	EVG-15, 5" x 12"	50 sheets/box	\$216.99	\$206.14	\$206.14	\$195.29

X-OMAT DBF Film

- Blue light sensitive

Panoramic

Blue Light Sensitive – Not Interleaved*

1990712	DF-76, 6" x 12"	50 sheets/box	\$216.99	\$206.14	\$206.14	\$195.29
---------	-----------------	---------------	----------	----------	----------	----------

CARESTREAM INTRAORAL FILM

Ultra-Speed Film

A D-speed film that generates the high-quality, detailed images you need to make accurate diagnoses.

Periapical Film in Super Poly-Soft Pack



DF-57 Size 2 (1-1/4" x 1")

1753664	130 2-film packets/box	\$129.79	\$123.30	\$123.30	\$116.81
---------	------------------------	----------	----------	----------	----------

DF-58 Size 2 (1-1/4" x 1")

1658194	150 1-film packets/box	\$107.49	\$102.12	\$102.12	\$96.74
---------	------------------------	----------	----------	----------	---------

DF-54 Size 0 (7/8" x 1")

1228840	100 1-film packets/box	\$71.29	\$67.73	\$67.73	\$64.16
---------	------------------------	---------	---------	---------	---------

Occlusal Film in Paper Packets:

DF-50 Size 4 (2-1/4" x 3")

1666163	25 1-film packets/box	\$128.49	\$122.07	\$122.07	\$115.64
---------	-----------------------	----------	----------	----------	----------

CARESTREAM INSIGHT FILM

The standard for intraoral film. With INSIGHT dental film, you'll get a true F-speed film that requires up to 60% less radiation per exposure.

Periapical Film in Super Poly-Soft Packets:

IP-22 Size 2 (1-1/4" x 1-5/8")

1798628	130 2-film packets/box	\$159.29	\$151.33	\$151.33	\$143.36
---------	------------------------	----------	----------	----------	----------

IP-21 Size 2 (1-1/4" x 1-5/8")

1163401	150 1-film packets/box	\$119.99	\$113.99	\$113.99	\$107.99
---------	------------------------	----------	----------	----------	----------

IP-02 Size 0 (7/8" x 1-3/8")

1280619	100 2-film packets/box	\$123.29	\$117.13	\$117.13	\$110.96
---------	------------------------	----------	----------	----------	----------

IP-01 Size 0 (7/8" x 1-3/8")

8675332	100 2-film packets/box	\$97.99	\$93.09	\$93.09	\$88.19
---------	------------------------	---------	---------	---------	---------

CARESTREAM DUPLICATING FILM

X-OMAT 2 Dental Duplicating Film

X-OMAT 2 Dental Duplicating film is available in multiple sizes for both extraoral and intraoral x-ray duplication. The film can be exposed using UV light from available duplicators or white light.

1215821	5" x 12"	50 sheets/box	\$173.49	\$164.82	\$164.82	\$156.14
1493691	15 x 30cm	50 sheets/box	\$247.29	\$234.93	\$234.93	\$222.56
1586460	1-1/4" x 1"	150 packets/box	\$117.49	\$111.62	\$111.62	\$105.74

X-RAY

CARESTREAM GBX CHEMICALS

Manual processing of Carestream Dental intraoral and extraoral dental films

Carestream Dental GBX chemicals are liquid chemical concentrates designed for manual processing of all Carestream Dental dental films (both intraoral and extraoral). Both developer and fixer are single-solution concentrates, which provide mixing convenience, improved chemical stability, and identical mixing and dilution procedures for tank replenishment solutions.

- Each 37oz bottle makes 5 liters of solution



Developer and Replenisher

5158597 † 37oz	6/case	\$204.49	\$194.27	\$194.27	\$184.04
----------------	--------	-----------------	-----------------	-----------------	-----------------

Fixer and Replenisher

5158605 † 37oz	6/case	\$170.79	\$162.25	\$162.25	\$153.71
----------------	--------	-----------------	-----------------	-----------------	-----------------

CARESTREAM RAPID ACCESS TWIN PACK DEVELOPER AND FIXER

With Rapid Access dental chemicals, ready-to-read radiographs are available for interpretation in 30 seconds, without leaving the patient's side.

- 30-second processing (15-second develop, 15-second fix)
- Room-temperature solutions (70-76°F, 21-24°C)
- Ready-to-use, no mixing convenience
- Compatible with both Ultra-speed and INSIGHT dental films
- Long shelf life when unopened

Contains:

- one 500ml bottle of Rapid Access Developer
- one 500ml bottle of Rapid Access Fixer

Rapid Access Twin Pack

8606899 17oz, 500ml	1 each/pack	\$23.99	\$22.79	\$22.79	\$21.59
---------------------	-------------	----------------	----------------	----------------	----------------



CARESTREAM READYMATIC DENTAL CHEMICALS

READYMATIC dental chemicals have been specially developed for use with INSIGHT and Ultra-speed dental films. The result: consistently high resolution images.



When our dental chemicals are used to process our dental film, you'll get:

- Outstanding image quality for accurate diagnosis
- Superior clarity and resolution
- Compatibility with all leading dental film processors
- Simplified workflow - ready-to-use formulas mean no starter solution, mixing or dilution

Superb Results with Minimal Hassle

READYMATIC dental chemicals are designed for use in all automatic roller-transport dental x-ray film processors. They deliver clear, consistent results for all intraoral and extraoral radiographs. Plus, they:

- Provide excellent film transportability and sensitometry
- Are specifically formulated to limit residue on rollers and in tanks
- Will keep for up to two years, unopened and stored between 5-30°C (40-85°F)
- Are available in convenient chem packs or as separate solutions

Developer and Replenisher

5285929 * 5 L	2 x 5 liters/case	\$66.49	\$63.17	\$63.17	\$59.84
---------------	-------------------	----------------	----------------	----------------	----------------

Fixer and Replenisher

5285937 * 5 L	2 x 5 liters/case	\$64.29	\$61.08	\$61.08	\$57.86
---------------	-------------------	----------------	----------------	----------------	----------------

Dental Chem Pack

Contains: one of each, Developer Replenisher and one Fixer and Replenisher.

8606873 * 5 L	2 x 5 liters/case	\$63.99	\$60.79	\$60.79	\$57.59
---------------	-------------------	----------------	----------------	----------------	----------------

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

See page 2 for Details


† A freight charge of \$5.50 will be added to the price of this product.

* A freight charge of \$8.00 will be added to the price of this product.

PROTECTIVE MOUNT (COIN) ENVELOPES

- Manilla

S6285	2-1/2" x 4-1/4"	500/box	\$27.49	\$26.12	\$26.12	\$24.74
-------	-----------------	---------	---------	---------	---------	---------



ACE ACETATE TRACING PAPER

- Matte finish
- 50 sheets per pad
- 8" x 10" sheets

0032986	Acetate Paper, .003"	50/pad	\$30.49	\$28.97	\$28.97	\$27.44
0032987	Acetate Paper, .005"	50/pad	\$39.49	\$37.52	\$37.52	\$35.54



PALMERO PAN-O-MARKER

- This permanent x-ray film marker writes on wet or dry film
- Will not wash off or rub off; three replacement writing tips included

83	Pan-O-Marker	each	\$16.29	\$15.48	\$15.48	\$14.66
----	--------------	------	---------	---------	---------	---------



FLOW DENTAL PERFECT POCKET TEAR-AWAY MOUNT

Mounts that allow you to custom fit your x-ray records to your needs

The first ever masked pocket tear-away mounts you can buy! With Flow's Tear Off Perfect Pocket Mounts you can use as many windows as you need, no more... no less.

The Perfect Pocket Tear Away Mounts are in a format that allows you to tear off two mounts at a time. And each one over one window section is a complete mount, with space for patient info and dates to be written in. The windows accommodate a 2 H film.

Perfect Pocket Tear Away Mounts come in a box of 80 mount sheets, 12 windows per sheet. That gives you a total of 960 windows per box. Each sheet is 11.5" long, 4.125" high, making storage of whole sheets easy in a standard file draw set to accommodate letter size pages. The Tear Offs have rounded corners to protect you from sharp edges and their rigid construction allows for easy x-ray insertion.

810500	Xray Mounts	960/box	\$58.79	\$55.85	\$55.85	\$52.91
--------	-------------	---------	---------	---------	---------	---------



X-RAY



WOLF X-RAY SENSIBLES™ UNIVERSAL SENSOR HOLDER

Use just one bite block for all of your vertical and horizontal X-rays and just one bite block per patient. Bite block locks down to form a custom fit around any size sensor, in either the vertical or horizontal position. After exposure, simply press a button to release the sensor.

The Sensibles™ System includes anterior and posterior aiming rings and positioning arms so there's no worry about cone-cutting. Your images are always in parallel alignment.

Sensibles™ comes in two sizes. Large is suited for most size 1 and 2 sensors in the vertical position. Medium is suited for most size 1 and 2 sensors in the horizontal position and size 0 sensors in either the vertical or horizontal positions.

- Color-coded sizes
- Economically priced
- Quick and easy to use
- Can be used with Nexus Sensors
- Versatile and autoclavable

Econo Starter Kit Includes:

- 2 Large & 2 Medium Bite Blocks
- 1 Anterior Positioning Arm
- 1 Anterior Aiming Ring
- 1 Posterior Positioning Arm
- 1 Posterior Aiming Ring
- 25 Deluxe Comfees Samples

A 40901	Sensibles™ Econo Starter Kit	each	\$40.49	\$38.47	\$38.47	\$36.44
----------------	------------------------------	------	---------	---------	---------	---------

Deluxe Starter Kit Includes:

- 4 Large & 8 Medium Bite Blocks
- 2 Anterior Positioning Arms
- 2 Anterior Aiming Rings
- 2 Posterior Positioning Arms
- 2 Posterior Aiming Rings
- 50 Deluxe Comfees Samples

B 40900	Sensibles™ Deluxe Starter Kit	each	\$100.29	\$95.28	\$95.28	\$90.26
----------------	-------------------------------	------	----------	---------	---------	---------



Sensibles™ Refill Kits and Accessories

40910	Large (two-tone) Bite Blocks	12/pkg	\$48.79	\$46.35	\$46.35	\$43.91
40909	Medium Bite Blocks	12/pkg	\$49.99	\$47.49	\$47.49	\$44.99
408570	Anterior Positioning Arm, Stainless Steel	each	\$28.79	\$27.35	\$27.35	\$25.91
408580	Posterior Positioning Arm, Stainless Steel	each	\$28.49	\$27.07	\$27.07	\$25.64
408600	Posterior Aiming Ring	each	\$7.79	\$7.40	\$7.40	\$7.01

WOLF X-RAY RAPiD RIGHT ANGLE POSITIONING DEVICE

For PSP and X-Ray Film

RAPiD bite blocks feature rounded corners and are flexible for easy loading and unloading. Tabs ensure your film (or PSP plates) won't fall out during X-rays. Bite blocks and aiming rings are color-coded so assembling the correct positioner is always quick and easy.

The RAPiD is a right angle positioning device for intraoral dental X-ray film that can be used for anterior, posterior and bite-wing dental X-ray exams. The RAPiD can be used with any brand and size (sizes 3, 2, 1 and 0) of dental film. The right angle positioning system assures correct film/cone alignment and prevents "cone cutting".

RAPiDs are autoclavable for over 100 uses. Note that you cannot chemiclave or dry heat sterilize or use the autoclave drying cycle.

RAPiD Paralleling Kit with Bite Wings Includes:

- 4 Aiming Rings: 2 Anterior/Bite-Wing, 2 Posterior
- 4 Indicator Arms: 1 Anterior, 2 Posterior, 1 Bite-Wing
- 16 Bite Blocks: 4 Anterior, 8 Posterior, 4 Bite-Wings


C 408560	RAPiD Paralleling Kit	each	\$105.79	\$100.50	\$100.50	\$95.21
-----------------	-----------------------	------	----------	----------	----------	---------



PINNACLE DISPOS-A-BITE®

The Dispos-a-Bite is a disposable bite block designed to replace the present permanent bite block found on your panoramic x-ray machine. Each Dispos-a-Bite is unit dose packaged, which enables the bite block to be opened in front of each patient. This assures the patient that the bite block has not been used previously and is intended for single patient use only.

7060	Fits Planmeca	100/box	\$61.49	\$58.42	\$58.42	\$55.34
PIN7040	Fits GE Panolipse and Gendex	100/box	\$61.49	\$58.42	\$58.42	\$55.34
7030	Fits Phillips, Den-Tal-Ez Panolite, Soredex	100/box	\$65.79	\$62.50	\$62.50	\$59.21
7010	Fits Siemens OP 10, Instrumentarium	100/box	\$65.79	\$62.50	\$62.50	\$59.21
7020	Fits Siemens OP 3, OP 5	100/box	\$65.79	\$62.50	\$62.50	\$59.21



CLIKRAY™ SENSOR COVERS

001-152	For ALL Size 2 Sensors, DEXIS and Carestream Dental Size 2	500/box	\$34.79	\$33.05	\$33.05	\$31.31
001-138	For ALL Size 0, 1 Sensors	500/box	\$34.79	\$33.05	\$33.05	\$31.31
001-145	For Carestream Dental Size 1 and Thin Size 2 Sensors such as Schick, Suni and Dent-X	500/box	\$34.79	\$33.05	\$33.05	\$31.31



PALMERO DIGITAL X-RAY SENSOR PROTECTORS

- Disposable plastic covers for digital x-ray sensors
- Prevents cross contamination, protecting patients and staff
- Sizes available to fit most systems

1970	1-3/8" x 8"	500/box	\$14.49	\$13.77	\$13.77	\$13.04
1971	1-5/8" x 8"	500/box	\$14.49	\$13.77	\$13.77	\$13.04
1915A	2-1/2" x 10"	500/box	\$14.49	\$13.77	\$13.77	\$13.04



PALMERO BITE BLOCK COVERS

- Clear, convenient plastic barriers for use with Siemens x-ray machine bite blocks
- Economical way to reduce cross-contamination, eliminating the need for disinfection or sterilization

1850	Regular 1" x 2"	1,000/box	\$25.79	\$24.50	\$24.50	\$23.21
1850S	Siemens/Large 1-1/2" x 3"	500/box	\$25.79	\$24.50	\$24.50	\$23.21



WOLF X-RAY COMFEE'S™ DELUXE™ BARRIER SLEEVES

Comfee's™ makes sensor imaging easier and quicker for the technician, safer and more comfortable for the patient. Comfee's Deluxe™ feature a super-soft textured vinyl material that makes loading the sensors effortless. The added strength of the material also protects against tears.

Feature:

- Rounded corners for better fit and added patient comfort
- Unique tear-away seams so you can easily peel back the Comfee instead of tugging on the sensor wire
- Easy to insert and remove

80210	Large - fits most size 2 sensors	500/box	\$85.29	\$81.03	\$81.03	\$76.76
-------	----------------------------------	---------	---------	---------	---------	---------



WOLF X-RAY COMFEE'S™ CABLE SAVER™ SENSOR CABLE PROTECTION

- Cable Saver™ protects your sensor cables from being stepped on, rolled over, twisted and tangled.
- Easy to install and easy to keep clean, Cable Savers will add years to the life of your sensor cables.
- Available in 8 ft lengths. Simply cut to size.

W080500	Cable Saver™	8 feet	\$36.99	\$35.14	\$35.14	\$33.29
---------	--------------	--------	---------	---------	---------	---------



DEFEND X-RAY HEAD SLEEVES

- These X-ray head sleeves are designed to fit over most X-ray heads.
- The barrier is made of clear plastic and is easy to apply.

BF6000	15" x 26" Head Sleeves	250/box	\$34.79	\$33.05	\$33.05	\$31.31
--------	------------------------	---------	---------	---------	---------	---------



X-RAY



ALSO AVAILABLE
IN LEAD-FREE



ALSO AVAILABLE
IN LEAD-FREE



ALSO AVAILABLE
IN LEAD-FREE



PALMERO CLING SHIELD™ VINYL X-RAY APRONS / PROTECTION

Colors Available: Teal, Taupe, Slate Blue, Royal, Mauve, Grey, Green, Forest Green, Brown Check, Blue, Black, White, Beige Leatherette, Beige, Steel Grey

Due to customization, some aprons may incur a 2-4 week manufacturer lead time

Adult Patient Apron

- Covers from over the shoulders to below the gonadal area
- Loops attached for easy hanging and storage

20	Adult Patient Apron, 22-1/4" x 24-1/4"	each	\$121.29	\$115.23	\$115.23	\$109.16
----	--	------	----------	----------	----------	----------

Neck Collar

- For the protection of the thyroid area
- Velcro closing
- One size fits all
- Use in conjunction with Adult Patient Apron

25	Neck Collar, 22-5/8" x 4-1/4"	each	\$57.29	\$54.43	\$54.43	\$51.56
----	-------------------------------	------	---------	---------	---------	---------

Adult Protectall Apron

- Neck collar permanently attached
- Covers from thyroid to below the gonadal area
- Velcro clasp on neck collar
- Loops attached for easy hanging and storage

24	Adult Protectall Apron, 22-1/4" x 25-1/2"	each	\$152.99	\$145.34	\$145.34	\$137.69
----	---	------	----------	----------	----------	----------

24U	LEAD-FREE Adult Protectall Apron, 22-1/4" x 25-1/2"	each	\$252.49	\$239.87	\$239.87	\$227.24
-----	--	------	----------	----------	----------	----------

Pano-Adult Dual Apron

- Covers in front from shoulders to below gonadal area
- Back is the same length as front
- Velcro clasp at shoulder

26	Pano-Adult Dual Apron, 22-1/2" x 26-1/2"	each	\$278.29	\$264.38	\$264.38	\$250.46
----	--	------	----------	----------	----------	----------

26U	LEAD-FREE Pano-Adult Dual Apron, 22-1/2" x 26-1/2"	each	\$452.29	\$429.68	\$429.68	\$407.06
-----	---	------	----------	----------	----------	----------

Petite/Child Patient Apron

- Covers from over shoulders to below the gonadal area
- Loops attached for easy hanging and storage

22	Petite/Child Patient Apron, 20" x 19-1/2"	each	\$103.29	\$98.13	\$98.13	\$92.96
----	---	------	----------	---------	---------	---------

Petite/Child Protectall Apron*

- Neck collar permanently attached
- Covers from thyroid to below the gonadal area
- Sized for petite adults or children
- Velcro clasp on neck collar
- Loops attached for easy hanging and storage

27	Petite/Child Protectall Apron*, 20" x 20"	each	\$119.29	\$113.33	\$113.33	\$107.36
----	---	------	----------	----------	----------	----------

27U	LEAD-FREE Petite/Child Protectall Apron*, 20" x 20"	each	\$210.29	\$199.78	\$199.78	\$189.26
-----	--	------	----------	----------	----------	----------

* Child's Protectall Apron not standard with Smiley Face design. Special designs (Smiley Face, Happy Tooth or Happy Bear) are available for an additional charge.

Pano-Petite/Child Dual Apron

- Covers in front from shoulders to below gonadal area
- Back is the same length as front
- Velcro clasp at shoulder

26C	Pano-Petite/Child Dual Apron, 19-7/8" x 19-1/2"	each	\$195.79	\$186.00	\$186.00	\$176.21
-----	---	------	----------	----------	----------	----------

X-RAY

WOLF X-RAY DENTAL PATIENT APRONS

- Designed for complete patient protection of the neck, throat, chest, shoulders and gonadal area. Single Velcro closure make these aprons fast and easy to use.
- Cleaning: Hand wash in mild detergent and cool water. Not machine washable.
- 0.3mm Lead
- Color: Beige

Adult Apron

7505023	Adult Apron, 24" x 27"	each	\$190.99	\$181.44	\$181.44	\$171.89
---------	------------------------	------	----------	----------	----------	----------

Adult Apron, with thyroid collar

7505123	Adult Apron with thyroid collar, 24" x 27"	each	\$170.49	\$161.97	\$161.97	\$153.44
---------	--	------	----------	----------	----------	----------

Panoramic Dual Apron

- Provides excellent front & back radiologic protection. Velcro shoulder closure for quick application.

7506123	Panoramic Dual Apron, 24" x 27"	each	\$189.79	\$180.30	\$180.30	\$170.81
---------	---------------------------------	------	----------	----------	----------	----------

Thyroid Collar

- Designed for thyroid protection with rear Velcro closure.

7505630	Thyroid Collar, 4-3/4" x 26"	each	\$47.29	\$44.93	\$44.93	\$42.56
---------	------------------------------	------	---------	---------	---------	---------



PALMERO APRON HANGERS

Wall Mount Hanger, Standard

- White-coated steel hanger for use with Palermo's adult and child Patient and Protectall aprons only
- Hardware included for hanging



28	15" x 2"	each	\$18.29	\$17.38	\$17.38	\$16.46
----	----------	------	---------	---------	---------	---------

Wall Mount Hanger, Chrome-Coated

- Sturdy chrome-coated steel hanger
- For use with Palermo's adult and child Patient and Protectall aprons only
- Wall mounting, hardware included



PAL32	15"	each	\$18.49	\$17.57	\$17.57	\$16.64
-------	-----	------	---------	---------	---------	---------

X-Ray Apron Hanger

- White-coated steel hanger
- Can be used with Palermo's adult and child Patient and Protectall style aprons as well as the Coat and Dual style aprons



33	20.875" x 10.25"	each	\$21.79	\$20.70	\$20.70	\$19.61
----	------------------	------	---------	---------	---------	---------

Deluxe Coat Panoramic Hanger

- Sturdy chrome-coated steel hanger
- Specially constructed for use with panoramic aprons such as Palermo's Dual and Coat aprons



30	21" x 9.25"	each	\$18.49	\$17.57	\$17.57	\$16.64
----	-------------	------	---------	---------	---------	---------

WOLF DENTAL APRON HANGER

- Hangs any protective apron on wall hook or in a conventional closet



16401	23-1/2" x 9"	each	\$34.49	\$32.77	\$32.77	\$31.04
-------	--------------	------	---------	---------	---------	---------

MISCELLANEOUS

**COVIDIEN ARGYLE SALEM SUMP™
DUAL LUMEN STOMACH TUBES**

- Double lumen PVC tube - one for suction drainage; one for sump vent
- X-ray opaque Sentinel Line® with Sentinel Eye®
- Integral funnel connector with removable 5-in-1 adapter
- Integral irrigation funnel for irrigating through sump vent lumen
- Funnel end of vent lumen can be used to cap 5-in-1 connector
- Sterile



264945	48", 14 Fr	each	\$2.79	\$2.65	\$2.65	\$2.51
264960	48", 16 Fr	each	\$2.79	\$2.65	\$2.65	\$2.51
264986	48", 18 Fr	each	\$2.79	\$2.65	\$2.65	\$2.51

**MEDLINE SKIN
MARKERS**

- Medical-grade, non-irritating gentian violet ink
- Constant ink flow for clear marking
- Sterile
- Includes ruler



Regular Tip

DYNJSM01	50/case	\$39.99	\$37.99	\$37.99	\$35.99
----------	---------	---------	---------	---------	---------

Regular Tip and 6 Labels

DYNJSM02	50/case	\$46.29	\$43.98	\$43.98	\$41.66
----------	---------	---------	---------	---------	---------

Fine Tip

DYNJSM03	50/case	\$46.79	\$44.45	\$44.45	\$42.11
----------	---------	---------	---------	---------	---------

Fine Tip and 6 Labels

DYNJSM04	50/case	\$70.29	\$66.78	\$66.78	\$63.26
----------	---------	---------	---------	---------	---------

**COVIDIEN DEVON®
SURGICAL MARKERS**

- Sterile, single-use markers
- Non-toxic, non-smearing, non-irritating gentian violet ink
- Ruler cap for measuring surgical site



31145926



31145819

Regular Tip with Ruler Cap

31145926	25/box	\$54.79	\$52.05	\$52.05	\$49.31
----------	--------	---------	---------	---------	---------

Dual Tip (Regular/X-Fine) with Ruler Cap

31145900	25/box	\$71.29	\$67.73	\$67.73	\$64.16
----------	--------	---------	---------	---------	---------

Regular Tip

31145819	25/box	\$30.99	\$29.44	\$29.44	\$27.89
----------	--------	---------	---------	---------	---------

Regular Tip with Flexi-Ruler™

31145785	25/box	\$35.29	\$33.53	\$33.53	\$31.76
----------	--------	---------	---------	---------	---------

**SEPTODONT
HAND CREAM®**

- Permeates the skin, replenishing the moisture lost from frequent washing
- Forms an effective barrier against the irritating properties of your gloves



3-1/3 oz Tube

0110300	each	\$11.29	\$10.73	\$10.73	\$10.16
---------	------	---------	---------	---------	---------

**MEDLINE
MEDCELL
ALKALINE
BATTERIES**



MPHBAAZ	AA	24/box	\$11.29	\$10.73	\$10.73	\$10.16
MPHB9VZ	9V	12/box	\$20.29	\$19.28	\$19.28	\$18.26
MPHBAAZ	AAA	24/box	\$13.49	\$12.82	\$12.82	\$12.14
MPHBCZ	C	12/box	\$18.29	\$17.38	\$17.38	\$16.46
MPHBDZ	D	12/box	\$23.99	\$22.79	\$22.79	\$21.59

PALMERO PUMP DISPENSERS

- Dispense the right amount of liquid by pressing lightly on the cap, the swing away cap springs back to prevent evaporation and spillage
- Unbreakable plastic
- 8 oz capacity



PAL47	Brown	each	\$13.49	\$12.82	\$12.82	\$12.82
46	Clear	each	\$13.49	\$12.82	\$12.82	\$12.82
PAL48	Clear, "alcohol" imprint	each	\$13.49	\$12.82	\$12.82	\$12.82

**GRAHAM-FIELD
UNLABELED
SUNDRY JARS**

- Unlabeled clear flint glass jar with overlapping stainless steel cover
- 7"x 4.25" diameter
- ⚠️ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov



3458	set of 6	\$39.79	\$37.80	\$37.80	\$35.81
------	----------	---------	---------	---------	---------

**GRAHAM-FIELD
LABELED
SUNDRY JARS**

- Clear flint glass jar with overlapping stainless steel cover
- 7" x 4.25" diameter
- Each jar permanently imprinted in blue: Applicators, Gauze, Tongue Depressors, Cotton and Bandages
- ⚠️ WARNING: Cancer and reproductive harm, www.p65warnings.ca.gov



3454	set of 5	\$46.99	\$44.64	\$44.64	\$42.29
------	----------	---------	---------	---------	---------

MISCELLANEOUS

BREWER SQUARE STEEL WASTE CANS*

- Steel square outer can and lid with removable plastic inner liner
- Step-on foot pedal opens lid
- Perfect for exam rooms, restrooms, and front office areas
- Meets OSHA standard on blood borne pathogens
- Hazardous waste label included
- Enamel finish



Small-20 Quart/5 Gal Capacity (11"D x 11.75"W x 17.25"H)

35269	white	each	\$113.29	\$107.63	\$107.63	\$101.96
35270	red	each	\$113.29	\$107.63	\$107.63	\$101.96
35271	beige	each	\$113.29	\$107.63	\$107.63	\$101.96

Large-32 Quart/8 Gal Capacity (11.5"D x 12.25"W x 18.25"H)

35266	white	each	\$128.79	\$122.35	\$122.35	\$115.91
35267	red	each	\$128.79	\$122.35	\$122.35	\$115.91
35268	beige	each	\$128.79	\$122.35	\$122.35	\$115.91

*Non-stock item, may incur a 2-4 week manufacturer lead time

Looking for waste can liners?
See page 98 in the Disposables section

BREWER PLASTIC UTILITY CART

- 3 ABS plastic shelves (21.75" x 16")
- 12" between shelves
- 4 non-locking, 3" swivel polyurethane ball-bearing casters
- Total capacity: 200 lbs
- Easy to assemble
- Dimensions with handle: 25.25"L x 17.5"W x 35"H
- Assembly required



63530	each	\$344.29	\$327.08	\$327.08	\$309.86
-------	------	----------	----------	----------	----------

BREWER STAINLESS STEEL ALL PURPOSE CART

- Non-corrosive stainless steel construction
- 3 stainless steel shelves (15.5"W x 24"L)
- Quietly rides on 3.5" rubber wheel casters
- Total weight capacity: 300 lbs
- Dimensions with handle: 28"L x 15.5"W x 33"H
- No assembly required



63500	each	\$429.29	\$407.83	\$407.83	\$386.36
-------	------	----------	----------	----------	----------

LAKESIDE HEALTHCARE UTILITY CARTS

- Stainless steel
- Non-stock, drop ship from manufacturer, please allow up to 2 weeks for delivery
- All-welded angle leg design provides durability
- Non-returnable



LS359 | 469 (Heavy Duty)



356



LS474 (Heavy Duty)



LS311



358

Part Number	Drawer Qty.	Carrying Capacity	Shelf size Depth x Width	Overall size Depth x Width x Height	Shelf Clearance	SS Top Thickness	SS Legs Thickness	4 swivel Casters	ACE SOUTHERN Tier Pricing			
									Regular Price	Member Price	Premier Price	Premier+ Price
LS359	-	200 lbs.	16-5/8" x 18-3/8"	16-3/4" x 18-3/4" x 30"	10-3/4"	22 gauge	20 gauge	3-1/2" dia.	\$414.79	\$394.05	\$394.05	\$373.31
469	-	400 lbs.	16-5/8" x 18-3/8"	16-3/4" x 18-3/4" x 30-7/8"	10-7/8"	18 gauge	16 gauge	4" dia.	\$563.99	\$535.79	\$535.79	\$507.59
356	1	200 lbs.	16-5/8" x 18-3/8"	16-3/4" x 18-3/4" x 30"	8-1/4"	22 gauge	20 gauge	3-1/2" dia.	\$794.49	\$754.77	\$754.77	\$715.04
LS474	2	400 lbs.	20" x 36"	20-1/8" x 36-3/8" x 35"	10"	20 gauge	16 gauge	4" dia.	\$1,781.99	\$1,692.89	\$1,692.89	\$1,603.79
LS311	-	*200 lbs.	15-1/2" x 24"	16-1/4" x 27-3/4" x 31-7/8"	11-3/4"	22 gauge	20 gauge	3-1/2" dia.	\$318.79	\$302.85	\$302.85	\$286.91
358	2	200 lbs.	16-5/8" x 18-3/8"	16-3/4" x 18-3/4" x 35"	17-1/8"	22 gauge	16 gauge	3-1/2" dia.	\$1,309.29	\$1,243.83	\$1,243.83	\$1,178.36

*Includes bumpers on legs and handles

Tier Pricing Details see page 2	Regular Price	Member Price	Premier Price	Premier+ Price
---------------------------------	---------------	--------------	---------------	----------------

MISCELLANEOUS



DYNAREX BASIC WHEELCHAIRS

- Steel frame powder-coated silver vein
- Black padded embossed vinyl upholstery
- Detachable arm with padded armrests (desk length)
- Dual axle provides easy transition of wheelchair standard height to hemi level
- Push to lock brakes
- Chart pocket on backrest provides additional convenience
- Equipped with tool-free adjustable swing-away footrests
- 24" wheels with solid rubber tires and composite hand rims

10222	18"	each	\$271.49	\$257.92	\$257.92	\$244.34
10226	20"	each	\$271.49	\$257.92	\$257.92	\$244.34



DYNAREX FOOT STOOL

- Heavy duty steel construction with silver vein powder coat finish provides in use stability and long term durability
- Slip resistant ribbed rubber platform and thick rubber feet provide grip and stability on floor surfaces

DY10550		each	\$33.99	\$32.29	\$32.29	\$30.59
---------	--	------	----------------	----------------	----------------	----------------



MABIS DMI FOOT STOOL

- Non-slip textured matting and reinforced rubber tips add stability; contains latex
- Dimensions: 9.5"H x 14"W x 10"D
- Weight capacity: 250 lbs


53919010000		each	\$53.49	\$50.82	\$50.82	\$48.14
-------------	--	------	----------------	----------------	----------------	----------------



BLICKMAN FOOT STOOL

- Black non-slip rubber treads and feet tips
- Dimensions: 8"H x 18"W x 12"D

7757MR		each	\$285.29	\$271.03	\$271.03	\$256.76
--------	--	------	-----------------	-----------------	-----------------	-----------------



MID CENTRAL MEDICAL STACKING/INTERLOCKING STEP STOOL

- Durable 16 gauge #304 stainless steel construction
- Slip proof rubber tread on top surface for safety
- Stools are designed to interlock and stack to create desired platform
- Tip resistant
- Dimensions: 12"D x 18"W x 5"H

MCM140		each	\$295.29	\$280.53	\$280.53	\$265.76
--------	--	------	-----------------	-----------------	-----------------	-----------------

MISCELLANEOUS

BREWER CENTURY SERIES EXAM STOOLS*

Brewer's Most Popular Seating Series

The Brewer Century Series of medical seating can be found in many offices and clinics, due to its flexibility and overall comfort.

- 5-leg base for overall support
- 4" poly-foam seat cushion for maximum comfort
- Pneumatic height adjustments
- Optional back rest
- 47 standard colors (See Stool Color Selections below)
- 250 lb weight capacity
- Height adjustment 17.75" - 24.5"
- 5-year warranty



Pneumatic Height Adjustment without Backrest

11001	each	\$363.79	\$345.60	\$345.60	\$327.41
-------	------	----------	----------	----------	----------

Pneumatic Height Adjustment with Backrest

11001B	each	\$486.79	\$461.95	\$461.95	\$438.10
--------	------	----------	----------	----------	----------

BREWER VALUE PLUS SERIES STOOLS*

The Brewer Value Plus Series of exam room seating offers affordability without compromising style and reliability.

- 5-leg base for overall support
- 4" poly-foam seat cushion for maximum comfort
- 15" wide seat cushion
- Optional adjustable contoured backrest available
- Available in pneumatic or spin lift options
- 47 standard colors (See Stool Color Selections below)
- 250 lb weight capacity
- Height adjustment 17" - 21.25"
- 5-year warranty



Pneumatic Height Adjustment without Backrest

22500	each	\$306.49	\$291.17	\$291.17	\$275.84
-------	------	----------	----------	----------	----------

Pneumatic Height Adjustment with Backrest

22500B	each	\$412.99	\$392.34	\$392.34	\$371.69
--------	------	----------	----------	----------	----------

Spinlift Height Adjustment without Backrest

22400	each	\$218.29	\$207.38	\$207.38	\$196.46
-------	------	----------	----------	----------	----------

Spinlift Height Adjustment with Backrest

22400B	each	\$374.49	\$355.77	\$355.77	\$337.04
--------	------	----------	----------	----------	----------

STOOL COLOR SELECTIONS

Black Plum CG 3767	Alabaster US 396	Strato Grey AMB 110	Black US 393	Dove US 413	Ash UR 001	Greystone US 425	Sand US 384	Mushroom LV 017	Gray Taupe TL 58	Clamshell PH 57	Putty US 403	Gunmetal US 392	Graphite SF 73	Chocolate HAM 122
Black HAM 114	Buckskin RU 840	Seafoam CG 3550	Teal LV 001	Sage HAM 105	Tradewind TL 60	Azure Blue LV 003	Emerald NEO 40	Forest Green HAM 124	Petal SF 88	Mauve AM 45	Spanish Tile HAM 106	Burgundy US 364	Slate/ Yacht Blue NEO 26	Wedgewood MS 68
Agate Blue LV 003	Royal US 354	Claret US 357	Blue Fog PR 46	Port NEO 62										

•Matte texture ••Shiny texture

Note: The stool item number will be the stool number as a prefix and the color number as the suffix. Example: ordering a Century Series 11001 stool in the color Putty would be #11001-US403.

*Due to customization options, some items may incur a 2-4 week manufacturer lead time.

INDEX

A	AlloOss Plus 72	Autoclave Pouch Racks 255	Blades, Surgical 185-187
Absorbable Hemostat 318	AlloOss Putty 71	Autoclave Pouches 258-259	Block, Allograft 71
Access Absorbable Membrane 77	Ambu Bag with Manometer 64	Autoclave Tape 257	Blood Collection Sets 92
Access Device, Vacutainer 219	Ambu Bags 64-65	Autoclave Tubing 259	Blood Collection Tubes, Vacutainer 219
ACE Blades 186	Ampule Breaker 22	Autoclave Wrap 260	Blood Control, IV Catheters with 217
ACE Mini Blades 187	Anesthesia Bags 44	Autoclaves 254	Blood Glucose Meters 245
ACE New Gauze 321	Anesthesia Breathing Circuits 44	Autografters 84	Blood Lancets 245
ACE Suture 296-298	Anesthesia Masks 46	Autografter Screen Assembly 84	Blood Pressure Devices and Units 234
Acetate Paper 335	Anesthetic Carpules/Cartridges 17-18	Automatic Defibrillators 68-69 236-237
Actcel Hemostatic Gauze 317	Anesthetic Cream 20	Automatic Valves 214	Blood Transfer Device 219
Adapter, Cannula 51	Anesthetic Gel 18-19	Avitene 31	Bloodstop™ Absorbent Hemostat 318
Adapter, Christmas Tree 51	Anesthetic Liquid 19-20		Bluetooth Stethoscopes 240
Adapter, CO ₂ Nasal Hood 52, 54-55	Anesthetics, Benzocaine Topical 18-19	B	Blunt Cannulas 210-211, 280-281
Adapter, High Volume Tubing 276	Anesthetics, Topical 18-20	Bacitracin Ointment 17, 329	Blunt Needles 210-211, 280-281
Adapter, Multi-Sample Luer 220	Antifog Masks 126-129	Backhaus Towel Clamps 168	Board, CPR 68
Adapter, Nasal Hood 52, 54-55	Antimicrobial Soaps 150	Bacterial/Viral Filters 45	Body Restraints, Velcro Straps 66, 222
Adapter, Nipple & Nut 51	Antiseptic Hand Gel 150-151	Bag Mask Resuscitators 64-65	Bolt Gauze 325
Adapter, Tubing 51, 276-277	Apical Elevators 188	Bags, Biohazard/Infectious Waste 154	Bondable Traction Hooks 181
Adapter, Vial 23	Applicator Brushes, Micro 93	Bags, Take Home 99	Bonding Eyewear 93
Adherent Wrap, Self 328	Applicator Tips, Micro 93	Band-Aids 327	Bonding Lens Shield 93
Adhesive Bandages 327	Applicators, Cotton-Tipped 97	Bandage Scissors 163	Bone Basin 85
Adhesive Discs 240	Applicators, Cotton-Tipped 97	Bandages, Adhesive 327	Bone Chisels 87, 170
Adhesive Dressing, IV 327	Apron Gowns 293	Banicide 251	Bone Collectors 84
Adhesive, Light Cure 93	Apron Hanger 339	Barrier, CPR 68	Bone Compactors 85
Adhesive, Liquid Dressing 327	Aprons, Patient 141	Barrier Film 147	Bone Crushers 85
Adhesive, Periodontal Tissue 316	Aprons, X-Ray 338-339	Barrier Milk 251	Bone Curettes 195
Adhesive Remover 327	Arch Bar Coil, Erich 180	Barrier Sleeves 146-148	Bone Cutting Burs 110
Adhesive Strips 315, 327	Arch Bar Wax 180	Basins, Emesis 199	Bone Expanders 91
Adhesive, Tape 328	Arch Bars 180	Bathroom Tissue 98	Bone Files 170
Adhesive, Topical Skin 315	Arm Restraints, Velcro Straps 66, 222	Batteries, Alkaline 340	Bone Grafting Condensers/Pluggers 86
Administration Sets 206-207	Armboards 221-222	Batteries, Replacement for Otoscope	Bone Grafting Holding Forceps 86
..... 209, 211, 212, 213	Armboard Pads 222	and Ophthalmoscope Handles 242-243	Bone Grafting Material 74-78
Adson Tissue Forceps 166	Armboard Strap 66, 222	Beacon, Osstell 104	Bone Grafting Material Spoons 86
Advance Plus Burs 113	Arthrocentesis Instrument 204	Beakers 253	Bone Grafting Syringe 85
Advance Plus Suture 299	Articulating Paper 120	Beaver Mini Blades 187	Bone Harvesting Instruments 84-85
AEDs 68-69	Articulating Paper Forceps 168	Berman Airways 42	Bone Mallet 170
Air Fresheners 153, 154	Articulator 204	Bibs, Oral Surgery 95	Bone Rongeurs 169
Air/Nitrogen Adapters 108	Aspirating Syringes 282	Bibs, Patient 94-95	Bone Scraper 85
Air/Water Syringe Tips 262-263	Aspirator Adapters, High-Volume 276	Bibs, Prophy 94-95	Bone Screws 81
Airways 32-35, 37-43	Aspirator Brushes 272	Biohazard Bags 154	Bone Spreaders 171
Airways, Berman 42	Aspirator Tips 268	Biohazard Infectious Waste Containers 287	Bone Syringe 85
Airways, Guedel 42	Aspirators 270-271	Biohazard Labels 24	Bone Tap 115
Airways, Nasopharyngeal 42-43	Atlas Monitor Replacement Parts 233	Biological Indicators 256-257	Bone Traps 84
Alcohol Dispenser 247, 340	Attest Biological Sterilizer System 256	Biological Indicators, Mail-in 256	Bone Wax 319
Alcohol Prep Pads 149	Atomization Devices, Mucosal 22	Biopsy Punches 182	Bottle Cap Assembly 265-266
Alconox Powdered Precision Cleaner 247	Atraumatic Extracting Forceps 156-159	Bipolar Jeweler Forceps 108	Bottle Caps 266
Alginate Impression Material 118-119	Autoclavable Bins 98, 154	Bite Block Covers 337	Bottles, Collection 265-266
Alginate Measuring Sets 122	Autoclavable Bowls 199	Bite Blocks 175	Bouffant Caps 130
Alginate Spatulas 120-121	Autoclavable Cleaning Brushes 247	Bite Registration Material 119-120	Bowl, Sponge 199
AlloGuide GF 76	Autoclavable Delrin Spools 180	Blade Flask 187	Bowls, Alginate Mixing 121
AlloOss 70-72	Autoclavable Bags 259	Blade Handles 184	Bowman Probes 172
AlloOss Block Graft 71	Autoclave Case 198	Blade Removal Systems 187	Bracket Holder 181
AlloOss Demineralized Bone 70-71	Autoclave Cleaner, Chamber Brite 255	Blades, Laryngoscope 26, 28-31	Bracket Holder Tweezers 181
AlloOss Particulate 70-71	Autoclave Door Gaskets 255	Blades, Periosteum 192-194	Bracket Tray Covers 96

- Breathing Bags 44
Breathing Circuits 44
Breathing Masks 51
Breathing Tubing & Adapters 45
Brushes, Aspirator 272
Brushes, Instrument Cleaning 198, 247
Brushes, Micro-Applicator 93
Brushes, Surgical Scrub 150
Bulb Irrigation Syringe 283
Bulbs, Laryngoscope 28, 30-31
Bulbs, Air/Water Syringe Tip 262
Bulbs, Inflation 237
Bur Block 116
Bur Brush 107
Bur Coolant Tube 105
Bur Guards 107
Bur Holder 116
Bur/Implant Fixture Holder 116, 198
Bur Kits 114-115, 117
Bur, Twist 111
Burs 110-117
Butterfly Infusion Sets 218-219
BV Filters 45
- C**
- C-Fold Towels 98
Cabinets, Narcotic 21
Calcium Sulfate 75
Calipers 172
Cancellous Allograft 70-71
Cancellous Xenograft 74
Cancellous/Cortical Allograft 70-71
Canisters, Suction 264-265
Cannula Adapter 51
Cannulas, Dual Nasal 49-50
Cannulas, Irrigation 263
Cannulas, Nasal 48-50
Capnographers 232
Caps, Bouffant 130
Caps, Surgeons 130-131
Carbon Dioxide Absorbent 47
Carbon Dioxide Detector 36
Carpject Syringe System Holder 23
Carpules 17-18
Carts, Crash 21
Carts, Medical 21
Carts, Utility 341
Cassettes, Instrument 198
Castroviejo Caliper 172
Castroviejo Needle Holders 164
Catania Cheek Retractor 177
Cath-Guide Airways 42
Catheter Dressing, IV 222
Catheter Forceps, Magill 168
Catheters with Blood Control, IV 217
Catheters, IV 214-217
Catheters, Suction 65, 276
Cautery Pens, Disposable 108
Cavicide 151-152
Caviwipes 151-152
Ceiling Suspension Units 203
Cement, Temporary 93
Cement Spatulas 200
Centrifuge 92
Chain and Grid Traction Refill 181
Chair Sleeves 146
Chamber Brite Autoclave Cleaner 255
Cheek & Lip Retractors 176-177
Chemical Indicator Strips 257
Chemical Integrators 257
Chest Pieces 240
Chisels 87, 91, 170
Chlorhexidine Oral Rinse 320
Chloroprene Exam Gloves 135
Chloroprene Surgical Gloves 137
Christmas Tree Adapter 51
Chronic Gut Suture 297, 299-300
..... 304-305, 310, 312
Cidex OPA Plus Solution 250
Circuits, Breathing 44
Circular Saws 111
Clamps, Backhaus Towel 168
Clean & Simple Ultrasonic
Cleaning Tablets 252
Cleaners, Evacuation System 278-279
Cleaners, Glutaraldehyde-Free
Instrument 250
Cleaners, Instrument 248-251
Cleaning Brushes, Instrument 198, 247
Cleaning Spray, Handpiece 105, 248-249
Cleaning Towelettes, Optical 130
Closed IV Catheter System 216
Coats, Lab 292-293
Coats and Gowns, Protective 139-141
..... 292-293
Codman Surgical Patties 323
Coecide™ 251
Cohesive Wrap 328
Coin Envelope, X-Ray Film 335
Cold Compresses 330
Collagen Foam 78
Collagen Membrane 77
Collagen Plug 78
Collagen Tape 78
Collagen, Nuoss 74
Collagen Wound Dressing Products 78
Collar, Thyroid 339
Collatene 318
Collection Bottles 265-266
Collectors, Bone 84
College Pliers 166
COMBITUBE Airway 41
Compactors/Condensers/Pluggers,
Bone 86
Compressed Air/Nitrogen Adapters 108
Cone Burs 111
Cone Masks 128
Cone Socket Mouth Mirrors
and Handles 179
ConFORM Collagen Membranes 77
Connective Tissue Graft 183
Connector, Gas Sampling Elbow 45
Connectors, Needle-Free 209-213
Contra Angle Burs 110-117
Contra Angle Handpieces 102
Contra Angle Tissue Punches 182-183
Controlled Substances
Online Ordering System 246
Coolant Tubes 105
Corrugated Breathing Tubes 45
Corrugated Tubing,
Scavenging System 52, 56-57, 59
Corn Suture Pliers 166
Cortical Allograft 70-71
Cortical/Cancellous Allograft 70-71
Cots, Finger 131
Cotton Balls 97
Cotton (College) Pliers 166
Cotton & Dressing Forceps 166
Cotton Rolls 325
Cotton-Tipped Applicators 97
Countersink Burs 80
Coupland Aspirators 270
Covers, Headrest 146-147
CPR Barrier 68
CPR Board 68
Cradle, Head 66
Crash Carts 21
Creams/Lotions/Ointments 16-17, 329
Crestotome Burs 111
Cricothyroidotomy Device 67
Crile Forceps 167
Crown Gripper 168
Cryogenic Freezer 72-73
CSOS 246
CSR Wrap 260
Cuff Inflator, Laryseal 40
Cuffs, Blood Pressure 225, 227
..... 229, 233-234, 236-237
Cups, Drinking 97
Cups, Instrument 199
Cups, Iodine 199
Cups, Medicine 97, 199
Cups, Waxed 97
Curettes, Bone 195
Curettes/Scalers 195
Curing Light, LED 93
Curved Tip Syringes 263, 282
Cutters, Wire 163
Cylinder Burs 110, 113
Cylinder Carts 63
Cylinder Wrench 63
- D**
- D5W 223
Dappen Dishes 200
DEA 222 Form Sample 2
Dean Scissors 160
Debaquey Atraumatic Tissue Forceps 167
Defibrillators, Automatic 68-69
Delivery Syringes 85
Demand Valves 63
Demineralized Allograft 71
Dental Carpules/Cartridges 17-18
Dental Implants 100
Dental Mirrors 178-179
Dental Needles 280
Dental Rolls 325
Dental Syringes, Aspirating 282
Denture Relining Material 118
Depressors, Tongue 97
Depth Gauges 172
Dermabond Advanced 315
Dermal Punches 182
Detergents, Enzymatic 249
Developer and Fixer 134
Diamond Saws and Burs 111
Digital Calipers 172
Digital Narcotic Safe 21
Digital Thermometers 241
Directional Flow Assemblies 59
Disinfectants, Surface 151-154
Disinfecting and Sterilization Solutions,
Instrument 248-251
Dispenser Alcohol Bottle 247, 340
Dispensers, Label 24-25
Dispenser, Pump 247, 340
Dispensing Gun and Tips 122
Dispensing Pins 23
Dispos-a-Traps 278
Donut Positioner 66, 222
Double Stick Discs 240
Drainage Products, Wound 329
Drains, Penrose 329
Drains, Wound 329
Drapes 143-145

- Hemostat, Absorbable 317
Hemostatic Agents 317
Hemostatic Gauze 317
Hibiclens Surgical Soap 148
High Temp Cauteries 108
High-Volume Aspirator Adapters 276
Holder, Vacutainer 219-220
Holding Chamber, Valved 22
Hoods, Nasal 52-57
Hose Assemblies, Oxygen/Nitrous 61
Hoses for Surgi Drill Systems 106
Hot Packs 330
Hot/Cold Packs 330
Huck Towels 144-145
HVAC Suction Tips 84
Hypodermic Needles (needle only) 281
- I**
- Ice Packs 330-331
ID Bands, Instrument 201
Identification Labels 24
i-gel Supraglottic Airway 41
I-Hex Implants 100
Immobilizer, Elbow 222
Impact Air 45 Degree Handpiece 107
Implant Drill 101-103
Implant Machine 101
Implant Packs 142
Implant Probes 173
Implant Prophy Kits 196
Implant Related Products 101-104
Implant Scalars 196
Implant Tubing 105
Implantmed Surgery Unit 101
Implants 100
Impression Accessories 120-123
Impression Material 118-120
Impression Material Spatulas 120-121
Impression Tray Cleaner 123
Impression Trays 122-123
Indicator Strips, Steam 257
Indicator Tape 257
Indicators, Biological 256-257
Infectious Waste Bags 154
Infuser Bag 223
Infusion Pump Sets 207, 213
Infusion Sets, Safety Winged 218-219
Infusion Sets, Winged 218-219
Inhaler Adapter, Meter Dose 22, 47
Injection Plugs 217-218
Injection Site, Intermittent 217-218
Injection Site, Needle-Free 209, 210-213
Inline Vacuum Kit 57
- Instant Hot and Cold Packs 130
Instrument Cassettes 198
Instrument Cleaners 248-251
Instrument Cleaners,
Glutaraldehyde-Free 248-251
Instrument Cleaning Brushes 198, 247
Instrument ID Bands and Guards 201
Instrument Lubricants 105, 107, 248-251
Instrument Marking Tape 201
Instrument Rings 201
Instrument Stands 202-203
Instrument Trays 198-199, 223
Instrument Wrap 260
Instruments 156-204
Integrators, Chemical 257
Integrators, Type 5 Chemical 257
Interchangeable Periostomes 192-194
Intermaxillary Rulers 174
Intraoral Film, X-Ray 333
Intraoral Mirrors 178
Intraosseous Access Device 67
Intravenous Fluids 223
Intubation Kit 37
Iodine Cup 199
Iodine Prep Pads 149
Iris Scissors 160
Irrigating Tips 262-263
Irrigation Clips 105
Irrigation Fluids 262
Irrigation Needles, Endodontic 280
Irrigation Pipe 105
Irrigation Solutions 262
Irrigation Syringe, Bulb 283
Irrigation Syringes, Curved Tip 263, 282
Irrigation Tubing 105
Isolation Gowns 139-140
ISQ, Osstell 101, 104
IV Administration Sets 206-207
..... 209, 211, 212, 213
IV Administration Sets,
Needle-Free 210-213
IV Administration Sets,
Prepierced 209-210
IV Armboard 221-222
IV Catheter Dressing 222
IV Catheters 214-217
IV Catheters with Blood Control 217
IV Catheters, Safety 216-217
IV Fluids 223
IV Flush Solution 218
IV Stabilization Dressing 222
IV Stands 203
IV Start Kits 220
- J**
- Jackets, Lab 292-293
Jackets, Scrub 291
Jar, Stainless Steel Dressing 199
Jar, Stainless Steel Forceps 199
Jars, Sundry 200, 340
- K**
- Kelly Forceps 167
Kelly Scissors 161
Kerlix Roll Gauze 323
Kerrison Rongeur 169
Key Filler, Sevoflurane 47
King Vision aBlade Video Laryngoscope 26
Kirkland Knives 194
Kirschner Burs 111
Kleenex 98
Knitted Vicryl Mesh 79
Knives, Periodontal 194
Knives, Plaster 120-121
- L**
- Lab Coats 292-293
Lab Jackets 292-293
Lactated Ringer's 223
Label Dispensers, Syringe 24, 25
Labels, Biohazard 24
Labels, Syringe 24, 25
Lacrimal Probes 172
Lancets, Blood 245
Laparotomy Sponges 322
Laryngeal Mask Airways 38-41
Laryngoscope Blades and Handles 26
..... 28-31
Laryngoscope Bulbs 28, 30-31
Laryngoscope, Video 26
Laryseal Cuff Inflator 40
Latex Exam Gloves 135-136
Latex Surgical Gloves 138
Latex-Free Exam Gloves 131-135
Latex-Free Surgical Gloves 137
Lead Aprons 338
LED Curing Light 93
Lidocaine Topicals 20
Ligament Elevators 193
Light Cure Adhesive 93
Light Handle Covers 145, 148
Light, LED Curing 93
Light Sleeves 148
Lights, Exam and Procedure 244
Lights, Pen 244
Limb Leads 227, 231
Liners, Biohazard 154
- Liners, Nasal Hood 55, 57
Liners, Waste Can 98
Lip and Cheek Retractors 176-177
Locking Pin 105
Log Book, Narcotic 22
Log Book, Sterilization 256
LOOK Suture 300-302
Lotion, Skin Moisturizing 340
Lotions/Ointments/Creams 16-17, 329
Low Waste Syringe 283
Low-Profile Nasal Hoods 52-53
Lubricant, Handpiece Spray 105
..... 107, 248, 249
Luxating Hybrids 194
Lysol 153
Lysol, Foam 153
- M**
- Magill Catheter Forceps 168
Magnetic Drapes 203
Magnification Calculator 175
Mail-in Biological Indicators 256
Mallet 170
Manual Resuscitator Bags 64-65
Manual Suction Systems 66, 266
Markers, Skin 183, 340
Mask Harness 45
Mask Liners 55, 57
Mask, Pocket Resuscitation 68
Masks, Anesthesia 46
Masks, Earloop 124-127
Masks, Laryngeal 38-41
Masks, Oxygen 51
Masks, Sensitive 127
Masks, Surgical 128-129
Maxillofacial Fracture Kits 180
Mayo Stand Covers 144
Mayo Stands 202
Measuring Instruments 172-175
Measuring Set, Impression Material 122
Medical Carts 21
Medication Labels 24-25
Medicine Cups 97, 199
Medium Concentration Mask 51
Membrane, Collagen 77
Membrane Forceps 86
Membrane Placement Instruments 86
Membrane Punch 83
Membranes, ConFORM 77
Membranes, RCM6 77
MemGuide 76
Mesh 79, 82-83
Meter Dose Inhaler Adapter 22, 47

INDEX

Metricide.....	250	Needles, Filter.....	280	Otoscope Specula.....	243	Pliers, College.....	166
Micro-Applicators.....	93	Needles, Hypodermic (needle only).....	281	Otosopes.....	242-243	Pliers, Corn Suture.....	166
Microsurgery Scissors.....	162-163	Needles, Irrigation.....	280	Outlet Station.....	61	Pluggers/Condensers/Compactors.....	86
Mineralized Allograft.....	70-71	Neoprene Surgical Gloves.....	137	Oval Burs.....	110	Plug, Nuoss Collagen.....	74
Mineralized Xenograft.....	74	Neosporin Ointment.....	17, 329	Overshield Eyewear.....	93, 130	Plugs, Collagen.....	74,78
Mini Blades.....	187	Nipple and Nut Adapter.....	51	Oxygen Aspirators.....	63	Plugs, Injection.....	217-218
Mirrors, Mouth.....	178-179	Nitrile Exam Gloves.....	131-134	Oxygen Assemblies.....	58-63	Plunger, Titanium.....	84
Mixing Bowls, Alginate.....	121	Nitrogen Adapters.....	108	Oxygen Carts.....	63	Pocket Mask.....	68
Mixing Pads.....	97, 120	Nitrogen Regulators.....	108	Oxygen Connectors.....	51	Polychloroprene Exam Gloves.....	134
Mixing Spatulas.....	120-121	Nitrous Hose Assembly.....	61	Oxygen Flowmeters.....	58-60, 62	Polyisoprene Surgical Gloves.....	137
Moisture Traps.....	225, 227, 232	Non-Glutaraldehyde Instrument Cleaners.....	248, 249	Oxygen Hose.....	61	Pop-Off Valve, Hank's.....	47
Molded Surgical Masks.....	128	Non-Rebreathing Masks.....	51	Oxygen Masks.....	51	Porcine Membrane.....	76
Molt Bone Curettes.....	195	Non-Woven Sponges.....	322	Oxygen Regulators.....	61-62	Porous Tape.....	328
Molt Mouth Gags.....	176	Normal Saline Bags/Bottles.....	223, 262	Oxygen Tubing.....	50	Portable Suction Units.....	66, 266-267
Monoscopes.....	240	Nuoss.....	74	P		Positioner, Patient.....	66, 222
Monitors, Capnography.....	224-232	Nuoss Collagen.....	74	Packing Strips.....	321	Positioning Devices, X-Ray.....	336
Monitors, Patient.....	224-228	Nuoss DS.....	74	Pants, Scrub.....	290-291	Posterior Aiming Ring.....	336
Monotex PTFE Suture.....	313	Nuoss Plug.....	74	Paper Cups.....	97	Posterior Positioning Arm.....	336
Mosquito Forceps, Halsted.....	167	Nuoss XC.....	74	Paper Towels.....	98	Pouches, Autoclave.....	258
Motors for Elcomed.....	101	Nurse Caps.....	130	Particulate Grafts.....	70-75	Povidone Iodine.....	149
Motors for Implantmed.....	101	O		Patient Aprons.....	141	Precordial Stethoscope.....	240
Mouth Gags, Replacement Tips.....	176	Octagon Implants.....	100	Patient Bibs.....	94-95	Pregnancy Test, Urine.....	245
Mouth Gags, Molt.....	176	Offset Osteotome Kit.....	89	Patient Gowns.....	141, 293	Prep Pads, Alcohol.....	149
Mouth Mirror Handles.....	179	Ointments/Creams/Lotions.....	16-17, 329	Patient Monitors.....	224-226, 228	Prep Pads, Povidone Iodine.....	149
Mouth Mirrors.....	178-179	Ointments, Topical Wound.....	16-17, 329	Patient Strap.....	66, 222	Prep and Scrub Solutions.....	148-149
Mouth Props.....	175	OK Autoclave Steam Indicator Strips.....	257	Patient Towel/Bib Holder.....	94-95	Prepierced IV Admin Sets.....	209, 210
Mucosal Atomization Devices.....	22	OPA Disinfecting Solutions.....	250	Pear Burs.....	111	Pressure Infuser Bag.....	223
Multi-Sample Luer Adapters.....	220	Ophthalmic Ointments and Solutions.....	16	Pen Lights.....	244	Probes/Explorers, Periodontal.....	173
N		Ophthalmoscope.....	242	Penrose Drains.....	329	Probes, Lacrimal.....	172
N95 Surgical Mask.....	128	Ophthalmoscope Handle Batteries.....	243	Pens, Cautery.....	108	Probes, Measuring.....	172-173
Napkin Holders.....	95	Ophthalmoscope Handles.....	242-243	Pens, Skin Marking.....	183, 340	Procedure Lights.....	244
Narcotic Cabinets.....	21	Ophthalmoscope Handles.....	242-243	PeriAcryl Tissue Adhesive.....	116	Procedure Masks.....	124-127
Narcotic Log Book.....	22	Optical Cleaning Towelettes.....	130	Peridex Oral Rinse.....	320	Prophy Angles.....	99
Narcotic Safes.....	21	OptiChamber Diamond VHC.....	22	Perio Brushes.....	99	Prophy Bibs.....	94-95
Nasal Cannulas.....	48-50	O.R. Towels.....	143-145	Perio Care.....	320	Prophy & Handpiece Lube.....	105, 107 248, 249
Nasal Cannulas, CO ₂ Sampling Only.....	49	Oral Medication Syringe.....	23, 282	Periodontal Dressing, Coe-Pak.....	320	Props, Mouth.....	175
Nasal Cannulas, Dual.....	49-50	Oral Rinse, Chlorhexidine.....	320	Periodontal Knives.....	194	Protective Coats and Gowns.....	139-141 292-293
Nasal Hood Adapters.....	61	Oral Suction Systems, Manual.....	66, 266	Periodontal Ligament Elevators.....	193	Protective Eyewear.....	93, 130
Nasal Hood Adapters, CO ₂	52, 54-55	Oral Surgery Burs.....	110-113	Periodontal Probes/Explorers.....	173	Protective Gowns and Coats.....	139-141 292-293
Nasal Hood Liners.....	55, 57	Oral Swabsticks.....	99	Periodontal Tissue Adhesive.....	316	Protective Sleeves.....	146-148
Nasal Hood Punch.....	55	Oraqix.....	19	Periosteal Elevators.....	191	Protractors.....	174
Nasal Hoods.....	52-57	Orasoothie® Sockit Gel.....	319	Periotomes.....	192-194	PRP Blood Collection Kit.....	92
Nasal Hoods, CO ₂ Monitoring.....	52, 54-55	Orban Knives.....	194	Periotomes, Interchangeable.....	192-194	PTFE Barrier Membrane.....	79
Nasal Hoods, Low-Profile.....	52-53	ORCA™ Foam.....	318	Pharmaceutical Waste Containers.....	25	PTFE Suture.....	296, 302, 313
Nasal Specula.....	176	ORCA™ Pain Relief Putty.....	319	Pharmaceuticals.....	5-18	Pull-Tight Seals.....	22
Nasopharyngeal Airways.....	42-43	Orthodontic Chains.....	181	Pharmaceuticals, Emergency.....	7-8	Pulse Oximeters.....	233, 235
Needle-Free Connectors.....	209-213	Orthotraction Pads.....	181	Pharmaceuticals, Schedule II.....	5	Punches, Biopsy.....	182
Needle-Free Injection Sites.....	209-213	Osseous Glider.....	85	Pharmaceuticals, Schedule III, IV, V.....	6-7	PUREVAC SC Evacuation Cleaner.....	278
Needle-Free IV Admin Sets.....	210-213	Osstell, Beacon, ISQ, SmartPegs.....	104	Pillowcases.....	96	Q	
Needle Holders.....	164-165	Osteocel Plus.....	73	Placental Membrane.....	76	Quicksmart Blade Removal System.....	187
Needle Recapping Device.....	286	Osteotomes.....	88-89, 171	Plain Gut Suture.....	297, 299-300 303, 305, 310, 312	QuickTrach.....	67
Needle Sheath Prop.....	286	Otoscope Ear Tips.....	243	Plastic Cups.....	97		
Needles, Blunt.....	210-211, 280-281	Otoscope Handle Batteries.....	243	Plastic Retractors.....	176		
Needles, Dental.....	280	Otoscope Handles.....	242-243				

S

Rack, Autoclave 255
 Rack, Blood Collection Tube 92
 Radiosurge Electrosurgery Unit 109
 Ramus Retractor 87
 Ratchet Extension 182
 Ratchets, Universal 91
 RCM6 Membrane 77
 READYMATIC Dental Chemicals 334
 Recapping Device 286
 Red Biohazard Bags 154
 Refrigerant, Skin 20
 Regulators 61-62
 Reline Material 118
 Resorbable Collagen 77
 Resuscitator Bag with Manometer 64-65
 Resuscitator Bags 64-65
 Retractors, Lip and Cheek 176-177
 Retractors, Ramus 87
 Ridge Expanding Kits 91
 Right Angle Driver Kit 117
 Rings, Instrument 201
 Rongeurs 169
 Root Forceps 157
 Root Tip Elevators 190
 Root Tip Picks 190
 Ruler, Intermaxillary 174
 Rulers 174
 Russian Tissue Forceps 167

S

Safes, Narcotic 21
 Safety Butterfly Infusion Sets 218-219
 Safety Catheters 216-217
 Safety Eyewear 93, 130
 Safety IV Administration Sets 209-213
 Safety IV Catheters 216-217
 Safety Lancets 245
 Safety Needles (needle only) 285
 Safety Scalpels 185
 Safety Syringe with Needle 286
 Safety Winged Infusion Sets 218-219
 Salem Sump Tubes 340
 Saline Bags/Bottles, Normal 223, 262
 Saliva Ejectors 268
 Sampling Lines, Gas 50
 Sani-Cloths 152
 Sani-Sleeve 147
 Scalars, Implant 196
 Scalars/Curettes 195-197
 Scalpel Blade Remover 187
 Scalpel Handles 184
 Scalpels, Safety 185
 Scavenger Vacuum Adapter 56
 Scavenging Cone Assembly 56-57

Schedule II Online Ordering System 246
 Schedule II Pharmaceuticals 5
 Schedule III, IV, V Pharmaceuticals 6-7
 Schrader Air Adapter 108
 Scissors 160-163
 Scissors, Dean 160
 Scissors, Wire 163
 Scraper, Bone 85
 Screen Assembly, Autograft 84
 Screws 81
 Scrub Brushes, Surgical 150
 Scrub and Prep Solutions 148-149
 Scrubs 288-291
 Scrubs Color Chart 294
 Scrubs Embroidery Information 295
 Scrubs Size Chart 294
 Scrubs Embroidery Information 295
 Seals, Pull-Tight 22
 Secure Fit Masks 125
 Self-Adherent Wrap 328
 Sensitive Masks 127
 Sensor Covers, X-Ray 337
 Sets, IV Administration 206-207
 209, 211, 212, 213
 Sevoflurane Key Filler 47
 Sharps Containers 287
 Sheets, Sterile Drape 143-144
 Shirts, Scrub 288-290
 Shoecovers 131
 Silhouette Flowmeter Connector 53
 Silicone Replacement Tips,
 Mouth Gag 176
 Silicone Spray 105
 Silk Suture 298-299
 301, 304, 308, 311-312
 Sinus Curette 90
 Sinus Lift Instruments 88-90
 Sinus Lift Kit 90
 Sinus Lift Osteotomes 88-89
 Sinus Lift Instruments 88-90
 Site Dilator Kit 88
 Size Chart, Scrubs 294
 Skin Adhesive, Topical 315
 Skin Closure Strips 315, 327
 Skin Markers 183, 340
 Skin Refrigerant 20
 Skin Staplers 315
 Sleeves, Protective 146-148
 SmartPegs, Osstell 104
 Soaking Tub 253
 Soaps, Hand 150
 Soaps/Scrubs, Surgical 148-149
 Socket Preservation Kit 78
 Sockit Oral Wound Dressing 319
 Soda Lime 47
 Sodium Chloride Bags/Bottles 223, 262

Sodium Chloride Flush 218
 Solidification System,
 Solidifier Canister Fluid 279
 Souffle Cups 97
 Spade Tip Elevators 188-189
 Spatula Osteotomes 171
 Spatulas, Impression Material 120-121
 Specula, Otolaryngology 243
 Sphygmomanometers 236-237
 Spinal Needles 282
 Split Sheets 144
 Sponge Bowl 199
 Sponges, Gauze 322
 Sponges, Laparotomy 322
 Sponges, Non-Woven 322
 Spool Wire 180
 Spore Check 255-256
 Spray Cap, Handpiece Lubricants 105
 Spray, Handpiece Cleaning 105, 248-249
 Spray Lubricant, Handpiece
 105, 107, 248, 249
 Spray Tubes 105
 Sprays, Disinfectant 151-154
 Sprays, Topical Anesthetic 19, 20
 Stacking Step Stool 342
 Stainless Steel Ware 198-199
 Stands, Mayo Instrument 202-203
 Staplers, Skin 315
 Steam Indicator Strips, Chemical 257
 Steam Sterilization Integrator 257
 Steiglitz Root Forceps 157
 Step Stools 342
 Steri-Strip Skin Closures 315, 327
 Sterile Water Bags/Bottles 223, 262
 Sterilization Indicator Tape 257
 Sterilization Pouches 258
 Sterilization Trays/Cassettes 105, 198-199
 Sterilization Tub 253
 Sterilization Tubing 259
 Sterilization Wrap 260
 Sterilizing and Disinfecting Solutions,
 Instrument 248-251
 Sterilizing Forceps 168
 Steri-Strips 315, 327
 Stethoscope, Precordial 240
 Stethoscopes 238-240
 Stockinets, Tubular 329
 Stools, Exam 343
 Stopcocks 214
 Strap, Armboard 66, 222
 Strep A Test 245
 Stretcher Sheets 144-145
 Stools, Exam 343
 Stopcocks 214
 Strap, Armboard 66, 222
 Strep A Test 245
 Stretcher Sheets 144-145
 Strips, Packing 321
 Stylets, Endotracheal Tube 32, 36
 Suction Assembly Units 265-266

Suction Canisters 264-265
 Suction Catheters 65, 276
 Suction Handles, Yankauer 269-270
 Suction Irrigator 262
 Suction Systems, Manual 66, 266
 Suction Tip Brushes 272
 Suction Tips 84, 271-272
 Suction Tube, Frazier Ferguson 272
 Suction Tubing 84, 274-275
 Suction Units, Handheld 66, 266
 Suction Units, Portable 66, 266-267
 Sundry Jars 200, 340
 Surface Barrier Film 147
 Surface Disinfectants/Wipes 151-154
 Surgeons Caps 130-131
 Surgeons Gloves 137-138
 Surgi Drill Systems 106
 Surgical Aspirators 268, 270-271
 Surgical Blades 186-187
 Surgical Caps 130-131
 Surgical Drapes 143-144
 Surgical Gloves 137-138
 Surgical Gowns 140-141
 Surgical Markers 340
 Surgical Masks 128-129
 Surgical Packs 142-143
 Surgical Patties, Codman 323
 Surgical Scrubs & Soaps 148-149
 Surgical Tape 328
 Surgicel 317
 Surgifoam 317
 Suspension Units, Ceiling 203
 Suture 296-313
 Suture Pliers 166
 Suture Spools 313
 Swabsticks, Oral 99
 Swabsticks, Povidone Iodine 149
 Sweetheart Retractors 177
 Synoss 75
 Synthetic Bone Grafting Material 75
 Syringe Label Dispensers 24, 25
 Syringe Labels 24, 25
 Syringe with Needle 284-286
 Syringe Only 283
 Syringe with Preattached
 Blunt Cannula 286
 Syringe Pump Sets 207, 213
 Syringe Sleeves 148
 Syringe Tip Cap 284
 Syringe Tips, Air/Water 262-263
 Syringes, Aspirating 282
 Syringes, Bone Graft Delivery 85
 Syringes, Bulb Irrigation 283
 Syringes, Curved Tip Irrigation 263, 282
 Syringe, Oral Medication 23, 282

INDEX

Syringes, Safety 286
 Syringes 283-286

T

Table Cover 96
 Table Paper 96
 Take Home Bags 99
 Tape, Autoclave 257
 Tape, Collagen 78
 Tape, Instrument Marking 201
 Tape, Surgical 328
 Tape, Waterproof 328
 Taper Burs 110, 113
 Taper Dome Burs 110, 113
 Tegaderm IV Dressing 222
 Telfa™ Ouchless Pads 323
 Temporary Cement 93
 Tension Slide Adjuster 57
 Tenting Screws 81
 Test, Strep A 245
 Test Tube Rack 92
 Test, Urine Pregnancy 245
 Thermometers 241
 Throat Packing 323
 Thumb Knob 182
 Tie-On Surgical Masks 128-129
 Tissue Adhesive, Periodontal 316
 Tissue, Bathroom 98
 Tissue, Facial 98
 Tissue Forceps 86, 166-167
 Tissue Harvester 183
 Tissue Punches 182-183
 Tissue Retractor 177
 Titanium Cutting Burs 111
 Titanium Foil 83
 Titanium Mesh 82-83
 TMJ Dysfunction Treatment 204
 Toilet Tissue 98
 Tongue Depressors 97
 Tongue Retractors 177
 Tongue Seizing Forceps 168
 Toothbrushes 99
 Toothbrushes with Xylitol, Single-Use 99
 Tooth Protector 36
 Topical Anesthetic Cream 20
 Topical Anesthetic Gel 18-19
 Topical Anesthetic Liquid 19-20
 Topical Anesthetic Patch 20
 Topical Anesthetic Sprays 19-20
 Topical Anesthetics 18-20
 Topical Ointments
 and Solutions 16-17, 329
 Topical Wound Ointments 16-17, 329
 Tourniquets 220

Towel Clamps, Backhaus 168
 Towels, Folded 98
 Towel Forceps 168
 Towels, O.R./Huck 143-145
 Towels, Patient 94-95
 Tracheostomy Device 67
 Tracing Paper 335
 Traction, Chain and Grid 181
 Transfer Device, Blood 219
 Transparent Tape 328
 Transpore Tape 328
 Traps, Bone 84
 Traps, Evacuation 277, 278
 Trash Bags 98
 Trash Can Liners 98
 Trash Cans 98, 341
 Tray Cleaner, Impression 123
 Tray Covers 96
 Tray Sleeves 147
 Trays, Impression 122-123
 Trays, Instrument 198-199, 223
 Trepine Burs 114, 115, 117
 Trepine Drill Kit 114, 117
 Tri-Cam Implants 100
 Triple Outlet Station 61
 truFIX Fixation Kit 80-81
 truScrews 81
 truTacks 81
 truTents 81
 Tube Clamps 105
 Tube Rack 92
 Tube Sleeve 147
 Tube Wire 180
 Tubing Adapters, High Volume 276
 Tubing, Irrigation 105
 Tubing, Oxygen 50
 Tubing, Sterilization 259
 Tubing, Suction 84, 274-275
 Tubing, Vacuum 56
 Tubular Stockinettes 329
 Tunics, Scrub 288-289
 Twist Burs 111
 Type 5 Chemical Integrators 257

U

U17 Utility Pick-Up Dressing Pliers 166
 Ultrasonic Cleaners 252-253
 Ultrasonic Cleaning Tablets,
 Clean & Simple 252
 Ultrasonic Cleaning Units 253
 Ultrasonic Solutions 248-251
 Umbilical Tape 313
 Underpads 145
 Uniforms 288-290

Unigrift 183
 Universal Drivers, Implant 117, 204
 University of Minnesota
 Cheek Retractors 177
 Urine Pregnancy Test 245
 Utility Carts 341
 Utility Drapes 143-144
 Utility Gloves 133, 247

V

V Vac Manual Suction System 66, 266
 Vacutainer Access Device 219
 Vacutainer Blood Collection Tubes 219
 Vacutainer Holder 219-220
 Vacutainer Needle 220
 Vacuum Hose 57
 Vacuum Tubing 56
 Valved Holding Chamber 22
 Valves, Demand 63
 Varijet Suction Irrigator 262
 Velcro Limb Leads 231
 Velcro Straps 66, 222
 Vein Finder 220
 Venoscope Transilluminator 220
 VHC, OptiChamber Diamond 22
 Vial Adapters 23
 Vicryl Mesh 79
 Video Laryngoscope 26
 Vital Signs Devices 225-226, 234

W

Waste Can Liners 98
 Waste Cans 98
 Waste Cans 341
 Water Cartridge 247
 Waste Container, Pharmaceutical 25
 Water, Sterile Bags/Bottles 223, 262
 Water Traps 225, 227, 232
 Waterline Treatment 247
 Waterproof Tape 328
 Wax, Bone 319
 Weider Cheek and Tongue Retractor 177
 Wenger Chest Pieces 240
 W&H Units and Accessories 101-105
 Wheelchairs 342
 Whitegaps Nasal Hood Liners 55
 Williams Probes 172
 Winged Infusion Sets 218-219
 Winter Arch Bars 180
 Wipes, Surface Disinfectant 151-152
 Wire Cutters 163
 Wire Cutting Scissors 163
 Wire Twisters 163
 Wire, Spool 180

Wire, Tube 180
 Woodson Filling Instrument 200
 Wrap, Cohesive 328
 Wrap, Self-Adherent 328
 Wrap, Sterilization 260

X

Xenograft Bone Particulate 74
 X-Ray Apron Hangers 339
 X-Ray Detectable Gauze 321
 X-Ray Film 333
 X-Ray Film Coin Envelopes 335
 X-Ray Film Marker 335
 X-Ray Film Mounts 335
 X-Ray Magnification Calculator 175
 X-Ray Processing Chemicals 334
 X-Ray Ruler 175
 X-Ray Sensor Covers 337
 X-Ray Sensor Protectors 337
 X-Ray Sleeves 148, 337
 X-Trac System 117
 Xylitol, Single-Use Toothbrush with 99

S

Y Connector 105
 Y Selector Valves 61
 Yankauer Suction Handles 269-270
 Young Tongue Seizing Forceps 168

Z

Zest Locator Abutments 100

TERMS

- **All invoices are net 30 days.**
- ACE SOUTHERN reserves the right to assess a late charge of \$25.00 on accounts past due.
- For inquiries on your account, please call our **Accounting Department** toll free at **1-800-455-8500** 8am - 7pm EST or contact them by email: **info@acesouthern.com**.
- Due to manufacturer policies, prices, discounts and/or product information are subject to change without notice.
- Items on long term backorder due to manufacturer supply are subject to cancellation.

PAYMENT

We accept checks, Visa, MasterCard, Discover, and American Express for payment.

STATEMENTS

Statements are mailed monthly. Please contact our Accounting Department Toll Free at **1-800-455-8500** 8am - 6pm EST or by email: **info@acesouthern.com**, if you have any questions regarding monthly statements.

CONTROLLED SUBSTANCES ORDERING

- **The Drug Enforcement Administration requires a copy of your current Federal DEA registration certificate to be on file prior to shipping controlled substances.**
- **Federal Law requires a properly completed DEA Form 222 for all Class II Pharmaceuticals.** Shipment of Class II Pharmaceuticals cannot be made until this form has been received. Improperly completed forms, or forms with alterations, will be returned. In lieu of using the DEA Form 222, orders may also be placed using our electronic 222 process.
- Regulations require us to limit the sale of Rx and controlled substances only to registered, licensed healthcare professionals. If you are a new customer or have recently moved, please furnish us with a copy of your updated state and federal registrations verifying your shipping address. Please note that all orders for controlled substances are subject to a due diligence review process. Schedule II controlled substances can be ordered electronically or by mail. For more information on our Controlled Substance Ordering System please visit <https://www.acesouthern.com/e222>; if you prefer to continue using Federal 222 Forms to order Schedule II controlled substances, mail the form to:

ACE SOUTHERN
One Southern Court
West Columbia, SC 29169

- ACE SOUTHERN restricts the sale or other transfer of medications to prisons/correctional facilities for use in lethal injections, based on our manufacturer agreements. The goods ACE SOUTHERN sells are intended to be used for their label-approved purposes or applicable standards of care, which do not include human lethal injection.

WARRANTIES

ACE SOUTHERN will pass through to the customer, at the time of sale, any transferable product warranties, indemnities and remedies provided to ACE SOUTHERN by the applicable manufacturer. EXCEPT AS OTHERWISE PROVIDED HEREIN, TO THE EXTENT PERMITTED BY LAW, ACE SOUTHERN PROVIDES NO WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT

LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE OR NON-INFRINGEMENT, AND THE CUSTOMER SHALL LOOK TO THE MANUFACTURER OF THE PRODUCT FOR ANY WARRANTY THEREON.

LIMITATION OF LIABILITY

THE CUSTOMER AGREES TO LOOK SOLELY TO THE MANUFACTURER OF THE PRODUCT FOR ANY CLAIM ARISING DUE TO LOSS, INJURY, DAMAGE OR DEATH RELATED TO THE USE OR SALE OF PRODUCTS. ACE SOUTHERN SHALL NOT BE LIABLE FOR INDIRECT, INCIDENTAL, PUNITIVE, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, LOST PROFITS AND LOSS OF GOODWILL, ARISING FROM OR RELATING TO ANY BREACH TO THIS AGREEMENT (OR ANY DUTY OF COMMON LAW, AND WHETHER OR NOT OCCASIONED BY THE NEGLIGENCE OF ACE SOUTHERN OR ITS AFFILIATES), REGARDLESS OF ANY NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.

SHIPPING

Order by 2:30pm EST and all in stock items will ship the same day (excluding weekends and federal holidays).

All orders are subject to shipping charges which may include:

- Minimum order fee
- Hazardous material and refrigerated product fees
- Excessive weight fee
- Express shipping charge
- Order handling fee

RETURN POLICY

If you are not satisfied with your purchase, contact our Customer Service department within 30 days of purchase at 800-624-5926. All returned products require a return merchandise authorization, RMA. Returns will not be accepted without this number.

- **Products must be in saleable condition.** Product that has been opened, used, written on or otherwise tampered with is not considered to be in saleable condition.
- Returns are subject to a restocking fee.

Please Note the Following Items are Non-Returnable:

- Refrigerated and hazardous items
- Equipment
- Narcotics and controlled substances
- Non-stock and drop-ship items
- Sterile solutions for injection and irrigation
- Products that are noted "non-returnable"

BACKORDERS

- Backorders older than six (6) months will be cancelled automatically.
- For questions regarding expired backorders, please contact Customer Service at **800.624.5926**.

DAMAGES/SHORTAGES

- Claims for missing or damaged items must be reported within 48 hours of receipt.
- Retain all packing material, boxes, and packing slips.
- Before reporting missing items, please examine the packing material. Occasionally product is lodged in the foam.

ARBITRATION

All Claims related to or arising under or relating to this Agreement are to be exclusively and finally determined by binding arbitration in the state of New York, or another location mutually agreeable to the parties. Any and all Claims must be arbitrated on an individual basis, and there shall be no right or authority for any Claims or disputes to be arbitrated on a class action or collective basis. For avoidance of doubt, each party irrevocably waives any right to: (i) have any Claim resolved in connection with any class action or collective action, or (ii) recover any damages or relief directly or indirectly as part of any class action or collective action. The arbitration shall be conducted on a confidential basis pursuant to the Commercial Arbitration Rules of the American Arbitration Association, or if applicable, under its Procedures for Large, Complex Commercial Disputes. Any decision or award as a result of any such arbitration proceeding shall be in writing and shall provide an explanation for all conclusions of law and fact and shall include the assessment of costs, expenses and reasonable attorneys' fees. Any such arbitration shall be conducted by an arbitrator experienced in the disputed subject matter and shall include a written record of the arbitration hearing. The parties reserve the right to object to any individual who shall be or has been at any time employed by or affiliated with a competing organization or entity. An award of arbitration may be confirmed in a court of competent jurisdiction. To the extent that any Claim or dispute is determined to not be subject to arbitration, all other Claims or disputes that would otherwise be subject to arbitration must be arbitrated. As used in this Agreement, "Claims" shall mean any and all liabilities, disputes and expenses whatsoever including, without limitation, claims, adversary proceedings (whether before a court, administrative agency or any other tribunal), damages (whether compensatory, multiple, exemplary or punitive), judgments, awards, penalties, settlements, investigations, costs, responses to subpoenas or other governmental directives and reasonable attorneys' fees and disbursements with respect to any claims that may be sustained, suffered or incurred by a Party hereto.

DISCOUNTS, REBATES, AND DISCLOSURES

Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion or discount program. You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.



INDEX

FORMING AN ALLIANCE

Rewarding Loyalty



ALLIANCE PROGRAM

Introducing the ACE SOUTHERN Alliance Program (ASAP)

When you create an alliance with **ACE SOUTHERN** you have a committed partner in your success. Part of that commitment is to reward you **for your loyalty with significant savings on every dollar you spend.**

As an ASAP Qualified Customer, you will benefit from savings and special promotions across all the products **ACE SOUTHERN** offers to help your practice prosper.

	Premier \$15K Annual Spend	Premier Plus \$25K Annual Spend*
5% savings on all Rx, Supplies and Equipment	●	
10% savings on all Instruments	●	●
15-20% savings on Bone Grafting Materials and Implants	●	●
10% savings on all Rx, Supplies and Equipment		●

*AAOMS/AAP members only

Become an ASAP Member by Enrolling • acesouthern.com/ASAP

I'm ready to buy.
Which Price Tier Column Do I Use?

Regular Price

Most customers fall into this category. It's not "suggested retail", but reflects the strong buying power of **ACE SOUTHERN**, which we pass along to you as **everyday low prices.**

Member Price

In partnership with the **AAOMS** and **AAP** societies, their members are given a special 5% discount on general merchandise, including Rx, supplies and equipment.

Premier Price

As shown above, **enrolled members** who have met the required annual spend of \$15,000 are qualified to receive the **Premier Discount.**

Premier+ Price

As a special reward for our **enrolled society members**, once they have met the required annual spend of \$25,000 they will receive the **Premier+ Discount.**

- 1 PHARMACEUTICALS
- 2 ANESTHESIA SUPPLIES
EMERGENCY RESPONSE
- 3 BONE GRAFTING
- 4 COSMETIC AND RESTORATIVE DENTISTRY
- 5 DISPOSABLES
- 6 HANDPIECES AND BURS
- 7 IMPRESSION MATERIAL
- 8 INFECTION CONTROL
- 9 INSTRUMENTS
- 10 IV ADMINISTRATION
- 11 PATIENT DIAGNOSTIC
- 12 STERILIZATION
- 13 SUCTION / IRRIGATION
- 14 SYRINGE AND NEEDLE
- 15 UNIFORMS
- 16 WOUND CLOSURE / SUTURE
- 17 WOUND MANAGEMENT
- 18 X-RAY
- 19 MISCELLANEOUS
- 20 INDEX

COLOR-CODED SECTION REFERENCE

MAKES FINDING WHAT YOU NEED A BREEZE.

Simply identify the section you want and its color, match it to the colors on the outside edges of the pages, and flip right to your section!

BE THE OFFICE HERO

Save Time and Money with...



REPLENISH

Our **REPLENISH** program allows you to set up scheduled deliveries for the products your practice uses on a regular basis. Ensure that your practice never runs out of essential items that keep your office running smoothly.

The ACE SOUTHERN Replenish Program, easy as 1-2-3

- 1 SHOP** from hundreds of eligible **Replenish** products
- 2 CHOOSE** your delivery frequency; every week, every month, **you're in control**
- 3 SAVE** an additional 3% on all **Replenish**-eligible products

Easy to Manage

Review, revise or skip orders at your convenience through your online account page. There, you can manage your Replenish items, change the quantities, change the frequency or update products. Prior to each delivery, we will send you an email showing the products that will be shipped. Put the Replenish program to work for you today and see how easy it is to save time and money.



For more details, scan the qr-code, or go to acesouthern.com/replenish



ACE SOUTHERN
Surgical Solutions

800.624.5926 | 800.441.3100

info@acesouthern.com

[ACESOUTHERN.COM](https://www.acesouthern.com) |   



PRESORTED
BOUND PRINTED
MATERIAL
US POSTAGE PAID
Columbia, SC
Permit No. 332

ACE SOUTHERN
One Southern Court
West Columbia, SC 29169

AN ALLIANCE

Built for Your Success

TRUST



ACE Surgical Supply Company and **Southern Anesthesia & Surgical, Inc.**
have come together as **ACE SOUTHERN**. With 90+ years of combined industry experience
ACE SOUTHERN will continue the legacy as a leading and trusted partner of dental surgical solutions.



800.624.5926 | 800.441.3100
info@acesouthern.com
acesouthern.com |   

