

HIGH PERFORMANCE FOR IMPROVED HEALING

Provides Strength and Permeability




collagen membrane

- Natural Collagen Structure Easily Conforms to Defect
- Dual Hydrophilic Surfaces to Support Osteogenesis
- Resists Tearing During Shaping and Tacking



perFORM COLLAGEN MEMBRANE

Activate regeneration! perFORM Collagen Membranes have dual permeable layers. The layers contain highly purified collagen fibers that act to support new bone growth. On the outer layer, densely packed collagen bundles prevent ingrowth of epithelial cells into the repair site. The perFORM membrane is fully resorbable with sufficient time for osteogenesis to occur and requires no second surgery to remove.



HIGH PERFORMANCE

Strength

- Easily cut wet or dry
- May be sutured or tacked without tearing

Hydrophilic


- Permeable dual layer design allows passage of fluids
- Loosely packed collagen fibers on the rough side support new bone growth
- Densely packed fibers on the smooth side impede tissue growth into the site

Conformable

- Easily conforms to the contours of the defect

perFORM is indicated for use in:

- Augmentation around implants placed in immediate extraction sockets
- Augmentation around implants placed in delayed extraction sockets
- Filling of bone defects after root resection or removal of retained teeth
- Guided Bone Regeneration (GBR) in dehiscence defects
- Guided Tissue Regeneration (GTR) procedures in intra-bony periodontal defects



perFORM				Regular	Member	Premier	Premier+
PERF1520	15 x 20mm	each		\$156.99	\$156.99	\$133.44	\$133.44
PERF2030	20 x 30mm	each		\$193.49	\$193.49	\$164.47	\$164.47
PERF3040	30 x 40mm	each		\$311.49	\$311.49	\$264.77	\$264.77

Shown Actual Size



Copyright © 2024 ACE SOUTHERN

perFORM is a registered trademark of ACE Surgical.
perFORM is manufactured for ACE Surgical.

All prices are subject to change without notice. Prices shown are domestic prices in U.S. dollars; does not include shipping or any duties. For discount details visit acesouthern.com/ASAP. And for the most up-to-date pricing, please contact ACE SOUTHERN or visit our website.