

DENTAL SURGICAL INSTRUMENTS

CATALOG





BE THE OFFICE HERO



The Replenish Program, easy as 1-2-3

Our **REPLENISH** program allows you to set up scheduled deliveries for the products your practice uses on a regular basis. Ensure that your practice never runs out of essential items that keep your office running smoothly.





SAVE an additional 3% on all **Replenish**-eligible products

For more details, scan the qr-code or visit us at...

acesouthern.com/lp/replenish



W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

ACE SOUTHERN

"NEED SOMETHING ELSE?"

In addition to ACE Surgical, our preferred brand of surgical instruments, ACE SOUTHERN also carries a wide selection from other top brands such as A-Titan, Hu-Friedy and Miltex.

Visit acesouthern.com today

TABLE OF CONTENTS

AUTOGEN	IOUS BONE COLLECTION	4
CLEANING	G AND STORAGE	8
CUTTING	AND LIFTING	10
HINGED I	NSTRUMENTS	16
MEASURI	NG INSTRUMENTS	22
MIRRORS	S	25
MOUTH P	ROPS AND RETRACTORS	26
SUCTION	AND IRRIGATION	28
SURGICA	L KITS	34
WIRE AND	WIRE ACCESSORIES	40
MISCELLA	ANFOLIS INSTRUMENTS	41

"WHERE'S THE PRICING?"

We want to ensure you have the most current pricing available, but understand sending multiple catalogs a year in order to achieve that has a negative impact on the environment. Therefore, in order to create less waste, we have eliminated pricing from this catalog. Following the link below, or scanning the QR code, will allow you to view current retail pricing. In order to view pricing specific to your account, you must log in to your account on our website. This will also provide real-time product availability information. Plus, ordering is now just a click away!

Simply select the page number or sku in the index and begin!

Visit acesouthern.com/catalog/instruments24 today



SOCIETY MEMBERS SAVE 5% ON RX, SUPPLIES AND EQUIPMENT





TRUST - CONVENIENCE - VALUE













ACE SOUTHERN's retail sales operation has integrated with Henry Schein's North America Distribution Group to create a unified sales team, specifically focused on the needs of oral surgeons and periodontists.

TRUSTED SOURCE FOR BONE REGENERATION PRODUCTS

ACE Surgical is one of the only dental companies to receive accreditation from the AATB® (American Association of Tissue Banks). This means ACE Surgical operates under the stringent guidelines set forth by the AATB to process and provide safe quality human tissue for ACE SOUTHERN's flagship alloOss® product portfolio.

Additionally, ACE SOUTHERN offers a complete portfolio of regenerative products including NuOss® xenograft, an array of resorbable membranes, and bone viable matrix products with growth factors.

ACESOUTHERN.COM

The online store at ACESOUTHERN.com is more than just easy online ordering. You can manage your account, set up scheduled deliveries with our Replenish Program, and receive exclusive online offers.

- Order anytime, 24/7/365
- Save favorites for reordering
- Order from history
- View order history
- Check shipping status
- Pay invoices

REPLENISH PROGRAM

Imagine your office running smoothly, never running out of the essentials used the most. This can be a reality with our Replenish Program. Set up scheduled deliveries for the products your practice uses the most and sit back, save, and be the hero of your office.

You can shop from hundreds of eligible products. Just look for this symbol on ACESOUTHERN.com. Once enrolled, you are in complete control - easily manage your Replenish items, change quantities, change frequency, and add or update products. And don't forget, Replenish eligible products receive an additional savings!

PRODUCT EXPERTS

Our legacy is founded by an oral surgeon and a pharmacist who put patients and practices first. ACE SOUTHERN has invested in industry expertise in regenerative medicine and pharmaceutical compliance. Our team is focused on making sure that we offer the best products at competitive prices. You will enjoy knowledgeable customer service associates and access to dedicated inside and outside sales representatives to provide support with a call, click or text.


ASSOCIATION SUPPORT

AAOMS and AAP members who purchase from ACE SOUTHERN receive a 5% upfront discount on supplies, pharmaceuticals, and equipment. In addition, a portion of each purchase made by AAOMS and AAP members with ACE SOUTHERN is paid to AAOMS and AAP to invest in programs that provide added member benefits. Our goal is to be a partner in the success of dental surgical practices nationwide.





AUTOGENOUS BONE COLLECTION





Simple method for obtaining autogenous bone

The design of the ACE AUTOgrafter allows the surgeon to filter and retain particles of bone which have been suctioned along with standard irrigation during surgery. The graft may be obtained from almost any bony site; examples include dental implant preparation sites, the mandible, maxilla, calvarium, etc. The collected graft material can be used either in a wet amorphous form, or a compressed osseous particulate. The drier particulate, which

is easier to handle, can be obtained by compressing the collected material in a screen with the titanium plunger. This forces out excess moisture through the screen.

NOTE: The ACE AUTOgrafter is not a routine suction device; therefore, bleeding should be controlled prior to use. The ilium should not be used as clots will clog this device.

tem lumber	Manufacturer Number		Description/Details	Quantity
Autografter				
092179	0034215		Autografter Start Up Kit, (includes: A, C, D and G)	1 kit
466760	0034212	A	ACE Autografter, for HVAC tip, autoclavable	each
212231	0034210	В	ACE Autografter, inline suction hose connection, autoclavable	each
\utoarafter	Accessories			
466733	0900024	C	HVAC Suction Tip, 4mm tip, 45° angle, 7" (18cm)	each
476175	0034211	D	Titanium Plunger, 4-1/4", 11cm	each
150072	9600449	Ε	Square Titanium Bowl, 12mm deep, 50 x 50mm	each
150073	1000216	F	Implant Curette with round cups, titanium, 5-1/2" (14cm) cups are 5 x 10mm & 6 x 20mm	each
092178	6462002	G	Screen Assembly, disposable, sterile, 12cc capacity	10/box
464526	0034238		Suction Connection Tubing, sterile, 3' (.9m)	20/case
464527	0034239		Suction Connection Tubing, sterile, 6' (1.8m)	20/case

 $\hbox{Autografter}^{\scriptscriptstyle{\text{TM}}} \hbox{ is a trademark of ACE Surgical Supply.}$



The Bone Crusher introduces an alternative to the Bone Mill. Insert the bone piece into the crusher and use the pestle in order to crush the bone. If necessary, you may additionally use a mallet.

Stainless steel

/s Item Number	Manufacturer Number	Description/Details	Quantity
1150092	9600050	20mm crushing surface	each



ACE BONE FILE

∕ Item Number	Manufacturer Number	Description/Details	Quantity
1463463	6400009	#21 Miller Bone File 7-1/4"	each



Autografted Bone Histology of Compressed Cortico/Cancellous Bone

NEED PRICES? Follow the link below, or scan the QR code, to view retail pricing.

AUTOGENOUS BONE COLLECTION

ACE SURGICAL BURS

Shank images (shown to the right) and drill tip images (shown below) have been enlarged to show detail. See descriptions for actual sizes.




CARBIDE	/ Item Number	Manufacturer Number		Head Material	Head Diameter	Bur Number	Total Length	Cutting Length	Features	Qty	6
POUND PURO	1463199	21441980	HP	Carbide	2.3 mm	8	44 mm	-	Round	10	•
ROUND BURS	1466699	21501980	HP	Carbide	2.3 mm	8	50 mm	-	Round	10	•
	1463426	21651980	HP	Carbide	2.3 mm	8	65 mm	-	Round	10	•
											_
702 TAPER BURS	1463200	61441902	HP	Carbide	1.6 mm	702	44 mm	5 mm	Tapered	10	•
<u>~</u>	1463008	61501902	HP	Carbide	1.6 mm	702	50 mm	5 mm	Tapered	10	
쁜	1463027	61651902	HP	Carbide	1.6 mm	702	65 mm	5 mm	Tapered	10	•
S											
703 TAPER BURS	1463221	61441909	HP	Carbide	2.1 mm	703	44 mm	5 mm	Tapered	10	•
Ш	1463209	61501909	HP	Carbide	2.1 mm	703	50 mm	5 mm	Tapered	10	•
A A	1463427	61651909	HP	Carbide	2.1 mm	703	65 mm	5 mm	Tapered	10	•
X											
1703L TAPER	1466739	63441903	HP	Carbide	2.1 mm	1703	44 mm	5 mm	Tapered	10	•
DOME BURS	1463195	63501903	HP	Carbide	2.1 mm	1703	50 mm	5 mm	Tapered	10	•
	1.400000	04441000	ш	0	0.1	17001	44	7	T	10	
1703L TAPER	1463029	64441903	HP	Carbide	2.1 mm	1703L	44 mm	7 mm	Tapered	10	
DOME BURS	1466883	64501903	HP	Carbide	2.1 mm	1703L	50 mm	7 mm	Tapered	10	•
	1466866	64651903	HP	Carbide	2.1 mm	1703L	65 mm	7 mm	Tapered	10	•
ROUND BURS	1466871	21440080	НР	Carbide	2.3 mm	8	44 mm	-	Round	10	
702 TAPER BURS	1466858	61440002	НР	Carbide	1.6 mm	702	44 mm	5 mm	Tapered	10	
		1									
703 TAPER BURS	1466860	61440003	HP	Carbide	2.1 mm	703	44 mm	5 mm	Tapered	10	
CYLINDER BUR	1466830	40500016	HP	Carbide	1.5 mm	-	50 mm	7.5 mm	-	10	•
OVAL BONE	1416845	504080501	HP	Carbide	4 mm	-	48 mm	8 mm	Cross Cut	1	•
CUTTING BURS	1463010	31680031	HP	Carbide	4 mm	-	70 mm	8 mm	Cross Cut	1	
	1463009	31440023	HP	Carbide	5.5 mm	-	44 mm	10 mm	Cross Cut	1	
	1463013	31500023	HP	Carbide	5.5 mm	-	50 mm	10 mm	Cross Cut	1	
	1466743	504080505	HP	Carbide	5.5 mm	-	70 mm	10 mm	Cross Cut	1	
	1466863	31440020S	SN	Carbide	4 mm	_	44 mm	8 mm	Cross Cut	1	Ť
	1466718	504080501S	SN	Carbide	4 mm	-	48 mm	8 mm	Cross Cut	1	
	1466855	504502ST	ST	Carbide	4 mm	_	54 mm	8 mm	Cross Cut	1	
	1466828	31440023ST	ST	Carbide	5.5 mm	_	44 mm	10 mm	Cross Cut	1	
	1466744	31500023ST	ST	Carbide	5.5 mm	_	50 mm	10 mm	Cross Cut	1	_
				Head	Head				01000 001		
SPECIALTY	/ Item Number	Manufacturer Number		Material	Diameter	Bur Number	Total Length	Cutting Length	Features	Qty	
LINDEMAN BURS	1463222	10440166	HP	SS	2 mm	-	44 mm	10 mm	-	1	•
LINDLINAN DONS	1463146	10600168	HP	SS	2 mm	-	60 mm	24 mm	-	1	•
								-			
PEAR BUR	1463212	00500013	HP	SS	6 mm	-	50 mm	12 mm	Contouring	1	
										-	
CONE BUR	1463428	00500015	HP	SS	4 mm	-	50 mm	14 mm	Contouring	1	
2 -STEP										-	
PILOT DRILL	1260097	454281001	CA	SS	0.81-1.6 mm	-	30 mm	-	-	1	•
	1467036	45440202	CA	SS	1.1 mm	_	29.8 mm	_	Twist	1	
PILOT TWIST DRILLS		502700045	CA	SS	0.45 mm	-	27 mm	-	Twist	1	•
DIAMOND	/> Item Number	Manufacturer Number	OA.	Head Material	Head Diameter	Bur	Total	Cutting Length	Features	Qty	
	IVUITIDEI	Mannaer				Number	Length		Features		
ROUND DIAMOND RUBS	4.400004	00440040									_
ROUND DIAMOND BURS	1469984	90440018	HP	Diamond	1.8 mm	6	44 mm	-	Round	1	
ROUND DIAMOND BURS	1469984 1468859	90440018 90440023	HP HP	Diamond Diamond	1.8 mm 2.3 mm	8	44 mm	-	Round	1	•

AUTOGENOUS BONE COLLECTION

ACE TREPHINE BURS - CONTRA-ANGLE

- For harvesting bone in autogenous bone graft placement
- May be used for retrieval of implants
- Laser marked score lines
- Manufactured from Titanium Alloy (Ti-6Al-4V ELI)
- Internally irrigated to minimize the risk of overheating the bone during drilling
- Can be used with the Drill Guide Paralleling System to eliminate drill chatter during graft or implant harvesting



/ ≥ Item	Manufacturer		I.D.	0.D.	Length	Cutting		
Number	Number		mm	mm	mm	Length	Markings	Quantity
1150104	04948001	A	2	2.8	35	18.5 mm	8,10,13,15,18 mm	each
1150105	04948002	В	2.75	3.5	35	18.5 mm	8,10,13,15,18 mm	each
1150106	04948201	C	3.75	4.6	35	18.5 mm	8,10,13,15,18 mm	each
1150107	04948202	D	4.25	5.10	35	18.5 mm	8,10,13,15,18 mm	each
1150108	04948401	E	5.25	6.10	32	15.5 mm	8,10,13,15 mm	each
1150109	04948501	F	6	6.9	32	15.7 mm	8,10,13,15 mm	each
1150110	04948502	G	7	7.9	32	15.8 mm	8,10,13,15 mm	each
1150111	04948701	Н	8	8.9	32	14 mm	8,10,13 mm	each

Suggested	Trephines	To Use	For Retrieve	Failed Implants	For Implants up to
					. or mipianto up tomi

3.3 mm	use 1150106	3.5–4.0 mm use 1150107	4.5–5.0 mm use 1150108	5.0–6.0 mm use 1150109
--------	-------------	-------------------------------	-------------------------------	-------------------------------



SAVE an additional 3% with scheduled deliveries on o marked items in the column. See inside front cover for details.

W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

acesouthern.com/catalog/instruments24/page/6

CLEANING AND STORAGE



ACE CLEANING BRUSHES

/> Item Number	Manufacturer Number		Description/Details	Quantity				
	Bur Brush - stainless steel Autoclavable							
1260234	00200106	A	Stainless Steel Bristles, 5-1/4" (13.5cm)	3/pack				
Cleaning Brush - plastic handle Not Autoclavable								
1463136	6039015	В	Nylon Bristles, 7" (18cm)	3/pack				
1463137	6039016	C	Stainless Steel Bristles, 7" (18cm)	3/pack				

ACE SQUARE TITANIUM BOWL

1150072	9600449	D	50 x 50mm, 12mm depth	Quantity each
/s Item Number	Manufacturer Number		Description/Details	Quantit





ACE INSTRUMENT TRAYS

Number	Manufacturer Number	Description/Details	Quantity
Mayo Tra	ys		
LustroiStainle	us, heavy-duty ins ss Steel	trument trays	
1464275	0531024	E 19"x12" x 5%"	each

Polymer (Plastic) Instrument Trays Tray with Finger Mat

 Autoclave up to 270° F 			
1396966 7651001		6" x 2-½" x ¾"	each
1469396 7651002	F	6-½" x 4", ¾"	each
1469397 7651011		10" x 6" x ¾"	each

ACE STAINLESS STEEL WARE

- Heavy Gauge, 18/8, 304 stainless steel throughout
- Seamless, crevice-free, sanitary construction
- Uniform flat bottoms
- Contoured for nursing efficiency and patient comfort
- Permanent, brilliant mirror finish

/> Item Number	Manufacturer Number		Description/Details	Quantity
Round Corn	er Instrument 1	Tray		
1466712	0531030	Α	8-%" x 5" x 2"	each

Forceps Jar - Heavy Base

1466887	0531027	В	4- ⁷ / ₁₆ " x 2"	each	1

Needle/Ointment Jar

1466849 0531008 C	1/4 qt.	3" x 2-½"	each
-------------------	---------	-----------	------

Dressing Jar - Recessed Handle

Sponge Bowls

1464362	0531009		1 qt. 5-¾" x 2-½"	each
1466872	0531010	Ε	1-3/8" qt. 6-1/2" x 3"	each

Medicine Cup, Graduated

1463233 0531005	F	2 oz., 4 tsp.	each
lodine or Oil Cup			
1464361 0531007	G	6 oz.	each

Tumbler

1464359	0531006	н	7 oz.	each

Emesis Basin

1464363	0531033	1	6" x 3-½" x 1-½"	each
1466715	0531034		8" x 3-¾" x 1-5/8"	each

CLEANING AND STORAGE















CUTTING AND LIFTING





BLADE HANDLES

/> Item Number	Manufacturer Number	Description/Details	Quantity
#3 Flat Bla	nde Handles		
 For surgion 	cal blades 10 through	15	
1463163	0030992	A 5" with metric markings	each
1463174	0030993	5"	each
	•Style Blade Hand cal blades 10 through 0030995		each
#7 Flat Bla	ade Handle		
 For surgi 	cal blades 10 through	ı 15	
1466790	0030997	C 6-½" straight	each

ACE STERILE SURGICAL BLADES 🥯

/> Item Number	Manufacturer Number	Description/Details	Quantity
Carbon St	eel		
1466975	0036112	#12	100/box
1463128	0036312	#12D	100/box
1463127	0036115	#15*	100/box
1463129	0036315	#15C	100/box
* the tin of t	the 15* hlade is sligh	atly different than a standard 15	

trie tip oi	me	10	blade is	Silgrilly	amerent	man a	stanuaru	IJ	

Stanness Steel								
	1463073	0037215	#15*	100/box				
	1463130	0037315	#15C	100/box				

^{*} the tip of the 15* blade is slightly different than a standard 15

ACE STERILE MINI BLADES 🥯

/s Item Number	Manufacturer Number	Description/Details	Quantity
Mini Stainl	ess Steel		
1466740	7866400	#64	10/pkg
1466800	7866400	#67	10/pkg
1463139	7866400	#69	10/pkg



SAVE an additional 3% with scheduled deliveries on marked items. See inside front cover for details.

ACE MUCOSA TISSUE PUNCHES

Reusable Contra Angle Punches

• 1" (25mm)


/s Item Number	Manufacturer Number	Description/Details	Quantit
1260101	9200002	2.0mm diameter	each
1150119	9200003	3.0mm diameter	each
1150120	9200004	4.0mm diameter	each
1150121	9200005	5.0mm diameter	each
1150122	9200006	6.0mm diameter	each

2mm	3mm	4mm	5mm	6mm	
•	•				

ACE PERIOTOME HANDLE & BLADES KIT

- Periotome is used to remove ligament attachment from the root of the tooth and the alveolar housing
- Single end handle can be tapped with a mallet or used with hand pressure
- Includes 6-Blade Storage Holder: stands up for easy access during procedures and folds flat for sterilization and storage

Manufacturer




Number	Number		Description/Details	Quantit
Periotome	Kit			
1466810	3441770		Periotome Kit	kit
Periotome	Kit Includes:	(ava	uilable separately)	
1467005	6400312S	Α	Surgical Mallet, 6-1/2" (16.5cm)	each
1476173	16901001	В	Single-End Handle only, 6" (15cm)	each
1466802	3441651	C	P3, 1.5mm, angled, 1.75" (4.5cm)	each
1463167	3441631	D	P1, 2mm, straight, 1.5" (4cm)	each
1466772	3441641	Е	P2, 1.5mm, straight, 1.5" (4cm)	each
1476174	3441691	F	2.5mm, straight, semi rigid, 1.5" (4cm)	each
1466822	3441661	G	Arrow Tip	each
			Blade Storage Holder, 6" (15.2cm), not sold individually	
			Instrument Tray, 175 x 140 x 24mm, not sold individually	

CUTTING AND LIFTING







In Use: Periotome is removing the ligament attachment from the root of the tooth and the alveolar housing.



CUTTING AND LIFTING CUTTING AND LIFTING



ACE SINUS LIFT KIT

/> Item Number	Manufacturer Number		Description/Details	Size	Quantity
Sinus Lift k	Cit				
1150080	0034194	Α	Sinus Lift Kit with Tray		each
Kit Include	s:				
1150091	7651050	В	Autoclavable Tray		each
1150081	0311088	C	Shepherd Sinus Lift Retractor	7" (18cm)	each
1150082	0311023	D	Freer Periosteal Elevator, strong curve	7-½" (19cm)	each
1150084	0311077	Е	#3 Clepper Tatum Antralplasty Elevator, DE	6-1/2" (16.5cm)	each
1150085	1000236	F	#1 Sinus Lift Curette, standard bend, DE, 3mm/5mm	7-½" (19cm)	each
1150086	1000237	G	#2 Sinus Lift Curette, severe bend, DE, 3mm/5mm	6-1/2" (16.5cm)	each
1150088	7140243	Н	Depth/Packing Gauge, DE, Titanium,	6-1/2" (16.5cm)	each
			markings at 8-10-13-15-18mm, 3.3/4mm ends		

ACE ELEVATOR KIT

/> Item Number

Kit Includes: 1150084

1466879

1466812

1466877

1476168

1466880

Antralplasty Elevator Kit for Sinus and Socket Lift Procedures

- Designed for exploring and walking antral floor in displaced root tip removal
- Instruments have dulled edges to prevent tearing of the membrane

Manufacturer Number

0311077

1000231


1000225

1000234

1000230

1000233

Antralplasty Elevator Kit 1466763 0000220



В

C

W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

For account specific pricing, please log into your account on our website.


acesouthern.com/catalog/instruments24/page/13



handle

CUTTING AND LIFTING

CUTTING AND LIFTING



ACE PERIO & ORAL SURGERY CURETTES

TC = Tungsten Carbide

Ti = Titaniui DE = Double					
/s Item Number	Manufacturer Number		Description/Details	Size	Quantity
Implant Cu	rette				
1150073	1000216	A	Ti, 5 x 10mm & 6 x 20mm, round cups	5-1/2"	each
Antral Cure	ette				
1466812	1000225	В	#1, double ended	6-1/2"	each
Lucas Cure	otto				
1466827	1000212	C	#88, double ended	7"	each
1466808	1000224	D	Serrated, angled, 3.5mm, 5 teeth	6-1/2"	each
1463164	1000226	Ε	Serrated, angled, 4.5mm tips, 7 teeth	6-1/2"	each
Miller Cure	ette				
1463168	1000203	F	#9, double ended	7-1/4"	each
1463169	1000204	G	#10, double ended	7-1/4"	each
1463170	1000205	Н	#11, double ended	7"	each
1466770	1000206	1	#12, double ended	7"	each
Molt Curet	to				
1463166	1000219	J	#2, single ended	6"	each
1463165	1000220	K	#4, single ended	6"	each

#2/4, double ended



1000221

1276209

SAVE an additional 3% with scheduled deliveries on \bullet marked items in the \bowtie column. See inside front cover for details.

ACE PERIO & ORAL SURGERY ELEVATORS

Number	Manufacturer Number		Description/Details	Size	Quantity
Elevators					
1466890	6300004	A	#301 Apexo	6"	each
1466870	6300008	В	#303 Apexo	6"	each
1466882	6300018	C	#8 Crane Pick	6"	each
1463423	6300011		#46R Heidbrink Elevator	6"	each
1466878	6300014		#34 Seldin	6"	each
1466813	6300140	D	#34S Seldin	6"	each

T-Bar Handle Elevators

1466706	6300063	Ε	#6X Potts, T-bar handle	6"	each
1466707	6300064	F	#7X Potts, T-bar handle	6"	each

Periosteal Elevators

1466881	0300080	G	#9A Allen	6-1/2"	each
1463420	0300139	Н	Buser	7"	each
1150084	0311077		#3 Clepper Tatum	6-1/2"	each
1150082	0311023	1	Freer, strong curve	7-1/2"	each
1463107	0300073	J	Molt #9, double ended	7"	each
1463447	0314042		Molt, double ended, SB	7-1/2"	each
1466756	0300072	K	Molt #9, double ended, titanium	7-1/2"	each
1463429	0300001		#23 Seldin, flat handle, double ended	7-¾"	each
1466890	0300004	L	#23 Seldin, octagonal handle, double ended	7-1/2"	each
1466758	0300065	M	Pritchard, double ended	7-¾"	each
1466875	0300016	N	#24G Raspatory straight, double ended	6-¾"	each

Root Tip Elevators / Picks

1463439	6315130	Q	Double-ended, serrated	7-¾"	each
1463116	6300118		#13/14 Heidbrink	7-1/2"	each
1463117	6300115	P	#3 Heidbrink	6-1/2"	each
1463114	6300114	0	#2 Heidbrink	6-1/2"	each
1466737	6300113		#1 Heidbrink	6-1/2"	each

NEED PRICES? Follow the link below, or scan the QR code, to view retail pricing.











W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

HINGED INSTRUMENTS



ACE MICROSURGICAL SCISSORS

/> Item Number	Manufacturer Number		Description/Details	Size	Quantity
Westcott					
1466970	4102042	Α	Tenotomy, curved, sharp	5-½"	each

ACE SURGICAL SCISSORS

/> Item Number	Manufacturer Number	Description/Details	Size	Quantity
Sullivan				
1466991	4000009	B #15 Sullivan, smooth	5-½"	each
_				
Dean		_		
1463400	4000001	C #9, smooth	6-1/2"	each
1463432	4015002	#9, serrated	6-1/2"	each
1467020	4005201	#9, serrated, tungsten carbide	6-3/4"	each
Iris				
1463419	4001019	D #17, straight, delicate,	4-1/2"	each
1466982	4002020	#18, curved, delicate	4-1/2"	each
Kelly				
•			0.4/11	. 1
1463113	4002005	E #1 curved, smooth	6-1/4"	each
1467004	4002006	#1 curved, serrated	6-1/4"	each
1467021	4005203	Curved, tungsten carbide	6-1/4"	each
Locklin				

ACE CASTROVIEJO NEEDLE HOLDERS

/ ⊱ Item	Manufacturer				
Number	Number		Description/Details	Size	Quantity
1467016	3105206	Α	Serrated, straight	5-½"	each
1466998	3105110		Serrated, straight, narrow, tungsten carbide	5-1/2"	each
1466980	3105010		Serrated, straight, broad, tungsten carbide	5-1/2"	each
1467008	3105016		Serrated, straight, broad, tungsten carbide	7"	each
1467010	3105018	В	Serrated, curved, broad, tungsten carbide	7"	each

ACE SURGICAL NEEDLE HOLDERS

/> Item Number	Manufacturer Number		Description/Details	Size	Quantity
Crile Wood	d				
1467014	3001012	C	Straight, serrated	6"	each
1463179	3005117		Straight, tungsten carbide	6"	each
1463178	3005118	D	Straight, serrated, tungsten carbide	7"	each
Mayo Heg	ar				
1467018	3001019		Straight, serrated	7"	each
1463112	3005109	Е	Straight, serrated, tungsten carbide	6 1/4"	each
1466972	3005110		Straight, serrated, tungsten carbide	7"	each
Webster					

HINGED INSTRUMENTS





HINGED INSTRUMENTS

ACE FORCEPS, PLIERS AND CLAMPS

	•				
/s Item Number	Manufacturer Number		Description/Details	Size	Quantit
Bone Graft	Holding Forceps	;			
1150077	9600026		Serrated, with open-end drill-relief tips	8"	each
1150078	9600027	A	Serrated, with open-end drill-relief tips	6"	each
Corn Sutur	e Pliers				
1276403	2000107	В	Angled	6"	each
Dressing F	orceps/Pliers				
1466973	2002145	C	#17 Dressing	6-1⁄4"	each
Hemostatio	Forceps				
1463435	2002202	D	#3 Halsted Mosquito, curved, serrated	4-¾"	each
1467033	2002204		Halsted Mosquito, curved, extra light tips	4-3/4"	each
1463120	2002221	Ε	Providence Kelly, curved	5-1/2"	each
1467011	2002223	F	Rankin Kelly, curved	6-1/4"	each
1463115	2002225	G	Crile, curved	5-1/2"	each
Tissue Fore	ceps				
1466983	2015115	J	Adson #41, serrated	4-3/4"	each
1466943	2000118		Adson #42, 1x2 teeth	4-¾"	each
1466941	2000119	K	Adson, 1x2 teeth	6"	each
1463444	2115037	L	Micro Adson, 1x2 teeth, straight	4-3/4"	each
(continued n	ext page)				



HINGED INSTRUMENTS

ACE FO	RCEPS, PLIE	RS AND CLAMPS		
/> Item Number	Manufacturer Number	Description/Details	Size	Quantity
Tissue Ford	ceps (continued)			
1463119	2000048 A	Debakey, 2mm wide	6"	each
1467009	2000049	Debakey, 2mm wide	7-3/4"	each
1463438	2000104 B	Gerald, 1x2 teeth, straight	7"	each
1463422	2015039	Gerald, straight, serrated	7"	each
Tissue Ford	eps - Russian			
1463393	2001108 C	Russian Tissue Forceps	6"	each
1466977	2001109	Russian Tissue Forceps	8"	each
Towel Forc	eps/Clamps			
1467032	2000305 D	Backhaus Towel Forceps	3-1/2"	each
1467031	2000306	Backhaus Towel Forceps	5"	each
Magill Cath	eter Forceps			
1463437	4419051 E	Magill Forceps, child	7-1/2"	each
1463173	4423050	Magill Forceps, adult	9-3/4"	each
Chalazion I	Lip Forceps			
1463188	0117148 F	Medium, 1" wide	3-1/2"	each
		Medium, 1" wide	3-1/2"	eac



W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

HINGED INSTRUMENTS

ACE Handle Reference

ACE EXTRACTION FORCEPS

Manufacturer Number		Description/Details	Size	Quantity
6300301	Α	#1 Standard, upper anterior	7"	each
6300453	В	#89 Upper Cowhorn	7"	each
6300452	С	#90 Upper Cowhorn	7"	each
6300444	D	#150 Upper Molar, smooth	7"	each
6315309	E	#150 Upper Molar, serrated	7"	each
6316473	F	#150S Upper	6"	each
6300345	G	#23 Lower Molar Forceps	6-1/2"	each
6300320	Н	#151 Lower Molar Forceps	6-1/2"	each
6316472	1	#151S Lower Molar Forceps	6"	each
6316480	J	#23C Cowhorn Lower Molar, pediatric	4-1/2"	each
	Number 6300301 6300453 6300452 6300444 6315309 6316473 6300345 6300320 6316472	Number 6300301 A 6300453 B 6300452 C 6300444 D 6315309 E 6316473 F 6300345 G 6300320 H 6316472 I	Number Description/Details 6300301 A #1 Standard, upper anterior 6300453 B #89 Upper Cowhorn 6300452 C #90 Upper Cowhorn 6300444 D #150 Upper Molar, smooth 6315309 E #150 Upper Molar, serrated 6316473 F #150S Upper 6300345 G #23 Lower Molar Forceps 6300320 H #151 Lower Molar Forceps 6316472 I #151S Lower Molar Forceps	Number Description/Details Size 6300301 A #1 Standard, upper anterior 7" 6300453 B #89 Upper Cowhorn 7" 6300452 C #90 Upper Cowhorn 7" 6300444 D #150 Upper Molar, smooth 7" 6315309 E #150 Upper Molar, serrated 7" 6316473 F #150S Upper 6" 6300345 G #23 Lower Molar Forceps 6-½" 6300320 H #151 Lower Molar Forceps 6-½" 6316472 I #151S Lower Molar Forceps 6"



















HINGED INSTRUMENTS

ACE EXTRACTION RONGEURS

/> Item Number	Manufacturer Number		Description/Details	Size	Quantity
1463396	1000001	A	#1 Rongeur	6-1/2"	each
1466765	1000043	В	Beyer, double action	7"	each
1463399	1000009	C	Standard Blumenthal, 30° Angle	6"	each
1466703	1000021	D	#4 Cleveland	5-1/2"	each
1463433	1000034	E	Friedman	5-1⁄2"	each
1463434	1010035		Mini Friedman	5-½"	each
1463110	1010036	F	Micro Friedman	5-1/2"	each







MEASURING INSTRUMENTS ACE CALIPERS Manufacturer Number /> Item Number Description/Details Quantity Castroviejo Calipers - Short Beak 1463186 0800015 A 3-¼" 0-20mm, straight each 3-1/4" 0-40mm, straight 1466749 0800017 each 0800018 3-½" 0-40mm, 45° angle 1152858 each Castroviejo Calipers - Long Beak 0800019 D 6-1/2" 0-40mm, straight 1150100 each 6-1/2" 0-40mm, 45° angle 1150099 0800020 each **Bone Implant Caliper** • 0-20mm F 3-1/2" (9cm) 1152847 0800050 each **Ridge Mapping Caliper** • 0-25mm **1150101** 0800059 **G** 6" (15cm) each





MIRRORS


MEASURING INSTRUMENTS











ACE RULERS

Number	Number Steel Ruler	Description/Details	Quantity
/> Item	Manufacturer	Description/Details	Quantity

A ½" x 6"

Plastic Ruler

- Standard and metric increments
- Cold sterilization only

1463995 0800011

1463996	0032937	В	7∕%" x 6"	10/pkg

Patient Opening Gauge

- Allows patient to monitor progress at home
- Cold sterilization only

X-Ray Magnification Calculator™

- Specially designed for dental implantology and orthodontics
- Eliminates the complex calculations for determining panoramic x-ray distortion
- Measurements can easily be obtained to determine the precise magnification of the x-ray and to accurately determine the amount of bone available for implant placement
- Can be used with the Delrin Bite Blocks with Stainless Steel Pin (2080023, 2080024 and 2080021) or Exacta Mouth Prop® with 6mm Marking Ball (2080022)
- Actual size: 3-34" x 9-1/2", 9.5cm x 24cm
- · Cold sterilization only

2060025 0031967 D ACE X-Ray Magnification Calculator™	each
---	------

ACE EXACTA® MOUTH PROPS

Manufacturer

For use with #2060025 X-Ray Magnification Calculator (above)

Exacta®	Mouth Prop with	Calculator Pin	
 Exacta 	mouth prop has a	20mm pin	
2080021	0100456	E 20mm calculator pin	each

Description/Details

Exacta® Mouth Prop with 6mm Marking Ball

· Exacta mouth prop has a 6mm spherical marking ball

080022 0100450	F	6mm marking ball	each
-----------------------	---	------------------	------

Exacta® Mouth Prop with 5mm and 10mm Pin

• The size of the pin on a panoramic x-ray is comparable to the original pin length. The amount of elongation is determined and the exact measurement extrapolated, (If a 10mm pin increases to 12mm, a 20% elongation is realized – the x-ray bone height is reduced by 20%)

2080023 0100350 G 5mm and 10mm pins each

Exacta® Mouth Prop with 10mm Pin

- Allows simple and exact determination of vertical bone height in dental implantology
- Autoclavable
 Contains a 10mm pin
- For single tooth or narrow edentulous area

2080024	0100448	Н	10mm pin	each
---------	---------	---	----------	------


ACE INTRAORAL PHOTO MIRRORS

Chromium Coated, Front Surface - Not shown actual size

- Material: float glass, thickness 1/8" (3mm) ground edges; absolutely flat surfaces even on edges (unlike metallic mirrors)
- Coating: chromium coating on both sides by sputtering vacuum metallization
- Reflection: reflection capacity = 65%; sharp images without any deformation and without double reflection (front surface)
- Sterilization: autoclavable up to 428°F (220°C)
- Can be used on both sides

each

/> Item Number	Manufacturer Number		Description/Details	Quantity
Buccal In	traoral			
1466861	0032989	A	Adult, 6"	each
Occlusal I	ntraoral			
1466886	0032971	В	Adult, 4"	each
1466816	5961009	С	Oversized, 5-1/4"	each



ACE DENTAL MIRRORS Mirrors and Handles are sold separately





• Simple Stem and Cone Socket mirrors are not identical and are not interchangeable

Simple Stem mirrors go with Simple Stem handles

· Cone Socket mirrors go with Cone Socket handles

Number	Manufacturer Number		Description/Details	Quantity
Simple Ste	em Mirror			
1463118	9700012	D	Size 4, plane, 1/4" (22mm) diameter	each
1463171	9700014	D	Size 5, plane, 15/16" (24mm) diameter	each
Cone Sock	et Mirror			
1466741	9700013	Ε	Size 4, plane, 1/4" (22mm) diameter	each
1466742	9700015	Е	Size 5, plane, 15/16" (24mm) diameter	each

Simple Stem Mirror Handle

1463451 9700030 F 5", stainless steel eac	1463451
---	---------

Cone Socket Mirror Handle

1463436 97000	1 G	4-1/2", stainless steel	each
----------------------	------------	-------------------------	------







W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

MOUTH PROPS AND RETRACTORS



ACE MOUTH GAGS

Autoclavable

Manufacturer Number		Description/Details	Quantity
th Gags			
0117016		Adult 5-1/2"	each
0116018	Α	Child 4-1/2"	each
ent Covers for M	/lolt l	Mouth Gags	
0117023	В	Child/adult, contains latex	2/pack
	Number th Gags 0117016 0116018 ent Covers for I	th Gags 0117016 0116018 A ent Covers for Molt I	Number Description/Details th Gags 0117016 Adult 5-½" 0116018 A Child 4-½" ent Covers for Molt Mouth Gags

ACE MCKESSON TYPE SILICONE MOUTH PROPS

Silicone mouth props support the mandible and relieve muscle tension by providing a soft, pacifying surface to bite on. Reduces patient gagging and fatigue. Ideal for children and adults. Finger channel for insertion. Deep corrugation for secure stabilization. Virtually indestructible silicone will provide years of reusable service.

• Autoclavable • Not shown actual size • Silicone — Latex-free

/s Item Number	Manufacturer Number		Description/Details	Quantity
Small Chil	d			
1463453	0116036		Tan	pair
1463431	0116039		Blue	pair
1463452	0116035	C	Black	pair
Child				
1463456	0116038		Tan	pair
1336121	0116040		Blue	pair
1463454	0116037	D	Black	pair
Adult				
1463446	0117042		Tan	pair
1463457	0117043		Blue	pair
1463445	0117041	Ε	Black	pair

MOUTH PROPS AND RETRACTORS

ACE PLASTIC LIP/CHEEK RETRACTORS 🥯

- Autoclavable
- Can be steam sterilized to 270°F (130°C)
- Dry heat not recommended

/> Item Number	Manufacturer Number		Description/Details	Quantity	
Double-En	ded				
1463172	0100030	A	Universal, 4" (40mm)	2/pkg	

ACE SELF-SPANNING PLASTIC LIP/CHEEK RETRACTORS 🥯

- Autoclavable
- Latex-free
- Can be steam sterilized to 250°F (121°C) for 3 to 3-1/2 minutes
- Dry heat not recommended

Number	Manufacturer Number		Description/Details	Quantity	
Self-Spani	ning				
1463459	0116032	В	Child, 35mm	each	ĺ
1463401	0117031		Adult, 50mm	each	
Self-Spani	ning with Holde	rs			
1466773	0116025		Child, 35mm	each	
1463183	0117025	C	Adult, 50mm	each	

SAVE an additional 3% with scheduled deliveries on marked items. See inside front cover for details.

ACE STAINLESS STEEL LIP/CHEEK RETRACTORS

 Autoclav 	able			
/s Item Number	Manufacturer Number		Description/Details	Quantity
Austin, #2				
1466848	0100158	D	3" (7.5cm)	each
Austin Rot	ractor, Universi	tu of	Minnesota	
1466869	0100160	E	3" (7.5cm)	each
Branemarl	k® Type			
1150075	6100390	F	6-¼" (16cm)	each
Columbia				
1466833	0100391	G	5" (13cm), double ended	each
Feldman				
1466845	0100394	Н	5-½" (14cm), double ended	each
University	of Minnesota, T	Tongı	ie/Cheek	
1463398	0100088	1	5-1/2" (14cm), double ended	each
Weider				
1463177	0116095	J	5" (13cm), child	each
1463109	0117094		5" (13cm), adult	each





Small Child C

SUCTION AND IRRIGATION

— 1/4" 0.D.

28

ACE ADAPTERS FOR HVAC SUCTION TO SUCTION TUBING

/ Item Number	Manufacturer Number		Description/Details	Quantity
rambor	Trainibo.		Boot provide Botalie	dadinity
Brass				
1466889	6592968	Α	0.D. = 1/4" x 7/16"	each
2080061	0032973	В	0.D. = 5/16" x 7/16	each
01-1-1	N1			
Stainless S	Steel			
1464508	0033008	C	0.D. = 5/16" x 3/8"	each
1101000			0.2. 0,10 X 0,0	

ACE SUCTION TUBING ADAPTER

Allows conversion from high-volume to low volume suction

/> Item Number	Manufacturer Number		Description/Details	Quantity
1464221	0900211	D	Tubing adapter, 5mm head	each

SUCTION AND IRRIGATION

ACE SUCTION CONNECTING TUBING

- Special female connectors will fit all suction instruments
- Economical
- Excellent quality
- Single use
- Latex-free
- Non conductive

∕ E Item Number	Manufacturer Number	Description/Details	Size	Quantit
Sterile - C	ear with Blue Ti	nt		
1464526	0034238	3' long (0.9m)	1/4" I.D. (6.4mm)	20/case
1464527	0034239	6' long (1.8m)	1/4" I.D. (6.4mm)	20/case
1464496	0034240	10' long (3.0m)	1/4" I.D. (6.4mm)	20/case
Hides flu	paque Blue with id from patient		14" LD (6.4mm)	20/0200
	•	6' long (1.8m) 10' long (3.0m)	¼" I.D. (6.4mm) ¼" I.D. (6.4mm)	20/case
Hides flu14644881464490	id from patient 0034241	6' long (1.8m) 10' long (3.0m)		
Hides flu14644881464490	0034241 0034242	6' long (1.8m) 10' long (3.0m)		
 Hides flu 1464488 1464490 Non-Steril 	id from patient 0034241 0034242 e - Clear with Blo	6' long (1.8m) 10' long (3.0m) ue Tint	1/4" I.D. (6.4mm)	20/case

**Adapter allows connecting of dental saliva evacuation unit to suction tubing



1464427 0034244

1260495 0034245

SAVE an additional 3% with scheduled deliveries on marked items. See inside front cover for details.

1/4" I.D. (6.4mm)

1/4" I.D. (6.4mm)

25/case

25/case

ACE VARIJET™ SUCTION/IRRIGATION SYSTEM

6' long (1.8m)

10' long (3.0m)

- Single-handed suction with controlled irrigation
- Flush and clean site with one convenient instrument
- Improved vision to the site by using fewer instruments


ltem nber	Manufacturer Number	Description/Details	Quantity	
30030	0900090	3mm, 8" (20cm)	each	0











SUCTION AND IRRIGATION ACE SUCTION TIPS Manufacturer Size Description/Details Number Number **Oral Surgery** Oral Surgery Suction Handle **1466791** 0900040 6" (15cm) each 1466798 0900041 Oral Suction Tip, size 0, 1mm 4" (10cm) each 1464529 0900042 Oral Suction Tip, size 1, 2mm 4" (10cm) each Oral Suction Tip, size 3, 3mm 1464525 0900043 4" (10cm) each 1466796 Oral Suction Tonsil Tip, size 3 0900044 5-¾" (14.5cm) each Yankauer 1464377 1128004 Yankaur Suction Tip, complete 11" (28cm) each 1464521 0900129 Yankaur Suction Tip, complete, HVAC end 11" (28cm) each 1466859 1102006 Yankaur Suction Tip, replacement tip only 1" (2.5cm) each **HVAC Tips** Includes adapter for **1466733 0900024 I** HVAC Suction Tip, 4mm tip, 45° angle additional connection 7" (18cm) each **Pelton Crane 1466889** 6834301 Suction Tip, 3mm, with flow control 6-1/4" (16cm) each (continued next page) SAVE an additional 3% with scheduled deliveries on o marked items in the column. See inside front cover for details. REPLENISH W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

SUCTION AND IRRIGATION

Stylet only

ACE SUCTION TIPS

Manufacturer Number		Description/Details	\odot	Size	Quantity
0900060		Suction Tip, 6 FR w/stylet, 30° angle		7" (18cm)	each
0900062	A	Suction Tip, 8 FR w/stylet, 30° angle		7" (18cm)	each
0911072		Suction Tip, 8 FR w/stylet, 30° angle, delica	ate tip	7" (18cm)	each
0900063		Suction Tip, 9 FR w/stylet, 30° angle	•	7" (18cm)	each
0900064		Suction Tip, 10 FR w/stylet, 30° angle		7" (18cm)	each
0911073		Suction Tip, 10 FR w/stylet, 30° angle, delice	cate tip	7" (18cm)	each
0900065		Suction Tip, 11 FR w/stylet, 30° angle		7" (18cm)	each
0900066		Suction Tip, 12 FR w/stylet, 30° angle	•	7" (18cm)	each
0911074		Suction Tip, 12 FR w/stylet, 30° angle, delic	cate tip	7" (18cm)	each
0900070	В	Replacement Stylet		7-¼" (18.5cm)	each
0900071		Replacement Stylet with Bulb		7-1/4" (18.5cm)	each
	0900060 0900062 0911072 0900063 0900064 0911073 0900065 0900066 0911074 0900070	0900060 0900062 A 0911072 0900063 0900064 0911073 0900065 0900066 0911074 0900070 B	Number Description/Details 0900060 Suction Tip, 6 FR w/stylet, 30° angle 0900062 A Suction Tip, 8 FR w/stylet, 30° angle 0911072 Suction Tip, 8 FR w/stylet, 30° angle, delicated 0900063 Suction Tip, 9 FR w/stylet, 30° angle 0900064 Suction Tip, 10 FR w/stylet, 30° angle 0911073 Suction Tip, 10 FR w/stylet, 30° angle, delicated 0900065 Suction Tip, 11 FR w/stylet, 30° angle 0900066 Suction Tip, 12 FR w/stylet, 30° angle 0911074 Suction Tip, 12 FR w/stylet, 30° angle, delicated 0900070 B Replacement Stylet	Number Description/Details 0900060 Suction Tip, 6 FR w/stylet, 30° angle 0900062 A Suction Tip, 8 FR w/stylet, 30° angle 0911072 Suction Tip, 8 FR w/stylet, 30° angle, delicate tip 0900063 Suction Tip, 9 FR w/stylet, 30° angle 0900064 Suction Tip, 10 FR w/stylet, 30° angle 0911073 Suction Tip, 10 FR w/stylet, 30° angle, delicate tip 0900065 Suction Tip, 11 FR w/stylet, 30° angle 0900066 Suction Tip, 12 FR w/stylet, 30° angle 0911074 Suction Tip, 12 FR w/stylet, 30° angle, delicate tip 0900070 B Replacement Stylet	Number Description/Details Size 0900060 Suction Tip, 6 FR w/stylet, 30° angle 7" (18cm) 0900062 A Suction Tip, 8 FR w/stylet, 30° angle 7" (18cm) 0911072 Suction Tip, 8 FR w/stylet, 30° angle, delicate tip 7" (18cm) 0900063 Suction Tip, 9 FR w/stylet, 30° angle 7" (18cm) 0900064 Suction Tip, 10 FR w/stylet, 30° angle 7" (18cm) 0911073 Suction Tip, 10 FR w/stylet, 30° angle, delicate tip 7" (18cm) 0900065 Suction Tip, 11 FR w/stylet, 30° angle 7" (18cm) 0900066 Suction Tip, 12 FR w/stylet, 30° angle, delicate tip 7" (18cm) 0911074 Suction Tip, 12 FR w/stylet, 30° angle, delicate tip 7" (18cm) 0900070 B Replacement Stylet 7-¼" (18.5cm)

McKesson Type

1464507	0900001		McKesson Suction Tip, 2mm		7" (18cm)	each
1464498	0900003	C	McKesson Suction Tip, 4mm	•	7" (18cm)	each
1416837	0900004	D	McKesson Suction Tip, 4mm, w/sideport	•	7" (18cm)	each
1464522	0900005		McKesson Suction Tip, 3mm		9" (23cm)	each
1416843	0900006		McKesson Suction Tip, 4mm		9" (23cm)	each
Makasas	Tuna Madif	:1				

464520	0900002	Ε	McKesson Modified Suction Tip, 3mm	•	7" (18cm)	each



acesouthern.com/catalog/instruments24/page/30

SAVE an additional 3% with scheduled deliveries on marked items in the Golumn. See inside front cover for details.














Tapered tip is designed to fit the irrigation cannulas.



SUCTION AND IRRIGATION





ACE AUTOCLAVABLE CLEANING BRUSHES FOR SUCTION TIPS

- Ideal for cleaning stainless steel cannulae, suction tubes, endoscopy scopes
- For use with or without detergent

Manufacturer Number	Description/ Details	Size	Quantity
0916117	Small	3mm bristle diameter	12/pkg
0917118	Medium	5.3mm bristle diameter	12/pkg
0918119	Large	11mm bristle diameter	12/pkg
0916118	Small	3mm bristle diameter	12/pkg
0917119	Medium	5.3mm bristle diameter	12/pkg
	0916117 0917118 0918119 0916118	Number Details 0916117 Small 0917118 Medium 0918119 Large 0916118 Small	Number Details Size 0916117 Small 3mm bristle diameter 0917118 Medium 5.3mm bristle diameter 0918119 Large 11mm bristle diameter 0916118 Small 3mm bristle diameter

ACE STAINLESS STEEL IRRIGATION CANNULAS 🥯

- · Luer slip fits any luer type syringe
- 45° angle

/> Item Number	Manufacturer Number	Description/Details	Quantity
1-1/4" (3.	5cm)		
1466709	0900159	12 Gauge 🕒	each
1466799	0900162	14 Gauge 👝	each
1463159	0900165	16 Gauge 👝	each
1463146	0900168	18 Gauge 🛑	each
1466769	0900171	20 Gauge	each
2" (5cm)			
1463157	0900160	12 Gauge	each
1466801	0900163	14 Gauge	each
1276484	0900166	16 Gauge	each
1463148	0900169	18 Gauge 🔵	each
3" (7.5cm)		
1466792	0900161	12 Gauge	each
1463158	0900164	14 Gauge	each
1463149	0900167	16 Gauge 🛑	each
1276012	0900170	18 Gauge	each
3" (7.5cm) Straight		
1466710	0900270	18 Gauge, straight	each



SAVE an additional 3% with scheduled deliveries on smarked items. See inside front cover for details.

W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

SUCTION AND IRRIGATION

ACE BULB SYRINGE

- Single-handed suction with controlled irrigation
- Flush and clean site with one convenient instrument
- Improved vision to the site by using fewer instruments
- Order the Tips and Bulb separately

Number	Number		Description/Details	Quantity
Autoclava	ble Stainles	s Ste	eel Syringe Tips	
1464501	0900102	A	Water Syringe Tip - 1.9mm opening, 4-3/4" (12cm)	each
1464505	0900103	В	Air Syringe Tin - 2 2mm opening 4-1/2" (11cm)	each

Bulbs for ACE Syringe Tips

/s Item Manufacturer

Steam sterilize at 270°F for no longer than 20 minutes

Otodini otornizo di 27 o	1 10	The length than 20 minuted		
1415333 0900108	C	Bulb only, black, 60 cc	10/box	

ACE SELF-FILLING IRRIGATION SYRINGE

- Includes one 2" angled cannula (replacements sold individually on page 32)
- Stainless Steel

/s Item Number	Manufacturer Number		Description/Details	Quantity
1466829	0900101	D	10cc volume, 1mm opening	each

ACE COOK-WAITE TYPE ASPIRATING SYRINGES

· Accepts 1.8ml standard dental cartridges

/s Item Number	Manufacturer Number		Description/Details	Quantity
1463106	6000001	E	5-1/2" (14cm)	each

ACE GRAFT DELIVERY SYRINGES WITH FUNNEL

Stainless Steel

/> Item Number	Manufacturer Number		Description/Details	Quantity
1463421	9800035		3.5mm, curved	each
1463460	9800045	F	4.5mm, curved	each
1467007	9800060		6.0mm, curved	each

ACE STERILE DISPOSABLE SYRINGES

- Carries material to surgical site Single use
- Can be used for delivery of any material to surgical site

/> Item Number	Manufacturer Number		Description/Details	Quantity
1367126	0032019	G	Sterile, 1cc capacity	10/pack



SAVE an additional 3% with scheduled deliveries on marked items. See inside front cover for details.



SURGICAL KITS

SURGICAL KITS

ACE SINUS LIFT OSTEOTOME KIT

- Concave tip, sharply honed cutting surfaceConcave tip can be used to carry and place bone grafting material into site
- All osteotomes are 6", 15cm Long



steotome	Implant OD
2.8mm	3.25-3.5mm
3.3mm	3.75-4.0mm
4.35mm	4.75mm
5.0mm	5.75mm

CONCAVE TIP DETAIL Sharply honed cutting surface 4mm from tip - Concaved Tip

/> Item Number	Manufacturer Number	Description/Details	Quantity
Sinus Lift Os	steotome Kit - (C	oncave Tips)	
1150149	4545000	Set of 11 Osteotomes with Tray	kit
Sinus Lift Os	steotome Kit Inclu	des:	
1151573	4545001	Step - 1.6mm/2.0mm	each
1151574	4545101	Step - 2.0mm/2.8mm	each
1151575	4545201	Step - 2.8mm/3.3mm	each
1151576	4546001	Step - 3.3mm/4.0mm	each
1151577	4546101	Step - 4.0mm/4.35mm	each
1476169	4546301	Step - 4.35mm/5.0mm	each
1151578	4545701	Final - 2.8mm	each
1151579	4545301	Final - 3.3mm	each
1151580	4545401	Final - 4.0mm	each
1151581	4546201	Final - 4.35mm	each
1151582	4546401	Final - 5.0mm	each
		Autoclavable Tray (not sold individually)	

ACE SITE DILATOR KIT

- Convex tip compresses bone, limiting the amount of bone that must be removed from the site
- All implant site dilators are 6", 15cm long



Osteotome	Implant OD	
2.8mm	3.25-3.5mm	
3.3mm	3.75-4.0mm	
4.35mm	4.75mm	
5.0mm	5.75mm	
OD = Outside Diameter		

/ ► Item Number	Manufacturer Number	Description/Details	Quantit
Implant Site	Dilator Osteotom	e Kit - (Convex Tips)	
1150150	4543000	Set of 11 Osteotomes with Tray	kit
Implant Site	Dilating Osteoton	ne Kit Includes:	
1151583	4543001	Step - 1.6mm/2.0mm	each
1151584	4543101	Step - 2.0mm/2.8mm	each
1151585	4543201	Step - 2.8mm/3.3mm	each
1151586	4544001	Step - 3.3mm/4.0mm	each
1151587	4544101	Step - 4.0mm/4.35mm	each
1151588	4544301	Step - 4.35mm/5.0mm	each
1151589	4543701	Final - 2.8mm	each
1151590	4543301	Final - 3.3mm	each
1151591	4543401	Final - 4.0mm	each
1151592	4544201	Final - 4.35mm	each
1151593	4544401	Final - 5.0mm	each
		Autoclavable Tray (not sold individually)	





SURGICAL KITS


SURGICAL KITS

ACE OFFSET OSTEOTOME KIT

- Offset advantage: force of the strike stays directly on target
- Concave tip, sharply honed cutting surface
- Concave tip can be used to carry and place bone grafting material into site
- All osteotomes are 6", 15cm Long

The Advantage of ACE Offset Osteotomes

ACE Osteotomes vs Other Offset Osteotomes



Osteotome	Implant OD	
2.8mm	3.25-3.5mm	
3.3mm	3.75-4.0mm	
4.35mm	4.75mm	
5.0mm	5.75mm	
OD = Outside Diameter		






Offset Osteotome Kit Includes:

- Step 1.6mm/2.0mm
- Step 2.0mm/2.8mm
- Step 2.8mm/3.3mm
- Step 3.3mm/4.0mm
- Final 2.8mm
- Final 3.3mm
- Final 4.0mm
- Autoclavable Tray

Looking For Individual Osteotomes? Call Us at 800-624-5926.

ACE BONE SPREADER OSTEOTOME KIT

- Single end handle can be tapped with a mallet
- Markings at 13mm, 16mm and 19mm
- D-contoured tips to avoid fracturing weaker boney wall
- Instruments and tray are Stainless Steel



Set of 4 w	ith mallet and tra	ay	
	eader Osteotor		
/> Item Number	Manufacturer Number	Description/Details	Quantity

Spreader Osteotome Kit Includes:

- 4.3mm 27mm scoop
- 5.45mm 26mm scoop
- 7.5mm 26mm scoop
- 10.0mm 27mm scoopSurgical mallet, 6-½" (16.5cm)
- Instrument tray, 175 x 140 x 24mm

Looking For Individual Bone Spreaders? Call Us at 800-624-5926.

Especially designed for taking Chin and Ramus monocortical block grafts.


NEED PRICES? Follow the link below, or scan the QR code, to view retail pricing.

For account specific pricing, please log into your account on our website.

In Use: Here the Bone Spreader Osteotome is in place, separating the block graft from the ramus on the mandible. The site was prepared and the crest burnished. The crest is then perforated to outline the graft using a .75mm round bur. The outline should be 4-5mm medial from the buccal plate. The length of the graft is determined using a caliper.

A 702L bur is used to connect the perforations to the length of the graft. Vertical anterior and posterior cuts are made with the 702L bur 10-12mm in length. At the vertical depth, the buccal plate is scored using a #8 round bur. The block graft is out-fractured using the bone spreader osteotome.







acesouthern.com/catalog/instruments24/page/36

CONCAVE TIP DETAIL

15mm

cutting surface

4mm from tip

- Concaved Tip

SURGICAL KITS

SURGICAL KITS



Ridge Expanding Kit Includes:

- 2.5mm Expander
- 3.1mm Expander
- 3.6mm Expander
- 4.3mm Expander
- Connection Driver Tool (Hex), for CA handpiece
- Universal Ratchet
- Driver Extender (Hex)
- Square Titanium Bowl
- Autoclavable Tray

Looking For Individual Expanders? Call Us at 800-624-5926.



Ridge Expanding Chisel Kit

kit 1469995 4544700 Set of 4 with Tray And Handle

Description/Details

Chisel Kit Includes:

- 3.0 mm Ridge Chisel
- 4.0 mm Ridge Chisel
- 5.0 mm Ridge Chisel
- 9.5 mm Ridge Chisel
- Chisel Handle (Uses 1.25mm Hex Driver)
- Autoclavable Tray

Looking For Individual Chisels? Call Us at 800-624-5926.

TUTORIAL: SINGLE STAGE RIDGE SPLIT





Volume 14, Number 4.



drill is introduced to

mark implant site.

Reference Article: Palti A. Ridge Splitting and Implant Techniques

for the Anterior Maxilla. Dental Implantology Update 2003;

(Option for Soft Bone) An initial crest osteotomy can be started with an ACE Mini Blade (1466740 p.10)

3



The procedures and techniques presented throughout this catalog are

intended to provide basic reference only and should not be substituted for


the complete surgical procedure being performed by a licensed practitioner. A complete understanding and execution of the techniques and procedures

Quantity

Vertical osteotomy cuts performed with Ridge Expanding Chisel through the buccal cortex in an oblique fashion followed by horizontal osteotomy at ridge crest.



can be used to prepare site is completed the expanded osteotomy using chosen implant for implant placement. protocols. Implant



presented within this catalog are the full responsibility of the licensed practitioner. These procedures and techniques are NOT TO BE SUBSTITUTED for the practitioner's complete knowledge, education, experience or for any instrument or product instruction for use.





WIRE AND WIRE ACCESSORIES



ACE WIRE CUTTERS, SCISSORS AND TWISTERS

/s Item Number	Manufacturer Number	Description/Details	Quantit
WCPU Wire	e Cutter		
1466836	0805192	A 5" (13cm) tungsten carbide (max 0.5)	each
 1 blade is 1466766 	s serrated 0805260	B 7" (18cm), tungsten carbide	each
1466766		B 7" (18cm), tungsten carbide	each
1466851	0800213	6" (15cm)	each
Wire Twist	er		

ACE SPOOL WIRE ON AUTOCLAVABLE DELRIN SPOOLS

- 8% nickel in spool line
- 4 oz. Spool
- 200 feet
- Stainless steel

/s Item	Manufacturer		
Number	Number	Description/Details	Quantity
1466771	0800087	28G, .012"	each
1464947	0800088	26G, .016"	each
1430266	0800089	25G, .018"	each
1464945	0800090	24G, .020"	each

ACE TUBE WIRE

- 8% nickel in tube wire
- Stainless steel
- 18" strands

/> Item Number	Manufacturer Number	Description/Details	Quantity
1464946	0800100	26G, .016"	50/tube
1466826	0800101	25G, .018"	50/tube
1464948	0800102	24G, .020"	50/tube

ACE ARCH BAR

- Meter roll (approximately 36")
- Not MRI compatible

Number	Manufacturer Number	Description/Details	Quantity
1463458	0800110	Arch Bar	each

MISCELLANEOUS INSTRUMENTS

Trocar for Superior Chamber

В

ACE TMJ DYSFUNCTION TREATMENT

- Levaging TM Joint Shepherd Arthrocentesis Instrument
- A small incision (1/8") is made into the skin; mosquito forceps are used to make a blunt dissection to the lateral ligament; the arthrocentesis instrument may then enter the capsule into the superior chamber.

Number	Manufacturer Number	Description/Details	Quantity
Shepherd	Arthrocentesis lı	nstrument	
1423719	0702042	18G x 2.5"	each

ACE MALLETS

el each

ACE WOODSON PLASTIC FILLING INSTRUMENT

- Double Ended
- Stainless steel

Number	Manufacturer Number		Description/Details	Quantity
1467000	6200097A	C	#1, 6-1/2"	each



W: acesouthern.com | P: 1.800.624.5926 | M-F 8am - 7pm EST

INDEX

Adapters, Tubing28	Forceps, Bone Graft Holding	18
Antralplasty Elevator Kit13	Forceps, Chalazion Lip	19
Arch Bar40	Forceps, Dressing	18
Autografter4	Forceps, Extraction	20
Blade Handle10	Forceps, Hemostatic	18
Bone Crusher5	Forceps, Magill Catheter	19
3one File5	Forceps, Tissue	18
Bowl, Titanium4	Forceps, Towel	19
Bur Brush, Stainless8	Forceps Jar, Heavy Base	9
Burs, Surgical6	Gauge, Depth	12, 23
Caliper, Bone Implant22	Gauge, Patient Opening	24
Caliper, Ridge Mapping22	Implant Curette	4, 14
Calipers, Castroviejo22	Implant Site Dilator, Osteotome	35
Cannulas, Irrigation32	Instrument Trays	8-9, 12
Castroviejo Calipers22	lodine / Oil Cup	9
Clamps, Towel19	irrigation, Cannulas	32
Cleaning Brush, Nylon & Stainless8	Kit, Osteotome Bone Spreader	37
Cleaning Brush, Suction Tips	Kit, Osteotome Offset	36
Curettes, Antral13-14	Kit, Osteotome Sinus Lift	34
Curettes, Lucas14	Kit, Osteotome Site Dilator	35
Curettes, Miller14	Kit, Ridge Expanding	38
Curettes, Molt14	Kit, Ridge Expanding Chisel	39
Curettes, Sinus Lift12	Lip / Cheek Retractors	27
Depth Gauges & Probe12, 23	Mallets	11, 41
Oressing Jar9	Medicine Cup, Graduated	9
Elevators12-13, 15	Microsurgical Scissors	16
Elevators, Clepper Tatum Antralplasty12-13	Mirrors, Intraoral & Exam	25
Elevators, Periosteal12, 15	Mouth Gags	26
Elevators, Root Tip15	Mouth Prop, Exacta®	24
Elevators, T-Bar Handle15	Mouth Props, McKesson Type	26
Emesis Basin9	Needle Holders	17
Exacta® Mouth Prop24	Needle / Ointment Jar	9
Extraction Forceps20	Osteotomes, Sinus Lift	34
Extraction Rongeurs21	Osteotomes, Site Dilator	35
Filling Instrument, Woodson41	Periotomes	11

Periotome Kit	11
Pliers, Corn Suture	18
Pliers, Dressing	18
Plunger, Titanium	4
Replacement Covers, Mouth Gags	26
Retractor, Shepherd Sinus Lift	12
Rongeurs, Extraction	21
Ruler, X-Ray Magnification Calculator	24
Rulers, Plastic & Stainless	24
Scissors	16, 40
Screen Assembly, Autografter	4
Sinus Lift Kit	12, 34
Sponge Bowls	9
Square Titanium Bowl	8
Stainless Steel Ware	9
Suction Tips	4, 30-31
Suction Tubing	4, 29
Surgical Blades	10
Surgical Burs	6
Syringe, Aspirating	33
Syringe, Bulb & Tips	33
Syringe, Irrigation	33
Syringes, Graft Delivery	33
Tissue Forceps	18
Tissue Punches, Mucosa	11
TMJ Shepherd Arthrocentesis	41
Trephine Burs	7
Tubing, Suction	4, 29
Tumbler	9
$VARIJET^{\intercal_M}\ Suction/Irrigation\ System.$	29
Wire, Spools & Tubes	40
Wire Cutters	40
Wire Twister	40

"WHERE'S THE PRICING?"

We want to ensure you have the most current pricing available, but understand sending multiple catalogs a year in order to achieve that has a negative impact on the environment. Therefore, in order to create less waste, we have eliminated pricing from this catalog. Following the link below, or scanning the QR code, will allow you to view current retail pricing. In order to view pricing specific to your account, you must log in to your account on our website. This will also provide real-time product availability information. Plus, ordering is now just a click away!

Simply select the page number or sku in the index and begin!

Visit acesouthern.com/catalog/instruments24 today



THE HENRY SCHEIN PRICE POLICY:

We endeavor to maintain prices for the duration of a catalog, but we reserve the right to make price adjustments in response to manufacturers' price increases or extraordinary circumstances. Prices are subject to change without notice.

Henry Schein, Inc. ("Henry Schein") and customer agree that the terms and conditions hereinafter set forth shall govern the relationship between Henry Schein and the customer to the extent that the parties do not have a written agreement in effect that conflicts with such terms and conditions. Customer acknowledges and accepts all such terms and conditions by placing an order for goods with Henry Schein, and upon Henry Schein's delivery of the order to the customer.

CHOOSE YOUR PAYMENT METHOD

Reduce the cost and administration of paving Henry Schein-Pay electronically (ACH Debit) or set up AutoPay, Please call Customer Service for details.

For your convenience, we provide several payment alternatives. Orders billed to your account may be paid by ACH Debit, Check by Phone, or Check. If you prefer you may use your Henry Schein Credit Card. American Express, Visa, MasterCard or Discover Card when placing your order. All sales are subject to our normal terms and conditions

Unless otherwise instructed, check payments must be mailed to: Henry Schein, Inc. • Dept. CH 10241 • Palatine, IL 60055-0241

All sales are subject to credit approval. Invoices are payable within agreed terms of sale. Henry Schein and its affiliates reserve the right to offer a dual pricing system for payment by check, Henry Schein-approved ACH and/or Henry Schein-approved EFT vs. any other payment method and, to the extent permitted by law, reserves the right to surcharge for payment by credit card and/or other methods incurring a fee to Henry Schein

Open Accounts Receivable:

All unpaid accounts receivable past due are subject to a 1.5% finance charge.

DELIVERY TERMS:

Delivery dates are not guaranteed. Purchase orders are subject to product availability. Customers' sole and exclusive remedy for failure to deliver shall be refund of monies paid for the products at issue. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Except as noted below, title passes at the time the shipment is loaded at the shipper's dock.

California: For all shipments of goods to customers located within California, title will pass upon receipt of goods by California customers. Continental U.S.: All orders will be subject to a handling charge. This charge includes freight, except for additional carrier charges related to snecial delivery services and hazardous material shipments. Special orders are subject to additional freight charges.

Alaska, Hawaii & Pacific Protectorates: Standard shipping methods provide direct, reduced cost, expedited air delivery service to all accounts in Alaska and Hawaii. Customers in the Pacific Protectorates are offered direct surface transport, or postal services for reliable delivery. No additional surcharges apply, except when special services are requested. Low-level hazardous items (dangerous goods in accepted quantities and Consumer Commodity ID 8000) are now available via UPS 2nd-day air.

Guam, Puerto Rico, U.S. Trust Territories & Virgin Islands:

 All orders will be subject to a handling charge. This charge includes freight through the United States Postal Service (USPS) • Special delivery orders and hazardous material shipments can be shipped via United Parcel Service (UPS) for an additional charge. No minimum order amount or weight applies. Speak to your International Representative for details.

Outside U.S. (50 States): If your order is being shipped outside the LLS (50 states) please refer to the International Terms & Conditions at henryschein.com. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Title passes at the time the shipment is loaded at the shipper's dock. Customer is responsible for compliance with any applicable import requirements.

RX PRODUCTS & CONTROLLED SUBSTANCES

· Regulations require us to limit the sale of Rx and controlled substances only to registered, licensed healthcare professionals, If you are a new customer or have recently moved, please furnish us with a copy of your updated state and federal registrations verifying your shipping address. Please note that all orders for controlled substances are subject to a due diligence review process. Schedule II controlled substances can be ordered electronically or by mail. For more information on our Controlled Substance Ordering System please visit www.henryschein.com/e222; if you prefer to continue using Federal 222 Forms to order Schedule II controlled substances, mail the form to: Henry Schein, Inc. 5315 West 74th Street • Indianapolis, IN 46268

 Henry Schein restricts the sale or other transfer of medications to prisons/correctional facilities for use in lethal injections, based on our manufacturer agreements. The goods Henry Schein sells are intended to be used for their label-approved purposes or applicable standards of care, which do not include human lethal injection.

REGULATORY REQUIREMENTS:

Local regulatory requirements may apply to use or installation of certain products. Be sure to understand and comply with any such requirements prior to purchase, use, or installation of products.

THE DRUG SUPPLY CHAIN SECURITY ACT (DSCSA):

(MN, DM, WH, M2) The Drug Supply Chain Security Act (DSCSA) information related to prescription drug products is available on our website www.HenrySchein.com/pedigree. If you have any problems accessing our website or would like to receive a copy of DSCSA documentation via fax, mail or email, please contact our customer service department.

RETURNS

WE CANNOT ACCEPT ANY RETURNS WITHOUT PRIOR AUTHORIZATION.

To arrange for a return, simply call our Customer Service department or contact your Sales Consultant

The following conditions must be complied with:

 All returns must be accompanied by a copy of your invoice and a reason for the return • Merchandise must be returned in its original container, unmarked, and properly packaged . Returned products must have been purchased within the previous thirty (30) days. Any returns past thirty (30) days are subject to a restocking fee • Shortages or errors in shipments must be reported within seven (7) days of invoice date to issue credit (if applicable) • Shipping charges will apply on all returns

The following special, customized, or government-regulated items are not returnable:

· Opened handpieces, small equipment, and custom-ordered equipment • Special order items (products that we do not ordinarily stock) • Personalized and imprinted items • Opened computer hardware and software • Hazardous/flammable materials • Expired products • Items that cannot be returned to the manufacturer • Any item marked non-returnable •Items required to be shipped and stored frozen •Any drop-shipped products

Prescription Drug Returns:

Please note that, in order to comply with Federal and State traceability requirements, prescription drugs may be returned providing that the following key elements are met:

- 1) Returns of prescription drugs will only be accepted if Henry Schein is notified within 30 calendar days of shipment date and valid return authorization is issued by Henry Schein.
- 2) The Prescription Drug Marketing Act requires any customer returning prescription drugs to complete and return a Prescription Drug Return Authorization form. Federal law requires that the healthcare entity returning prescription drugs document that the product was kept under proper storage and handling conditions while in their possession and during the return of the product. To get a copy of the form and proper return authorization, please contact Customer
- 3) In addition, traceability regulations require that the healthcare entity returning prescription drugs certifies that the product being returned is the same exact product purchased from HSI.
- 4) Henry Schein will not issue credit for any returned prescription drugs which return was not authorized as provided herein, have been tampered with or where the labeling has been altered in any way.

INSTITUTIONAL. GOVERNMENT & CORPORATE ACCOUNTS:

Terms of Sale follow the same guidelines unless denoted differently in a contract. When applicable, freight is calculated by the weight. There is no minimum order; however, there is a nominal handling charge on orders under \$125. Some offers and promotions outlined in this catalog may not apply. Requests for bids and proposals may

Henry Schein, Inc. • Bid Department-E255 135 Duryea Rd. • Melville, NY 11747

INTERNATIONAL ORDERS:

We proudly serve healthcare professionals, governments, and dealers throughout the world. To place orders or for inquiries on export terms and conditions please contact the International Department (USA) at: export@henryschein.com

ARBITRATION: All Claims related to or arising under or relating to this Agreement (except for Claims for non-payment of amounts due to

TERMS AND CONDITIONS

Henry Schein for the sale of products) are to be exclusively and finally determined by binding arbitration in the state of New York, or another location mutually agreeable to the parties. Any and all Claims must be arbitrated on an individual basis, and there shall be no right or authority for any Claims or disputes to be arbitrated on a class action or collective basis. For avoidance of doubt, each party irrevocably waives any right to: (i) have any Claim resolved in connection with any class action or collective action, or (ii) recover any damages or relief directly or indirectly as part of any class action or collective action. The arbitration shall be conducted on a confidential basis pursuant to the Commercial Arbitration Rules of the American Arbitration Association, or if applicable, under its Procedures for Large, Complex Commercial Disputes. Any decision or award as a result of any such arbitration. proceeding shall be in writing and shall provide an explanation for all conclusions of law and fact and shall include the assessment of costs. expenses and reasonable attorneys' fees. Any such arbitration shall be conducted by an arbitrator experienced in the disputed subject matter and shall include a written record of the arbitration hearing. The parties reserve the right to object to any individual who shall be or has been at any time employed by or affiliated with a competing organization or entity. An award of arbitration may be confirmed in a court of competent jurisdiction. To the extent that any Claim or dispute is determined to not be subject to arbitration, all other Claims or disputes that would otherwise be subject to arbitration must be arbitrated. As used in this Agreement, "Claims" shall mean any and all liabilities, disputes and expenses whatsoever including, without limitation, claims, adversary proceedings (whether before a court, administrative agency or any other tribunal), damages (whether compensatory, multiple, exemplary or punitive), judgments, awards penalties settlements investigations costs responses to subpoenas or other governmental directives and reasonable attorneys' fees and disbursements with respect to any claims that may be sustained, suffered or incurred by a Party hereto.

WARRANTIES: Henry Schein will pass through to the customer, at the time of sale, any transferable product warranties, indemnities and remedies provided to Henry Schein by the applicable manufacturer. EXCEPT AS OTHERWISE PROVIDED HEREIN, TO THE EXTENT PERMITTED BY LAW, HENRY SCHEIN PROVIDES NO WARRANTIES OF ANY KIND EXPRESS OR IMPLIED INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE OR NON-INFRINGEMENT AND THE CUSTOMER SHALL LOOK TO THE MANUFACTURER OF THE PRODUCT FOR ANY WARRANTY THEREON.

LIMITATION OF LIABILITY: THE CUSTOMER AGREES TO LOOK SOLELY TO THE MANUFACTURER OF THE PRODUCT FOR ANY CLAIM ARISING DUE TO LOSS, INJURY, DAMAGE OR DEATH RELATED TO THE USE OR SALE OF PRODUCTS. HENRY SCHEIN SHALL NOT BE LIABLE FOR INDIRECT, INCIDENTAL, PUNITIVE, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, LOST PROFITS AND LOSS OF GOODWILL ARISING FROM OR RELATING TO ANY BREACH TO THIS AGREEMENT (OR ANY DUTY OF COMMON I AW. AND WHETHER OR NOT OCCASIONED BY THE NEGLIGENCE OF HENRY SCHEIN OR ITS AFFILIATES), REGARDLESS OF ANY NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.

DISCOUNTS, REBATES AND DISCLOSURES: Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion or discount program. You must fully and accurately report to Medicare. Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion or discount program, Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements

Terms of Sale may change without notice. For current information, go to www.henryschein.com/legalterms



Melville, NY 11747



Dental Surgical Instruments Catalog

