



ACE SOUTHERN

A HENRY SCHEIN® BRAND

KEEP YOUR PLAN IN PLACE

Reliable, Predictable, Pliable and Fully Resorbable.



conFORM

resorbable collagen membrane

- Fully resorbable with a predictable resorption time of 3-4 months¹
- Easy to handle, conforms to the surgical site without re-shaping and does not adhere to itself
- Guides healing and regeneration of bone and surrounding tissue
- Unique fiber orientation provides high tensile strength for *in situ* stability of membrane
- Cell occlusive, designed to retard and/or prevent epithelium down growth

¹ Reference on back



conFORM RESORBABLE COLLAGEN MEMBRANE

conFORM is a resorbable collagen membrane matrix made from Type I Bovine Collagen. Its pliable characteristics and predictable resorption time of 3-4 months make it an ideal membrane for guiding the healing of bone and surrounding tissue. When hydrated, conFORM will naturally adapt to surgical defects and maintain a high tensile strength, which makes it the ideal membrane for suturing.

With its in vivo stability, permeability and mechanical strength, conFORM is a reliable choice for various clinical applications. When used in conjunction with a bone grafting particulate, conFORM can help mitigate the migration of the particulate from the surgical site. Its resorbability eliminates the need for surgical removal. conFORM is packaged sterile, non-pyrogenic, and for single use.

Features

- White, non-friable membrane matrix engineered from highly purified Type I Collagen
- Fully resorbable with a predictable resorption time of 3-4 months ¹
- Guides healing and regeneration of bone and surrounding tissue
- Unique fiber orientation provides high tensile strength for *in situ* stability of membrane
- Cell occlusive, designed to retard and/or prevent epithelium down growth

Advantages


- Conforms to surgical site without re-shaping
- Easy to handle, does not adhere to itself
- Sufficient resorption time for sustained function in bone repair
- Can easily be trimmed during surgery for a final fit
- May be placed with either side against the tissue or bone

Typical Uses

- Ridge augmentation
- Bone regeneration around implants
- Sinus lift
- Periodontal defects
- Apical root resection
- Extraction sites

conFORM

1099383	15 x 20mm	each
1096584	20 x 30mm	each
1098764	30 x 40mm	each



Copyright © Henry Schein, Inc.
Henry Schein is a registered trademark of HS TM, LLC.

For the most up-to-date pricing, please contact ACE SOUTHERN or visit our website at www.acesouthern.com.

Footnote References:

- ¹ D.Yuen, J. Ulreich, G. Zuclich, H-B. Lin and S-T Li "Prediction of In Vivo Stability of a Resorbable, Reconstituted Type I Collagen Membrane by In Vitro Methods." 2000 Society for Biomaterials - Sixth World Congress.