



Accessories


incl. in scope of delivery



BKHEVSN40SWM
Heater Module



VSN40SBT
Bake-out tent



VSN40SFLIGHTCASE
Flightcase

CTD50
Liquid Nitrogen Dewar CF50
LN2 filling volume: ~0.5lt

VG40VAT (UNF)
Gatevalve with UNF threads

ADTSVSNWMXY
Wobblestick Fine
Alignment Adapterplate

WM40
Single shaft wobblestick
available linear travel:
~ 700mm

LSN1
Controller for
NEG activation

LSA3
Battery powered
ion pump controller

PUSAESNEXTorr®D200-5
Ion combination pump 200l/s H2





AFCF40KF40L25
CF - KF Adapter

VSN40SRKKF
Buffer Volume KF40

Grabber tools for various
Sample Holder systems
-Flagstyle Sample plates
-Dovetail Cry Shuttles
-Leica Blocks
-Atom Probe Pucks - available directly from CAMECA

Ra 0.8 (✓) $\begin{matrix} +0.2 \\ +0.1 \end{matrix}$ $\begin{matrix} -0.1 \\ -0.2 \end{matrix}$

Material: various		nominal ≤30mm		DIN ISO 2768-f-H	
		dimensions >30mm		DIN ISO 2768-m-K excluding Sym./Runout	
		Sym./Runout>30mm		DIN ISO 2768-H	
Ultra High Vacuum Cryo Transfer Suitcase - Long throw UHVCTSL		formed part dimensions		DIN ISO 2768-f-H	
		Scale		Drawn	L.Jeromin
		1 : 4 (1:15) (1:20)		Date	12.09.2023
				Changed	23.12.2025 LJ
Ferrovac <i>CH-8050 Zurich</i>		THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF FERROVAC. ANY REPRODUCTION IN PART OR AS WHOLE WITHOUT OUR PERMISSION IS PROHIBITED.		UHVCTSL A3 sheet 1/2	
		Rev. B			



Material: various

Ultra High Vacuum Cryo Transfer
Suitcase - Long throw
UHVCTSL

Ferrovac
CH-8050 Zurich

THE INFORMATION CONTAINED IN THIS
DRAWING IS THE SOLE PROPERTY OF
FERROVAC. ANY REPRODUCTION
IN PART OR AS WHOLE WITHOUT OUR
PERMISSION IS PROHIBITED.

nominal <30mm	DIN ISO 2768-f-H	
dimensions >30mm	DIN ISO 2768-m-K excluding Sym./Runout	
Sym./Runout >30mm	DIN ISO 2768-H	
formed part dimensions	DIN ISO 2768-f-H	
Scale	Drawn	L.Jeromin
1:5	Date	12.09.2023
	Changed	23.12.2025 LJ
Rev.	B	UHVCTSL A3 sheet 2/2