

Semi Rigid Ureteroscope

strauss[®]
SURGICAL

straussurgical.com

ATRAUMATIC DISTAL TIP

Specially designed distal tip allows for a smooth insert into the ureteric orifice offering a clear view of the ureter.

NO RE-FOCUSING

Due to special engineering, there is no need for camera re-focusing and no loss of picture while the scope is bent.

STRICT STANDARDS

In compliance with the rigorous protocols established by the World Health Organization (WHO), FDA and Health Canada.




- ✓ Simple lithotripsy possible: straight channel for rigid probes
- ✓ Outstanding performance (picture quality, insertion, irrigation)
- ✓ Safe and Reliable (robust, laser welded, easy to clean)
- ✓ Durability extends life time, resulting in low cost of ownership
- ✓ Easy to use due to unique click-lock mechanism
- ✓ Smooth insertion into ureteric orifice by special distal tip

TECHNICAL DATA

Item Number	Working	Working Channel	Working Length	Direction of View	Field of View	Image Bundle
STS-UR0-6075-310-HA	Distal 7.5 Fr. to 9.8 Fr.	1x5.0 Fr. or 2x3.0 Fr.	310 mm	8° to working channel	85°	30,000 px
STS-UR0-6075-425-HA	Distal 7.5 Fr. to 9.8 Fr.	1x5.0 Fr. or 2x3.0 Fr.	425 mm	8° to working channel	85°	30,000 px
STS-UR0-6060-310-HA	Distal 6.0 Fr. to 7.5 Fr.	1x4.0 Fr. or 2x2.4Fr.	310 mm	8° to working channel	100°	17,000 px
STS-UR0-6060-425-HA	Distal 6.0 Fr. to 7.5 Fr.	1 x4.0 Fr. or 2x2.4Fr.	425 mm	8° to working channel	100°	17,000 px

ACCESSORIES

Description	Working Channel	Field of View
Instrument Bridge	STS-UR0-0617-25LC1 or LC2	Yes
Three-slot membrane; pack of 10	STS-ACC-6000-MEMD-OOO	Yes
Cleaning Brush for STS-UR0-6060	STS-ACC-6060-520N-BRU	Sold individually
Light guide adapter for Karl Storz	STS-ACC-0000-LGUA	Yes
Light guide adapter for Richard Wolf	STS-ACC-0000-LGUA-RWF	Yes
EMS LC2 adapter; for all URS mentioned above	STS-ACC-6000-EMSA-LC2	Sold individually
EMS LC1 adapter; for all URS mentioned above	STS-ACC-6000-EMSA-LCI	Sold individually
Cleaning brushes for STS-UR0-6075	STS-ACC-6075-520N-BRU	Sold individually
Storage basket	STS-ACC-5600-TRYA-OOO	Sold individually



DETAILS

- > Highly optimized distal tip design
- > Made in Germany
- > 1 year warranty
- > Autoclavable: 18 min. at 134°C
- > Laser welded
- > Semi-rigid working length: 15° bending in all directions at full lifetime
- > Instrument bridge with lubrication-free stop cocks
- > Health Canada and FDA Certificated

strauss
SURGICAL